



Alcohol Testing Program

INSTRUMENT PROCESSING SHEET

Agency Vaporaiso PD

S/N 80-000742

Date In 10/4/16

Date Out 10/4/16

Ship P/U H/D CMI EE

Intake Performed By <u>PWS</u> <input type="checkbox"/> Registration <input checked="" type="checkbox"/> Annual <input type="checkbox"/> Return from CMI <input type="checkbox"/> Return from Enforcement Electronics <input type="checkbox"/> Other _____ Visual Inspection: <input checked="" type="checkbox"/> Case <input checked="" type="checkbox"/> Handle <input checked="" type="checkbox"/> Dry Gas Holder <input checked="" type="checkbox"/> Feet <input checked="" type="checkbox"/> Keyboard/Plug <input checked="" type="checkbox"/> Back/Plugs <input checked="" type="checkbox"/> Screws tight <input checked="" type="checkbox"/> Breath Hose <input checked="" type="checkbox"/> Other Equipment: <input type="checkbox"/> Power cord <input type="checkbox"/> Printer Cable <input type="checkbox"/> Other _____ Notes: _____ _____ _____		Quality Checks Performed By <u>PWS</u> <input checked="" type="checkbox"/> Breath Tube Screen <input checked="" type="checkbox"/> Replace O-Rings <input checked="" type="checkbox"/> Instrument Set Up Verified <input checked="" type="checkbox"/> R-Value <u>215</u> <input checked="" type="checkbox"/> Flow Verification (L/s) Flow Column # <u>ADP102</u> 32mm <u>148</u> (.139 - .169) 36mm <u>167</u> (.156 - .190) 53mm <u>238</u> (.228 - .278) 103mm <u>500</u> (.447 - .547) <input checked="" type="checkbox"/> Barometric Pressure Check Gauge ID # <u>28427</u> <input checked="" type="checkbox"/> Stability Checks		Flow Calibration Performed By _____ <input checked="" type="checkbox"/> Flow Calibration N/A <input type="checkbox"/> Flow Calibration Complete Flow Column # _____ <input type="checkbox"/> 5L/min - 17mm <input type="checkbox"/> 15L/min - <u>40mm</u> <input type="checkbox"/> 30L/min - 103mm <input type="checkbox"/> R-Value _____ <input type="checkbox"/> Post Calibration Verification (L/s) Flow Column # _____ 32mm _____ (.139 - .169) 36mm _____ (.156 - .190) 53mm _____ (.228 - .278) 103mm _____ (.447 - .547)														
		<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td>0.05</td> <td>SD1018</td> <td>201507A 7/14/17</td> </tr> <tr> <td>0.08</td> <td>SD1011</td> <td>201601F 1/26/18</td> </tr> <tr> <td>0.20</td> <td>SD1025</td> <td>201604C 4/5/18</td> </tr> <tr> <td>0.08 DGS</td> <td>N/A</td> <td>AG612405 5/3/18</td> </tr> </tbody> </table>		Simulator	Serial #	Lot #/Exp	0.05	SD1018	201507A 7/14/17	0.08	SD1011	201601F 1/26/18	0.20	SD1025	201604C 4/5/18	0.08 DGS	N/A	AG612405 5/3/18
Simulator	Serial #	Lot #/Exp																
0.05	SD1018	201507A 7/14/17																
0.08	SD1011	201601F 1/26/18																
0.20	SD1025	201604C 4/5/18																
0.08 DGS	N/A	AG612405 5/3/18																
		Maintenance Performed By _____ <input type="checkbox"/> Battery Replacement <input type="checkbox"/> Dry Gas Regulator Replacement <input type="checkbox"/> Breath Tube Replacement <input type="checkbox"/> Other _____																
		Suggested Service _____ _____																

RECEIVED
OCT 06 2016
FDLE
Alcohol Testing Program

Optical Bench Calibration Performed By <u>PWS</u> <input type="checkbox"/> Optical Bench Calibration N/A <input checked="" type="checkbox"/> Optical Bench Calibration Complete Barometric Pressure Gauge <u>1012</u> ID # <u>26932</u>			
Simulator	Serial Number	Lot Number	Expiration
0.000	SD1016	N/A	N/A
0.040	SD1024	16101	2/2/18
0.100	SD1013	16001	5/8/18
0.200	SD1012	16103	6/14/18
0.400	G6621	16102	3/22/18
0.080 DGS	N/A	1561508042	7/5/17
<input checked="" type="checkbox"/> Post Calibration Stability Checks			
Simulator	Serial Number	Lot Number	Expiration
0.05	SD1018	201507A	7/14/17
0.08	SD1011	201601F	1/26/18
0.20	SD1025	201604C	4/5/18
0.08 DGS	N/A	AG612405	5/3/18

Department Inspection Performed By <u>PWS</u> <input checked="" type="checkbox"/> Barometric Pressure <u>1014</u> Gauge ID# <u>28427</u> <u>1011</u> Instrument Mouth Alcohol Solution Lot # <u>2016-A</u> Acetone Stock Solution Lot # <u>2016-B</u>	
Simulator	Serial Number
0.00	SD1019
Interferent	SD1021
0.05	SD1018
0.08	SD1011
0.20	SD1025

Attachments <input checked="" type="checkbox"/> Form 41 <input checked="" type="checkbox"/> Pre-Stability Tests <input type="checkbox"/> Flow Calibration <input checked="" type="checkbox"/> Optical Bench Cal <input checked="" type="checkbox"/> Post-Stability Tests <input type="checkbox"/> Other _____	
--	--

Notes: QC:SP

 Instrument Complies with Chapter 11D-8, FAC
 Instrument Does Not Comply with Chapter 11D-8, FAC
 Return to/Place into Evidentiary Use
 Remain Out of Evidentiary Use
 Conduct an Agency Inspection Before Evidentiary Use

Brett Kirkland
Quality Control Review

10/6/16
Date

Stability Tests - Valparaiso PD # 80-000742 10/4/16
 Pre-Calibration

VALPARAISO PD
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-000742
 10/04/2016
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:17
Control Test	0.049	10:18
Air Blank	0.000	10:19
Control Test	0.049	10:19
Air Blank	0.000	10:20
Control Test	0.050	10:20
Air Blank	0.000	10:21
Control Test Stats		
Average	0.0493	
Std Dev	0.0006	
Rel Std Dev(%)	1.1703	

[Signature]
 Operator's Signature

VALPARAISO PD
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-000742
 10/04/2016
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:22
Control Test	0.079	10:23
Air Blank	0.000	10:23
Control Test	0.079	10:24
Air Blank	0.000	10:24
Control Test	0.079	10:25
Air Blank	0.000	10:25
Control Test Stats		
Average	0.0790	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

[Signature]
 Operator's Signature

VALPARAISO PD
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-000742
 10/04/2016
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:27
Control Test	0.193	10:27
Air Blank	0.000	10:28
Control Test	0.195	10:28
Air Blank	0.000	10:29
Control Test	0.195	10:30
Air Blank	0.000	10:30
Control Test Stats		
Average	0.1943	
Std Dev	0.0012	
Rel Std Dev(%)	0.5942	

[Signature]
 Operator's Signature

VALPARAISO PD
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-000742
 10/04/2016
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:31
Control Test	0.077	10:31
Air Blank	0.000	10:32
Control Test	0.077	10:32
Air Blank	0.000	10:33
Control Test	0.077	10:33
Air Blank	0.000	10:33
Control Test Stats		
Average	0.0770	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

[Signature]
 Operator's Signature

BK
 SP

Optical Bench
Calibration

V-Paraiso PD

#80-000742

10/4/16

Soil Value = 0.040 g/210L ***
Fit Value = 0.1905 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12645, Sum Io = 13305
<<<<< CHANNEL 1 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.7990 (-0.0240)
Sample #2 = 0.7980 (0.0160)
Sample #3 = 0.7730 (0.0040)
Sample #4 = 0.7830 (0.0310)
Avg % Abs = 0.7847 (0.0170)
STD DEV = 0.0126 (0.0135)
REL STD DEV = 1.604 (79.575)

<<<<< CHANNEL 2 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 1.4690 (-0.0040)
Sample #2 = 1.4700 (-0.0050)
Sample #3 = 1.4630 (-0.0030)
Sample #4 = 1.4880 (-0.0080)
Avg % Abs = 1.4737 (-0.0053)
STD DEV = 0.0129 (0.0025)
REL STD DEV = 0.875 (47.186)

Soil Value = 0.100 g/210L ***
Fit Value = 0.4762 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12639, Sum Io = 13303
<<<<< CHANNEL 1 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 1.8650 (-0.0260)
Sample #2 = 1.8690 (0.0060)
Sample #3 = 1.8680 (0.0010)
Sample #4 = 1.8590 (0.0390)
Avg % Abs = 1.8653 (0.0153)
STD DEV = 0.0055 (0.0206)
REL STD DEV = 0.295 (134.660)

<<<<< CHANNEL 2 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 3.5970 (-0.0090)
Sample #2 = 3.6010 (-0.0160)
Sample #3 = 3.5980 (-0.0180)
Sample #4 = 3.5910 (-0.0090)
Avg % Abs = 3.5967 (-0.0143)
STD DEV = 0.0051 (0.0047)
REL STD DEV = 0.143 (32.971)

Auto Calibration
Max Power Res Value = 12
Auto Range Res Value = 8

Soil Value = 0.000 g/210L ***
Fit Value = 0.0000 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12654, Sum Io = 13305
<<<<< CHANNEL 1 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.1030 (-0.0170)
Sample #2 = 0.0830 (0.0160)
Sample #3 = 0.0610 (0.0200)
Sample #4 = 0.0670 (0.0270)
Avg % Abs = 0.0703 (0.0210)
STD DEV = 0.0114 (0.0056)
REL STD DEV = 16.169 (26.513)

<<<<< CHANNEL 2 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.0480 (-0.0110)
Sample #2 = 0.0830 (-0.0240)
Sample #3 = 0.0910 (-0.0340)
Sample #4 = 0.0710 (-0.0220)
Avg % Abs = 0.0817 (-0.0267)
STD DEV = 0.0101 (0.0064)
REL STD DEV = 12.326 (24.105)

**** AUTO CAL DATA ****
<<<<< CHANNEL 1 >>>>>
Soil Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.070
Std Dev = 0.01 Rel Std Dev = 16.17
Soil Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 0.785
Std Dev = 0.01 Rel Std Dev = 1.60
Soil Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 1.865
Std Dev = 0.01 Rel Std Dev = 0.30
Soil Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 3.561
Std Dev = 0.01 Rel Std Dev = 0.22
Soil Val = 1.9048 mg/l or 0.400 g/210L
% Abs = 6.832
Std Dev = 0.03 Rel Std Dev = 0.42
Zero Order Coef = -114.71
First Order Coef = 2483.41
Second Order Coef = 46.68
Standard Deviation = 68.450943

<<<<< CHANNEL 2 >>>>>
Soil Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.082
Std Dev = 0.01 Rel Std Dev = 12.33
Soil Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 1.474
Std Dev = 0.01 Rel Std Dev = 0.88
Soil Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 3.597
Std Dev = 0.01 Rel Std Dev = 0.14
Soil Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 6.951
Std Dev = 0.00 Rel Std Dev = 0.03
Soil Val = 1.9048 mg/l or 0.400 g/210L
% Abs = 12.772
Std Dev = 0.05 Rel Std Dev = 0.42
Zero Order Coef = -44.47
First Order Coef = 1255.94
Second Order Coef = 18.61
Standard Deviation = 57.210785

Soil Value = 0.080 g/210L ***
Fit Value = 0.3810 mg/l %%%
Samples Taken = 4, Discarded = 1
**** CHANNEL 1

Sample #1 = 3226.00
Sample #2 = 3307.00
Sample #3 = 3275.00
Sample #4 = 3407.00
Average Result = 3329.6667
STD DEV = 68.8573
REL STD DEV = 2.068

**** CHANNEL 2

Sample #1 = 3397.00
Sample #2 = 3447.00
Sample #3 = 3425.00
Sample #4 = 3433.00
Average Result = 3435.0000
STD DEV = 11.1355
REL STD DEV = 0.324

Dry Gas H2O Adjust Results *****
Barometric Pressure = 1011
3 um H2O Adjust (mg/l*10,000) = 480
9 um H2O Adjust (mg/l*10,000) = 374
**** AUTO CAL PASS

SP

ASK

Stability Tests - Valparaiso PD #80-000742 10/4/16
 Post-Calibration

VALPARAISO PD
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-000742
 10/04/2016
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:32
Control Test	0.049	11:32
Air Blank	0.000	11:33
Control Test	0.049	11:34
Air Blank	0.000	11:34
Control Test	0.049	11:35
Air Blank	0.000	11:35
Control Test Stats		
Average	0.0490	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

[Signature]
 Operator's Signature

VALPARAISO PD
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-000742
 10/04/2016
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:37
Control Test	0.060	11:38
Air Blank	0.000	11:39
Control Test	0.060	11:39
Air Blank	0.000	11:40
Control Test	0.060	11:40
Air Blank	0.000	11:41
Control Test Stats		
Average	0.0600	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

[Signature]
 Operator's Signature

VALPARAISO PD
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-000742
 10/04/2016
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:42
Control Test	0.196	11:43
Air Blank	0.000	11:44
Control Test	0.198	11:44
Air Blank	0.000	11:45
Control Test	0.199	11:45
Air Blank	0.000	11:46
Control Test Stats		
Average	0.1977	
Std Dev	0.0015	
Rel Std Dev(%)	0.7728	

[Signature]
 Operator's Signature

VALPARAISO PD
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-000742
 10/04/2016
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:51
Control Test	0.083	11:51
Air Blank	0.000	11:51
Control Test	0.082	11:52
Air Blank	0.000	11:52
Control Test	0.082	11:53
Air Blank	0.000	11:53
Control Test Stats		
Average	0.0823	
Std Dev	0.0006	
Rel Std Dev(%)	0.7012	

TGS

[Signature]
 Operator's Signature

SP BK