



Alcohol Testing Program

INSTRUMENT PROCESSING SHEET

Agency Franklin CSO

S/N 80-000952

Date In 9/22/16

Date Out 9/22/16

Ship P/U H/D CMI EE

Intake	Performed By <u>PWS</u>	Quality Checks	Performed By	Flow Calibration	Performed By																								
<input type="checkbox"/> Registration <input checked="" type="checkbox"/> Annual <input type="checkbox"/> Return from CMI <input type="checkbox"/> Return from Enforcement Electronics <input type="checkbox"/> Other _____ Visual Inspection: <input checked="" type="checkbox"/> Case <u>OK</u> <input checked="" type="checkbox"/> Handle <u>OK</u> <input checked="" type="checkbox"/> Dry Gas Holder <u>OK</u> <input checked="" type="checkbox"/> Feet <u>OK</u> <input checked="" type="checkbox"/> Keyboard/Plug <u>OK</u> <input checked="" type="checkbox"/> Back/Plugs <u>OK</u> <input checked="" type="checkbox"/> Screws tight <u>OK</u> <input checked="" type="checkbox"/> Breath Hose <u>OK</u> Other Equipment: <input type="checkbox"/> Power cord <input type="checkbox"/> Printer Cable <input type="checkbox"/> Other _____ Notes: _____ _____ _____		<input checked="" type="checkbox"/> Breath Tube Screen <input checked="" type="checkbox"/> Replace O-Rings <input checked="" type="checkbox"/> Instrument Set Up Verified <input checked="" type="checkbox"/> R-Value <u>222</u> <input checked="" type="checkbox"/> Flow Verification (L/s) Flow Column # <u>ATP102</u> 32mm <u>148</u> (.139 - .169) 36mm <u>167</u> (.156 - .190) 53mm <u>230</u> (.228 - .278) 103mm <u>511</u> (.447 - .547) <input checked="" type="checkbox"/> Barometric Pressure Check Gauge ID # <u>28427</u> <input checked="" type="checkbox"/> Stability Checks		<input checked="" type="checkbox"/> Flow Calibration N/A <input type="checkbox"/> Flow Calibration Complete Flow Column # _____ <input type="checkbox"/> 5L/min - 17mm <input type="checkbox"/> 15L/min - 53mm <input type="checkbox"/> 30L/min - 103mm <input type="checkbox"/> R-Value _____ <input type="checkbox"/> Post Calibration Verification (L/s) Flow Column # _____ 32mm _____ (.139 - .169) 36mm _____ (.156 - .190) 53mm _____ (.228 - .278) 103mm _____ (.447 - .547)																									
		<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td>0.05</td> <td>SD1018</td> <td>201507A 7/14/17</td> </tr> <tr> <td>0.08</td> <td>SD1011</td> <td>201601F 1/26/18</td> </tr> <tr> <td>0.20</td> <td>SD1025</td> <td>201604C 4/5/18</td> </tr> <tr> <td>0.08 DGS</td> <td>N/A</td> <td>AG12405 5/3/18</td> </tr> </tbody> </table>	Simulator	Serial #	Lot #/Exp	0.05	SD1018	201507A 7/14/17	0.08	SD1011	201601F 1/26/18	0.20	SD1025	201604C 4/5/18	0.08 DGS	N/A	AG12405 5/3/18	<table border="1"> <thead> <tr> <th>Maintenance</th> <th>Performed By</th> </tr> </thead> <tbody> <tr> <td><input type="checkbox"/> Battery Replacement</td> <td></td> </tr> <tr> <td><input type="checkbox"/> Dry Gas Regulator Replacement</td> <td></td> </tr> <tr> <td><input type="checkbox"/> Breath Tube Replacement</td> <td></td> </tr> <tr> <td><input type="checkbox"/> Other _____</td> <td></td> </tr> </tbody> </table>		Maintenance	Performed By	<input type="checkbox"/> Battery Replacement		<input type="checkbox"/> Dry Gas Regulator Replacement		<input type="checkbox"/> Breath Tube Replacement		<input type="checkbox"/> Other _____	
Simulator	Serial #	Lot #/Exp																											
0.05	SD1018	201507A 7/14/17																											
0.08	SD1011	201601F 1/26/18																											
0.20	SD1025	201604C 4/5/18																											
0.08 DGS	N/A	AG12405 5/3/18																											
Maintenance	Performed By																												
<input type="checkbox"/> Battery Replacement																													
<input type="checkbox"/> Dry Gas Regulator Replacement																													
<input type="checkbox"/> Breath Tube Replacement																													
<input type="checkbox"/> Other _____																													

RECEIVED
SEP 26 2016
FDLE
Alcohol Testing Program

Optical Bench Calibration	Performed By <u>PWS</u>	Department Inspection	Performed By <u>PWS</u>																																								
<input type="checkbox"/> Optical Bench Calibration N/A <input checked="" type="checkbox"/> Optical Bench Calibration Complete Barometric Pressure Gauge <u>1015</u> ID# <u>26932</u>		<input checked="" type="checkbox"/> Barometric Pressure <u>1016</u> Gauge ID# <u>28427</u> <u>1014</u> Instrument Mouth Alcohol Solution Lot # <u>2016-A</u> Acetone Stock Solution Lot # <u>2016-B</u>																																									
<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td>SD1016</td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td>SD1024</td> <td>16101</td> <td>2/2/18</td> </tr> <tr> <td>0.100</td> <td>DR1279</td> <td>16001</td> <td>5/8/18</td> </tr> <tr> <td>0.200</td> <td>DR3856</td> <td>16103</td> <td>6/14/18</td> </tr> <tr> <td>0.400</td> <td>SD1013</td> <td>16102</td> <td>3/22/18</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>15615080A2</td> <td>7/5/17</td> </tr> </tbody> </table>	Simulator	Serial Number	Lot Number	Expiration	0.000	SD1016	N/A	N/A	0.040	SD1024	16101	2/2/18	0.100	DR1279	16001	5/8/18	0.200	DR3856	16103	6/14/18	0.400	SD1013	16102	3/22/18	0.080 DGS	N/A	15615080A2	7/5/17		<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.00</td> <td>SD1019</td> </tr> <tr> <td>Interferent</td> <td>SD1021</td> </tr> <tr> <td>0.05</td> <td>SD1018</td> </tr> <tr> <td>0.08</td> <td>SD1011</td> </tr> <tr> <td>0.20</td> <td>SD1025</td> </tr> </tbody> </table>	Simulator	Serial Number	0.00	SD1019	Interferent	SD1021	0.05	SD1018	0.08	SD1011	0.20	SD1025	
Simulator	Serial Number	Lot Number	Expiration																																								
0.000	SD1016	N/A	N/A																																								
0.040	SD1024	16101	2/2/18																																								
0.100	DR1279	16001	5/8/18																																								
0.200	DR3856	16103	6/14/18																																								
0.400	SD1013	16102	3/22/18																																								
0.080 DGS	N/A	15615080A2	7/5/17																																								
Simulator	Serial Number																																										
0.00	SD1019																																										
Interferent	SD1021																																										
0.05	SD1018																																										
0.08	SD1011																																										
0.20	SD1025																																										
<input checked="" type="checkbox"/> Post Calibration Stability Checks		<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.05</td> <td>SD1018</td> <td>201507A</td> <td>7/14/17</td> </tr> <tr> <td>0.08</td> <td>SD1011</td> <td>201601F</td> <td>1/26/18</td> </tr> <tr> <td>0.20</td> <td>SD1025</td> <td>201604C</td> <td>4/5/18</td> </tr> <tr> <td>0.08 DGS</td> <td>N/A</td> <td>AG12405</td> <td>5/3/18</td> </tr> </tbody> </table>		Simulator	Serial Number	Lot Number	Expiration	0.05	SD1018	201507A	7/14/17	0.08	SD1011	201601F	1/26/18	0.20	SD1025	201604C	4/5/18	0.08 DGS	N/A	AG12405	5/3/18																				
Simulator	Serial Number	Lot Number	Expiration																																								
0.05	SD1018	201507A	7/14/17																																								
0.08	SD1011	201601F	1/26/18																																								
0.20	SD1025	201604C	4/5/18																																								
0.08 DGS	N/A	AG12405	5/3/18																																								
<table border="1"> <thead> <tr> <th>Attachments</th> <th></th> </tr> </thead> <tbody> <tr> <td><input checked="" type="checkbox"/> Form 41</td> <td><input checked="" type="checkbox"/> Optical Bench Cal</td> </tr> <tr> <td><input checked="" type="checkbox"/> Pre-Stability Tests</td> <td><input checked="" type="checkbox"/> Post-Stability Tests</td> </tr> <tr> <td><input type="checkbox"/> Flow Calibration</td> <td><input type="checkbox"/> Other _____</td> </tr> </tbody> </table>		Attachments		<input checked="" type="checkbox"/> Form 41	<input checked="" type="checkbox"/> Optical Bench Cal	<input checked="" type="checkbox"/> Pre-Stability Tests	<input checked="" type="checkbox"/> Post-Stability Tests	<input type="checkbox"/> Flow Calibration	<input type="checkbox"/> Other _____																																		
Attachments																																											
<input checked="" type="checkbox"/> Form 41	<input checked="" type="checkbox"/> Optical Bench Cal																																										
<input checked="" type="checkbox"/> Pre-Stability Tests	<input checked="" type="checkbox"/> Post-Stability Tests																																										
<input type="checkbox"/> Flow Calibration	<input type="checkbox"/> Other _____																																										

Notes: OK 9/20/16

<input checked="" type="checkbox"/> Instrument Complies with Chapter 11D-8, FAC
<input type="checkbox"/> Instrument Does Not Comply with Chapter 11D-8, FAC
<input checked="" type="checkbox"/> Return to/Place into Evidentiary Use
<input type="checkbox"/> Remain Out of Evidentiary Use
<input checked="" type="checkbox"/> Conduct an Agency Inspection Before Evidentiary Use

Brett Hultstrand

Quality Control Review

9/26/16
Date

Stability Tests - Franklin CSO #80-000952 9/22/16
 Pre-Calibration

FRANKLIN COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-000952
 09/22/2016
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:40
Control Test	0.047	11:41
Air Blank	0.000	11:42
Control Test	0.049	11:42
Air Blank	0.000	11:43
Control Test	0.049	11:44
Air Blank	0.000	11:44
Control Test Stats		
Average	0.0483	
Std Dev	0.0012	
Rel Std Dev(%)	2.3890	

PWS

Operator's Signature

FRANKLIN COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-000952
 09/22/2016
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:45
Control Test	0.078	11:46
Air Blank	0.000	11:46
Control Test	0.078	11:47
Air Blank	0.000	11:48
Control Test	0.079	11:48
Air Blank	0.000	11:49
Control Test Stats		
Average	0.0783	
Std Dev	0.0006	
Rel Std Dev(%)	0.7370	

PWS

Operator's Signature

FRANKLIN COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-000952
 09/22/2016
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:36
Control Test	0.195	11:36
Air Blank	0.000	11:37
Control Test	0.196	11:38
Air Blank	0.000	11:38
Control Test	0.194	11:39
Air Blank	0.000	11:39
Control Test Stats		
Average	0.1950	
Std Dev	0.0010	
Rel Std Dev(%)	0.5128	

PWS

Operator's Signature

FRANKLIN COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-000952
 09/22/2016
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:53
Control Test	0.082	11:54
Air Blank	0.000	11:54
Control Test	0.083	11:54
Air Blank	0.000	11:55
Control Test	0.082	11:55
Air Blank	0.000	11:56
Control Test Stats		
Average	0.0823	
Std Dev	0.0006	
Rel Std Dev(%)	0.7012	

DAS

BK

PWS

Operator's Signature

Optical Bench
 Calibration
 #80-000952
 Franklin Co
 9/22/16

Sol Value = 0.1000 mg/l ***
 Fit Value = 0.1000 mg/l ****
 Samples Taken = 4, Discarded = 1
 Sum to = 12827, Sum to = 13126
 <<<< CHANNEL 1 >>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 0.9100 (-0.0300)
 Sample #2 = 0.9500 (-0.0340)
 Sample #3 = 0.9190 (-0.0000)
 Sample #4 = 0.8850 (-0.0270)
 Avg % Abs = 0.9180 (-0.0023)
 STD DEV = 0.0325 (0.0306)
 REL STD DEV = 3.542 (1310.008)

<<<< CHANNEL 2 >>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 1.5280 (-0.0460)
 Sample #2 = 1.5890 (-0.0380)
 Sample #3 = 1.5190 (-0.0090)
 Sample #4 = 1.5270 (-0.0020)
 Avg % Abs = 1.5450 (-0.0163)
 STD DEV = 0.0383 (0.0191)
 REL STD DEV = 2.480 (116.862)

FRANKLIN COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Model 8000
 SN 80-000952
 12:09:27

Auto Calibration
 Max Power Res Value = 54
 Auto Range Res Value = 38

Sol Value = 0.0000 mg/l ***
 Fit Value = 0.0000 mg/l ****
 Samples Taken = 4, Discarded = 1
 Sum to = 12851, Sum to = 13131
 <<<< CHANNEL 1 >>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 0.1790 (0.0090)
 Sample #2 = 0.1540 (0.0790)
 Sample #3 = 0.1590 (0.1210)
 Sample #4 = 0.1730 (0.1530)
 Avg % Abs = 0.1620 (0.1177)
 STD DEV = 0.0098 (0.0371)
 REL STD DEV = 6.080 (31.540)

<<<< CHANNEL 2 >>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 0.1560 (-0.0270)
 Sample #2 = 0.1470 (-0.0030)
 Sample #3 = 0.1670 (0.0150)
 Sample #4 = 0.1770 (0.0000)
 Avg % Abs = 0.1637 (0.0040)
 STD DEV = 0.0153 (0.0096)
 REL STD DEV = 9.333 (241.051)

Sol Value = 0.0000 mg/l ***
 Fit Value = 0.0000 mg/l ****
 Samples Taken = 4, Discarded = 1
 Sum to = 12820, Sum to = 13123
 <<<< CHANNEL 1 >>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 3.7210 (-0.0010)
 Sample #2 = 3.7350 (0.0160)
 Sample #3 = 3.7760 (0.0070)
 Sample #4 = 3.7480 (0.0270)
 Avg % Abs = 3.7530 (0.0167)
 STD DEV = 0.0210 (0.0100)
 REL STD DEV = 0.558 (60.100)

<<<< CHANNEL 2 >>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 6.7680 (-0.0070)
 Sample #2 = 6.7920 (0.0120)
 Sample #3 = 6.8130 (0.0190)
 Sample #4 = 6.8030 (0.0140)
 Avg % Abs = 6.8027 (0.0150)
 STD DEV = 0.0105 (0.0036)
 REL STD DEV = 0.154 (24.037)

Sol Value = 0.4000 g/210L ***
 Fit Value = 1.9048 mg/l ****
 Samples Taken = 4, Discarded = 1
 Sum to = 12819, Sum to = 13123
 <<<< CHANNEL 1 >>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 7.1100 (-0.0190)
 Sample #2 = 7.0780 (0.0150)
 Sample #3 = 7.0660 (0.0350)
 Sample #4 = 7.0530 (0.0240)
 Avg % Abs = 7.0637 (0.0247)
 STD DEV = 0.0129 (0.0100)
 REL STD DEV = 0.183 (40.608)

<<<< CHANNEL 2 >>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 12.7230 (-0.0210)
 Sample #2 = 12.6910 (0.0420)
 Sample #3 = 12.6870 (0.0510)
 Sample #4 = 12.7280 (0.0380)
 Avg % Abs = 12.6993 (0.0437)
 STD DEV = 0.0180 (0.0167)
 REL STD DEV = 0.142 (15.248)

***** AUTO CAL DATA *****
 <<<< CHANNEL 1 >>>>
 Sol Val = 0.0000 mg/l or 0.000 g/210L
 % Abs = 0.162
 Std Dev = 0.01 Rel Std Dev = 6.08
 Sol Val = 0.1905 mg/l or 0.040 g/210L
 % Abs = 0.918
 Std Dev = 0.03 Rel Std Dev = 3.54
 Sol Val = 0.4762 mg/l or 0.100 g/210L
 % Abs = 1.984
 Std Dev = 0.01 Rel Std Dev = 0.61
 Sol Val = 0.9524 mg/l or 0.200 g/210L
 % Abs = 3.753
 Std Dev = 0.02 Rel Std Dev = 0.56
 Sol Val = 1.9048 mg/l or 0.400 g/210L
 % Abs = 7.064
 Std Dev = 0.01 Rel Std Dev = 0.18
 Zero Order Coef = -421.15
 First Order Coef = 2536.07
 Second Order Coef = 31.13
 Standard Deviation = 21.615335

<<<< CHANNEL 2 >>>>
 Sol Val = 0.0000 mg/l or 0.000 g/210L
 % Abs = 0.164
 Std Dev = 0.02 Rel Std Dev = 9.33
 Sol Val = 0.1905 mg/l or 0.040 g/210L
 % Abs = 1.545
 Std Dev = 0.04 Rel Std Dev = 2.48
 Sol Val = 0.4762 mg/l or 0.100 g/210L
 % Abs = 3.547
 Std Dev = 0.01 Rel Std Dev = 0.24
 Sol Val = 0.9524 mg/l or 0.200 g/210L
 % Abs = 6.803
 Std Dev = 0.01 Rel Std Dev = 0.15
 Sol Val = 1.9048 mg/l or 0.400 g/210L
 % Abs = 12.699
 Std Dev = 0.02 Rel Std Dev = 0.14
 Zero Order Coef = -203.92
 First Order Coef = 1349.98
 Second Order Coef = 13.74
 Standard Deviation = 25.300144

Sol Value = 0.080 g/210L ***
 Fit Value = 0.3810 mg/l ****
 Samples Taken = 4, Discarded = 1
 ***** CHANNEL 1 *****
 Sample #1 = 2694.00
 Sample #2 = 2665.00
 Sample #3 = 2735.00
 Sample #4 = 2740.00
 Average Result = 2713.3333
 STD DEV = 41.9325
 REL STD DEV = 1.545
 ***** CHANNEL 2 *****
 Sample #1 = 3094.00
 Sample #2 = 3065.00
 Sample #3 = 3059.00
 Sample #4 = 3108.00
 Average Result = 3078.6667
 STD DEV = 25.8908
 REL STD DEV = 0.841

 Dry Gas H2O Adjust Results *****
 Barometric Pressure = 1014
 3 um H2O Adjust (mg/l * 10,000) = 1096
 9 um H2O Adjust (mg/l * 10,000) = 731
 ***** AUTO CAL PRESS *****

<<<< CHANNEL 2 >>>>
 Sol Val = 0.0000 mg/l or 0.000 g/210L
 % Abs = 0.164
 Std Dev = 0.02 Rel Std Dev = 9.33
 Sol Val = 0.1905 mg/l or 0.040 g/210L
 % Abs = 1.545
 Std Dev = 0.04 Rel Std Dev = 2.48
 Sol Val = 0.4762 mg/l or 0.100 g/210L
 % Abs = 3.547
 Std Dev = 0.01 Rel Std Dev = 0.24
 Sol Val = 0.9524 mg/l or 0.200 g/210L
 % Abs = 6.803
 Std Dev = 0.01 Rel Std Dev = 0.15
 Sol Val = 1.9048 mg/l or 0.400 g/210L
 % Abs = 12.699
 Std Dev = 0.02 Rel Std Dev = 0.14
 Zero Order Coef = -203.92
 First Order Coef = 1349.98
 Second Order Coef = 13.74
 Standard Deviation = 25.300144

ASK

Stability Tests
Post-Calibration

Franklin CSD #80-000952 9/22/16

FRANKLIN COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000952
09/22/2016
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:52
Control Test	0.049	12:53
Air Blank	0.000	12:53
Control Test	0.050	12:54
Air Blank	0.000	12:55
Control Test	0.050	12:55
Air Blank	0.000	12:56
Control Test Stats		
Average	0.0497	
Std Dev	0.0006	
Rel Std Dev(%)	1.1625	

[Signature]

Operator's Signature

FRANKLIN COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000952
09/22/2016
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:06
Control Test	0.079	13:07
Air Blank	0.000	13:07
Control Test	0.080	13:08
Air Blank	0.000	13:09
Control Test	0.080	13:09
Air Blank	0.000	13:10
Control Test Stats		
Average	0.0797	
Std Dev	0.0006	
Rel Std Dev(%)	0.7247	

[Signature]

Operator's Signature

FRANKLIN COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000952
09/22/2016
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:02
Control Test	0.198	13:02
Air Blank	0.000	13:03
Control Test	0.199	13:04
Air Blank	0.000	13:04
Control Test	0.198	13:05
Air Blank	0.000	13:05
Control Test Stats		
Average	0.1983	
Std Dev	0.0006	
Rel Std Dev(%)	0.2911	

[Signature]

Operator's Signature

FRANKLIN COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000952
09/22/2016
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:11
Control Test	0.081	13:11
Air Blank	0.000	13:12
Control Test	0.080	13:12
Air Blank	0.000	13:13
Control Test	0.080	13:13
Air Blank	0.000	13:13
Control Test Stats		
Average	0.0803	
Std Dev	0.0006	
Rel Std Dev(%)	0.7187	

[Signature]

Operator's Signature

DKS BK