

INSTRUMENT PROCESSING SHEET

Agency Haines City PD S/N 80-001007
 Date In 8/1/16 Date Out 8/10/16 Ship P/U H/D CMI EE

Intake <input type="checkbox"/> Registration <input checked="" type="checkbox"/> Annual <input type="checkbox"/> Return from CMI <input type="checkbox"/> Return from Enforcement Electronics <input type="checkbox"/> Other _____ Visual Inspection: <input checked="" type="checkbox"/> Case <input checked="" type="checkbox"/> Handle <input checked="" type="checkbox"/> Dry Gas Holder <input checked="" type="checkbox"/> Feet <input checked="" type="checkbox"/> Keyboard/Plug <input checked="" type="checkbox"/> Back/Plugs <input checked="" type="checkbox"/> Screws tight <input checked="" type="checkbox"/> Breath Hose Other Equipment: <input type="checkbox"/> Power cord <input type="checkbox"/> Printer Cable <input type="checkbox"/> Other _____ Notes: _____ _____ _____	Quality Checks Performed By <u>RWS</u> <input checked="" type="checkbox"/> Breath Tube Screen <input checked="" type="checkbox"/> Replace O-Rings <input checked="" type="checkbox"/> Instrument Set Up Verified <input checked="" type="checkbox"/> R-Value <u>188</u> <input checked="" type="checkbox"/> Flow Verification (L/s) Flow Column # <u>ATP102</u> 32mm <u>0.152</u> (.139 - .169) 36mm <u>0.167</u> (.156 - .190) 53mm <u>0.238</u> (.228 - .278) 103mm <u>0.500</u> (.447 - .547) <input checked="" type="checkbox"/> Barometric Pressure Check Gauge ID # <u>28427</u> <input checked="" type="checkbox"/> Stability Checks <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td>0.05</td> <td>SD1018</td> <td>201507A 7/14/17</td> </tr> <tr> <td>0.08</td> <td>SD1011</td> <td>201601F 1/26/18</td> </tr> <tr> <td>0.20</td> <td>SD1025</td> <td>201505A 5/12/17</td> </tr> <tr> <td>0.08 DGS</td> <td>N/A</td> <td>AG612405 5/3/18</td> </tr> </tbody> </table>	Simulator	Serial #	Lot #/Exp	0.05	SD1018	201507A 7/14/17	0.08	SD1011	201601F 1/26/18	0.20	SD1025	201505A 5/12/17	0.08 DGS	N/A	AG612405 5/3/18	Flow Calibration Performed By <u>RWS</u> <input checked="" type="checkbox"/> Flow Calibration N/A <input type="checkbox"/> Flow Calibration Complete Flow Column # _____ <input type="checkbox"/> 5L/min - 17mm <input type="checkbox"/> 15L/min - 17mm <input type="checkbox"/> 30L/min - 103mm <input type="checkbox"/> R-Value _____ <input type="checkbox"/> Post Calibration Verification (L/s) Flow Column # _____ 32mm _____ (.139 - .169) 36mm _____ (.156 - .190) 53mm _____ (.228 - .278) 103mm _____ (.447 - .547) Maintenance Performed By _____ <input type="checkbox"/> Battery Replacement <input type="checkbox"/> Dry Gas Regulator Replacement <input type="checkbox"/> Breath Tube Replacement <input type="checkbox"/> Other _____ Suggested Service _____ _____
Simulator	Serial #	Lot #/Exp															
0.05	SD1018	201507A 7/14/17															
0.08	SD1011	201601F 1/26/18															
0.20	SD1025	201505A 5/12/17															
0.08 DGS	N/A	AG612405 5/3/18															

RECEIVED
 AUG 11 2016
 FDLE
 Alcohol Testing Program

Optical Bench Calibration Performed By <u>RWS</u> <input type="checkbox"/> Optical Bench Calibration N/A <input checked="" type="checkbox"/> Optical Bench Calibration Complete Barometric Pressure Gauge <u>1016</u> ID# <u>26932</u> <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td>DRI275</td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td>G2882</td> <td>16101</td> <td>2/2/18</td> </tr> <tr> <td>0.100</td> <td>G2078</td> <td>16001</td> <td>5/8/18</td> </tr> <tr> <td>0.200</td> <td>G2408</td> <td>15104</td> <td>5/27/17</td> </tr> <tr> <td>0.400</td> <td>G5358</td> <td>16102</td> <td>3/22/18</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>DB415080A1</td> <td>3/5/17</td> </tr> </tbody> </table> <input checked="" type="checkbox"/> Post Calibration Stability Checks <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.05</td> <td>SD1018</td> <td>201507A</td> <td>7/14/17</td> </tr> <tr> <td>0.08</td> <td>SD1011</td> <td>201601F</td> <td>1/26/18</td> </tr> <tr> <td>0.20</td> <td>SD1025</td> <td>201505A</td> <td>5/12/17</td> </tr> <tr> <td>0.08 DGS</td> <td>N/A</td> <td>AG612405</td> <td>5/3/18</td> </tr> </tbody> </table>	Simulator	Serial Number	Lot Number	Expiration	0.000	DRI275	N/A	N/A	0.040	G2882	16101	2/2/18	0.100	G2078	16001	5/8/18	0.200	G2408	15104	5/27/17	0.400	G5358	16102	3/22/18	0.080 DGS	N/A	DB415080A1	3/5/17	Simulator	Serial Number	Lot Number	Expiration	0.05	SD1018	201507A	7/14/17	0.08	SD1011	201601F	1/26/18	0.20	SD1025	201505A	5/12/17	0.08 DGS	N/A	AG612405	5/3/18	Department Inspection Performed By <u>RWS</u> <input checked="" type="checkbox"/> Barometric Pressure <u>1017</u> Gauge ID# <u>28427</u> <u>1014</u> Instrument Mouth Alcohol Solution Lot # <u>2015-B A</u> Acetone Stock Solution Lot # <u>2016-B</u> <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.00</td> <td>G2879</td> </tr> <tr> <td>Interferent</td> <td>G8144</td> </tr> <tr> <td>0.05</td> <td>G11739</td> </tr> <tr> <td>0.08</td> <td>G8149</td> </tr> <tr> <td>0.20</td> <td>G11621</td> </tr> </tbody> </table> Attachments <input checked="" type="checkbox"/> Form 41 <input checked="" type="checkbox"/> Optical Bench Cal <input checked="" type="checkbox"/> Pre-Stability Tests <input checked="" type="checkbox"/> Post-Stability Tests <input type="checkbox"/> Flow Calibration <input type="checkbox"/> Other _____	Simulator	Serial Number	0.00	G2879	Interferent	G8144	0.05	G11739	0.08	G8149	0.20	G11621
Simulator	Serial Number	Lot Number	Expiration																																																										
0.000	DRI275	N/A	N/A																																																										
0.040	G2882	16101	2/2/18																																																										
0.100	G2078	16001	5/8/18																																																										
0.200	G2408	15104	5/27/17																																																										
0.400	G5358	16102	3/22/18																																																										
0.080 DGS	N/A	DB415080A1	3/5/17																																																										
Simulator	Serial Number	Lot Number	Expiration																																																										
0.05	SD1018	201507A	7/14/17																																																										
0.08	SD1011	201601F	1/26/18																																																										
0.20	SD1025	201505A	5/12/17																																																										
0.08 DGS	N/A	AG612405	5/3/18																																																										
Simulator	Serial Number																																																												
0.00	G2879																																																												
Interferent	G8144																																																												
0.05	G11739																																																												
0.08	G8149																																																												
0.20	G11621																																																												

Notes: Optical bench calibration completed to bring values closer to nominal. RWS
QC-RWS
Brett Kirkland
 Quality Control Review

Instrument Complies with Chapter 11D-8, FAC
 Instrument Does Not Comply with Chapter 11D-8, FAC
 Return to/Place into Evidentiary Use
 Remain Out of Evidentiary Use
 Conduct an Agency Inspection Before Evidentiary Use

8/11/16
Date

Pre-Calibration Stability Checks

#80-001007 Haines City P.D. 8/1/16 RNS

RNS

HAINES CITY PD
Intoxilyzer - Alcohol Analyzer
Model: 8000 SN 80-001007
08/01/2016
Software: 8100.27

HAINES CITY PD
Intoxilyzer - Alcohol Analyzer
Model: 8000 SN 80-001007
08/01/2016
Software: 8100.27

HAINES CITY PD
Intoxilyzer - Alcohol Analyzer
Model: 8000 SN 80-001007
08/01/2016
Software: 8100.27

HAINES CITY PD
Intoxilyzer - Alcohol Analyzer
Model: 8000 SN 80-001007
08/01/2016
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	17:55
Control Test	0.050	17:56
Air Blank	0.000	17:56
Control Test	0.049	17:57
Air Blank	0.000	17:57
Control Test	0.049	17:58
Air Blank	0.000	17:59
Control Test Stats		
Average	0.0493	
Std Dev	0.0006	
Rel Std Dev(%)	1.1703	

Test	g/210L	Time
Air Blank	0.000	18:00
Control Test	0.080	18:01
Air Blank	0.000	18:01
Control Test	0.079	18:02
Air Blank	0.000	18:03
Control Test	0.080	18:03
Air Blank	0.000	18:04
Control Test Stats		
Average	0.0797	
Std Dev	0.0006	
Rel Std Dev(%)	0.7247	

Test	g/210L	Time
Air Blank	0.000	18:06
Control Test	0.197	18:07
Air Blank	0.000	18:08
Control Test	0.195	18:08
Air Blank	0.000	18:09
Control Test	0.195	18:10
Air Blank	0.000	18:10
Control Test Stats		
Average	0.1957	
Std Dev	0.0012	
Rel Std Dev(%)	0.5901	

Test	g/210L	Time
Air Blank	0.000	18:12
Control Test	0.077	18:13
Air Blank	0.000	18:13
Control Test	0.078	18:14
Air Blank	0.000	18:14
Control Test	0.078	18:14
Air Blank	0.000	18:15
Control Test Stats		
Average	0.0777	
Std Dev	0.0006	
Rel Std Dev(%)	0.7434	

RNS BK

RNS

Operator's Signature

RNS

Operator's Signature

RNS

Operator's Signature

RNS

Operator's Signature

Optical Bench Calibration Data #80-001007 Haines City P.D. 8/1/16 DCS

HAINES CITY PD
 Intoxilizer - Alcohol Analyzer
 Model 8000
 SN 80-001007
 19:01:25

Auto Calibration
 Max Power Res Value = 21
 Auto Range Res Value = 18

**** AUTO CAL DATA ****
 <<<< CHANNEL 2 >>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 1.4950 (-0.0340)
 Sample #2 = 1.4770 (-0.0290)
 Sample #3 = 1.4740 (-0.0130)
 Sample #4 = 1.4570 (-0.0070)
 Avg % Abs = 1.4693 (-0.0163)
 STD DEV = 0.0108 (0.0114)
 REL STD DEV = 0.734 (69.627)

**** AUTO CAL DATA ****
 <<<< CHANNEL 1 >>>>
 Sol Val = 0.000 mg/l or 0.000 g/210L
 % Abs = 0.105
 Std Dev = 0.03 Rel Std Dev = 26.69
 Sol Val = 0.1905 mg/l or 0.040 g/210L
 % Abs = 0.804
 Std Dev = 0.01 Rel Std Dev = 0.76
 Sol Val = 0.4762 mg/l or 0.100 g/210L
 % Abs = 1.888
 Std Dev = 0.01 Rel Std Dev = 0.53
 Sol Val = 0.9524 mg/l or 0.200 g/210L
 % Abs = 3.686
 Std Dev = 0.02 Rel Std Dev = 0.44
 Sol Val = 1.9048 mg/l or 0.400 g/210L
 % Abs = 6.939
 Std Dev = 0.04 Rel Std Dev = 0.65
 Zero Order Coef = -202.25
 First Order Coef = 2528.12
 Second Order Coef = 35.13
 Standard Deviation = 64.000046

Sol Value = 0.000 g/210L ***
 Fit Value = 0.0000 mg/l %%%
 Samples Taken = 4, Discarded = 1
 Sum Io = 12980, Sum Io = 12511
 <<<< CHANNEL 1 >>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 0.1610 (-0.0340)
 Sample #2 = 0.0990 (0.0080)
 Sample #3 = 0.1350 (0.0260)
 Sample #4 = 0.0800 (0.0500)
 Avg % Abs = 0.1047 (0.0280)
 STD DEV = 0.0279 (0.0211)
 REL STD DEV = 26.689 (75.255)

Sol Value = 0.100 g/210L ***
 Fit Value = 0.4762 mg/l %%%
 Samples Taken = 4, Discarded = 1
 Sum Io = 12974, Sum Io = 12497
 <<<< CHANNEL 1 >>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 1.8870 (-0.0150)
 Sample #2 = 1.8770 (0.0030)
 Sample #3 = 1.8890 (0.0110)
 Sample #4 = 1.8970 (0.0190)
 Avg % Abs = 1.8877 (0.0110)
 STD DEV = 0.0101 (0.0080)
 REL STD DEV = 0.533 (72.727)

Sol Value = 0.400 g/210L ***
 Fit Value = 1.9048 mg/l %%%
 Samples Taken = 4, Discarded = 1
 Sum Io = 12965, Sum Io = 12502
 <<<< CHANNEL 1 >>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 7.1020 (-0.0020)
 Sample #2 = 6.9710 (0.1570)
 Sample #3 = 6.9590 (0.1640)
 Sample #4 = 6.8880 (0.2160)
 Avg % Abs = 6.9393 (0.1757)
 STD DEV = 0.0449 (0.0265)
 REL STD DEV = 0.646 (115.086)

<<<< CHANNEL 2 >>>>
 Sol Val = 0.000 mg/l or 0.000 g/210L
 % Abs = 0.073
 Std Dev = 0.02 Rel Std Dev = 20.89
 Sol Val = 0.1905 mg/l or 0.040 g/210L
 % Abs = 1.469
 Std Dev = 0.01 Rel Std Dev = 0.73
 Sol Val = 0.4762 mg/l or 0.100 g/210L
 % Abs = 3.606
 Std Dev = 0.01 Rel Std Dev = 0.37
 Sol Val = 0.9524 mg/l or 0.200 g/210L
 % Abs = 7.068
 Std Dev = 0.03 Rel Std Dev = 0.42
 Sol Val = 1.9048 mg/l or 0.400 g/210L
 % Abs = 13.174
 Std Dev = 0.05 Rel Std Dev = 0.35
 Zero Order Coef = -32.82
 First Order Coef = 1262.40
 Second Order Coef = 14.03
 Standard Deviation = 60.178032

<<<< CHANNEL 2 >>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 0.1630 (-0.0230)
 Sample #2 = 0.0590 (0.0780)
 Sample #3 = 0.0890 (0.0590)
 Sample #4 = 0.0700 (0.0730)
 Avg % Abs = 0.0727 (0.0780)
 STD DEV = 0.0152 (0.0098)
 REL STD DEV = 20.885 (14.070)

<<<< CHANNEL 2 >>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 3.6170 (-0.0330)
 Sample #2 = 3.5910 (0.0030)
 Sample #3 = 3.6120 (0.0040)
 Sample #4 = 3.6160 (0.0080)
 Avg % Abs = 3.6063 (0.0050)
 STD DEV = 0.0134 (0.0026)
 REL STD DEV = 0.372 (52.915)

<<<< CHANNEL 2 >>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 13.3400 (-0.0260)
 Sample #2 = 13.2190 (0.1460)
 Sample #3 = 13.1760 (0.1620)
 Sample #4 = 13.1260 (0.1910)
 Avg % Abs = 13.1737 (0.1663)
 STD DEV = 0.0465 (0.0228)
 REL STD DEV = 0.353 (13.714)

<<<< CHANNEL 2 >>>>
 Sol Val = 0.000 mg/l or 0.000 g/210L
 % Abs = 0.073
 Std Dev = 0.02 Rel Std Dev = 20.89
 Sol Val = 0.1905 mg/l or 0.040 g/210L
 % Abs = 1.469
 Std Dev = 0.01 Rel Std Dev = 0.73
 Sol Val = 0.4762 mg/l or 0.100 g/210L
 % Abs = 3.606
 Std Dev = 0.01 Rel Std Dev = 0.37
 Sol Val = 0.9524 mg/l or 0.200 g/210L
 % Abs = 7.068
 Std Dev = 0.03 Rel Std Dev = 0.42
 Sol Val = 1.9048 mg/l or 0.400 g/210L
 % Abs = 13.174
 Std Dev = 0.05 Rel Std Dev = 0.35
 Zero Order Coef = -32.82
 First Order Coef = 1262.40
 Second Order Coef = 14.03
 Standard Deviation = 60.178032

Sol Value = 0.040 g/210L ***
 Fit Value = 0.1905 mg/l %%%
 Samples Taken = 4, Discarded = 1
 Sum Io = 12977, Sum Io = 12498
 <<<< CHANNEL 1 >>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 0.8330 (-0.0240)
 Sample #2 = 0.8050 (-0.0040)
 Sample #3 = 0.7970 (0.0210)
 Sample #4 = 0.8090 (0.0040)
 Avg % Abs = 0.8037 (0.0070)
 STD DEV = 0.0061 (0.0128)
 REL STD DEV = 0.760 (182.388)

Sol Value = 0.200 g/210L ***
 Fit Value = 0.9524 mg/l %%%
 Samples Taken = 4, Discarded = 1
 Sum Io = 12974, Sum Io = 12496
 <<<< CHANNEL 1 >>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 3.7360 (-0.0120)
 Sample #2 = 3.6770 (0.0380)
 Sample #3 = 3.6770 (0.0680)
 Sample #4 = 3.7150 (0.0400)
 Avg % Abs = 3.6863 (0.0487)
 STD DEV = 0.0162 (0.0168)
 REL STD DEV = 0.439 (34.465)

Sol Value = 0.400 g/210L ***
 Fit Value = 1.9048 mg/l %%%
 Samples Taken = 4, Discarded = 1
 Sum Io = 12965, Sum Io = 12502
 <<<< CHANNEL 1 >>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 7.1020 (-0.0020)
 Sample #2 = 6.9710 (0.1570)
 Sample #3 = 6.9590 (0.1640)
 Sample #4 = 6.8880 (0.2160)
 Avg % Abs = 6.9393 (0.1757)
 STD DEV = 0.0449 (0.0265)
 REL STD DEV = 0.646 (115.086)

<<<< CHANNEL 2 >>>>
 Sol Val = 0.000 mg/l or 0.000 g/210L
 % Abs = 0.073
 Std Dev = 0.02 Rel Std Dev = 20.89
 Sol Val = 0.1905 mg/l or 0.040 g/210L
 % Abs = 1.469
 Std Dev = 0.01 Rel Std Dev = 0.73
 Sol Val = 0.4762 mg/l or 0.100 g/210L
 % Abs = 3.606
 Std Dev = 0.01 Rel Std Dev = 0.37
 Sol Val = 0.9524 mg/l or 0.200 g/210L
 % Abs = 7.068
 Std Dev = 0.03 Rel Std Dev = 0.42
 Sol Val = 1.9048 mg/l or 0.400 g/210L
 % Abs = 13.174
 Std Dev = 0.05 Rel Std Dev = 0.35
 Zero Order Coef = -32.82
 First Order Coef = 1262.40
 Second Order Coef = 14.03
 Standard Deviation = 60.178032

Solution Stats Quadratic Fit Chan 1

Act	Fit	Residual
g/210L	g/210L	g/210L
0.000	0.001	-0.0013
0.040	0.039	0.0011
0.100	0.099	0.0014
0.200	0.201	-0.0015
0.400	0.400	0.0003

Solution Stats Quadratic Fit Chan 2

Act	Fit	Residual
g/210L	g/210L	g/210L
0.000	0.001	-0.0012
0.040	0.039	0.0011
0.100	0.099	0.0013
0.200	0.201	-0.0014
0.400	0.400	0.0003

Sol Value = 0.080 g/210L ***
 Fit Value = 0.3810 mg/l %%%
 Samples Taken = 4, Discarded = 1
 ***** CHANNEL 1 *****
 Sample #1 = 3194.00
 Sample #2 = 3121.00
 Sample #3 = 3097.00
 Sample #4 = 3049.00
 Average Result = 3089.0000
 STD DEV = 36.6606
 REL STD DEV = 1.187
 ***** CHANNEL 2 *****
 Sample #1 = 3497.00
 Sample #2 = 3497.00
 Sample #3 = 3534.00
 Sample #4 = 3521.00
 Average Result = 3517.3333
 STD DEV = 18.7705
 REL STD DEV = 0.534
 ***** CHANNEL 1 *****
 Dry Gas H2O Adjust Results *****
 Barometric Pressure = 1016
 3 um H2O Adjust (mg/l*10,000) = 720
 9 um H2O Adjust (mg/l*10,000) = 292
 **** AUTO CAL PASS

gms

Post-Cal
Stability Checks # 80-001007 Haines City P.D. 8/1/16 ~~DAB~~

DAB

HAINES CITY PD
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001007
08/01/2016
Software: 8100.27

Test	9/210L	Time
Air Blank	0.000	20:17
Control Test	0.050	20:18
Air Blank	0.000	20:19
Control Test	0.050	20:19
Air Blank	0.000	20:20
Control Test	0.050	20:21
Air Blank	0.000	20:21
Control Test Stats		
Average	0.0500	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

HAINES CITY PD
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001007
08/01/2016
Software: 8100.27

Test	9/210L	Time
Air Blank	0.000	20:30
Control Test	0.060	20:31
Air Blank	0.000	20:31
Control Test	0.062	20:32
Air Blank	0.000	20:32
Control Test	0.062	20:33
Air Blank	0.000	20:34
Control Test Stats		
Average	0.0613	
Std Dev	0.0012	
Rel Std Dev(%)	1.4197	

HAINES CITY PD
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001007
08/01/2016
Software: 8100.27

Test	9/210L	Time
Air Blank	0.000	20:44
Control Test	0.201	20:45
Air Blank	0.000	20:45
Control Test	0.198	20:46
Air Blank	0.000	20:47
Control Test	0.198	20:47
Air Blank	0.000	20:48
Control Test Stats		
Average	0.1990	
Std Dev	0.0017	
Rel Std Dev(%)	0.8704	

HAINES CITY PD
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001007
08/01/2016
Software: 8100.27

Test	9/210L	Time
Air Blank	0.000	19:51
Control Test	0.061	19:51
Air Blank	0.000	19:52
Control Test	0.060	19:52
Air Blank	0.000	19:52
Control Test	0.061	19:53
Air Blank	0.000	19:53
Control Test Stats		
Average	0.0607	
Std Dev	0.0006	
Rel Std Dev(%)	0.7157	

QNS
13K

DAB

Operator's Signature

DAB

Operator's Signature

DAB

Operator's Signature

DAB

Operator's Signature