



Alcohol Testing Program

INSTRUMENT PROCESSING SHEET

Agency DeLand PD S/N 80-001136
Date In 2/4/16 Date Out 2/14/16 3/16/16 4/26/16

Intake Performed By [Signature]
Registration
Annual
Return from CMI
Return from Enforcement Electronics
Other
Visual Inspection:
Case OK Handle OK
Dry Gas Holder OK Feet OK
Keyboard/Plug OK Back/Plugs OK
Screws tight OK Breath Hose OK
Other Equipment:
Power cord
Printer Cable
Other Static Bag

Quality Checks Performed By [Signature]
Breath Tube Screen
Replace O-Rings
Instrument Set Up Verified
R-Value 206
Flow Verification (L/s)
Flow Column # AT1105
32mm 152 (.139 -.169)
36mm 174 (.156 -.190)
53mm 251 (.228 -.278)
103mm 511 (.447 -.547)
Barometric Pressure Check
Gauge ID # 28427
Stability Checks
Table with Simulator, Serial #, Lot #/Exp

Flow Calibration Performed By [Signature]
Flow Calibration N/A
Flow Calibration Complete
Flow Column #
5L/min - 17mm
15L/min - 53mm
30L/min - 103mm
R-Value
Post Calibration Verification (L/s)
Flow Column #
32mm (.139 -.169)
36mm (.156 -.190)
53mm (.228 -.278)
103mm (.447 -.547)
Maintenance Performed By
Battery Replacement
Dry Gas Regulator Replacement
Breath Tube Replacement
Other
Suggested Service

RECEIVED APR 26 2016 FDLE Alcohol Testing Program

Optical Bench Calibration Performed By [Signature]
Optical Bench Calibration N/A
Optical Bench Calibration Complete
Barometric Pressure Gauge 1014 ID# 26932
Table with Simulator, Serial Number, Lot Number, Expiration
Post Calibration Stability Checks
Table with Simulator, Serial Number, Lot Number, Expiration

Department Inspection Performed By [Signature]
Barometric Pressure ID# 28427 Gauge 1015/1014 Instrument 1017/1014
Mouth Alcohol Solution Lot # 2015-A
Acetone Stock Solution Lot # 2015-B
Table with Simulator, Serial Number
Attachments
Form 41
Pre-Stability Tests
Flow Calibration
Optical Bench Cal
Post-Stability Tests
Other G.M.C.

Notes: QC-18K
QC-205

Instrument Complies with Chapter 11D-8, FAC
Instrument Does Not Comply with Chapter 11D-8, FAC
Return to/Place into Evidentiary Use
Remain Out of Evidentiary Use
Conduct an Agency Inspection Before Evidentiary Use

Brett Wickland

Quality Control Review

4/26/16

Date

80-001136  
 Stability Checks  
 2/2/16

INTOXILYZER 8000  
 Instrument Initialization  
 08:08 02/08/2016

DELAND PD  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-001136  
 02/08/2016  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:52
Control Test	0.050	08:53
Air Blank	0.000	08:53
Control Test	0.051	08:54
Air Blank	0.000	08:54
Control Test	0.051	08:55
Air Blank	0.000	08:56
Control Test Stats		
Average	0.0507	
Std Dev	0.0006	
Rel Std Dev(%)	1.1395	

Operator's Signature

DELAND PD  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-001136  
 02/08/2016  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:57
Control Test	0.080	08:57
Air Blank	0.000	08:58
Control Test	0.081	08:58
Air Blank	0.000	08:59
Control Test	0.081	08:59
Air Blank	0.000	08:59
Control Test Stats		
Average	0.0807	
Std Dev	0.0006	
Rel Std Dev(%)	0.7157	

Operator's Signature

DELAND PD  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-001136  
 02/08/2016  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:21
Control Test	0.206	09:21
Air Blank	0.000	09:22
Control Test	0.206	09:23
Air Blank	0.000	09:23
Control Test	0.206	09:23
Air Blank	0.000	09:24
Control Test	0.206	09:24
Air Blank	0.000	09:25
Control Test Stats		
Average	0.2060	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

Operator's Signature

BK

DGS

Operator's Signature

INTOXILYZER 8000  
Instrument Initialization  
08:12 04/26/2016

DELAND PD  
Intoxilyzer - Alcohol Analyzer  
Model 8000  
SN 80-001136  
08:36:31  
04/26/2016

Auto Calibration  
Max Power Res Value = 14  
Auto Range Res Value = 10

Sol Value = 0.000 g/210L \*\*\*  
Fit Value = 0.000 mg/l \*\*\*\*  
Samples Taken = 4, Discarded = 1  
Sum Io = 12683, Sum Io = 14130

\*\*\*\* CHANNEL 1 \*\*\*\*

Sample % Abs (% Abs Ref)  
Sample #1 = 0.0950 (0.0070)  
Sample #2 = 0.1020 (0.0080)  
Sample #3 = 0.0970 (0.0040)  
Sample #4 = 0.0950 (0.0670)  
Avg % Abs = 0.0980 (0.0397)  
STD DEV = 0.0036 (0.0297)  
REL STD DEV = 3.679 (74.969)

\*\*\*\* CHANNEL 2 \*\*\*\*

Sample % Abs (% Abs Ref)  
Sample #1 = 0.0970 (0.0020)  
Sample #2 = 0.0850 (0.0140)  
Sample #3 = 0.0920 (0.0140)  
Sample #4 = 0.0850 (0.0300)  
Avg % Abs = 0.0873 (0.0193)  
STD DEV = 0.0040 (0.0992)  
REL STD DEV = 4.526 (47.783)

Sol Value = 0.000 g/210L \*\*\*  
Fit Value = 0.000 mg/l \*\*\*\*  
Samples Taken = 4, Discarded = 1  
Sum Io = 12643, Sum Io = 14103

\*\*\*\* CHANNEL 1 \*\*\*\*

Sample % Abs (% Abs Ref)  
Sample #1 = 0.8030 (0.0160)  
Sample #2 = 0.8250 (0.0130)  
Sample #3 = 0.8200 (0.0340)  
Sample #4 = 0.8160 (0.0540)  
Avg % Abs = 0.8203 (0.0337)  
STD DEV = 0.0045 (0.0205)  
REL STD DEV = 0.550 (60.897)

\*\*\*\* CHANNEL 2 \*\*\*\*

Sample % Abs (% Abs Ref)  
Sample #1 = 1.5280 (0.0000)  
Sample #2 = 1.5310 (0.0250)  
Sample #3 = 1.5390 (0.0320)  
Sample #4 = 1.5130 (0.0560)  
Avg % Abs = 1.5277 (0.0377)  
STD DEV = 0.0133 (0.0163)  
REL STD DEV = 0.872 (43.164)

Sol Value = 0.100 g/210L \*\*\*  
Fit Value = 0.4762 mg/l \*\*\*\*  
Samples Taken = 4, Discarded = 1  
Sum Io = 12635, Sum Io = 14086

\*\*\*\* CHANNEL 1 \*\*\*\*

Sample % Abs (% Abs Ref)  
Sample #1 = 1.9230 (-0.0110)  
Sample #2 = 1.9280 (0.0260)  
Sample #3 = 1.9080 (0.0460)  
Sample #4 = 1.8960 (0.0500)  
Avg % Abs = 1.9107 (0.0440)  
STD DEV = 0.0162 (0.0171)  
REL STD DEV = 0.846 (38.836)

\*\*\*\* CHANNEL 2 \*\*\*\*

Sample % Abs (% Abs Ref)  
Sample #1 = 3.7470 (-0.0270)  
Sample #2 = 3.7030 (0.0350)  
Sample #3 = 3.6900 (0.0440)  
Sample #4 = 3.6960 (0.0450)  
Avg % Abs = 3.6963 (0.0413)  
STD DEV = 0.0065 (0.0055)  
REL STD DEV = 0.176 (13.325)

\*\*\*\* AUTO CAL DATA \*\*\*\*  
\*\*\*\* CHANNEL 1 \*\*\*\*

Sol Val = 0.0000 mg/l or 0.000 g/210L  
% Abs = 0.098  
Std Dev = 0.00 Rel Std Dev = 3.68

Sol Val = 0.1905 mg/l or 0.040 g/210L  
% Abs = 0.820  
Std Dev = 0.00 Rel Std Dev = 0.55

Sol Val = 0.4762 mg/l or 0.100 g/210L  
% Abs = 1.911  
Std Dev = 0.02 Rel Std Dev = 0.85

Sol Val = 0.9524 mg/l or 0.200 g/210L  
% Abs = 3.683  
Std Dev = 0.02 Rel Std Dev = 0.44

Sol Val = 1.9048 mg/l or 0.400 g/210L  
% Abs = 6.940  
Std Dev = 0.01 Rel Std Dev = 0.17

Zero Order Coef = -214.48  
First Order Coef = 2517.12  
Second Order Coef = 37.06  
Standard Deviation = 32.666877

\*\*\*\* CHANNEL 2 \*\*\*\*

Sol Val = 0.0000 mg/l or 0.000 g/210L  
% Abs = 0.087  
Std Dev = 0.00 Rel Std Dev = 4.63

Sol Val = 0.1905 mg/l or 0.040 g/210L  
% Abs = 1.528  
Std Dev = 0.01 Rel Std Dev = 0.87

Sol Val = 0.4762 mg/l or 0.100 g/210L  
% Abs = 3.696  
Std Dev = 0.01 Rel Std Dev = 0.18

Sol Val = 0.9524 mg/l or 0.200 g/210L  
% Abs = 7.112  
Std Dev = 0.01 Rel Std Dev = 0.20

Sol Val = 1.9048 mg/l or 0.400 g/210L  
% Abs = 13.188  
Std Dev = 0.02 Rel Std Dev = 0.15

Zero Order Coef = -71.09  
First Order Coef = 1242.91  
Second Order Coef = 15.63  
Standard Deviation = 35.572418

\*\*\*\* CHANNEL 1 \*\*\*\*

Sample % Abs (% Abs Ref)  
Sample #1 = 3314.00  
Sample #2 = 3333.00  
Sample #3 = 3358.00  
Sample #4 = 3335.00  
Average Result = 3342.0000  
STD DEV = 13.8924  
REL STD DEV = 0.416

\*\*\*\* CHANNEL 2 \*\*\*\*

Sample % Abs (% Abs Ref)  
Sample #1 = 13.3190 (-0.0050)  
Sample #2 = 13.1900 (0.1780)  
Sample #3 = 13.1670 (0.1940)  
Sample #4 = 13.2060 (0.1930)  
Avg % Abs = 13.1877 (0.1683)  
STD DEV = 0.0196 (0.0190)  
REL STD DEV = 0.149 (4.755)

80-001136  
Optical Bench Calibration  
4/26/16

SK

Solution Stats Quadratic Fit Chan 2

Act Fit Residual  
g/210L g/210L g/210L  
0.000 0.001 -0.0008  
0.040 0.039 0.0009  
0.100 0.099 0.0005  
0.200 0.201 -0.0008  
0.400 0.400 0.0002

Sol Value = 0.080 g/210L \*\*\*

Fit Value = 0.3810 mg/l \*\*\*\*

Samples Taken = 4, Discarded = 1

\*\*\*\* CHANNEL 1 \*\*\*\*

Sample #1 = 2956.00  
Sample #2 = 2937.00  
Sample #3 = 3019.00  
Sample #4 = 2995.00  
Average Result = 2983.6667  
STD DEV = 42.1584  
REL STD DEV = 1.413

\*\*\*\* CHANNEL 2 \*\*\*\*

Sample #1 = 3314.00  
Sample #2 = 3333.00  
Sample #3 = 3358.00  
Sample #4 = 3335.00  
Average Result = 3342.0000  
STD DEV = 13.8924  
REL STD DEV = 0.416

\*\*\*\* CHANNEL 1 \*\*\*\*

Dry Gas H2O Adjust Results \*\*\*\*\*  
Barometric Pressure = 1014  
3 um H2O Adjust (mg/l\*10,000) = 826  
9 um H2O Adjust (mg/l\*10,000) = 467

\*\*\*\* AUTO CAL PASS \*\*\*\*

Solution Stats Quadratic Fit Chan 1

Act Fit Residual  
g/210L g/210L g/210L  
0.000 0.001 -0.0007  
0.040 0.039 0.0005  
0.100 0.099 0.0005  
0.200 0.201 -0.0008  
0.400 0.400 0.0002

80-001136  
 Post of Bench Cal. Stability checks  
 4/24/16

DELAND PD  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000  
 04/26/2016  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:29
Control Test	0.080	09:29
Air Blank	0.000	09:30
Control Test	0.080	09:30
Air Blank	0.000	09:31
Control Test	0.080	09:31
Air Blank	0.000	09:31
Control Test Stats		
Average	0.0800	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

DGS.

Operator's Signature

DELAND PD  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000  
 04/26/2016  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:33
Control Test	0.051	09:34
Air Blank	0.000	09:35
Control Test	0.050	09:35
Air Blank	0.000	09:36
Control Test	0.051	09:37
Air Blank	0.000	09:37
Control Test Stats		
Average	0.0507	
Std Dev	0.0006	
Rel Std Dev(%)	1.1395	

Operator's Signature

DELAND PD  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000  
 04/26/2016  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:38
Control Test	0.080	09:38
Air Blank	0.000	09:39
Control Test	0.081	09:39
Air Blank	0.000	09:40
Control Test	0.080	09:40
Air Blank	0.000	09:40
Control Test Stats		
Average	0.0803	
Std Dev	0.0006	
Rel Std Dev(%)	0.7187	

Operator's Signature

DELAND PD  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000  
 04/26/2016  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:41
Control Test	0.201	09:42
Air Blank	0.000	09:42
Control Test	0.200	09:43
Air Blank	0.000	09:44
Control Test	0.200	09:44
Air Blank	0.000	09:45
Control Test Stats		
Average	0.2003	
Std Dev	0.0006	
Rel Std Dev(%)	0.2882	

PAS  
 ASK

Operator's Signature