



Alcohol Testing Program

INSTRUMENT PROCESSING SHEET

Agency Daytona Beach PD s/N 80-005377
Date In 3/9/16 Date Out 3/26/16
Ship P/U H/D CMI EE

Intake Performed By TB
Registration
Annual
Return from CMI
Return from Enforcement Electronics
Other
Visual Inspection:
Case Handle
Dry Gas Holder Feet
Keyboard/Plug Back/Plugs
Screws tight Breath Hose
Other Equipment:
Power cord
Printer Cable
Other Static Bag

Quality Checks Performed By RB
Breath Tube Screen
Replace O-Rings
Instrument Set Up Verified
R-Value 190
Flow Verification (L/s)
Flow Column # A7P103
32mm 149 (.139 -.169)
36mm 167 (.156 -.190)
53mm 244 (.228 -.278)
103mm 506 (.447 -.547)
Barometric Pressure Check
Gauge ID # 28427
Stability Checks
Table with Simulator, Serial #, Lot #/Exp

Flow Calibration Performed By
Flow Calibration N/A
Flow Calibration Complete
Flow Column #
5L/min - 17mm
15L/min - 53mm
30L/min - 103mm
R-Value
Post Calibration Verification (L/s)
Flow Column #
32mm (.139 -.169)
36mm (.156 -.190)
53mm (.228 -.278)
103mm (.447 -.547)
Maintenance Performed By
Battery Replacement
Dry Gas Regulator Replacement
Breath Tube Replacement
Other
Suggested Service

RECEIVED
MAR 28 2016
FDLE
Alcohol Testing Program

Optical Bench Calibration Performed By
Optical Bench Calibration N/A
Optical Bench Calibration Complete
Barometric Pressure Gauge ID #
Table with Simulator, Serial Number, Lot Number, Expiration
Post Calibration Stability Checks
Table with Simulator, Serial Number, Lot Number, Expiration

Department Inspection Performed By RB
Barometric Pressure 1016 Gauge
ID# 28427 1017 Instrument
Mouth Alcohol Solution Lot # 2016-A
Acetone Stock Solution Lot # 2016-R
Table with Simulator, Serial Number
Attachments
Form 41
Pre-Stability Tests
Flow Calibration
Optical Bench Cal
Post-Stability Tests
Other Form 40

Notes:
QA/QC PLAN

Instrument Complies with Chapter 11D-8, FAC
Instrument Does Not Comply with Chapter 11D-8, FAC
Return to/Place into Evidentiary Use
Remain Out of Evidentiary Use
Conduct an Agency Inspection Before Evidentiary Use

Brett Highland

Quality Control Review

3/28/16

Date

80-005537  
 Stability Checks  
 3/26/16

13K

INTOXILYZER 8000  
 Instrument Initialization  
 12:04 03/26/2016

DRAN

DAYTONA BEACH PC  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-005537  
 03/26/2016  
 Software: 8100.27

DAYTONA BEACH PC  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-005537  
 03/26/2016  
 Software: 8100.27

DAYTONA BEACH PC  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-005537  
 03/26/2016  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:45
Control Test	0.196	13:46
Air Blank	0.000	13:46
Control Test	0.197	13:47
Air Blank	0.000	13:47
Control Test	0.197	13:48
Air Blank	0.000	13:49
Control Test Stats		
Average	0.1967	
Std Dev	0.0006	
Rel Std Dev(%)	0.2936	

Test	g/210L	Time
Air Blank	0.090	13:23
Control Test	0.079	13:24
Air Blank	0.000	13:24
Control Test	0.079	13:24
Air Blank	0.000	13:25
Control Test	0.079	13:25
Air Blank	0.090	13:26
Control Test Stats		
Average	0.0790	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

Test	g/210L	Time
Air Blank	0.000	13:02
Control Test	0.047	13:03
Air Blank	0.000	13:04
Control Test	0.048	13:04
Air Blank	0.000	13:05
Control Test	0.048	13:05
Air Blank	0.000	13:06
Control Test Stats		
Average	0.0477	
Std Dev	0.0005	
Rel Std Dev(%)	1.2112	

Test	g/210L	Time
Air Blank	0.000	12:40
Control Test	0.078	12:40
Air Blank	0.000	12:40
Control Test	0.078	12:41
Air Blank	0.000	12:41
Control Test	0.079	12:42
Air Blank	0.000	12:42
Control Test Stats		
Average	0.0783	
Std Dev	0.0005	
Rel Std Dev(%)	0.7370	

DGS

Operator's Signature

Operator's Signature

Operator's Signature

Operator's Signature