



Alcohol Testing Program

INSTRUMENT PROCESSING SHEET

Agency Homestead Police Department S/N 80-006199  
Date In 11/4/2016 Date Out 11/7/2016  Ship  P/U  H/D  CMI  EE

**Intake** Performed By DERR **Quality Checks** Performed By DERR **Flow Calibration** Performed By \_\_\_\_\_

Registration  
 Annual  
 Return from CMI  
 Return from Enforcement Electronics  
 Other \_\_\_\_\_

Visual Inspection:  
OK Case OK Handle  
OK Dry Gas Holder OK Feet  
OK Keyboard/Plug OK Back/Plugs  
OK Screws tight OK Breath Hose

Other Equipment:  
 Power cord  
 Printer Cable  
 Other Anti Static Bag

Notes: \_\_\_\_\_

Breath Tube Screen  
 Replace O-Rings  
 Instrument Set Up Verified  
 R-Value 198  
 Flow Verification (L/s)  
 Flow Column # ATP.101  
 32mm 156 (.139 - .169)  
 36mm 171 (.156 - .190)  
 53mm 238 (.228 - .278)  
 103mm 484 (.447 - .547)  
 Barometric Pressure Check  
 Gauge ID # 281623  
 Stability Checks

Simulator	Serial #	Lot #/Exp
0.05	<u>503967</u>	<u>201507A</u> <u>07/14/2017</u>
0.08	<u>503968</u>	<u>201601F</u> <u>01/26/2018</u>
0.20	<u>503969</u>	<u>201505A</u> <u>05/12/2017</u>
0.08 DGS	N/A	<u>A6600504</u> <u>01/05/2018</u>

Flow Calibration N/A  
 Flow Calibration Complete  
 Flow Column # \_\_\_\_\_  
 5L/min - 17mm  
 15L/min - 53mm  
 30L/min - 103mm  
 R-Value \_\_\_\_\_  
 Post Calibration Verification (L/s)  
 Flow Column # \_\_\_\_\_  
 32mm \_\_\_\_\_ (.139 - .169)  
 36mm \_\_\_\_\_ (.156 - .190)  
 53mm \_\_\_\_\_ (.228 - .278)  
 103mm \_\_\_\_\_ (.447 - .547)

**Maintenance** Performed By \_\_\_\_\_  
 Battery Replacement  
 Dry Gas Regulator Replacement  
 Breath Tube Replacement  
 Other \_\_\_\_\_

**Suggested Service** \_\_\_\_\_

RECEIVED  
NOV 16 2016  
FDLE  
Alcohol Testing Program

**Optical Bench Calibration** Performed By DERR

Optical Bench Calibration N/A  
 Optical Bench Calibration Complete  
 Barometric Pressure Gauge 1020 ID # 28199

Simulator	Serial Number	Lot Number	Expiration
0.000	<u>2235</u>	N/A	N/A
0.040	<u>2236</u>	<u>16101</u>	<u>02/02/2018</u>
0.100	<u>2237</u>	<u>15001</u>	<u>05/20/2017</u>
0.200	<u>2238</u>	<u>15104</u>	<u>05/27/2017</u>
0.400	<u>2239</u>	<u>16102</u>	<u>03/22/2018</u>
0.080 DGS	N/A	<u>03415080A</u>	<u>03/05/2017</u>

Post Calibration Stability Checks

Simulator	Serial Number	Lot Number	Expiration
0.05	<u>SD3967</u>	<u>201507A</u>	<u>07/14/2017</u>
0.08	<u>SD3968</u>	<u>201601F</u>	<u>01/26/2018</u>
0.20	<u>SD3969</u>	<u>201505A</u>	<u>05/12/2017</u>
0.08 DGS	N/A	<u>A6600504</u>	<u>01/05/2018</u>

**Department Inspection** Performed By DERR

Barometric Pressure 1019 Gauge  
 ID# 68639 1019 Instrument

Mouth Alcohol Solution Lot # 2016-A  
 Acetone Stock Solution Lot # 2016-B

Simulator	Serial Number
0.00	<u>SD3965</u>
Interferent	<u>SD3966</u>
0.05	<u>SD3967</u>
0.08	<u>SD3968</u>
0.20	<u>SD3969</u>

**Attachments**

Form 41  
 Pre-Stability Tests  
 Flow Calibration  
 Optical Bench Cal  
 Post-Stability Tests  
 Other \_\_\_\_\_

Notes:  E-MAILED 11/7/2016  APPROVED

QA/OC AKASPM 11/15/16

Brett K. Kurland

Instrument Complies with Chapter 11D-8, FAC  
 Instrument Does Not Comply with Chapter 11D-8, FAC  
 Return to/Place into Evidentiary Use  
 Remain Out of Evidentiary Use  
 Conduct an Agency Inspection Before Evidentiary Use

Quality Control Review

Date 11/16/16

<b>TYPE OF TEST</b>	<b>SERIAL NUMBER</b>	<b>AGENCY</b>	<b>DATE</b>	<b>PERFORMED BY</b>
Post Stabilities	80-006199	Homestead Police Department	11/07/2016	<i>[Signature]</i>

*[Handwritten initials]*

<b>0.05g/210L</b>	<b>0.08g/210L</b>	<b>0.20g/210L</b>	<b>DGS 0.08g/210L</b>
-------------------	-------------------	-------------------	-----------------------

<p>0.047 to 0.053 <input checked="" type="checkbox"/></p> <p>HOMESTEAD PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-006199 11/07/2016 Software: 8100.27</p> <p>Test g/210L Time</p> <p>Air Blank 0.000 10:33 Control Test 0.050 10:33 Air Blank 0.000 10:34 Control Test 0.050 10:35 Air Blank 0.000 10:35 Control Test 0.050 10:36 Air Blank 0.000 10:36 Control Test Stats Average 0.0500 Std Dev 0.0000 Rel Std Dev(%) 0.0000</p> <p>Operator's Signature <i>[Signature]</i></p>	<p>0.077 to 0.083 <input checked="" type="checkbox"/></p> <p>HOMESTEAD PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-006199 11/07/2016 Software: 8100.27</p> <p>Test g/210L Time</p> <p>Air Blank 0.000 10:37 Control Test 0.080 10:36 Air Blank 0.000 10:39 Control Test 0.079 10:39 Air Blank 0.000 10:40 Control Test 0.080 10:41 Air Blank 0.000 10:41 Control Test Stats Average 0.0797 Std Dev 0.0006 Rel Std Dev(%) 0.7247</p> <p>Operator's Signature <i>[Signature]</i></p>	<p>0.194 to 0.206 <input checked="" type="checkbox"/></p> <p>HOMESTEAD PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-006199 11/07/2016 Software: 8100.27</p> <p>Test g/210L Time</p> <p>Air Blank 0.000 10:42 Control Test 0.201 10:43 Air Blank 0.000 10:43 Control Test 0.201 10:44 Air Blank 0.000 10:45 Control Test 0.200 10:45 Air Blank 0.000 10:46 Control Test Stats Average 0.2007 Std Dev 0.0006 Rel Std Dev(%) 0.2877</p> <p>Operator's Signature <i>[Signature]</i></p>	<p>0.077 to 0.083 <input checked="" type="checkbox"/></p> <p>HOMESTEAD PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-006199 11/07/2016 Software: 8100.27</p> <p>Test g/210L Time</p> <p>Air Blank 0.000 10:48 Control Test 0.082 10:48 Air Blank 0.000 10:49 Control Test 0.081 10:49 Air Blank 0.000 10:49 Control Test 0.081 10:50 Air Blank 0.000 10:50 Control Test Stats Average 0.0813 Std Dev 0.0006 Rel Std Dev(%) 0.7099</p> <p>Operator's Signature <i>[Signature]</i></p>
---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

*[Handwritten initials]*

HOMESTEAD PD  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000  
 11/07/2016  
 SN 80-006199  
 09:46:05

Auto Calibration  
 Max Power Res Value = 88  
 Auto Range Res Value = 54

Soil Value = 0.009 g/210L \*\*\*  
 Fit value = 0.0000 mg/l %  
 Samples Taken = 4, Discarded = 1  
 Sum Io = 12497, Sum Io = 12757

Channel 1 >>>>  
 Sample % Abs (% Abs Ref)  
 Sample #1 = 0.1180 (-0.02210)  
 Sample #2 = 0.1060 (-0.0140)  
 Sample #3 = 0.0710 (-0.0410)  
 Sample #4 = 0.1110 (-0.0330)  
 Avg % Abs = 0.0960 (-0.0293)  
 STD DEV = 0.0218 (-0.0139)  
 REL STD DEV = 22.713 (47.279)

Soil Value = 0.040 g/210L \*\*\*  
 Fit value = 0.1905 mg/l %  
 Samples Taken = 4, Discarded = 1  
 Sum Io = 12489, Sum Io = 12757

Channel 2 >>>>  
 Sample % Abs (% Abs Ref)  
 Sample #1 = 1.5130 (-0.0020)  
 Sample #2 = 1.4830 (-0.0000)  
 Sample #3 = 1.5240 (-0.0040)  
 Sample #4 = 1.5130 (-0.0010)  
 Avg % Abs = 1.5067 (-0.0010)  
 STD DEV = 0.0212 (-0.0026)  
 REL STD DEV = 1.408 (264.575)

Soil Value = 0.100 g/210L \*\*\*  
 Fit value = 0.4762 mg/l %  
 Samples Taken = 4, Discarded = 1  
 Sum Io = 12483, Sum Io = 12754

Channel 1 >>>>  
 Sample % Abs (% Abs Ref)  
 Sample #1 = 1.8470 (-0.0210)  
 Sample #2 = 1.8530 (-0.0030)  
 Sample #3 = 1.8480 (-0.0020)  
 Sample #4 = 1.8910 (-0.0090)  
 Avg % Abs = 1.8607 (-0.0027)  
 STD DEV = 0.0178 (-0.0060)  
 REL STD DEV = 0.956 (226.039)

Channel 2 >>>>  
 Sample % Abs (% Abs Ref)  
 Sample #1 = 3.5560 (-0.0230)  
 Sample #2 = 3.5630 (-0.0270)  
 Sample #3 = 3.5080 (-0.0010)  
 Sample #4 = 3.5530 (-0.0030)  
 Avg % Abs = 3.5413 (-0.0077)  
 STD DEV = 0.0293 (-0.0168)  
 REL STD DEV = 0.827 (218.778)

Soil Value = 0.200 g/210L \*\*\*  
 Fit value = 0.9524 mg/l %  
 Samples Taken = 4, Discarded = 1  
 Sum Io = 12479, Sum Io = 12754

Channel 1 >>>>  
 Sample % Abs (% Abs Ref)  
 Sample #1 = 3.5250 (-0.0080)  
 Sample #2 = 3.5360 (-0.0290)  
 Sample #3 = 3.5510 (-0.0310)  
 Sample #4 = 3.5770 (-0.0370)  
 Avg % Abs = 3.5547 (-0.0323)  
 STD DEV = 0.0207 (-0.0042)  
 REL STD DEV = 0.594 (12.876)

Channel 2 >>>>  
 Sample % Abs (% Abs Ref)  
 Sample #1 = 6.7370 (-0.0100)  
 Sample #2 = 6.7440 (-0.0300)  
 Sample #3 = 6.7320 (-0.0570)  
 Sample #4 = 6.7500 (-0.0390)  
 Avg % Abs = 6.7420 (-0.0420)  
 STD DEV = 0.0092 (-0.0137)  
 REL STD DEV = 0.136 (32.733)

Soil Value = 0.400 g/210L \*\*\*  
 Fit value = 1.9048 mg/l %  
 Samples Taken = 4, Discarded = 1  
 Sum Io = 12475, Sum Io = 12751

Channel 1 >>>>  
 Sample % Abs (% Abs Ref)  
 Sample #1 = 6.7630 (-0.0060)  
 Sample #2 = 6.8080 (-0.0210)  
 Sample #3 = 6.7610 (-0.0280)  
 Sample #4 = 6.7960 (-0.0140)  
 Avg % Abs = 6.7883 (-0.0210)  
 STD DEV = 0.0244 (-0.0070)  
 REL STD DEV = 0.361 (33.333)

Channel 2 >>>>  
 Sample % Abs (% Abs Ref)  
 Sample #1 = 12.5000 (-0.0100)  
 Sample #2 = 12.5850 (-0.0120)  
 Sample #3 = 12.5130 (-0.0420)  
 Sample #4 = 12.5590 (-0.0200)  
 Avg % Abs = 12.5523 (-0.0247)  
 STD DEV = 0.0365 (-0.0155)  
 REL STD DEV = 0.290 (62.979)

Soil Value = 0.9524 mg/l %  
 Samples Taken = 4, Discarded = 1  
 Sum Io = 12479, Sum Io = 12754

Channel 1 >>>>  
 Sample % Abs (% Abs Ref)  
 Sample #1 = 3.5250 (-0.0080)  
 Sample #2 = 3.5360 (-0.0290)  
 Sample #3 = 3.5510 (-0.0310)  
 Sample #4 = 3.5770 (-0.0370)  
 Avg % Abs = 3.5547 (-0.0323)  
 STD DEV = 0.0207 (-0.0042)  
 REL STD DEV = 0.594 (12.876)

Optical Calibration	
SN: 80-006199	1
Agency: Homestead P.D.	
Date: 11/07/2016	
Quadratic Fit: +/-0.002g/210L	
By: <i>[Signature]</i>	

\*\*\*\*\* AUTO CRL DATA \*\*\*\*\*  
 Channel 1 >>>>  
 Soil Val = 0.0000 mg/l or 0.000 g/210L  
 % Abs = 0.096  
 Std Dev = 0.02 Rel Std Dev = 22.76  
 Soil Val = 0.1905 mg/l or 0.040 g/210L  
 % Abs = 0.758  
 Std Dev = 0.01 Rel Std Dev = 1.66  
 Soil Val = 0.4762 mg/l or 0.100 g/210L  
 % Abs = 1.861  
 Std Dev = 0.02 Rel Std Dev = 0.96  
 Soil Val = 0.9524 mg/l or 0.200 g/210L  
 % Abs = 3.555  
 Std Dev = 0.02 Rel Std Dev = 0.58  
 Soil Val = 1.9048 mg/l or 0.400 g/210L  
 % Abs = 6.788  
 Std Dev = 0.02 Rel Std Dev = 0.36  
 Zero Order Coef = -197.87  
 First Order Coef = 2828.35  
 Second Order Coef = 30.38  
 Standard Deviation = 56.896935

Channel 2 >>>>  
 Soil Val = 0.0000 mg/l or 0.000 g/210L  
 % Abs = 0.182  
 Std Dev = 0.02 Rel Std Dev = 8.79  
 Soil Val = 0.1905 mg/l or 0.040 g/210L  
 % Abs = 1.507  
 Std Dev = 0.02 Rel Std Dev = 1.41  
 Soil Val = 0.4762 mg/l or 0.100 g/210L  
 % Abs = 3.541  
 Std Dev = 0.03 Rel Std Dev = 0.83  
 Soil Val = 0.9524 mg/l or 0.200 g/210L  
 % Abs = 6.742  
 Std Dev = 0.01 Rel Std Dev = 0.14  
 Soil Val = 1.9048 mg/l or 0.400 g/210L  
 % Abs = 12.552  
 Std Dev = 0.04 Rel Std Dev = 0.29  
 Zero Order Coef = -207.43  
 First Order Coef = 1347.62  
 Second Order Coef = 14.86  
 Standard Deviation = 34.311481

Channel 1 >>>>  
 Soil Val = 0.0000 mg/l or 0.000 g/210L  
 % Abs = 0.000  
 Std Dev = 0.000  
 Soil Val = 0.0000 mg/l or 0.000 g/210L  
 % Abs = 0.000  
 Std Dev = 0.000  
 Soil Val = 0.0000 mg/l or 0.000 g/210L  
 % Abs = 0.000  
 Std Dev = 0.000  
 Soil Val = 0.0000 mg/l or 0.000 g/210L  
 % Abs = 0.000  
 Std Dev = 0.000

Solution Stats Quadratic Fit Chan 1:

Act	Fit	Residual
g/210L	g/210L	g/210L
0.000	0.001	-0.0008
0.040	0.039	0.0018
0.100	0.100	0.0002
0.200	0.201	-0.0006
0.400	0.400	0.0002

Soil Value = 0.080 g/210L \*\*\*  
 Fit value = 0.3810 mg/l %  
 Samples Taken = 4, Discarded = 1

Channel 1 >>>>  
 Sample #1 = 3294.00  
 Sample #2 = 3299.00  
 Sample #3 = 3223.00  
 Sample #4 = 3368.00  
 Average Result = 3296.6667  
 STD DEV = 72.5282  
 REL STD DEV = 2.200

Channel 2 >>>>  
 Sample #1 = 3340.00  
 Sample #2 = 3256.00  
 Sample #3 = 3316.00  
 Sample #4 = 3316.00  
 Average Result = 3309.3333  
 STD DEV = 11.5470  
 REL STD DEV = 0.349

\*\*\*\*\*  
 Dry Gas H2O Adjust Results \*\*\*\*\*  
 Barometric Pressure = 1020  
 3 um H2O Adjust (mg/l\*10.000) = 513  
 9 um H2O Adjust (mg/l\*10.000) = 500  
 \*\*\*\*\* AUTO CRL PASS \*\*\*\*\*

*[Signature]*  
 SK

TYPE OF TEST	SERIAL NUMBER	AGENCY	DATE	PERFORMED BY
Stabilities	80-006199	Homestead Police Department	11/07/2016	<i>[Signature]</i>

0.05g/210L <input checked="" type="checkbox"/>	0.08g/210L <input checked="" type="checkbox"/>	0.20g/210L <input checked="" type="checkbox"/>	DGS 0.08g/210L <input checked="" type="checkbox"/>
0.047 to 0.053 <input checked="" type="checkbox"/>	0.077 to 0.083 <input checked="" type="checkbox"/>	0.194 to 0.206 <input checked="" type="checkbox"/>	0.077 to 0.083 <input checked="" type="checkbox"/>

HOMESTEAD PD  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-006199  
11/07/2016  
Software: 8100.27

Test g/210L Time

Air Blank 0.000 08:19  
Control Test 0.049 08:19  
Air Blank 0.000 08:20  
Control Test 0.049 08:21  
Air Blank 0.000 08:21  
Control Test 0.048 08:22  
Air Blank 0.000 08:22

Control Test Stats  
Average 0.0487  
Std Dev 0.0006  
Rel Std Dev(%) 1.1853

HOMESTEAD PD  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-006199  
11/07/2016  
Software: 8100.27

Test g/210L Time

Air Blank 0.000 08:23  
Control Test 0.080 08:24  
Air Blank 0.000 08:25  
Control Test 0.079 08:25  
Air Blank 0.000 08:26  
Control Test 0.078 08:26  
Air Blank 0.000 08:27

Control Test Stats  
Average 0.0787  
Std Dev 0.0012  
Rel Std Dev(%) 1.4678

HOMESTEAD PD  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-006199  
11/07/2016  
Software: 8100.27

Test g/210L Time

Air Blank 0.000 08:28  
Control Test 0.198 08:29  
Air Blank 0.000 08:29  
Control Test 0.198 08:30  
Air Blank 0.000 08:30  
Control Test 0.198 08:31  
Air Blank 0.000 08:32

Control Test Stats  
Average 0.1986  
Std Dev 0.0008  
Rel Std Dev(%) 0.0000

HOMESTEAD PD  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-006199  
11/07/2016  
Software: 8100.27

Test g/210L Time

Air Blank 0.000 08:33  
Control Test 0.082 08:33  
Air Blank 0.000 08:34  
Control Test 0.084 08:34  
Air Blank 0.000 08:34  
Control Test 0.084 08:35  
Air Blank 0.000 08:35

Control Test Stats  
Average 0.0833  
Std Dev 0.0012  
Rel Std Dev(%) 1.3855

*DEER*  
Operator's Signature

*[Signature]*  
Operator's Signature

*[Signature]*  
Operator's Signature

*[Signature]*  
Operator's Signature

*PK*  
*CPM*