

INSTRUMENT PROCESSING SHEET

Agency Maratee County SO S/N 80-006236
 Date In 7/12/16 Date Out 7/19/16 Ship P/U H/D CMI EE

Intake Performed By DP

Registration
 Annual
 Return from CMI
 Return from Enforcement Electronics
 Other _____

Visual Inspection:
OK Case OK Handle
OK Dry Gas Holder OK Feet
OK Keyboard/Plug OK Back/Plugs
OK Screws tight OK Breath Hose

Other Equipment:
 Power cord
 Printer Cable
 Other Static Bag

Notes: _____

Quality Checks Performed By PWS

Breath Tube Screen
 Replace O-Rings
 Instrument Set Up Verified
 R-Value 241
 Flow Verification (L/s)
 Flow Column # ATP102
 32mm 160 (.139 - .169)
 36mm 179 (.156 - .190)
 53mm 246 (.228 - .278)
 103mm 500 (.447 - .547)

Barometric Pressure Check
 Gauge ID # 28427

Stability Checks

Simulator	Serial #	Lot #/Exp
0.05	SD1018	201507A 7/14/17
0.08	SD1011	201601F 1/26/18
0.20	SD1025	201505A 5/12/17
0.08 DGS	N/A	AG431502 11/11/16

Flow Calibration Performed By _____

Flow Calibration N/A
 Flow Calibration Complete
 Flow Column # _____
 5L/min - 47mm
 15L/min - 53mm
 30L/min - 103mm

R-Value _____
 Post Calibration Verification (L/s)
 Flow Column # _____
 32mm _____ (.139 - .169)
 36mm _____ (.156 - .190)
 53mm _____ (.228 - .278)
 103mm _____ (.447 - .547)

Maintenance Performed By _____

Battery Replacement
 Dry Gas Regulator Replacement
 Breath Tube Replacement
 Other _____

Suggested Service _____

RECEIVED
 JUL 21 2016
 FDLE
 Alcohol Testing Program

Optical Bench Calibration Performed By PWS

Optical Bench Calibration N/A
 Optical Bench Calibration Complete
 Barometric Pressure Gauge 1020 ID# 26932

Simulator	Serial Number	Lot Number	Expiration
0.000	DR1275	N/A	N/A
0.040	G2882	16101	2/2/18
0.100	G2078	16001	5/8/18
0.200	G2408	15104	5/27/17
0.400	G5358	16102	3/22/18
0.080 DGS	N/A	03415080A1	3/5/17

Post Calibration Stability Checks

Simulator	Serial Number	Lot Number	Expiration
0.05	SD1018	201507A	7/14/17
0.08	SD1011	201601F	1/26/18
0.20	SD1025	201505A	5/12/17
0.08 DGS	N/A	AG431502	11/11/16

Department Inspection Performed By PWS

Barometric Pressure 1020 Gauge
 ID# 28427 1018 Instrument

Mouth Alcohol Solution Lot # 2015-A
 Acetone Stock Solution Lot # 2016-B

Simulator	Serial Number
0.00	SD1019
Interferent	SD1021
0.05	SD1018
0.08	SD1011
0.20	SD1025

Attachments

Form 41
 Pre-Stability Tests
 Flow Calibration

Optical Bench Cal
 Post-Stability Tests
 Other _____

Notes: All stability checks within tolerance. Calibrated optical bench to bring values closer to nominal. (PWS)

DA/DC OK WSM

Brett Kirkland
 Quality Control Review

Instrument Complies with Chapter 11D-8, FAC
 Instrument Does Not Comply with Chapter 11D-8, FAC

Return to/Place into Evidentiary Use
 Remain Out of Evidentiary Use

Conduct an Agency Inspection Before Evidentiary Use

7/21/16
 Date

Stability Tests
Pre-Calibration

- Manatee CSO

80-006236

7/14/16

MANATEE COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006236
07/19/2016
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:57
Control Test	0.047	09:58
Air Blank	0.000	09:58
Control Test	0.047	09:59
Air Blank	0.000	09:59
Control Test	0.047	10:00
Air Blank	0.000	10:01
Control Test Stats		
Average	0.0470	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

QMS

QMS

Operator's Signature

MANATEE COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006236
07/19/2016
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:11
Control Test	0.078	10:12
Air Blank	0.000	10:12
Control Test	0.078	10:13
Air Blank	0.000	10:14
Control Test	0.078	10:14
Air Blank	0.000	10:15
Control Test Stats		
Average	0.0780	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

QMS

Operator's Signature

MANATEE COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006236
07/19/2016
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:02
Control Test	0.199	10:03
Air Blank	0.000	10:04
Control Test	0.198	10:04
Air Blank	0.000	10:05
Control Test	0.199	10:05
Air Blank	0.000	10:06
Control Test Stats		
Average	0.1987	
Std Dev	0.0006	
Rel Std Dev(%)	0.2906	

QMS

Operator's Signature

MANATEE COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006236
07/19/2016
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:07
Control Test	0.080	10:08
Air Blank	0.000	10:08
Control Test	0.080	10:09
Air Blank	0.000	10:09
Control Test	0.081	10:09
Air Blank	0.000	10:10
Control Test Stats		
Average	0.0803	
Std Dev	0.0006	
Rel Std Dev(%)	0.7187	

QMS

BSK

QMS

Operator's Signature

Optical Bench
 Calibration
 Manatee CSO
 # 30-006234
 7/19/14

Sol Value = 0.040 g/210L ***
 Fit value = 0.1905 mg/l %%%
 Samples Taken = 4, Discarded = 1
 Sum Io = 12513, Sum Io = 12894
 <<<<< CHANNEL 1 >>>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 0.7190 (0.0000)
 Sample #2 = 0.7610 (0.0060)
 Sample #3 = 0.7390 (0.0440)
 Sample #4 = 0.7760 (0.0310)
 Avg % Abs = 0.7587 (0.0270)
 STD DEV = 0.0186 (0.0193)
 REL STD DEV = 2.453 (71.530)

<<<<< CHANNEL 2 >>>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 1.4710 (0.0020)
 Sample #2 = 1.5020 (-0.0960)
 Sample #3 = 1.4880 (0.0100)
 Sample #4 = 1.5020 (0.0170)
 Avg % Abs = 1.4973 (0.0070)
 STD DEV = 0.0081 (0.0118)
 REL STD DEV = 0.540 (168.426)

Sol Value = 0.100 g/210L ***
 Fit value = 0.4762 mg/l %%%
 Samples Taken = 4, Discarded = 1
 Sum Io = 12595, Sum Io = 12890
 <<<<< CHANNEL 1 >>>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 1.7140 (0.0050)
 Sample #2 = 1.7540 (0.0130)
 Sample #3 = 1.7390 (0.0440)
 Sample #4 = 1.7690 (0.0300)
 Avg % Abs = 1.7540 (0.0290)
 STD DEV = 0.0150 (0.0155)
 REL STD DEV = 0.855 (53.532)

<<<<< CHANNEL 2 >>>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 3.4500 (0.0100)
 Sample #2 = 3.4990 (-0.0010)
 Sample #3 = 3.4590 (0.0290)
 Sample #4 = 3.5090 (0.0070)
 Avg % Abs = 3.4923 (0.0117)
 STD DEV = 0.0208 (0.0155)
 REL STD DEV = 0.596 (133.156)

Sol Value = 0.200 g/210L ***
 Fit value = 0.9524 mg/l %%%
 Samples Taken = 4, Discarded = 1
 Sum Io = 12501, Sum Io = 12889
 <<<<< CHANNEL 1 >>>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 3.4020 (-0.0030)
 Sample #2 = 3.3580 (0.0540)
 Sample #3 = 3.3800 (0.0590)
 Sample #4 = 3.4200 (0.0370)
 Avg % Abs = 3.3860 (0.0580)
 STD DEV = 0.0314 (0.0115)
 REL STD DEV = 0.928 (23.065)

<<<<< CHANNEL 2 >>>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 6.7300 (-0.0100)
 Sample #2 = 6.7110 (0.0560)
 Sample #3 = 6.7490 (0.0520)
 Sample #4 = 6.7820 (0.0450)
 Avg % Abs = 6.7473 (0.0510)
 STD DEV = 0.0355 (0.0056)
 REL STD DEV = 0.527 (10.917)

Sol Value = 0.400 g/210L ***
 Fit value = 1.9048 mg/l %%%
 Samples Taken = 4, Discarded = 1
 Sum Io = 12497, Sum Io = 12888
 <<<<< CHANNEL 1 >>>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 6.3600 (-0.0110)
 Sample #2 = 6.3730 (0.0570)
 Sample #3 = 6.3720 (0.0690)
 Sample #4 = 6.3730 (0.0640)
 Avg % Abs = 6.3727 (0.0700)
 STD DEV = 0.0006 (0.0135)
 REL STD DEV = 0.009 (19.325)

<<<<< CHANNEL 2 >>>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 12.4110 (0.0030)
 Sample #2 = 12.3930 (0.1070)
 Sample #3 = 12.4230 (0.1100)
 Sample #4 = 12.4200 (0.1260)
 Avg % Abs = 12.4120 (0.1143)
 STD DEV = 0.0165 (0.0102)
 REL STD DEV = 0.133 (8.934)

**** AUTO CAL DATA ****
 <<<<< CHANNEL 1 >>>>>
 Sol Val = 0.0000 mg/l or 0.000 g/210L
 % Abs = 0.099
 Std Dev = 0.02 Rel Std Dev = 21.88
 Sol Val = 0.1905 mg/l or 0.040 g/210L
 % Abs = 0.759
 Std Dev = 0.02 Rel Std Dev = 2.45
 Sol Val = 0.4762 mg/l or 0.100 g/210L
 % Abs = 1.754
 Std Dev = 0.02 Rel Std Dev = 0.86
 Sol Val = 0.9524 mg/l or 0.200 g/210L
 % Abs = 3.386
 Std Dev = 0.03 Rel Std Dev = 0.93
 Sol Val = 1.9048 mg/l or 0.400 g/210L
 % Abs = 6.373
 Std Dev = 0.00 Rel Std Dev = 0.01
 Zero Order Coef = -236.28
 First Order Coef = 2749.49
 Second Order Coef = 43.17
 Standard Deviation = 39.393608

<<<<< CHANNEL 2 >>>>>
 Sol Val = 0.0000 mg/l or 0.000 g/210L
 % Abs = 0.168
 Std Dev = 0.01 Rel Std Dev = 7.94
 Sol Val = 0.1905 mg/l or 0.040 g/210L
 % Abs = 1.497
 Std Dev = 0.01 Rel Std Dev = 0.54
 Sol Val = 0.4762 mg/l or 0.100 g/210L
 % Abs = 3.492
 Std Dev = 0.02 Rel Std Dev = 0.60
 Sol Val = 0.9524 mg/l or 0.200 g/210L
 % Abs = 6.747
 Std Dev = 0.04 Rel Std Dev = 0.53
 Sol Val = 1.9048 mg/l or 0.400 g/210L
 % Abs = 12.412
 Std Dev = 0.02 Rel Std Dev = 0.13
 Zero Order Coef = -166.66
 First Order Coef = 1328.85
 Second Order Coef = 17.55
 Standard Deviation = 63.841362

Solution Stats Quadratic Fit Chan 1

Act	Fit	Residual
g/210L	g/210L	g/210L
0.000	0.001	-0.0008
0.040	0.039	0.0006
0.100	0.099	0.0009
0.200	0.201	-0.0009
0.400	0.400	0.0002

Solution Stats Quadratic Fit Chan 2

Act	Fit	Residual
g/210L	g/210L	g/210L
0.000	0.001	-0.0012
0.040	0.039	0.0009
0.100	0.098	0.0015
0.200	0.202	-0.0016
0.400	0.400	0.0004

Sol Value = 0.080 g/210L ***
 Fit value = 0.3810 mg/l %%%
 Samples Taken = 4, Discarded = 1
 ***** CHANNEL 1 *****
 Sample #1 = 2996.00
 Sample #2 = 3040.00
 Sample #3 = 2906.00
 Sample #4 = 2994.00
 Average Result = 2980.0000
 STD DEV = 68.0882
 REL STD DEV = 2.285
 ***** CHANNEL 2 *****
 Sample #1 = 3177.00
 Sample #2 = 3213.00
 Sample #3 = 3179.00
 Sample #4 = 3216.00
 Average Result = 3202.6667
 STD DEV = 20.5508
 REL STD DEV = 0.642

 Dry Gas H2O Adjust Results *****
 Barometric Pressure = 1019
 3 um H2O Adjust (mg/l*10,000) = 829
 9 um H2O Adjust (mg/l*10,000) = 607
 ***** AUTO CAL PASS *****

ASK

Stability Tests - Maratee CSO # 80-006236 7/19/16
 Post-Calibration

ABK

MANATEE COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-006236
 07/19/2016
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:14
Control Test	0.048	13:15
Air Blank	0.000	13:15
Control Test	0.050	13:16
Air Blank	0.000	13:17
Control Test	0.050	13:17
Air Blank	0.000	13:18
Control Test Stats		
Average	0.0493	
Std Dev	0.0012	
Rel Std Dev(%)	2.3406	

BM

[Signature]

Operator's Signature

MANATEE COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-006236
 07/19/2016
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:19
Control Test	0.081	13:20
Air Blank	0.000	13:20
Control Test	0.081	13:21
Air Blank	0.000	13:22
Control Test	0.080	13:22
Air Blank	0.000	13:23
Control Test Stats		
Average	0.0807	
Std Dev	0.0006	
Rel Std Dev(%)	0.7157	

[Signature]

Operator's Signature

MANATEE COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-006236
 07/19/2016
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:24
Control Test	0.201	13:25
Air Blank	0.000	13:25
Control Test	0.202	13:26
Air Blank	0.000	13:27
Control Test	0.202	13:27
Air Blank	0.000	13:28
Control Test Stats		
Average	0.2017	
Std Dev	0.0006	
Rel Std Dev(%)	0.2863	

[Signature]

Operator's Signature

MANATEE COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-006236
 07/19/2016
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:29
Control Test	0.081	13:30
Air Blank	0.000	13:30
Control Test	0.081	13:31
Air Blank	0.000	13:31
Control Test	0.081	13:32
Air Blank	0.000	13:32
Control Test Stats		
Average	0.0810	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

DS

[Signature]

Operator's Signature