



Alcohol Testing Program

INSTRUMENT PROCESSING SHEET

Agency Clay County SC

S/N 80-006237

Date In 2/1/16

Date Out 2/4/2016

Ship P/U H/D CMI EE

Intake Performed By [Signature]

Registration
 Annual
 Return from CMI
 Return from Enforcement Electronics
 Other _____

Visual Inspection:
OK Case OK Handle
OK Dry Gas Holder OK Feet
OK Keyboard/Plug OK Back/Plugs
OK Screws tight OK Breath Hose

Other Equipment:
 Power cord
 Printer Cable
 Other Steric Bag

Notes: _____

Quality Checks Performed By gg

Breath Tube Screen
 Replace O-Rings
 Instrument Set Up Verified
 R-Value 239
 Flow Verification (L/s)
 Flow Column # APP103
 32mm .160 (.139 - .169)
 36mm .179 (.156 - .190)
 53mm .242 (.228 - .278)
 103mm .500 (.447 - .547)

Barometric Pressure Check
 Gauge ID # 28427

Stability Checks

| Simulator | Serial # | Lot #/Exp |
|-----------|----------|------------------------|
| 0.05 | SD1018 | 201507A 07/14/2017 |
| 0.08 | SD1011 | 201502G 02/21/2017 |
| 0.20 | SD1025 | 201505A 05/12/2017 |
| 0.08 DGS | N/A | AG511701 01/27/2017 |

Flow Calibration Performed By _____

Flow Calibration N/A
 Flow Calibration Complete
 Flow Column # _____
 5L/min - 17mm
 15L/min - 4/53mm
 30L/min - 103mm

R-Value _____
 Post Calibration Verification (L/s)
 Flow Column # _____
 32mm _____ (.139 - .169)
 36mm _____ (.156 - .190)
 53mm _____ (.228 - .278)
 103mm _____ (.447 - .547)

Maintenance Performed By _____

Battery Replacement
 Dry Gas Regulator Replacement
 Breath Tube Replacement
 Other _____

Suggested Service

RECEIVED
FEB 05 2016
FDLE
Alcohol Testing Program

Optical Bench Calibration Performed By _____

Optical Bench Calibration N/A
 Optical Bench Calibration Complete

Barometric Pressure Gauge ID # _____

| Simulator | Serial Number | Lot Number | Expiration |
|-----------|---------------|------------|------------|
| 0.000 | | N/A | N/A |
| 0.040 | | | |
| 0.100 | | | |
| 0.200 | | | |
| 0.400 | | | |
| 0.080 DGS | N/A | | |

Post Calibration Stability Checks

| Simulator | Serial Number | Lot Number | Expiration |
|-----------|---------------|------------|------------|
| 0.05 | | | |
| 0.08 | | | |
| 0.20 | | | |
| 0.08 DGS | N/A | | |

Department Inspection Performed By gg

Barometric Pressure Gauge ID# 28427 1018
 Instrument ID# 1019

Mouth Alcohol Solution Lot # 2015-A
 Acetone Stock Solution Lot # 2015-B

| Simulator | Serial Number |
|-------------|---------------|
| 0.00 | SD1022 |
| Interferent | SD1021 |
| 0.05 | SD1018 |
| 0.08 | SD1011 |
| 0.20 | SD1025 |

Attachments

Form 41
 Pre-Stability Tests
 Flow Calibration
 Optical Bench Cal
 Post-Stability Tests
 Other _____

Notes: QC-BK

Instrument Complies with Chapter 11D-8, FAC
 Instrument Does Not Comply with Chapter 11D-8, FAC

Return to/Place into Evidentiary Use
 Remain Out of Evidentiary Use

Conduct an Agency Inspection Before Evidentiary Use

2/5/16
Date

Patrick Murphy
Quality Control Review

Stability Checks- Instrument #80-006237 Clay County SO 2/4/2016 29

POAM

CLAY COUNTY SO
Intoxilizer - Alcohol Analyzer
Model 8000 SN 80-006237
02/04/2016
Software: 8100.27

| Test | 9/210L | Time |
|--------------------|--------|-------|
| Air Blank | 0.000 | 10:29 |
| Control Test | 0.048 | 10:30 |
| Air Blank | 0.000 | 10:30 |
| Control Test | 0.049 | 10:31 |
| Air Blank | 0.000 | 10:31 |
| Control Test | 0.049 | 10:32 |
| Air Blank | 0.000 | 10:33 |
| Control Test Stats | | |
| Average | 0.0487 | |
| Std Dev | 0.0006 | |
| Rel Std Dev(%) | 1.1863 | |

Operator's Signature

CLAY COUNTY SO
Intoxilizer - Alcohol Analyzer
Model 8000 SN 80-006237
02/04/2016
Software: 8100.27

| Test | 9/210L | Time |
|--------------------|--------|-------|
| Air Blank | 0.000 | 10:34 |
| Control Test | 0.079 | 10:34 |
| Air Blank | 0.000 | 10:35 |
| Control Test | 0.079 | 10:36 |
| Air Blank | 0.000 | 10:36 |
| Control Test | 0.078 | 10:37 |
| Air Blank | 0.000 | 10:38 |
| Control Test Stats | | |
| Average | 0.0787 | |
| Std Dev | 0.0006 | |
| Rel Std Dev(%) | 0.7339 | |

Operator's Signature

CLAY COUNTY SO
Intoxilizer - Alcohol Analyzer
Model 8000 SN 80-006237
02/04/2016
Software: 8100.27

| Test | 9/210L | Time |
|--------------------|--------|-------|
| Air Blank | 0.000 | 10:39 |
| Control Test | 0.205 | 10:40 |
| Air Blank | 0.000 | 10:40 |
| Control Test | 0.204 | 10:41 |
| Air Blank | 0.000 | 10:41 |
| Control Test | 0.204 | 10:42 |
| Air Blank | 0.000 | 10:43 |
| Control Test Stats | | |
| Average | 0.2043 | |
| Std Dev | 0.0006 | |
| Rel Std Dev(%) | 0.2826 | |

Operator's Signature

CLAY COUNTY SO
Intoxilizer - Alcohol Analyzer
Model 8000 SN 80-006237
02/04/2016
Software: 8100.27

| Test | 9/210L | Time |
|--------------------|--------|-------|
| Air Blank | 0.000 | 10:44 |
| Control Test | 0.080 | 10:44 |
| Air Blank | 0.000 | 10:45 |
| Control Test | 0.080 | 10:45 |
| Air Blank | 0.000 | 10:46 |
| Control Test | 0.080 | 10:46 |
| Air Blank | 0.000 | 10:47 |
| Control Test Stats | | |
| Average | 0.0800 | |
| Std Dev | 0.0000 | |
| Rel Std Dev(%) | 0.0000 | |

Operator's Signature

BK

D65