

INSTRUMENT PROCESSING SHEET

Agency Kennedy Space Center S/N 80-006621
 Date In 3/18/16 Date Out 3/22/16 Ship P/U H/D CMI EE

Intake Performed By [Signature]

Registration
 Annual
 Return from CMI
 Return from Enforcement Electronics
 Other _____

Visual Inspection:
 Case Handle
 Dry Gas Holder Feet
 Keyboard/Plug Back/Plugs
 Screws tight Breath Hose

Other Equipment:
 Power cord
 Printer Cable
 Other Static Bag

Notes: Two AC power cords, one DC power cord

Quality Checks Performed By [Signature]

Breath Tube Screen
 Replace O-Rings
 Instrument Set Up Verified
 R-Value 218
 Flow Verification (L/s)
 Flow Column # 97P103
 32mm 150 (.139 - .169)
 36mm 174 (.156 - .190)
 53mm 249 (.228 - .278)
 103mm 511 (.447 - .547)

Barometric Pressure Check
 Gauge ID # 28427

Stability Checks

Simulator	Serial #	Lot #/Exp
0.05	SD1018	201507A 7/14/17
0.08	SD1011	201601F 1/26/18
0.20	G4444	201505A 5/12/17
0.08 DGS	N/A	AG6005UP 1/9/15

Flow Calibration Performed By _____

Flow Calibration N/A
 Flow Calibration Complete
 Flow Column # _____
 5L/min - 17mm
 15L/min - 53mm
 30L/min - 103mm
 R-Value _____
 Post Calibration Verification (L/s)
 Flow Column # _____
 32mm _____ (.139 - .169)
 36mm _____ (.156 - .190)
 53mm _____ (.228 - .278)
 103mm _____ (.447 - .547)

Maintenance Performed By _____

Battery Replacement
 Dry Gas Regulator Replacement
 Breath Tube Replacement
 Other _____

Suggested Service

RECEIVED
 MAR 28 2016
 FDLE
 Alcohol Testing Program

Optical Bench Calibration Performed By _____

Optical Bench Calibration N/A
 Optical Bench Calibration Complete

Barometric Pressure Gauge ID # _____

Simulator	Serial Number	Lot Number	Expiration
0.000		N/A	N/A
0.040			
0.100			
0.200			
0.400			
0.080 DGS	N/A		

Post Calibration Stability Checks

Simulator	Serial Number	Lot Number	Expiration
0.05			
0.08			
0.20			
0.08 DGS	N/A		

Department Inspection Performed By [Signature]

Barometric Pressure 1016 Gauge
 ID# 28427 1017 Instrument

Mouth Alcohol Solution Lot # 2016-A
 Acetone Stock Solution Lot # 2016-B

Simulator	Serial Number
0.00	SD1022
Interferent	SD1021
0.05	SD1018
0.08	SD1011
0.20	G4444

Attachments

Form 41
 Pre-Stability Tests
 Flow Calibration

Optical Bench Cal
 Post-Stability Tests
 Other _____

Notes: QA/DC GPM

Instrument Complies with Chapter 11D-8, FAC
 Instrument Does Not Comply with Chapter 11D-8, FAC

Return to/Place into Evidentiary Use
 Remain Out of Evidentiary Use

Conduct an Agency Inspection Before Evidentiary Use

Brett Hokland
 Quality Control Review

3/28/16
 Date

80-006621
 Steady Clear
 3/26/16

INTOXILYZER 8000
 Instrument Initialization
 11:05 03/26/2016

BGM

BK

CMI INC.
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-006621
 03/26/2016
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:24
Control Test	0.078	12:24
Air Blank	0.000	12:25
Control Test	0.078	12:25
Air Blank	0.000	12:25
Control Test	0.078	12:25
Air Blank	0.000	12:26
Control Test	0.078	12:26
Air Blank	0.000	12:26
Control Test	0.000	12:26
Average	0.0780	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

CMI INC.
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-006621
 03/26/2016
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:24
Control Test	0.078	12:24
Air Blank	0.000	12:25
Control Test	0.078	12:25
Air Blank	0.000	12:25
Control Test	0.078	12:25
Air Blank	0.000	12:26
Control Test	0.078	12:26
Air Blank	0.000	12:26
Control Test	0.000	12:26
Average	0.0780	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

CMI INC.
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-006621
 03/26/2016
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:07
Control Test	0.047	12:08
Air Blank	0.000	12:09
Control Test	0.048	12:09
Air Blank	0.000	12:10
Control Test	0.048	12:10
Air Blank	0.000	12:11
Control Test	0.000	12:11
Average	0.0477	
Std Dev	0.0006	
Rel Std Dev(%)	1.2112	

CMI INC.
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-006621
 03/26/2016
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:40
Control Test	0.077	11:40
Air Blank	0.000	11:40
Control Test	0.077	11:41
Air Blank	0.000	11:41
Control Test	0.078	11:42
Air Blank	0.000	11:42
Control Test	0.000	11:42
Average	0.0773	
Std Dev	0.0006	
Rel Std Dev(%)	0.7466	

DS

Operator's Signature

Operator's Signature

Operator's Signature

Operator's Signature

Kennedy Space Center.

Showing accurate
agency names

INTOXILYZER 8000
Instrument Initialization
07:38 04/26/2016

BLK

KENNEDY SPACE CENTER
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006621
04/26/2016
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:04
Control Test	0.047	08:05
Air Blank	0.000	08:05
Control Test	0.048	08:06
Air Blank	0.000	08:06
Control Test	0.048	08:07
Air Blank	0.000	08:08
Control Test Stats		
Average	0.0477	
Std Dev	0.0006	
Rel Std Dev(%)	1.2112	

Operator's Signature