



INSTRUMENT PROCESSING SHEET

Agency FHP Tampa S/N 80-006632
Date In 7/18/16 Date Out 11/30/16 Ship P/U H/D CMI EE

Intake Performed By [Signature]

Registration
 Annual
 Return from CMI
 Return from Enforcement Electronics
 Other _____

Visual Inspection:
 Case Handle
 Dry Gas Holder Feet
 Keyboard/Plug Back/Plugs
 Screws tight Breath Hose

Other Equipment:
 Power cord
 Printer Cable
 Other Static Bag

Notes: _____

Quality Checks Performed By [Signature]

Breath Tube Screen
 Replace O-Rings
 Instrument Set Up Verified
 R-Value 207
 Flow Verification (L/s)
 Flow Column # ATP102
 32mm 0.164 (.139 - .169)
 36mm 0.179 (.156 - .190)
 53mm 0.250 (.228 - .278)
 103mm 0.500 (.447 - .547)

Barometric Pressure Check
 Gauge ID # 28427

Stability Checks

Simulator	Serial #	Lot #/Exp
0.05	SD1018	201507A 7/14/17
0.08	SD1011	201601F 1/26/18
0.20	SD1025	201604C 4/5/18
0.08 DGS	N/A	AG612405 5/3/18

Flow Calibration Performed By [Signature]

Flow Calibration N/A
 Flow Calibration Complete
 Flow Column # _____
 5L/min - 100mm
 15L/min - 53mm
 30L/min - 103mm

R-Value _____
 Post Calibration Verification (L/s)
 Flow Column # _____
 32mm _____ (.139 - .169)
 36mm _____ (.156 - .190)
 53mm _____ (.228 - .278)
 103mm _____ (.447 - .547)

Maintenance Performed By [Signature]

Battery Replacement
 Dry Gas Regulator Replacement
 Breath Tube Replacement
 Other _____

Suggested Service

RECEIVED
DEC 01 2016
FDLE
Alcohol Testing Program

Optical Bench Calibration Performed By [Signature]

Optical Bench Calibration N/A
 Optical Bench Calibration Complete
 Barometric Pressure Gauge 1012 ID # 20932 x2 (Both Sept Cal)

Simulator	Serial Number	Lot Number	Expiration
0.000	SD1016	N/A	N/A
0.040	SD1024	16101	2/2/18
0.100	DR1279	16001	5/8/18
0.200	DR3856	16103	6/14/18
0.400	SD1013	16102	3/22/18
0.080 DGS	N/A	15615080A2	7/5/17

Post Calibration Stability Checks

Simulator	Serial Number	Lot Number	Expiration
0.05	SD1018	201507A	7/14/17
0.08	SD1011	201601F	1/26/18
0.20	SD1025	201604C	4/5/18
0.08 DGS	N/A	AG612405	5/3/18

Department Inspection Performed By [Signature]

Barometric Pressure ID# 28427 Gauge 1011 Instrument 1009

Mouth Alcohol Solution Lot # 2015-A
 Acetone Stock Solution Lot # 2016-B

Simulator	Serial Number
0.00	SD1019
Interferent	SD1021
0.05	DR2035
0.08	SD1011
0.20	SD1025

Attachments

Form 41
 Pre-Stability Tests
 Flow Calibration
 Optical Bench Cal x2 x 3
 Post-Stability Tests x2
 Other Form 40 x2
 Diagnostic Check

Notes: Quality checks performed after regulator replacement @WS
Performed optical bench calibrations to bring values closer to nominal @S
See IPS pg 2 for 3rd calibration documentation @S
PA/OC OK SP 11/30
Brett Kirkland

Instrument Complies with Chapter 11D-8, FAC
 Instrument Does Not Comply with Chapter 11D-8, FAC
 Return to/Place into Evidentiary Use
 Remain Out of Evidentiary Use
 Conduct an Agency Inspection Before Evidentiary Use

Date 12/1/16

Quality Control Review

Date

Agency FL Highway Patrol

S/N 80-006632

Date In _____

Date Out _____

Ship P/U H/D CMI EE

<p>Intake Performed By _____</p> <p><input checked="" type="checkbox"/> Registration <input type="checkbox"/> Annual <input type="checkbox"/> Return from CMI <input type="checkbox"/> Return from Enforcement <input type="checkbox"/> Electronics <input type="checkbox"/> Other _____</p> <p>Visual Inspection: _____ Case _____ Handle _____ Dry Gas Holder _____ Feet _____ Keyboard/Plug _____ Back/Plugs _____ Screws tight _____ Breath Hose</p> <p>Other Equipment: <input type="checkbox"/> Power cord <input type="checkbox"/> Printer Cable <input type="checkbox"/> Other _____</p> <p>Notes: _____ _____ _____</p>	<p>Quality Checks Performed By _____</p> <p><input type="checkbox"/> Breath Tube Screen <input type="checkbox"/> Replace O-Rings <input type="checkbox"/> Instrument Set Up Verified <input type="checkbox"/> R-Value _____ <input type="checkbox"/> Flow Verification (L/s) Flow Column # _____ 32mm _____ (.139 - .169) 36mm _____ (.156 - .190) 53mm _____ (.228 - .278) 103mm _____ (.447 - .547)</p> <p><input type="checkbox"/> Barometric Pressure Check Gauge ID # _____</p> <p><input type="checkbox"/> Stability Checks</p> <table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td>0.05</td> <td></td> <td></td> </tr> <tr> <td>0.08</td> <td></td> <td></td> </tr> <tr> <td>0.20</td> <td></td> <td></td> </tr> <tr> <td>0.08 DGS</td> <td>N/A</td> <td></td> </tr> </tbody> </table>	Simulator	Serial #	Lot #/Exp	0.05			0.08			0.20			0.08 DGS	N/A		<p>Flow Calibration Performed By _____</p> <p><input type="checkbox"/> Flow Calibration N/A <input type="checkbox"/> Flow Calibration Complete Flow Column # _____ <input type="checkbox"/> 5L/min - 17mm <input type="checkbox"/> 15L/min - 53mm <input type="checkbox"/> 30L/min - 103mm</p> <p><input type="checkbox"/> R-Value _____ <input type="checkbox"/> Post Calibration Verification (L/s) Flow Column # _____ 32mm _____ (.139 - .169) 36mm _____ (.156 - .190) 53mm _____ (.228 - .278) 103mm _____ (.447 - .547)</p> <p>Maintenance Performed By _____</p> <p><input type="checkbox"/> Battery Replacement <input type="checkbox"/> Dry Gas Regulator Replacement <input type="checkbox"/> Breath Tube Replacement <input type="checkbox"/> Other _____</p> <p>Suggested Service</p> <p>_____</p> <p>_____</p>
Simulator	Serial #	Lot #/Exp															
0.05																	
0.08																	
0.20																	
0.08 DGS	N/A																

Optical Bench Calibration Performed By SP

Optical Bench Calibration N/A #2 (11/30/16)
 Optical Bench Calibration Complete

Barometric Pressure Gauge 1010 ID # 26932

Simulator	Serial Number	Lot Number	Expiration
0.000	G4444	N/A	N/A
0.040	SD1024	116101	2/2/18
0.100	SD1013	116001	5/8/18
0.200	G2403	116103	4/14/18
0.400	G6621	116102	3/22/18
0.080 DGS	N/A	03415080A1	3/5/17

Post Calibration Stability Checks

Simulator	Serial Number	Lot Number	Expiration
0.05	DR2035	201507A	7/14/18
0.08	SD1011	201601F	1/26/18
0.20	SD1025	201604C	4/5/18
0.08 DGS	N/A	AG626605	9/22/18

Department Inspection Performed By _____

Barometric Pressure _____ Gauge
 ID# _____ Instrument

Mouth Alcohol Solution Lot # _____
 Acetone Stock Solution Lot # _____

Simulator	Serial Number
0.00	
Interferent	
0.05	
0.08	
0.20	

Attachments

Form 41 Optical Bench Cal
 Pre-Stability Tests Post-Stability Tests
 Flow Calibration Other _____

Notes: SP

Instrument Complies with Chapter 11D-8, FAC
 Instrument Does Not Comply with Chapter 11D-8, FAC

Return to/Place into Evidentiary Use
 Remain Out of Evidentiary Use
 Conduct an Agency Inspection Before Evidentiary Use

Pre-Cal
Stability Checks

#80-000632 FL Highway Patrol 9/21/16 *QDS*

QDS

SP
BK

FL HIGHWAY PATROL
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006632
09/21/2016
Software: 8100.27

Test	9/21/16	Time
Air Blank	0.000	16:38
Control Test	0.047	16:38
Air Blank	0.000	16:39
Control Test	0.047	16:40
Air Blank	0.000	16:40
Control Test	0.047	16:41
Air Blank	0.000	16:41
Control Test Stats		
Average	0.0470	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

QDS
Operator's Signature

FL HIGHWAY PATROL
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006632
09/21/2016
Software: 8100.27

Test	9/21/16	Time
Air Blank	0.000	16:43
Control Test	0.074	16:44
Air Blank	0.000	16:45
Control Test	0.076	16:45
Air Blank	0.000	16:46
Control Test	0.076	16:47
Air Blank	0.000	16:47
Control Test Stats		
Average	0.0753	
Std Dev	0.0012	
Rel Std Dev(%)	1.5328	

QDS
Operator's Signature

FL HIGHWAY PATROL
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006632
09/21/2016
Software: 8100.27

Test	9/21/16	Time
Air Blank	0.000	16:49
Control Test	0.194	16:50
Air Blank	0.000	16:51
Control Test	0.193	16:51
Air Blank	0.000	16:52
Control Test	0.194	16:53
Air Blank	0.000	16:53
Control Test Stats		
Average	0.1937	
Std Dev	0.0006	
Rel Std Dev(%)	0.2981	

QDS
Operator's Signature

FL HIGHWAY PATROL
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006632
09/21/2016
Software: 8100.27

Test	9/21/16	Time
Air Blank	0.000	16:55
Control Test	0.081	16:55
Air Blank	0.000	16:56
Control Test	0.081	16:56
Air Blank	0.000	16:56
Control Test	0.081	16:57
Air Blank	0.000	16:57
Control Test Stats		
Average	0.0810	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

QDS
Operator's Signature

Optical Bench Calibration Data # 80-006632 FL Highway Patrol 9/23/16 RAS

SJP

FL HIGHWAY PATROL
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006632
09/22/2016 15:56:06

Auto Calibration
Max Power Res Value = 90
Auto Range Res Value = 56

Sol Value = 0.000 g/210L ***
Fit value = 0.0000 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12551, Sum Io = 12912

<<<< CHANNEL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = -0.0230 (0.0000)
Sample #2 = -0.0120 (0.0000)
Sample #3 = 0.0430 (-0.0200)
Sample #4 = -0.0030 (0.0100)
Avg % Abs = 0.0093 (-0.0033)
STD DEV = 0.0295 (0.0153)
REL STD DEV = 316.087 (458.258)

<<<< CHANNEL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.0570 (-0.0040)
Sample #2 = 0.0940 (-0.0120)
Sample #3 = 0.1450 (-0.0370)
Sample #4 = 0.1080 (-0.0160)
Avg % Abs = 0.1157 (-0.0217)
STD DEV = 0.0264 (0.0134)
REL STD DEV = 22.781 (61.979)

Sol Value = 0.040 g/210L ***
Fit value = 0.1905 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12548, Sum Io = 12916

<<<< CHANNEL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.7370 (-0.0200)
Sample #2 = 0.7170 (-0.0080)
Sample #3 = 0.7620 (-0.0050)
Sample #4 = 0.6980 (0.0270)
Avg % Abs = 0.7257 (0.0047)
STD DEV = 0.0329 (0.0194)
REL STD DEV = 4.529 (415.700)

<<<< CHANNEL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 1.4800 (0.0000)
Sample #2 = 1.4920 (0.0100)
Sample #3 = 1.5000 (0.0090)
Sample #4 = 1.4940 (0.0200)
Avg % Abs = 1.4953 (0.0130)
STD DEV = 0.0042 (0.0061)
REL STD DEV = 0.278 (46.790)

Sol Value = 0.040 g/210L ***
Fit value = 0.1905 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12546, Sum Io = 12912

<<<< CHANNEL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.7630 (-0.0220)
Sample #2 = 0.7630 (-0.0220)
Sample #3 = 0.6950 (0.0140)
Sample #4 = 0.7450 (-0.0040)
Avg % Abs = 0.7343 (-0.0040)
STD DEV = 0.0352 (0.0180)
REL STD DEV = 4.798 (450.080)

<<<< CHANNEL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 1.5060 (-0.0120)
Sample #2 = 1.5340 (-0.0240)
Sample #3 = 1.5040 (-0.0100)
Sample #4 = 1.5250 (0.0000)
Avg % Abs = 1.5210 (-0.0113)
STD DEV = 0.0154 (0.0121)
REL STD DEV = 1.012 (106.371)

**** AUTO CAL FAIL

2nd Calibration Data #80-060632 FL Highway Patrol 9/22/16 DWS

FL HIGHWAY PATROL
Intoxilyzer - Alcohol Analyzer
Model 8000
09/22/2016
SN 80-06632
16:51:50
Auto Calibration
Max Power Res Value = 89
Auto Range Res Value = 56

<<<< CHANNEL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 1.4860 (-0.0230)
Sample #2 = 1.5140 (-0.0200)
Sample #3 = 1.4900 (0.0060)
Sample #4 = 1.5090 (0.0000)
Avg % Abs = 1.5043 (-0.0047)
STD DEV = 0.0127 (0.0136)
REL STD DEV = 0.842 (291.723)

Sol Value = 0.100 g/210L ***
Fit Value = 0.4762 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12551, Sum Io = 12914
<<<< CHANNEL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 1.8110 (-0.0030)
Sample #2 = 1.8260 (0.0140)
Sample #3 = 1.7900 (0.0500)
Sample #4 = 1.8250 (0.0150)
Avg % Abs = 1.8137 (0.0263)
STD DEV = 0.0205 (0.0205)
REL STD DEV = 1.130 (77.856)

<<<< CHANNEL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 3.5220 (-0.0070)
Sample #2 = 3.5770 (-0.0030)
Sample #3 = 3.5780 (0.0010)
Sample #4 = 3.5810 (0.0010)
Avg % Abs = 3.5787 (-0.0003)
STD DEV = 0.0021 (0.0023)
REL STD DEV = 0.058 (692.820)

Sol Value = 0.200 g/210L ***
Fit Value = 0.9524 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12546, Sum Io = 12912
<<<< CHANNEL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 3.5190 (-0.0360)
Sample #2 = 3.5070 (0.0080)
Sample #3 = 3.5010 (0.0140)
Sample #4 = 3.5180 (0.0280)
Avg % Abs = 3.5087 (0.0167)
STD DEV = 0.0086 (0.0103)
REL STD DEV = 0.246 (61.579)

***** AUTO CAL DATA *****
<<<< CHANNEL 1 >>>>
Sol Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.023
Std Dev = 0.03 Rel Std Dev = 146.99
Sol Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 0.739
Std Dev = 0.01 Rel Std Dev = 1.07
Sol Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 1.814
Std Dev = 0.02 Rel Std Dev = 1.13
Sol Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 3.509
Std Dev = 0.01 Rel Std Dev = 0.25
Sol Val = 1.9048 mg/l or 0.400 g/210L
% Abs = 6.744
Std Dev = 0.03 Rel Std Dev = 0.45
Zero Order Coef = -58.96
First Order Coef = 2612.35
Second Order Coef = 32.75
Standard Deviation = 16.15897

<<<< CHANNEL 2 >>>>
Sol Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.097
Std Dev = 0.01 Rel Std Dev = 12.15
Sol Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 1.504
Std Dev = 0.01 Rel Std Dev = 0.84
Sol Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 3.579
Std Dev = 0.00 Rel Std Dev = 0.06
Sol Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 6.857
Std Dev = 0.01 Rel Std Dev = 0.22
Sol Val = 1.9048 mg/l or 0.400 g/210L
% Abs = 12.759
Std Dev = 0.02 Rel Std Dev = 0.16
Zero Order Coef = -105.72
First Order Coef = 1298.89
Second Order Coef = 15.82
Standard Deviation = 20.052130

<<<< CHANNEL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 12.6280 (0.0030)
Sample #2 = 12.7380 (0.0170)
Sample #3 = 12.7610 (0.0170)
Sample #4 = 12.7780 (0.0160)
Avg % Abs = 12.7590 (0.0167)
STD DEV = 0.0201 (0.0066)
REL STD DEV = 0.157 (3.464)

Sol Value = 0.040 g/210L ***
Fit Value = 0.1915 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12551, Sum Io = 12914
<<<< CHANNEL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.7460 (-0.0290)
Sample #2 = 0.7480 (-0.0080)
Sample #3 = 0.7360 (0.0120)
Sample #4 = 0.7330 (0.0230)
Avg % Abs = 0.7390 (0.0090)
STD DEV = 0.0079 (0.0157)
REL STD DEV = 1.074 (174.625)

Solution Stats Quadratic Fit Chan 1
Act Fit Residual
g/210L g/210L g/210L
0.000 0.000 -0.0000
0.040 0.040 0.0003
0.100 0.101 -0.0005
0.200 0.200 0.0003
0.400 0.400 -0.0000

Solution Stats Quadratic Fit Chan 2
Act Fit Residual
g/210L g/210L g/210L
0.000 0.000 -0.0004
0.040 0.040 0.0004
0.100 0.100 0.0004
0.200 0.200 -0.0004
0.400 0.400 0.0001

Sol Value = 0.080 g/210L ***
Fit Value = 0.3810 mg/l %%%
Samples Taken = 4, Discarded = 1
***** CHANNEL 1
Sample #1 = 3449.00
Sample #2 = 3357.00
Sample #3 = 3306.00
Sample #4 = 3316.00
Average Result = 3326.3333
STD DEV = 27.0247
REL STD DEV = 0.812
***** CHANNEL 2
Sample #1 = 3353.00
Sample #2 = 3327.00
Sample #3 = 3325.00
Sample #4 = 3334.00
Average Result = 3328.6667
STD DEV = 4.7258
REL STD DEV = 0.142

Dry Gas H2O Adjust Results *****
Barometric Pressure = 1012
3 um H2O Adjust (mg/l*10,000) = 483
9 um H2O Adjust (mg/l*10,000) = 481
***** AUTO CAL PASS

SP
RSK

Post-Car Stability Checks #80-006632 FZ Highway Patrol 9/22/16 DS

DS

FL HIGHWAY PATROL
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006632
09/22/2016
Software: 8100.27

FL HIGHWAY PATROL
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006632
09/22/2016
Software: 8100.27

FL HIGHWAY PATROL
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006632
09/22/2016
Software: 8100.27

FL HIGHWAY PATROL
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006632
09/22/2016
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	18:34
Control Test	0.082	18:34
Air Blank	0.000	18:35
Control Test	0.082	18:35
Air Blank	0.000	18:36
Control Test	0.081	18:36
Air Blank	0.000	18:37
Control Test Stats		
Average	0.0817	
Std Dev	0.0006	
Rel Std Dev(%)	0.7070	

Test	g/210L	Time
Air Blank	0.000	18:29
Control Test	0.194	18:29
Air Blank	0.000	18:30
Control Test	0.196	18:31
Air Blank	0.000	18:31
Control Test	0.196	18:32
Air Blank	0.000	18:32
Control Test Stats		
Average	0.1953	
Std Dev	0.0012	
Rel Std Dev(%)	0.5911	

Test	g/210L	Time
Air Blank	0.000	18:24
Control Test	0.078	18:25
Air Blank	0.000	18:25
Control Test	0.079	18:26
Air Blank	0.000	18:26
Control Test	0.080	18:27
Air Blank	0.000	18:28
Control Test Stats		
Average	0.0790	
Std Dev	0.0010	
Rel Std Dev(%)	1.2658	

Test	g/210L	Time
Air Blank	0.000	18:19
Control Test	0.048	18:20
Air Blank	0.000	18:20
Control Test	0.049	18:21
Air Blank	0.000	18:22
Control Test	0.050	18:22
Air Blank	0.000	18:23
Control Test Stats		
Average	0.0498	
Std Dev	0.0010	
Rel Std Dev(%)	2.0408	

SP
BSK

RMS
Operator's Signature

RMS
Operator's Signature

RMS
Operator's Signature

RMS
Operator's Signature

Optical Bench Calibration Data # 80-006632 FL Highway Patrol 11/30/16

FL HIGHWAY PATROL
Intoxilyzer - Alcohol Analyzer
Model 8000
11/30/2016

Auto Calibration
Max Power Res Value = 90
Auto Range Res Value = 56
Sol Value = 0.000 g/210L ***
Fit value = 0.000 mg/l ****

Samples Taken = 4, Discarded = 1
Sum Io = 12638, Sum Io = 12904
Sol Value = 0.000 mg/l ****
Fit value = 0.000 mg/l ****

Channel 2 Data
Sample % Abs (% Abs Ref)
Sample #1 = 1.5150 (-0.0200)
Sample #2 = 1.4940 (-0.0050)
Sample #3 = 1.5050 (-0.0100)
Sample #4 = 1.5220 (-0.0020)
Avg % Abs = 1.5070 (-0.0083)
STD DEV = 0.0141 (0.0085)
REL STD DEV = 0.936 (102.055)

Channel 2 Data
Sol Value = 0.100 g/210L ***
Fit value = 0.4762 mg/l ****
Samples Taken = 4, Discarded = 1
Sum Io = 12625, Sum Io = 12901

Channel 2 Data
Sample % Abs (% Abs Ref)
Sample #1 = 1.7680 (-0.0100)
Sample #2 = 1.7680 (-0.0160)
Sample #3 = 1.7930 (-0.0030)
Sample #4 = 1.7950 (0.0000)
Avg % Abs = 1.7853 (-0.0063)
STD DEV = 0.0150 (0.0085)
REL STD DEV = 0.843 (134.288)

Channel 2 Data
Sample % Abs (% Abs Ref)
Sample #1 = 3.5360 (-0.0100)
Sample #2 = 3.5370 (-0.0050)
Sample #3 = 3.5350 (-0.0010)
Sample #4 = 3.5400 (0.0020)
Avg % Abs = 3.5400 (-0.0017)
STD DEV = 0.0070 (0.0040)
REL STD DEV = 0.198 (242.487)

Channel 2 Data
Sol Value = 0.200 g/210L ***
Fit value = 0.9524 mg/l ****
Samples Taken = 4, Discarded = 1
Sum Io = 12626, Sum Io = 12902

Channel 2 Data
Sample % Abs (% Abs Ref)
Sample #1 = 6.9290 (0.0000)
Sample #2 = 6.9380 (0.0060)
Sample #3 = 6.9400 (0.0120)
Sample #4 = 6.9520 (0.0170)
Avg % Abs = 6.9433 (0.0117)
STD DEV = 0.0076 (0.0055)
REL STD DEV = 0.109 (47.208)

Channel 2 Data
Sol Value = 0.400 g/210L ***
Fit value = 1.9048 mg/l ****
Samples Taken = 4, Discarded = 1
Sum Io = 12619, Sum Io = 12900

Channel 2 Data
Sample % Abs (% Abs Ref)
Sample #1 = 6.7210 (-0.0170)
Sample #2 = 6.7480 (0.0030)
Sample #3 = 6.7710 (-0.0040)
Sample #4 = 6.7350 (0.0160)
Avg % Abs = 6.7513 (0.0050)
STD DEV = 0.0182 (0.0101)
REL STD DEV = 0.270 (202.978)

Channel 2 Data
Sample % Abs (% Abs Ref)
Sample #1 = 12.7060 (0.0010)
Sample #2 = 12.7170 (0.0340)
Sample #3 = 12.7490 (0.0250)
Sample #4 = 12.7160 (0.0430)
Avg % Abs = 12.7273 (0.0340)
STD DEV = 0.0188 (0.0090)
REL STD DEV = 0.147 (26.471)

Channel 2 Data
Sol Value = 0.200 g/210L ***
Fit value = 0.9524 mg/l ****
Samples Taken = 4, Discarded = 1
Sum Io = 12626, Sum Io = 12902

Channel 1 Data
Sol Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.001
Std Dev = 0.02 Rel Std Dev = 1757.84
Sol Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 0.740
Std Dev = 0.01 Rel Std Dev = 1.76
Sol Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 1.785
Std Dev = 0.02 Rel Std Dev = 0.84
Sol Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 3.596
Std Dev = 0.01 Rel Std Dev = 0.28
Sol Val = 1.9048 mg/l or 0.400 g/210L
% Abs = 6.751
Std Dev = 0.02 Rel Std Dev = 0.27
Zero Order Coef = 38.12
First Order Coef = 2493.80
Second Order Coef = 47.28
Standard Deviation = 79.386086

Channel 2 Data
Sol Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.098
Std Dev = 0.01 Rel Std Dev = 8.28
Sol Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 1.507
Std Dev = 0.01 Rel Std Dev = 0.94
Sol Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 3.540
Std Dev = 0.01 Rel Std Dev = 0.20
Sol Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 6.943
Std Dev = 0.01 Rel Std Dev = 0.11
Sol Val = 1.9048 mg/l or 0.400 g/210L
% Abs = 12.727
Std Dev = 0.02 Rel Std Dev = 0.15
Zero Order Coef = -69.99
First Order Coef = 1269.57
Second Order Coef = 18.14
Standard Deviation = 79.251427

Channel 1 Data
Sol Value = 0.080 g/210L ***
Fit value = 0.3810 mg/l ****
Samples Taken = 4, Discarded = 1

Channel 1 Data
Sample #1 = 3307.00
Sample #2 = 3316.00
Sample #3 = 3286.00
Sample #4 = 3307.00
Average Result = 3303.0000
STD DEV = 15.3948
REL STD DEV = 0.466

Channel 2 Data
Sample #1 = 3277.00
Sample #2 = 3315.00
Sample #3 = 3323.00
Sample #4 = 3335.00
Average Result = 3324.3333
STD DEV = 10.0664
REL STD DEV = 0.303

Solution Stats Quadratic Fit Chan 1
Act g/210L
Fit g/210L
Residual g/210L

Solution Stats Quadratic Fit Chan 2
Act g/210L
Fit g/210L
Residual g/210L

Sol Value = 0.080 g/210L ***
Fit value = 0.3810 mg/l ****
Samples Taken = 4, Discarded = 1

Sol Value = 0.200 g/210L ***
Fit value = 0.9524 mg/l ****
Samples Taken = 4, Discarded = 1

Dry Gas H2O Adjust Results *****
Barometric Pressure = 1009
3 um H2O Adjust (mg/l*10,000) = 506
9 um H2O Adjust (mg/l*10,000) = 485
**** AUTO CAL PASS

ABK
SP

FL HIGHWAY PATROL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-006632
 11/30/2016 12:23:03

Auto Calibration

pg 1 of 2

	<<<<<	3um	>>>>>	<<<<<	9um	>>>>>

Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1						
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)	
Sample #1	-0.0020	(-0.0130)		0.0880	(-0.0110)	
Sample #2	0.0280	(-0.0110)		0.1050	(-0.0200)	
Sample #3	-0.0160	(0.0180)		0.0890	(-0.0150)	
Sample #4	-0.0080	(0.0230)		0.0990	(-0.0220)	
Avg % Abs	0.0013	(0.0100)		0.0977	(-0.0190)	
STD DEV	0.0234	(0.0184)		0.0081	(0.0036)	
REL STD DEV	1757.839	(183.576)		8.276	(18.977)	

Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1						
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)	
Sample #1	0.7480	(-0.0200)		1.5150	(-0.0200)	
Sample #2	0.7400	(-0.0100)		1.4940	(-0.0050)	
Sample #3	0.7270	(-0.0020)		1.5050	(-0.0180)	
Sample #4	0.7530	(0.0060)		1.5220	(-0.0020)	
Avg % Abs	0.7400	(-0.0020)		1.5070	(-0.0083)	
STD DEV	0.0130	(0.0080)		0.0141	(0.0085)	
REL STD DEV	1.757	(400.000)		0.936	(102.059)	

Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1						
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)	
Sample #1	1.7680	(-0.0180)		3.5360	(-0.0100)	
Sample #2	1.7680	(-0.0160)		3.5370	(-0.0060)	
Sample #3	1.7930	(-0.0030)		3.5350	(-0.0010)	
Sample #4	1.7950	(0.0000)		3.5480	(0.0020)	
Avg % Abs	1.7853	(-0.0063)		3.5400	(-0.0017)	
STD DEV	0.0150	(0.0085)		0.0070	(0.0040)	
REL STD DEV	0.843	(134.288)		0.198	(242.487)	

Solution = 0.200 g/210L or 0.9524 mg/l, Samples = 4, Discarded = 1						
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)	
Sample #1	3.5700	(-0.0140)		6.9290	(0.0000)	
Sample #2	3.6050	(-0.0020)		6.9380	(0.0060)	
Sample #3	3.5850	(0.0100)		6.9400	(0.0120)	
Sample #4	3.5970	(0.0140)		6.9520	(0.0170)	
Avg % Abs	3.5957	(0.0073)		6.9433	(0.0117)	
STD DEV	0.0101	(0.0083)		0.0076	(0.0055)	
REL STD DEV	0.280	(113.545)		0.109	(47.208)	

Solution = 0.400 g/210L or 1.9048 mg/l, Samples = 4, Discarded = 1						
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)	
Sample #1	6.7210	(-0.0170)		12.7060	(0.0010)	
Sample #2	6.7480	(0.0030)		12.7170	(0.0340)	
Sample #3	6.7710	(-0.0040)		12.7490	(0.0250)	
Sample #4	6.7350	(0.0160)		12.7160	(0.0430)	
Avg % Abs	6.7513	(0.0050)		12.7273	(0.0340)	
STD DEV	0.0182	(0.0101)		0.0188	(0.0090)	
REL STD DEV	0.270	(202.978)		0.147	(26.471)	

SP

FL HIGHWAY PATROL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-006632
 11/30/2016 12:23:03

Auto Calibration

pg 2 of 2

```

<<<<<      3um      >>>>>
-----
Zero Order Coef    38.12
First Order Coef   2493.80
Second Order Coef  47.28
  
```

```

<<<<<      9um      >>>>>
-----
                -69.99
                1269.57
                18.14
  
```

```

-----
Act      Fit      Residual
(g/210L) (g/210L) (g/210L)
0.000    0.001    -0.0009
0.040    0.040    -0.0001
0.100    0.097     0.0025
0.200    0.202    -0.0019
0.400    0.400     0.0004
  
```

```

-----
Act      Fit      Residual
(g/210L) (g/210L) (g/210L)
0.000    0.001    -0.0011
0.040    0.040     0.0004
0.100    0.098     0.0023-
0.200    0.202    -0.0020-
0.400    0.400     0.0004
  
```

```

<<<<<      3um      >>>>>
-----
Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1
Sample
  
```

```

<<<<<      9um      >>>>>
-----
  
```

```

Sample #1          3307.00          3277.00
Sample #2          3316.00          3315.00
Sample #3          3286.00          3323.00
Sample #4          3307.00          3335.00
Avg                3303.0000        3324.3333
STD DEV            15.3948           10.0664
REL STD DEV        0.466            0.303
H2O adjust (mg/l*10k) 506          485
  
```

Barometric Pressure = 1009

SP BK

*****CALIBRATION SUCCESSFUL*****

Post-Calibration Stability Checks # 80-200632 FL Highway Patrol 11/30/16

DES

FL HIGHWAY PATROL
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006632
11/30/2016
Software: 8100.27

Test	9/210L	Time
Air Blank	0.000	13:25
Control Test	0.049	13:26
Air Blank	0.000	13:27
Control Test	0.049	13:27
Air Blank	0.000	13:28
Control Test	0.050	13:29
Air Blank	0.000	13:29
Control Test Stats		
Average	0.0493	
Std Dev	0.0006	
Rel Std Dev(%)	1.1703	

[Signature]
Operator's Signature

FL HIGHWAY PATROL
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006632
11/30/2016
Software: 8100.27

Test	9/210L	Time
Air Blank	0.000	13:30
Control Test	0.078	13:31
Air Blank	0.000	13:31
Control Test	0.078	13:32
Air Blank	0.000	13:33
Control Test	0.079	13:33
Air Blank	0.000	13:34
Control Test Stats		
Average	0.0783	
Std Dev	0.0006	
Rel Std Dev(%)	0.7370	

[Signature]
Operator's Signature

FL HIGHWAY PATROL
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006632
11/30/2016
Software: 8100.27

Test	9/210L	Time
Air Blank	0.000	13:35
Control Test	0.194	13:36
Air Blank	0.000	13:36
Control Test	0.195	13:37
Air Blank	0.000	13:37
Control Test	0.196	13:38
Air Blank	0.000	13:39
Control Test Stats		
Average	0.1950	
Std Dev	0.0010	
Rel Std Dev(%)	0.5128	

[Signature]
Operator's Signature

FL HIGHWAY PATROL
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006632
11/30/2016
Software: 8100.27

Test	9/210L	Time
Air Blank	0.000	13:39
Control Test	0.080	13:40
Air Blank	0.000	13:40
Control Test	0.080	13:41
Air Blank	0.000	13:41
Control Test	0.080	13:41
Air Blank	0.000	13:42
Control Test Stats		
Average	0.0800	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

SP
BK

[Signature]
Operator's Signature

INSTRUMENT PROCESSING SHEET

Agency Florida Highway Patrol S/N 8D-006632

Date In 5/26/2016 Date Out 5/25/2016 Ship P/U H/D CMI EE

Intake Performed By DELL

Registration
 Annual
 Return from CMI
 Return from Enforcement
 Electronics
 Other _____

Visual Inspection:
OK Case OK Handle
OK Dry Gas Holder OK Feet
OK Keyboard/Plug OK Back/Plugs
OK Screws tight OK Breath Hose

Other Equipment:
 Power cord
 Printer Cable
 Other ANT STATIC BAG

Notes: AC POWER CORD

Quality Checks Performed By DELL

Breath Tube Screen
 Replace O-Rings
 Instrument Set Up Verified
 R-Value 213
 Flow Verification (L/s)
 Flow Column # ATP 105
 32mm 152 (.139 - .169)
 36mm 171 (.156 - .190)
 53mm 242 (.228 - .278)
 103mm 515 (.447 - .547)

Barometric Pressure Check
 Gauge ID # 28427

Stability Checks

Simulator	Serial #	Lot #/Exp
0.05	G2403	201507A 07/14/2017
0.08	SD3964	201601F 01/26/2018
0.20	G4444	201505A 05/12/2017
0.08 DGS	N/A	A66025D4 01/05/2015

Flow Calibration Performed By _____

Flow Calibration **RECEIVED**
 Flow Calibration Complete
 Flow Column # _____
 5L/min - 17mm
 15L/min - 53mm
 30L/min - 103mm

R-Value _____
 Post Calibration Verification (L/s)
 Flow Column # _____
 32mm _____ (.139 - .169)
 36mm _____ (.156 - .190)
 53mm _____ (.228 - .278)
 103mm _____ (.447 - .547)

Maintenance Performed By _____

Battery Replacement
 Dry Gas Regulator Replacement
 Breath Tube Replacement
 Other _____

Suggested Service

Optical Bench Calibration Performed By _____

Optical Bench Calibration N/A
 Optical Bench Calibration Complete

Barometric Pressure Gauge ID # _____

Simulator	Serial Number	Lot Number	Expiration
0.000		N/A	N/A
0.040			
0.100			
0.200			
0.400			
0.080 DGS	N/A		

Post Calibration Stability Checks

Simulator	Serial Number	Lot Number	Expiration
0.05			
0.08			
0.20			
0.08 DGS	N/A		

Department Inspection Performed By DELL

Barometric Pressure 1020 Gauge
 ID# 26932 1020 Instrument

Mouth Alcohol Solution Lot # 2016A
 Acetone Stock Solution Lot # 2016D

Simulator	Serial Number
0.00	G2880
Interferent	G2834
0.05	G2403
0.08	SD3964
0.20	G4444

Attachments

Form 41 Optical Bench Cal
 Pre-Stability Tests Post-Stability Tests
 Flow Calibration Other _____

Notes: OK PWS

Instrument Complies with Chapter 11D-8, FAC
 Instrument Does Not Comply with Chapter 11D-8, FAC

Return to/Place into Evidentiary Use
 Remain Out of Evidentiary Use

Conduct an Agency Inspection Before Evidentiary Use

Brett Hurdler
 Quality Control Review

5/26/16
 Date

TYPE OF TEST	SERIAL NUMBER	AGENCY	DATE	PERFORMED BY
Stabilities	80-0066302	Florida Highway Patrol	05/24/2016	<i>PELL</i>

SK

0.05g/210L 0.047 to 0.053 <input checked="" type="checkbox"/>	0.08g/210L 0.077 to 0.083 <input checked="" type="checkbox"/>	0.20g/210L 0.194 to 0.206 <input checked="" type="checkbox"/>	DGS 0.08g/210L 0.077 to 0.083 <input checked="" type="checkbox"/>
<p>FHP Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-006632 05/25/2016 Software: 8100.27</p> <p>Test g/210L Time</p> <p>Air Blank 0.000 10:56 Control Test 0.047 10:56 Air Blank 0.000 10:57 Control Test 0.048 10:58 Air Blank 0.000 10:58 Control Test 0.048 10:59 Air Blank 0.000 10:59</p> <p>Control Test Stats Average 0.0477 Std Dev 0.0006 Rel Std Dev(%) 1.2112</p> <p>Operator's Signature <i>PELL</i></p>	<p>FHP Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-006632 05/25/2016 Software: 8100.27</p> <p>Test g/210L Time</p> <p>Air Blank 0.000 11:01 Control Test 0.077 11:01 Air Blank 0.000 11:02 Control Test 0.078 11:03 Air Blank 0.000 11:03 Control Test 0.077 11:04 Air Blank 0.000 11:04</p> <p>Control Test Stats Average 0.0773 Std Dev 0.0006 Rel Std Dev(%) 0.7466</p> <p>Operator's Signature <i>PELL</i></p>	<p>FHP Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-006632 05/25/2016 Software: 8100.27</p> <p>Test g/210L Time</p> <p>Air Blank 0.000 11:07 Control Test 0.196 11:07 Air Blank 0.000 11:08 Control Test 0.198 11:09 Air Blank 0.000 11:09 Control Test 0.198 11:10 Air Blank 0.000 11:10</p> <p>Control Test Stats Average 0.1973 Std Dev 0.0012 Rel Std Dev(%) 0.5852</p> <p>Operator's Signature <i>PELL</i></p>	<p>FHP Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-006632 05/25/2016 Software: 8100.27</p> <p>Test g/210L Time</p> <p>Air Blank 0.000 11:12 Control Test 0.079 11:13 Air Blank 0.000 11:13 Control Test 0.079 11:14 Air Blank 0.000 11:14 Control Test 0.079 11:14 Air Blank 0.000 11:15</p> <p>Control Test Stats Average 0.0790 Std Dev 0.0000 Rel Std Dev(%) 0.0000</p> <p>Operator's Signature <i>PELL</i></p>

SK

TYPE OF TEST	SERIAL NUMBER	AGENCY	DATE	PERFORMED BY
Stabilities	80-006632	Florida Highway Patrol	05/24/2016	<i>Will</i>

0.05g/210L 0.047 to 0.053 <input checked="" type="checkbox"/>	0.08g/210L 0.077 to 0.083 <input checked="" type="checkbox"/>	0.20g/210L 0.194 to 0.206 <input checked="" type="checkbox"/>	DGS 0.08g/210L 0.077 to 0.083 <input checked="" type="checkbox"/>
--	--	--	--

FHP
Intoxilizer - Alcohol Analyzer
Model 8000 SN 80-006632
05/24/2016
Software: 8100.27

Test g/210L Time

Air Blank 0.000 10:37
Control Test 0.048 10:38
Air Blank 0.000 10:39
Control Test 0.049 10:39
Air Blank 0.000 10:40
Control Test 0.049 10:41
Air Blank 0.000 10:41
Control Test 0.000 10:41
Control Test Stats
Average 0.0487
Std Dev 0.0006
Rel Std Dev(x) 1.1863

Operator's Signature *Will*

FHP
Intoxilizer - Alcohol Analyzer
Model 8000 SN 80-006632
05/24/2016
Software: 8100.27

Test g/210L Time

Air Blank 0.000 10:52
Control Test 0.072 10:53
Air Blank 0.000 10:54
Control Test 0.073 10:54
Air Blank 0.000 10:55
Control Test 0.074 10:56
Air Blank 0.000 10:56
Control Test 0.000 10:56
Control Test Stats
Average 0.0730
Std Dev 0.0010
Rel Std Dev(x) 1.3699

Operator's Signature *Will*

Fubbing Simulator

FHP
Intoxilizer - Alcohol Analyzer
Model 8000 SN 80-006632
05/24/2016
Software: 8100.27

Test g/210L Time

Air Blank 0.000 10:48
Control Test 0.198 10:48
Air Blank 0.000 10:49
Control Test 0.198 10:50
Air Blank 0.000 10:50
Control Test 0.198 10:51
Air Blank 0.000 10:51
Control Test 0.000 10:51
Control Test Stats
Average 0.1980
Std Dev 0.0000
Rel Std Dev(x) 0.0000

Operator's Signature *Will*

FHP
Intoxilizer - Alcohol Analyzer
Model 8000 SN 80-006632
05/24/2016
Software: 8100.27

Test g/210L Time

Air Blank 0.000 11:03
Control Test 0.080 11:03
Air Blank 0.000 11:03
Control Test 0.080 11:04
Air Blank 0.000 11:04
Control Test 0.079 11:05
Air Blank 0.000 11:05
Control Test 0.000 11:05
Control Test Stats
Average 0.0797
Std Dev 0.0006
Rel Std Dev(x) 0.7247

Operator's Signature *Will*

Will