



INSTRUMENT PROCESSING SHEET

Agency Pinellas County SOS/N 80-000889

Florida Department of Law Enforcement

Date In 01/12/2018 DI Completion Date 01/24/2018 Ship P/U H/D CMI EE

Intake Performed By <u>TG</u> <input checked="" type="checkbox"/> Annual <input type="checkbox"/> Registration <input type="checkbox"/> Return from CMI / EE Visual Inspection: <input checked="" type="checkbox"/> Case <input checked="" type="checkbox"/> Handle <input checked="" type="checkbox"/> Keyboard <input checked="" type="checkbox"/> Dry Gas Shelf <input checked="" type="checkbox"/> Feet <input checked="" type="checkbox"/> Breath Tube <input checked="" type="checkbox"/> Ports <input checked="" type="checkbox"/> Screws Tight Other Equipment/ Accessories: <input type="checkbox"/> Power cord <input type="checkbox"/> Printer Cable <input type="checkbox"/> Static Bag <input type="checkbox"/> 12V DC Cable Notes: _____ _____ _____	Quality Checks Performed By <u>SDS</u> <input checked="" type="checkbox"/> Breath Tube Screen <input checked="" type="checkbox"/> Replace External O-Rings <input checked="" type="checkbox"/> Instrument Set Up Verified <input checked="" type="checkbox"/> R-Value <u>214</u> <input checked="" type="checkbox"/> Flow Verification (L/s) Flow Column # <u>ATP 105</u> 32 mm <u>0.144</u> (.139 - .169) 36 mm <u>0.160</u> (.156 - .190) 53 mm <u>0.218</u> (.228 - .278) 103 mm <u>0.447</u> (.447 - .547) <input checked="" type="checkbox"/> Barometric Pressure Check Gauge ID # <u>28427</u> <input checked="" type="checkbox"/> Stability Checks	Flow Calibration Performed By <u>SDS</u> Flow Column # <u>ATP 102</u> <input checked="" type="checkbox"/> 5L/min - 17mm <input checked="" type="checkbox"/> 15L/min - 53mm <input checked="" type="checkbox"/> 30L/min - 103mm <input checked="" type="checkbox"/> R-Value <u>214</u> <input checked="" type="checkbox"/> Post Calibration Verification (L/s) Flow Column # <u>ATP 105</u> 32 mm <u>0.152</u> (.139 - .169) 36 mm <u>0.171</u> (.156 - .190) 53 mm <u>0.259</u> (.228 - .278) 103 mm <u>0.542</u> (.447 - .547)
---	---	---

Final Release Date
FDLE
JAN 26 2018
Alcohol Testing Program

Simulator	Serial #	Lot #/Exp
0.050	G2835	201707D 07/25/2019
0.080	SD1013	201707E 07/25/2019
0.200	SD1025	201707C 07/24/2019
0.080 DGS	N/A	AG708807 03/29/2019

Maintenance Performed By SDS
 Battery Replacement
 Dry Gas Regulator Replacement
 Breath Tube Replacement
 Other _____

Temperature Checks Performed By SDS
 Lab Temp °C 21.9
 External Digital Therm. ID#: 300505
 34°C +/- .2 Serial #: G2835
 34°C +/- .2 Serial #: G2882
 34°C +/- .2 Serial #: SD1025

Calibration Adjustment Performed By SDS
 Barometric Pressure Gauge 1025 ID # 28662

Simulator	Serial Number	Lot Number	Expiration
0.000	SD1016	N/A	N/A
0.040	SD1018	16320	10/21/2018
0.100	G2834	17280	09/11/2019
0.200	SD1011	17090	02/24/2019
0.300	DR1275	17140	05/15/2019
0.080 DGS	N/A	17817080A2	08/05/2019

Post Calibration Adjustment Stability Checks

Simulator	Serial Number	Lot Number	Expiration
0.050	G2835	201707D	07/25/2019
0.080	SD1013/G2882	201707E	07/25/2019
0.200	SD1025	201707C	07/24/2019
0.080 DGS	N/A	AG708807	03/29/2019

Department Inspection Performed By SDS
 Barometric Pressure ID# 28427
 Gauge 1022 Instrument 1022
 Mouth Alcohol Solution Lot # 2016-C
 Acetone Stock Solution Lot # 2017-A

Simulator	Serial Number
0.000	G2880
Interferent	G2840
0.050	G2835
0.080	G2882
0.200	SD1025

Attachments

<input checked="" type="checkbox"/> Form 41	<input checked="" type="checkbox"/> Post-Stability Checks
<input checked="" type="checkbox"/> Stability Checks	<input checked="" type="checkbox"/> Flow Calibration
<input checked="" type="checkbox"/> Calibration Certificate	<input checked="" type="checkbox"/> Form 40
<input checked="" type="checkbox"/> Calibration Adjustment	<input type="checkbox"/> Other _____

Notes/Suggested Service: Performed optical bench calibration adjustment to bring values closer to nominal. DMB 1/24/28

Instrument Complies with Chapter 11D-8, FAC
 Instrument Does Not Comply with Chapter 11D-8, FAC

Return to/Place into Evidentiary Use
 Remain Out of Evidentiary Use

Conduct an Agency Inspection Before Evidentiary Use

SP 1/25/18 JJ [Signature] 1/24/18
 Tech Review / Date Admin Review / Date

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: PINELLAS COUNTY SO
Time of Inspection: 15:33

Date of Inspection: 01/24/2018

Serial Number: 80-000889
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:201707D Exp: 07/25/2019	0.08g/210L Test (g/210L) Lot#:201707E Exp: 07/25/2019	0.20g/210L Test (g/210L) Lot#:201707C Exp: 07/24/2019	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG708807 Exp: 03/29/2019
0.000	0.049	0.081	0.202	0.078
0.000	0.050	0.083	0.204	0.079
0.000	0.051	0.081	0.206	0.080
0.000	0.050	0.082	0.207	0.078
0.000	0.050	0.082	0.206	0.078
0.000	0.050	0.083	0.205	0.078
0.000	0.051	0.082	0.207	0.079
0.000	0.050	0.083	0.205	0.079
0.000	0.052	0.083	0.205	0.079
0.000	0.050	0.082	0.206	0.078

Standard Deviations	0.0008	0.0007	0.0014	0.0006
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0008 Number of Simulators Used: 5

Remarks:

SP

The above instrument complies () does not comply () with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Danielle M Bell

DANIELLE M BELL

Signature and Printed Name

01/24/2018
Date

1/26/18
JD

Stability Checks # 80-000889 Pinellas County S.O. 1/22/18 ~~SP~~

DES

PINELAS COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Mode: 8000 SN 80-000889
 01/22/2018
 Software: 8100.27

Test	g/21.0L	Time
Air Blank	0.000	16:09
Control Test	0.051	16:09
Air Blank	0.000	16:09
Control Test	0.052	16:10
Air Blank	0.000	16:11
Control Test	0.050	16:11
Air Blank	0.000	16:12
Control Test Stats		
Average	0.0510	
Std Dev	0.0010	
Rel. Std Dev(%)	1.9508	

[Signature]

Operator's Signature

1/26/18
 SP

PINELAS COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Mode: 8000 SN 80-000889
 01/22/2018
 Software: 8100.27

Test	g/21.0L	Time
Air Blank	0.000	16:13
Control Test	0.082	16:13
Air Blank	0.000	16:14
Control Test	0.084	16:15
Air Blank	0.000	16:15
Control Test	0.083	16:15
Air Blank	0.000	16:17
Control Test Stats		
Average	0.0830	
Std Dev	0.0010	
Rel. Std Dev(%)	1.2048	

[Signature]

Operator's Signature

PINELAS COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Mode: 8000 SN 80-000889
 01/22/2018
 Software: 8100.27

Test	g/21.0L	Time
Air Blank	0.000	16:17
Control Test	0.203	16:19
Air Blank	0.000	16:19
Control Test	0.203	16:19
Air Blank	0.000	16:20
Control Test	0.202	16:20
Air Blank	0.000	16:21
Control Test Stats		
Average	0.2027	
Std Dev	0.0006	
Rel. Std Dev(%)	0.2849	

[Signature]

Operator's Signature

PINELAS COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Mode: 8000 SN 80-000889
 01/22/2018
 Software: 8100.27

Test	g/21.0L	Time
Air Blank	0.000	16:22
Control Test	0.084	16:23
Air Blank	0.000	16:23
Control Test	0.084	16:23
Air Blank	0.000	16:24
Control Test	0.084	16:24
Air Blank	0.000	16:25
Control Test Stats		
Average	0.0840	
Std Dev	0.0000	
Rel. Std Dev(%)	0.0000	

SP

[Signature]

Operator's Signature



Florida Department of Law Enforcement
 Alcohol Testing Program
 2729 Fort Knox Blvd.
 Bldg. 2, Suite 1300
 Tallahassee, FL 32308

Calibration Certificate

This is to certify the calibration of Intoxilyzer 8000 serial number 80-000889, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-000889</u>	UNCERTAINTY* ±
Owning Agency:	<u>PINELLAS COUNTY SO</u>	0.05 g/ 210 L
Calibration Date:	<u>01/24/2018</u>	0.08 g/ 210 L
Calibration Time:	<u>15:33</u>	0.20 g/ 210 L
		0.080 g/ 210 L Dry Gas Control
		0.004
		0.005
		0.008
		0.005

All results are reported in g/ 210 L.

Bias is limited by Calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the Target Alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99% level of confidence (k=3).

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS) ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Thermometer temperatures are checked with NIST traceable Eutechnics 4400 digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the uses of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

01/24/2018

Date

DANIELLE M BELL,
 Department Inspector

FDLE/ATP Form 69 January 2018
 Issuing Authority: Alcohol Testing Program

SP

1/24/18
 SP

Optical Bench Calibration Adjustment Data #80-000889 Pinellas County S.O. 10/24/18

PINE...S COUNTY SO
Integrator - Alcohol Analyzer
Mode: 8000
01/21/2018 10:38:14
SN 80-000889

Auto Calibration
Max Power Res Value = 61
Auto Range Res Value = 53
Sol Value = 0.000 g/210L ***
Fit Value = 0.0000 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Lo = 12742, Sum Hi = 13939

Sample % Abs (% Abs Ref)
Sample #1 = 1.6010 (-0.0140)
Sample #2 = 1.5330 (0.0210)
Sample #3 = 1.5450 (0.0260)
Sample #4 = 1.4860 (0.0350)
Avg % Abs = 1.5213 (0.0273)
STD DEV = 0.0312 (0.0071)
REL STD DEV = 2.050 (25.956)

Sample % Abs (% Abs Ref)
Sample #1 = 2.0020 (-0.0240)
Sample #2 = 1.9660 (0.0610)
Sample #3 = 1.9450 (0.0570)
Sample #4 = 1.9510 (0.0510)
Avg % Abs = 1.9540 (0.0563)
STD DEV = 0.0118 (0.0050)
REL STD DEV = 0.554 (8.935)

Sample % Abs (% Abs Ref)
Sample #1 = 0.1820 (-0.0240)
Sample #2 = 0.1580 (0.0060)
Sample #3 = 0.1640 (0.0080)
Sample #4 = 0.1700 (0.0150)
Avg % Abs = 0.1640 (0.0097)
STD DEV = 0.0060 (0.0047)
REL STD DEV = 3.659 (48.888)

Auto Calibration
Max Power Res Value = 61
Auto Range Res Value = 53
Sol Value = 0.000 g/210L ***
Fit Value = 0.0000 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Lo = 12742, Sum Hi = 13939

Sample % Abs (% Abs Ref)
Sample #1 = 1.6010 (-0.0140)
Sample #2 = 1.5330 (0.0210)
Sample #3 = 1.5450 (0.0260)
Sample #4 = 1.4860 (0.0350)
Avg % Abs = 1.5213 (0.0273)
STD DEV = 0.0312 (0.0071)
REL STD DEV = 2.050 (25.956)

Sample % Abs (% Abs Ref)
Sample #1 = 2.0020 (-0.0240)
Sample #2 = 1.9660 (0.0610)
Sample #3 = 1.9450 (0.0570)
Sample #4 = 1.9510 (0.0510)
Avg % Abs = 1.9540 (0.0563)
STD DEV = 0.0118 (0.0050)
REL STD DEV = 0.554 (8.935)

Sample % Abs (% Abs Ref)
Sample #1 = 0.1820 (-0.0240)
Sample #2 = 0.1580 (0.0060)
Sample #3 = 0.1640 (0.0080)
Sample #4 = 0.1700 (0.0150)
Avg % Abs = 0.1640 (0.0097)
STD DEV = 0.0060 (0.0047)
REL STD DEV = 3.659 (48.888)

Sample % Abs (% Abs Ref)
Sample #1 = 0.1820 (-0.0240)
Sample #2 = 0.1580 (0.0060)
Sample #3 = 0.1640 (0.0080)
Sample #4 = 0.1700 (0.0150)
Avg % Abs = 0.1640 (0.0097)
STD DEV = 0.0060 (0.0047)
REL STD DEV = 3.659 (48.888)

Auto Calibration
Max Power Res Value = 61
Auto Range Res Value = 53
Sol Value = 0.000 g/210L ***
Fit Value = 0.0000 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Lo = 12742, Sum Hi = 13939

Sample % Abs (% Abs Ref)
Sample #1 = 1.6010 (-0.0140)
Sample #2 = 1.5330 (0.0210)
Sample #3 = 1.5450 (0.0260)
Sample #4 = 1.4860 (0.0350)
Avg % Abs = 1.5213 (0.0273)
STD DEV = 0.0312 (0.0071)
REL STD DEV = 2.050 (25.956)

Sample % Abs (% Abs Ref)
Sample #1 = 2.0020 (-0.0240)
Sample #2 = 1.9660 (0.0610)
Sample #3 = 1.9450 (0.0570)
Sample #4 = 1.9510 (0.0510)
Avg % Abs = 1.9540 (0.0563)
STD DEV = 0.0118 (0.0050)
REL STD DEV = 0.554 (8.935)

Sample % Abs (% Abs Ref)
Sample #1 = 0.1820 (-0.0240)
Sample #2 = 0.1580 (0.0060)
Sample #3 = 0.1640 (0.0080)
Sample #4 = 0.1700 (0.0150)
Avg % Abs = 0.1640 (0.0097)
STD DEV = 0.0060 (0.0047)
REL STD DEV = 3.659 (48.888)

Sample % Abs (% Abs Ref)
Sample #1 = 0.1820 (-0.0240)
Sample #2 = 0.1580 (0.0060)
Sample #3 = 0.1640 (0.0080)
Sample #4 = 0.1700 (0.0150)
Avg % Abs = 0.1640 (0.0097)
STD DEV = 0.0060 (0.0047)
REL STD DEV = 3.659 (48.888)

Auto Calibration
Max Power Res Value = 61
Auto Range Res Value = 53
Sol Value = 0.000 g/210L ***
Fit Value = 0.0000 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Lo = 12742, Sum Hi = 13939

Sample % Abs (% Abs Ref)
Sample #1 = 1.6010 (-0.0140)
Sample #2 = 1.5330 (0.0210)
Sample #3 = 1.5450 (0.0260)
Sample #4 = 1.4860 (0.0350)
Avg % Abs = 1.5213 (0.0273)
STD DEV = 0.0312 (0.0071)
REL STD DEV = 2.050 (25.956)

Sample % Abs (% Abs Ref)
Sample #1 = 2.0020 (-0.0240)
Sample #2 = 1.9660 (0.0610)
Sample #3 = 1.9450 (0.0570)
Sample #4 = 1.9510 (0.0510)
Avg % Abs = 1.9540 (0.0563)
STD DEV = 0.0118 (0.0050)
REL STD DEV = 0.554 (8.935)

Sample % Abs (% Abs Ref)
Sample #1 = 0.1820 (-0.0240)
Sample #2 = 0.1580 (0.0060)
Sample #3 = 0.1640 (0.0080)
Sample #4 = 0.1700 (0.0150)
Avg % Abs = 0.1640 (0.0097)
STD DEV = 0.0060 (0.0047)
REL STD DEV = 3.659 (48.888)

Sample % Abs (% Abs Ref)
Sample #1 = 0.1820 (-0.0240)
Sample #2 = 0.1580 (0.0060)
Sample #3 = 0.1640 (0.0080)
Sample #4 = 0.1700 (0.0150)
Avg % Abs = 0.1640 (0.0097)
STD DEV = 0.0060 (0.0047)
REL STD DEV = 3.659 (48.888)

Auto Calibration
Max Power Res Value = 61
Auto Range Res Value = 53
Sol Value = 0.000 g/210L ***
Fit Value = 0.0000 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Lo = 12742, Sum Hi = 13939

Sample % Abs (% Abs Ref)
Sample #1 = 1.6010 (-0.0140)
Sample #2 = 1.5330 (0.0210)
Sample #3 = 1.5450 (0.0260)
Sample #4 = 1.4860 (0.0350)
Avg % Abs = 1.5213 (0.0273)
STD DEV = 0.0312 (0.0071)
REL STD DEV = 2.050 (25.956)

Sample % Abs (% Abs Ref)
Sample #1 = 2.0020 (-0.0240)
Sample #2 = 1.9660 (0.0610)
Sample #3 = 1.9450 (0.0570)
Sample #4 = 1.9510 (0.0510)
Avg % Abs = 1.9540 (0.0563)
STD DEV = 0.0118 (0.0050)
REL STD DEV = 0.554 (8.935)

Sample % Abs (% Abs Ref)
Sample #1 = 0.1820 (-0.0240)
Sample #2 = 0.1580 (0.0060)
Sample #3 = 0.1640 (0.0080)
Sample #4 = 0.1700 (0.0150)
Avg % Abs = 0.1640 (0.0097)
STD DEV = 0.0060 (0.0047)
REL STD DEV = 3.659 (48.888)

Sample % Abs (% Abs Ref)
Sample #1 = 0.1820 (-0.0240)
Sample #2 = 0.1580 (0.0060)
Sample #3 = 0.1640 (0.0080)
Sample #4 = 0.1700 (0.0150)
Avg % Abs = 0.1640 (0.0097)
STD DEV = 0.0060 (0.0047)
REL STD DEV = 3.659 (48.888)

***** AUTO CAL DATA *****

Channel 1
Sol Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.164
Std Dev = 0.01 Rel Std Dev = 3.66
Sol Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 0.918
Std Dev = 0.03 Rel Std Dev = 2.90
Sol Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 1.954
Std Dev = 0.01 Rel Std Dev = 0.55
Sol Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 3.664
Std Dev = 0.01 Rel Std Dev = 0.40
Sol Val = 1.4286 mg/l or 0.300 g/210L
% Abs = 5.303
Std Dev = 0.01 Rel Std Dev = 0.20
Zero Order Coef = -451.47
First Order Coef = 2587.07
Second Order Coef = 36.33
Standard Deviation = 29.552658

***** CHANNEL 1 *****

Channel 2
Sol Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.119
Std Dev = 0.07 Rel Std Dev = 57.76
Sol Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 1.521
Std Dev = 0.03 Rel Std Dev = 2.05
Sol Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 3.521
Std Dev = 0.01 Rel Std Dev = 0.41
Sol Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 6.783
Std Dev = 0.01 Rel Std Dev = 0.15
Sol Val = 1.4286 mg/l or 0.300 g/210L
% Abs = 9.817
Std Dev = 0.02 Rel Std Dev = 0.25
Zero Order Coef = -157.28
First Order Coef = 1342.68
Second Order Coef = 13.01
Standard Deviation = 20.794878

***** CHANNEL 2 *****

Channel 2
Sol Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.119
Std Dev = 0.07 Rel Std Dev = 57.76
Sol Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 1.521
Std Dev = 0.03 Rel Std Dev = 2.05
Sol Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 3.521
Std Dev = 0.01 Rel Std Dev = 0.41
Sol Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 6.783
Std Dev = 0.01 Rel Std Dev = 0.15
Sol Val = 1.4286 mg/l or 0.300 g/210L
% Abs = 9.817
Std Dev = 0.02 Rel Std Dev = 0.25
Zero Order Coef = -157.28
First Order Coef = 1342.68
Second Order Coef = 13.01
Standard Deviation = 20.794878

***** CHANNEL 2 *****

Channel 2
Sample % Abs (% Abs Ref)
Sample #1 = 6.9640 (-0.0410)
Sample #2 = 6.7870 (0.1310)
Sample #3 = 6.7910 (0.1650)
Sample #4 = 6.7720 (0.1670)
Avg % Abs = 6.7833 (0.1543)
STD DEV = 0.0100 (0.0202)
REL STD DEV = 0.148 (13.109)

***** CHANNEL 2 *****

Channel 2
Sample % Abs (% Abs Ref)
Sample #1 = 5.3870 (0.0000)
Sample #2 = 5.2950 (0.1320)
Sample #3 = 5.3000 (0.1270)
Sample #4 = 5.3150 (0.1510)
Avg % Abs = 5.3033 (0.1367)
STD DEV = 0.0104 (0.0127)
REL STD DEV = 0.196 (9.265)

***** CHANNEL 2 *****

Channel 2
Sample % Abs (% Abs Ref)
Sample #1 = 10.0230 (-0.0330)
Sample #2 = 9.8440 (0.1930)
Sample #3 = 9.7950 (0.2380)
Sample #4 = 9.8110 (0.2510)
Avg % Abs = 9.8167 (0.2273)
STD DEV = 0.0250 (0.0304)
REL STD DEV = 0.255 (13.388)

***** CHANNEL 2 *****

Channel 2
Sample % Abs (% Abs Ref)
Sample #1 = 1.6010 (-0.0140)
Sample #2 = 1.5330 (0.0210)
Sample #3 = 1.5450 (0.0260)
Sample #4 = 1.4860 (0.0350)
Avg % Abs = 1.5213 (0.0273)
STD DEV = 0.0312 (0.0071)
REL STD DEV = 2.050 (25.956)

***** CHANNEL 2 *****

Channel 2
Sample % Abs (% Abs Ref)
Sample #1 = 2.0020 (-0.0240)
Sample #2 = 1.9660 (0.0610)
Sample #3 = 1.9450 (0.0570)
Sample #4 = 1.9510 (0.0510)
Avg % Abs = 1.9540 (0.0563)
STD DEV = 0.0118 (0.0050)
REL STD DEV = 0.554 (8.935)

***** CHANNEL 2 *****

Channel 2
Sample % Abs (% Abs Ref)
Sample #1 = 3.7590 (-0.0240)
Sample #2 = 3.6790 (0.0430)
Sample #3 = 3.6500 (0.1040)
Sample #4 = 3.6640 (0.0790)
Avg % Abs = 3.6643 (0.0753)
STD DEV = 0.0145 (0.0307)
REL STD DEV = 0.396 (40.706)

***** CHANNEL 2 *****

Channel 2
Sample % Abs (% Abs Ref)
Sample #1 = 0.1690 (-0.0070)
Sample #2 = 0.1740 (-0.0090)
Sample #3 = 0.1400 (0.0070)
Sample #4 = 0.1420 (0.0690)
Avg % Abs = 0.1187 (0.0223)
STD DEV = 0.0695 (0.0412)
REL STD DEV = 57.756 (184.472)

***** CHANNEL 2 *****

Channel 2
Sample % Abs (% Abs Ref)
Sample #1 = 0.1690 (-0.0070)
Sample #2 = 0.1740 (-0.0090)
Sample #3 = 0.1400 (0.0070)
Sample #4 = 0.1420 (0.0690)
Avg % Abs = 0.1187 (0.0223)
STD DEV = 0.0695 (0.0412)
REL STD DEV = 57.756 (184.472)

***** CHANNEL 2 *****

Channel 2
Sample % Abs (% Abs Ref)
Sample #1 = 0.9550 (-0.0370)
Sample #2 = 0.8870 (-0.0060)
Sample #3 = 0.9340 (0.0000)
Sample #4 = 0.9320 (0.0020)
Avg % Abs = 0.9177 (-0.0020)
STD DEV = 0.0266 (0.0053)
REL STD DEV = 2.896 (264.575)

Dry Gas H2O Adjust Results *****
Barometric Pressure = 1024
3 um H2O Adjust (mg/l*10,000) = 931
9 um H2O Adjust (mg/l*10,000) = 354
***** AUTO CAL PASS

SP

10/24/18

Post-Calibration Adjust Stability Checks

#80-000889 Pinellas County S.D. 1/24/18 ~~0089~~

DS

PINELLAS COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model: 8000
01/24/2018
Software: 8100.27

PINELLAS COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model: 8000
01/24/2018
Software: 8100.27

PINELLAS COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model: 8000
01/24/2018
Software: 8100.27

PINELLAS COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model: 8000
01/24/2018
Software: 8100.27

PINELLAS COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model: 8000
01/24/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:21
Control Test	0.052	12:31
Air Blank	0.000	12:3
Control Test	0.051	12:3
Air Blank	0.000	12:3
Control Test	0.050	12:3
Air Blank	0.000	12:3
Control Test Stats		
Average	0.0510	
Std Dev	0.0010	
Rel Std Dev(%)	1.9608	

Test	g/210L	Time
Air Blank	0.000	12:34
Control Test	0.084	12:35
Air Blank	0.000	12:36
Control Test	0.083	12:36
Air Blank	0.000	12:37
Control Test	0.082	12:38
Air Blank	0.000	12:38
Control Test Stats		
Average	0.0830	
Std Dev	0.0010	
Rel Std Dev(%)	1.2048	

Test	g/210L	Time
Air Blank	0.000	13:47
Control Test	0.080	13:48
Air Blank	0.000	13:48
Control Test	0.081	13:49
Air Blank	0.000	13:50
Control Test	0.082	13:51
Air Blank	0.000	13:51
Control Test Stats		
Average	0.0810	
Std Dev	0.0010	
Rel Std Dev(%)	1.2346	

Test	g/210L	Time
Air Blank	0.000	12:44
Control Test	0.202	12:45
Air Blank	0.000	12:45
Control Test	0.205	12:46
Air Blank	0.000	12:47
Control Test	0.202	12:47
Air Blank	0.000	12:48
Control Test Stats		
Average	0.2030	
Std Dev	0.0017	
Rel Std Dev(%)	0.8532	

Test	g/210L	Time
Air Blank	0.000	12:24
Control Test	0.079	12:24
Air Blank	0.000	12:25
Control Test	0.079	12:25
Air Blank	0.000	12:26
Control Test	0.079	12:26
Air Blank	0.000	12:26
Control Test Stats		
Average	0.0790	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

SID1013
Changed sims

G2882

RMS
Operator's Signature

RMS
Operator's Signature

RMS
Operator's Signature

RMS
Operator's Signature

RMS
Operator's Signature

1/26/18
SP

SP

PINELLAS COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000889
01/22/2018
Software: 8100.27

Flow Rate Calibration*****

1: Rate (Liters/min) = 5

SQRT(Diff)) = 6.633

2: Rate (Liters/min) = 15

SQRT(Diff)) = 10.535

3: Rate (Liters/min) = 30

SQRT(Diff)) = 17.887

Dependent Data Scale Factor = 100000 L/min

Independent Data Scale Factor = 256

Rounded Slope = 859

Rounded Intercept = -901692

Correlation = 0.99819

Flow Calibration Adjust Data

80-000889

Pinellas County Sheriff's Office

01/22/18



SP

1/24/18
SO

Florida Department of Law Enforcement Alcohol Testing Program

AGENCY INSPECTION REPORT - INTOXILYZER 8000

Agency: PINELLAS COUNTY SO
Time of Inspection: 15:45

Date of Inspection: 01/22/2018

Serial Number: 80-000889
Software: 8100.27

Check or Test	YES	NO
Date and/or Time Adjusted		No
Diagnostic Check (Pre-Inspection): OK		No
Alcohol Free Subject Test: 0.000		No
Mouth Alcohol Test: Slope Not Met		No
Interferent Detect Test: Interferent Detect		No
Diagnostic Check (Post-Inspection): OK		No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#: _____ Exp: _____	0.08g/210L Test (g/210L) Lot#: _____ Exp: _____	0.20g/210L Test (g/210L) Lot#: _____ Exp: _____	0.08 g/210L Dry Gas Std Test (g/210L) Lot#: _____ Exp: _____

Number of Simulators Used: _____

Remarks:
BYPASSED AI TO OPERATE INSTRUMENT. NOT A COMPLIANCE CHECK.

SP

N/A Compliance Not Determined

The above instrument complies () does not comply () with Chapter 11D-8, FAC.

I certify that I hold a valid Florida Department of Law Enforcement Agency Inspector Permit and that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

[Handwritten Signature]

DANIELLE M BELL

Signature and Printed Name

01/22/2018
Date

1/22/18
JTB