



INSTRUMENT PROCESSING SHEET

Agency FFWCCS/N 80-000903

Florida Department of Law Enforcement

Date In 10/30/2018DI Completion Date 11/16/18 Ship P/U H/D CMI EE

Intake Performed By <u>SP</u> <input checked="" type="checkbox"/> Annual <input type="checkbox"/> Registration <input type="checkbox"/> Return from CMI / EE Visual Inspection: <input checked="" type="checkbox"/> Case <input checked="" type="checkbox"/> Handle <input checked="" type="checkbox"/> Keyboard <input checked="" type="checkbox"/> Dry Gas Shelf <input checked="" type="checkbox"/> Feet <input checked="" type="checkbox"/> Breath Tube <input checked="" type="checkbox"/> Ports <input checked="" type="checkbox"/> Screws Tight Other Equipment/ Accessories: <input type="checkbox"/> Power cord <input type="checkbox"/> Printer Cable <input checked="" type="checkbox"/> Static Bag <input type="checkbox"/> 12V DC Cable Notes: _____ _____ _____	Quality Checks Performed By <u>SP</u> <input checked="" type="checkbox"/> Breath Tube Screen <input checked="" type="checkbox"/> Replace External O-Rings <input checked="" type="checkbox"/> Instrument Set Up Verified <input checked="" type="checkbox"/> R-Value <u>155</u> <input checked="" type="checkbox"/> Flow Verification (L/s) Flow Column # <u>ATP105</u> 32 mm <u>.144</u> (.139 - .169) 36 mm <u>.164</u> (.156 - .190) 53 mm <u>.234</u> (.228 - .278) 103 mm <u>.511</u> (.447 - .547) <input checked="" type="checkbox"/> Barometric Pressure Check Gauge ID # <u>26932</u> <input checked="" type="checkbox"/> Stability Checks	Flow Calibration Performed By _____ Flow Column # _____ <input type="checkbox"/> 5L/min - 17mm <input type="checkbox"/> 15L/min - 53mm <input type="checkbox"/> 30L/min - 103mm <input type="checkbox"/> R-Value _____ <input type="checkbox"/> Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)																																								
Final Release Date <p style="text-align: center;">FDLE</p> <p style="text-align: center;">NOV 16 2018</p> <p style="text-align: center;">Alcohol Testing Program</p>	<table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td><u>SD1012</u></td> <td>201707D 07/25/2019</td> </tr> <tr> <td>0.080</td> <td><u>DR1279</u></td> <td>201707E 07/25/2019</td> </tr> <tr> <td>0.200</td> <td><u>DR3856</u></td> <td>201707C 07/24/2019</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td><u>AG805702</u> <u>2-26-20</u></td> </tr> </tbody> </table>	Simulator	Serial #	Lot #/Exp	0.050	<u>SD1012</u>	201707D 07/25/2019	0.080	<u>DR1279</u>	201707E 07/25/2019	0.200	<u>DR3856</u>	201707C 07/24/2019	0.080 DGS	N/A	<u>AG805702</u> <u>2-26-20</u>	Maintenance Performed By _____ <input type="checkbox"/> Battery Replacement <input type="checkbox"/> Dry Gas Regulator Replacement <input type="checkbox"/> Breath Tube Replacement <input type="checkbox"/> Other _____ Temperature Checks Performed By <u>SP</u> <input checked="" type="checkbox"/> Lab Temp °C <u>22.1</u> External Digital Therm. ID#: <u>300948</u> <input checked="" type="checkbox"/> 34°C +/- .2 Serial #: <u>SD1012</u> <input checked="" type="checkbox"/> 34°C +/- .2 Serial #: <u>DR1279</u> <input checked="" type="checkbox"/> 34°C +/- .2 Serial #: <u>DR3856</u>																									
Simulator	Serial #	Lot #/Exp																																								
0.050	<u>SD1012</u>	201707D 07/25/2019																																								
0.080	<u>DR1279</u>	201707E 07/25/2019																																								
0.200	<u>DR3856</u>	201707C 07/24/2019																																								
0.080 DGS	N/A	<u>AG805702</u> <u>2-26-20</u>																																								
Calibration Adjustment Performed By <u>SP</u> Barometric Pressure Gauge <u>1021</u> ID # <u>28662</u> <table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td><u>G2834</u></td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td><u>SD1022</u></td> <td><u>17410</u></td> <td><u>12-6-19</u></td> </tr> <tr> <td>0.100</td> <td><u>SD3964</u></td> <td><u>18070</u></td> <td><u>2-26-20</u></td> </tr> <tr> <td>0.200</td> <td><u>SD1025</u></td> <td><u>17340</u></td> <td><u>10-9-19</u></td> </tr> <tr> <td>0.300</td> <td><u>SD1024</u></td> <td><u>18110</u></td> <td><u>4-2-20</u></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td><u>1781708DA2</u></td> <td><u>8-5-19</u></td> </tr> </tbody> </table> <input checked="" type="checkbox"/> Post Calibration Adjustment Stability Checks	Simulator	Serial Number	Lot Number	Expiration	0.000	<u>G2834</u>	N/A	N/A	0.040	<u>SD1022</u>	<u>17410</u>	<u>12-6-19</u>	0.100	<u>SD3964</u>	<u>18070</u>	<u>2-26-20</u>	0.200	<u>SD1025</u>	<u>17340</u>	<u>10-9-19</u>	0.300	<u>SD1024</u>	<u>18110</u>	<u>4-2-20</u>	0.080 DGS	N/A	<u>1781708DA2</u>	<u>8-5-19</u>	Department Inspection Performed By <u>SP</u> Barometric Pressure ID# <u>26932</u> Gauge <u>1022</u> Instrument <u>1023</u> Mouth Alcohol Solution Lot # <u>2018-B</u> Acetone Stock Solution Lot # <u>2018-A</u> <table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td><u>G2408</u></td> </tr> <tr> <td>Interferent</td> <td><u>G2882</u></td> </tr> <tr> <td>0.050</td> <td><u>SD1012</u></td> </tr> <tr> <td>0.080</td> <td><u>DR1279</u></td> </tr> <tr> <td>0.200</td> <td><u>DR3856</u></td> </tr> </tbody> </table>	Simulator	Serial Number	0.000	<u>G2408</u>	Interferent	<u>G2882</u>	0.050	<u>SD1012</u>	0.080	<u>DR1279</u>	0.200	<u>DR3856</u>	Attachments <input checked="" type="checkbox"/> Form 41 <input checked="" type="checkbox"/> Stability Checks <input checked="" type="checkbox"/> Calibration Certificate <input checked="" type="checkbox"/> Calibration Adjustment <input checked="" type="checkbox"/> Post-Stability Checks <input type="checkbox"/> Flow Calibration <input type="checkbox"/> Form 40 <input type="checkbox"/> Other _____
Simulator	Serial Number	Lot Number	Expiration																																							
0.000	<u>G2834</u>	N/A	N/A																																							
0.040	<u>SD1022</u>	<u>17410</u>	<u>12-6-19</u>																																							
0.100	<u>SD3964</u>	<u>18070</u>	<u>2-26-20</u>																																							
0.200	<u>SD1025</u>	<u>17340</u>	<u>10-9-19</u>																																							
0.300	<u>SD1024</u>	<u>18110</u>	<u>4-2-20</u>																																							
0.080 DGS	N/A	<u>1781708DA2</u>	<u>8-5-19</u>																																							
Simulator	Serial Number																																									
0.000	<u>G2408</u>																																									
Interferent	<u>G2882</u>																																									
0.050	<u>SD1012</u>																																									
0.080	<u>DR1279</u>																																									
0.200	<u>DR3856</u>																																									
Notes/Suggested Service: _____ _____ _____ _____ _____	<input checked="" type="checkbox"/> Instrument Complies with Chapter 11D-8, FAC <input type="checkbox"/> Instrument Does Not Comply with Chapter 11D-8, FAC <input checked="" type="checkbox"/> Return to/Place into Evidentiary Use <input type="checkbox"/> Remain Out of Evidentiary Use <input checked="" type="checkbox"/> Conduct an Agency Inspection Before Evidentiary Use <u>J. Dehn 11/16/18</u> <u>Brett Kirkland 11/16/18</u> Tech Review / Date Admin Review / Date																																									

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: FFWCC SOUTH REGION
Time of Inspection: 11:50

Date of Inspection: 11/16/2018

Serial Number: 80-000903
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:201707D Exp: 07/25/2019	0.08g/210L Test (g/210L) Lot#:201707E Exp: 07/25/2019	0.20g/210L Test (g/210L) Lot#:201707C Exp: 07/24/2019	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG805701 Exp: 02/26/2020
0.000	0.051	0.082	0.200	0.079
0.000	0.052	0.082	0.202	0.080
0.000	0.051	0.083	0.202	0.080
0.000	0.052	0.082	0.202	0.080
0.000	0.052	0.083	0.203	0.079
0.000	0.051	0.083	0.202	0.079
0.000	0.051	0.082	0.202	0.079
0.000	0.051	0.082	0.203	0.080
0.000	0.052	0.083	0.202	0.080
0.000	0.052	0.082	0.202	0.079

Standard Deviations	0.0005	0.0005	0.0008	0.0005
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0005 Number of Simulators Used: 5

Remarks:

The above instrument complies () does not comply () with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Shayla Platt

Signature and Printed Name SHAYLA D PLATT

11/16/2018
Date

11/16/18
JD
BK

STABILITY CHECKS # 80-000903

FFWCC SOUTH REGION
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000903
11/09/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:59
Control Test	0.050	09:59
Air Blank	0.000	10:00
Control Test	0.050	10:00
Air Blank	0.000	10:01
Control Test	0.050	10:02
Air Blank	0.000	10:02
Control Test Stats		
Average	0.0500	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

SJP

Operator's Signature

FFWCC SOUTH REGION
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000903
11/09/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:43
Control Test	0.081	09:43
Air Blank	0.000	09:44
Control Test	0.080	09:44
Air Blank	0.000	09:45
Control Test	0.081	09:46
Air Blank	0.000	09:46
Control Test Stats		
Average	0.0807	
Std Dev	0.0006	
Rel Std Dev(%)	0.7157	

SJP

Operator's Signature

FFWCC SOUTH REGION
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000903
11/09/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:37
Control Test	0.200	09:38
Air Blank	0.000	09:39
Control Test	0.199	09:39
Air Blank	0.000	09:40
Control Test	0.200	09:41
Air Blank	0.000	09:41
Control Test Stats		
Average	0.1997	
Std Dev	0.0006	
Rel Std Dev(%)	0.2892	

SJP

Operator's Signature

FFWCC SOUTH REGION
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000903
11/09/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:49
Control Test	0.076	09:49
Air Blank	0.000	09:50
Control Test	0.076	09:50
Air Blank	0.000	09:51
Control Test	0.076	09:51
Air Blank	0.000	09:52
Control Test Stats		
Average	0.0760	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

DAS

SJP

Operator's Signature

11/16/18
JA
BK



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
2729 Fort Knox Blvd.
Bldg. 2, Suite 1300
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-000903, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-000903</u>	UNCERTAINTY* ±	
Owning Agency:	<u>FFWCC SOUTH REGION</u>	0.050 g/ 210 L	0.004
Calibration Date:	<u>11/16/2018</u>	0.080 g/ 210 L	0.005
Calibration Time:	<u>11:50</u>	0.200 g/ 210 L	0.008
		0.080 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.
*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence (k=3).

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Thermometer temperatures are checked with NIST traceable Eutechnics 4400 digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the uses of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

11/16/2018

Date

Shayla Platt

SHAYLA D PLATT,
Department Inspector

FDLE/ATP Form 69 July 2018

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality

*11/16/18
JD
BSK*

FFLACC SOUTH REGION

Intoxilyzer - Alcohol Analyzer
Model: 8000
11/14/2018

Auto Calibration
Max Power Res Value = 38
Auto Range Res Value = 31

<<<<< CHANNEL 2 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 1.5720 (-0.0120)
Sample #2 = 1.5400 (0.0110)
Sample #3 = 1.5530 (0.0120)
Sample #4 = 1.5300 (0.0140)
Avg % Abs = 1.5410 (0.0123)
STD DEV = 0.0115 (0.0015)
REL STD DEV = 0.748 (12.385)

Soil Value = 0.100 g/210L ***
Fit Value = 0.4762 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12738, Sum Io = 14316

<<<<< CHANNEL 1 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 1.8920 (-0.0320)
Sample #2 = 1.8470 (0.0210)
Sample #3 = 1.8710 (0.0100)
Sample #4 = 1.8110 (0.0430)
Avg % Abs = 1.8430 (0.0247)
STD DEV = 0.0302 (0.0168)
REL STD DEV = 1.639 (68.119)

<<<<< CHANNEL 2 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 3.6000 (-0.0110)
Sample #2 = 3.6040 (0.0300)
Sample #3 = 3.6410 (0.0090)
Sample #4 = 3.5910 (0.0310)
Avg % Abs = 3.6120 (0.0233)
STD DEV = 0.0259 (0.0124)
REL STD DEV = 0.718 (53.242)

Soil Value = 0.200 g/210L ***
Fit Value = 0.9524 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12734, Sum Io = 14303

<<<<< CHANNEL 1 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 3.5410 (-0.0140)
Sample #2 = 3.5550 (0.0020)
Sample #3 = 3.5300 (0.0420)
Sample #4 = 3.5610 (0.0480)
Avg % Abs = 3.5487 (0.0307)
STD DEV = 0.0164 (0.0250)
REL STD DEV = 0.463 (81.543)

<<<<< CHANNEL 2 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 6.7220 (0.0050)
Sample #2 = 6.7830 (0.0260)
Sample #3 = 6.7470 (0.0410)
Sample #4 = 6.7790 (0.0440)
Avg % Abs = 6.7697 (0.0370)
STD DEV = 0.0197 (0.0096)
REL STD DEV = 0.291 (26.164)

Soil Value = 0.300 g/210L ***
Fit Value = 1.4286 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12728, Sum Io = 14297

<<<<< CHANNEL 1 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 5.2180 (-0.0170)
Sample #2 = 5.2280 (0.0200)
Sample #3 = 5.2570 (0.0190)
Sample #4 = 5.2680 (0.0190)
Avg % Abs = 5.2510 (0.0193)
STD DEV = 0.0207 (0.0006)
REL STD DEV = 0.394 (2.986)

<<<<< CHANNEL 2 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 9.8200 (-0.0190)
Sample #2 = 9.8470 (0.0030)
Sample #3 = 9.8940 (-0.0040)
Sample #4 = 9.8750 (0.0260)
Avg % Abs = 9.8720 (0.0083)
STD DEV = 0.0236 (0.0157)
REL STD DEV = 0.239 (188.340)

***** AUTO CAL DATA *****

<<<<< CHANNEL 1 >>>>>
Sol Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.045
Std Dev = 0.02 Rel Std Dev = 53.93
Sol Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 0.755
Std Dev = 0.02 Rel Std Dev = 2.26
Sol Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 1.843
Std Dev = 0.03 Rel Std Dev = 1.64
Sol Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 3.549
Std Dev = 0.02 Rel Std Dev = 0.46
Sol Val = 1.4286 mg/l or 0.300 g/210L
% Abs = 5.251
Std Dev = 0.02 Rel Std Dev = 0.39
Zero Order Coef = -116.73
First Order Coef = 2630.28
Second Order Coef = 21.78
Standard Deviation = 29.503632

<<<<< CHANNEL 2 >>>>>
Sol Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.154
Std Dev = 0.01 Rel Std Dev = 7.30
Sol Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 1.541
Std Dev = 0.01 Rel Std Dev = 0.75
Sol Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 3.612
Std Dev = 0.03 Rel Std Dev = 0.72
Sol Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 6.770
Std Dev = 0.02 Rel Std Dev = 0.29
Sol Val = 1.4286 mg/l or 0.300 g/210L
% Abs = 9.872
Std Dev = 0.02 Rel Std Dev = 0.24
Zero Order Coef = -219.33
First Order Coef = 1348.83
Second Order Coef = 12.35
Standard Deviation = 36.391979

CAL ADJUSTMENT SP
#80-000903

11/16/18
JA
BK

Solution Stats Quadratic Fit Chan 1
Act Fit Residual
g/210L g/210L g/210L
0.000 0.000 -0.0000
0.040 0.040 0.0005
0.100 0.101 -0.0009
0.200 0.199 0.0007
0.300 0.300 -0.0002

Solution Stats Quadratic Fit Chan 2
Act Fit Residual
g/210L g/210L g/210L
0.000 -0.000 0.0001
0.040 0.040 0.0003
0.100 0.101 -0.0011
0.200 0.199 0.0010
0.300 0.300 -0.0003

Sol Value = 0.080 g/210L ***
Fit Value = 0.3810 mg/l %%%
Samples Taken = 4, Discarded = 1
***** CHANNEL 1
Sample #1 = 3430.00
Sample #2 = 3537.00
Sample #3 = 3435.00
Sample #4 = 3424.00
Average Result = 3465.3333
STD DEV = 62.3084
REL STD DEV = 1.798

***** CHANNEL 2
Sample #1 = 3294.00
Sample #2 = 3327.00
Sample #3 = 3291.00
Sample #4 = 3271.00
Average Result = 3296.3333
STD DEV = 28.3784
REL STD DEV = 0.861

Dry Gas H2O Adjust Results *****
Barometric Pressure = 1020
3 um H2O Adjust (mg/l*10,000) = 344
9 um H2O Adjust (mg/l*10,000) = 513
**** AUTO CAL PASS

POST CAL ADJUST
STABILITY CHECKS # 80-000903

FFWCC SOUTH REGION
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000903
11/15/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	14:38
Control Test	0.050	14:38
Air Blank	0.000	14:39
Control Test	0.049	14:40
Air Blank	0.000	14:40
Control Test	0.049	14:41
Air Blank	0.000	14:41
Control Test Stats		
Average	0.0493	
Std Dev	0.0006	
Rel Std Dev(%)	1.1703	

SP

Operator's Signature

FFWCC SOUTH REGION
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000903
11/15/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	14:32
Control Test	0.080	14:33
Air Blank	0.000	14:34
Control Test	0.081	14:34
Air Blank	0.000	14:35
Control Test	0.081	14:36
Air Blank	0.000	14:36
Control Test Stats		
Average	0.0807	
Std Dev	0.0006	
Rel Std Dev(%)	0.7157	

SP

Operator's Signature

FFWCC SOUTH REGION
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000903
11/15/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	14:42
Control Test	0.198	14:43
Air Blank	0.000	14:43
Control Test	0.198	14:44
Air Blank	0.000	14:45
Control Test	0.200	14:45
Air Blank	0.000	14:46
Control Test Stats		
Average	0.1987	
Std Dev	0.0012	
Rel Std Dev(%)	0.5812	

SP

Operator's Signature

FFWCC SOUTH REGION
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000903
11/15/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	14:29
Control Test	0.079	14:29
Air Blank	0.000	14:30
Control Test	0.079	14:30
Air Blank	0.000	14:31
Control Test	0.079	14:31
Air Blank	0.000	14:31
Control Test Stats		
Average	0.0790	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

DGS

SP

Operator's Signature

11/16/18
JD
BSK