



INSTRUMENT PROCESSING SHEET

Agency Walton County SOS/N 80-001312Florida Department of
Law EnforcementDate In 06/04/2018DI Completion Date 6/6/2018 Ship P/U H/D CMI EE

Intake Performed By <u>JR</u> <input checked="" type="checkbox"/> Annual <input type="checkbox"/> Registration <input type="checkbox"/> Return from CMI / EE Visual Inspection: <input checked="" type="checkbox"/> Case <input checked="" type="checkbox"/> Handle <input checked="" type="checkbox"/> Keyboard <input checked="" type="checkbox"/> Dry Gas Shelf <input checked="" type="checkbox"/> Feet <input checked="" type="checkbox"/> Breath Tube <input checked="" type="checkbox"/> Ports <input checked="" type="checkbox"/> Screws Tight Other Equipment/ Accessories: <input type="checkbox"/> Power cord <input type="checkbox"/> Printer Cable <input type="checkbox"/> Static Bag <input type="checkbox"/> 12V DC Cable Notes: _____ _____ _____	Quality Checks Performed By <u>QPM</u> <input checked="" type="checkbox"/> Breath Tube Screen <input checked="" type="checkbox"/> Replace External O-Rings <input checked="" type="checkbox"/> Instrument Set Up Verified <input checked="" type="checkbox"/> R-Value <u>1.02</u> <input checked="" type="checkbox"/> Flow Verification (L/s) Flow Column # <u>ATP103</u> 32 mm <u>.144</u> (.139 - .169) 36 mm <u>.164</u> (.156 - .190) 53 mm <u>.230</u> (.228 - .278) 103 mm <u>.484</u> (.447 - .547) <input checked="" type="checkbox"/> Barometric Pressure Check Gauge ID # <u>28427</u> <input checked="" type="checkbox"/> Stability Checks	Flow Calibration Performed By _____ Flow Column # _____ <input type="checkbox"/> 5L/min - 17mm <input type="checkbox"/> 15L/min - 53mm <input type="checkbox"/> 30L/min - 103mm <input type="checkbox"/> R-Value _____ <input type="checkbox"/> Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)
---	---	--

Final Release Date FDLE JUN 08 2018 Alcohol Testing Program	<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td><u>SD1018</u></td> <td><u>201707D</u> <u>07/25/2019</u></td> </tr> <tr> <td>0.080</td> <td><u>SD3962</u></td> <td><u>201707E</u> <u>07/25/2019</u></td> </tr> <tr> <td>0.200</td> <td><u>G2078</u></td> <td><u>201707C</u> <u>07/24/2019</u></td> </tr> <tr> <td>0.080 DGS</td> <td><u>N/A</u></td> <td><u>AG805702</u> <u>2/26/20</u></td> </tr> </tbody> </table>	Simulator	Serial #	Lot #/Exp	0.050	<u>SD1018</u>	<u>201707D</u> <u>07/25/2019</u>	0.080	<u>SD3962</u>	<u>201707E</u> <u>07/25/2019</u>	0.200	<u>G2078</u>	<u>201707C</u> <u>07/24/2019</u>	0.080 DGS	<u>N/A</u>	<u>AG805702</u> <u>2/26/20</u>	Maintenance Performed By _____ <input type="checkbox"/> Battery Replacement <input type="checkbox"/> Dry Gas Regulator Replacement <input type="checkbox"/> Breath Tube Replacement <input type="checkbox"/> Other _____ Temperature Checks Performed By <u>QPM</u> <input checked="" type="checkbox"/> Lab Temp °C <u>21.2</u> External Digital Therm. ID#: <u>300503</u> <input checked="" type="checkbox"/> 34°C +- .2 Serial #: <u>SD1021</u> <input checked="" type="checkbox"/> 34°C +- .2 Serial #: <u>DR1275</u> <input checked="" type="checkbox"/> 34°C +- .2 Serial #: <u>SD1019</u>
Simulator	Serial #	Lot #/Exp															
0.050	<u>SD1018</u>	<u>201707D</u> <u>07/25/2019</u>															
0.080	<u>SD3962</u>	<u>201707E</u> <u>07/25/2019</u>															
0.200	<u>G2078</u>	<u>201707C</u> <u>07/24/2019</u>															
0.080 DGS	<u>N/A</u>	<u>AG805702</u> <u>2/26/20</u>															

Calibration Adjustment Performed By <u>QPM</u> Barometric Pressure Gauge <u>1010</u> ID # <u>28662</u> <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td><u>G8144</u></td> <td><u>N/A</u></td> <td><u>N/A</u></td> </tr> <tr> <td>0.040</td> <td><u>G2403</u></td> <td><u>17410</u></td> <td><u>12/6/19</u></td> </tr> <tr> <td>0.100</td> <td><u>G2879</u></td> <td><u>18070</u></td> <td><u>2/26/20</u></td> </tr> <tr> <td>0.200</td> <td><u>G3709</u></td> <td><u>17340</u></td> <td><u>10/9/19</u></td> </tr> <tr> <td>0.300</td> <td><u>G8149</u></td> <td><u>18110</u></td> <td><u>4/2/20</u></td> </tr> <tr> <td>0.080 DGS</td> <td><u>N/A</u></td> <td><u>22817080A3</u></td> <td><u>10/5/19</u></td> </tr> </tbody> </table> <input checked="" type="checkbox"/> Post Calibration Adjustment Stability Checks <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td><u>SD1021</u></td> <td><u>201707D</u></td> <td><u>7/25/19</u></td> </tr> <tr> <td>0.080</td> <td><u>DR1275</u></td> <td><u>201707E</u></td> <td><u>7/25/19</u></td> </tr> <tr> <td>0.200</td> <td><u>SD1019</u></td> <td><u>201707C</u></td> <td><u>7/24/19</u></td> </tr> <tr> <td>0.080 DGS</td> <td><u>N/A</u></td> <td><u>AG805702</u></td> <td><u>2/26/20</u></td> </tr> </tbody> </table>	Simulator	Serial Number	Lot Number	Expiration	0.000	<u>G8144</u>	<u>N/A</u>	<u>N/A</u>	0.040	<u>G2403</u>	<u>17410</u>	<u>12/6/19</u>	0.100	<u>G2879</u>	<u>18070</u>	<u>2/26/20</u>	0.200	<u>G3709</u>	<u>17340</u>	<u>10/9/19</u>	0.300	<u>G8149</u>	<u>18110</u>	<u>4/2/20</u>	0.080 DGS	<u>N/A</u>	<u>22817080A3</u>	<u>10/5/19</u>	Simulator	Serial Number	Lot Number	Expiration	0.050	<u>SD1021</u>	<u>201707D</u>	<u>7/25/19</u>	0.080	<u>DR1275</u>	<u>201707E</u>	<u>7/25/19</u>	0.200	<u>SD1019</u>	<u>201707C</u>	<u>7/24/19</u>	0.080 DGS	<u>N/A</u>	<u>AG805702</u>	<u>2/26/20</u>	Department Inspection Performed By <u>QPM</u> Barometric Pressure ID# <u>28427</u> Gauge <u>1011</u> Instrument <u>1010</u> Mouth Alcohol Solution Lot # <u>2017-B</u> Acetone Stock Solution Lot # <u>2018-A</u> <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td><u>G11621</u></td> </tr> <tr> <td>Interferent</td> <td><u>DR3855</u></td> </tr> <tr> <td>0.050</td> <td><u>SD1021</u></td> </tr> <tr> <td>0.080</td> <td><u>DR1275</u></td> </tr> <tr> <td>0.200</td> <td><u>SD1019</u></td> </tr> </tbody> </table>	Simulator	Serial Number	0.000	<u>G11621</u>	Interferent	<u>DR3855</u>	0.050	<u>SD1021</u>	0.080	<u>DR1275</u>	0.200	<u>SD1019</u>
Simulator	Serial Number	Lot Number	Expiration																																																										
0.000	<u>G8144</u>	<u>N/A</u>	<u>N/A</u>																																																										
0.040	<u>G2403</u>	<u>17410</u>	<u>12/6/19</u>																																																										
0.100	<u>G2879</u>	<u>18070</u>	<u>2/26/20</u>																																																										
0.200	<u>G3709</u>	<u>17340</u>	<u>10/9/19</u>																																																										
0.300	<u>G8149</u>	<u>18110</u>	<u>4/2/20</u>																																																										
0.080 DGS	<u>N/A</u>	<u>22817080A3</u>	<u>10/5/19</u>																																																										
Simulator	Serial Number	Lot Number	Expiration																																																										
0.050	<u>SD1021</u>	<u>201707D</u>	<u>7/25/19</u>																																																										
0.080	<u>DR1275</u>	<u>201707E</u>	<u>7/25/19</u>																																																										
0.200	<u>SD1019</u>	<u>201707C</u>	<u>7/24/19</u>																																																										
0.080 DGS	<u>N/A</u>	<u>AG805702</u>	<u>2/26/20</u>																																																										
Simulator	Serial Number																																																												
0.000	<u>G11621</u>																																																												
Interferent	<u>DR3855</u>																																																												
0.050	<u>SD1021</u>																																																												
0.080	<u>DR1275</u>																																																												
0.200	<u>SD1019</u>																																																												

Notes/Suggested Service: _____ _____ _____ _____ _____	<input checked="" type="checkbox"/> Instrument Complies with Chapter 11D-8, FAC <input type="checkbox"/> Instrument Does Not Comply with Chapter 11D-8, FAC <input checked="" type="checkbox"/> Return to/Place into Evidentiary Use <input type="checkbox"/> Remain Out of Evidentiary Use <input checked="" type="checkbox"/> Conduct an Agency Inspection Before Evidentiary Use <u>SP 6/8/18</u> <u>J. Zaban 6/8/18</u> Tech Review / Date Admin Review / Date
--	--

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: WALTON COUNTY SO
Time of Inspection: 11:06

Date of Inspection: 06/06/2018

Serial Number: 80-001312
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:201707D Exp: 07/25/2019	0.08g/210L Test (g/210L) Lot#:201707E Exp: 07/25/2019	0.20g/210L Test (g/210L) Lot#:201707C Exp: 07/24/2019	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG805702 Exp: 02/26/2020
0.000	0.048	0.081	0.199	0.079
0.000	0.049	0.081	0.200	0.079
0.000	0.050	0.081	0.200	0.079
0.000	0.049	0.082	0.200	0.079
0.000	0.050	0.081	0.200	0.079
0.000	0.049	0.081	0.200	0.079
0.000	0.049	0.081	0.200	0.079
0.000	0.049	0.082	0.200	0.079
0.000	0.050	0.082	0.199	0.079
0.000	0.050	0.082	0.200	0.079

Standard Deviations	0.0006	0.0005	0.0004	0.0000
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0003 Number of Simulators Used: 5

Remarks:

SP

The above instrument complies () does not comply () with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Patrick J Murphy

PATRICK J MURPHY

Signature and Printed Name

06/06/2018
Date

6/8/18
J2

WALTON COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001312
06/05/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:19
Control Test	0.049	12:20
Air Blank	0.000	12:21
Control Test	0.050	12:21
Air Blank	0.000	12:22
Control Test	0.050	12:23
Air Blank	0.000	12:23
Control Test Stats		
Average	0.0497	
Std Dev	0.0006	
Rel Std Dev(%)	1.1625	

P Murphy
Operator's Signature

WALTON COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001312
06/05/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:25
Control Test	0.082	12:26
Air Blank	0.000	12:26
Control Test	0.082	12:27
Air Blank	0.000	12:27
Control Test	0.082	12:28
Air Blank	0.000	12:29
Control Test Stats		
Average	0.0820	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

P Murphy
Operator's Signature

WALTON COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001312
06/05/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:30
Control Test	0.198	12:31
Air Blank	0.000	12:31
Control Test	0.202	12:32
Air Blank	0.000	12:33
Control Test	0.202	12:33
Air Blank	0.000	12:34
Control Test Stats		
Average	0.2007	
Std Dev	0.0023	
Rel Std Dev(%)	1.1509	

P Murphy
Operator's Signature

WALTON COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001312
06/05/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:35
Control Test	0.078	12:35
Air Blank	0.000	12:36
Control Test	0.078	12:36
Air Blank	0.000	12:37
Control Test	0.078	12:37
Air Blank	0.000	12:38
Control Test Stats		
Average	0.0780	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

DGS

P Murphy
Operator's Signature

SP
6/8/18
JD



Florida Department of Law Enforcement
 Alcohol Testing Program
 2729 Fort Knox Blvd.
 Bldg. 2, Suite 1300
 Tallahassee, FL 32308

Calibration Certificate

This is to certify the calibration of Intoxilyzer 8000 serial number 80-001312, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-001312</u>	UNCERTAINTY* ±	
Owning Agency:	<u>WALTON COUNTY SO</u>	0.050 g/ 210 L	0.004
Calibration Date:	<u>06/06/2018</u>	0.080 g/ 210 L	0.005
Calibration Time:	<u>11:06</u>	0.200 g/ 210 L	0.008
		0.080 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.
 Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.
 *Uncertainty is based on fleet-wide data and is expressed to a 99% level of confidence (k=3).

TRACEABILITY INFORMATION
 This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Thermometer temperatures are checked with NIST traceable Eutechnics 4400 digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the uses of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

SP

06/06/2018

Date

Patrick J Murphy

PATRICK J MURPHY,
 Department Inspector

FDLE/ATP Form 69 March 2018
 Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality

6/8/18
JD

WALTON COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001312
06/06/2018 07:58:29

Auto Calibration
Max Power Res Value = 33
Auto Range Res Value = 7

Sol Value = 0.000 g/210L ***
Fit value = 0.0000 mg/l %%%
Samples Taken = 4, Discarded = 1
3um lo = 12659, 9um lo = 14060

<<<< CHANNEL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.1450 (-0.0030)
Sample #2 = 0.1340 (0.0570)
Sample #3 = 0.1580 (0.1070)
Sample #4 = 0.1390 (0.1450)
Avg % Abs = 0.1437 (0.1030)
STD DEU = 0.0127 (0.0441)
REL STD DEU = 8.814 (42.851)

<<<< CHANNEL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.1560 (0.0000)
Sample #2 = 0.1320 (0.0280)
Sample #3 = 0.1320 (0.0380)
Sample #4 = 0.1630 (0.0210)
Avg % Abs = 0.1423 (0.0290)
STD DEU = 0.0179 (0.0085)
REL STD DEU = 12.575 (29.462)

Sol Value = 0.040 g/210L ***
Fit value = 0.1905 mg/l %%%
Samples Taken = 4, Discarded = 1
3um lo = 12630, 9um lo = 14048

<<<< CHANNEL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.8710 (-0.0160)
Sample #2 = 0.8640 (0.0350)
Sample #3 = 0.8760 (0.0440)
Sample #4 = 0.8620 (0.0900)
Avg % Abs = 0.8673 (0.0563)
STD DEU = 0.0076 (0.0295)
REL STD DEU = 0.873 (52.369)

<<<< CHANNEL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 1.5840 (-0.0210)
Sample #2 = 1.5730 (-0.0070)
Sample #3 = 1.5800 (-0.0120)
Sample #4 = 1.5790 (0.0010)
Avg % Abs = 1.5773 (-0.0060)
STD DEU = 0.0038 (0.0066)
REL STD DEU = 0.240 (109.291)

Sol Value = 0.100 g/210L ***
Fit value = 0.4762 mg/l %%%
Samples Taken = 4, Discarded = 1
3um lo = 12614, 9um lo = 14044

<<<< CHANNEL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 1.9020 (-0.0140)
Sample #2 = 1.9320 (0.0360)
Sample #3 = 1.9390 (0.0420)
Sample #4 = 1.9300 (0.0690)
Avg % Abs = 1.9337 (0.0490)
STD DEU = 0.0047 (0.0176)
REL STD DEU = 0.244 (35.874)

<<<< CHANNEL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 3.6540 (-0.0120)
Sample #2 = 3.6790 (0.0020)
Sample #3 = 3.6830 (0.0000)
Sample #4 = 3.6710 (0.0190)
Avg % Abs = 3.6777 (0.0070)
STD DEU = 0.0061 (0.0104)
REL STD DEU = 0.166 (149.147)

Sol Value = 0.200 g/210L ***
Fit value = 0.9524 mg/l %%%
Samples Taken = 4, Discarded = 1
3um lo = 12599, 9um lo = 14038

<<<< CHANNEL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 3.5470 (-0.0110)
Sample #2 = 3.5710 (0.0000)
Sample #3 = 3.5750 (0.0360)
Sample #4 = 3.5690 (0.0420)
Avg % Abs = 3.5717 (0.0260)
STD DEU = 0.0031 (0.0227)
REL STD DEU = 0.086 (87.368)

<<<< CHANNEL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 6.8310 (0.0070)
Sample #2 = 6.8490 (0.0030)
Sample #3 = 6.8870 (0.0150)
Sample #4 = 6.8600 (0.0210)
Avg % Abs = 6.8653 (0.0130)
STD DEU = 0.0196 (0.0092)
REL STD DEU = 0.285 (70.501)

Sol Value = 0.300 g/210L ***
Fit value = 1.4286 mg/l %%%
Samples Taken = 4, Discarded = 1
3um lo = 12589, 9um lo = 14034

<<<< CHANNEL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 5.2450 (-0.0270)
Sample #2 = 5.2600 (-0.0100)
Sample #3 = 5.2580 (0.0240)
Sample #4 = 5.2710 (0.0350)
Avg % Abs = 5.2630 (0.0163)
STD DEU = 0.0070 (0.0235)
REL STD DEU = 0.133 (143.628)

<<<< CHANNEL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 9.9650 (-0.0070)
Sample #2 = 10.0080 (-0.0020)
Sample #3 = 10.0190 (0.0150)
Sample #4 = 10.0480 (0.0100)
Avg % Abs = 10.0250 (0.0077)
STD DEU = 0.0207 (0.0087)
REL STD DEU = 0.206 (113.960)

***** AUTO CAL DATA *****
<<<< CHANNEL 1 >>>>
Sol Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.144
Std Dev = 0.01 Rel Std Dev = 8.81
Sol Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 0.867
Std Dev = 0.01 Rel Std Dev = 0.87
Sol Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 1.934
Std Dev = 0.00 Rel Std Dev = 0.24
Sol Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 3.572
Std Dev = 0.00 Rel Std Dev = 0.09
Sol Val = 1.4286 mg/l or 0.300 g/210L
% Abs = 5.263
Std Dev = 0.01 Rel Std Dev = 0.13
Zero Order Coef = -418.17
First Order Coef = 2673.11
Second Order Coef = 24.00
Standard Deviation = 63.570698

SP
6/8/18
J2

<<<< CHANNEL 2 >>>>

Sol Val = 0.0000 mg/l or 0.000 g/210L
 % Abs = 0.142
 Std Dev = 0.02 Rel Std Dev = 12.57
 Sol Val = 0.1905 mg/l or 0.040 g/210L
 % Abs = 1.577
 Std Dev = 0.00 Rel Std Dev = 0.24
 Sol Val = 0.4762 mg/l or 0.100 g/210L
 % Abs = 3.678
 Std Dev = 0.01 Rel Std Dev = 0.17
 Sol Val = 0.9524 mg/l or 0.200 g/210L
 % Abs = 6.865
 Std Dev = 0.02 Rel Std Dev = 0.28
 Sol Val = 1.4286 mg/l or 0.300 g/210L
 % Abs = 10.025
 Std Dev = 0.02 Rel Std Dev = 0.21
 Zero Order Coef = -208.38
 First Order Coef = 1321.38
 Second Order Coef = 12.63
 Standard Deviation = 46.798744

Solution Stats Quadratic Fit Chan 1

Act	Fit	Residual
g/210L	g/210L	g/210L
0.000	-0.001	0.0007
0.040	0.040	-0.0003
0.100	0.102	-0.0016
0.200	0.198	0.0019
0.300	0.301	-0.0006

Solution Stats Quadratic Fit Chan 2

Act	Fit	Residual
g/210L	g/210L	g/210L
0.000	-0.000	0.0004
0.040	0.040	-0.0001
0.100	0.101	-0.0013
0.200	0.199	0.0014
0.300	0.300	-0.0005

Sol Value = 0.080 g/210L ***
 Fit value = 0.3810 mg/l ***
 Samples Taken = 4, Discarded = 1
 ***** CHANNEL 1
 Sample #1 = 2913.00
 Sample #2 = 2908.00
 Sample #3 = 2867.00
 Sample #4 = 2842.00
 Average Result = 2872.3333
 STD DEV = 33.3217
 REL STD DEV = 1.160

 ***** CHANNEL 2
 Sample #1 = 3297.00
 Sample #2 = 3299.00
 Sample #3 = 3332.00
 Sample #4 = 3316.00
 Average Result = 3315.6667
 STD DEV = 16.5025
 REL STD DEV = 0.498

 Dry Gas H2O Adjust Results *****
 Barometric Pressure = 1010
 3 um H2O Adjust (mg/l*10,000) = 937
 9 um H2O Adjust (mg/l*10,000) = 494
 **** AUTO CAL PASS

WALTON COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-001312
 06/06/2018
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:40
Control Test	0.049	08:40
Air Blank	0.000	08:41
Control Test	0.049	08:41
Air Blank	0.000	08:42
Control Test	0.049	08:43
Air Blank	0.000	08:43
Control Test Stats		
Average	0.0490	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

P. Murphy
 Operator's Signature

WALTON COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-001312
 06/06/2018
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:51
Control Test	0.197	08:52
Air Blank	0.000	08:52
Control Test	0.197	08:53
Air Blank	0.000	08:54
Control Test	0.198	08:54
Air Blank	0.000	08:55
Control Test Stats		
Average	0.1973	
Std Dev	0.0006	
Rel Std Dev(%)	0.2926	

P. Murphy
 Operator's Signature

80-001312 Post stability

WALTON COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-001312
 06/06/2018
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:45
Control Test	0.079	08:45
Air Blank	0.000	08:46
Control Test	0.079	08:47
Air Blank	0.000	08:47
Control Test	0.080	08:48
Air Blank	0.000	08:48
Control Test Stats		
Average	0.0793	
Std Dev	0.0006	
Rel Std Dev(%)	0.7277	

P. Murphy
 Operator's Signature

WALTON COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-001312
 06/06/2018
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:56
Control Test	0.079	08:57
Air Blank	0.000	08:57
Control Test	0.079	08:57
Air Blank	0.000	08:58
Control Test	0.080	08:58
Air Blank	0.000	08:59
Control Test Stats		
Average	0.0793	
Std Dev	0.0006	
Rel Std Dev(%)	0.7277	

DGS

P. Murphy
 Operator's Signature

SP
 6/8/18
 JLB