



INSTRUMENT PROCESSING SHEET

Agency Holmes Beach PDS/N 80-006692

Florida Department of Law Enforcement

Date In 05/25/2018DI Completion Date 6/7/18 Ship P/U H/D CMI EE

Intake Performed By <u>JA</u>	Quality Checks Performed By <u>SP</u>	Flow Calibration Performed By _____
<input checked="" type="checkbox"/> Annual <input type="checkbox"/> Registration <input type="checkbox"/> Return from CMI / EE Visual Inspection: <input checked="" type="checkbox"/> Case <input checked="" type="checkbox"/> Handle <input checked="" type="checkbox"/> Keyboard <input checked="" type="checkbox"/> Dry Gas Shelf <input checked="" type="checkbox"/> Feet <input checked="" type="checkbox"/> Breath Tube <input checked="" type="checkbox"/> Ports <input checked="" type="checkbox"/> Screws Tight Other Equipment/ Accessories: <input type="checkbox"/> Power cord <input type="checkbox"/> Printer Cable <input type="checkbox"/> Static Bag <input type="checkbox"/> 12V DC Cable Notes: <u>Printer cover came off during shipping. It was reinstalled when the box was opened.</u>	<input checked="" type="checkbox"/> Breath Tube Screen <input checked="" type="checkbox"/> Replace External O-Rings <input checked="" type="checkbox"/> Instrument Set Up Verified <input checked="" type="checkbox"/> R-Value <u>228</u> <input checked="" type="checkbox"/> Flow Verification (L/s) Flow Column # <u>ATP105</u> 32 mm <u>.156</u> (.139 - .169) 36 mm <u>.171</u> (.156 - .190) 53 mm <u>.242</u> (.228 - .278) 103 mm <u>.511</u> (.447 - .547) <input checked="" type="checkbox"/> Barometric Pressure Check Gauge ID # <u>28662</u> <input checked="" type="checkbox"/> Stability Checks	Flow Column # _____ <input type="checkbox"/> 5L/min - 17mm <input type="checkbox"/> 15L/min - 53mm <input type="checkbox"/> 30L/min - 103mm <input type="checkbox"/> R-Value _____ <input type="checkbox"/> Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)

Final Release Date	Maintenance Performed By _____	Temperature Checks Performed By <u>SP</u>
FDLE JUN 07 2018 Alcohol Testing Program	<input type="checkbox"/> Battery Replacement <input type="checkbox"/> Dry Gas Regulator Replacement <input type="checkbox"/> Breath Tube Replacement <input type="checkbox"/> Other _____	<input checked="" type="checkbox"/> Lab Temp °C <u>21.3</u> External Digital Therm. ID#: <u>300505</u> <input checked="" type="checkbox"/> 34°C +/- .2 Serial #: <u>SD1018</u> <input checked="" type="checkbox"/> 34°C +/- .2 Serial #: <u>SD3962</u> <input checked="" type="checkbox"/> 34°C +/- .2 Serial #: <u>G2078</u>

Calibration Adjustment Performed By <u>SP</u>																												
Barometric Pressure Gauge <u>1009</u> ID # <u>28427</u>																												
<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td><u>G8144</u></td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td><u>G2403</u></td> <td><u>17410</u></td> <td><u>12-6-19</u></td> </tr> <tr> <td>0.100</td> <td><u>G2879</u></td> <td><u>18070</u></td> <td><u>2-26-20</u></td> </tr> <tr> <td>0.200</td> <td><u>G3709</u></td> <td><u>17340</u></td> <td><u>10-9-19</u></td> </tr> <tr> <td>0.300</td> <td><u>G8149</u></td> <td><u>18110</u></td> <td><u>4-2-20</u></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td><u>1781708A2</u></td> <td><u>8/5/19</u></td> </tr> </tbody> </table>	Simulator	Serial Number	Lot Number	Expiration	0.000	<u>G8144</u>	N/A	N/A	0.040	<u>G2403</u>	<u>17410</u>	<u>12-6-19</u>	0.100	<u>G2879</u>	<u>18070</u>	<u>2-26-20</u>	0.200	<u>G3709</u>	<u>17340</u>	<u>10-9-19</u>	0.300	<u>G8149</u>	<u>18110</u>	<u>4-2-20</u>	0.080 DGS	N/A	<u>1781708A2</u>	<u>8/5/19</u>
Simulator	Serial Number	Lot Number	Expiration																									
0.000	<u>G8144</u>	N/A	N/A																									
0.040	<u>G2403</u>	<u>17410</u>	<u>12-6-19</u>																									
0.100	<u>G2879</u>	<u>18070</u>	<u>2-26-20</u>																									
0.200	<u>G3709</u>	<u>17340</u>	<u>10-9-19</u>																									
0.300	<u>G8149</u>	<u>18110</u>	<u>4-2-20</u>																									
0.080 DGS	N/A	<u>1781708A2</u>	<u>8/5/19</u>																									
<input checked="" type="checkbox"/> Post Calibration Adjustment Stability Checks																												
<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td><u>SD1021</u></td> <td><u>201707D</u></td> <td><u>7-25-19</u></td> </tr> <tr> <td>0.080</td> <td><u>DR1275</u></td> <td><u>201707E</u></td> <td><u>7-25-19</u></td> </tr> <tr> <td>0.200</td> <td><u>SD1019</u></td> <td><u>201707C</u></td> <td><u>7-24-19</u></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td><u>AG805702</u></td> <td><u>2-26-20</u></td> </tr> </tbody> </table>	Simulator	Serial Number	Lot Number	Expiration	0.050	<u>SD1021</u>	<u>201707D</u>	<u>7-25-19</u>	0.080	<u>DR1275</u>	<u>201707E</u>	<u>7-25-19</u>	0.200	<u>SD1019</u>	<u>201707C</u>	<u>7-24-19</u>	0.080 DGS	N/A	<u>AG805702</u>	<u>2-26-20</u>								
Simulator	Serial Number	Lot Number	Expiration																									
0.050	<u>SD1021</u>	<u>201707D</u>	<u>7-25-19</u>																									
0.080	<u>DR1275</u>	<u>201707E</u>	<u>7-25-19</u>																									
0.200	<u>SD1019</u>	<u>201707C</u>	<u>7-24-19</u>																									
0.080 DGS	N/A	<u>AG805702</u>	<u>2-26-20</u>																									

Department Inspection Performed By <u>SP</u>												
Barometric Pressure ID# <u>28662</u>												
Gauge <u>1014</u> Instrument <u>1013</u>												
Mouth Alcohol Solution Lot # <u>2016-C</u>												
Acetone Stock Solution Lot # <u>2018-A</u>												
<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td><u>G4444</u></td> </tr> <tr> <td>Interferent</td> <td><u>G6621</u></td> </tr> <tr> <td>0.050</td> <td><u>SD1018</u></td> </tr> <tr> <td>0.080</td> <td><u>SD3962</u></td> </tr> <tr> <td>0.200</td> <td><u>G2078</u></td> </tr> </tbody> </table>	Simulator	Serial Number	0.000	<u>G4444</u>	Interferent	<u>G6621</u>	0.050	<u>SD1018</u>	0.080	<u>SD3962</u>	0.200	<u>G2078</u>
Simulator	Serial Number											
0.000	<u>G4444</u>											
Interferent	<u>G6621</u>											
0.050	<u>SD1018</u>											
0.080	<u>SD3962</u>											
0.200	<u>G2078</u>											
Attachments												
<input checked="" type="checkbox"/> Form 41 <input checked="" type="checkbox"/> Stability Checks <input checked="" type="checkbox"/> Calibration Certificate <input checked="" type="checkbox"/> Calibration Adjustment	<input checked="" type="checkbox"/> Post-Stability Checks <input type="checkbox"/> Flow Calibration <input type="checkbox"/> Form 40 <input type="checkbox"/> Other _____											

Notes/Suggested Service: _____

Instrument Complies with Chapter 11D-8, FAC

Instrument Does Not Comply with Chapter 11D-8, FAC

Return to/Place into Evidentiary Use

Remain Out of Evidentiary Use

Conduct an Agency Inspection Before Evidentiary Use

6/7/18 JJ Drake 6/7/18
 Tech Review / Date Admin Review / Date

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: HOLMES BEACH PD
Time of Inspection: 10:21

Date of Inspection: 06/07/2018

Serial Number: 80-006692
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:201707D Exp: 07/25/2019	0.08g/210L Test (g/210L) Lot#:201707E Exp: 07/25/2019	0.20g/210L Test (g/210L) Lot#:201707C Exp: 07/24/2019	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG805702 Exp: 02/26/2020
0.000	0.048	0.080	0.198	0.078
0.000	0.049	0.081	0.200	0.078
0.000	0.049	0.081	0.200	0.078
0.000	0.049	0.081	0.201	0.078
0.000	0.049	0.081	0.200	0.078
0.000	0.050	0.081	0.201	0.078
0.000	0.050	0.082	0.201	0.078
0.000	0.050	0.081	0.201	0.078
0.000	0.050	0.082	0.201	0.078
0.000	0.050	0.082	0.200	0.078

Standard Deviations	0.0006	0.0006	0.0009	0.0000
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0005 Number of Simulators Used: 5

Remarks:

The above instrument complies () does not comply () with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Shayla Platt

SHAYLA D PLATT

Signature and Printed Name

06/07/2018
Date

*6/7/18
JIA*

STABILITY CHECKS #80-006692

HOLMES BEACH PD
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006692
06/01/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:48
Control Test	0.046	08:49
Air Blank	0.000	08:49
Control Test	0.046	08:50
Air Blank	0.000	08:50
Control Test	0.047	08:51
Air Blank	0.000	08:51
Control Test Stats		
Average	0.0463	
Std Dev	0.0006	
Rel Std Dev(%)	1.2461	

SJP

Operator's Signature

8/2/18
SJP

HOLMES BEACH PD
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006692
06/01/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:01
Control Test	0.079	09:02
Air Blank	0.000	09:02
Control Test	0.079	09:03
Air Blank	0.000	09:04
Control Test	0.079	09:04
Air Blank	0.000	09:05
Control Test Stats		
Average	0.0790	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

SJP

Operator's Signature

HOLMES BEACH PD
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006692
06/01/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:53
Control Test	0.193	08:53
Air Blank	0.000	08:54
Control Test	0.195	08:54
Air Blank	0.000	08:55
Control Test	0.195	08:56
Air Blank	0.000	08:56
Control Test Stats		
Average	0.1943	
Std Dev	0.0012	
Rel Std Dev(%)	0.5942	

SJP

Operator's Signature

HOLMES BEACH PD
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006692
06/01/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:57
Control Test	0.077	08:58
Air Blank	0.000	08:58
Control Test	0.078	08:58
Air Blank	0.000	08:59
Control Test	0.076	08:59
Air Blank	0.000	09:00
Control Test Stats		
Average	0.0770	
Std Dev	0.0010	
Rel Std Dev(%)	1.2987	

DGS

SJP

Operator's Signature



Florida Department of Law Enforcement
 Alcohol Testing Program
 2729 Fort Knox Blvd.
 Bldg. 2, Suite 1300
 Tallahassee, FL 32308

Calibration Certificate

This is to certify the calibration of Intoxilyzer 8000 serial number 80-006692, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-006692</u>	UNCERTAINTY* ±
Owning Agency:	<u>HOLMES BEACH PD</u>	0.050 g/ 210 L 0.004
Calibration Date:	<u>06/07/2018</u>	0.080 g/ 210 L 0.005
Calibration Time:	<u>10:21</u>	0.200 g/ 210 L 0.008
		0.080 g/ 210 L Dry Gas Control 0.005


All results are reported in g/ 210 L.
 Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.
 *Uncertainty is based on fleet-wide data and is expressed to a 99% level of confidence (k=3).

TRACEABILITY INFORMATION
 This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Thermometer temperatures are checked with NIST traceable Eutechnics 4400 digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the uses of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.


06/07/2018 SHAYLA D PLATT,
 Date Department Inspector

Handwritten initials: 6/7/18 JP

HOLMES BENCH PD
Intoxilyzer - Alconco Analyzer
Model 8000
46/06/2018
10:46:19
SN 80-006692

Auto Calibration
Max Power Res Value = 53
Auto Range Res Value = 30

Sol Value = 0.000 g/210L ***
Fit value = 0.0000 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12716, Sum Io = 12996
<<<<< CHANNEL 1 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.0840 (0.0070)
Sample #2 = 0.0710 (0.0660)
Sample #3 = 0.0870 (0.1330)
Sample #4 = 0.0830 (0.1760)
Avg % Abs = 0.0803 (0.1317)
STD DEV = 0.0083 (0.0451)
REL STD DEV = 10.365 (34.188)

<<<<< CHANNEL 2 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.1510 (0.0020)
Sample #2 = 0.1460 (0.0300)
Sample #3 = 0.1310 (0.0500)
Sample #4 = 0.1490 (0.0660)
Avg % Abs = 0.1420 (0.0487)
STD DEV = 0.0096 (0.0190)
REL STD DEV = 6.791 (37.062)

Sol Value = 0.040 g/210L ***
Fit value = 0.1905 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12688, Sum Io = 12984
<<<<< CHANNEL 1 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.7640 (0.0080)
Sample #2 = 0.7550 (0.0480)
Sample #3 = 0.7650 (0.0700)
Sample #4 = 0.7720 (0.0940)
Avg % Abs = 0.7640 (0.0707)
STD DEV = 0.0085 (0.0230)
REL STD DEV = 1.116 (32.557)

6/7/18
JD

<<<<< CHANNEL 2 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 1.5540 (0.0090)
Sample #2 = 1.5740 (0.0200)
Sample #3 = 1.5640 (0.0220)
Sample #4 = 1.5950 (0.0300)
Avg % Abs = 1.5777 (0.0240)
STD DEV = 0.0158 (0.0053)
REL STD DEV = 1.003 (22.048)

Sol Value = 0.100 g/210L ***
Fit value = 0.4762 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12644, Sum Io = 12959
<<<<< CHANNEL 1 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 1.7880 (0.0050)
Sample #2 = 1.7800 (0.0230)
Sample #3 = 1.7870 (0.0320)
Sample #4 = 1.7750 (0.0510)
Avg % Abs = 1.7807 (0.0353)
STD DEV = 0.0060 (0.0143)
REL STD DEV = 0.339 (40.456)

<<<<< CHANNEL 2 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 3.6010 (-0.0010)
Sample #2 = 3.6390 (-0.0050)
Sample #3 = 3.6470 (0.0020)
Sample #4 = 3.6410 (0.0010)
Avg % Abs = 3.6423 (-0.0007)
STD DEV = 0.0042 (0.0038)
REL STD DEV = 0.114 (567.891)

Sol Value = 0.200 g/210L ***
Fit value = 0.9524 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12625, Sum Io = 12949
<<<<< CHANNEL 1 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 3.3740 (0.0000)
Sample #2 = 3.4060 (0.0150)
Sample #3 = 3.3870 (0.0340)
Sample #4 = 3.3900 (0.0470)
Avg % Abs = 3.3943 (0.0320)
STD DEV = 0.0102 (0.0161)
REL STD DEV = 0.301 (50.292)

<<<<< CHANNEL 2 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 6.8070 (-0.0120)
Sample #2 = 6.8420 (0.0000)
Sample #3 = 6.8710 (-0.0060)
Sample #4 = 6.8550 (-0.0060)
Avg % Abs = 6.8560 (-0.0040)
STD DEV = 0.0145 (0.0035)
REL STD DEV = 0.212 (86.603)

Sol Value = 0.300 g/210L ***
Fit value = 1.4286 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12619, Sum Io = 12950
<<<<< CHANNEL 1 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 4.9690 (0.0070)
Sample #2 = 4.9960 (0.0280)
Sample #3 = 5.0110 (0.0360)
Sample #4 = 5.0000 (0.0470)
Avg % Abs = 5.0023 (0.0370)
STD DEV = 0.0078 (0.0095)
REL STD DEV = 0.155 (25.782)

<<<<< CHANNEL 2 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 9.9150 (0.0000)
Sample #2 = 9.9630 (0.0130)
Sample #3 = 9.9820 (0.0220)
Sample #4 = 9.9760 (0.0160)
Avg % Abs = 9.9737 (0.0170)
STD DEV = 0.0097 (0.0046)
REL STD DEV = 0.097 (26.956)

***** AUTO CAL DATA *****
<<<<< CHANNEL 1 >>>>>
Sol Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.080
Std Dev = 0.01 Rel Std Dev = 10.37
Sol Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 0.764
Std Dev = 0.01 Rel Std Dev = 1.12
Sol Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 1.781
Std Dev = 0.01 Rel Std Dev = 0.34
Sol Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 3.394
Std Dev = 0.01 Rel Std Dev = 0.30
Sol Val = 1.4286 mg/l or 0.300 g/210L
% Abs = 5.002
Std Dev = 0.01 Rel Std Dev = 0.16
Zero Order Coef = -233.88
First Order Coef = 2779.77
Second Order Coef = 25.00
Standard Deviation = 25.080271

<<<<< CHANNEL 2 >>>>>
Sol Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.142
Std Dev = 0.01 Rel Std Dev = 6.79
Sol Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 1.578
Std Dev = 0.02 Rel Std Dev = 1.00
Sol Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 3.642
Std Dev = 0.00 Rel Std Dev = 0.11
Sol Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 6.856
Std Dev = 0.01 Rel Std Dev = 0.21
Sol Val = 1.4286 mg/l or 0.300 g/210L
% Abs = 9.974
Std Dev = 0.01 Rel Std Dev = 0.10
Zero Order Coef = -205.05
First Order Coef = 1324.10
Second Order Coef = 13.05
Standard Deviation = 26.64592

CAL ADJUSTMENT
80-006692 SP

Solution Stats Quadratic Fit Chan 1
Act Fit Residual
g/210L g/210L g/210L
0.000 -0.000 0.0002
0.040 0.040 0.0000
0.100 0.101 -0.0007
0.200 0.199 0.0007
0.300 0.300 -0.0002

Solution Stats Quadratic Fit Chan 2
Act Fit Residual
g/210L g/210L g/210L
0.000 -0.000 0.0004
0.040 0.040 -0.0002
0.100 0.101 -0.0006
0.200 0.199 0.0008
0.300 0.300 -0.0003

Sol Value = 0.080 g/210L ***
Fit value = 0.3810 mg/l %%%
Samples Taken = 4, Discarded = 1
***** CHANNEL 1 *****
Sample #1 = 3311.00
Sample #2 = 3311.00
Sample #3 = 3326.00
Sample #4 = 3318.00
Average Result = 3318.3333
STD DEV = 7.5056
REL STD DEV = 0.226

***** CHANNEL 2 *****
Sample #1 = 3342.00
Sample #2 = 3342.00
Sample #3 = 3355.00
Sample #4 = 3341.00
Average Result = 3346.0000
STD DEV = 7.8102
REL STD DEV = 0.233

Dry Gas H20 Adjust Results *****
Barometric Pressure = 1010
3 um H20 Adjust (mg/*10,000) = 491
9 um H20 Adjust (mg/*10,000) = 463
**** AUTO CAL PASS

POST CAL ADJUST STABILITY CHECKS #80-006692

HOLMES BEACH PD
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006692
06/06/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:55
Control Test	0.048	12:56
Air Blank	0.000	12:56
Control Test	0.049	12:57
Air Blank	0.000	12:58
Control Test	0.049	12:58
Air Blank	0.000	12:59
Control Test Stats		
Average	0.0487	
Std Dev	0.0006	
Rel Std Dev(%)	1.1863	

SP

Operator's Signature

6/7/18
JD

HOLMES BEACH PD
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006692
06/06/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:17
Control Test	0.080	13:17
Air Blank	0.000	13:18
Control Test	0.080	13:19
Air Blank	0.000	13:19
Control Test	0.080	13:20
Air Blank	0.000	13:20
Control Test Stats		
Average	0.0800	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

SP

Operator's Signature

HOLMES BEACH PD
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006692
06/06/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:46
Control Test	0.193	12:47
Air Blank	0.000	12:48
Control Test	0.197	12:48
Air Blank	0.000	12:49
Control Test	0.199	12:50
Air Blank	0.000	12:50
Control Test Stats		
Average	0.1953	
Std Dev	0.0031	
Rel Std Dev(%)	1.5561	

TIGHTENED SIM 1st
Head Between
and 2nd sample SP

SP

Operator's Signature

HOLMES BEACH PD
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006692
06/06/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:51
Control Test	0.199	12:51
Air Blank	0.000	12:52
Control Test	0.198	12:53
Air Blank	0.000	12:53
Control Test	0.198	12:54
Air Blank	0.000	12:54
Control Test Stats		
Average	0.1983	
Std Dev	0.0006	
Rel Std Dev(%)	0.2911	

SP

Operator's Signature

HOLMES BEACH PD
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006692
06/06/2018
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:43
Control Test	0.079	12:44
Air Blank	0.000	12:44
Control Test	0.079	12:44
Air Blank	0.000	12:45
Control Test	0.079	12:45
Air Blank	0.000	12:46
Control Test Stats		
Average	0.0790	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

DCS

SP

Operator's Signature