



INSTRUMENT PROCESSING SHEET

Agency Marion County

S/N 80-007171

Florida Department of Law Enforcement

Date In 3/2/2020

DI Completion Date 3/16/2020

Ship P/U H/D CMI EE

<b>Intake</b> Performed By <u>DP</u> <input checked="" type="checkbox"/> Annual <input type="checkbox"/> Registration <input type="checkbox"/> Return from CMI / EE Visual Inspection: <input checked="" type="checkbox"/> Case <input checked="" type="checkbox"/> Handle <input checked="" type="checkbox"/> Keyboard <input checked="" type="checkbox"/> Dry Gas Shelf <input checked="" type="checkbox"/> Feet <input checked="" type="checkbox"/> Breath Tube <input checked="" type="checkbox"/> Ports <input checked="" type="checkbox"/> Screws Tight Other Equipment/ Accessories: <input checked="" type="checkbox"/> Power cord <input type="checkbox"/> Printer Cable <input checked="" type="checkbox"/> Static Bag <input type="checkbox"/> 12V DC Cable Notes: _____	<b>Quality Checks</b> Performed By <u>SP</u> <input checked="" type="checkbox"/> Breath Tube Screen <input checked="" type="checkbox"/> Replace External O-Rings <input checked="" type="checkbox"/> Instrument Set Up Verified <input checked="" type="checkbox"/> R-Value <u>199</u> <input checked="" type="checkbox"/> Flow Verification (L/s) Flow Column # <u>ATP102</u> 32 mm <u>.160</u> (.139 - .169) 36 mm <u>.175</u> (.156 - .190) 53 mm <u>.242</u> (.228 - .278) 103 mm <u>.515</u> (.447 - .547) <input checked="" type="checkbox"/> Barometric Pressure Check Gauge ID # <u>30793</u> <input checked="" type="checkbox"/> Stability Checks	<b>Flow Calibration</b> Performed By _____ Flow Column # _____ <input type="checkbox"/> 5L/min - 17mm <input type="checkbox"/> 15L/min - 53mm <input type="checkbox"/> 30L/min - 103mm <input type="checkbox"/> R-Value _____ <input type="checkbox"/> Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)
---	---	--

**Final Release Date**

**FDLE**

**MAR 16 2020**

**Alcohol Testing Program**

Simulator	Serial #	Lot #/Exp
0.050	SD1021	201905A 5-14-21
0.080	SD3902	201905B 5-14-21
0.200	SD3933	201904D 4-30-21
0.080 DGS	N/A	AG931603 11-12-21

**Maintenance** Performed By \_\_\_\_\_

Battery Replacement  
 Dry Gas Regulator Replacement  
 Breath Tube Replacement  
 Other \_\_\_\_\_

**Temperature Checks** Performed By SP

Lab Temp °C 22.1  
 External Digital Therm. ID#: 300502  
 34°C +-2 Serial #: SD1021  
 34°C +-2 Serial #: DR1279 DR1275  
 34°C +-2 Serial #: SD1011

**Calibration Adjustment** Performed By SP

Barometric Pressure Gauge 1020 ID # 28421

Simulator	Serial Number	Lot Number	Expiration
0.000	G8144	N/A	N/A
0.040	G2403	19080	3-4-21
0.100	G2879	19160	7-9-21
0.200	G3709	19040	1-29-21
0.300	G8149	19010	1-3-21
0.080 DGS	N/A	08819080A1	6-5-21

Post Calibration Adjustment Stability Checks

Simulator	Serial Number	Lot Number	Expiration
0.050	SD1021	201905A	5-14-21
0.080	DR1275	201905B	5-14-21
0.200	SD1011	201904D	4-30-21
0.080 DGS	N/A	AG931603	11-12-21

**Department Inspection** Performed By SP

Barometric Pressure ID# 20932  
 Gauge 1023 Instrument 1023  
 Mouth Alcohol Solution Lot # 2019-B  
 Acetone Stock Solution Lot # 2019-A

Simulator	Serial Number
0.000	G8144
Interferent	DR3855
0.050	SD1021
0.080	DR1279 DR1275 SP
0.200	SD1011

**Attachments**

Form 41  Post-Stability Checks  
 Stability Checks  Flow Calibration  
 Calibration Certificate  Form 40  
 Calibration Adjustment  Other \_\_\_\_\_

Notes/Suggested Service: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

Instrument Complies with Chapter 11D-8, FAC  
 Instrument Does Not Comply with Chapter 11D-8, FAC

Return to/Place into Evidentiary Use  
 Remain Out of Evidentiary Use

Conduct an Agency Inspection Before Evidentiary Use

SPM 3/16/20 Scott Kirkland 3/16/2020  
 Tech Review / Date Admin Review / Date



# Calibration Certificate

Florida Department of Law Enforcement  
Alcohol Testing Program  
2729 Fort Knox Blvd.  
Bldg. 2, Suite 1300  
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-007171, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-007171</u>	UNCERTAINTY* ±	
Owning Agency:	<u>MARION COUNTY SO</u>	0.050 g/ 210 L	0.004
Calibration Date:	<u>03/16/2020</u>	0.080 g/ 210 L	0.005
Calibration Time:	<u>12:25</u>	0.200 g/ 210 L	0.007
		0.080 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.  
\*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence (k=3).

### TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Thermometer temperatures are checked with NIST traceable Eutechnics 4400 digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the uses of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

03/16/2020 Date  
*Shayla Platt*

**SHAYLA D PLATT,**  
**Department Inspector**

FDLE/ATP Form 69 January 2020  
Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality

*ADM TK  
3/16/2020*

# Florida Department of Law Enforcement Alcohol Testing Program

## DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: MARION COUNTY SO  
Time of Inspection: 12:25

Date of Inspection: 03/16/2020

Serial Number: 80-007171  
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:201905A Exp: 05/14/2021	0.08g/210L Test (g/210L) Lot#:201905B Exp: 05/14/2021	0.20g/210L Test (g/210L) Lot#:201904D Exp: 04/30/2021	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG931603 Exp: 11/12/2021
0.000	0.048	0.078	0.198	0.080
0.000	0.048	0.078	0.198	0.080
0.000	0.049	0.079	0.199	0.080
0.000	0.049	0.079	0.199	0.080
0.000	0.049	0.079	0.199	0.079
0.000	0.049	0.079	0.199	0.079
0.000	0.049	0.080	0.199	0.079
0.000	0.048	0.080	0.199	0.080
0.000	0.049	0.080	0.199	0.080
0.000	0.049	0.080	0.198	0.079

Standard Deviations	0.0004	0.0007	0.0004	0.0005
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0005 Number of Simulators Used: 5

Remarks:

JPM  
TSK  
3/16/2020

The above instrument complies (  ) does not comply (  ) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Shayla Platt

SHAYLA D PLATT

Signature and Printed Name

03/16/2020  
Date

# Stability Checks # 80-007171

MARION COUNTY SO  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-007171  
 03/11/2020  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:28
Control Test	0.046	10:29
Air Blank	0.000	10:29
Control Test	0.046	10:30
Air Blank	0.000	10:30
Control Test	0.046	10:31
Air Blank	0.000	10:32
Control Test Stats		
Average	0.0460	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

*SP*  
 -----  
 Operator's Signature

MARION COUNTY SO  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-007171  
 03/11/2020  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:23
Control Test	0.076	10:24
Air Blank	0.000	10:24
Control Test	0.075	10:25
Air Blank	0.000	10:25
Control Test	0.076	10:26
Air Blank	0.000	10:27
Control Test Stats		
Average	0.0757	
Std Dev	0.0006	
Rel Std Dev(%)	0.7630	

*SP*  
 -----  
 Operator's Signature

MARION COUNTY SO  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-007171  
 03/11/2020  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:18
Control Test	0.193	10:19
Air Blank	0.000	10:20
Control Test	0.194	10:20
Air Blank	0.000	10:21
Control Test	0.194	10:22
Air Blank	0.000	10:22
Control Test Stats		
Average	0.1937	
Std Dev	0.0006	
Rel Std Dev(%)	0.2981	

*SP*  
 -----  
 Operator's Signature

MARION COUNTY SO  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-007171  
 03/11/2020  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:36
Control Test	0.077	10:36
Air Blank	0.000	10:37
Control Test	0.078	10:37
Air Blank	0.000	10:38
Control Test	0.078	10:38
Air Blank	0.000	10:38
Control Test Stats		
Average	0.0777	
Std Dev	0.0006	
Rel Std Dev(%)	0.7434	

DGS

*SP*  
 -----  
 Operator's Signature

*OBM*  
*SK*  
*3/16/2020*

MARLON COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000  
03/13/2020

Auto Calibration

Max Power Res Value = 113  
Auto Range Res Value = 89

Sol Value = 0.000 g/210L \*\*\*  
Fit value = 0.0000 mg/l %%%  
Samples Taken = 4, Discarded = 1  
Sum Io = 12584, Sum Io = 13118  
<<<<< CHANNEL 1 >>>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = -0.0160 (-0.0150)  
Sample #2 = -0.0130 (-0.0100)  
Sample #3 = -0.0150 (0.0000)  
Sample #4 = -0.0270 (0.0020)  
Avg % Abs = -0.0183 (0.0007)  
STD DEV = 0.0076 (0.0110)  
REL STD DEV = 41.301 (1652.271)

<<<<< CHANNEL 2 >>>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 1.3950 (-0.0070)  
Sample #2 = 1.3990 (-0.0030)  
Sample #3 = 1.3980 (0.0000)  
Sample #4 = 1.4020 (-0.0060)  
Avg % Abs = 1.3963 (-0.0030)  
STD DEV = 0.0074 (0.0030)  
REL STD DEV = 0.528 (100.000)

Sol Value = 0.100 g/210L \*\*\*  
Fit value = 0.4762 mg/l %%%  
Samples Taken = 4, Discarded = 1  
Sum Io = 12563, Sum Io = 13108  
<<<<< CHANNEL 1 >>>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 1.7920 (0.0060)  
Sample #2 = 1.8140 (0.0000)  
Sample #3 = 1.8000 (0.0300)  
Sample #4 = 1.7800 (0.0340)  
Avg % Abs = 1.7980 (0.0213)  
STD DEV = 0.0171 (0.0186)  
REL STD DEV = 0.950 (87.108)

<<<<< CHANNEL 2 >>>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 3.3260 (0.0020)  
Sample #2 = 3.3070 (0.0190)  
Sample #3 = 3.3290 (0.0170)  
Sample #4 = 3.3180 (0.0150)  
Avg % Abs = 3.3180 (0.0170)  
STD DEV = 0.0110 (0.0020)  
REL STD DEV = 0.332 (11.765)

Sol Value = 0.200 g/210L \*\*\*  
Fit value = 0.9524 mg/l %%%  
Samples Taken = 4, Discarded = 1  
Sum Io = 12554, Sum Io = 13103  
<<<<< CHANNEL 1 >>>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 3.5030 (-0.0150)  
Sample #2 = 3.5130 (-0.0040)  
Sample #3 = 3.5520 (-0.0080)  
Sample #4 = 3.5190 (0.0030)  
Avg % Abs = 3.5280 (-0.0030)  
STD DEV = 0.0210 (0.0056)  
REL STD DEV = 0.595 (185.592)

Sol Value = 0.040 g/210L \*\*\*  
Fit value = 0.1905 mg/l %%%  
Samples Taken = 4, Discarded = 1  
Sum Io = 12567, Sum Io = 13108  
<<<<< CHANNEL 1 >>>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 0.7460 (-0.0140)  
Sample #2 = 0.7830 (0.0070)  
Sample #3 = 0.7320 (0.0080)  
Sample #4 = 0.7460 (-0.0040)  
Avg % Abs = 0.7270 (0.0037)  
STD DEV = 0.0219 (0.0067)  
REL STD DEV = 3.017 (181.591)

<<<<< CHANNEL 2 >>>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 6.3440 (-0.0050)  
Sample #2 = 6.3630 (-0.0060)  
Sample #3 = 6.4070 (-0.0070)  
Sample #4 = 6.4070 (-0.0040)  
Avg % Abs = 6.3923 (-0.0057)  
STD DEV = 0.0254 (0.0015)  
REL STD DEV = 0.397 (26.956)

Sol Value = 0.300 g/210L \*\*\*  
Fit value = 1.4286 mg/l %%%  
Samples Taken = 4, Discarded = 1  
Sum Io = 12554, Sum Io = 13105  
<<<<< CHANNEL 1 >>>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 5.1950 (-0.0020)  
Sample #2 = 5.2000 (0.0250)  
Sample #3 = 5.1890 (0.0200)  
Sample #4 = 5.2130 (0.0120)  
Avg % Abs = 5.2007 (0.0190)  
STD DEV = 0.0120 (0.0066)  
REL STD DEV = 0.231 (34.513)

<<<<< CHANNEL 2 >>>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 9.2640 (-0.0010)  
Sample #2 = 9.2740 (0.0370)  
Sample #3 = 9.2770 (0.0320)  
Sample #4 = 9.2640 (0.0300)  
Avg % Abs = 9.2717 (0.0330)  
STD DEV = 0.0068 (0.0036)  
REL STD DEV = 0.073 (10.926)

Sol Value = 0.000 mg/l or 0.000 g/210L  
Fit value = 0.3810 mg/l %%%  
Samples Taken = 4, Discarded = 1  
<<<<< CHANNEL 1 >>>>>  
Sample #1 = 3567.00  
Sample #2 = 3523.00  
Sample #3 = 3560.00  
Sample #4 = 3553.00  
Average Result = 3545.3333  
STD DEV = 19.6554  
REL STD DEV = 0.554  
<<<<< CHANNEL 2 >>>>>  
Sample #1 = 3470.00  
Sample #2 = 3443.00  
Sample #3 = 3462.00  
Sample #4 = 3460.00  
Average Result = 3455.0000  
STD DEV = 10.4403  
REL STD DEV = 0.302  
\*\*\*\*\*  
Dry Gas H2O Adjust Results \*\*\*\*\*  
Barometric Pressure = 1020  
3 um H2O Adjust (mg/l x 10,000) = 264  
9 um H2O Adjust (mg/l x 10,000) = 354  
\*\*\*\*\* AUTO CAL PASS

CAL ADJUSTMENT  
# 80-007171 SP

BBM  
3/16/2020

\*\*\*\*\* AUTO CAL DATA \*\*\*\*\*  
<<<<< CHANNEL 1 >>>>>  
Sol Val = 0.0000 mg/l or 0.000 g/210L  
% Abs = -0.018  
Std Dev = 0.01 Rel Std Dev = 41.30  
Sol Val = 0.1905 mg/l or 0.040 g/210L  
% Abs = 0.727  
Std Dev = 0.02 Rel Std Dev = 3.02  
Sol Val = 0.4762 mg/l or 0.100 g/210L  
% Abs = 1.798  
Std Dev = 0.02 Rel Std Dev = 0.95  
Sol Val = 0.9524 mg/l or 0.200 g/210L  
% Abs = 3.528  
Std Dev = 0.02 Rel Std Dev = 0.60  
Sol Val = 1.4286 mg/l or 0.300 g/210L  
% Abs = 5.201  
Std Dev = 0.01 Rel Std Dev = 0.23  
Zero Order Coef = 35.19  
First Order Coef = 2571.01  
Second Order Coef = 32.68  
Standard Deviation = 11.992397

<<<<< CHANNEL 2 >>>>>  
Sol Val = 0.0000 mg/l or 0.000 g/210L  
% Abs = 0.053  
Std Dev = 0.01 Rel Std Dev = 11.45  
Sol Val = 0.1905 mg/l or 0.040 g/210L  
% Abs = 1.396  
Std Dev = 0.01 Rel Std Dev = 0.53  
Sol Val = 0.4762 mg/l or 0.100 g/210L  
% Abs = 3.318  
Std Dev = 0.01 Rel Std Dev = 0.33  
Sol Val = 0.9524 mg/l or 0.200 g/210L  
% Abs = 6.392  
Std Dev = 0.03 Rel Std Dev = 0.40  
Sol Val = 1.4286 mg/l or 0.300 g/210L  
% Abs = 9.272  
Std Dev = 0.01 Rel Std Dev = 0.07  
Zero Order Coef = -75.37  
First Order Coef = 1401.60  
Second Order Coef = 15.85  
Standard Deviation = 8.502403

<<<<< CHANNEL 1 >>>>>  
Sample #1 = 3567.00  
Sample #2 = 3523.00  
Sample #3 = 3560.00  
Sample #4 = 3553.00  
Average Result = 3545.3333  
STD DEV = 19.6554  
REL STD DEV = 0.554  
\*\*\*\*\*  
<<<<< CHANNEL 2 >>>>>  
Sample #1 = 3470.00  
Sample #2 = 3443.00  
Sample #3 = 3462.00  
Sample #4 = 3460.00  
Average Result = 3455.0000  
STD DEV = 10.4403  
REL STD DEV = 0.302  
\*\*\*\*\*  
Dry Gas H2O Adjust Results \*\*\*\*\*  
Barometric Pressure = 1020  
3 um H2O Adjust (mg/l x 10,000) = 264  
9 um H2O Adjust (mg/l x 10,000) = 354  
\*\*\*\*\* AUTO CAL PASS

Solution Stats Quadratic Fit Chan 1  
Act Fit Residual  
g/210L g/210L g/210L  
0.000 -0.000 0.0003  
0.040 0.040 -0.0004  
0.100 0.100 -0.0000  
0.200 0.200 0.0002  
0.300 0.300 -0.0001

Solution Stats Quadratic Fit Chan 2  
Act Fit Residual  
g/210L g/210L g/210L  
0.000 -0.000 0.0000  
0.040 0.040 -0.0002  
0.100 0.100 0.0003  
0.200 0.200 -0.0002  
0.300 0.300 0.0001

Sol Value = 0.080 g/210L \*\*\*  
Fit value = 0.3810 mg/l %%%  
Samples Taken = 4, Discarded = 1  
\*\*\*\*\* CHANNEL 1

Sol Value = 0.080 g/210L \*\*\*  
Fit value = 0.3810 mg/l %%%  
Samples Taken = 4, Discarded = 1  
\*\*\*\*\* CHANNEL 1  
Sample #1 = 3567.00  
Sample #2 = 3523.00  
Sample #3 = 3560.00  
Sample #4 = 3553.00  
Average Result = 3545.3333  
STD DEV = 19.6554  
REL STD DEV = 0.554  
\*\*\*\*\*  
\*\*\*\*\* CHANNEL 2

Sol Value = 0.080 g/210L \*\*\*  
Fit value = 0.3810 mg/l %%%  
Samples Taken = 4, Discarded = 1  
\*\*\*\*\* CHANNEL 1  
Sample #1 = 3567.00  
Sample #2 = 3523.00  
Sample #3 = 3560.00  
Sample #4 = 3553.00  
Average Result = 3545.3333  
STD DEV = 19.6554  
REL STD DEV = 0.554  
\*\*\*\*\*  
\*\*\*\*\* CHANNEL 2  
Sample #1 = 3470.00  
Sample #2 = 3443.00  
Sample #3 = 3462.00  
Sample #4 = 3460.00  
Average Result = 3455.0000  
STD DEV = 10.4403  
REL STD DEV = 0.302  
\*\*\*\*\*  
Dry Gas H2O Adjust Results \*\*\*\*\*  
Barometric Pressure = 1020  
3 um H2O Adjust (mg/l x 10,000) = 264  
9 um H2O Adjust (mg/l x 10,000) = 354  
\*\*\*\*\* AUTO CAL PASS

# Post Cal - Adjust Stability Checks # 80-007171

MARION COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-007171  
03/16/2020  
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:49
Control Test	0.049	09:50
Air Blank	0.000	09:50
Control Test	0.049	09:51
Air Blank	0.000	09:52
Control Test	0.049	09:52
Air Blank	0.000	09:53
Control Test Stats		
Average	0.0490	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

SP

Operator's Signature

MARION COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-007171  
03/16/2020  
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:54
Control Test	0.079	09:55
Air Blank	0.000	09:55
Control Test	0.079	09:56
Air Blank	0.000	09:56
Control Test	0.079	09:57
Air Blank	0.000	09:58
Control Test Stats		
Average	0.0790	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

SP

Operator's Signature

MARION COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-007171  
03/16/2020  
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:45
Control Test	0.197	09:46
Air Blank	0.000	09:46
Control Test	0.197	09:47
Air Blank	0.000	09:47
Control Test	0.198	09:48
Air Blank	0.000	09:49
Control Test Stats		
Average	0.1973	
Std Dev	0.0006	
Rel Std Dev(%)	0.2926	

SP

Operator's Signature

MARION COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-007171  
03/16/2020  
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:06
Control Test	0.081	10:06
Air Blank	0.000	10:07
Control Test	0.080	10:07
Air Blank	0.000	10:08
Control Test	0.080	10:08
Air Blank	0.000	10:08
Control Test Stats		
Average	0.0803	
Std Dev	0.0006	
Rel Std Dev(%)	0.7187	

DAS

SP

Operator's Signature

DDM

TSK

3/16/2020

INTOXYLIZER 80-007172

PLEASE ADJUST PRINTER AS IT WILL NOT FEED CORRECTLY THROUGH THE COVER.

*WRB*