

INSTRUMENT PROCESSING SHEET

Agency St. Petersburg PD

S/N 80-001653

lorida Dep aw Enforce				DI Completion	Date 10/17/2022	Ship	<u>0</u> P/U			DEE
Intake	ByT	DG Qu	ality Che	cks By TDG	Date 9/26/2022	Flow Calib	oration B	y	Date_	
Annual	1	< t 📄	Breath T	ube Screen		Flow Colu	mn #			
🛛 Registrati	on		Replace External O-Rings				/min – 17	mm		
🔳 Return fro	om CMI / EE		Instrument Set Up Verified				L/min – 5	3mm		
			R-Value	210		30	L/min – 1	03mm		
Visual Inspec			low Ver	ification (L/s)		R-Value	3			
Case Handle Keyboard Dry Gas Shelf			w Colum	n # ATP104		D Post Ca			ion (L/s)	
			2 mm 0	.148	(.139169)					
Feet	Breath Tub		6 mm 0	.167	(.156190)	Flow Colu 32 mm			(.139	169
Ports	Screws Tig	^{nt} 5	3 mm 0	.234	(.228278)					
Other Equip	ment/ Accessories:	10	3 mm 0	.503	(.447547)					
Power co	rd 🛛 🖬 Printer Cal			ric Pressure Ch		103 mm				
Static Bag	🗖 12V DC Ca		uge ID #_							
Notes:			Stability							
			nulator		Lot #/Exp	Malat			0	
_		_			The second s	Maintena			Ву	
		0.0	050	MP6286	202201C	Battery				
				101 0200	01/11/2024	Dry Gas		50 C		
		0.0	080	MP6287	202201D	Breath			It	
				IVIP0207	01/18/2024	Other_				
		0.2	200	MP6288	202201E					
					01/18/2024					
		0.0	80 DGS	N/A	00521080A2					
4					02/05/2023					
Calibration /	Adjustment		a second	ByTDG	Department Inspec	tion		97 2 C	By TD	G
	Pressure Gauge 10	13 11	# 2866		Barometric Pressure	and the second se	9			
Simulator		Lot #	E	Expiration	Gauge <u>1013</u>	In	strument	1012		
0.000	MP5099	N/A		N/A	Mouth Alcohol Solu					
0.040	MP5096	2107) 0	3/01/2023	Acetone Stock Solut	tion Lot # 2	021-C			
0.100	MP5098	2138		9/13/2023	Simulator		Serial N	lumber	1.212	
0.200	MP5100	2051		2/03/2022	0.000 🐒			MP62	284	
0.300	MP5101	21420		0/20/2023	Interferent		2	MP62		
0.080 DGS	N/A				0.050			MP62		
		AG1159		6/08/2023	0.080			MP62		
Post Calib	ration Adjustment	Stability Che	cks				1	MP62	.00	
Simulator	Serial #	Lot #	E	xpiration	Attachments	Section and		Carlo Core		
0.050	MP6286	202201	C 0	1/11/2024	Form 41			-Stability		
0.080	MP6287	202201	D 0	1/18/2024	Stability Checks			v Calibra		
0.200	MP6288	202201		1/18/2024	Calibration Cert		Forr	n 40 (x2)		
0.080 DGS	- N/A	0052108	22-22	2/05/2023	Calibration Adju	ustment	Oth Oth	er		
cal adjust t	ested Service: <u>Perf</u> o bring values clo ed due to Hurricar	ser to nom	nal. Pro		 Instrument Cor Instrument Doc Return to/Place Remain Out of Conduct an Age 	es Not Comp e into Evide Evidentiary	oly with C ntiary Us Use	Chapter : e	l1D-8, F/	
	(<u>*</u>)			<u> </u>	DERR Digitally signed Date: 2022.10.11 10:04:07 - 04'00'	by DERR	el Soto	Digitally	signed by Israe 22.10.18 10:45:	el Soto
	2		I		Tech Review / Da	ite	Admin	Review	/ Date	

FDLE/ATP Form 48 January 2022 Issuing Authority: Alcohol Testing Program PRINTED COPIES UNCONTROLLED For Internal ATP Use ONLY

Florida Department of Law Enforcement Alcohol Testing Program

AGENCY INSPECTION REPORT - INTOXILYZER 8000

Agency: ST PETERSBURG PD Time of Inspection: 11:56

Date of Inspection: 09/26/2022

Serial Number: 80-001653 Software: 8100.27

Check or Test	YES	NO
Date and/or Time Adjusted		No
Diagnostic Check (Pre-Inspection): OK		No
Alcohol Free Subject Test: 0.000		No
Mouth Alcohol Test: Slope Not Met		No
Interferent Detect Test: Interferent Detect		No
Diagnostic Check (Post-Inspection): OK	а. С	No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#: Exp:	0.08g/210L Test (g/210L) Lot#: Exp:	0.20g/210L Test (g/210L) Lot#: Exp:	0.08 g/210L Dry Gas Std Test (g/210L) Lot#: Exp:

Number of Simulators Used:

Remarks:

BYPASSED AI. COMPLIANCE NOT DETERMINED.

The above instrument complies (X) does not comply (

) with Chapter 11D-8, FAC.

I certify that I hold a valid Florida Department of Law Enforcement Agency Inspector Permit and that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Jayor Kutalan	TAYLOR D GUTSCHOW	
and Contract	Signature and Printed Name	
	09/26/2022	

09/26/2022 Date

Florida Department of Law Enforcement Alcohol Testing Program

AGENCY INSPECTION REPORT - INTOXILYZER 8000

Agency: ST PETERSBURG PD Time of Inspection: 09:24

Date of Inspection: 10/17/2022

Serial Number: 80-001653 Software: 8100.27

10/17/22

Check or Test	YES	NO
Date and/or Time Adjusted		No
Diagnostic Check (Pre-Inspection): OK	х.	No
Alcohol Free Subject Test: 0.000		No
Mouth Alcohol Test: Slope Not Met		No
Interferent Detect Test: Interferent Detect		No
Diagnostic Check (Post-Inspection): OK		No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#: Exp:	0.08g/210L Test (g/210L) Lot#: Exp:	0.20g/210L Test (g/210L) Lot#: Exp:	0.08 g/210L Dry Gas Std Test (g/210L) Lot#: Exp:
12		Z [*]		

Number of Simulators Used:

Remarks:

BYPASSED AI TO OPERATE. AI NOT CONDUCTED. Processing delay due to Hurricane Ian.

Not determined m6-10/17/2022

\$

The above instrument complies (X) does not comply (

) with Chapter 11D-8, FAC.

I	certify	that	I	hold	a	valid	Florida	Departmen	nt of	Law	Enforcement	Agency	Inspector	Permit	and	that	I
pe	rformed	this	ing	pecti	on	in acc	cordance	with the	provis	sions	of Chapter TA	11D-8, I	FAC.				
		ybr	X	hu	th	The	~				ТА	YLOR D C	GUTSCHOW				

Signature and Printed Name

10/17/2022 Date

Type of Test	Serial Number	Agency	Date	36		Perform	ned By
Stabilities	80-00 1653	St. Petersburg PD	09/	24	2022	TDG	MG
2							

0.05g/210L	0.08g/210L	0.20g/210L	DGS 0.08g/210L
0.047 to 0.053	0.077 to 0.083	0.194 to 0.206	0.077 to 0.083 ✓ ≤0.003 of Wet ✓
ST PETERSBURG PD intoxilyzer - Alconol Analyzer Model 8000 SM 80-001653 09/25/2022	ST PETERSBURG PO intoxilyzar - 41contmätyzen Model 8000 SN 60-001653 09/26/2022 Software: 8100.27	ST FETERSSURG FD Intoxilyzer - Piconol Analyzar Model 8000 SN 80-001653 09/26/2022 Software: 8100.27	ST PETERSSURG PD Intoxilyzer - Piconol Analyzer Model 8000 SN 80-001653 09/26/2022 Software: €100.27
Software: 8190.27	Test g/210L Time	Test g/210L Time	Test g/210_ Time
Test g/210L Time Air Blank 0.000 13:18 Control Test 0.048 13:19 Air Blank 0.000 13:19 Control Test 0.048 13:20 Control Test 0.048 13:20 Control Test 0.049 13:20 Control Test 0.049 13:21 Air Blank 0.000 13:22 Control Test 1.049 13:22 Control Test 0.0480 13:22 Control Test Stats Auerage 0.0480 Std Deu 9.1090 Rel Std Deu(%)	Air Blank 0.000 13:32 Control Test 0.077 13:33 Air Blank 0.000 13:33 Control Test 0.077 13:34 Air Blank 0.000 13:33 Control Test 0.077 13:34 Air Blank 0.000 13:34 Control Test 0.077 13:35 Air Blank 0.000 13:35 Control Test 0.077 13:35 Air Blank 0.000 13:35 Control Test 0.0770 13:36 Control Test Stats Auenage 0.0000 Rel Std Deu(%) 0.0000 13:36	Air Blank 1.000 13:39 Control Test 0.197 13:40 Aic Blank 0.000 13:40 Control Test 0.196 13:41 Air Blank 0.000 13:42 Control Test 0.196 13:42 Control Test 0.196 13:42 Control Test 0.196 13:43 Control Test 0.196 13:43 Control Test 0.196 13:43 Control Test 0.200 13:43 Control Test 0.200 13:43 Control Test 0.2005 Rel Std Deu(%) 0.2941	Air Blank 0.000 13:23 Control Test 0.079 13:24 Air Blank 0.000 13:24 Control Test 0.080 13:24 Air Blank 0.100 13:25 Control Test 0.080 13:25 Control Test 0.080 13:25 Control Test 0.080 13:26 Control Test Stats Auerage 0.0797 Std Dev 0.085 Rel Std Dev(%) 0.7247
Openatoris Signature	MU . Openator 's Signature	Cpenetorije Styreture	Coenatoris Signature
		8	

ST PETERSBURG PD Intoxilyzer - Alconol Analyzer Model 8000 SN 80-001653 10/17/2022 05:27:10	<pre></pre>	<pre></pre>	<pre>xxxxx AUTO CAL DATA ***** <<<<< CHANNEL 1 >>>> Sol Val = 0.0000 mg/l or 0.000 g/210L % Abs = 0.129 Std Deu = 0.01 Rel Std Deu = 5.29 Sol Val = 0.1905 mg/l or 0.040 g/210L % Abs = 0.840 Std Deu = 0.02 Rel Std Deu = 1.86 Sol Val = 0.4762 mg/l or 0.100 g/210L % Abs = 1.945</pre>	Solution Stats Quadratic Fit Chan 2 Act Fit Residual g/210L g/210L g/210L 0.000 0.000 -0.0001 0.040 0.040 0.0001 0.100 0.100 0.0001 0.200 0.200 -0.0001 0.300 0.300 0.0001 Sol Ualue = 0.080 g/210L ***
Auto Calibration Max Power Res Ualue = 99 Auto Range Res Ualue = 82 Sol Ualue = 0.000 g/210L *** Fit ualue = 0.000 mg/1 %%% Samples TaKen = 4, Discarded = 1 Jum Io = 12562, 9um Io = 13896	501 Ualue = 0.100 g/210L *** Fit Ualue = 0.4762 mg/l %%% Samples Taken = 4, Discarded = 1 3um Io = 12527, 9um Io = 13877 <<<<<< CHANNEL 1 >>>>> Sample % Abs (% Abs Ref) Sample #1 = 1.9490 (-0.0100) Sample #2 = 1.9350 (-0.0120) Sample #3 = 1.9430 (0.0060) Sample #4 = 1.9580 (0.0260) Aug % Abs = 1.9453 (0.0067) STD DEU = 0.0117 (0.0190) REL STD DEU = 0.600 (285.132)	Sol Ualue = 0.300 g/210L *** Fit ualue = 1.4286 mg/1 %%% Samples Taken = 4, Discarded = 1 3um Io = 12506, 9um Io = 13866	Std Deu = 0.01 Rel Std Deu = 0.60 Sol Ual = 0.9524 mg/l or 0.200 g/210L % Abs = 3.711 Std Deu = 0.01 Rel Std Deu = 0.17 Sol Ual = 1.4286 mg/l or 0.300 g/210L % Abs = 5.403 Std Deu = 0.02 Rel Std Deu = 0.34 Zero Order Coef = -300.35 First Order Coef = 25.98 Standard Deviation = 27.839069	Sol UBIUE = 0.080 G/2:01 *** Fit ualue = 0.3810 mg/1 %%% Samples Taken = 4, Discarded = 1 ***** CHANNEL 1 Sample #1 = 2890.00 Sample #2 = 2929.00 Sample #4 = 2880.00 Auerage Result = 2896.6667 STD DEU = 28.0060 *REL STD DEU = 0.967 ************************************
Aug % Abs = 0.1610 (0.0060) STD DEU = 0.0056 (0.0066) REL STD DEU = 3.458 (109.291) Sol Ualue = 0.040 g/210L *** Fit ualue = 0.1905 mg/l %%% Samples Taken = 4, Discarded = 1	Sol Value = 0.200 g/210L *** Fit value = 0.9524 mg/l %%% Samples Taken = 4, Discarded = 1 3um Io = 12515, 9um Io = 13870		Sol Ual = 1.4286 mg/l or 0.300 g/210L % Abs = 9.900 Std Dev = 0.01 Rel Std Dev = 0.10 Zero Order Coef = -210.88 First Order Coef = 1337.50 Second Order Coef = 12.77 Standard Deviation = 5.179947	
3um Io = 12543, 9um Io = 13886 <<<<<<>>CHANNEL I >>>>> Sample & Abs (% Abs Ref) Sample #1 = 0.8730 (-0.0100) Sample #2 = 0.8230 (0.0380) Sample #3 = 0.8420 (0.0500) Sample #4 = 0.8540 (0.0540) Sample #4 = 0.8397 (0.0473) STD DEU = 0.0156 (0.0083) REL STD DEU = 1.862 (17.592)	Sample #2 = 3.7090 (0.0000) Sample #3 = 3.7180 (0.0050) Sample #4 = 3.7060 (0.0330) Aug % Abs = 3.7110 (0.0127) STD DEU = 0.0062 (0.0178) REL STD DEU = 0.168 (140.414)	Optical Calibration SN: 80-00 [453 Agency: SJ. (decsburg P() Date: 10 17 2022 Quadratic Fit: +/- 0.002g/210L √ By: TDG M	Solution Stats Quadratic Fit Chan 1 Act Fit Residual g/210L g/210L g/210L 0.000 0.001 -0.0006 0.040 0.039 0.0009 0.100 0.100 -0.0001 0.200 0.200 -0.0004 0.300 0.300 0.0002	

Type of Test	Serial Number	Agency	Date	Performed By
Stabilities (Post-GI)	80-001653	St. Petersburg PD	10 17 2022	TDG MG

0.05g/210L	0.08g/210L	0.20g/210L	DGS 0.08g/210L
0.047 to 0.053 🗸	0.077 to 0.083 🗸	0.194 to 0.206 🖌	0.077 to 0.083 🗸 ≤0.003 of Wet 🗸
ST PETERSBURG PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001653 10/17/2022 Software: 8100.27	ST PETERSBURG PD Intoxilyzer - Alconol Analyzer Model 8000 SN 80-001653 10/17/2022 Software: 8100.27	ST PETERSBURG PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001653 10/17/2022 Software: 8100.27	ST PETERSBURG PD Intoxilyzer - Alconol Analyzer Model 8000 SN 80-001653 10/17/2022 Software: 8100.27
Test g/210L Time	Test g/210L Time	Test g/210L Time	Test g/210L Time
Air Blank 0.000 10:34 Control Test 0.048 10:35 Air Blank 0.000 10:35 Control Test 0.048 10:36 Air Blank 0.000 10:37 Control Test 0.048 10:37 Control Test 0.048 10:37 Control Test 0.048 10:38 Control Test 0.048 10:38 Control Test Stats Auerage 0.0480 Std Deu 0.0000 Rel Std Deu(%) 0.0000	Air Blank 0.000 10;42 Control Test 0.078 10:43 Air Blank 0.000 10:43 Control Test 0.078 10:44 Air Blank 0.000 10:45 Control Test 0.078 10:45 Control Test 0.078 10:45 Control Test 0.078 10:45 Air Blank 0.000 10:46 Control-Test Stds 10:46 Control-Test 0.0780 10:46 Control-Test Std Deu 0.0000 Rel Std Deu(%) 0.0000 10:46	Air Blank 0.000 10:51 Control Test 0.199 10:51 Air Blank 0.000 10:52 Control Test 0.198 10:53 Air Blank 0.000 10:53 Control Test 0.198 10:54 Air Blank 0.000 10:54 Control Test 0.198 10:54 Control Test Stats Auerage 0.1983 Std Deu 0.0006 Rel Std Deu(%) 0.2911	Air Blank 0.000 10:56 Control Test 0.079 10:56 Air Blank 0.000 10:57 Air Blank 0.000 10:57 Control Test 0.079 10:57 Control Test 0.079 10:58 Air Blank 0.000 10:58 Control Test 0.079 10:58 Air Blank 0.000 10:58 Control Test Stats Auerage 0.0790 Std Deu 0.0000 Rel Std Deu(%) 0.0000
Operator's Signature	Operator's Signature	Operator's Signature	Operator's Signature

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: ST PETERSBURG PD Time of Inspection: 14:02

Date of Inspection: 10/17/2022

Serial Number: 80-001653 Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202201C Exp: 01/11/2024	0.08g/210L Test (g/210L) Lot#:202201D Exp: 01/18/2024	0.20g/210L Test (g/210L) Lot#:202201E Exp: 01/18/2024	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:00521080A2 Exp: 02/05/2023		
0.000	0.048	0.078	0.198	0.079		
0.000	0.048	0.078	0.198	0.079		
0.000	0.048	0.078	0.198	0.079		
0.000	0.048	0.078	0.198	0.079		
0.000	0.049	0.078	0.198	0.079		
0.000	0.048	0.078	0.198	0.079		
0.000	0.048	0.078	0.198	0.079		
0.000	0.049	0.078	0.198	0.079		
0.000	0.049	0.079	0.198	0.079		
0.000	0.049	0.078	0.198	0.079		
3			\$			
Standard Deviations	0.0005	0.0003	0.0000	0.0000		

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0002 Number of Simulators Used: 5

Remarks:

The above instrument complies (X) does not comply (

) with Chapter 11D-8, FAC.

I	certify	that	I performed	this	inspection	in	accordance	with	the	provisions	of	Chapter	11D-8,	FAC.
	1	\cap	$(I \cap$											
	. 0	14	V hours .	2						TAYLOR I	D GI	JTSCHOW		

Signature and Printed Name

10/17/2022 Date

X http:

Jay G



Calibration Certificate

Florida Department of Law Enforcement Alcohol Testing Program 4700 Terminal Drive, Suite 1 Ft. Myers, FL 33907

This is to certify the calibration of Intoxilyzer 8000 serial number <u>80-001653</u>, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	80-001653	UNCERTAINTY* ±	
Owning Agency:	ST PETERSBURG PD	0.050 g/ 210 L	0.004
Calibration Date:	10/17/2022	0.080 g/ 210 L	0.004
Calibration Time:	14:02	0.200 g/ 210 L	0.007
	1/15	0.080 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration. *Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence (k=3). The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

10/17/2022

Date

TAYLOR D GUTSCHOW. **Department Inspector**

FDLE/ATP Form 69 December 2021 Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality

Page 1 of 1