

0.05g/210L 0.047 to 0.053	0.08g/210L 0.077 to 0.083	0.20g/210L 0.194 to 0.206	DGS 0.08g/210L 0.077 to 0.083																																																																																																																																																
<p> FORT MYERS PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-005468 09/19/2023 Software: 8100.27 </p> <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>14:44</td></tr> <tr><td>Control Test</td><td>0.049</td><td>14:44</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:45</td></tr> <tr><td>Control Test</td><td>0.049</td><td>14:46</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:46</td></tr> <tr><td>Control Test</td><td>0.049</td><td>14:47</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:47</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0490</td><td></td></tr> <tr><td>Std Dev</td><td>0.0000</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr> </tbody> </table> <p>Operator's Signature _____</p>	Test	g/210L	Time	Air Blank	0.000	14:44	Control Test	0.049	14:44	Air Blank	0.000	14:45	Control Test	0.049	14:46	Air Blank	0.000	14:46	Control Test	0.049	14:47	Air Blank	0.000	14:47	Control Test Stats			Average	0.0490		Std Dev	0.0000		Rel Std Dev(%)	0.0000		<p> FORT MYERS PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-005468 09/19/2023 Software: 8100.27 </p> <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>14:50</td></tr> <tr><td>Control Test</td><td>0.076</td><td>14:50</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:51</td></tr> <tr><td>Control Test</td><td>0.077</td><td>14:52</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:52</td></tr> <tr><td>Control Test</td><td>0.078</td><td>14:53</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:54</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0770</td><td></td></tr> <tr><td>Std Dev</td><td>0.0010</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>1.2987</td><td></td></tr> </tbody> </table> <p>Operator's Signature _____</p>	Test	g/210L	Time	Air Blank	0.000	14:50	Control Test	0.076	14:50	Air Blank	0.000	14:51	Control Test	0.077	14:52	Air Blank	0.000	14:52	Control Test	0.078	14:53	Air Blank	0.000	14:54	Control Test Stats			Average	0.0770		Std Dev	0.0010		Rel Std Dev(%)	1.2987		<p> FORT MYERS PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-005468 09/19/2023 Software: 8100.27 </p> <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>14:56</td></tr> <tr><td>Control Test</td><td>0.199</td><td>14:57</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:58</td></tr> <tr><td>Control Test</td><td>0.198</td><td>14:58</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:59</td></tr> <tr><td>Control Test</td><td>0.198</td><td>15:00</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:00</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.1983</td><td></td></tr> <tr><td>Std Dev</td><td>0.0016</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.2911</td><td></td></tr> </tbody> </table> <p>Operator's Signature _____</p>	Test	g/210L	Time	Air Blank	0.000	14:56	Control Test	0.199	14:57	Air Blank	0.000	14:58	Control Test	0.198	14:58	Air Blank	0.000	14:59	Control Test	0.198	15:00	Air Blank	0.000	15:00	Control Test Stats			Average	0.1983		Std Dev	0.0016		Rel Std Dev(%)	0.2911		<p> FORT MYERS PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-005468 09/19/2023 Software: 8100.27 </p> <p>OKS</p> <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>15:01</td></tr> <tr><td>Control Test</td><td>0.080</td><td>15:01</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:02</td></tr> <tr><td>Control Test</td><td>0.079</td><td>15:02</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:03</td></tr> <tr><td>Control Test</td><td>0.079</td><td>15:03</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:04</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0793</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.7277</td><td></td></tr> </tbody> </table> <p>Operator's Signature _____</p>	Test	g/210L	Time	Air Blank	0.000	15:01	Control Test	0.080	15:01	Air Blank	0.000	15:02	Control Test	0.079	15:02	Air Blank	0.000	15:03	Control Test	0.079	15:03	Air Blank	0.000	15:04	Control Test Stats			Average	0.0793		Std Dev	0.0006		Rel Std Dev(%)	0.7277	
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	14:44																																																																																																																																																	
Control Test	0.049	14:44																																																																																																																																																	
Air Blank	0.000	14:45																																																																																																																																																	
Control Test	0.049	14:46																																																																																																																																																	
Air Blank	0.000	14:46																																																																																																																																																	
Control Test	0.049	14:47																																																																																																																																																	
Air Blank	0.000	14:47																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0490																																																																																																																																																		
Std Dev	0.0000																																																																																																																																																		
Rel Std Dev(%)	0.0000																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	14:50																																																																																																																																																	
Control Test	0.076	14:50																																																																																																																																																	
Air Blank	0.000	14:51																																																																																																																																																	
Control Test	0.077	14:52																																																																																																																																																	
Air Blank	0.000	14:52																																																																																																																																																	
Control Test	0.078	14:53																																																																																																																																																	
Air Blank	0.000	14:54																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0770																																																																																																																																																		
Std Dev	0.0010																																																																																																																																																		
Rel Std Dev(%)	1.2987																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	14:56																																																																																																																																																	
Control Test	0.199	14:57																																																																																																																																																	
Air Blank	0.000	14:58																																																																																																																																																	
Control Test	0.198	14:58																																																																																																																																																	
Air Blank	0.000	14:59																																																																																																																																																	
Control Test	0.198	15:00																																																																																																																																																	
Air Blank	0.000	15:00																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.1983																																																																																																																																																		
Std Dev	0.0016																																																																																																																																																		
Rel Std Dev(%)	0.2911																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	15:01																																																																																																																																																	
Control Test	0.080	15:01																																																																																																																																																	
Air Blank	0.000	15:02																																																																																																																																																	
Control Test	0.079	15:02																																																																																																																																																	
Air Blank	0.000	15:03																																																																																																																																																	
Control Test	0.079	15:03																																																																																																																																																	
Air Blank	0.000	15:04																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0793																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	0.7277																																																																																																																																																		

FORT MYERS PD
Intoxilyzer - Alcohol Analyzer
Model 8000
09/21/2023
SN 80-005468
10:14:43

Auto Calibration
Max Power Res Value = 23
Auto Range Res Value = 10

<<<< CHANNEL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 1.5950 (-0.0070)
Sample #2 = 1.5710 (0.0020)
Sample #3 = 1.6120 (-0.0090)
Sample #4 = 1.5730 (0.0240)
Avg % Abs = 1.5853 (0.0057)
STD DEV = 0.0231 (0.0168)
REL STD DEV = 1.458 (296.520)

Sol Value = 0.100 g/210L ***
Fit Value = 0.4762 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12523, Sum Io = 13717

<<<< CHANNEL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 1.7140 (0.0020)
Sample #2 = 1.7490 (0.0230)
Sample #3 = 1.7760 (0.0040)
Sample #4 = 1.7710 (0.0380)
Avg % Abs = 1.7653 (0.0217)
STD DEV = 0.0144 (0.0170)
REL STD DEV = 0.814 (78.642)

Sol Value = 0.000 g/210L ***
Fit Value = 0.0000 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12544, Sum Io = 13729

<<<< CHANNEL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.0490 (-0.0100)
Sample #2 = 0.0900 (0.0150)
Sample #3 = 0.0900 (0.0340)
Sample #4 = 0.0750 (0.0520)
Avg % Abs = 0.0850 (0.0337)
STD DEV = 0.0087 (0.0185)
REL STD DEV = 10.189 (54.957)

<<<< CHANNEL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.1430 (-0.0130)
Sample #2 = 0.1790 (-0.0080)
Sample #3 = 0.1790 (0.0140)
Sample #4 = 0.1750 (0.0210)
Avg % Abs = 0.1777 (0.0090)
STD DEV = 0.0023 (0.0151)
REL STD DEV = 1.300 (168.142)

Sol Value = 0.040 g/210L ***
Fit Value = 0.1905 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12532, Sum Io = 13722

<<<< CHANNEL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.7740 (-0.0080)
Sample #2 = 0.7470 (0.0110)
Sample #3 = 0.7800 (0.0140)
Sample #4 = 0.7550 (0.0420)
Avg % Abs = 0.7607 (0.0223)
STD DEV = 0.0172 (0.0171)
REL STD DEV = 2.263 (76.557)

<<<< CHANNEL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 3.6000 (-0.0060)
Sample #2 = 3.6260 (0.0140)
Sample #3 = 3.6290 (0.0160)
Sample #4 = 3.6120 (0.0330)
Avg % Abs = 3.6223 (0.0210)
STD DEV = 0.0091 (0.0104)
REL STD DEV = 0.250 (49.716)

Sol Value = 0.200 g/210L ***
Fit Value = 0.9524 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12514, Sum Io = 13711

<<<< CHANNEL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 3.4220 (-0.0100)
Sample #2 = 3.3760 (0.0280)
Sample #3 = 3.3900 (0.0400)
Sample #4 = 3.3810 (0.0470)
Avg % Abs = 3.3823 (0.0383)
STD DEV = 0.0071 (0.0096)
REL STD DEV = 0.210 (25.067)

***** AUTO CAL DATA *****

<<<< CHANNEL 1 >>>>
Sol Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.085
Std Dev = 0.01 Rel Std Dev = 10.19
Sol Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 0.761
Std Dev = 0.02 Rel Std Dev = 2.26
Sol Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 1.765
Std Dev = 0.01 Rel Std Dev = 0.81
Sol Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 3.382
Std Dev = 0.01 Rel Std Dev = 0.21
Sol Val = 1.4286 mg/l or 0.300 g/210L
% Abs = 4.995
Std Dev = 0.03 Rel Std Dev = 0.56
Zero Order Coef = -248.73
First Order Coef = 2819.36
Second Order Coef = 18.45
Standard Deviation = 18.571564

<<<< CHANNEL 2 >>>>
Sol Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.178
Std Dev = 0.00 Rel Std Dev = 1.30
Sol Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 1.585
Std Dev = 0.02 Rel Std Dev = 1.46
Sol Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 3.622
Std Dev = 0.01 Rel Std Dev = 0.25
Sol Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 6.881
Std Dev = 0.01 Rel Std Dev = 0.10
Sol Val = 1.4286 mg/l or 0.300 g/210L
% Abs = 9.951
Std Dev = 0.02 Rel Std Dev = 0.17
Zero Order Coef = -235.41
First Order Coef = 1330.61
Second Order Coef = 12.90
Standard Deviation = 5.816043

Solution Stats Quadratic Fit Chan 2
Act Fit Residual
g/210L g/210L g/210L
0.000 0.000 -0.0000
0.040 0.040 -0.0000
0.100 0.100 0.0002
0.200 0.200 -0.0002
0.300 0.300 0.0001
Sol Value = 0.080 g/210L ***
Fit Value = 0.3810 mg/l %%%
Samples Taken = 4, Discarded = 1
***** CHANNEL 1 *****
Sample #1 = 3483.00
Sample #2 = 3435.00
Sample #3 = 3427.00
Sample #4 = 3442.00
Average Result = 3434.6667
STD DEV = 7.5056
REL STD DEV = 0.219

***** CHANNEL 2 *****
Sample #1 = 3338.00
Sample #2 = 3284.00
Sample #3 = 3318.00
Sample #4 = 3298.00
Average Result = 3300.0000
STD DEV = 17.0880
REL STD DEV = 0.518

Dry Gas H2O Adjust Results *****
Barometric Pressure = 1016
3 um H2O Adjust (mg/l x 10,000) = 375
9 um H2O Adjust (mg/l x 10,000) = 509
***** AUTO CAL PASS *****

Optical Calibration

SN: 80-005468
Agency: Fort Myers PD
Date: 09/21/2023
Quadratic Fit: +/- 0.002g/210L ✓
By: TDG MB

Solution Stats Quadratic-Fit Chan 1
Act Fit Residual
g/210L g/210L g/210L
0.000 -0.000 0.0002
0.040 0.040 -0.0000
0.100 0.101 -0.0005
0.200 0.199 0.0005
0.300 0.300 -0.0002

Post-Cal Stability Checks

0.05g/210L 0.047 to 0.053	0.08g/210L 0.077 to 0.083	0.20g/210L 0.194 to 0.206	DGS 0.08g/210L 0.077 to 0.083	≤0.003 of Wet
<p>FORT MYERS PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-005468 09/21/2023 Software: 8100.27</p> <p>Test g/210L Time</p> <p>Air Blank 0.000 11:23 Control Test 0.048 11:24 Air Blank 0.000 11:24 Control Test 0.048 11:25 Air Blank 0.000 11:25 Control Test 0.048 11:26 Air Blank 0.000 11:27</p> <p>Control Test Stats Average 0.0480 Std Dev 0.0000 Rel Std Dev(%) 0.0000</p> <p>Operator's Signature <i>MC</i></p>	<p>FORT MYERS PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-005468 09/21/2023 Software: 8100.27</p> <p>Test g/210L Time</p> <p>Air Blank 0.000 11:28 Control Test 0.077 11:29 Air Blank 0.000 11:29 Control Test 0.078 11:30 Air Blank 0.000 11:30 Control Test 0.077 11:31 Air Blank 0.000 11:32</p> <p>Control Test Stats Average 0.0773 Std Dev 0.0006 Rel Std Dev(%) 0.7466</p> <p>Operator's Signature <i>MC</i></p>	<p>FORT MYERS PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-005468 09/21/2023 Software: 8100.27</p> <p>Test g/210L Time</p> <p>Air Blank 0.000 11:38 Control Test 0.199 11:39 Air Blank 0.000 11:39 Control Test 0.200 11:40 Air Blank 0.000 11:40 Control Test 0.199 11:41 Air Blank 0.000 11:41</p> <p>Control Test Stats Average 0.1993 Std Dev 0.0006 Rel Std Dev(%) 0.2896</p> <p>Operator's Signature <i>MC</i></p>	<p><i>DGS</i></p> <p>FORT MYERS PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-005468 09/21/2023 Software: 8100.27</p> <p>Test g/210L Time</p> <p>Air Blank 0.000 11:18 Control Test 0.078 11:18 Air Blank 0.000 11:19 Control Test 0.078 11:19 Air Blank 0.000 11:20 Control Test 0.078 11:20 Air Blank 0.000 11:20</p> <p>Control Test Stats Average 0.0780 Std Dev 0.0000 Rel Std Dev(%) 0.0000</p> <p>Operator's Signature <i>MC</i></p>	<input checked="" type="checkbox"/>

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: FORT MYERS PD
Time of Inspection: 14:41

Date of Inspection: 09/21/2023

Serial Number: 80-005468
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202201C Exp: 01/11/2024	0.08g/210L Test (g/210L) Lot#:202201D Exp: 01/18/2024	0.20g/210L Test (g/210L) Lot#:202201E Exp: 01/18/2024	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG223802 Exp: 08/26/2024
0.000	0.048	0.078	0.200	0.078
0.000	0.048	0.077	0.199	0.078
0.000	0.048	0.077	0.199	0.078
0.000	0.048	0.077	0.199	0.078
0.000	0.048	0.077	0.200	0.078
0.000	0.048	0.078	0.200	0.078
0.000	0.049	0.078	0.199	0.078
0.000	0.048	0.077	0.199	0.078
0.000	0.048	0.078	0.199	0.077
0.000	0.049	0.077	0.199	0.078


Standard Deviations	0.0004	0.0005	0.0004	0.0003
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0004 Number of Simulators Used: 5

Remarks:

The above instrument complies (X) does not comply () with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



TAYLOR D GUTSCHOW

Signature and Printed Name

09/21/2023
Date



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
4700 Terminal Drive, Suite 1
Ft. Myers, FL 33907

This is to certify the calibration of Intoxilyzer 8000 serial number 80-005468, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	80-005468	UNCERTAINTY* \pm
Owning Agency:	FORT MYERS PD	0.050 g/ 210 L 0.004
Calibration Date:	09/21/2023	0.080 g/ 210 L 0.004
Calibration Time:	14:41	0.200 g/ 210 L 0.007
		0.080 g/ 210 L Dry Gas Control 0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards. This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

09/21/2023

Date

Taylor D Gutenschow
Department Inspector

FDLE/ATP Form 69 December 2021

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality