



Agency Clearwater PD

S/N 80-001740

Date In 3/1/17

Date Out 3/8/17

☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Optical Bench Calibration		Performed By <u>QWS</u>		Department Inspection		Performed By <u>QWS</u>	
☐ Optical Bench Calibration N/A ☒ Optical Bench Calibration Complete				☒ Barometric Pressure ID# <u>28421</u>		<u>1021</u> Gauge <u>1022</u> Instrument	
Barometric Pressure Gauge <u>1028</u> ID # <u>28427</u>				Mouth Alcohol Solution Lot # <u>2D16-C</u> Acetone Stock Solution Lot # <u>2D17-A</u>			
Simulator	Serial Number	Lot Number	Expiration	Simulator		Serial Number	
0.000	<u>SD1016</u>	<u>N/A</u>	<u>N/A</u>	0.00		<u>G2879</u>	
0.040	<u>SD1024</u>	<u>16101</u>	<u>2/2/18</u>	Interferent		<u>G8144</u>	
0.100	<u>SD1022</u>	<u>16001</u>	<u>5/8/18</u>	0.05		<u>G11739</u>	
0.200	<u>DR3855</u>	<u>16103</u>	<u>6/14/18</u>	0.08		<u>G8149</u>	
0.400	<u>G2407</u>	<u>16102</u>	<u>3/22/18</u>	0.20		<u>G11621</u>	
0.080 DGS	<u>N/A</u>	<u>D3415080A1</u>	<u>3/5/17</u>				
☒ Post Calibration Stability Checks				Attachments ☒ Form 41 ☒ Pre-Stability Tests ☐ Flow Calibration ☒ Optical Bench Cal ☒ Post-Stability Tests ☐ Other _____			
Simulator	Serial Number	Lot Number	Expiration				
0.05	<u>DR2035</u>	<u>201603D</u>	<u>3/8/18</u>				
0.08	<u>SD1011</u>	<u>201601F</u>	<u>1/26/18</u>				
0.20	<u>SD1025</u>	<u>201604C</u>	<u>4/5/18</u>				
0.08 DGS	<u>N/A</u>	<u>AG626605</u>	<u>9/22/18</u>				
Notes: <u>Performed optical bench calibration adjustment to bring values closer to nominal. QWS</u> <u>PAAC OK SP</u> <u>P # 11/11/18</u>				☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use			

Quality Control Review

Date _____

Revised November 2012

Pre-Cal

Stability Checks # 80-001740 Clearwater P.D. 3/3/17 *ARS*

DES
CLEARWATER PD
Intoxilyzer - Alcohol Analyzer
Model 8000
03/03/2017
Software: 8100.27
SN 80-001740

Test	g/210L	Time
Air Blank	0.000	12:57
Control Test	0.079	12:57
Air Blank	0.000	12:58
Control Test	0.078	12:58
Air Blank	0.000	12:59
Control Test	0.079	12:59
Air Blank	0.000	13:00
Control Test Stats		
Average	0.0787	
Std Dev	0.0006	
Rel Std Dev(%)	0.7339	

131K

SP

ARS
Operator's Signature

CLEARWATER PD
Intoxilyzer - Alcohol Analyzer
Model 8000
03/03/2017
Software: 8100.27
SN 80-001740

Test	g/210L	Time
Air Blank	0.000	12:52
Control Test	0.199	12:52
Air Blank	0.000	12:53
Control Test	0.200	12:53
Air Blank	0.000	12:54
Control Test	0.199	12:55
Air Blank	0.000	12:55
Control Test Stats		
Average	0.1993	
Std Dev	0.0006	
Rel Std Dev(%)	0.2896	

ARS
Operator's Signature

CLEARWATER PD
Intoxilyzer - Alcohol Analyzer
Model 8000
03/03/2017
Software: 8100.27
SN 80-001740

Test	g/210L	Time
Air Blank	0.000	12:47
Control Test	0.082	12:47
Air Blank	0.000	12:48
Control Test	0.082	12:49
Air Blank	0.000	12:49
Control Test	0.082	12:50
Air Blank	0.000	12:50
Control Test Stats		
Average	0.0820	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

ARS
Operator's Signature

CLEARWATER PD
Intoxilyzer - Alcohol Analyzer
Model 8000
03/03/2017
Software: 8100.27
SN 80-001740

Test	g/210L	Time
Air Blank	0.000	12:41
Control Test	0.051	12:42
Air Blank	0.000	12:42
Control Test	0.050	12:43
Air Blank	0.000	12:43
Control Test	0.051	12:44
Air Blank	0.000	12:45
Control Test Stats		
Average	0.0507	
Std Dev	0.0006	
Rel Std Dev(%)	1.1395	

ARS
Operator's Signature

Optical bench calibration data #80-001740 Clearwater P.D. 3/3/17

CLEARWATER PD

Intoxilizer - Alcohol Analyzer

SN 80-001740

14:20:20

03/03/2017

Auto Calibration

Max Power Res Value = 45

Auto Range Res Value = 23

Sol Value = 0.000 g/210L ***

Fit value = 0.0000 mg/l %%%

Samples Taken = 4, Discarded = 1

Sum Io = 12607, Sum Io = 13853

Sample % Abs (% Abs Ref)

Sample #1 = 0.1560 (-0.0170)

Sample #2 = 0.1340 (0.0180)

Sample #3 = 0.1390 (0.0130)

Sample #4 = 0.1290 (0.0550)

Avg % Abs = 0.1340 (0.0287)

STD DEV = 0.0050 (0.0229)

REL STD DEV = 3.731 (80.030)

Sample % Abs (% Abs Ref)

Sample #1 = 0.1500 (-0.0140)

Sample #2 = 0.1630 (-0.0050)

Sample #3 = 0.1630 (-0.0120)

Sample #4 = 0.1510 (0.0020)

Avg % Abs = 0.1590 (-0.0050)

STD DEV = 0.0069 (0.0070)

REL STD DEV = 4.357 (140.000)

Sol Value = 0.040 g/210L ***

Fit value = 0.1905 mg/l %%%

Samples Taken = 4, Discarded = 1

Sum Io = 12597, Sum Io = 13850

Sample % Abs (% Abs Ref)

Sample #1 = 0.8190 (-0.0180)

Sample #2 = 0.8280 (0.0070)

Sample #3 = 0.8250 (0.0160)

Sample #4 = 0.8350 (0.0030)

Avg % Abs = 0.8293 (0.0087)

STD DEV = 0.0051 (0.0067)

REL STD DEV = 0.619 (76.827)

Channel 2 >>>>

Sample % Abs (% Abs Ref)

Sample #1 = 1.5340 (-0.0140)

Sample #2 = 1.5450 (0.0000)

Sample #3 = 1.5420 (0.0020)

Sample #4 = 1.5260 (0.0040)

Avg % Abs = 1.5377 (0.0020)

STD DEV = 0.0102 (0.0020)

REL STD DEV = 0.664 (100.000)

Sol Value = 0.100 g/210L ***

Fit value = 0.4762 mg/l %%%

Samples Taken = 4, Discarded = 1

Sum Io = 12594, Sum Io = 13849

Sample % Abs (% Abs Ref)

Sample #1 = 1.8950 (-0.0220)

Sample #2 = 1.8730 (0.0000)

Sample #3 = 1.9060 (-0.0010)

Sample #4 = 1.8650 (0.0080)

Avg % Abs = 1.8813 (0.0023)

STD DEV = 0.0217 (0.0049)

REL STD DEV = 1.155 (211.409)

Channel 2 >>>>

Sample % Abs (% Abs Ref)

Sample #1 = 3.5830 (-0.0140)

Sample #2 = 3.5980 (-0.0070)

Sample #3 = 3.5980 (0.0000)

Sample #4 = 3.5930 (0.0020)

Avg % Abs = 3.5897 (-0.0017)

STD DEV = 0.0029 (0.0047)

REL STD DEV = 0.080 (283.549)

Sol Value = 0.200 g/210L ***

Fit value = 0.9524 mg/l %%%

Samples Taken = 4, Discarded = 1

Sum Io = 12590, Sum Io = 13848

Sample % Abs (% Abs Ref)

Sample #1 = 3.5660 (-0.0240)

Sample #2 = 3.5360 (-0.0080)

Sample #3 = 3.5520 (0.0060)

Sample #4 = 3.5480 (-0.0060)

Avg % Abs = 3.5453 (-0.0027)

STD DEV = 0.0083 (0.0076)

REL STD DEV = 0.235 (283.945)

Channel 2 >>>>

Sample % Abs (% Abs Ref)

Sample #1 = 6.8810 (-0.0100)

Sample #2 = 6.8690 (0.0020)

Sample #3 = 6.8600 (0.0140)

Sample #4 = 6.8610 (0.0060)

Avg % Abs = 6.8633 (0.0073)

STD DEV = 0.0049 (0.0061)

REL STD DEV = 0.072 (83.320)

Sol Value = 0.400 g/210L ***

Fit value = 1.9048 mg/l %%%

Samples Taken = 4, Discarded = 1

Sum Io = 12588, Sum Io = 13849

Sample % Abs (% Abs Ref)

Sample #1 = 6.7070 (-0.0190)

Sample #2 = 6.7550 (-0.0060)

Sample #3 = 6.7710 (0.0100)

Sample #4 = 6.7610 (0.0170)

Avg % Abs = 6.7623 (0.0070)

STD DEV = 0.0081 (0.0118)

REL STD DEV = 0.120 (168.426)

Channel 2 >>>>

Sample % Abs (% Abs Ref)

Sample #1 = 12.8130 (0.0030)

Sample #2 = 12.8610 (0.0280)

Sample #3 = 12.8710 (0.0290)

Sample #4 = 12.8710 (0.0390)

Avg % Abs = 12.8677 (0.0320)

STD DEV = 0.0058 (0.0061)

REL STD DEV = 0.045 (19.009)

Channel 1 >>>>

Sol Val = 0.0000 mg/l or 0.000 g/210L

% Abs = 0.134

Std Dev = 0.00 Rel Std Dev = 3.73

Sol Val = 0.1905 mg/l or 0.040 g/210L

% Abs = 0.829

Std Dev = 0.01 Rel Std Dev = 0.62

Sol Val = 0.4762 mg/l or 0.100 g/210L

% Abs = 1.881

Std Dev = 0.02 Rel Std Dev = 1.16

Sol Val = 0.9524 mg/l or 0.200 g/210L

% Abs = 3.545

Std Dev = 0.01 Rel Std Dev = 0.23

Sol Val = 1.9048 mg/l or 0.400 g/210L

% Abs = 6.762

Std Dev = 0.01 Rel Std Dev = 0.12

Zero Order Coef = -359.02

First Order Coef = 2686.64

Second Order Coef = 27.14

Standard Deviation = 19.043007

Channel 2 >>>>

Sol Val = 0.0000 mg/l or 0.000 g/210L

% Abs = 0.159

Std Dev = 0.01 Rel Std Dev = 4.36

Sol Val = 0.1905 mg/l or 0.040 g/210L

% Abs = 1.538

Std Dev = 0.01 Rel Std Dev = 0.66

Sol Val = 0.4762 mg/l or 0.100 g/210L

% Abs = 3.590

Std Dev = 0.00 Rel Std Dev = 0.08

Sol Val = 0.9524 mg/l or 0.200 g/210L

% Abs = 6.863

Std Dev = 0.00 Rel Std Dev = 0.07

Sol Val = 1.9048 mg/l or 0.400 g/210L

% Abs = 12.868

Std Dev = 0.01 Rel Std Dev = 0.04

Zero Order Coef = -192.58

First Order Coef = 1330.83

Second Order Coef = 12.75

Standard Deviation = 18.191330

Solution Stats Quadratic Fit Chan 1

Act Fit Residual

g/210L g/210L g/210L

0.000 0.000 -0.0000

0.040 0.040 -0.0004

0.100 0.101 -0.0006

0.200 0.200 0.0003

0.400 0.400 -0.0001

Solution Stats Quadratic Fit Chan 2

Act Fit Residual

g/210L g/210L g/210L

0.000 0.000 -0.0004

0.040 0.040 0.0004

0.100 0.100 0.0003

0.200 0.200 -0.0004

0.400 0.400 0.0001

Sol Value = 0.080 g/210L ***

Fit value = 0.3810 mg/l %%%

Samples Taken = 4, Discarded = 1

Channel 1

Sample #1 = 3031.00

Sample #2 = 3086.00

Sample #3 = 3144.00

Sample #4 = 3036.00

Average Result = 3088.6667

STD DEV = 54.0494

REL STD DEV = 1.750

Channel 2

Sample #1 = 3376.00

Sample #2 = 3377.00

Sample #3 = 3383.00

Sample #4 = 3364.00

Average Result = 3374.6667

STD DEV = 9.7125

REL STD DEV = 0.288

Channel 1

Dry Gas H2O Adjust Results *****

Barometric Pressure = 1027

3 um H2O Adjust (mg/l*10,000) = 721

9 um H2O Adjust (mg/l*10,000) = 435

***** AUTO CAL PASS

Post-Cal Stability Checks #80-001740

Clearwater P.D.

3/3/17

DCS

CLEARWATER PD
Intoxilyzer - Alcohol Analyzer
Model 8000
03/03/2017
Software: 8100.27
SN 80-001740

Test	g/210L	Time
Air Blank	0.000	20:15
Control Test	0.050	20:16
Air Blank	0.000	20:17
Control Test	0.050	20:17
Air Blank	0.000	20:18
Control Test	0.049	20:19
Air Blank	0.000	20:19
Control Test Stats		
Average	0.0497	
Std Dev	0.0006	
Rel Std Dev(%)	1.1625	

CLEARWATER PD
Intoxilyzer - Alcohol Analyzer
Model 8000
03/03/2017
Software: 8100.27
SN 80-001740

Test	g/210L	Time
Air Blank	0.000	20:20
Control Test	0.080	20:21
Air Blank	0.000	20:21
Control Test	0.081	20:22
Air Blank	0.000	20:22
Control Test	0.081	20:23
Air Blank	0.000	20:24
Control Test Stats		
Average	0.0807	
Std Dev	0.0006	
Rel Std Dev(%)	0.7157	

CLEARWATER PD
Intoxilyzer - Alcohol Analyzer
Model 8000
03/03/2017
Software: 8100.27
SN 80-001740

Test	g/210L	Time
Air Blank	0.000	20:32
Control Test	0.197	20:33
Air Blank	0.000	20:33
Control Test	0.198	20:34
Air Blank	0.000	20:35
Control Test	0.198	20:35
Air Blank	0.000	20:36
Control Test Stats		
Average	0.1977	
Std Dev	0.0006	
Rel Std Dev(%)	0.2921	

CLEARWATER PD
Intoxilyzer - Alcohol Analyzer
Model 8000
03/03/2017
Software: 8100.27
SN 80-001740

Test	g/210L	Time
Air Blank	0.000	20:25
Control Test	0.082	20:25
Air Blank	0.000	20:26
Control Test	0.082	20:26
Air Blank	0.000	20:26
Control Test	0.082	20:27
Air Blank	0.000	20:27
Control Test Stats		
Average	0.0820	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

13K

Operator's Signature

Operator's Signature

Operator's Signature

Operator's Signature

SP



Alcohol Testing Program

INSTRUMENT PROCESSING SHEET

Agency

Clearwater PD

S/N

80-001740

Date In

12/28/16

Date Out

1/6/17

☐ Ship☐ P/U☐ H/D☐ CMI☒ EE

Intake Performed By <u>TD</u>		Quality Checks Performed By <u>EDB</u>		Flow Calibration Performed By <u>EDB</u>																
☐ Registration ☒ Annual ☐ Return from CMI ☐ Return from Enforcement Electronics ☐ Other _____ Visual Inspection: <u>OK</u> Case <u>OK</u> Handle <u>OK</u> Dry Gas Holder <u>OK</u> Feet <u>OK</u> Keyboard/Plug <u>OK</u> Back/Plugs <u>OK</u> Screws tight <u>OK</u> Breath Hose Other Equipment: ☐ Power cord ☐ Printer Cable ☐ Other _____ Notes: _____ _____ _____ _____		☒ Breath Tube Screen ☒ Replace O-Rings ☒ Instrument Set Up Verified ☐ R-Value _____ ☐ Flow Verification (L/s) Flow Column # _____ 32mm _____ (.139 - .169) 36mm _____ (.156 - .190) 53mm _____ (.228 - .278) 103mm _____ (.447 - .547) ☐ Barometric Pressure Check Gauge ID # _____ ☐ Stability Checks <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td>0.05</td> <td></td> <td></td> </tr> <tr> <td>0.08</td> <td></td> <td></td> </tr> <tr> <td>0.20</td> <td></td> <td></td> </tr> <tr> <td>0.08 DGS</td> <td>N/A</td> <td></td> </tr> </tbody> </table>		Simulator	Serial #	Lot #/Exp	0.05			0.08			0.20			0.08 DGS	N/A		☐ Flow Calibration N/A ☐ Flow Calibration Complete Flow Column # _____ ☐ 5L/min - 47mm ☐ 15L/min - 53mm ☐ 30L/min - 103mm ☐ R-Value _____ ☐ Post Calibration Verification (L/s) Flow Column # _____ 32mm _____ (.139 - .169) 36mm _____ (.156 - .190) 53mm _____ (.228 - .278) 103mm _____ (.447 - .547) Maintenance Performed By _____ ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other _____ Suggested Service _____ _____ _____	
Simulator	Serial #	Lot #/Exp																		
0.05																				
0.08																				
0.20																				
0.08 DGS	N/A																			

Optical Bench Calibration Performed By _____		Department Inspection Performed By _____																																																													
☐ Optical Bench Calibration N/A ☐ Optical Bench Calibration Complete Barometric Pressure Gauge _____ ID # _____ <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td></td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.100</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.400</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table> ☐ Post Calibration Stability Checks <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.05</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.08</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.20</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.08 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table>		Simulator	Serial Number	Lot Number	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.400				0.080 DGS	N/A			Simulator	Serial Number	Lot Number	Expiration	0.05				0.08				0.20				0.08 DGS	N/A			☐ Barometric Pressure _____ Gauge ID# _____ Instrument Mouth Alcohol Solution Lot # _____ Acetone Stock Solution Lot # _____ <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.00</td> <td></td> </tr> <tr> <td>Interferent</td> <td></td> </tr> <tr> <td>0.05</td> <td></td> </tr> <tr> <td>0.08</td> <td></td> </tr> <tr> <td>0.20</td> <td></td> </tr> </tbody> </table> Attachments ☐ Form 41 ☐ Pre-Stability Tests ☐ Flow Calibration ☐ Optical Bench Cal ☐ Post-Stability Tests ☐ Other _____		Simulator	Serial Number	0.00		Interferent		0.05		0.08		0.20	
Simulator	Serial Number	Lot Number	Expiration																																																												
0.000		N/A	N/A																																																												
0.040																																																															
0.100																																																															
0.200																																																															
0.400																																																															
0.080 DGS	N/A																																																														
Simulator	Serial Number	Lot Number	Expiration																																																												
0.05																																																															
0.08																																																															
0.20																																																															
0.08 DGS	N/A																																																														
Simulator	Serial Number																																																														
0.00																																																															
Interferent																																																															
0.05																																																															
0.08																																																															
0.20																																																															

Notes: Modern Test on Diagnostic checks
failng. Coupled w/ Results from last AI -
recommend repair.

- ☐ Instrument Complies with Chapter 11D-8, FAC
☒ Instrument Does Not Comply with Chapter 11D-8, FAC
☐ Return to/Place into Evidentiary Use
☒ Remain Out of Evidentiary Use
☐ Conduct an Agency Inspection Before Evidentiary Use

Quality Control Review

Date

Revised November 2012