



# INSTRUMENT PROCESSING SHEET

Agency Gadsden County Sheriff's OfficeS/N 80-000959

Florida Department of Law Enforcement

Date In 09/23/2020 DI Completion Date 09-24-2020 Ship  P/U  H/D  CMI  EE9-25-2020

<b>Intake</b> Performed By <u>RAW</u> <input checked="" type="checkbox"/> Annual <input checked="" type="checkbox"/> Registration <input checked="" type="checkbox"/> Return from CMI / EE Visual Inspection: <input checked="" type="checkbox"/> Case <input checked="" type="checkbox"/> Handle <input checked="" type="checkbox"/> Keyboard <input checked="" type="checkbox"/> Dry Gas Shelf <input checked="" type="checkbox"/> Feet <input checked="" type="checkbox"/> Breath Tube <input checked="" type="checkbox"/> Ports <input checked="" type="checkbox"/> Screws Tight Other Equipment/ Accessories: <input type="checkbox"/> Power cord <input type="checkbox"/> Printer Cable <input checked="" type="checkbox"/> Static Bag <input type="checkbox"/> 12V DC Cable Notes: _____ _____ _____		<b>Quality Checks</b> Performed By <u>IS</u> <input checked="" type="checkbox"/> Breath Tube Screen <input checked="" type="checkbox"/> Replace External O-Rings <input checked="" type="checkbox"/> Instrument Set Up Verified <input checked="" type="checkbox"/> R-Value <u>202</u> <input checked="" type="checkbox"/> Flow Verification (L/s) Flow Column # <u>ATP-105</u> 32 mm <u>0.160</u> (.139 - .169) 36 mm <u>0.175</u> (.156 - .190) 53 mm <u>0.246</u> (.228 - .278) 103 mm <u>0.507</u> (.447 - .547) <input checked="" type="checkbox"/> Barometric Pressure Check Gauge ID # <u>28427</u> <input checked="" type="checkbox"/> Stability Checks		<b>Flow Calibration</b> Performed By _____ Flow Column # _____ <input type="checkbox"/> 5L/min - 17mm <input type="checkbox"/> 15L/min - 53mm <input type="checkbox"/> 30L/min - 103mm <input type="checkbox"/> R-Value _____ <input type="checkbox"/> Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)																																																											
<b>Final Release Date</b> FDLE Alcohol Testing Program Digitally signed by FDLE Alcohol Testing Program Date: 2020.09.25 11:19:59 -04'00'		<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td>MP5088</td> <td>201905A 05-14-2021</td> </tr> <tr> <td>0.080</td> <td>MP5089</td> <td>201905B 05-18-2021</td> </tr> <tr> <td>0.200</td> <td>MP5090</td> <td>201904D 04-30-2021</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>AG931603 11-12-2021</td> </tr> </tbody> </table>		Simulator	Serial #	Lot #/Exp	0.050	MP5088	201905A 05-14-2021	0.080	MP5089	201905B 05-18-2021	0.200	MP5090	201904D 04-30-2021	0.080 DGS	N/A	AG931603 11-12-2021	<b>Maintenance</b> Performed By _____ <input type="checkbox"/> Battery Replacement <input type="checkbox"/> Dry Gas Regulator Replacement <input type="checkbox"/> Breath Tube Replacement <input type="checkbox"/> Other _____																																												
Simulator	Serial #	Lot #/Exp																																																													
0.050	MP5088	201905A 05-14-2021																																																													
0.080	MP5089	201905B 05-18-2021																																																													
0.200	MP5090	201904D 04-30-2021																																																													
0.080 DGS	N/A	AG931603 11-12-2021																																																													
<b>Calibration Adjustment</b> Performed By _____ Barometric Pressure Gauge ID # _____ <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td></td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.100</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.300</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table> <input type="checkbox"/> Post Calibration Adjustment Stability Checks <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table>		Simulator	Serial Number	Lot Number	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			Simulator	Serial Number	Lot Number	Expiration	0.050				0.080				0.200				0.080 DGS	N/A			<b>Department Inspection</b> Performed By <u>IS</u> Barometric Pressure ID# <u>30793</u> Gauge <u>1012</u> Instrument <u>1011</u> Mouth Alcohol Solution Lot # <u>2019-B</u> Acetone Stock Solution Lot # <u>2019-A</u> <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td>MP5086</td> </tr> <tr> <td>Interferent</td> <td>MP5087</td> </tr> <tr> <td>0.050</td> <td>MP5088</td> </tr> <tr> <td>0.080</td> <td>MP5089</td> </tr> <tr> <td>0.200</td> <td>MP5090</td> </tr> </tbody> </table>		Simulator	Serial Number	0.000	MP5086	Interferent	MP5087	0.050	MP5088	0.080	MP5089	0.200	MP5090
Simulator	Serial Number	Lot Number	Expiration																																																												
0.000		N/A	N/A																																																												
0.040																																																															
0.100																																																															
0.200																																																															
0.300																																																															
0.080 DGS	N/A																																																														
Simulator	Serial Number	Lot Number	Expiration																																																												
0.050																																																															
0.080																																																															
0.200																																																															
0.080 DGS	N/A																																																														
Simulator	Serial Number																																																														
0.000	MP5086																																																														
Interferent	MP5087																																																														
0.050	MP5088																																																														
0.080	MP5089																																																														
0.200	MP5090																																																														
Notes/Suggested Service: <u>Admin Review: Switched "Ship" to "P/U", corrected 0.08ARS to 805-14-2021</u> <u>9-25-2020</u>		<input checked="" type="checkbox"/> Instrument Complies with Chapter 11D-8, FAC <input type="checkbox"/> Instrument Does Not Comply with Chapter 11D-8, FAC <input checked="" type="checkbox"/> Return to/Place into Evidentiary Use <input type="checkbox"/> Remain Out of Evidentiary Use <input checked="" type="checkbox"/> Conduct an Agency Inspection Before Evidentiary Use																																																													
_____ _____ _____		<table border="1"> <tr> <td> <u>Michael D. Hargrave</u>            2020.09.25            10:05:13 -04'00'         </td> <td> <u>CR</u>            2020.09.25            11:17:01 -04'00'         </td> </tr> <tr> <td>Tech Review / Date</td> <td>Admin Review / Date</td> </tr> </table>		<u>Michael D. Hargrave</u> 2020.09.25 10:05:13 -04'00'	<u>CR</u> 2020.09.25 11:17:01 -04'00'	Tech Review / Date	Admin Review / Date																																																								
<u>Michael D. Hargrave</u> 2020.09.25 10:05:13 -04'00'	<u>CR</u> 2020.09.25 11:17:01 -04'00'																																																														
Tech Review / Date	Admin Review / Date																																																														

# Florida Department of Law Enforcement Alcohol Testing Program

## DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: GADSDEN COUNTY SO  
Time of Inspection: 14:20

Date of Inspection: 09/24/2020

Serial Number: 80-000959  
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:201905A Exp: 05/14/2021	0.08g/210L Test (g/210L) Lot#:201905B Exp: 05/14/2021	0.20g/210L Test (g/210L) Lot#:201904D Exp: 04/30/2021	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG011102 Exp: 04/20/2022
0.000	0.049	0.081	0.203	0.081
0.000	0.049	0.081	0.203	0.081
0.000	0.049	0.081	0.204	0.080
0.000	0.050	0.081	0.203	0.080
0.000	0.049	0.081	0.203	0.080
0.000	0.049	0.081	0.203	0.079
0.000	0.049	0.081	0.203	0.080
0.000	0.049	0.081	0.204	0.080
0.000	0.050	0.081	0.203	0.080
0.000	0.049	0.081	0.203	0.080

Standard Deviations	0.0004	0.0000	0.0004	0.0005
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0003 Number of Simulators Used: 5

Remarks:

The above instrument complies (  ) does not comply (  ) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

*Israel Soto*

ISRAEL SOTO

Signature and Printed Name

09/24/2020  
Date

*MX*

*SP*  
2020.0  
9.25  
11:17:24  
-04'00'

# Stability Checks

GADSDEN COUNTY SO  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-000959  
 09/24/2020  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:12
Control Test	0.050	12:13
Air Blank	0.000	12:13
Control Test	0.050	12:14
Air Blank	0.000	12:14
Control Test	0.049	12:15
Air Blank	0.000	12:16
Control Test Stats		
Average	0.0497	
Std Dev	0.0006	
Rel Std Dev(%)	1.1625	

  
 -----  
 Operator's Signature

GADSDEN COUNTY SO  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-000959  
 09/24/2020  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:20
Control Test	0.081	12:21
Air Blank	0.000	12:22
Control Test	0.080	12:22
Air Blank	0.000	12:23
Control Test	0.081	12:24
Air Blank	0.000	12:24
Control Test Stats		
Average	0.0807	
Std Dev	0.0006	
Rel Std Dev(%)	0.7157	

wet

  
 -----  
 Operator's Signature

GADSDEN COUNTY SO  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-000959  
 09/24/2020  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:26
Control Test	0.203	12:26
Air Blank	0.000	12:27
Control Test	0.203	12:28
Air Blank	0.000	12:28
Control Test	0.202	12:29
Air Blank	0.000	12:29
Control Test Stats		
Average	0.2027	
Std Dev	0.0006	
Rel Std Dev(%)	0.2849	


  
 -----  
 Operator's Signature

GADSDEN COUNTY SO  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-000959  
 09/24/2020  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:30
Control Test	0.082	12:30
Air Blank	0.000	12:31
Control Test	0.082	12:31
Air Blank	0.000	12:32
Control Test	0.083	12:32
Air Blank	0.000	12:33
Control Test Stats		
Average	0.0823	
Std Dev	0.0006	
Rel Std Dev(%)	0.7012	

Dry

  
 -----  
 Operator's Signature

MK  
  
 2020.09.  
 25  
 11:18:08  
 -04'00'



# Calibration Certificate

Florida Department of Law Enforcement  
Alcohol Testing Program  
2729 Fort Knox Blvd.  
Bldg. 2, Suite 1300  
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-000959, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-000959</u>	UNCERTAINTY* ±	
Owning Agency:	<u>GADSDEN COUNTY SO</u>	0.050 g/210 L	0.004
Calibration Date:	<u>09/24/2020</u>	0.080 g/210 L	0.005
Calibration Time:	<u>14:20</u>	0.200 g/210 L	0.007
		0.080 g/210 L Dry Gas Control	0.005

All results are reported in g/210 L.  
Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.  
\*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence (K=3).  
The instrument results before and after any adjustment are found in the associated pre and post stability checks.

### TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/IEC 17025 Standards.  
Simulator temperatures are traceable to NIST. Thermometer temperatures are checked with NIST traceable Eutechnics 4400 digital thermometers calibrated by Precision Metrology in accordance with ISO/IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the uses of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/IEC 17025 standards.  
This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

FDLE/ATP Form 69 April 2020  
Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality

09/24/2020 Date *Israel Soto*  
ISRAEL SOTO,  
Department Inspector

2020.09.  
25  
11:18:25  
-04'00'

**FDLE** REQUEST FOR REGISTRATION  
Florida Department of  
Law Enforcement

MAKE AND MODEL OF INSTRUMENT: Intoxilyzer 8000

SERIAL NUMBER: 80-000959

OWNING AGENCY: Gadsden County Sheriff Office

DATE OF DEPARTMENT INSPECTION: 09-24-2020

AGENCY INSPECTOR: Fredrick Hening

ADDRESS: 339 E Jefferson St

CITY, STATE, ZIP: Quincy, FL, 32351

TELEPHONE NUMBER: (850) 627-9233

FAX NUMBER: \_\_\_\_\_

EMAIL ADDRESS (if available): FredrickHening@flhsmv.gov

***For Program Office Use Only:***

- Registration Issued 2020.09.2 *BK* 5 1230:56 0400
- Instrument Added to Evidentiary Instrument Database
- Instrument Added to Monthly Statistics Database 2020.09.25
- Contact Information Added to Instrument Database -0400

*MX*