



INSTRUMENT PROCESSING SHEET

Agency Perry PDS/N 80-005337

Florida Department of Law Enforcement

Date In 02/28/2020 DI Completion Date 3/16/20 Ship P/U H/D CMI EE

Intake Performed By <u>RAW</u> <input checked="" type="checkbox"/> Annual <input type="checkbox"/> Registration <input checked="" type="checkbox"/> Return from CMI / EE Visual Inspection: <input checked="" type="checkbox"/> Case <input checked="" type="checkbox"/> Handle <input checked="" type="checkbox"/> Keyboard <input checked="" type="checkbox"/> Dry Gas Shelf <input checked="" type="checkbox"/> Feet <input checked="" type="checkbox"/> Breath Tube <input checked="" type="checkbox"/> Ports <input checked="" type="checkbox"/> Screws Tight Other Equipment/ Accessories: <input checked="" type="checkbox"/> Power cord <input type="checkbox"/> Printer Cable <input checked="" type="checkbox"/> Static Bag <input type="checkbox"/> 12V DC Cable Notes: _____	Quality Checks Performed By <u>JS</u> <input checked="" type="checkbox"/> Breath Tube Screen <input checked="" type="checkbox"/> Replace External O-Rings <input checked="" type="checkbox"/> Instrument Set Up Verified <input checked="" type="checkbox"/> R-Value <u>209</u> <input checked="" type="checkbox"/> Flow Verification (L/s) Flow Column # <u>ATP-103</u> 32 mm <u>0.160</u> (.139 - .169) 36 mm <u>0.171</u> (.156 - .190) 53 mm <u>0.238</u> (.228 - .278) 103 mm <u>0.492</u> (.447 - .547) <input checked="" type="checkbox"/> Barometric Pressure Check Gauge ID # <u>28427</u> <input checked="" type="checkbox"/> Stability Checks	Flow Calibration Performed By _____ Flow Column # _____ <input type="checkbox"/> 5L/min - 17mm <input type="checkbox"/> 15L/min - 53mm <input type="checkbox"/> 30L/min - 103mm <input type="checkbox"/> R-Value _____ <input type="checkbox"/> Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)
---	--	--

Final Release Date

FDLE FDLE

~~MAR 16 2020~~ ~~16 2020~~ DC

Alcohol Testing Testing Program

Simulator	Serial #	Lot #/Exp
0.050	SD1018	201905A 05-17-2021
0.080	SD3962	201905B 05-17-2021
0.200	G2078	201904D 04-30-2021
0.080 DGS	N/A	AG931603 11-12-2021

Maintenance Performed By _____

Battery Replacement
 Dry Gas Regulator Replacement
 Breath Tube Replacement
 Other _____

Temperature Checks Performed By SP

Lab Temp °C 22.1
 External Digital Therm. ID#: 300502
 34°C +-2 Serial #: SD1021
 34°C +-2 Serial #: DR1275
 34°C +-2 Serial #: SD1011

Calibration Adjustment Performed By SP

Barometric Pressure Gauge 1019 ID # 30793

Simulator	Serial Number	Lot Number	Expiration
0.000	G8144	N/A	N/A
0.040	G2403	19080	3-4-21
0.100	G2879	19160	7-9-21
0.200	SP G3709	19040	1-29-21
0.300	G8149	19010	1-3-21
0.080 DGS	N/A	08819080A1	6-5-21

Post Calibration Adjustment Stability Checks

Simulator	Serial Number	Lot Number	Expiration
0.050	SD1021	201905A	5-14-21
0.080	DR1275	201905B	5-14-21
0.200	SD1011	201904D	4-30-21
0.080 DGS	N/A	AG931603	11-12-21

Department Inspection Performed By SP

Barometric Pressure ID# 28421
 Gauge 1024 Instrument 1022
 Mouth Alcohol Solution Lot # 2019-B
 Acetone Stock Solution Lot # 2019-A

Simulator	Serial Number
0.000	G8144
Interferent	DR3855 SP DR3855
0.050	SD1021
0.080	DR1275
0.200	SD1011

Notes/Suggested Service: _____

Attachments

Form 41 Post-Stability Checks
 Stability Checks Flow Calibration
 Calibration Certificate Form 40
 Calibration Adjustment Other _____

Instrument Complies with Chapter 11D-8, FAC
 Instrument Does Not Comply with Chapter 11D-8, FAC

Return to/Place into Evidentiary Use
 Remain Out of Evidentiary Use

Conduct an Agency Inspection Before Evidentiary Use

SPM 3/16/20 Brett Kivland 3/16/2020
 Tech Review / Date Admin Review / Date



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
2729 Fort Knox Blvd.
Bldg. 2, Suite 1300
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-005337, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-005337</u>	UNCERTAINTY* ±	
Owning Agency:	<u>PERRY POLICE DEPT</u>	0.050 g/ 210 L	0.004
Calibration Date:	<u>03/16/2020</u>	0.080 g/ 210 L	0.005
Calibration Time:	<u>12:50</u>	0.200 g/ 210 L	0.007
		0.080 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration. *Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence (k=3).

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Thermometer temperatures are checked with NIST traceable Eutechnics 4400 digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the uses of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

03/16/2020

Date

Shayla D Platt

SHAYLA D PLATT,
Department Inspector

FDLE/ATP Form 69 January 2020

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality

DDM BK 3/16/2020

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: PERRY POLICE DEPT
Time of Inspection: 12:50

Date of Inspection: 03/16/2020

Serial Number: 80-005337
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:201905A Exp: 05/14/2021	0.08g/210L Test (g/210L) Lot#:201905B Exp: 05/14/2021	0.20g/210L Test (g/210L) Lot#:201904D Exp: 04/30/2021	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG931603 Exp: 11/12/2021
0.000	0.048	0.078	0.199	0.080
0.000	0.048	0.078	0.200	0.080
0.000	0.048	0.078	0.200	0.080
0.000	0.048	0.078	0.200	0.079
0.000	0.048	0.078	0.201	0.079
0.000	0.048	0.079	0.201	0.079
0.000	0.049	0.079	0.201	0.079
0.000	0.048	0.079	0.200	0.079
0.000	0.048	0.079	0.200	0.079
0.000	0.049	0.079	0.200	0.079

Standard Deviations	0.0004	0.0005	0.0006	0.0004
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0004 Number of Simulators Used: 5

Remarks:

BAM
BK
3/16/2020

The above instrument complies () does not comply () with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Shayla Platt

SHAYLA D PLATT

Signature and Printed Name

03/16/2020
Date

Stability Checks

PERRY POLICE DEPT
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-005337
 02/28/2020
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	14:01
Control Test	0.042	14:01
Air Blank	0.000	14:02
Control Test	0.043	14:02
Air Blank	0.000	14:03
Control Test	0.042	14:04
Air Blank	0.000	14:04
Control Test Stats		
Average	0.0423	
Std Dev	0.0006	
Rel Std Dev(%)	1.3638	

SS

Operator's Signature

PERRY POLICE DEPT
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-005337
 02/28/2020
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	14:06
Control Test	0.079	14:06
Air Blank	0.000	14:07
Control Test	0.078	14:07
Air Blank	0.000	14:08
Control Test	0.079	14:09
Air Blank	0.000	14:09
Control Test Stats		
Average	0.0787	
Std Dev	0.0006	
Rel Std Dev(%)	0.7339	

Wet

SS

Operator's Signature

PERRY POLICE DEPT
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-005337
 02/28/2020
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	14:11
Control Test	0.198	14:11
Air Blank	0.000	14:12
Control Test	0.199	14:13
Air Blank	0.000	14:13
Control Test	0.199	14:14
Air Blank	0.000	14:14
Control Test Stats		
Average	0.1987	
Std Dev	0.0006	
Rel Std Dev(%)	0.2906	

SS

Operator's Signature

PERRY POLICE DEPT
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-005337
 02/28/2020
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	14:17
Control Test	0.080	14:17
Air Blank	0.000	14:18
Control Test	0.079	14:18
Air Blank	0.000	14:19
Control Test	0.079	14:19
Air Blank	0.000	14:19
Control Test Stats		
Average	0.0793	
Std Dev	0.0006	
Rel Std Dev(%)	0.7277	

Dry

SS

Operator's Signature

GMM
 TBK
 3/16/2020

PERRY POLICE DEPT
Intoxilyzer - Alcohol Analyzer
Model 8000
03/13/2020
SN 80-005537
09:22:11

Auto Calibration
Max Power Res Value = 45
Auto Range Res Value = 28

Sol Value = 0.000 g/210L ***
Fit Value = 0.0000 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12540, Sum Io = 13430
<<<<< CHANNEL 1 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.0530 (-0.0140)
Sample #2 = 0.0690 (-0.0140)
Sample #3 = 0.0030 (0.0280)
Sample #4 = 0.0530 (0.0180)
Avg % Abs = 0.0417 (0.0107)
STD DEV = 0.0344 (0.0219)
REL STD DEV = 82.629 (205.661)

<<<<< CHANNEL 2 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.1210 (-0.0150)
Sample #2 = 0.1220 (-0.0060)
Sample #3 = 0.0990 (0.0020)
Sample #4 = 0.1260 (-0.0070)
Avg % Abs = 0.1157 (-0.0037)
STD DEV = 0.0146 (0.0049)
REL STD DEV = 12.598 (134.533)

Sol Value = 0.040 g/210L ***
Fit Value = 0.1905 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12534, Sum Io = 13429
<<<<< CHANNEL 1 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.7650 (-0.0160)
Sample #2 = 0.8050 (-0.0190)
Sample #3 = 0.7500 (0.0070)
Sample #4 = 0.7890 (0.0000)
Avg % Abs = 0.7813 (-0.0040)
STD DEV = 0.0283 (0.0135)
REL STD DEV = 3.621 (336.341)

<<<<< CHANNEL 2 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 6.7270 (-0.0170)
Sample #2 = 6.6850 (0.0470)
Sample #3 = 6.7070 (0.0370)
Sample #4 = 6.6800 (0.0620)
Avg % Abs = 6.6907 (0.0487)
STD DEV = 0.0144 (0.0126)
REL STD DEV = 0.215 (25.856)

Sol Value = 0.300 g/210L ***
Fit Value = 1.4286 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12530, Sum Io = 13427
<<<<< CHANNEL 1 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 5.2880 (-0.0050)
Sample #2 = 5.2510 (0.0640)
Sample #3 = 5.2330 (0.0710)
Sample #4 = 5.2950 (0.0440)
Avg % Abs = 5.2597 (0.0597)
STD DEV = 0.0319 (0.0140)
REL STD DEV = 0.506 (23.484)

<<<<< CHANNEL 2 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 9.6890 (0.0000)
Sample #2 = 9.6320 (0.1020)
Sample #3 = 9.6130 (0.1150)
Sample #4 = 9.6680 (0.1000)
Avg % Abs = 9.6377 (0.1057)
STD DEV = 0.0279 (0.0081)
REL STD DEV = 0.290 (7.708)

Sol Value = 0.200 g/210L ***
Fit Value = 0.9524 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12534, Sum Io = 13427
<<<<< CHANNEL 1 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 3.6430 (-0.0210)
Sample #2 = 3.5650 (0.0410)
Sample #3 = 3.5940 (0.0280)
Sample #4 = 3.5670 (0.0710)
Avg % Abs = 3.5753 (0.0467)
STD DEV = 0.0162 (0.0221)
REL STD DEV = 0.453 (47.256)

***** AUTO CAL DATA *****
<<<<< CHANNEL 1 >>>>>
Sol Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.042
Std Dev = 0.03 Rel Std Dev = 82.63
Sol Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 0.781
Std Dev = 0.03 Rel Std Dev = 3.62
Sol Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 1.850
Std Dev = 0.04 Rel Std Dev = 2.04
Sol Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 3.575
Std Dev = 0.02 Rel Std Dev = 0.45
Sol Val = 1.4286 mg/l or 0.300 g/210L
% Abs = 5.260
Std Dev = 0.03 Rel Std Dev = 0.61
Zero Order Coef = -121.56
First Order Coef = 2592.23
Second Order Coef = 28.17
Standard Deviation = 14.497541

<<<<< CHANNEL 2 >>>>>
Sol Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.116
Std Dev = 0.01 Rel Std Dev = 12.60
Sol Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 1.532
Std Dev = 0.02 Rel Std Dev = 1.28
Sol Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 3.544
Std Dev = 0.05 Rel Std Dev = 1.45
Sol Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 6.691
Std Dev = 0.01 Rel Std Dev = 0.21
Sol Val = 1.4286 mg/l or 0.300 g/210L
% Abs = 9.638
Std Dev = 0.03 Rel Std Dev = 0.29
Zero Order Coef = -159.93
First Order Coef = 1325.90
Second Order Coef = 17.98
Standard Deviation = 6.782514

Sol Value = 0.080 g/210L ***
Fit Value = 0.3810 mg/l %%%
Samples Taken = 4, Discarded = 1
***** CHANNEL 1 *****
Sample #1 = 3436.00
Sample #2 = 3418.00
Sample #3 = 3411.00
Sample #4 = 3411.00
Average Result = 3413.3333
STD DEV = 4.0415
REL STD DEV = 0.118

***** CHANNEL 2 *****
Sample #1 = 3216.00
Sample #2 = 3202.00
Sample #3 = 3221.00
Sample #4 = 3234.00
Average Result = 3219.0000
STD DEV = 16.0935
REL STD DEV = 0.500

Dry Gas H2O Adjust Results *****
Barometric Pressure = 1019
3 um H2O Adjust (mg/l*10.000) = 386
9 um H2O Adjust (mg/l*10.000) = 590
***** AUTO CAL PASS *****

Sol Value = 0.040 g/210L ***
Fit Value = 0.1905 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12534, Sum Io = 13429
<<<<< CHANNEL 1 >>>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.7650 (-0.0160)
Sample #2 = 0.8050 (-0.0190)
Sample #3 = 0.7500 (0.0070)
Sample #4 = 0.7890 (0.0000)
Avg % Abs = 0.7813 (-0.0040)
STD DEV = 0.0283 (0.0135)
REL STD DEV = 3.621 (336.341)

Poem BK 3/16/2020

Post Cal-Adjust Stability Checks #80-005337

PERRY POLICE DEPT
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005337
03/16/2020
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:43
Control Test	0.049	09:44
Air Blank	0.000	09:44
Control Test	0.048	09:45
Air Blank	0.000	09:45
Control Test	0.049	09:46
Air Blank	0.000	09:47
Control Test Stats		
Average	0.0487	
Std Dev	0.0006	
Rel Std Dev(%)	1.1863	

SP

Operator's Signature

PERRY POLICE DEPT
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005337
03/16/2020
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:48
Control Test	0.079	09:49
Air Blank	0.000	09:50
Control Test	0.079	09:50
Air Blank	0.000	09:51
Control Test	0.079	09:51
Air Blank	0.000	09:52
Control Test Stats		
Average	0.0790	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

SP

Operator's Signature

PERRY POLICE DEPT
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005337
03/16/2020
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:53
Control Test	0.200	09:54
Air Blank	0.000	09:54
Control Test	0.200	09:55
Air Blank	0.000	09:56
Control Test	0.200	09:56
Air Blank	0.000	09:57
Control Test Stats		
Average	0.2000	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

SP

Operator's Signature

PERRY POLICE DEPT
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005337
03/16/2020
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:59
Control Test	0.081	09:59
Air Blank	0.000	09:59
Control Test	0.080	10:00
Air Blank	0.000	10:00
Control Test	0.081	10:00
Air Blank	0.000	10:01
Control Test Stats		
Average	0.0807	
Std Dev	0.0006	
Rel Std Dev(%)	0.7157	

SP

SP

Operator's Signature

PPM 1SK 3/16/2020

Return Material Authorization

Ship to: CMI, Inc.

Enforcement Electronics

Shipment to repair facility authorized by: Lin Gray on 12/17/2019

Items Returned: Instrument Supplies Other Describe: _____

Instrument Model: 8000 Serial Number: 80-005337

Bill To Address:
Lin Gray
Perry Police Department
211 S. Washington Street
Perry, FL 32347

Ship to Address:
FDLE Off-Site Mail Facility
c/o Florida Department of Law Enforcemen
Alcohol Testing Program
813B Lake Bradford Road
Tallahassee FL 32304

Reason for Return:

Voltage across flow sensor at rest (R value) is below 100. Flow sensor needs replacement

Please choose one of the following options:

- 1. I _____, authorize all repairs.
- 2. I _____, authorize repairs up to \$ _____.
- 3. I require an estimate **BEFORE** any repairs will be authorized and/ or conducted.

Please contact: Name: Lin Gray

Phone #: 850-838-6596 Email: lin.gray@perrypolice.net

ATP Contact Name: Patrick Murphy ATP Email: patrickmurphy@fdle.state.fl.us



INSTRUMENT PROCESSING SHEET

Agency Perry PD

S/N 80-005337

Florida Department of Law Enforcement

Date In 12/9/2019 DI Completion Date _____

Ship P/U H/D CMI EE

Intake Performed By <u>DP</u> <input checked="" type="checkbox"/> Annual <input type="checkbox"/> Registration <input type="checkbox"/> Return from CMI / EE Visual Inspection: <input checked="" type="checkbox"/> Case <input checked="" type="checkbox"/> Handle <input checked="" type="checkbox"/> Keyboard <input checked="" type="checkbox"/> Dry Gas Shelf <input checked="" type="checkbox"/> Feet <input checked="" type="checkbox"/> Breath Tube <input checked="" type="checkbox"/> Ports <input checked="" type="checkbox"/> Screws Tight Other Equipment/ Accessories: <input type="checkbox"/> Power cord <input type="checkbox"/> Printer Cable <input type="checkbox"/> Static Bag <input type="checkbox"/> 12V DC Cable Notes: _____ _____ _____	Quality Checks Performed By _____ <input checked="" type="checkbox"/> Breath Tube Screen <input checked="" type="checkbox"/> Replace External O-Rings <input checked="" type="checkbox"/> Instrument Set Up Verified <input type="checkbox"/> R-Value <u>90</u> <input type="checkbox"/> Flow Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547) <input type="checkbox"/> Barometric Pressure Check Gauge ID # _____ <input type="checkbox"/> Stability Checks <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td></td> <td></td> </tr> <tr> <td>0.080</td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> </tr> </tbody> </table>	Simulator	Serial #	Lot #/Exp	0.050			0.080			0.200			0.080 DGS	N/A		Flow Calibration Performed By _____ Flow Column # _____ <input type="checkbox"/> 5L/min - 17mm <input type="checkbox"/> 15L/min - 53mm <input type="checkbox"/> 30L/min - 103mm <input type="checkbox"/> R-Value _____ <input type="checkbox"/> Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547) Maintenance Performed By _____ <input type="checkbox"/> Battery Replacement <input type="checkbox"/> Dry Gas Regulator Replacement <input type="checkbox"/> Breath Tube Replacement <input type="checkbox"/> Other _____ Temperature Checks Performed By _____ <input type="checkbox"/> Lab Temp °C _____ External Digital Therm. ID#: _____ <input type="checkbox"/> 34°C +-2 Serial #: _____ <input type="checkbox"/> 34°C +-2 Serial #: _____ <input type="checkbox"/> 34°C +-2 Serial #: _____
Simulator	Serial #	Lot #/Exp															
0.050																	
0.080																	
0.200																	
0.080 DGS	N/A																
Final Release Date _____ _____																	

Calibration Adjustment Performed By _____ Barometric Pressure Gauge _____ ID # _____ <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td></td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.100</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.300</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table> <input type="checkbox"/> Post Calibration Adjustment Stability Checks <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table>	Simulator	Serial Number	Lot Number	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			Simulator	Serial Number	Lot Number	Expiration	0.050				0.080				0.200				0.080 DGS	N/A			Department Inspection Performed By _____ Barometric Pressure ID# _____ Gauge _____ Instrument _____ Mouth Alcohol Solution Lot # _____ Acetone Stock Solution Lot # _____ <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td></td> </tr> <tr> <td>Interferent</td> <td></td> </tr> <tr> <td>0.050</td> <td></td> </tr> <tr> <td>0.080</td> <td></td> </tr> <tr> <td>0.200</td> <td></td> </tr> </tbody> </table> Attachments <input type="checkbox"/> Form 41 <input type="checkbox"/> Post-Stability Checks <input type="checkbox"/> Stability Checks <input type="checkbox"/> Flow Calibration <input type="checkbox"/> Calibration Certificate <input type="checkbox"/> Form 40 <input type="checkbox"/> Calibration Adjustment <input type="checkbox"/> Other _____	Simulator	Serial Number	0.000		Interferent		0.050		0.080		0.200	
Simulator	Serial Number	Lot Number	Expiration																																																										
0.000		N/A	N/A																																																										
0.040																																																													
0.100																																																													
0.200																																																													
0.300																																																													
0.080 DGS	N/A																																																												
Simulator	Serial Number	Lot Number	Expiration																																																										
0.050																																																													
0.080																																																													
0.200																																																													
0.080 DGS	N/A																																																												
Simulator	Serial Number																																																												
0.000																																																													
Interferent																																																													
0.050																																																													
0.080																																																													
0.200																																																													

Notes/Suggested Service: _____ _____ _____ _____ _____ _____	<input type="checkbox"/> Instrument Complies with Chapter 11D-8, FAC <input type="checkbox"/> Instrument Does Not Comply with Chapter 11D-8, FAC <input type="checkbox"/> Return to/Place into Evidentiary Use <input type="checkbox"/> Remain Out of Evidentiary Use <input type="checkbox"/> Conduct an Agency Inspection Before Evidentiary Use _____ Tech Review / Date _____ Admin Review / Date _____
---	---

Florida Department of Law Enforcement Alcohol Testing Program

AGENCY INSPECTION REPORT - INTOXILYZER 8000

Agency: PERRY POLICE DEPT
Time of Inspection: 09:59

Date of Inspection: 12/09/2019

Serial Number: 80-005337
Software: 8100.27

Check or Test	YES	NO
Date and/or Time Adjusted		No
Diagnostic Check (Pre-Inspection): OK		No
Alcohol Free Subject Test: 0.000		No
Mouth Alcohol Test: Slope Not Met		No
Interferent Detect Test: Interferent Detect		No
Diagnostic Check (Post-Inspection): OK		No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#: _____ Exp: _____	0.08g/210L Test (g/210L) Lot#: _____ Exp: _____	0.20g/210L Test (g/210L) Lot#: _____ Exp: _____	0.08 g/210L Dry Gas Std Test (g/210L) Lot#: _____ Exp: _____

Number of Simulators Used: _____

Remarks:

SKIPPED AI TO OPERATE INSTRUMENT

The above instrument complies () does not comply () with Chapter 11D-8, FAC.

I certify that I hold a valid Florida Department of Law Enforcement Agency Inspector Permit and that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Signature and Printed Name

12/09/2019
Date