

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: SANTA ROSA COUNTY SO
Time of Inspection: 13:26

Date of Inspection: 08/06/2021

Serial Number: 80-001303
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202010A Exp: 10/05/2022	0.08g/210L Test (g/210L) Lot#:202010B Exp: 10/05/2022	0.20g/210L Test (g/210L) Lot#:202010D Exp: 10/06/2022	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG011102 Exp: 04/20/2022
0.000	0.049	0.079	0.201	0.078
0.000	0.049	0.079	0.201	0.078
0.000	0.049	0.079	0.201	0.078
0.000	0.050	0.079	0.201	0.078
0.000	0.049	0.079	0.201	0.078
0.000	0.049	0.080	0.202	0.078
0.000	0.050	0.080	0.202	0.078
0.000	0.050	0.081	0.202	0.078
0.000	0.051	0.080	0.202	0.078
0.000	0.050	0.080	0.203	0.078


Standard Deviations	0.0006	0.0006	0.0006	0.0000
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0004 Number of Simulators Used: 5

Remarks:

The above instrument complies () does not comply () with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



ISRAEL SOTO

Signature and Printed Name

08/06/2021
Date

stability checks

SANTA ROSA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001303
08/06/2021
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:48
Control Test	0.049	11:48
Air Blank	0.000	11:49
Control Test	0.048	11:50
Air Blank	0.000	11:50
Control Test	0.048	11:51
Air Blank	0.000	11:52
Control Test Stats		
Average	0.0483	
Std Dev	0.0006	
Rel Std Dev(%)	1.1945	



Operator's Signature

SANTA ROSA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001303
08/06/2021
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:53
Control Test	0.079	11:53
Air Blank	0.000	11:54
Control Test	0.078	11:55
Air Blank	0.000	11:55
Control Test	0.079	11:56
Air Blank	0.000	11:57
Control Test Stats		
Average	0.0787	
Std Dev	0.0006	
Rel Std Dev(%)	0.7339	

wet



Operator's Signature

SANTA ROSA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001303
08/06/2021
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:57
Control Test	0.201	11:58
Air Blank	0.000	11:59
Control Test	0.199	11:59
Air Blank	0.000	12:00
Control Test	0.199	12:01
Air Blank	0.000	12:01
Control Test Stats		
Average	0.1997	
Std Dev	0.0012	
Rel Std Dev(%)	0.5783	



Operator's Signature

SANTA ROSA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001303
08/06/2021
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:02
Control Test	0.078	12:03
Air Blank	0.000	12:03
Control Test	0.078	12:04
Air Blank	0.000	12:04
Control Test	0.079	12:04
Air Blank	0.000	12:05
Control Test Stats		
Average	0.0783	
Std Dev	0.0006	
Rel Std Dev(%)	0.7370	

Dry



Operator's Signature



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
2729 Fort Knox Blvd.
Bldg. 2, Suite 1300
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-001303, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-001303</u>	UNCERTAINTY* ±	
Owning Agency:	<u>SANTA ROSA COUNTY SO</u>	0.050 g/ 210 L	0.005
Calibration Date:	<u>08/06/2021</u>	0.080 g/ 210 L	0.004
Calibration Time:	<u>13:26</u>	0.200 g/ 210 L	0.007
		0.080 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence (k=3).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Thermometer temperatures are checked with NIST traceable Eutechnics 4400 digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the uses of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

Israel Soto
Digitally signed by Israel Soto
Date: 2021.08.06 14:40:12
-04'00'

08/06/2021

Date

ISRAEL SOTO,

Department Inspector

FDLE/ATP Form 69 January 2021

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality



INSTRUMENT PROCESSING SHEET

Agency Santa Rosa County Sheriff's Office

S/N 80-001303

Florida Department of Law Enforcement

Date In 11/8/2021

DI Completion Date 11/9/2021

Ship P/U H/D CMI EE

Intake, Quality Checks, Flow Calibration, Maintenance, DI Temp. Checks, Notes

Calibration Adjustment, Barometric Pressure Gauge, Post Calibration Adjustment Stability Checks

Department Inspection, Barometric Pressure, Mouth Alcohol Solution, Acetone Stock Solution, Attachments

Notes/Suggested Service: Values on stability checks nominal, performed extra stability checks...

Instrument Complies with Chapter 11D-8, FAC, Return to/Place into Evidentiary Use, Conduct an Agency Inspection Before Evidentiary Use

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: SANTA ROSA COUNTY SO
Time of Inspection: 10:44

Date of Inspection: 11/09/2021

Serial Number: 80-001303
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202010A Exp: 10/05/2022	0.08g/210L Test (g/210L) Lot#:202010B Exp: 10/05/2022	0.20g/210L Test (g/210L) Lot#:202010D Exp: 10/06/2022	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG113403 Exp: 05/14/2023
0.000	0.049	0.079	0.199	0.080
0.000	0.049	0.079	0.198	0.080
0.000	0.049	0.079	0.199	0.080
0.000	0.050	0.080	0.199	0.080
0.000	0.050	0.080	0.199	0.080
0.000	0.050	0.079	0.199	0.080
0.000	0.050	0.080	0.199	0.079
0.000	0.050	0.080	0.199	0.079
0.000	0.051	0.080	0.199	0.079
0.000	0.050	0.080	0.199	0.079


Standard Deviations	0.0006	0.0005	0.0003	0.0005
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0004 Number of Simulators Used: 5

Remarks:

The above instrument complies () does not comply () with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



ISRAEL SOTO

Signature and Printed Name

11/09/2021
Date

Stability Checks

SANTA ROSA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001303
11/08/2021
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:12
Control Test	0.051	13:13
Air Blank	0.000	13:14
Control Test	0.050	13:14
Air Blank	0.000	13:15
Control Test	0.050	13:16
Air Blank	0.000	13:16
Control Test Stats		
Average	0.0503	
Std Dev	0.0006	
Rel Std Dev(%)	1.1471	

SANTA ROSA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001303
11/08/2021
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:17
Control Test	0.080	13:18
Air Blank	0.000	13:18
Control Test	0.079	13:19
Air Blank	0.000	13:20
Control Test	0.079	13:20
Air Blank	0.000	13:21
Control Test Stats		
Average	0.0793	
Std Dev	0.0006	
Rel Std Dev(%)	0.7277	


SANTA ROSA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001303
11/08/2021
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:23
Control Test	0.202	13:24
Air Blank	0.000	13:24
Control Test	0.201	13:25
Air Blank	0.000	13:25
Control Test	0.201	13:26
Air Blank	0.000	13:27
Control Test Stats		
Average	0.2013	
Std Dev	0.0006	
Rel Std Dev(%)	0.2868	

wet



Operator's Signature



Operator's Signature




Operator's Signature

SANTA ROSA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001303
11/08/2021
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:28
Control Test	0.080	13:29
Air Blank	0.000	13:29
Control Test	0.080	13:30
Air Blank	0.000	13:30
Control Test	0.080	13:30
Air Blank	0.000	13:31
Control Test Stats		
Average	0.0800	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

Dry



Operator's Signature



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
2729 Fort Knox Blvd.
Bldg. 2, Suite 1300
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-001303, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-001303</u>	UNCERTAINTY* ±	
Owning Agency:	<u>SANTA ROSA COUNTY SO</u>	0.050 g/ 210 L	0.005
Calibration Date:	<u>11/09/2021</u>	0.080 g/ 210 L	0.004
Calibration Time:	<u>10:44</u>	0.200 g/ 210 L	0.007
		0.080 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence (k=3).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Thermometer temperatures are checked with NIST traceable Eutechnics 4400 digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the uses of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

Israel Soto
Digitally signed by Israel Soto
Date: 2021.11.09 11:56:51
-05'00'

11/09/2021

Date

ISRAEL SOTO,

Department Inspector

FDLE/ATP Form 69 January 2021

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality

Optical Bench Calibration

Adjust

SANTA ROSA COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-001303
 11/09/2021 07:54:29

Auto Calibration
 Max Power Res Value = 41
 Auto Range Res Value = 27

Sol Value = 0.000 g/210L ***
 Fit value = 0.0000 mg/l %%%
 Samples Taken = 4, Discarded = 1
 Sum lo = 12726, Sum hi = 13418

Channel 1 Data:
 Sample % Abs (% Abs Ref)
 Sample #1 = 0.2300 (-0.0180)
 Sample #2 = 0.1860 (0.0650)
 Sample #3 = 0.1720 (0.1260)
 Sample #4 = 0.1750 (0.1700)
 Avg % Abs = 0.1777 (0.1203)
 STD DEV = 0.0074 (0.0527)
 REL STD DEV = 4.149 (43.819)

Channel 2 Data:
 Sample % Abs (% Abs Ref)
 Sample #1 = 0.1550 (-0.0070)
 Sample #2 = 0.1500 (0.0050)
 Sample #3 = 0.1380 (0.0350)
 Sample #4 = 0.1550 (0.0480)
 Avg % Abs = 0.1477 (0.0297)
 STD DEV = 0.0087 (0.0215)
 REL STD DEV = 5.917 (72.478)

Sol Value = 0.040 g/210L ***
 Fit value = 0.1905 mg/l %%%
 Samples Taken = 4, Discarded = 1
 Sum lo = 12685, Sum hi = 13400

Channel 1 Data:
 Sample % Abs (% Abs Ref)
 Sample #1 = 0.9390 (0.0050)
 Sample #2 = 0.9330 (0.0440)
 Sample #3 = 0.9180 (0.0830)
 Sample #4 = 0.9140 (0.1100)
 Avg % Abs = 0.9217 (0.0790)
 STD DEV = 0.0100 (0.0332)
 REL STD DEV = 1.087 (42.002)

Channel 2 Data:
 Sample % Abs (% Abs Ref)
 Sample #1 = 1.5900 (0.0040)
 Sample #2 = 1.5720 (0.0290)
 Sample #3 = 1.5840 (0.0470)
 Sample #4 = 1.5800 (0.0560)
 Avg % Abs = 1.5787 (0.0440)
 STD DEV = 0.0061 (0.0137)
 REL STD DEV = 0.387 (31.245)

Sol Value = 0.100 g/210L ***
 Fit value = 0.4762 mg/l %%%
 Samples Taken = 4, Discarded = 1
 Sum lo = 12558, Sum hi = 13384

Channel 1 Data:
 Sample % Abs (% Abs Ref)
 Sample #1 = 2.0290 (-0.0180)
 Sample #2 = 2.0010 (0.0600)
 Sample #3 = 1.9940 (0.0670)
 Sample #4 = 2.0170 (0.0890)
 Avg % Abs = 2.0040 (0.0729)
 STD DEV = 0.0118 (0.0151)
 REL STD DEV = 0.588 (21.018)

Channel 2 Data:
 Sample % Abs (% Abs Ref)
 Sample #1 = 3.6630 (-0.0050)
 Sample #2 = 3.6610 (0.0370)
 Sample #3 = 3.6400 (0.0360)
 Sample #4 = 3.6580 (0.0490)
 Avg % Abs = 3.6530 (0.0407)
 STD DEV = 0.0114 (0.0072)
 REL STD DEV = 0.311 (17.789)

Sol Value = 0.200 g/210L ***
 Fit value = 0.9524 mg/l %%%
 Samples Taken = 4, Discarded = 1
 Sum lo = 12626, Sum hi = 13358

Channel 1 Data:
 Sample % Abs (% Abs Ref)
 Sample #1 = 3.8220 (-0.0130)
 Sample #2 = 3.7970 (0.0360)
 Sample #3 = 3.7990 (0.0780)
 Sample #4 = 3.7720 (0.0930)
 Avg % Abs = 3.7893 (0.0690)
 STD DEV = 0.0150 (0.0295)
 REL STD DEV = 0.397 (42.821)

Channel 2 Data:
 Sample % Abs (% Abs Ref)
 Sample #1 = 7.0180 (-0.0170)
 Sample #2 = 6.9470 (0.0490)
 Sample #3 = 6.9800 (0.0590)
 Sample #4 = 6.9510 (0.0730)
 Avg % Abs = 6.9593 (0.0603)
 STD DEV = 0.0180 (0.0121)
 REL STD DEV = 0.259 (19.981)

Sol Value = 0.300 g/210L ***
 Fit value = 1.4286 mg/l %%%
 Samples Taken = 4, Discarded = 1
 Sum lo = 12610, Sum hi = 13357

Channel 1 Data:
 Sample % Abs (% Abs Ref)
 Sample #1 = 5.4980 (-0.0110)
 Sample #2 = 5.5330 (0.0040)
 Sample #3 = 5.4850 (0.0820)
 Sample #4 = 5.4890 (0.0520)
 Avg % Abs = 5.5023 (0.0493)
 STD DEV = 0.0266 (0.0405)
 REL STD DEV = 0.484 (82.122)

Channel 2 Data:
 Sample % Abs (% Abs Ref)
 Sample #1 = 10.1080 (-0.0170)
 Sample #2 = 10.0910 (0.0230)
 Sample #3 = 10.0820 (0.0370)
 Sample #4 = 10.0750 (0.0390)
 Avg % Abs = 10.0827 (0.0330)
 STD DEV = 0.0080 (0.0087)
 REL STD DEV = 0.080 (26.418)

***** AUTO CAL DATA *****
 Channel 1 Data:
 Sol Val = 0.0000 mg/l or 0.000 g/210L
 % Abs = 0.178
 Std Dev = 0.01 Rel Std Dev = 4.15
 Sol Val = 0.1905 mg/l or 0.040 g/210L
 % Abs = 0.922
 Std Dev = 0.01 Rel Std Dev = 1.09
 Sol Val = 0.4762 mg/l or 0.100 g/210L
 % Abs = 2.004
 Std Dev = 0.01 Rel Std Dev = 0.59
 Sol Val = 0.9524 mg/l or 0.200 g/210L
 % Abs = 3.789
 Std Dev = 0.02 Rel Std Dev = 0.40
 Sol Val = 1.4286 mg/l or 0.300 g/210L
 % Abs = 5.502
 Std Dev = 0.03 Rel Std Dev = 0.48
 Zero Order Coef = -452.91
 First Order Coef = 2544.95
 Second Order Coef = 24.14
 Standard Deviation = 12.174147

Channel 2 Data:
 Sol Val = 0.0000 mg/l or 0.000 g/210L
 % Abs = 0.148
 Std Dev = 0.01 Rel Std Dev = 5.92
 Sol Val = 0.1905 mg/l or 0.040 g/210L
 % Abs = 1.579
 Std Dev = 0.01 Rel Std Dev = 0.39
 Sol Val = 0.4762 mg/l or 0.100 g/210L
 % Abs = 3.653
 Std Dev = 0.01 Rel Std Dev = 0.31
 Sol Val = 0.9524 mg/l or 0.200 g/210L
 % Abs = 6.959
 Std Dev = 0.02 Rel Std Dev = 0.26
 Sol Val = 1.4286 mg/l or 0.300 g/210L
 % Abs = 10.083
 Std Dev = 0.01 Rel Std Dev = 0.08
 Zero Order Coef = -192.99
 First Order Coef = 1309.56
 Second Order Coef = 12.52
 Standard Deviation = 2.782640

Optical Bench
Calibration Adjust
continued

sd

Post stability checks

Solution Stats Quadratic Fit Chan 1

Act	Fit	Residual
g/210L	g/210L	g/210L
0.000	0.000	-0.0000
0.040	0.040	-0.0002
0.100	0.100	0.0004
0.200	0.200	-0.0003
0.300	0.300	0.0001

Solution Stats Quadratic Fit Chan 2

Act	Fit	Residual
g/210L	g/210L	g/210L
0.000	0.000	-0.0000
0.040	0.040	-0.0000
0.100	0.100	0.0001
0.200	0.200	-0.0001
0.300	0.300	0.0000

Sol Value = 0.080 g/210L ***
 Fit value = 0.3810 mg/l %%%
 Samples Taken = 4, Discarded = 1
 ***** CHANNEL 1
 Sample #1 = 2902.00
 Sample #2 = 2866.00
 Sample #3 = 2833.00
 Sample #4 = 2868.00
 Average Result = 2855.6667
 STD DEV = 19.6554
 REL STD DEV = 0.688

 ***** CHANNEL 2
 Sample #1 = 3285.00
 Sample #2 = 3279.00
 Sample #3 = 3279.00
 Sample #4 = 3268.00
 Average Result = 3275.3333
 STD DEV = 6.3509
 REL STD DEV = 0.194

 Dry Gas H2O Adjust Results *****
 Barometric Pressure = 1024
 3 um H2O Adjust (mg/l*10,000) = 954
 9 um H2O Adjust (mg/l*10,000) = 534
 **** AUTO CAL PASS

SANTA ROSA COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-001303
 11/09/2021
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:34
Control Test	0.049	08:35
Air Blank	0.000	08:35
Control Test	0.048	08:36
Air Blank	0.000	08:37
Control Test	0.048	08:37
Air Blank	0.000	08:38
Control Test Stats		
Average	0.0483	
Std Dev	0.0006	
Rel Std Dev(%)	1.1945	

sd

 Operator's Signature

SANTA ROSA COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-001303
 11/09/2021
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:43
Control Test	0.198	08:44
Air Blank	0.000	08:44
Control Test	0.198	08:45
Air Blank	0.000	08:46
Control Test	0.198	08:46
Air Blank	0.000	08:47
Control Test Stats		
Average	0.1980	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

sd

 Operator's Signature

SANTA ROSA COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-001303
 11/09/2021
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:39
Control Test	0.079	08:39
Air Blank	0.000	08:40
Control Test	0.078	08:41
Air Blank	0.000	08:41
Control Test	0.078	08:42
Air Blank	0.000	08:42
Control Test Stats		
Average	0.0783	
Std Dev	0.0006	
Rel Std Dev(%)	0.7370	

wet

sd

 Operator's Signature

SANTA ROSA COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-001303
 11/09/2021
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:48
Control Test	0.081	08:48
Air Blank	0.000	08:49
Control Test	0.081	08:49
Air Blank	0.000	08:50
Control Test	0.081	08:50
Air Blank	0.000	08:50
Control Test Stats		
Average	0.0810	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

Dry

sd

 Operator's Signature

SANTA ROSA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001303
11/08/2021
Software: 8100.27

Extra stability checks

Test	g/210L	Time
Air Blank	0.000	13:32
Control Test	0.081	13:33
Air Blank	0.000	13:34
Control Test	0.080	13:34
Air Blank	0.000	13:35
Control Test	0.080	13:36
Air Blank	0.000	13:36
Control Test	0.079	13:37
Air Blank	0.000	13:37
Control Test	0.080	13:38
Air Blank	0.000	13:39
Control Test	0.080	13:39
Air Blank	0.000	13:40
Control Test	0.080	13:41
Air Blank	0.000	13:41
Control Test	0.079	13:42
Air Blank	0.000	13:42
Control Test	0.079	13:43
Air Blank	0.000	13:44
Control Test	0.079	13:44
Air Blank	0.000	13:45
Control Test	0.080	13:46
Air Blank	0.000	13:46
Control Test	0.079	13:47
Air Blank	0.000	13:47
Control Test	0.079	13:48
Air Blank	0.000	13:49
Control Test	0.080	13:49
Air Blank	0.000	13:50
Control Test	0.080	13:51
Air Blank	0.000	13:51
Control Test	0.080	13:52
Air Blank	0.000	13:52
Control Test	0.080	13:53
Air Blank	0.000	13:54
Control Test	0.080	13:54
Air Blank	0.000	13:55
Control Test	0.079	13:56
Air Blank	0.000	13:56
Control Test	0.080	13:57
Air Blank	0.000	13:57
Control Test Stats		
Average	0.0797	
Std Dev	0.0006	
Rel Std Dev(%)	0.7167	

wet


Operator's Signature