



INSTRUMENT PROCESSING SHEET

Agency Hillsborough CSO

S/N 80-000808

Florida Department of
Law Enforcement

Date In 3/11/2022 DI Completion Date 3/17/2022

☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake	By TDG	Quality Checks	By TDG	Date <u>3/14/2022</u>	Flow Calibration	By	Date																																								
☒ Annual ☐ Registration ☐ Return from CMI / EE Visual Inspection: ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/ Accessories: ☐ Power cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes: _____ _____ _____ _____ _____ _____ _____ _____		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value <u>187</u> ☒ Flow Verification (L/s) Flow Column # <u>ATP104</u> 32 mm <u>0.152</u> (.139 - .169) 36 mm <u>0.167</u> (.156 - .190) 53 mm <u>0.242</u> (.228 - .278) 103 mm <u>0.500</u> (.447 - .547) ☒ Barometric Pressure Check Gauge ID # <u>68639</u> ☒ Stability Checks			Flow Column # _____ ☐ 5L/min - 17mm ☐ 15L/min - 53mm ☐ 30L/min - 103mm ☐ R-Value _____ ☐ Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)																																										
		<table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td>MP6286</td> <td>202010A 10/05/2022</td> </tr> <tr> <td>0.080</td> <td>MP6287</td> <td>202010B 10/05/2022</td> </tr> <tr> <td>0.200</td> <td>MP6288</td> <td>202010D 10/06/2022</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>AG115904 06/08/2023</td> </tr> </tbody> </table>			Simulator	Serial #	Lot #/Exp	0.050	MP6286	202010A 10/05/2022	0.080	MP6287	202010B 10/05/2022	0.200	MP6288	202010D 10/06/2022	0.080 DGS	N/A	AG115904 06/08/2023	Maintenance By _____ ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other _____ _____ _____ _____ _____																											
Simulator	Serial #	Lot #/Exp																																													
0.050	MP6286	202010A 10/05/2022																																													
0.080	MP6287	202010B 10/05/2022																																													
0.200	MP6288	202010D 10/06/2022																																													
0.080 DGS	N/A	AG115904 06/08/2023																																													
Calibration Adjustment By TDG _____ Barometric Pressure Gauge <u>1018</u> ID # <u>28663</u>				Department Inspection By TDG _____ Barometric Pressure ID# <u>28199</u> Gauge <u>1017</u> Instrument <u>1018</u> Mouth Alcohol Solution Lot # <u>2021-B</u> Acetone Stock Solution Lot # <u>2021-C</u>																																											
<table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td>MP5095</td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td>MP5096</td> <td>21070</td> <td>03/01/2023</td> </tr> <tr> <td>0.100</td> <td>MP5098</td> <td>21080</td> <td>03/08/2023</td> </tr> <tr> <td>0.200</td> <td>MP5100</td> <td>20510</td> <td>12/03/2022</td> </tr> <tr> <td>0.300</td> <td>MP5101</td> <td>21030</td> <td>02/02/2023</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>00521080A2</td> <td>02/05/2023</td> </tr> </tbody> </table>				Simulator	Serial #	Lot #	Expiration	0.000	MP5095	N/A	N/A	0.040	MP5096	21070	03/01/2023	0.100	MP5098	21080	03/08/2023	0.200	MP5100	20510	12/03/2022	0.300	MP5101	21030	02/02/2023	0.080 DGS	N/A	00521080A2	02/05/2023	<table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td>MP6284</td> </tr> <tr> <td>Interferent</td> <td>MP6285</td> </tr> <tr> <td>0.050</td> <td>MP6286</td> </tr> <tr> <td>0.080</td> <td>MP6287</td> </tr> <tr> <td>0.200</td> <td>MP6288</td> </tr> </tbody> </table>				Simulator	Serial Number	0.000	MP6284	Interferent	MP6285	0.050	MP6286	0.080	MP6287	0.200	MP6288
Simulator	Serial #	Lot #	Expiration																																												
0.000	MP5095	N/A	N/A																																												
0.040	MP5096	21070	03/01/2023																																												
0.100	MP5098	21080	03/08/2023																																												
0.200	MP5100	20510	12/03/2022																																												
0.300	MP5101	21030	02/02/2023																																												
0.080 DGS	N/A	00521080A2	02/05/2023																																												
Simulator	Serial Number																																														
0.000	MP6284																																														
Interferent	MP6285																																														
0.050	MP6286																																														
0.080	MP6287																																														
0.200	MP6288																																														
☒ Post Calibration Adjustment Stability Checks				Attachments ☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☒ Calibration Adjustment ☒ Post-Stability Checks ☐ Flow Calibration ☐ Form 40 ☒ Other <u>Extra Stabilities</u>																																											
Notes/Suggested Service: <u>During the 0.05 Stability Test, a value jumped out of range. Performed extra Stabilities and all passed. During the 0.20 Stability Test, another value jumped. Checked the source by flicking the inlet and checking the noise levels while running the pump. The noise levels were higher than suggested by CMI. Spoke to the AI. Instrument will be returned and either deregistered or kept out of service until sent to repair. (TDG 3/14/2022)</u>				☐ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☐ Return to/Place into Evidentiary Use ☒ Remain Out of Evidentiary Use ☐ Conduct an Agency Inspection Before Evidentiary Use																																											
David Eliezer Reyes-Rivera Digitally signed by David Eliezer Reyes-Rivera Date: 2022.03.23 14:44:37 -0400				 <u>29</u> <u>13:41:38</u> <u>04/00</u>																																											
Tech Review / Date				Admin Review / Date																																											

Type of Test	Serial Number	Agency	Date	Performed By
Stabilities	80-000808	Hillsborough CSO	03/14/2022	TDG <i>MC</i>

0.05g/210L 0.047 to 0.053	0.08g/210L 0.077 to 0.083	0.20g/210L 0.194 to 0.206	0.08g/210L 0.077 to 0.083	0.20g/210L 0.194 to 0.206	0.08g/210L 0.077 to 0.083	0.20g/210L 0.194 to 0.206
HILLSBOROUGH CO SO Intoxilyzer - Alcohol Analyzer Model: 8000 03/14/2022 Software: 8100.27 SN 80-000808 Test g/210L Time Air Blank: 0.000 13:55 Control Test 0.052 13:56 Air Blank: 0.000 13:57 Control Test 0.058 13:57 Air Blank: 0.000 13:58 Control Test 0.051 13:59 Air Blank: 0.000 13:59 Control Test Stats Average 0.0537 Std Dev 0.0038 Rel. Std Dev(%) 7.0545 <i>MC</i> Operator's Signature	HILLSBOROUGH CO SO Intoxilyzer - Alcohol Analyzer Model: 8000 03/14/2022 Software: 8100.27 SN 80-000808 Test g/210L Time Air Blank: 0.000 14:02 Control Test 0.082 14:03 Air Blank: 0.000 14:04 Control Test 0.081 14:04 Air Blank: 0.000 14:05 Control Test 0.081 14:06 Air Blank: 0.000 14:06 Control Test Stats Average 0.0813 Std Dev 0.0006 Rel. Std Dev(%) 0.7199 <i>MC</i> Operator's Signature	HILLSBOROUGH CO SO Intoxilyzer - Alcohol Analyzer Model: 8000 03/14/2022 Software: 8100.27 SN 80-000808 Test g/210L Time Air Blank: 0.000 14:45 Control Test 0.206 14:46 Air Blank: 0.000 14:46 Control Test 0.201 14:47 Air Blank: 0.000 14:48 Control Test 0.201 14:48 Air Blank: 0.000 14:49 Control Test Stats Average 0.2027 Std Dev 0.0029 Rel. Std Dev(%) 1.4244 <i>MC</i> Operator's Signature	HILLSBOROUGH CO SO Intoxilyzer - Alcohol Analyzer Model: 8000 03/14/2022 Software: 8100.27 SN 80-000808 Test g/210L Time Air Blank: 0.000 13:45 Control Test 0.081 13:45 Air Blank: 0.000 13:46 Control Test 0.081 13:46 Air Blank: 0.000 13:47 Control Test 0.081 13:47 Air Blank: 0.000 13:48 Control Test Stats Average 0.0810 Std Dev 0.0000 Rel. Std Dev(%) 0.0000 <i>MC</i> Operator's Signature	HILLSBOROUGH CO SO Intoxilyzer - Alcohol Analyzer Model: 8000 03/14/2022 Software: 8100.27 SN 80-000808 Test g/210L Time Air Blank: 0.000 13:45 Control Test 0.081 13:45 Air Blank: 0.000 13:46 Control Test 0.081 13:46 Air Blank: 0.000 13:47 Control Test 0.081 13:47 Air Blank: 0.000 13:48 Control Test Stats Average 0.0810 Std Dev 0.0000 Rel. Std Dev(%) 0.0000 <i>MC</i> Operator's Signature	HILLSBOROUGH CO SO Intoxilyzer - Alcohol Analyzer Model: 8000 03/14/2022 Software: 8100.27 SN 80-000808 Test g/210L Time Air Blank: 0.000 13:45 Control Test 0.081 13:45 Air Blank: 0.000 13:46 Control Test 0.081 13:46 Air Blank: 0.000 13:47 Control Test 0.081 13:47 Air Blank: 0.000 13:48 Control Test Stats Average 0.0810 Std Dev 0.0000 Rel. Std Dev(%) 0.0000 <i>MC</i> Operator's Signature	HILLSBOROUGH CO SO Intoxilyzer - Alcohol Analyzer Model: 8000 03/14/2022 Software: 8100.27 SN 80-000808 Test g/210L Time Air Blank: 0.000 13:45 Control Test 0.081 13:45 Air Blank: 0.000 13:46 Control Test 0.081 13:46 Air Blank: 0.000 13:47 Control Test 0.081 13:47 Air Blank: 0.000 13:48 Control Test Stats Average 0.0810 Std Dev 0.0000 Rel. Std Dev(%) 0.0000 <i>MC</i> Operator's Signature

Comments:

80-000808

3/14/2022

MG

HILLSBOROUGH CO, SO

Intoxilyzer - Alcohol Analyzer

Model 8000

SN 80-000808

03/14/2022

Software: 8100.27

Extra 0.05 Stabilities

Test	g/210L	Time
Air Blank	0.000	14:12
Control Test	0.052	14:12
Air Blank	0.000	14:13
Control Test	0.050	14:14
Air Blank	0.000	14:14
Control Test	0.051	14:15
Air Blank	0.000	14:16
Control Test	0.051	14:16
Air Blank	0.000	14:17
Control Test	0.051	14:18
Air Blank	0.000	14:18
Control Test	0.051	14:19
Air Blank	0.000	14:19
Control Test	0.050	14:20
Air Blank	0.000	14:21
Control Test	0.051	14:21
Air Blank	0.000	14:22
Control Test	0.051	14:23
Air Blank	0.000	14:23
Control Test	0.051	14:24
Air Blank	0.000	14:25
Control Test	0.050	14:25
Air Blank	0.000	14:26
Control Test	0.051	14:27
Air Blank	0.000	14:27
Control Test	0.050	14:28
Air Blank	0.000	14:28
Control Test	0.050	14:29
Air Blank	0.000	14:30
Control Test	0.050	14:31
Air Blank	0.000	14:31
Control Test	0.051	14:32
Air Blank	0.000	14:32
Control Test	0.051	14:33
Air Blank	0.000	14:34
Control Test	0.051	14:34
Air Blank	0.000	14:35
Control Test	0.052	14:36
Air Blank	0.000	14:36
Control Test	0.051	14:37
Air Blank	0.000	14:37
Control Test Stats		
Average	0.0508	
Std Dev	0.0006	
Rel Std Dev(%)	1.2118	

MG

Operator's Signature

```
***** AUTO CAL DATA *****
<<<< CHANNEL 1 >>>>
Sol Val = 0.000 % of 0.000
% Abs = 0.1105
Std Dev = 0.01 Rel Std Dev =
Sol Val = 0.195 % of 0.000
% Abs = 0.781
Std Dev = 0.01 Rel Std Dev =
Sol Val = 0.4762 % of 0.100
```

Sol Val = 0.1905 mg/l or 0.040
% Abs = 0.781
Std Dev = 0.01 Rel Std Dev =
Sol Val = 0.4762 mg/l or 0.100

Std Dev = 0.01 Rel Std Dev =
Sol. Val = 0.9524 mg/L or 0.201
% Abs = 3.489
Std Dev = 0.01 Rel Std Dev =
Sol. Val = 1.4296 mg/L or 0.330

```

% Ads = 5.103
Std Dev = 0.01 Rel Std Dev =
Zero Order Coef = -267.97
First Order Coef = 2742.85
Second Order Coef = 20.92
Standard Deviation = 24.838474
-----
<<<< CHANNEL 2 >>>>

```

Sample	Concentration (mg/L)	Std Dev	% RSD	Recovery (%)
1	0.000	0.000	0.179	100.0
2	0.040	0.005	1.250	100.0
3	0.100	0.015	1.557	100.0
4	0.200	0.032	1.600	100.0
5	0.400	0.072	1.800	100.0
6	0.800	0.144	1.800	100.0
7	1.600	0.288	1.800	100.0
8	3.200	0.576	1.800	100.0
9	6.400	1.152	1.800	100.0
10	12.800	2.304	1.800	100.0
11	25.600	4.608	1.800	100.0
12	51.200	9.216	1.800	100.0
13	102.400	18.432	1.800	100.0
14	204.800	36.864	1.800	100.0
15	409.600	73.728	1.800	100.0
16	819.200	147.456	1.800	100.0
17	1638.400	294.912	1.800	100.0
18	3276.800	589.824	1.800	100.0
19	6553.600	1179.648	1.800	100.0
20	13107.200	2359.296	1.800	100.0
21	26214.400	4718.592	1.800	100.0
22	52428.800	9437.184	1.800	100.0
23	104857.600	18874.368	1.800	100.0
24	209715.200	37748.736	1.800	100.0
25	419430.400	75497.472	1.800	100.0
26	838860.800	150994.944	1.800	100.0
27	1677721.600	301989.888	1.800	100.0
28	3355443.200	603979.776	1.800	100.0
29	6710886.400	1207959.552	1.800	100.0
30	13421772.800	2415919.104	1.800	100.0
31	26843545.600	4831838.208	1.800	100.0
32	53687091.200	9663676.416	1.800	100.0
33	107374182.400	19327352.832	1.800	100.0
34	214748364.800	38654705.664	1.800	100.0
35	429496729.600	77309411.328	1.800	100.0
36	858993459.200	154618822.656	1.800	100.0
37	1717986918.400	309237645.312	1.800	100.0
38	3435973836.800	618475290.624	1.800	100.0
39	6871947673.600	1236950581.248	1.800	100.0
40	13743895347.200	2473901162.496	1.800	100.0
41	27487790694.400	4947802324.992	1.800	100.0
42	54975581388.800	9895604649.984	1.800	100.0
43	109951162777.600	19791209299.968	1.800	100.0
44	219902325555.200	39582418599.936	1.800	100.0
45	439804651110.400	79164837199.872	1.800	100.0
46	879609302220.800	158329674399.744	1.800	100.0
47	1759218604441.600	316659348799.488	1.800	100.0
48	3518437208883.200	633318697598.976	1.800	100.0
49	7036874417766.400	1266637395197.952	1.800	100.0
50	14073748835532.800	2533274790395.904	1.800	100.0
51	28147497671065.600	5066549580791.808	1.800	100.0
52	56294995342131.200	10133099161583.616	1.800	100.0
53	112589990684262.400	20266198323167.232	1.800	100.0
54	225179981368524.800	40532396646334.464	1.800	100.0
55	450359962737049.600	81064793292668.928	1.800	100.0
56	900719925474099.200	162129586585337.856	1.800	100.0
57	1801439850948198.400	324259173170675.712		

```

      & A05 = 9.931
Std Dev = 0.01 Rel Std Dev =
Zero Order Coef = -239.24
First Order Coef = 1354.05
Second Order Coef = 10.93
Standard Deviation = 6.94811

```

```

First Order Coef = 235.24
Second Order Coef = 1354.05
Third Order Coef = 10.93
Standard Deviation = 6.94811

```

Solution State Quadratic Fit Coefficients			
Act	Fit	Residual	
g210L	g210L	g210L	
0.000	0.000	-0.0004	
0.040	0.040	0.0004	
0.100	0.100	0.0005	
0.200	0.201	-0.0007	

Note #1

Auto Calibration

Max Power Res Value = 93

Auto Range Res Value = 71

Sol Value = 0.000 g/200 L ***

Fit Value = 0.0000 mg/l % % % %

Note #2

Sample # 4 = 0.1050 (0.0110)

Avg % Abs = 0.1050 (-0.0073)

STD DEV = 0.0052 (0.0251)

TLC 03/17/2022

Type of Test	Serial Number	Agency	Date	Performed By
Stabilities (Post-Cal)	80-000808	Hillsborough CSO	03/17/2022	TDG ML

0.05g/210L	0.08g/210L	0.20g/210L	DGS 0.08g/210L
0.047 to 0.053	0.077 to 0.083	0.194 to 0.206	0.077 to 0.083
HILLSBOROUGH CO SO Intoxilyzer - Alcohol Analyzer Model: 8000 03/17/2022 Software: 8100.27 SN 80-000808 Test g/210L Time Air Blank 0.000 11:11 Control Test 0.051 11:11 Air Blank 0.000 11:12 Control Test 0.050 11:13 Air Blank 0.000 11:13 Control Test 0.049 11:14 Air Blank 0.000 11:15 Control Test Stats Average 0.0500 Std Dev 0.0010 Rel Std Dev(%) 2.0000 Operator's Signature ML	HILLSBOROUGH CO SO Intoxilyzer - Alcohol Analyzer Model: 8000 03/17/2022 Software: 8100.27 SN 80-000808 Test g/210L Time Air Blank 0.000 11:18 Control Test 0.081 11:19 Air Blank 0.000 11:19 Control Test 0.080 11:20 Air Blank 0.000 11:20 Control Test 0.079 11:21 Air Blank 0.000 11:22 Control Test Stats Average 0.0800 Std Dev 0.0010 Rel Std Dev(%) 1.2500 Operator's Signature ML	HILLSBOROUGH CO SO Intoxilyzer - Alcohol Analyzer Model: 8000 03/17/2022 Software: 8100.27 SN 80-000808 Test g/210L Time Air Blank 0.000 11:29 Control Test 0.201 11:30 Air Blank 0.000 11:30 Control Test 0.199 11:31 Air Blank 0.000 11:32 Control Test 0.199 11:32 Air Blank 0.000 11:33 Control Test Stats Average 0.1997 Std Dev 0.0012 Rel Std Dev(%) 0.5783 Operator's Signature ML	HILLSBOROUGH CO SO Intoxilyzer - Alcohol Analyzer Model: 8000 03/17/2022 Software: 8100.27 SN 80-000808 Test g/210L Time Air Blank 0.000 11:04 Control Test 0.079 11:05 Air Blank 0.000 11:05 Control Test 0.080 11:05 Air Blank 0.000 11:06 Control Test 0.080 11:06 Air Blank 0.000 11:07 Control Test Stats Average 0.0797 Std Dev 0.0006 Rel Std Dev(%) 0.7247 Operator's Signature ML

Comments:

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: HILLSBOROUGH CO SO
Time of Inspection: 14:51

Date of Inspection: 03/17/2022

Serial Number: 80-000808
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202201C Exp: 01/11/2024	0.08g/210L Test (g/210L) Lot#:202201D Exp: 01/18/2024	0.20g/210L Test (g/210L) Lot#:202201E Exp: 01/18/2024	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG115904 Exp: 06/08/2023
0.000	0.050	0.080	0.201	0.080
0.000	0.049	0.080	0.200	0.079
0.000	0.050	0.080	0.201	0.079
0.000	0.050	0.080	0.201	0.079
0.000	0.050	0.080	0.200	0.080
0.000	0.050	0.080	0.201	0.080
0.000	0.050	0.080	0.202	0.081
0.000	0.050	0.080	0.201	0.080
0.000	0.050	0.080	0.201	0.080
0.000	0.049	0.080	0.200	0.079


Standard Deviations	0.0004	0.0000	0.0006	0.0006
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0004 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.


I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



TAYLOR D GUTSCHOW

Signature and Printed Name

03/17/2022
Date



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
4700 Terminal Drive, Suite 1
Ft. Myers, FL 33907

This is to certify the calibration of Intoxilyzer 8000 serial number 80-000808, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-000808</u>	UNCERTAINTY* $\pm$
Owning Agency:	<u>HILLSBOROUGH CO SO</u>	0.050 g/ 210 L 0.004
Calibration Date:	<u>03/17/2022</u>	0.080 g/ 210 L 0.004
Calibration Time:	<u>14:51</u>	0.200 g/ 210 L 0.007
		0.080 g/ 210 L Dry Gas Control 0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards. This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.



03/17/2022

Date

TAYLOR D GUTSCHOW,

Department Inspector

FDLE/ATP Form 69 December 2021

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality