

#### **INSTRUMENT PROCESSING SHEET**

Agency Jacksonville Sheriff's Office S/N 80-001083

Florida Department of Law Enforcement

Date In 11/22/2022 DI Completion Date 11/29/2022 ■Ship □P/U □H/D □CMI □EE

Intake	By_PN	١	Quality C	heck	s By PN	Date_11/29/2022	Flow Calibi	ration By	Date
■ Annual			Breath	Tub	e Screen		Flow Colum	nn #	
☐ Registration	l		■ Replac	e Ext	ernal O-Rin	gs	□ 5L/r	min – 17mm	
Return from	CMI / EE		Instrun	nent	Set Up Veri	fied	☐ 15L,	/min – 53mm	
Visual Inspection	on:		R-Value	e <u>13</u>	32		□ 30L,	/min – 103mm	
Case	Handle		Flow V	erific	cation (L/s)		☐ R-Value		
	■ Case				# ATP 105			ibration Verifica	ation (L/s)
•	·		32 mm	.16	7	(.139169)		nn #	
Feet	■ Breath Tube		36 mm .175			(.156190)	32 mm		(.139169)
■ Ports	Screws Tigh	It	53 mm	.25	0	(.228278)			(.156190)
	ent/ Accessories:		103 mm	.51	5	(.447547)			(.228278)
	Printer Cab				Pressure Ch		103 mm		(.447547)
Static Bag	12V DC Cab	le	Gauge ID	# <u>28</u>	421				
Notes:			■ Stabilit	y Ch	ecks				
			Simulato	r	Serial #	Lot #/Exp	Maintenan	ice	Ву
			0.050			202201C	☐ Battery	Replacement	
					MP5088	01/11/2024		Regulator Repla	
			0.080			202201D		ube Replaceme	ent
					MP5089	01/18/2024	☐ Other _		
			0.200	1		202201E			
					MP5084	01/18/2024			
			0.080 DG	20	N/A	AG229803	·		
			0.080 DC	,,	IN/ A	10/25/2024	<del></del>		
Calibration Adj	iustmont			D	/PN	Department Inspec	tion		By PN
	essure Gauge <u>101</u>	18	ID#30		/	Barometric Pressur			<u> Бу </u>
	Serial #	Lot #		1	piration	Gauge 1017		trument 1016	
0.000	MP5091		N/A	LV	N/A	Mouth Alcohol Solu			
0.040	MP5082		1070	03/	01/2023	Acetone Stock Solu			
0.100				<b>—</b>	1	Simulator	1011 EOC# 20	Serial Number	
	MP5083		1380		13/2023	0.000			5086
0.200	MP6297		0510		03/2022	Interferent			5087
0.300	MP5085	2	1420	10/	20/2023	0.050		MP5088	
0.080 DGS	N/A	0812	21080A1	05/	05/2023	0.080			5089
Post Calibrat	tion Adjustment S	Stability	/ Checks			0.200		MP	5084
Simulator S	erial #	Lot#		Exp	iration	Attachments			
0.050	MP5088	20:	2201C	01/	11/2024	Form 41		Post-Stabil	ity Checks
0.080	MP5089		2201D		18/2024	■ Stability Check	5	☐ Flow Calibr	ation
0.200	MP5084	20:	2201E		18/2024	Calibration Cer	tificate	☐ Form 40	
0.080 DGS	N/A		229803		25/2024	Calibration Adj	ustment	☐ Other	
						- In at more and Ca		Sharatari 11D O	FAC
Notes/Suggest	ed Service:					Instrument Co			
						☐ Instrument Do			11D-8, FAC
					<del></del>	■ Return to/Plac ■ Remain Out of		•	
						Conduct an Ag	ency Inspecti	on Before Evid	entiary Use
					[	Taylor Digitally signed by Gutschow Date: 2022.11.30 -05'00'	3:02:08 Israe		rally signed by Israel Soto : 2022.12.01 15:18:11 -05'00'
						Tech Review / D	ate	Admin Review	w / Date

Type of Test	Serial Number	Agency	Date	Performed By
Stability Checks	80-001083	Jacksonville Sheriff's Office	11/29/2022	PN

0.05g/210L 0.08g/210L		0.20g/210L	DGS 0.08g/210L
0.047 to 0.053	0.077 to 0.083	0.194 to 0.206	0.077 to 0.083 ≤0.003 of Wet
ICKSONVILLE SO  Itoxilyzer - Alcoho! Analyzer  Idel 8000 SN 80-001083  /29/2022  If ware: 8100.27  Ist g/210L Time  In Blank 0.000 10:10  Intro! Test 0.047 10:11  In Blank 0.000 10:12  In Blank 0.000 10:13  Intro! Test 0.047 10:13  Intro! Test 0.048 10:13  In Blank 0.000 10:14  Intro! Test Stats  Verage 0.0473  It Dev 0.0006  It Std Dev(%) 1.2198	JACKSONUILLE SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001083 11/29/2022 Software: 8100.27  Test g/210L Time	JACKSONVILLE SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-601083 11/29/2022 Software: 8100.27  Test g/210L Time	JACKSONUILLE SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001083 11/29/2022 Software: 8100.27  Test g/210L Time
uperator's Signature	Operator's Signature	Operator's signature	Operator's Signature

# Florida Department of Law Enforcement Alcohol Testing Program

### DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: JACKSONVILLE SO Time of Inspection: 14:00

Date of Inspection: 11/29/2022

Serial Number: 80-001083

Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check			Date and/or Time Adjusted		
(Pre-Inspection): OK	Yes				No
Minimum Sample Volume			Barometric Pressure Sensor		
Check: OK	Yes		Check: OK	Yes	
Alcohol Free Subject			Mouth Alcohol Test:		
Test: 0.000	Yes		Slope Not Met	Yes	
Interferent Detect Test:			Diagnostic Check		
Interferent Detect	Yes		(Post-Inspection): OK	Yes	

Alcohol Free	0.05g/210L Test	0.08g/210L Test	0.20g/210L Test	0.08 g/210L
Test (g/210L)	(g/210L) Lot#:202201C Exp: 01/11/2024	(g/210L) Lot#:202201D Exp: 01/18/2024	(g/210L) Lot#:202201E Exp: 01/18/2024	Dry Gas Std Test (g/210L) Lot#:AG229803 Exp: 10/25/2024
0.000	0.049	0.078	0.199	0.080
0.000	0.049	0.078	0.199	0.079
0.000	0.049	0.078	0.199	0.079
0.000	0.049	0.078	0.199	0.080
0.000	0.049	0.078	0.199	0.080
0.000	0.048	0.078	0.199	0.080
0.000	0.049	0.078	0.200	0.080
0.000	0.049	0.078	0.199	0.080
0.000	0.049	0.078	0.199	0.080
0.000	0.049	0.078	0.199	0.080
Standard Deviations	0.0003	0.0000	0.0003	0.0004

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0002 Number of Simulators Used: 5

The above instrument complies ( X ) does not comply ( ) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

PHIL NICODEMO

Signature and Printed Name

11/29/2022 Date JACKSONVILLE SO Intoxilyzer - Alcohol Analyzer

Model 8000 SN 80-001083 11/29/2022 10:59:38 OPTICAL BENCH CACIBRATION
ADTUSTMENT
PLGS



Auto Calibration

pg 1 of 2

	<<<<	3um >>>>	<<<<	9um >>>>
Solution = 0 Sample #1 Sample #2 Sample #3 Sample #4 Avg % Abs STD DEV REL STD DEV	.000 g/210L % Abs 0.1710 0.1510 0.1510 0.1470 0.1497 0.0023 1.543	or 0.0000 mg/l, (% Abs Ref) (-0.0280) (0.0310) (0.0000) (0.0350) (0.0220) (0.0192) (87.078)	Samples = 4, % Abs 0.1360 0.1260 0.1110 0.1210 0.1193 0.0076 6.400	Discarded = 1 (% Abs Ref) (-0.0150) (0.0030) (-0.0010) (0.0080) (0.0033) (0.0045) (135.278)
Solution = 0 Sample Sample #1 Sample #2 Sample #3 Sample #4 Avg % Abs STD DEV REL STD DEV	.040 g/210L % Abs 0.8650 0.8980 0.8900 0.8740 0.8873 0.0122 1.377	or 0.1905 mg/l, (% Abs Ref) (-0.0050) (-0.0060) (0.0160) (0.0020) (0.0040) (0.0111) (278.388)	Samples = 4. % Abs 1.5130 1.5310 1.5280 1.5480 1.5357 0.0108 0.702	Discarded = 1 (% Abs Ref) (-0.0100) (-0.0200) (-0.0020) (-0.0280) (-0.0167) (0.0133) (79.900)
Solution = 0. Sample Sample #1 Sample #2 Sample #3 Sample #4 Avg % Abs STD DEV REL STD DEV	% Abs 2.0120 1.9770 2.0030 2.0010 1.9937 0.0145 0.726	or 0.4762 mg/l, (% Abs Ref) (-0.0190) (0.0030) (0.0090) (-0.0130) (-0.0003) (0.0114) (3411.742)	Samples = 4, % Abs 3.5960 3.5940 3.6110 3.5900 3.5983 0.0112 0.310	
Solution = 0. Sample Sample #1 Sample #2 Sample #3 Sample #4 Avg % Abs STD DEV REL STD DEV	200 g/210L % Abs 3.8400 3.8030 3.8220 3.8070 3.8107 0.0100 0.263	or 0.9524 mg/l, (% Abs Ref) (-0.0140) (0.0310) (0.0200) (0.0110) (0.0207) (0.0100) (48.468)	% Abs 6.8720 6.8790 6.8570 6.8480 6.8613 0.0159 0.232	(% Abs Ref) (-0.0140) (0.0020) (0.0010) (-0.0010) (0.0007) (0.0015) (229.129)
Sample #1 Sample #2 Sample #3	300 g/210L % Abs	or 1.4286 mg/l, (% Abs Ref) (-0.0110) (-0.0240) (-0.0200) (0.0060) (-0.0127) (0.0163) (128.598)	% Abs 9.9810 9.9730 9.9720 9.9750 9.9733 0.0015	Discarded = 1 (% Abs Ref) (-0.0100) (0.0000) (-0.0120) (0.0140)

JACKSONVILLE SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001083
11/29/2022 10:59:38



Auto Calibration

pg 2 of 2

	<<<< 3u	.m >>>>	<<<<	9um	>>>>
Zero Order Co First Order of Second Order	Coef 2501.5	•		0.41 5.57 4	
Act (g/210L) 0.000 0.040 0.100 0.200 0.300	Fit (g/210L) 0.000 0.040 0.099 0.201 0.300	Residual (g/210L) -0.0004 0.0005 -0.0007 0.0003	Act (g/210L) 0.000 0.040 0.100 0.200 0.300	Fit (g/210L -0.000 0.040 0.100 0.200 0.300	

	<<<<	3um	>>>>	<<<<	9um	>>>>
Sample	).080 g/210L	or 0.3	810 mg/l,	Samples = $4$ ,	Discard	led = 1
Sample #1		2858	.00		323	8.00
Sample #2		2908	.00		328	3.00
Sample #3		2795	.00		323	2.00
Sample #4		2920	.00		329	5.00
Avg		2874	.3333		327	0.0000
STD DEV		68.9	662		33.	4515
REL STD DEV		2.39	9		1.0	123
H2O adjust (	mg/l*10k)	935			539	)

Barometric Pressure = 1017

Type of Test	Serial Number	Agency	Date	Performed By
Post Calibration Adjustment Stability Checks	80-001083	Jacksonville Sheriff's Office	11/29/2022	PN

0.05g/210L		0.08g/210L		0.20g/210L		DGS 0.0	8g/210L
0.047 to 0.053	X	0.077 to 0.083	X	0.194 to 0.206	X	0.077 to 0.083 X	
ACKSONUILLE SO ntoxilyzer - Alcohol Analyzer odel 8000 SN 80-0 1/29/2022 oftware: 8100.27 est g/210L	Time  12:12	JACKSONUILLE SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001083 11/29/2022 Software: 8100.27 Test g/210L Ti	me ,	JACKSONUILLE SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001083 11/29/2022 Software: 8100.27 Test g/210L Time		JACKSONUILLE SO Intoxilyzer - Alcohol Model 8000 11/29/2022 Software: 8100.27 Test 9/2	SN 80-001083
Control Test 0.049 Air Blank 0.000 Control Test 0.048 Air Blank 0.000 Control Test 0.049 Air Blank 0.000 Control Test 0.049 Air Blank 0.000 Control Test Stats Average 0.0487 Std Dev 0.0006 Rel Std Dev(%) 1.1863	12:12 12:13 12:14 12:14 12:15 12:15	Control Test 0.077 12: Air Blank 0.000 12: Control Test 0.077 12: Air Blank 0.000 12: Control Test 0.078 12: Control Test 0.078 12: Control Test Stats Average 0.0773 13: Std Dev 0.0006 Rel Std Dev(%) 0.7466	17 18 19	Control Test 0.000 12:22  Air Blank 0.000 12:22  Control Test 0.199 • 12:23  Air Blank 0.000 12:24  Control Test 0.198 • 12:24  Air Blank 0.000 12:24  Air Blank 0.000 12:25  Control Test Stats  Auerage 0.1987  Std Deu 0.0006  Rel Std Deu(%) 0.2906		Control Test 0.00 Air Blank 0.01 Control Test 0.00 Air Blank 0.0 Control Test 0.00 Air Blank 0.0 Control Test 5.00	12:26 12:27 79 12:27 100 12:27 12:27 79 12:28 12:28 12:28
Plittuse) Operator's Signature		Parturds Operator's Signature	-	<b>Phol Mos</b> Operator's Signature		Phil W Operator's Si	gnature



## **Calibration Certificate**

Florida Department of Law Enforcement Alcohol Testing Program 2331 Phillips Road. Suite B1032 Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-001083, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	80-001083	UNCERTAINTY* ±	
Owning Agency:	JACKSONVILLE SO	0.050 g/ 210 L	0.004
Calibration Date:	11/29/2022	0.080 g/ 210 L	0.004
Calibration Time:	14:00	0.200 g/ 210 L	0.007
	0	0.080 g/210 L Dry Gas Control	0.005

All results are reported in g/210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration. \*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence (k=3).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

#### TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/IEC 17025 standards. Digitally signed by Phil

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

11/29/2022

Phil Nicodemo Date: 2022.11.29 14:23:58

-05'00'

Date

PHIL NICODEMO, Department Inspector

FDLE/ATP Form 69 March 2022

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality

Page 1 of 1