

INSTRUMENT PROCESSING SHEET

Agency Florida Department of Law Enforcement S/N 80-000227 DI Completion Date 02/06/2024 □Ship □P/U □H/D □CMI □EE Florida Department of Date In <u>n</u>/a Law Enforcement Date 02/06/2024 Quality Checks By TDG Date 02/06/2024 Intake By TDG Flow Calibration By Date Annual ■ Breath Tube Screen Flow Column # ☐ 5L/min – 17mm ■ Registration ■ Replace External O-Rings ☐ Return from CMI / EE ■ Instrument Set Up Verified ☐ 15L/min – 53mm ■ R-Value 225 □ 30L/min – 103mm Visual Inspection: ■ Flow Verification (L/s) □ R-Value Case ■ Handle Flow Column # ATP106 ☐ Post Calibration Verification (L/s) Dry Gas Shelf Keyboard Flow Column #_____ 32 mm 0.156 (.139 - .169)■ Feet ■ Breath Tube 36 mm 0.171 32 mm _____(.139 - .169) (.156 - .190)Ports Screws Tight 53 mm 0.242 (.228 - .278) 36 mm _____(.156 - .190) Other Equipment/ Accessories: 103 mm 0.511 (.447 - .547) 53 mm _____ (.228 - .278) Power cord ☐ Printer Cable ■ Barometric Pressure Check 103 mm _____ (.447 - .547) Static Bag ☐ 12V DC Cable Gauge ID # 26932 Stability Checks Notes: ____ Simulator Serial # Lot #/Exp Maintenance By____ _Date_ ☐ Battery Replacement 0.050 202303K MP5094 ☐ Dry Gas Regulator Replacement 03/29/2025 ☐ Breath Tube Replacement 0.080 202303L MP5095 ☐ Other ____ 03/29/2025 0.200 202304C MP5096 04/05/2025 0.080 DGS N/A 01923080A3 02/05/2025 By TDG **Calibration Adjustment** Department Inspection By Barometric Pressure ID# 26932 Barometric Pressure Gauge ID# Expiration ` Gauge 1009 Instrument 1009 Simulator | Serial # Lot# 0.000 N/A N/A Mouth Alcohol Solution Lot # 2023-A 0.040 Acetone Stock Solution Lot # 2022-B 0.100 Simulator Serial Number 0.000 MP5092 0.200 Interferent MP5093 0.300 0.050 MP5094 0.080 DGS N/A 0.080 MP5095 0.200 MP5096 ☐ Post Calibration Adjustment Stability Checks Attachments Simulator | Serial # Lot# Expiration 0.050 Form 41 ☐ Post-Stability Checks Stability Checks 0.080 ☐ Flow Calibration 0.200 Calibration Certificate Form 40 ☐ Calibration Adjustment ☐ Other __ 0.080 DGS N/A Notes/Suggested Service: Instrument does not print to the Instrument Complies with Chapter 11D-8, FAC internal printer. Repeated the DGS stability test to print ☐ Instrument Does Not Comply with Chapter 11D-8, FAC the results from external printer. All subsequent stability Return to/Place into Evidentiary Use checks printed from external printer. (TDG) ☐ Remain Out of Evidentiary Use Conduct an Agency Inspection Before Evidentiary Use Instrument will be retained at FDLE until requested by Shayla Digitally signed by Shayla Platt an agency. (TDG) Date: 2024.03.06 Phil Nicodemo Date: 2024.03.15 08:22:21-04:00' Tech Review / Date Admin Review / Date

Florida Department of Law Enforcement **Alcohol Testing Program**

AGENCY INSPECTION REPORT - INTOXILYZER 8000

Agency: FDLE

Time of Inspection: 09:58

Date of Inspection: 02/06/2024

Serial Number: 80-000227

Software: 8100.27

Check or Test	YES	NO
Date and/or Time Adjusted		No
Diagnostic Check (Pre-Inspection): OK		No
Alcohol Free Subject Test: 0.000		No
Mouth Alcohol Test: Slope Not Met		No
Interferent Detect Test: Interferent Detect	2	No
Diagnostic Check (Post-Inspection): OK	1 X ₂ 8	No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#: Exp:	0.08g/210L Test (g/210L) Lot#: Exp:	0.20g/210L Test (g/210L) Lot#: Exp:	0.08 g/210L Dry Gas Std Test (g/210L) Lot#: Exp:
		" 1 .		
S		× × × × × × × × × × × × × × × × × × ×		
		53		

	· · · · · · · · · · · · · · · · · · ·	
9	X	
Number of Simulators Used:	· N	

Remarks: AI NOT CONDUCTED. BYPASSED TO BRING OUT OF DISABLED MODE.

Not determined 02/06/2024 The above instrument complies (X) does not comply () with Chapter 11D-8, FAC.

I certify that I hold a valid Florida Department of Law Enforcement Agency Inspector Permit and that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Signature and Printed Name

02/06/2024 Date

Intoxilyzer - Alcohol Analyzer

Model 8000

SN 80-000227

02/06/2024

Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:59
Control Test	0.051	11:00
Air Blank	0.000	11:00
Control Test	0.051	11:01
Air Blank	0.000	11:01
Control Test	0.050	11:02
Air Blank	0.000	11:03
Control Test Stats	5	
Average	0.0507	
Std Dev	0.0006	
Rel Std Dev(%)	1.1395	

Mb

Operator's Signature

Intoxilyzer - Alcohol Analyzer

Model 8000

SN 80-000227

02/06/2024

Software: 8100.27

Net

Test	g/210L	Time
Air Blank	0.000	11:06
Control Test	0.080	11:07
Air Blank	0.000	11:07
Control Test	0.080	11:08
Air Blank	0.000	11:08
Control Test	0.079	11:09
Air Blank	0.000	11:09
Control Test Stat	cs	
Average	0.0797	
Std Dev	0.0006	
Rel Std Dev(%)	0.7247	

Operator's Signature

3

ř.

Intoxilyzer - Alcohol Analyzer

Model 8000

SN 80-000227

02/06/2024

Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:22
Control Test	0.198	11:22
Air Blank	0.000	11:23
Control Test	0.198	11:23
Air Blank	0.000	11:24
Control Test	0.198	11:25
Air Blank	0.000	11:25
Control Test Stat	cs	
Average	0.1980	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

MG

Operator's Signature

Intoxilyzer - Alcohol Analyzer

Model 8000

SN 80-000227

02/06/2024

Software: 8100.27

065

Test	g/210L	Time
Air Blank	0.000	10:56
Control Test	0.081	10:56
Air Blank	0.000	10:56
Control Test	0.081	10:57
Air Blank	0.000	10:57
Control Test	0.080	10:58
Air Blank	0.000	10:58
Control Test Stats		
Average	0.0807	
Std Dev	0.0006	
Rel Std Dev(%)	0.7157	

Operator's Signature

Florida Department of Law Enforcement **Alcohol Testing Program**

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: FDLE

Alcohol Free

Serial Number: 80-000227

0.08 g/210L

0.20g/210L Test

Time of Inspection: 13:53

0.05g/210L Test

Date of Inspection: 02/06/2024

Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	1
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

0.08g/210L Test

Test (g/210L)	(g/210L) Lot#:202303K Exp: 03/29/2025	(g/210L) Lot#:202303L Exp: 03/29/2025	(g/210L) Lot#:202304C Exp: 04/05/2025	Dry Gas Std Test* (g/210L) Lot#:01923080A3 Exp: 02/05/2025
0.000	0.050	0.080	0.196	0.081
0.000	0.050	0.080	0.198	0.080
0.000	0.050	0.080	0.198	0.080
0.000	0.050	0.080	0.197	0.080
0.000	0.051	0.080	0.197	0.080
0.000	0.050	0.081	0.198	0.080
0.000	0.050	0.080	0.197	0.080
0.000	0.050	0.080	0.197	0.080
0.000	0.050	0.081	0.197	0.080
0.000	0.050	0.080	0.197	0.079
Standard Deviations	0 0003	0 0004	0.0006	0.0004

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0004 Number of Simulators Used: 5 Remarks:

The above instrument complies (X) does not comply () with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Signature and Printed Name

02/06/2024 Date



Calibration Certificate

Florida Department of Law Enforcement Alcohol Testing Program 4700 Terminal Drive, Suite 1 Ft. Myers, FL 33907

This is to certify the calibration of Intoxilyzer 8000 serial number 80-000227, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

UNCERTAINTY* ±	0.050 g/210 L	0.080 g/210 L	- 0.200 g/210 L
80-000227	FDLE	02/06/2024	<u>13:53</u>
Serial Number:	Owning Agency:	on Date:	Calibration Time:

All results are reported in g/210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration. *Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence (k=3).

0.005

0.080 g/210 L Dry Gas Control

0.004

0.004

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

FRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards. Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards. This document shall not be reproduced except in full,

without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

02/06/2024

FAYLOR D GUTSCHOW Department Inspector

Service • Integrity • Respect • Quality

Issuing Authority: Alcohol Testing Program

FDLE/ATP Form 69 December 2021

Page 1 of