



# INSTRUMENT PROCESSING SHEET

Agency Marianna PDS/N 80-001305Florida Department of  
Law EnforcementDate In 09/05/2024DI Completion Date 09/18/2024 Ship  P/U  H/D  CMI  EE

<b>Intake</b> By <u>ALL *</u> Date <u>09/05/2024 *</u> <input checked="" type="checkbox"/> Annual *      * DA <u>9/24/24</u> <input type="checkbox"/> Registration <input type="checkbox"/> Return from CMI / EE  Visual Inspection: <input checked="" type="checkbox"/> Case* <input checked="" type="checkbox"/> Handle * <input checked="" type="checkbox"/> Keyboard* <input checked="" type="checkbox"/> Dry Gas Shelf * <input checked="" type="checkbox"/> Feet * <input checked="" type="checkbox"/> Breath Tube * <input checked="" type="checkbox"/> Ports * <input checked="" type="checkbox"/> Screws Tight *  Other Equipment/ Accessories: <input type="checkbox"/> Power cord <input type="checkbox"/> Printer Cable <input type="checkbox"/> Static Bag <input type="checkbox"/> 12V DC Cable  Notes: <u>dropoff with DGS</u>       	<b>Quality Checks</b> By <u>DA</u> Date <u>9/18/2024</u> <input checked="" type="checkbox"/> Breath Tube Screen <input checked="" type="checkbox"/> Replace External O-Rings <input checked="" type="checkbox"/> Instrument Set Up Verified <input checked="" type="checkbox"/> R-Value <u>213</u> <input checked="" type="checkbox"/> Flow Verification (L/s) Flow Column # <u>ATP102</u> 32 mm <u>0.152</u> (.139 - .169) 36 mm <u>0.171</u> (.156 - .190) 53 mm <u>0.242</u> (.228 - .278) 103 mm <u>0.511</u> (.447 - .547) <input checked="" type="checkbox"/> Barometric Pressure Check Gauge ID # <u>28663</u> <input checked="" type="checkbox"/> Stability Checks <table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td rowspan="2">0.050</td> <td rowspan="2">MP5088</td> <td>202303K</td> </tr> <tr> <td>03/29/2025</td> </tr> <tr> <td rowspan="2">0.080</td> <td rowspan="2">MP5089</td> <td>202303L</td> </tr> <tr> <td>03/29/2025</td> </tr> <tr> <td rowspan="2">0.200</td> <td rowspan="2">MP5090</td> <td>202304C</td> </tr> <tr> <td>04/05/2025</td> </tr> <tr> <td rowspan="2">0.080 DGS</td> <td rowspan="2">N/A</td> <td>06723080A5</td> </tr> <tr> <td>04/05/2025</td> </tr> </tbody> </table>	Simulator	Serial #	Lot #/Exp	0.050	MP5088	202303K	03/29/2025	0.080	MP5089	202303L	03/29/2025	0.200	MP5090	202304C	04/05/2025	0.080 DGS	N/A	06723080A5	04/05/2025	<b>Flow Calibration</b> By _____ Date _____ Flow Column # _____ <input type="checkbox"/> 5L/min – 17mm <input type="checkbox"/> 15L/min – 53mm <input type="checkbox"/> 30L/min – 103mm <input type="checkbox"/> R-Value _____ <input type="checkbox"/> Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)  <b>Maintenance</b> By _____ Date _____ <input type="checkbox"/> Battery Replacement <input type="checkbox"/> Dry Gas Regulator Replacement <input type="checkbox"/> Breath Tube Replacement <input type="checkbox"/> Other _____    
Simulator	Serial #	Lot #/Exp																			
0.050	MP5088	202303K																			
		03/29/2025																			
0.080	MP5089	202303L																			
		03/29/2025																			
0.200	MP5090	202304C																			
		04/05/2025																			
0.080 DGS	N/A	06723080A5																			
		04/05/2025																			

<b>Calibration Adjustment</b> By _____ Barometric Pressure Gauge _____ ID # _____ <table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td></td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.100</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.300</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table> <input type="checkbox"/> Post Calibration Adjustment Stability Checks <table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table>	Simulator	Serial #	Lot #	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			Simulator	Serial #	Lot #	Expiration	0.050				0.080				0.200				0.080 DGS	N/A			<b>Department Inspection</b> By <u>DA</u> Barometric Pressure ID# <u>28663</u> Gauge <u>1011</u> Instrument <u>1010</u> Mouth Alcohol Solution Lot # <u>2024-A</u> Acetone Stock Solution Lot # <u>2023-B</u> <table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td>MP5086</td> </tr> <tr> <td>Interferent</td> <td>MP5087</td> </tr> <tr> <td>0.050</td> <td>MP5088</td> </tr> <tr> <td>0.080</td> <td>MP5089</td> </tr> <tr> <td>0.200</td> <td>MP5090</td> </tr> </tbody> </table> <b>Attachments</b> <input checked="" type="checkbox"/> Form 41 <input checked="" type="checkbox"/> Stability Checks <input checked="" type="checkbox"/> Calibration Certificate <input type="checkbox"/> Calibration Adjustment <input type="checkbox"/> Post-Stability Checks <input type="checkbox"/> Flow Calibration <input type="checkbox"/> Form 40 <input type="checkbox"/> Other _____  <input checked="" type="checkbox"/> Instrument Complies with Chapter 11D-8, FAC <input type="checkbox"/> Instrument Does Not Comply with Chapter 11D-8, FAC <input checked="" type="checkbox"/> Return to/Place into Evidentiary Use <input type="checkbox"/> Remain Out of Evidentiary Use <input checked="" type="checkbox"/> Conduct an Agency Inspection Before Evidentiary Use 	Simulator	Serial Number	0.000	MP5086	Interferent	MP5087	0.050	MP5088	0.080	MP5089	0.200	MP5090
Simulator	Serial #	Lot #	Expiration																																																										
0.000		N/A	N/A																																																										
0.040																																																													
0.100																																																													
0.200																																																													
0.300																																																													
0.080 DGS	N/A																																																												
Simulator	Serial #	Lot #	Expiration																																																										
0.050																																																													
0.080																																																													
0.200																																																													
0.080 DGS	N/A																																																												
Simulator	Serial Number																																																												
0.000	MP5086																																																												
Interferent	MP5087																																																												
0.050	MP5088																																																												
0.080	MP5089																																																												
0.200	MP5090																																																												

Notes/Suggested Service: <u>Repeated minimum sample volume check because alcohol was detected. DA 9/18/2024</u>  * Tech review: Intake was performed on 09/05/2024 by ALL. Transcribed into file from hard copy by DA on 9/24/24.    	<table border="1" style="width:100%; border-collapse: collapse;"> <tr> <td style="width:50%;">                 Phil Nicodemo  <small>Digitally signed by Phil Nicodemo Date: 2024.09.24 13:00:35 -04'00'</small> </td> <td style="width:50%; text-align: right;">                 Shayla Platt                  Date: 2024.09.24                  14:08:33 -04'00'             </td> </tr> <tr> <td style="text-align: center;">Tech Review / Date</td> <td style="text-align: center;">Admin Review / Date</td> </tr> </table>	Phil Nicodemo <small>Digitally signed by Phil Nicodemo Date: 2024.09.24 13:00:35 -04'00'</small>	Shayla Platt Date: 2024.09.24 14:08:33 -04'00'	Tech Review / Date	Admin Review / Date
Phil Nicodemo <small>Digitally signed by Phil Nicodemo Date: 2024.09.24 13:00:35 -04'00'</small>	Shayla Platt Date: 2024.09.24 14:08:33 -04'00'				
Tech Review / Date	Admin Review / Date				

# Florida Department of Law Enforcement Alcohol Testing Program

## DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: MARIANNA PD

Serial Number: 80-001305

Time of Inspection: 15:27

Date of Inspection: 09/18/2024

Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#: 202303K Exp: 03/29/2025	0.08g/210L Test (g/210L) Lot#: 202303L Exp: 03/29/2025	0.20g/210L Test (g/210L) Lot#: 202304C Exp: 04/05/2025	0.08 g/210L Dry Gas Std Test (g/210L) Lot#: 06723080A5 Exp: 04/05/2025
0.000	0.048	0.077	0.195	0.080
0.000	0.048	0.078	0.196	0.081
0.000	0.048	0.077	0.196	0.080
0.000	0.048	0.077	0.195	0.080
0.000	0.048	0.078	0.195	0.080
0.000	0.048	0.078	0.195	0.080
0.000	0.048	0.078	0.195	0.080
0.000	0.049	0.078	0.195	0.081
0.000	0.049	0.078	0.195	0.080
0.000	0.048	0.078	0.196	0.080


Standard Deviations	0.0004	0.0004	0.0004	0.0004
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0004      Number of Simulators Used: 5

**Remarks:**  
 REPEATED MIN SAMPLE VOLUME CHECK ALCOHOL DETECTED  
 \*Residual mouth alcohol remained from previous instrument. Allowed mouth alcohol to dissipate and repeated minimum sample volume check. DA 9/23/24

The above instrument complies (  ) does not comply (      ) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



DESTINEE N ARMSTRONG

Signature and Printed Name


09/18/2024

Date

Stability Checks 80-001305 9/18/24 DA

MARIANNA PD  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-001305  
09/18/2024  
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:48
Control Test	0.049	12:48
Air Blank	0.000	12:49
Control Test	0.049	12:49
Air Blank	0.000	12:50
Control Test	0.048	12:51
Air Blank	0.000	12:51
Control Test Stats		
Average	0.0487	
Std Dev	0.0006	
Rel Std Dev(%)	1.1863	




Operator's Signature

wet

MARIANNA PD  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-001305  
09/18/2024  
Software: 8100.27


Test	g/210L	Time
-----		
Air Blank	0.000	12:53
Control Test	0.078	12:54
Air Blank	0.000	12:54
Control Test	0.078	12:55
Air Blank	0.000	12:55
Control Test	0.078	12:56
Air Blank	0.000	12:57
Control Test Stats		
Average	0.0780	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	



Operator's Signature

MARIANNA PD  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-001305  
09/18/2024  
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:59
Control Test	0.197	12:59
Air Blank	0.000	13:00
Control Test	0.197	13:01
Air Blank	0.000	13:01
Control Test	0.196	13:02
Air Blank	0.000	13:02
Control Test Stats		
Average	0.1967	
Std Dev	0.0006	
Rel Std Dev(%)	0.2936	




---

Operator's Signature

DGS

MARIANNA PD  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-001305  
09/18/2024  
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:06
Control Test	0.081	13:06
Air Blank	0.000	13:07
Control Test	0.081	13:07
Air Blank	0.000	13:07
Control Test	0.081	13:08
Air Blank	0.000	13:08
Control Test Stats		
Average	0.0810	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	



---

Operator's Signature



# Calibration Certificate

Florida Department of Law Enforcement  
Alcohol Testing Program  
2331 Phillips Road.  
Suite B1032  
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-001305, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-001305</u>	UNCERTAINTY* $\pm$	
Owning Agency:	<u>MARIANNA PD</u>	0.050 g/ 210 L	0.004
Calibration Date:	<u>09/18/2024</u>	0.080 g/ 210 L	0.004
Calibration Time:	<u>15:27</u>	0.200 g/ 210 L	0.007
		0.080 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within  $\pm 0.005$  or 5%, whichever is greater, of the target alcohol concentration.

\*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence (k=3).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

### TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

Destinee  
Armstrong  
Digitally signed by Destinee  
Armstrong  
Date: 2024.09.23 11:59:17  
-04'00'

09/18/2024

Date

DESTINEE N ARMSTRONG,  
Department Inspector

FDLE/ATP Form 69 March 2022

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality