



## Agency Madison County SO

S/N 80-001307

Date In 05/14/2024      DI Completion Date 6/19/24

☐ Ship   ☒ P/U   ☐ H/D   ☐ CMI   ☐ EE

## Calibration Adjustment

By SP

## Barometric Pressure Gauge 1010 ID # 28663

### ■ Post Calibration Adjustment Stability Checks

Notes/Suggested Service:

Performed cal adjustment due to barometric pressure greater than 1% difference. Second DI completed with DGS lot different from cal adjustment. SP

## Department Inspection

By SP

Barometric Pressure ID# 28662

|                 |                             |
|-----------------|-----------------------------|
| Gauge 1010/1016 | <u>Instrument 1011/1017</u> |
|-----------------|-----------------------------|

Mouth Alcohol Solution Lot # 2024-A

Acetone Stock Solution Lot # 2023-B

## Attachments

- |                                                                                                                                                                                                                                  |                                                                                                                                                                                         |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input checked="" type="checkbox"/> Form 41<br><input checked="" type="checkbox"/> Stability Checks<br><input checked="" type="checkbox"/> Calibration Certificate<br><input checked="" type="checkbox"/> Calibration Adjustment | <input checked="" type="checkbox"/> Post-Stability Checks<br><input type="checkbox"/> Flow Calibration<br><input checked="" type="checkbox"/> Form 40<br><input type="checkbox"/> Other |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

- ☒ Instrument Complies with Chapter 11D-8, FAC
  - ☐ Instrument Does Not Comply with Chapter 11D-8, FAC
  - ☒ Return to/Place into Evidentiary Use
  - ☐ Remain Out of Evidentiary Use
  - ☒ Conduct an Agency Inspection Before Evidentiary Use

Taylor  
Gutschow

Digitally signed by Taylor Gutschow  
Date: 2024.06.20 11:56:55  
+04'00'

Phil Nicodemo

Digitally signed by Phil  
Nicodemo  
Date: 2024.06.27 13:26:52 -04'00'


Tech Review / Date

Admin Review / Date

Stability Checks 80-001307 5/28/24 DA

MADISON COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-001307  
05/28/2024  
Software: 8100.27

| Test               | g/210L | Time  |
|--------------------|--------|-------|
| Air Blank          | 0.000  | 11:12 |
| Control Test       | 0.049  | 11:12 |
| Air Blank          | 0.000  | 11:13 |
| Control Test       | 0.049  | 11:14 |
| Air Blank          | 0.000  | 11:14 |
| Control Test       | 0.049  | 11:15 |
| Air Blank          | 0.000  | 11:15 |
| Control Test Stats |        |       |
| Average            | 0.0490 |       |
| Std Dev            | 0.0000 |       |
| Rel Std Dev(%)     | 0.0000 |       |



Operator's Signature

MADISON COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-001307  
05/28/2024  
Software: 8100.27

| Test               | g/210L | Time  |
|--------------------|--------|-------|
| Air Blank          | 0.000  | 11:17 |
| Control Test       | 0.079  | 11:18 |
| Air Blank          | 0.000  | 11:18 |
| Control Test       | 0.078  | 11:19 |
| Air Blank          | 0.000  | 11:19 |
| Control Test       | 0.078  | 11:20 |
| Air Blank          | 0.000  | 11:20 |
| Control Test Stats |        |       |
| Average            | 0.0783 |       |
| Std Dev            | 0.0006 |       |
| Rel Std Dev(%)     | 0.7370 |       |

  
Operator's Signature

MADISON COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-001307  
05/28/2024  
Software: 8100.27

| Test               | g/210L | Time  |
|--------------------|--------|-------|
| Air Blank          | 0.000  | 11:22 |
| Control Test       | 0.197  | 11:23 |
| Air Blank          | 0.000  | 11:24 |
| Control Test       | 0.197  | 11:24 |
| Air Blank          | 0.000  | 11:25 |
| Control Test       | 0.196  | 11:25 |
| Air Blank          | 0.000  | 11:26 |
| Control Test Stats |        |       |
| Average            | 0.1967 |       |
| Std Dev            | 0.0006 |       |
| Rel Std Dev(%)     | 0.2936 |       |



Operator's Signature

MADISON COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-001307  
05/28/2024  
Software: 8100.27

| Test               | g/210L | Time  |
|--------------------|--------|-------|
| Air Blank          | 0.000  | 11:28 |
| Control Test       | 0.080  | 11:28 |
| Air Blank          | 0.000  | 11:29 |
| Control Test       | 0.080  | 11:29 |
| Air Blank          | 0.000  | 11:30 |
| Control Test       | 0.080  | 11:30 |
| Air Blank          | 0.000  | 11:30 |
| Control Test Stats |        |       |
| Average            | 0.0800 |       |
| Std Dev            | 0.0000 |       |
| Rel Std Dev(%)     | 0.0000 |       |

DGS



Operator's Signature

# Florida Department of Law Enforcement Alcohol Testing Program

## AGENCY INSPECTION REPORT - INTOXILYZER 8000

Agency: MADISON COUNTY SO  
Time of Inspection: 11:09

Date of Inspection: 05/28/2024

Serial Number: 80-001307  
Software: 8100.27

| Check or Test                               | YES | NO |
|---------------------------------------------|-----|----|
| Date and/or Time Adjusted                   |     | No |
| Diagnostic Check (Pre-Inspection): OK       |     | No |
| Alcohol Free Subject Test: 0.000            |     | No |
| Mouth Alcohol Test: Slope Not Met           |     | No |
| Interferent Detect Test: Interferent Detect |     | No |
| Diagnostic Check (Post-Inspection): OK      |     | No |

|                            |                                           |                                           |                                           |                                                        |
|----------------------------|-------------------------------------------|-------------------------------------------|-------------------------------------------|--------------------------------------------------------|
| Alcohol Free Test (g/210L) | 0.05g/210L Test (g/210L)<br>Lot#:<br>Exp: | 0.08g/210L Test (g/210L)<br>Lot#:<br>Exp: | 0.20g/210L Test (g/210L)<br>Lot#:<br>Exp: | 0.08 g/210L Dry Gas Std Test (g/210L)<br>Lot#:<br>Exp: |
|                            |                                           |                                           |                                           |                                                        |
|                            |                                           |                                           |                                           |                                                        |
|                            |                                           |                                           |                                           |                                                        |

Number of Simulators Used: \_\_\_\_\_


Remarks:

BYPASS AI TO OPERATE INSTRUMENT COMPLIANCE NOT DETERMINED

*Compliance not determined DA 5/28/24*

~~The above instrument complies ( ☒ ) does not comply ( ☐ ) with Chapter 11D-8, FAC.~~

I certify that I hold a valid Florida Department of Law Enforcement Agency Inspector Permit and that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



DESTINEE N ARMSTRONG

Signature and Printed Name

05/28/2024  
Date

MAJISON COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000  
06/06/2024  
10:40:18  
SN 80-001307

Auto Calibration  
Max Power Res Value = 41  
Auto Range Res Value = 36

-----  
<<<< CHANNEL 2 >>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 1.5320 (-0.0080)  
Sample #2 = 1.5300 (-0.0030)  
Sample #3 = 1.5240 (-0.0020)  
Sample #4 = 1.5230 (-0.0060)  
Avg % Abs = 1.5257 (0.0003)  
STD DEV = 0.0038 (0.0049)  
REL STD DEV = 0.248 (1479.865)

-----  
Sol Value = 0.000 g/210L \*\*\*  
Fit value = 0.0000 mg/l %%%  
Samples Taken = 4, Discarded = 1  
Sum Io = 12858, Sum Io = 12766  
-----  
<<<< CHANNEL 1 >>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 0.0950 (0.0061)  
Sample #2 = 0.1120 (0.0220)  
Sample #3 = 0.1290 (0.0151)  
Sample #4 = 0.1020 (0.0531)  
Avg % Abs = 0.1143 (0.0300)  
STD DEV = 0.0137 (1.0202)  
REL STD DEV = 11.939 (67.412)

-----  
<<<< CHANNEL 2 >>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 0.1020 (-0.0010)  
Sample #2 = 0.1180 (-0.0120)  
Sample #3 = 0.1240 (-0.0070)  
Sample #4 = 0.1200 (-0.0060)  
Avg % Abs = 0.1207 (-0.0083)  
STD DEV = 0.0031 (0.0032)  
REL STD DEV = 2.532 (38.575)

-----  
Sol Value = 0.040 g/210L \*\*\*  
Fit value = 0.1905 mg/l %%%  
Samples Taken = 4, Discarded = 1  
Sum Io = 12845, Sum Io = 12763  
-----  
<<<< CHANNEL 1 >>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 0.8750 (-0.0420)  
Sample #2 = 0.8410 (-0.0240)  
Sample #3 = 0.8650 (-0.0150)  
Sample #4 = 0.8560 (0.0000)  
Avg % Abs = 0.8540 (-0.0130)  
STD DEV = 0.0121 (0.0121)  
REL STD DEV = 1.420 (93.264)

-----  
<<<< CHANNEL 2 >>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 6.8310 (-0.0050)  
Sample #2 = 6.8210 (0.0000)  
Sample #3 = 6.8070 (0.0120)  
Sample #4 = 6.8110 (0.0160)  
Avg % Abs = 6.8130 (0.0093)  
STD DEV = 0.0072 (0.0083)  
REL STD DEV = 0.116 (89.214)

-----  
Sol Value = 0.300 g/210L \*\*\*  
Fit value = 1.4286 mg/l %%%  
Samples Taken = 4, Discarded = 1  
Sum Io = 12824, Sum Io = 12751  
-----  
<<<< CHANNEL 1 >>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 5.3830 (-0.0190)  
Sample #2 = 5.3720 (-0.0180)  
Sample #3 = 5.3770 (0.0110)  
Sample #4 = 5.3880 (-0.0030)  
Avg % Abs = 5.3790 (-0.0033)  
STD DEV = 0.0082 (0.0145)  
REL STD DEV = 0.152 (435.086)

-----  
<<<< CHANNEL 2 >>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 9.8770 (0.0120)  
Sample #2 = 9.8810 (0.0160)  
Sample #3 = 9.8970 (0.0150)  
Sample #4 = 9.9180 (0.0020)  
Avg % Abs = 9.8987 (0.0110)  
STD DEV = 0.0186 (0.0078)  
REL STD DEV = 0.187 (71.002)

-----  
Sol Value = 0.200 g/210L \*\*\*  
Fit value = 0.9524 mg/l %%%  
Samples Taken = 4, Discarded = 1  
Sum Io = 12835, Sum Io = 12758  
-----  
<<<< CHANNEL 1 >>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 3.6900 (-0.0160)  
Sample #2 = 3.6690 (-0.0030)  
Sample #3 = 3.6760 (-0.0020)  
Sample #4 = 3.6670 (0.0180)  
Avg % Abs = 3.6707 (0.0043)  
STD DEV = 0.0047 (0.0118)  
REL STD DEV = 0.129 (273.375)

-----  
\*\*\*\*\* AUTO CAL DATA \*\*\*\*\*  
<<<< CHANNEL 1 >>>>  
Sol Val = 0.0000 mg/l or 0.000 g/210L  
% Abs = 0.114  
Std Dev = 0.01 Rel Std Dev = 11.94  
Sol Val = 0.1905 mg/l or 0.040 g/210L  
% Abs = 0.854  
Std Dev = 0.01 Rel Std Dev = 1.42  
Sol Val = 0.4762 mg/l or 0.100 g/210L  
% Abs = 1.944  
Std Dev = 0.01 Rel Std Dev = 0.29  
Sol Val = 0.9524 mg/l or 0.200 g/210L  
% Abs = 3.671  
Std Dev = 0.00 Rel Std Dev = 0.13  
Sol Val = 1.4286 mg/l or 0.300 g/210L  
% Abs = 5.379  
Std Dev = 0.01 Rel Std Dev = 0.15  
Zero Order Coef = -306.24  
First Order Coef = 2570.30  
Second Order Coef = 26.87  
Standard Deviation = 23.947496

-----  
<<<< CHANNEL 2 >>>>  
Sol Val = 0.0000 mg/l or 0.000 g/210L  
% Abs = 0.121  
Std Dev = 0.00 Rel Std Dev = 2.53  
Sol Val = 0.1905 mg/l or 0.040 g/210L  
% Abs = 1.526  
Std Dev = 0.00 Rel Std Dev = 0.25  
Sol Val = 0.4762 mg/l or 0.100 g/210L  
% Abs = 3.583  
Std Dev = 0.01 Rel Std Dev = 0.37  
Sol Val = 0.9524 mg/l or 0.200 g/210L  
% Abs = 6.813  
Std Dev = 0.01 Rel Std Dev = 0.11  
Sol Val = 1.4286 mg/l or 0.300 g/210L  
% Abs = 9.899  
Std Dev = 0.02 Rel Std Dev = 0.19  
Zero Order Coef = -161.83  
First Order Coef = 1331.55  
Second Order Coef = 12.96  
Standard Deviation = 9.571328

CAL ADJUSTMENT  
80-001307 SP  
6/6/24

-----  
Solution Stats Quadratic Fit Chan 1  
Act Fit Residual  
g/210L g/210L g/210L  
-0.000 -0.000 0.0003  
0.040 0.040 -0.0001  
0.100 0.101 -0.0006  
0.200 0.199 0.0007  
0.300 0.300 -0.0002

-----  
Solution Stats Quadratic Fit Chan 2  
Act Fit Residual  
g/210L g/210L g/210L  
0.000 -0.000 0.0000  
0.040 0.040 0.0001  
0.100 0.100 -0.0003  
0.200 0.200 0.0003  
0.300 0.300 -0.0001

-----  
Sol Value = 0.080 g/210L \*\*\*  
Fit value = 0.3810 mg/l %%%  
Samples Taken = 4, Discarded = 1

-----  
\*\*\*\*\* CHANNEL 1 \*\*\*\*\*  
Sample #1 = 3043.00  
Sample #2 = 3077.00  
Sample #3 = 3109.00  
Sample #4 = 3018.00  
Average Result = 3068.0000  
STD DEV = 45.1628  
REL STD DEV = 1.505  
\*\*\*\*\* CHANNEL 2 \*\*\*\*\*  
Sample #1 = 3354.00  
Sample #2 = 3382.00  
Sample #3 = 3362.00  
Sample #4 = 3373.00  
Average Result = 3372.3333  
STD DEV = 11.0167  
REL STD DEV = 0.297

-----  
\*\*\*\*\*  
Dry Gas H2O Adjust Results \*\*\*\*\*  
Barometric Pressure = 1010  
3 in H2O Adjust (mg/l\*10,000) = 741  
9 in H2O Adjust (mg/l\*10,000) = 437  
\*\*\*\* AUTO CAL PASS

# Post Cal Adjust Stability Checks #80-001307

MADISON COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-001307  
06/06/2024  
Software: 8100.27

| Test               | g/210L | Time  |
|--------------------|--------|-------|
| Air Blank          | 0.000  | 12:05 |
| Control Test       | 0.048  | 12:05 |
| Air Blank          | 0.000  | 12:06 |
| Control Test       | 0.047  | 12:06 |
| Air Blank          | 0.000  | 12:07 |
| Control Test       | 0.047  | 12:08 |
| Air Blank          | 0.000  | 12:08 |
| Control Test Stats |        |       |
| Average            | 0.0473 |       |
| Std Dev            | 0.0006 |       |
| Rel Std Dev(%)     | 1.2198 |       |

SP  
Operator's Signature

MADISON COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-001307  
06/06/2024  
Software: 8100.27

| Test               | g/210L | Time  |
|--------------------|--------|-------|
| Air Blank          | 0.000  | 12:14 |
| Control Test       | 0.079  | 12:15 |
| Air Blank          | 0.000  | 12:15 |
| Control Test       | 0.078  | 12:16 |
| Air Blank          | 0.000  | 12:16 |
| Control Test       | 0.079  | 12:17 |
| Air Blank          | 0.000  | 12:18 |
| Control Test Stats |        |       |
| Average            | 0.0787 |       |
| Std Dev            | 0.0006 |       |
| Rel Std Dev(%)     | 0.7339 |       |

SP  
Operator's Signature

MADISON COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-001307  
06/06/2024  
Software: 8100.27

| Test               | g/210L | Time  |
|--------------------|--------|-------|
| Air Blank          | 0.000  | 12:09 |
| Control Test       | 0.197  | 12:10 |
| Air Blank          | 0.000  | 12:11 |
| Control Test       | 0.198  | 12:11 |
| Air Blank          | 0.000  | 12:12 |
| Control Test       | 0.197  | 12:12 |
| Air Blank          | 0.000  | 12:13 |
| Control Test Stats |        |       |
| Average            | 0.1973 |       |
| Std Dev            | 0.0006 |       |
| Rel Std Dev(%)     | 0.2926 |       |

SP  
Operator's Signature

MADISON COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-001307  
06/06/2024  
Software: 8100.27

| Test               | g/210L | Time  |
|--------------------|--------|-------|
| Air Blank          | 0.000  | 12:27 |
| Control Test       | 0.081  | 12:28 |
| Air Blank          | 0.000  | 12:28 |
| Control Test       | 0.081  | 12:28 |
| Air Blank          | 0.000  | 12:29 |
| Control Test       | 0.081  | 12:29 |
| Air Blank          | 0.000  | 12:30 |
| Control Test Stats |        |       |
| Average            | 0.0810 |       |
| Std Dev            | 0.0000 |       |
| Rel Std Dev(%)     | 0.0000 |       |

DGS

SP  
Operator's Signature

# Florida Department of Law Enforcement Alcohol Testing Program

## DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: MADISON COUNTY SO  
Time of Inspection: 10:22

Date of Inspection: 06/14/2024

Serial Number: 80-001307  
Software: 8100.27

| Check or Test                                  | YES | NO | Check or Test                             | YES | NO |
|------------------------------------------------|-----|----|-------------------------------------------|-----|----|
| Diagnostic Check<br>(Pre-Inspection): OK       | Yes |    | Date and/or Time Adjusted                 |     | No |
| Minimum Sample Volume<br>Check: OK             | Yes |    | Barometric Pressure Sensor<br>Check: OK   | Yes |    |
| Alcohol Free Subject<br>Test: 0.000            | Yes |    | Mouth Alcohol Test:<br>Slope Not Met      | Yes |    |
| Interferent Detect Test:<br>Interferent Detect | Yes |    | Diagnostic Check<br>(Post-Inspection): OK | Yes |    |

| Alcohol Free<br>Test<br>(g/210L) | 0.05g/210L Test<br>(g/210L)<br>Lot#:202303K<br>Exp: 03/29/2025 | 0.08g/210L Test<br>(g/210L)<br>Lot#:202303L<br>Exp: 03/29/2025 | 0.20g/210L Test<br>(g/210L)<br>Lot#:202304C<br>Exp: 04/05/2024 | 0.08 g/210L<br>Dry Gas Std Test<br>(g/210L)<br>Lot#:AG310901<br>Exp: 04/19/2025 |
|----------------------------------|----------------------------------------------------------------|----------------------------------------------------------------|----------------------------------------------------------------|---------------------------------------------------------------------------------|
| 0.000                            | 0.050                                                          | 0.079                                                          | 0.197                                                          | 0.080                                                                           |
| 0.000                            | 0.050                                                          | 0.079                                                          | 0.198                                                          | 0.079                                                                           |
| 0.000                            | 0.050                                                          | 0.079                                                          | 0.197                                                          | 0.079                                                                           |
| 0.000                            | 0.050                                                          | 0.080                                                          | 0.197                                                          | 0.080                                                                           |
| 0.000                            | 0.050                                                          | 0.079                                                          | 0.198                                                          | 0.079                                                                           |
| 0.000                            | 0.050                                                          | 0.079                                                          | 0.197                                                          | 0.080                                                                           |
| 0.000                            | 0.050                                                          | 0.079                                                          | 0.197                                                          | 0.079                                                                           |
| 0.000                            | 0.050                                                          | 0.080                                                          | 0.198                                                          | 0.080                                                                           |
| 0.000                            | 0.050                                                          | 0.079                                                          | 0.197                                                          | 0.080                                                                           |
| 0.000                            | 0.050                                                          | 0.080                                                          | 0.198                                                          | 0.079                                                                           |
| Standard Deviations              | 0.0000                                                         | 0.0004                                                         | 0.0005                                                         | 0.0005                                                                          |

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0003 Number of Simulators Used: 5

Remarks:

The above instrument complies ( ☒ ) does not comply ( ☐ ) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

*Shayla Platt*

SHAYLA D PLATT

Signature and Printed Name

06/14/2024  
Date

# Florida Department of Law Enforcement Alcohol Testing Program

## DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: MADISON COUNTY SO  
Time of Inspection: 15:48

Date of Inspection: 06/19/2024

Serial Number: 80-001307  
Software: 8100.27

| Check or Test                                  | YES | NO | Check or Test                             | YES | NO |
|------------------------------------------------|-----|----|-------------------------------------------|-----|----|
| Diagnostic Check<br>(Pre-Inspection): OK       | Yes |    | Date and/or Time Adjusted                 |     | No |
| Minimum Sample Volume<br>Check: OK             | Yes |    | Barometric Pressure Sensor<br>Check: OK   | Yes |    |
| Alcohol Free Subject<br>Test: 0.000            | Yes |    | Mouth Alcohol Test:<br>Slope Not Met      | Yes |    |
| Interferent Detect Test:<br>Interferent Detect | Yes |    | Diagnostic Check<br>(Post-Inspection): OK | Yes |    |

| Alcohol Free<br>Test<br>(g/210L) | 0.05g/210L Test<br>(g/210L)<br>Lot#:202303K<br>Exp: 03/29/2025 | 0.08g/210L Test<br>(g/210L)<br>Lot#:202303L<br>Exp: 03/29/2025 | 0.20g/210L Test<br>(g/210L)<br>Lot#:202304C<br>Exp: 04/05/2025 | 0.08 g/210L<br>Dry Gas Std Test<br>(g/210L)<br>Lot#:06723080A5<br>Exp: 04/05/2025 |
|----------------------------------|----------------------------------------------------------------|----------------------------------------------------------------|----------------------------------------------------------------|-----------------------------------------------------------------------------------|
| 0.000                            | 0.049                                                          | 0.079                                                          | 0.199                                                          | 0.079                                                                             |
| 0.000                            | 0.049                                                          | 0.079                                                          | 0.198                                                          | 0.079                                                                             |
| 0.000                            | 0.049                                                          | 0.079                                                          | 0.199                                                          | 0.080                                                                             |
| 0.000                            | 0.050                                                          | 0.079                                                          | 0.199                                                          | 0.079                                                                             |
| 0.000                            | 0.050                                                          | 0.079                                                          | 0.198                                                          | 0.079                                                                             |
| 0.000                            | 0.050                                                          | 0.079                                                          | 0.198                                                          | 0.079                                                                             |
| 0.000                            | 0.050                                                          | 0.079                                                          | 0.198                                                          | 0.079                                                                             |
| 0.000                            | 0.050                                                          | 0.079                                                          | 0.198                                                          | 0.079                                                                             |
| 0.000                            | 0.050                                                          | 0.079                                                          | 0.198                                                          | 0.079                                                                             |
| 0.000                            | 0.050                                                          | 0.079                                                          | 0.199                                                          | 0.080                                                                             |


|                     |        |        |        |        |
|---------------------|--------|--------|--------|--------|
| Standard Deviations | 0.0004 | 0.0000 | 0.0005 | 0.0004 |
|---------------------|--------|--------|--------|--------|

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0003 Number of Simulators Used: 5

Remarks:

The above instrument complies ( ☒ ) does not comply ( ☐ ) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



SHAYLA D PLATT

Signature and Printed Name

06/19/2024  
Date



# Calibration Certificate

Florida Department of Law Enforcement  
Alcohol Testing Program  
2331 Phillips Road.  
Suite B1032  
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-001307, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

|                   |                          |                                |       |
|-------------------|--------------------------|--------------------------------|-------|
| Serial Number:    | <u>80-001307</u>         | UNCERTAINTY* $\pm$             |       |
| Owning Agency:    | <u>MADISON COUNTY SO</u> | 0.050 g/ 210 L                 | 0.004 |
| Calibration Date: | <u>06/19/2024</u>        | 0.080 g/ 210 L                 | 0.004 |
| Calibration Time: | <u>15:48</u>             | 0.200 g/ 210 L                 | 0.007 |
|                   |                          | 0.080 g/ 210 L Dry Gas Control | 0.005 |

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within  $\pm 0.005$  or 5%, whichever is greater, of the target alcohol concentration.

\*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ( $k=3$ ).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

## TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

**Shayla**

Digitally signed by  
Shayla Platt  
Date: 2024.06.20  
11:42:43 -04'00'

**06/19/2024**

**Platt**

Date

**SHAYLA D PLATT,**

**Department Inspector**

FDLE/ATP Form 69 March 2022

Issuing Authority: Alcohol Testing Program

*Service • Integrity • Respect • Quality*