



INSTRUMENT PROCESSING SHEET

Agency Homestead PDS/N 80-006199Florida Department of
Law EnforcementDate In 02/15/2024 DI Completion Date 02/16/2024☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake	By TDG	Date	Quality Checks	By TDG	Date	Flow Calibration	By TDG	Date																																																											
☒ Annual ☐ Registration ☒ Return from CMI / EE Visual Inspection: ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/ Accessories: ☐ Power cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes: _____ _____ _____ _____ _____ _____ _____ _____ _____ _____		<u>02/15/2024</u>	☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value <u>189</u> ☒ Flow Verification (L/s) Flow Column # <u>ATP106</u> 32 mm <u>0.132*</u> (.139 - .169) 36 mm <u>0.152*</u> (.156 - .190) 53 mm <u>0.214*</u> (.228 - .278) 103 mm <u>0.480</u> (.447 - .547) ☒ Barometric Pressure Check Gauge ID # <u>26932</u> ☒ Stability Checks <table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #/Exp</th></tr></thead><tbody><tr><td>0.050</td><td>MP5094</td><td>202303K 03/29/2025</td></tr><tr><td>0.080</td><td>MP5095</td><td>202303L 03/29/2025</td></tr><tr><td>0.200</td><td>MP5096</td><td>202304C 04/05/2025</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>01923080A3 02/05/2025</td></tr></tbody></table>	Simulator	Serial #	Lot #/Exp	0.050	MP5094	202303K 03/29/2025	0.080	MP5095	202303L 03/29/2025	0.200	MP5096	202304C 04/05/2025	0.080 DGS	N/A	01923080A3 02/05/2025		<u>02/15/2024</u>	Flow Column # <u>ATP101 (x3)</u> ☒ 5L/min - 17mm ☒ 15L/min - 53mm ☒ 30L/min - 103mm ☒ R-Value <u>191 / 192 / 193</u> ☒ Post Calibration Verification (L/s) Flow Column # <u>ATP106 (x3)</u> 32 mm <u>0.136* / 0.132* / 0.136*</u> (.139 - .169) 36 mm <u>0.156 / 0.144* / 0.144*</u> (.156 - .190) 53 mm <u>0.230 / 0.214* / 0.222*</u> (.228 - .278) 103 mm <u>0.496 / 0.507 / 0.500</u> (.447 - .547) Maintenance By _____ Date _____ ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other _____ _____ _____ _____ _____ _____		<u>02/15/2024</u>																																												
Simulator	Serial #	Lot #/Exp																																																																	
0.050	MP5094	202303K 03/29/2025																																																																	
0.080	MP5095	202303L 03/29/2025																																																																	
0.200	MP5096	202304C 04/05/2025																																																																	
0.080 DGS	N/A	01923080A3 02/05/2025																																																																	
Calibration Adjustment By _____			Department Inspection By <u>TDG</u>																																																																
Barometric Pressure Gauge ID # _____ <table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.000</td><td></td><td>N/A</td><td>N/A</td></tr><tr><td>0.040</td><td></td><td></td><td></td></tr><tr><td>0.100</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.300</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table> ☐ Post Calibration Adjustment Stability Checks <table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.050</td><td></td><td></td><td></td></tr><tr><td>0.080</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table>			Simulator	Serial #	Lot #	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			Simulator	Serial #	Lot #	Expiration	0.050				0.080				0.200				0.080 DGS	N/A			Barometric Pressure ID# <u>26932 (x2)</u> Gauge <u>1019 / 1016</u> Instrument <u>1020 / 1017</u> Mouth Alcohol Solution Lot # <u>2023-A (x2)</u> Acetone Stock Solution Lot # <u>2022-B (x2)</u> <table border="1"><thead><tr><th>Simulator</th><th>Serial Number</th></tr></thead><tbody><tr><td>0.000</td><td>MP5092 (x2)</td></tr><tr><td>Interferent</td><td>MP5093 (x2)</td></tr><tr><td>0.050</td><td>MP5094 (x2)</td></tr><tr><td>0.080</td><td>MP5095 (x2)</td></tr><tr><td>0.200</td><td>MP5096 (x2)</td></tr></tbody></table> Attachments <table border="1"><tbody><tr><td>☒ Form 41 (x2) ☒ Stability Checks ☒ Calibration Certificate (x2) ☐ Calibration Adjustment</td><td>☐ Post-Stability Checks ☒ Flow Calibrations ☐ Form 40 ☐ Other _____</td></tr></tbody></table>			Simulator	Serial Number	0.000	MP5092 (x2)	Interferent	MP5093 (x2)	0.050	MP5094 (x2)	0.080	MP5095 (x2)	0.200	MP5096 (x2)	☒ Form 41 (x2) ☒ Stability Checks ☒ Calibration Certificate (x2) ☐ Calibration Adjustment	☐ Post-Stability Checks ☒ Flow Calibrations ☐ Form 40 ☐ Other _____
Simulator	Serial #	Lot #	Expiration																																																																
0.000		N/A	N/A																																																																
0.040																																																																			
0.100																																																																			
0.200																																																																			
0.300																																																																			
0.080 DGS	N/A																																																																		
Simulator	Serial #	Lot #	Expiration																																																																
0.050																																																																			
0.080																																																																			
0.200																																																																			
0.080 DGS	N/A																																																																		
Simulator	Serial Number																																																																		
0.000	MP5092 (x2)																																																																		
Interferent	MP5093 (x2)																																																																		
0.050	MP5094 (x2)																																																																		
0.080	MP5095 (x2)																																																																		
0.200	MP5096 (x2)																																																																		
☒ Form 41 (x2) ☒ Stability Checks ☒ Calibration Certificate (x2) ☐ Calibration Adjustment	☐ Post-Stability Checks ☒ Flow Calibrations ☐ Form 40 ☐ Other _____																																																																		
Notes/Suggested Service: <u>*Flow values were outside the nominal range. (TDG)</u> _____ <u>After the Department Inspection, a flow cal adjust was performed using a different flow column (ATP104). The post-cal R-value was 194, and the post-cal verification passed (see below) using flow column ATP106. (TDG)</u> <u>0.144 (32 mm) / 0.160 (36 mm) / 0.230 (53 mm) / 0.507 (103 mm)</u>			☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use Phil Nicodemo Digitally signed by Phil Nicodemo Date: 2024.02.19 11:14:36 -05'00' Shayla Platt by Shayla Platt Date: 2024.02.20 13:49:14 -05'00' Tech Review / Date Admin Review / Date																																																																

Flow Cal Adjusts

MG

#1

HOMESTEAD PD
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006199
02/15/2024
Software: 8100.27

Flow Rate Calibration*****

1: Rate (Liters/min) = 5
SQRT(Diff)) = 7.547
2: Rate (Liters/min) = 15
SQRT(Diff)) = 11.957
3: Rate (Liters/min) = 30
SQRT(Diff)) = 21.141
Dependent Data Scale Factor = 100000 L/min
Independent Data Scale Factor = 256
Rounded Slope = 706
Rounded Intercept = -782609
Correlation = 0.99638

#2

HOMESTEAD PD
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006199
02/15/2024
Software: 8100.27

Flow Rate Calibration*****

1: Rate (Liters/min) = 5
SQRT(Diff)) = 8.000
2: Rate (Liters/min) = 15
SQRT(Diff)) = 12.203
3: Rate (Liters/min) = 30
SQRT(Diff)) = 21.234
Dependent Data Scale Factor = 100000 L/min
Independent Data Scale Factor = 256
Rounded Slope = 724
Rounded Intercept = -892467
Correlation = 0.99571

#3

HOMESTEAD PD
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006199
02/15/2024
Software: 8100.27

Flow Rate Calibration*****

1: Rate (Liters/min) = 5
SQRT(Diff)) = 7.809
2: Rate (Liters/min) = 15
SQRT(Diff)) = 12.203
3: Rate (Liters/min) = 30
SQRT(Diff)) = 21.211
Dependent Data Scale Factor = 100000 L/min
Independent Data Scale Factor = 256
Rounded Slope = 717
Rounded Intercept = -855682
Correlation = 0.99670

Stability Checks

0.05g/210L	0.08g/210L	0.20g/210L	DGS 0.08g/210L
0.047 to 0.053	0.077 to 0.083	0.194 to 0.206	0.077 to 0.083
✓ HOMESTEAD PD Intoxilyzer - Alcohol Analyzer Model 8000 02/15/2024 Software: 8100.27 SN 80-006199 Test g/210L Time Air Blank 0.000 14:02 Control Test 0.049 14:02 Air Blank 0.000 14:03 Control Test 0.049 14:03 Air Blank 0.000 14:04 Control Test 0.049 14:05 Air Blank 0.000 14:05 Control Test Stats Average 0.0490 Std Dev 0.0000 Rel Std Dev(%) 0.0000	✓ HOMESTEAD PD Intoxilyzer - Alcohol Analyzer Model 8000 02/15/2024 Software: 8100.27 SN 80-006199 Test g/210L Time Air Blank 0.000 14:08 Control Test 0.079 14:08 Air Blank 0.000 14:09 Control Test 0.078 14:09 Air Blank 0.000 14:10 Control Test 0.080 14:11 Air Blank 0.000 14:11 Control Test Stats Average 0.0790 Std Dev 0.0010 Rel Std Dev(%) 1.2658	✓ HOMESTEAD PD Intoxilyzer - Alcohol Analyzer Model 8000 02/15/2024 Software: 8100.27 SN 80-006199 Test g/210L Time Air Blank 0.000 14:14 Control Test 0.197 14:15 Air Blank 0.000 14:15 Control Test 0.198 14:16 Air Blank 0.000 14:17 Control Test 0.197 14:17 Air Blank 0.000 14:18 Control Test Stats Average 0.1973 Std Dev 0.0016 Rel Std Dev(%) 0.2956	✓ HOMESTEAD PD Intoxilyzer - Alcohol Analyzer Model 8000 02/15/2024 Software: 8100.27 SN 80-006199 Test g/210L Time Air Blank 0.000 13:56 Control Test 0.079 13:57 Air Blank 0.000 13:57 Control Test 0.079 13:57 Air Blank 0.000 13:58 Control Test 0.080 13:58 Air Blank 0.000 13:59 Control Test Stats Average 0.0793 Std Dev 0.0006 Rel Std Dev(%) 0.7277
Operator's Signature <i>W</i>	Operator's Signature <i>W</i>	Operator's Signature <i>W</i>	Operator's Signature <i>W</i>

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: HOMESTEAD PD
Time of Inspection: 12:17

Date of Inspection: 02/16/2024

Serial Number: 80-006199
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202303K Exp: 03/29/2025	0.08g/210L Test (g/210L) Lot#:202303L Exp: 03/29/2025	0.20g/210L Test (g/210L) Lot#:202304C Exp: 04/05/2025	0.08 g/210L Dry Gas Std Test* (g/210L) Lot#:01923080A3 Exp: 02/05/2025
0.000	0.049	0.078	0.196	0.080
0.000	0.051	0.078	0.197	0.081
0.000	0.049	0.079	0.197	0.079
0.000	0.050	0.078	0.197	0.079
0.000	0.050	0.078	0.197	0.080
0.000	0.050	0.079	0.197	0.078
0.000	0.050	0.078	0.197	0.079
0.000	0.050	0.079	0.197	0.079
0.000	0.049	0.079	0.197	0.078
0.000	0.050	0.079	0.197	0.079

Standard Deviations	0.0006	0.0005	0.0003	0.0009
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0005 Number of Simulators Used: 5


Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Taylor D Gutschow TAYLOR D GUTSCHOW
Signature and Printed Name

02/16/2024
Date



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
4700 Terminal Drive, Suite 1
Ft. Myers, FL 33907

This is to certify the calibration of Intoxilyzer 8000 serial number 80-006199, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-006199</u>	UNCERTAINTY* $\pm$
Owning Agency:	<u>HOMESTEAD PD</u>	0.050 g/ 210 L 0.004
Calibration Date:	<u>02/16/2024</u>	0.080 g/ 210 L 0.004
Calibration Time:	<u>12:17</u>	0.200 g/ 210 L 0.007
		0.080 g/ 210 L Dry Gas Control 0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence (k=3).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards. This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

02/16/2024

Date

TAYLOR D GUTSCHOW,
Department Inspector

FDLE/ATP Form 69 December 2021

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality

Flow Cal Adjust Conducted After
Department Inspection

2/16/2024

mg

Used ATP104

#4

HOMESTEAD PD
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-006199
02/16/2024
Software: 8100.27

Flow Rate Calibration*****

1: Rate (Liters/min) = 5

SQRT(Diff) = 7.348

2: Rate (Liters/min) = 15

* SQRT(Diff) = 11.703

3: Rate (Liters/min) = 30

SQRT(Diff) = 21.000

Dependent Data Scale Factor = 100000 L/min

Independent Data Scale Factor = 256

Rounded Slope = 702

Rounded Intercept = -732248

Correlation = 0.99586

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: HOMESTEAD PD
Time of Inspection: 15:27

Date of Inspection: 02/16/2024

Serial Number: 80-006199
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202303K Exp: 03/29/2025	0.08g/210L Test (g/210L) Lot#:202303L Exp: 03/29/2025	0.20g/210L Test (g/210L) Lot#:202304C Exp: 04/05/2025	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:01923080A3 Exp: 02/05/2025
0.000	0.048	0.079	0.196	0.079
0.000	0.049	0.078	0.197	0.078
0.000	0.048	0.079	0.198	0.078
0.000	0.049	0.078	0.197	0.078
0.000	0.049	0.079	0.197	0.078
0.000	0.049	0.079	0.198	0.078
0.000	0.049	0.079	0.197	0.078
0.000	0.049	0.079	0.197	0.078
0.000	0.049	0.079	0.198	0.077
0.000	0.050	0.079	0.197	0.079


Standard Deviations	0.0005	0.0004	0.0006	0.0005
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0005 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



TAYLOR D GUTSCHOW

Signature and Printed Name

02/16/2024
Date



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
4700 Terminal Drive, Suite 1
Ft. Myers, FL 33907

This is to certify the calibration of Intoxilyzer 8000 serial number 80-006199, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-006199</u>	UNCERTAINTY* $\pm$
Owning Agency:	<u>HOMESTEAD PD</u>	0.050 g/ 210 L 0.004
Calibration Date:	<u>02/16/2024</u>	0.080 g/ 210 L 0.004
Calibration Time:	<u>15:27</u>	0.200 g/ 210 L 0.007
		0.080 g/ 210 L Dry Gas Control 0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards. This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

02/16/2024

Date



TAYLOR D GUTSCHOW,

Department Inspector

FDLE/ATP Form 69 December 2021

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality