



INSTRUMENT PROCESSING SHEET

Agency Bradford County SOS/N 80-001290Florida Department of
Law EnforcementDate In 06/04/2024 DI Completion Date _____☐ Ship ☐ P/U ☐ H/D ☒ CMI ☐ EE

Intake By <u>ALL</u> Date <u>06/04/2024</u> ☒ Annual ☐ Registration ☐ Return from CMI / EE Visual Inspection: ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/ Accessories: ☒ Power cord ☐ Printer Cable ☐ Static Bag ☐ 12V DC Cable Notes: _____ _____ _____ _____ _____ _____ _____ _____ _____ _____	Quality Checks By <u>ALL</u> Date <u>06/11/2024</u> ☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value <u>173</u> ☒ Flow Verification (L/s) Flow Column # <u>ATP-102</u> 32 mm <u>.152</u> (.139 - .169) 36 mm <u>.171</u> (.156 - .190) 53 mm <u>.238</u> (.228 - .278) 103 mm <u>.492</u> (.447 - .547) ☒ Barometric Pressure Check Gauge ID # <u>28427</u> ☒ Stability Checks <table border="1" style="width:100%; border-collapse: collapse;"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #/Exp</th></tr></thead><tbody><tr><td>0.050</td><td>MP5088</td><td>202303K 03/29/2025</td></tr><tr><td>0.080</td><td>MP6296</td><td>202303L 03/29/2025</td></tr><tr><td>0.200</td><td>MP5090</td><td>202304C 04/05/2025</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>AG310901 04/19/2025</td></tr></tbody></table>	Simulator	Serial #	Lot #/Exp	0.050	MP5088	202303K 03/29/2025	0.080	MP6296	202303L 03/29/2025	0.200	MP5090	202304C 04/05/2025	0.080 DGS	N/A	AG310901 04/19/2025	Flow Calibration By _____ Date _____ Flow Column # _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value _____ ☐ Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547) Maintenance By _____ Date _____ ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other _____ _____ _____ _____ _____ _____																																	
Simulator	Serial #	Lot #/Exp																																																
0.050	MP5088	202303K 03/29/2025																																																
0.080	MP6296	202303L 03/29/2025																																																
0.200	MP5090	202304C 04/05/2025																																																
0.080 DGS	N/A	AG310901 04/19/2025																																																
Calibration Adjustment By <u>SP</u> Barometric Pressure Gauge <u>1011</u> ID # <u>28421</u> <table border="1" style="width:100%; border-collapse: collapse;"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.000</td><td>MP6294</td><td>N/A</td><td>N/A</td></tr><tr><td>0.040</td><td>MP6295</td><td>23400</td><td>10-24-25</td></tr><tr><td>0.100</td><td>MP6296</td><td>23450</td><td>12-5-25</td></tr><tr><td>0.200</td><td>MP6297</td><td>24080</td><td>2-13-26</td></tr><tr><td>0.300</td><td>MP6298</td><td>23410</td><td>11-1-25</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>06723080A5</td><td>4-5-25</td></tr></tbody></table> ☒ Post Calibration Adjustment Stability Checks <table border="1" style="width:100%; border-collapse: collapse;"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.050</td><td>MP5088</td><td>202303K</td><td>3-29-25</td></tr><tr><td>0.080</td><td>MP5089</td><td>202303L</td><td>3-29-25</td></tr><tr><td>0.200</td><td>MP5090</td><td>202304C</td><td>4-5-25</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>AG310901</td><td>4-19-25</td></tr></tbody></table> Notes/Suggested Service: <u>Received inrerferent detect on first .050 stability check, ensured connections were secure, ensured space around instrument was clear repeated check, results nominal. 06/11/2024 ALL</u> <u>First DGS stability checks outside tolerenace, ensured connections were secure, repeated check with results outside tolerance. 06/11/2024 ALL</u> _____			Simulator	Serial #	Lot #	Expiration	0.000	MP6294	N/A	N/A	0.040	MP6295	23400	10-24-25	0.100	MP6296	23450	12-5-25	0.200	MP6297	24080	2-13-26	0.300	MP6298	23410	11-1-25	0.080 DGS	N/A	06723080A5	4-5-25	Simulator	Serial #	Lot #	Expiration	0.050	MP5088	202303K	3-29-25	0.080	MP5089	202303L	3-29-25	0.200	MP5090	202304C	4-5-25	0.080 DGS	N/A	AG310901	4-19-25
Simulator	Serial #	Lot #	Expiration																																															
0.000	MP6294	N/A	N/A																																															
0.040	MP6295	23400	10-24-25																																															
0.100	MP6296	23450	12-5-25																																															
0.200	MP6297	24080	2-13-26																																															
0.300	MP6298	23410	11-1-25																																															
0.080 DGS	N/A	06723080A5	4-5-25																																															
Simulator	Serial #	Lot #	Expiration																																															
0.050	MP5088	202303K	3-29-25																																															
0.080	MP5089	202303L	3-29-25																																															
0.200	MP5090	202304C	4-5-25																																															
0.080 DGS	N/A	AG310901	4-19-25																																															
Department Inspection By <u>SP</u> Barometric Pressure ID# <u>28662</u> Gauge <u>1010</u> Instrument <u>1011</u> Mouth Alcohol Solution Lot # <u>2024-A</u> Acetone Stock Solution Lot # <u>2023-B</u> <table border="1" style="width:100%; border-collapse: collapse;"><thead><tr><th>Simulator</th><th>Serial Number</th></tr></thead><tbody><tr><td>0.000</td><td>MP5086</td></tr><tr><td>Interferent</td><td>MP5087</td></tr><tr><td>0.050</td><td>MP5088</td></tr><tr><td>0.080</td><td>MP5089</td></tr><tr><td>0.200</td><td>MP5090</td></tr></tbody></table> Attachments <table border="1" style="width:100%; border-collapse: collapse;"><tbody><tr><td>☒ Form 41 ☒ Stability Checks ☐ Calibration Certificate ☒ Calibration Adjustment</td><td>☒ Post-Stability Checks ☐ Flow Calibration ☒ Form 40 ☒ Other <u>Form 51</u></td></tr></tbody></table> ☐ Instrument Complies with Chapter 11D-8, FAC ☒ Instrument Does Not Comply with Chapter 11D-8, FAC ☐ Return to/Place into Evidentiary Use ☒ Remain Out of Evidentiary Use ☐ Conduct an Agency Inspection Before Evidentiary Use Tech Review / Date _____ Admin Review / Date _____			Simulator	Serial Number	0.000	MP5086	Interferent	MP5087	0.050	MP5088	0.080	MP5089	0.200	MP5090	☒ Form 41 ☒ Stability Checks ☐ Calibration Certificate ☒ Calibration Adjustment	☒ Post-Stability Checks ☐ Flow Calibration ☒ Form 40 ☒ Other <u>Form 51</u>																																		
Simulator	Serial Number																																																	
0.000	MP5086																																																	
Interferent	MP5087																																																	
0.050	MP5088																																																	
0.080	MP5089																																																	
0.200	MP5090																																																	
☒ Form 41 ☒ Stability Checks ☐ Calibration Certificate ☒ Calibration Adjustment	☒ Post-Stability Checks ☐ Flow Calibration ☒ Form 40 ☒ Other <u>Form 51</u>																																																	

Florida Department of Law Enforcement

Alcohol Testing Program

AGENCY INSPECTION REPORT - INTOXILYZER 8000

Agency: BRADFORD COUNTY SO

Serial Number: 80-001290

Time of Inspection: 12:48

Date of Inspection: 06/11/2024

Software: 8100.27

Check or Test	YES	NO
Date and/or Time Adjusted		No
Diagnostic Check (Pre-Inspection): OK	Yes	
Alcohol Free Subject Test: 0.000		No
Mouth Alcohol Test: Slope Not Met		No
Interferent Detect Test: Interferent Detect		No
Diagnostic Check (Post-Inspection): OK		No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#: Exp:	0.08g/210L Test (g/210L) Lot#: Exp:	0.20g/210L Test (g/210L) Lot#: Exp:	0.08 g/210L Dry Gas Std Test (g/210L) Lot #: Exp:


Number of Simulators Used: 0

Remarks:

BYPASS AI FOR INSPECTION COMPLIANCE NOT DETERMINED

The above instrument complies (X) does not comply () with Chapter 11D-8, FAC.

I certify that I hold a valid Florida Department of Law Enforcement Agency Inspector Permit and that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



ALFONSO L LOWRY

Signature and Printed Name

06/11/2024

Date

Stability checks 80-001290 06/11/24

BRADFORD COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001290
06/11/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:02
Control Test	INT*	13:02
Air Blank	0.000	13:03
Control Test	INT*	13:04
Air Blank	0.000	13:04
Control Test	0.052	13:05
Air Blank	0.000	13:05
Control Test Stats		
Average	0.0533	
Std Dev	0.0023	
Rel Std Dev(%)	4.3301	

*Interferent Detect

AZ

Operator's Signature

BRADFORD COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001290
06/11/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:08
Control Test	0.051	13:09
Air Blank	0.000	13:09
Control Test	0.051	13:10
Air Blank	0.000	13:10
Control Test	0.052	13:11
Air Blank	0.000	13:12
Control Test Stats		
Average	0.0513	
Std Dev	0.0006	
Rel Std Dev(%)	1.1247	

AZ

Operator's Signature

BRADFORD COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001290
06/11/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:13
Control Test	0.081	13:13
Air Blank	0.000	13:14
Control Test	0.081	13:14
Air Blank	0.000	13:15
Control Test	0.080	13:16
Air Blank	0.000	13:16
Control Test Stats		
Average	0.0807	
Std Dev	0.0006	
Rel Std Dev(%)	0.7157	

Wet

AZ

Operator's Signature

BRADFORD COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001290
06/11/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:18
Control Test	0.197	13:18
Air Blank	0.000	13:19
Control Test	0.196	13:20
Air Blank	0.000	13:20
Control Test	0.197	13:21
Air Blank	0.000	13:22
Control Test Stats		
Average	0.1967	
Std Dev	0.0006	
Rel Std Dev(%)	0.2936	

AZ

Operator's Signature

BRADFORD COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001290
06/11/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:25
Control Test	0.045	13:25
Air Blank	0.000	13:26
Control Test	0.044	13:26
Air Blank	0.000	13:26
Control Test	0.041	13:27
Air Blank	0.000	13:27
Control Test Stats		
Average	0.0433	
Std Dev	0.0021	
Rel Std Dev(%)	4.8038	

065

AZ

Operator's Signature

BRADFORD COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001290
06/11/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:29
Control Test	0.042	13:29
Air Blank	0.000	13:30
Control Test	0.042	13:30
Air Blank	0.000	13:31
Control Test	0.042	13:31
Air Blank	0.000	13:32
Control Test Stats		
Average	0.0420	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

065

AZ

Operator's Signature

BRADFORD COUNTY, SC
Intoxilyzer - Alcohol Analyzer
Model 8000
06/26/2024 08:12:48
SN 80-001290

Auto Calibration
Max Power Res Value = 34
Auto Range Res Value = 22

Sol Value = 0.000 g/210L ***
Fit Value = 0.0000 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12825, Sum Io = 14226

Channel 1
Sample % Abs (% Abs Ref)
Sample #1 = 0.2770 (0.0030)
Sample #2 = 0.2470 (0.0410)
Sample #3 = 0.2530 (0.0510)
Sample #4 = 0.2420 (0.0790)
Avg % Abs = 0.2473 (0.0570)
STD DEV = 0.0055 (0.0197)
REL STD DEV = 2.227 (34.557)

Channel 2

Sample % Abs (% Abs Ref)
Sample #1 = 0.2340 (0.0040)
Sample #2 = 0.1970 (0.0340)
Sample #3 = 0.1950 (0.0290)
Sample #4 = 0.1990 (0.0560)
Avg % Abs = 0.1970 (0.0397)
STD DEV = 0.0020 (0.0144)
REL STD DEV = 1.015 (36.213)

Sol Value = 0.040 g/210L ***
Fit Value = 0.1905 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12819, Sum Io = 14219

Channel 1
Sample % Abs (% Abs Ref)
Sample #1 = 1.0120 (0.0120)
Sample #2 = 0.9660 (0.0460)
Sample #3 = 0.9550 (0.0660)
Sample #4 = 0.9590 (0.0320)
Avg % Abs = 0.9577 (0.0466)
STD DEV = 0.0025 (0.0171)
REL STD DEV = 0.263 (35.600)

Channel 2
Sample % Abs (% Abs Ref)
Sample #1 = 1.6130 (0.0380)
Sample #2 = 1.5680 (0.0420)
Sample #3 = 1.5940 (0.0440)
Sample #4 = 1.5750 (0.0320)
Avg % Abs = 1.5790 (0.0393)
STD DEV = 0.0135 (0.0064)
REL STD DEV = 0.852 (15.345)

Sol Value = 0.100 g/210L ***
Fit Value = 0.4762 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12815, Sum Io = 14213

Channel 1
Sample % Abs (% Abs Ref)
Sample #1 = 2.0720 (0.0180)
Sample #2 = 2.0480 (0.0510)
Sample #3 = 2.0420 (0.0640)
Sample #4 = 2.0620 (0.0370)
Avg % Abs = 2.0507 (0.0507)
STD DEV = 0.0103 (0.0135)
REL STD DEV = 0.500 (26.651)

Channel 2

Sample % Abs (% Abs Ref)
Sample #1 = 3.6350 (0.0170)
Sample #2 = 3.5950 (0.0530)
Sample #3 = 3.5830 (0.0560)
Sample #4 = 3.6270 (0.0280)
Avg % Abs = 3.6017 (0.0457)
STD DEV = 0.0227 (0.0154)
REL STD DEV = 0.532 (33.564)

Sol Value = 0.200 g/210L ***
Fit Value = 0.9524 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12811, Sum Io = 14210

Channel 1
Sample % Abs (% Abs Ref)
Sample #1 = 3.7910 (0.0290)
Sample #2 = 3.7460 (0.0340)
Sample #3 = 3.7460 (0.0340)
Sample #4 = 3.7540 (0.0390)
Avg % Abs = 3.7487 (0.0357)
STD DEV = 0.0046 (0.0029)
REL STD DEV = 0.123 (8.094)

Channel 2

Sample % Abs (% Abs Ref)
Sample #1 = 6.8350 (0.0050)
Sample #2 = 6.7320 (0.1010)
Sample #3 = 6.7390 (0.0910)
Sample #4 = 6.7560 (0.0770)
Avg % Abs = 6.7423 (0.0897)
STD DEV = 0.0123 (0.0121)
REL STD DEV = 0.183 (13.445)

Sol Value = 0.300 g/210L ***
Fit Value = 1.4286 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum Io = 12802, Sum Io = 14191

Channel 1
Sample % Abs (% Abs Ref)
Sample #1 = 5.4670 (0.0440)
Sample #2 = 5.4540 (0.0460)
Sample #3 = 5.4250 (0.0440)
Sample #4 = 5.4370 (0.0300)
Avg % Abs = 5.4220 (0.0380)
STD DEV = 0.0167 (0.0072)
REL STD DEV = 0.308 (18.977)

Channel 2

Sample % Abs (% Abs Ref)
Sample #1 = 9.8460 (0.0280)
Sample #2 = 9.7760 (0.0700)
Sample #3 = 9.8170 (0.0570)
Sample #4 = 9.8170 (0.0500)
Avg % Abs = 9.8033 (0.0590)
STD DEV = 0.0237 (0.0101)
REL STD DEV = 0.241 (17.202)

AUTO CAL DATA *****

Channel 1

Sol Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.247
Std Dev = 0.01 Rel Std Dev = 2.23
Sol Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 0.958
Std Dev = 0.00 Rel Std Dev = 0.26
Sol Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 2.051
Std Dev = 0.01 Rel Std Dev = 0.50
Sol Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 3.749
Std Dev = 0.00 Rel Std Dev = 0.12
Sol Val = 1.4286 mg/l or 0.300 g/210L
% Abs = 5.422
Std Dev = 0.02 Rel Std Dev = 0.31
Zero Order Coef = -634.16
First Order Coef = 2589.88
Second Order Coef = 30.10
Standard Deviation = 29.604048

Channel 2

Sol Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.197
Std Dev = 0.00 Rel Std Dev = 1.02
Sol Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 1.579
Std Dev = 0.01 Rel Std Dev = 0.85
Sol Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 3.602
Std Dev = 0.02 Rel Std Dev = 0.63
Sol Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 6.742
Std Dev = 0.01 Rel Std Dev = 0.18
Sol Val = 1.4286 mg/l or 0.300 g/210L
% Abs = 9.803
Std Dev = 0.02 Rel Std Dev = 0.24
Zero Order Coef = -280.11
First Order Coef = 1364.23
Second Order Coef = 12.53
Standard Deviation = 26.179705

CAL ADJUSTMENT *sp*

Dry Gas H2O Adjust Results *****
Barometric Pressure = 1011
3 um H2O Adjust (mg/l*10,000) = 2528
9 um H2O Adjust (mg/l*10,000) = 1879
**** AUTO CAL PASS

Solution Stats Quadratic Fit Chan 1
Act Fit Residual
g/210L g/210L g/210L
0.000 0.000 -0.0002
0.040 0.039 0.0007
0.100 0.101 -0.0009
0.200 0.199 0.0006
0.300 0.300 -0.0002

Solution Stats Quadratic Fit Chan 2
Act Fit Residual
g/210L g/210L g/210L
0.000 -0.000 0.0002
0.040 0.040 -0.0000
0.100 0.101 -0.0007
0.200 0.199 0.0008
0.300 0.300 -0.0003

Sol Value = 0.080 g/210L ***
Fit Value = 0.3810 mg/l %%%
Samples Taken = 4, Discarded = 1
***** CHANNEL 1
Sample #1 = 1701.00
Sample #2 = 1275.00
Sample #3 = 1282.00
Sample #4 = 987.00
Average Result = 1181.3334
STD DEV = 168.3340
REL STD DEV = 14.249

***** CHANNEL 2
Sample #1 = 2096.00
Sample #2 = 2029.00
Sample #3 = 1855.00
Sample #4 = 1907.00
Average Result = 1930.3334
STD DEV = 89.3159
REL STD DEV = 4.627

Dry Gas H2O Adjust Results *****
Barometric Pressure = 1011
3 um H2O Adjust (mg/l*10,000) = 2528
9 um H2O Adjust (mg/l*10,000) = 1879
**** AUTO CAL PASS

Post Cal Adjust Stability Checks

BRADFORD COUNTY SO
Intoxilyzer - Alcotest Analyzer
Model 8000 SN 80-001290
06/26/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:03
Control Test	0.051	10:04
Air Blank	0.000	10:04
Control Test	0.050	10:05
Air Blank	0.000	10:06
Control Test	0.050	10:06
Air Blank	0.000	10:07
Control Test Stats		
Average	0.0503	
Std Dev	0.0006	
Rel Std Dev(%)	1.1471	

SP

Operator's Signature

BRADFORD COUNTY SO
Intoxilyzer - Alcotest Analyzer
Model 8000 SN 80-001290
06/26/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:58
Control Test	0.079	09:59
Air Blank	0.000	10:00
Control Test	0.079	10:00
Air Blank	0.000	10:01
Control Test	0.079	10:01
Air Blank	0.000	10:02
Control Test Stats		
Average	0.0790	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

SP

Operator's Signature

BRADFORD COUNTY SO
Intoxilyzer - Alcotest Analyzer
Model 8000 SN 80-001290
06/26/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:54
Control Test	0.199	09:55
Air Blank	0.000	09:55
Control Test	0.197	09:56
Air Blank	0.000	09:56
Control Test	0.199	09:57
Air Blank	0.000	09:58
Control Test Stats		
Average	0.1983	
Std Dev	0.0112	
Rel Std Dev(%)	0.5622	

SP

Operator's Signature

BRADFORD COUNTY SO
Intoxilyzer - Alcotest Analyzer
Model 8000 SN 80-001290
06/26/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:08
Control Test	0.079	10:09
Air Blank	0.000	10:09
Control Test	0.078	10:10
Air Blank	0.000	10:10
Control Test	0.078	10:10
Air Blank	0.000	10:11
Control Test Stats		
Average	0.0783	
Std Dev	0.0006	
Rel Std Dev(%)	0.7370	

SP

Operator's Signature

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: BRADFORD COUNTY SO
Time of Inspection: 12:25

Date of Inspection: 06/26/2024

Serial Number: 80-001290
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK		No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202303K Exp: 03/29/2025	0.08g/210L Test (g/210L) Lot#:202303L Exp: 03/29/2025	0.20g/210L Test (g/210L) Lot#:202304C Exp: 04/05/2025	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG310901 Exp: 04/19/2025
0.000	0.050	0.081	0.200	0.079 / 0.080
0.000	0.051	0.081	0.201	0.084 / 0.073
0.000	0.051	0.081	0.201	0.087 / 0.066
0.000	0.051	0.081	0.201	0.088 / 0.065
0.000	0.051	0.081	0.201	0.087 / 0.065
0.000	0.050	0.081	0.201	0.088 / 0.067
0.000	0.051	0.081	0.200	0.089 / 0.069
0.000	0.051	0.081	0.201	0.078 / 0.070
0.000	0.051	0.081	0.201	0.078 / 0.070
0.000	0.051	0.081	0.201	0.079 / 0.071
Standard Deviations	0.0004	0.0000	0.0004	0.0046 / 0.0045


Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0013 Number of Simulators Used: 5

Remarks:

08: Control Outside Tolerance, Control Outside Tolerance.

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



SHAYLA D PLATT

Signature and Printed Name

06/26/2024
Date

Return Material Authorization

Ship to:

☒ CMI, Inc.

☐ Enforcement Electronics

Shipment to repair facility authorized by: Jayson Levy on _____

Items Returned: Instrument ☒ Supplies ☐ Other ☐ Describe: _____

Instrument Model: Intoxilyzer 8000 Serial Number: 80-001290

Bill To Address:

Bradford County So

Ship to Address:

FDLE- Tallahassee

Reason for Return:

Instrument fails Dry Gas standard tests. See attached.

Please choose one of the following options:

- ☐ 1. I _____, authorize all repairs.
- ☐ 2. I _____, authorize repairs up to \$_____.
- ☒ 3. I require an estimate **BEFORE** any repairs will be authorized and/ or conducted.

Please contact: Name: Jayson Levy

Phone #: 352-339-0143 Email: _____

ATP Contact Name: Shayla Platt ATP Email: _____



INSTRUMENT PROCESSING SHEET

Agency Bradford County SOS/N 80-001290Florida Department of
Law EnforcementDate In 1/29/2024DI Completion Date 1/31/2024☐ Ship ☒ P/U ☐ H/D ☐ CMI ☐ EE

Intake	By BS	Date	Quality Checks	By BS	Date	Flow Calibration	By	Date																												
☒ Annual ☐ Registration ☐ Return from CMI / EE Visual Inspection: ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/ Accessories: ☒ Power cord ☐ Printer Cable ☐ Static Bag ☐ 12V DC Cable Notes: <u>Dropped off without box</u>		<u>1/29/2024</u>	☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value <u>171</u> ☒ Flow Verification (L/s) Flow Column # <u>ATP-102</u> 32 mm <u>0.160</u> (.139 - .169) 36 mm <u>0.175</u> (.156 - .190) 53 mm <u>0.246</u> (.228 - .278) 103 mm <u>0.507</u> (.447 - .547) ☒ Barometric Pressure Check Gauge ID # <u>28663</u> ☒ Stability Checks		<u>1/31/2024</u>	Flow Column # _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value _____ ☐ Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)																														
			<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #/Exp</th></tr></thead><tbody><tr><td>0.050</td><td>MP5088</td><td>202303K 3/29/2025</td></tr><tr><td>0.080</td><td>MP5089</td><td>202303L 3/29/2025</td></tr><tr><td>0.200</td><td>MP5090</td><td>202304C 4/5/2025</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>AG310901 4/19/2025</td></tr></tbody></table>	Simulator	Serial #	Lot #/Exp	0.050	MP5088	202303K 3/29/2025	0.080	MP5089	202303L 3/29/2025	0.200	MP5090	202304C 4/5/2025	0.080 DGS	N/A	AG310901 4/19/2025																		
Simulator	Serial #	Lot #/Exp																																		
0.050	MP5088	202303K 3/29/2025																																		
0.080	MP5089	202303L 3/29/2025																																		
0.200	MP5090	202304C 4/5/2025																																		
0.080 DGS	N/A	AG310901 4/19/2025																																		
Calibration Adjustment			Department Inspection																																	
Barometric Pressure Gauge _____ ID # _____			Barometric Pressure ID# <u>28663</u>																																	
<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.000</td><td></td><td>N/A</td><td>N/A</td></tr><tr><td>0.040</td><td></td><td></td><td></td></tr><tr><td>0.100</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.300</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table>			Simulator	Serial #	Lot #	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			Gauge <u>1018</u> Instrument <u>1013</u> Mouth Alcohol Solution Lot # <u>2022-A</u> Acetone Stock Solution Lot # <u>2023-B</u>					
Simulator	Serial #	Lot #	Expiration																																	
0.000		N/A	N/A																																	
0.040																																				
0.100																																				
0.200																																				
0.300																																				
0.080 DGS	N/A																																			
☐ Post Calibration Adjustment Stability Checks			<table border="1"><thead><tr><th>Simulator</th><th>Serial Number</th></tr></thead><tbody><tr><td>0.000</td><td>MP5086</td></tr><tr><td>Interferent</td><td>MP5087</td></tr><tr><td>0.050</td><td>MP5088</td></tr><tr><td>0.080</td><td>MP5089</td></tr><tr><td>0.200</td><td>MP5090</td></tr></tbody></table>						Simulator	Serial Number	0.000	MP5086	Interferent	MP5087	0.050	MP5088	0.080	MP5089	0.200	MP5090																
Simulator	Serial Number																																			
0.000	MP5086																																			
Interferent	MP5087																																			
0.050	MP5088																																			
0.080	MP5089																																			
0.200	MP5090																																			
<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.050</td><td></td><td></td><td></td></tr><tr><td>0.080</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table>			Simulator	Serial #	Lot #	Expiration	0.050				0.080				0.200				0.080 DGS	N/A			Attachments ☒ Form 41 ☐ Post-Stability Checks ☒ Stability Checks ☐ Flow Calibration ☒ Calibration Certificate ☒ Form 40 ☐ Calibration Adjustment ☐ Other _____													
Simulator	Serial #	Lot #	Expiration																																	
0.050																																				
0.080																																				
0.200																																				
0.080 DGS	N/A																																			
Notes/Suggested Service: _____ _____ _____ _____ _____ _____			☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use																																	
			Phil Nicodemus Digitally signed by Phil Nicodemus Date: 2024.01.31 14:08:10 -0500 Tech Review / Date _____																																	
			Shayla Platt Digitally signed by Shayla Platt Date: 2024.01.31 15:04:12 -0500 Admin Review / Date _____																																	

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: BRADFORD COUNTY SO
Time of Inspection: 11:54

Date of Inspection: 01/31/2024

Serial Number: 80-001290
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202303K Exp: 03/29/2025	0.08g/210L Test (g/210L) Lot#:202303L Exp: 03/29/2025	0.20g/210L Test (g/210L) Lot#:202304C Exp: 04/05/2025	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG310901 Exp: 04/19/2025
0.000	0.050	0.078	0.196	0.079
0.000	0.050	0.079	0.198	0.078
0.000	0.049	0.079	0.198	0.079
0.000	0.050	0.079	0.197	0.078
0.000	0.050	0.079	0.198	0.078
0.000	0.050	0.079	0.198	0.078
0.000	0.050	0.079	0.198	0.078
0.000	0.050	0.079	0.198	0.078
0.000	0.050	0.079	0.197	0.078
0.000	0.050	0.080	0.197	0.079

Standard Deviations	0.0003	0.0004	0.0007	0.0004
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0004 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



BENJAMIN W SIDDOWAY
Signature and Printed Name

01/31/2024
Date

Florida Department of Law Enforcement Alcohol Testing Program

AGENCY INSPECTION REPORT - INTOXILYZER 8000

Agency: BRADFORD COUNTY SO

Time of Inspection: 09:29

Date of Inspection: 01/31/2024

Serial Number: 80-001290

Software: 8100.27

Check or Test	YES	NO
Date and/or Time Adjusted		No
Diagnostic Check (Pre-Inspection): OK		No
Alcohol Free Subject Test: 0.000		No
Mouth Alcohol Test: Slope Not Met		No
Interferent Detect Test: Interferent Detect		No
Diagnostic Check (Post-Inspection): OK		No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#: Exp:	0.08g/210L Test (g/210L) Lot#: Exp:	0.20g/210L Test (g/210L) Lot#: Exp:	0.08 g/210L Dry Gas Std Test (g/210L) Lot#: Exp:

Number of Simulators Used: _____

Remarks:

COMPLIANCE NOT DETERMINED, BYPASS INSPECTION TO OPERATE

N/A compliance not determined

The above instrument ~~complies (X)~~ does not comply () with Chapter 11D-8, FAC.

I certify that I hold a valid Florida Department of Law Enforcement Agency Inspector Permit and that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Benjamin S. Siddoway

BENJAMIN W SIDDOWAY

Signature and Printed Name


01/31/2024

Date

Stability Checks


BRADFORD COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001290
01/31/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:56
Control Test	0.049	09:57
Air Blank	0.000	09:57
Control Test	0.049	09:58
Air Blank	0.000	09:59
Control Test	0.049	09:59
Air Blank	0.000	10:00
Control Test		
Control Test Stats		
Average	0.0490	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	


Operator's Signature


BRADFORD COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001290
01/31/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:02
Control Test	0.079	10:03
Air Blank	0.000	10:03
Control Test	0.078	10:04
Air Blank	0.000	10:04
Control Test	0.079	10:05
Air Blank	0.000	10:06
Control Test		
Control Test Stats		
Average	0.0787	
Std Dev	0.0006	
Rel Std Dev(%)	0.7339	


Operator's Signature

BRADFORD COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001290
01/31/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:07
Control Test	0.198	10:08
Air Blank	0.000	10:08
Control Test	0.197	10:09
Air Blank	0.000	10:10
Control Test	0.198	10:10
Air Blank	0.000	10:11
Control Test		
Control Test Stats		
Average	0.1977	
Std Dev	0.0006	
Rel Std Dev(%)	0.2921	


Operator's Signature

BRADFORD COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001290
01/31/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:12
Control Test	0.079	10:12
Air Blank	0.000	10:13
Control Test	0.080	10:13
Air Blank	0.000	10:14
Control Test	0.079	10:14
Air Blank	0.000	10:14
Control Test		
Control Test Stats		
Average	0.0793	
Std Dev	0.0006	
Rel Std Dev(%)	0.7277	

DGS


Operator's Signature



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
2331 Phillips Road.
Suite B1032
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-001290 , manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-001290</u>	UNCERTAINTY* $\pm$
Owning Agency:	<u>BRADFORD COUNTY SO</u>	0.050 g/ 210 L 0.004
Calibration Date:	<u>01/31/2024</u>	0.080 g/ 210 L 0.004
Calibration Time:	<u>11:54</u>	0.200 g/ 210 L 0.007
		0.080 g/ 210 L Dry Gas Control 0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.
*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

Benjamin
Siddoway
Digitally signed by
Benjamin Siddoway
Date: 2024.01.31
12:19:13 -05'00'

01/31/2024

Date

BENJAMIN W SIDDOWAY,
Department Inspector

FDLE/ATP Form 69 March 2022

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality