

Alcohol Testing Program - Instrument Processing Sheet

Agency: PINELLAS PARK PD

Instrument Serial Number: 80-001732



Date In: 7/30/2026

DI Completion Date: 8/10/26

☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake By: <u>WKP</u> Date: <u>7/30/2026</u>		Quality Checks By: <u>KTS</u> Date: <u>8/4/2026</u>		Flow Adjustment By: _____																																															
☒ Annual ☐ Dropped Off ☐ Registration ☒ Return from CMI / EE ☐ Training Instrument Visual Inspection ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/Accessories ☐ Power Cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes:		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value: <u>217</u> ☒ Flow Verification (L/s) Flow Column #: <u>ATP103</u> 32 mm <u>0.160</u> (.139-.169) 36 mm <u>0.171</u> (.156-.190) 53 mm <u>0.242</u> (.228-.278) 103 mm <u>0.507</u> (.447-.547) ☒ Barometric Pressure Check Gauge ID #: <u>28427</u> Gauge: <u>1010</u> Instrument: <u>1013</u> ☒ Stability Checks		Flow Column #: _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value: _____ ☐ Post Adjustment Verification (L/S) Flow Column #: _____ 32 mm _____ (.139-.169) 36 mm _____ (.156-.190) 53 mm _____ (.228-.278) 103 mm _____ (.447-.547)																																															
		<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot#/Exp</th></tr></thead><tbody><tr><td>0.050</td><td>MP5088</td><td>202603A 3/3/2028</td></tr><tr><td>0.080</td><td>MP5089</td><td>202603B 3/3/2028</td></tr><tr><td>0.200</td><td>MP5090</td><td>202603D 3/4/2028</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>23825080A1 10/5/2027</td></tr></tbody></table>		Simulator	Serial #	Lot#/Exp	0.050	MP5088	202603A 3/3/2028	0.080	MP5089	202603B 3/3/2028	0.200	MP5090	202603D 3/4/2028	0.080 DGS	N/A	23825080A1 10/5/2027	<table border="1"><thead><tr><th>Maintenance</th><th>By:</th><th>Date:</th></tr></thead><tbody><tr><td>☐ Battery Replacement</td><td></td><td></td></tr><tr><td>☐ Dry Gas Regulator Replacement</td><td></td><td></td></tr><tr><td>☐ Tank Sensor Tare</td><td></td><td></td></tr><tr><td>☐ Breath Tube Replacement</td><td></td><td></td></tr><tr><td>☐ Other:</td><td></td><td></td></tr></tbody></table>		Maintenance	By:	Date:	☐ Battery Replacement			☐ Dry Gas Regulator Replacement			☐ Tank Sensor Tare			☐ Breath Tube Replacement			☐ Other:															
Simulator	Serial #	Lot#/Exp																																																	
0.050	MP5088	202603A 3/3/2028																																																	
0.080	MP5089	202603B 3/3/2028																																																	
0.200	MP5090	202603D 3/4/2028																																																	
0.080 DGS	N/A	23825080A1 10/5/2027																																																	
Maintenance	By:	Date:																																																	
☐ Battery Replacement																																																			
☐ Dry Gas Regulator Replacement																																																			
☐ Tank Sensor Tare																																																			
☐ Breath Tube Replacement																																																			
☐ Other:																																																			
Optical Bench Adjustment		By: <u>IAS</u>		Department Inspection																																															
Barometric Pressure Gauge: <u>1018</u> ID#: <u>34420</u>		Barometric Pressure ID#: <u>28427</u>		By: <u>KTS</u>																																															
<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.000</td><td>MP6294</td><td>N/A</td><td>N/A</td></tr><tr><td>0.040</td><td>MP6295</td><td>25090</td><td>3/11/27</td></tr><tr><td>0.100</td><td>MP6296</td><td>24390</td><td>10/29/26</td></tr><tr><td>0.200</td><td>MP6297</td><td>25020</td><td>1/14/27</td></tr><tr><td>0.300</td><td>MP6298</td><td>24430</td><td>12/10/26</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>AG526603</td><td>9/23/27</td></tr></tbody></table>		Simulator	Serial #	Lot #	Expiration	0.000	MP6294	N/A	N/A	0.040	MP6295	25090	3/11/27	0.100	MP6296	24390	10/29/26	0.200	MP6297	25020	1/14/27	0.300	MP6298	24430	12/10/26	0.080 DGS	N/A	AG526603	9/23/27	<table border="1"><thead><tr><th>Gauge: <u>1014</u></th><th>Instrument: <u>1017</u></th></tr></thead><tbody><tr><td>Mouth Alcohol Solution Lot #: <u>2025-D</u></td><td>Exp: <u>9/25/27</u></td></tr><tr><td>Acetone Stock Solution Lot #: <u>2025-B</u></td><td>Exp: <u>9/22/27</u></td></tr></tbody></table>		Gauge: <u>1014</u>	Instrument: <u>1017</u>	Mouth Alcohol Solution Lot #: <u>2025-D</u>	Exp: <u>9/25/27</u>	Acetone Stock Solution Lot #: <u>2025-B</u>	Exp: <u>9/22/27</u>	<table border="1"><thead><tr><th>Simulator</th><th>Serial Number</th></tr></thead><tbody><tr><td>0.000</td><td>MP5086</td></tr><tr><td>Interferent</td><td>MP5087</td></tr><tr><td>0.050</td><td>MP5088</td></tr><tr><td>0.080</td><td>MP5089</td></tr><tr><td>0.200</td><td>MP5090</td></tr></tbody></table>		Simulator	Serial Number	0.000	MP5086	Interferent	MP5087	0.050	MP5088	0.080	MP5089	0.200	MP5090
Simulator	Serial #	Lot #	Expiration																																																
0.000	MP6294	N/A	N/A																																																
0.040	MP6295	25090	3/11/27																																																
0.100	MP6296	24390	10/29/26																																																
0.200	MP6297	25020	1/14/27																																																
0.300	MP6298	24430	12/10/26																																																
0.080 DGS	N/A	AG526603	9/23/27																																																
Gauge: <u>1014</u>	Instrument: <u>1017</u>																																																		
Mouth Alcohol Solution Lot #: <u>2025-D</u>	Exp: <u>9/25/27</u>																																																		
Acetone Stock Solution Lot #: <u>2025-B</u>	Exp: <u>9/22/27</u>																																																		
Simulator	Serial Number																																																		
0.000	MP5086																																																		
Interferent	MP5087																																																		
0.050	MP5088																																																		
0.080	MP5089																																																		
0.200	MP5090																																																		
☒ Post Optical Bench Adjustment Stability Checks		☒ Form 41		☒ Post-Stability Checks																																															
<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.050</td><td>MP5088</td><td>202603A</td><td>3/3/28</td></tr><tr><td>0.080</td><td>MP5089</td><td>202603B</td><td>3/3/2028</td></tr><tr><td>0.200</td><td>MP5090</td><td>202603D</td><td>3/4/28</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>23725080A1</td><td>10/5/27</td></tr></tbody></table>		Simulator	Serial #	Lot #	Expiration	0.050	MP5088	202603A	3/3/28	0.080	MP5089	202603B	3/3/2028	0.200	MP5090	202603D	3/4/28	0.080 DGS	N/A	23725080A1	10/5/27	☒ Stability Checks		☐ Flow Adjustment																											
Simulator	Serial #	Lot #	Expiration																																																
0.050	MP5088	202603A	3/3/28																																																
0.080	MP5089	202603B	3/3/2028																																																
0.200	MP5090	202603D	3/4/28																																																
0.080 DGS	N/A	23725080A1	10/5/27																																																
Gauge ID #: <u>28427</u> Gauge: <u>1014</u> Instrument: <u>1017</u>		☒ Calibration Certificate		☐ Form 40																																															
		☒ Optical Bench Adjustment		☒ Other: Calibration Certificate Calculations																																															
Notes/Suggested Service:		☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC																																																	
		☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use																																																	
		☒ Conduct an Agency Inspection Before Evidentiary Use																																																	
		S. LeAndra Higginbotham Date: 2026.08.14 10:51:15 -04'00'		Wen-Chi Pierson Date: 2026.08.18 09:26:28 -04'00'																																															
		Tech Review		Admin Review																																															

Stability Checks

0.050 g/210L	0.080 g/210L	0.200 g/210L	DGS 0.080 g/210L																																																																																																																																																
✓ 0.047 to 0.053 g/210L	✓ 0.077 to 0.083 g/210L	✓ 0.194 to 0.206 g/210L	✓ 0.077 to 0.083 g/210L 0.003 g/210L of Wet																																																																																																																																																
PINELLAS PARK PD Intoxilyzer - Alcotest Analyzer Model 8000 08/04/2026 Software: 8100.27 SN 80-001732 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>09:14</td></tr> <tr><td>Control Test</td><td>0.049</td><td>09:14</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:15</td></tr> <tr><td>Control Test</td><td>0.050</td><td>09:16</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:16</td></tr> <tr><td>Control Test</td><td>0.050</td><td>09:17</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:17</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0497</td><td></td></tr> <tr><td>Std Dev</td><td>0.0005</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>1.1625</td><td></td></tr> </tbody> </table>	Test	g/210L	Time	Air Blank	0.000	09:14	Control Test	0.049	09:14	Air Blank	0.000	09:15	Control Test	0.050	09:16	Air Blank	0.000	09:16	Control Test	0.050	09:17	Air Blank	0.000	09:17	Control Test Stats			Average	0.0497		Std Dev	0.0005		Rel Std Dev(%)	1.1625		PINELLAS PARK PD Intoxilyzer - Alcotest Analyzer Model 8000 08/04/2026 Software: 8100.27 SN 80-001732 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>09:20</td></tr> <tr><td>Control Test</td><td>0.079</td><td>09:21</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:21</td></tr> <tr><td>Control Test</td><td>0.080</td><td>09:22</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:22</td></tr> <tr><td>Control Test</td><td>0.078</td><td>09:23</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:24</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0790</td><td></td></tr> <tr><td>Std Dev</td><td>0.0010</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>1.2658</td><td></td></tr> </tbody> </table>	Test	g/210L	Time	Air Blank	0.000	09:20	Control Test	0.079	09:21	Air Blank	0.000	09:21	Control Test	0.080	09:22	Air Blank	0.000	09:22	Control Test	0.078	09:23	Air Blank	0.000	09:24	Control Test Stats			Average	0.0790		Std Dev	0.0010		Rel Std Dev(%)	1.2658		PINELLAS PARK PD Intoxilyzer - Alcotest Analyzer Model 8000 08/04/2026 Software: 8100.27 SN 80-001732 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>09:25</td></tr> <tr><td>Control Test</td><td>0.201</td><td>09:25</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:26</td></tr> <tr><td>Control Test</td><td>0.201</td><td>09:27</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:27</td></tr> <tr><td>Control Test</td><td>0.201</td><td>09:28</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:28</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.2010</td><td></td></tr> <tr><td>Std Dev</td><td>0.0000</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr> </tbody> </table>	Test	g/210L	Time	Air Blank	0.000	09:25	Control Test	0.201	09:25	Air Blank	0.000	09:26	Control Test	0.201	09:27	Air Blank	0.000	09:27	Control Test	0.201	09:28	Air Blank	0.000	09:28	Control Test Stats			Average	0.2010		Std Dev	0.0000		Rel Std Dev(%)	0.0000		PINELLAS PARK PD Intoxilyzer - Alcotest Analyzer Model 8000 08/04/2026 Software: 8100.27 SN 80-001732 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>08:59</td></tr> <tr><td>Control Test</td><td>0.083</td><td>09:00</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:00</td></tr> <tr><td>Control Test</td><td>0.083</td><td>09:00</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:01</td></tr> <tr><td>Control Test</td><td>0.082</td><td>09:01</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:02</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0827</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.6984</td><td></td></tr> </tbody> </table>	Test	g/210L	Time	Air Blank	0.000	08:59	Control Test	0.083	09:00	Air Blank	0.000	09:00	Control Test	0.083	09:00	Air Blank	0.000	09:01	Control Test	0.082	09:01	Air Blank	0.000	09:02	Control Test Stats			Average	0.0827		Std Dev	0.0006		Rel Std Dev(%)	0.6984	
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	09:14																																																																																																																																																	
Control Test	0.049	09:14																																																																																																																																																	
Air Blank	0.000	09:15																																																																																																																																																	
Control Test	0.050	09:16																																																																																																																																																	
Air Blank	0.000	09:16																																																																																																																																																	
Control Test	0.050	09:17																																																																																																																																																	
Air Blank	0.000	09:17																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0497																																																																																																																																																		
Std Dev	0.0005																																																																																																																																																		
Rel Std Dev(%)	1.1625																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	09:20																																																																																																																																																	
Control Test	0.079	09:21																																																																																																																																																	
Air Blank	0.000	09:21																																																																																																																																																	
Control Test	0.080	09:22																																																																																																																																																	
Air Blank	0.000	09:22																																																																																																																																																	
Control Test	0.078	09:23																																																																																																																																																	
Air Blank	0.000	09:24																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0790																																																																																																																																																		
Std Dev	0.0010																																																																																																																																																		
Rel Std Dev(%)	1.2658																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	09:25																																																																																																																																																	
Control Test	0.201	09:25																																																																																																																																																	
Air Blank	0.000	09:26																																																																																																																																																	
Control Test	0.201	09:27																																																																																																																																																	
Air Blank	0.000	09:27																																																																																																																																																	
Control Test	0.201	09:28																																																																																																																																																	
Air Blank	0.000	09:28																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.2010																																																																																																																																																		
Std Dev	0.0000																																																																																																																																																		
Rel Std Dev(%)	0.0000																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	08:59																																																																																																																																																	
Control Test	0.083	09:00																																																																																																																																																	
Air Blank	0.000	09:00																																																																																																																																																	
Control Test	0.083	09:00																																																																																																																																																	
Air Blank	0.000	09:01																																																																																																																																																	
Control Test	0.082	09:01																																																																																																																																																	
Air Blank	0.000	09:02																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0827																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	0.6984																																																																																																																																																		
 Operator's Signature	 Operator's Signature	 Operator's Signature	 Operator's Signature																																																																																																																																																

PINELLAS PARK PD
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-001732
 08/06/2026 08:55:26

Auto Calibration

pg 1 of 2

<<<<< 3um >>>>>			<<<<< 9um >>>>>		
Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1					
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)	
Sample #1	0.0940	(-0.0160)	0.1440	(-0.0040)	
Sample #2	0.0940	(0.0000)	0.1350	(0.0050)	
Sample #3	0.0720	(0.0300)	0.1440	(0.0040)	
Sample #4	0.0500	(0.0440)	0.1550	(-0.0070)	
Avg % Abs	0.0720	(0.0247)	0.1447	(0.0007)	
STD DEV	0.0220	(0.0225)	0.0100	(0.0067)	
REL STD DEV	30.556	(91.134)	6.924	(998.749)	
Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1					
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)	
Sample #1	0.8220	(-0.0100)	1.5990	(-0.0170)	
Sample #2	0.8330	(-0.0060)	1.6010	(-0.0240)	
Sample #3	0.7940	(0.0230)	1.5660	(-0.0090)	
Sample #4	0.8370	(0.0140)	1.5960	(-0.0110)	
Avg % Abs	0.8213	(0.0103)	1.5877	(-0.0147)	
STD DEV	0.0238	(0.0148)	0.0189	(0.0081)	
REL STD DEV	2.892	(143.648)	1.192	(55.531)	
Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1					
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)	
Sample #1	1.8800	(-0.0110)	3.6440	(-0.0170)	
Sample #2	1.8790	(-0.0020)	3.6560	(-0.0060)	
Sample #3	1.8730	(0.0060)	3.6290	(0.0120)	
Sample #4	1.9210	(-0.0040)	3.6940	(-0.0220)	
Avg % Abs	1.8910	(-0.0000)	3.6597	(-0.0053)	
STD DEV	0.0262	(0.0053)	0.0327	(0.0170)	
REL STD DEV	1.383	(6818049536.000)	0.892	(318.934)	
Solution = 0.200 g/210L or 0.9524 mg/l, Samples = 4, Discarded = 1					
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)	
Sample #1	3.6290	(0.0010)	6.9480	(0.0150)	
Sample #2	3.6430	(0.0030)	6.9370	(0.0450)	
Sample #3	3.6570	(0.0120)	6.9320	(0.0310)	
Sample #4	3.6310	(0.0460)	6.9330	(0.0380)	
Avg % Abs	3.6437	(0.0203)	6.9340	(0.0380)	
STD DEV	0.0130	(0.0227)	0.0026	(0.0070)	
REL STD DEV	0.357	(111.536)	0.038	(18.421)	
Solution = 0.300 g/210L or 1.4286 mg/l, Samples = 4, Discarded = 1					
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)	
Sample #1	5.3820	(-0.0140)	10.1550	(-0.0100)	
Sample #2	5.3420	(0.0100)	10.1030	(0.0280)	
Sample #3	5.3230	(0.0310)	10.0590	(0.0550)	
Sample #4	5.3220	(0.0460)	10.1020	(0.0400)	
Avg % Abs	5.3290	(0.0290)	10.0880	(0.0410)	
STD DEV	0.0113	(0.0181)	0.0251	(0.0135)	
REL STD DEV	0.211	(62.356)	0.249	(32.995)	

IAS

PINELLAS PARK PD
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-001732
 08/06/2026 08:55:26

Auto Calibration

pg 2 of 2

<<<< 3um >>>>

 Zero Order Coef -195.34
 First Order Coef 2558.98
 Second Order Coef 29.76

<<<< 9um >>>>

 -203.84
 1317.29
 11.88

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	-0.000	0.0002
0.040	0.040	-0.0005
0.100	0.100	0.0002
0.200	0.200	-0.0000
0.300	0.300	-0.0000

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	-0.000	0.0003
0.040	0.040	-0.0003
0.100	0.100	-0.0003
0.200	0.200	0.0005
0.300	0.300	-0.0002

	<<<< 3um >>>> -----	<<<< 9um >>>> -----
Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1		
Sample		
Sample #1	3326.00	3376.00
Sample #2	3242.00	3354.00
Sample #3	3282.00	3344.00
Sample #4	3243.00	3312.00
Avg	3255.6667	3336.6667
STD DEV	22.8108	21.9393
REL STD DEV	0.701	0.658
H2O adjust (mg/l*10k)	554	473

Barometric Pressure = 1018

*****CALIBRATION SUCCESSFUL***** IAS

Post-Optical

80-001732

Stability Checks

0.050 g/210L	0.080 g/210L	0.200 g/210L	DGS 0.080 g/210L
0.047 to 0.053 g/210L	0.077 to 0.083 g/210L	0.194 to 0.206 g/210L	0.077 to 0.083 g/210L ≤0.003 g/210L of Wet
☒ PINELLAS PARK PD Intoxilyzer - Alcohol Analyzer Model: 8000 08/16/2026 Software: 8100.27 SN 80-001732 Test g/210L Time Air Blank 0.000 11:49 Control Test 0.050 11:50 Air Blank 0.000 11:51 Control Test 0.050 11:51 Air Blank 0.000 11:52 Control Test 0.051 11:52 Air Blank 0.000 11:53 Control Test Stats Average 0.0503 Std Dev 0.0016 Rel Std Dev(%) 1.1471 <i>Isabella Snygar</i> Operator's Signature	☒ PINELLAS PARK PD Intoxilyzer - Alcohol Analyzer Model: 8000 08/16/2026 Software: 8100.27 SN 80-001732 Test g/210L Time Air Blank 0.000 12:14 Control Test 0.078 12:14 Air Blank 0.000 12:15 Control Test 0.079 12:15 Air Blank 0.000 12:16 Control Test 0.079 12:17 Air Blank 0.000 12:17 Control Test Stats Average 0.0787 Std Dev 0.0016 Rel Std Dev(%) 0.7339 <i>Isabella Snygar</i> Operator's Signature	☒ PINELLAS PARK PD Intoxilyzer - Alcohol Analyzer Model: 8000 08/16/2026 Software: 8100.27 SN 80-001732 Test g/210L Time Air Blank 0.000 12:20 Control Test 0.197 12:21 Air Blank 0.000 12:21 Control Test 0.197 12:22 Air Blank 0.000 12:23 Control Test 0.197 12:23 Air Blank 0.000 12:24 Control Test Stats Average 0.1970 Std Dev 0.0000 Rel Std Dev(%) 0.0000 <i>Isabella Snygar</i> Operator's Signature	☒ PINELLAS PARK PD Intoxilyzer - Alcohol Analyzer Model: 8000 08/16/2026 Software: 8100.27 SN 80-001732 Test g/210L Time Air Blank 0.010 12:35 Control Test 0.080 12:37 Air Blank 0.010 12:37 Control Test 0.079 12:38 Air Blank 0.010 12:38 Control Test 0.079 12:38 Air Blank 0.010 12:39 Control Test Stats Average 0.0793 Std Dev 0.0136 Rel Std Dev(%) 0.7277 <i>Isabella Snygar</i> Operator's Signature

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: PINELLAS PARK PD
Time of Inspection: 09:16

Date of Inspection: 08/10/2026

Serial Number: 80-001732
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202603A Exp: 03/03/2028	0.08g/210L Test (g/210L) Lot#:202603B Exp: 03/03/2028	0.20g/210L Test (g/210L) Lot#:202603D Exp: 03/04/2028	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:23825080A1 Exp: 10/05/2027
0.000	0.049	0.078	0.198	0.081
0.000	0.050	0.078	0.198	0.081
0.000	0.050	0.078	0.199	0.081
0.000	0.050	0.078	0.199	0.081
0.000	0.051	0.078	0.199	0.080
0.000	0.050	0.078	0.198	0.081
0.000	0.050	0.077	0.198	0.080
0.000	0.050	0.078	0.199	0.080
0.000	0.050	0.078	0.198	0.081
0.000	0.050	0.078	0.199	0.081


Standard Deviations	0.0004	0.0003	0.0005	0.0004
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0004 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



KATIE T SPEARIN
Signature and Printed Name

08/10/2026
Date



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
2331 Phillips Road
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-001732, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-001732</u>	Average of y $\pm$ Uncertainty*
Owning Agency:	<u>PINELLAS PARK PD</u>	0.050 g/ 210 L 0.004
Calibration Date:	<u>08/10/2026</u>	0.078 g/ 210 L 0.004
Calibration Time:	<u>09:16</u>	0.199 g/ 210 L 0.008
		0.081 g/ 210 L Dry Gas Control 0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence (k=3).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

Kaitlyn Spearin
Digitally signed by Kaitlyn Spearin
Date: 2026.08.11 08:53:58 -0400

08/10/2026

Date

KATIE T SPEARIN,

Department Inspector

FDLE/ATP Form 69 August 2026
Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality

Alcohol Testing Program - Instrument Processing Sheet

Agency: PINELLAS PARK PD Instrument Serial Number: 80-001732

Date In: 4/20/2026 DI Completion Date: N/A ☒ Ship ☐ P/U ☐ H/D ☒ CMI ☐ EE

Intake By: <u>IAS</u> Date: <u>4/20/2026</u> ☒ Annual ☐ Dropped Off ☐ Registration ☐ Return from CMI / EE ☐ Training Instrument Visual Inspection ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/Accessories ☐ Power Cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes:	Quality Checks By: <u>KTS</u> Date: <u>4/22/2026</u> ☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☐ R-Value: _____ ☐ Flow Verification (L/s) Flow Column #: _____ 32 mm _____ (.139-.169) 36 mm _____ (.156-.190) 53 mm _____ (.228-.278) 103 mm _____ (.447-.547) ☐ Barometric Pressure Check Gauge ID #: _____ Gauge: _____ Instrument: _____ ☐ Stability Checks	Flow Adjustment By: _____ Flow Column #: _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value: _____ ☐ Post Adjustment Verification (L/S) Flow Column #: _____ 32 mm _____ (.139-.169) 36 mm _____ (.156-.190) 53 mm _____ (.228-.278) 103 mm _____ (.447-.547)																																								
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 15%;">Simulator</th> <th style="width: 20%;">Serial #</th> <th style="width: 15%;">Lot#/Exp</th> <th style="width: 50%;">Maintenance By: _____ Date: _____</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td></td> <td></td> <td rowspan="5"> ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Tank Sensor Tare ☐ Breath Tube Replacement ☐ Other: _____ </td> </tr> <tr> <td>0.080</td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> </tr> <tr> <td></td> <td></td> <td></td> </tr> </tbody> </table>			Simulator	Serial #	Lot#/Exp	Maintenance By: _____ Date: _____	0.050			☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Tank Sensor Tare ☐ Breath Tube Replacement ☐ Other: _____	0.080			0.200			0.080 DGS	N/A																								
Simulator	Serial #	Lot#/Exp	Maintenance By: _____ Date: _____																																							
0.050			☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Tank Sensor Tare ☐ Breath Tube Replacement ☐ Other: _____																																							
0.080																																										
0.200																																										
0.080 DGS	N/A																																									
Optical Bench Adjustment By: _____ Barometric Pressure Gauge: _____ ID#: _____		Department Inspection By: _____ Barometric Pressure ID#: _____																																								
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td></td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.100</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.300</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table>	Simulator	Serial #	Lot #	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			Gauge: _____ Instrument: _____ Mouth Alcohol Solution Lot #: _____ Exp: _____ Acetone Stock Solution Lot #: _____ Exp: _____ <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 50%;">Simulator</th> <th style="width: 50%;">Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td></td> </tr> <tr> <td>Interferent</td> <td></td> </tr> <tr> <td>0.050</td> <td></td> </tr> <tr> <td>0.080</td> <td></td> </tr> <tr> <td>0.200</td> <td></td> </tr> </tbody> </table>		Simulator	Serial Number	0.000		Interferent		0.050		0.080		0.200	
Simulator	Serial #	Lot #	Expiration																																							
0.000		N/A	N/A																																							
0.040																																										
0.100																																										
0.200																																										
0.300																																										
0.080 DGS	N/A																																									
Simulator	Serial Number																																									
0.000																																										
Interferent																																										
0.050																																										
0.080																																										
0.200																																										
☐ Post Optical Bench Adjustment Stability Checks <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table>		Simulator	Serial #	Lot #	Expiration	0.050				0.080				0.200				0.080 DGS	N/A			Attachments ☐ Form 41 ☐ Stability Checks ☐ Calibration Certificate ☐ Optical Bench Adjustment ☐ Post-Stability Checks ☐ Flow Adjustment ☐ Form 40 ☒ Other: Return Material Authorization																				
Simulator	Serial #	Lot #	Expiration																																							
0.050																																										
0.080																																										
0.200																																										
0.080 DGS	N/A																																									
Gauge ID #: _____ Gauge: _____ Instrument: _____		☒ Instrument Complies with Chapter 11D-8, FAC ☒ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☒ Remain Out of Evidentiary Use ☐ Conduct an Agency Inspection Before Evidentiary Use																																								
Notes/Suggested Service: Instrument will not go into ready mode due to DSP fail. Could not conduct remaining quality checks. KTS 4/22/26 Tech Review: Corrected checkmarked boxes for instrument compliance.KTS 5/7/26		LeAndra Higginbotham Tech Review Digitally signed by LeAndra Higginbotham Date: 2026.05.08 10:00:23 -04'00' Admin Review																																								

Return Material Authorization

Ship to:☒ CMI, Inc.☐ Enforcement ElectronicsShipment to repair facility authorized by: Sterling Adams on 4/22/26Items Returned: ☒ Instrument ☐ Supplies ☐ Other Describe: _____Instrument Model: Intoxilyzer 8000 Serial Number: 80-001732

Bill To Address:
Sterling Adams

Pinellas Park Police Department

7700 59th St N, Pinellas Park, FL 33781

Ship to Address:
FL Dept of Law Enforcement
Alcohol Testing Program

813 B Lake Bradford Rd, FL 32304

Reason for Return:
Instrument will not go into ready mode due to DSP fail.

I require an estimate **BEFORE** any repairs will be authorized and/or conducted

Please contact: Sterling Adams Phone #: 727-946-3695Email: SAdams@Pinellas-park.comATP Contact Name: Katie Spearin ATP Email: KaitlynSpearin@fdle.state.fl.us