

Alcohol Testing Program - Instrument Processing Sheet

Agency: ORANGE COUNTY SO

Instrument Serial Number: 80-001420

Date In: 6/10/2026

DI Completion Date: 6/12/2026

☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake By: <u>WKP</u> Date: <u>6/10/2026</u>		Quality Checks By: <u>IAS</u> Date: <u>6/11/2026</u>		Flow Adjustment By: _____																													
☒ Annual ☐ Dropped Off ☐ Registration ☐ Return from CMI / EE ☐ Training Instrument Visual Inspection ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/Accessories ☐ Power Cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes:		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value: <u>213</u> ☒ Flow Verification (L/s) Flow Column #: <u>ATP102</u> 32 mm <u>.156</u> (.139-.169) 36 mm <u>.171</u> (.156-.190) 53 mm <u>.246</u> (.228-.278) 103 mm <u>.496</u> (.447-.547) ☒ Barometric Pressure Check Gauge ID #: <u>28427</u> Gauge: <u>1011</u> Instrument: <u>1011</u> ☒ Stability Checks		Flow Column #: _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value: _____ ☐ Post Adjustment Verification (L/S) Flow Column #: _____ 32 mm _____ (.139-.169) 36 mm _____ (.156-.190) 53 mm _____ (.228-.278) 103 mm _____ (.447-.547)																													
		<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot#/Exp</th> </tr> </thead> <tbody> <tr> <td rowspan="2">0.050</td> <td rowspan="2">MP6291</td> <td>202603A</td> </tr> <tr> <td>3/3/2028</td> </tr> <tr> <td rowspan="2">0.080</td> <td rowspan="2">MP5089</td> <td>202603B</td> </tr> <tr> <td>3/3/2028</td> </tr> <tr> <td rowspan="2">0.200</td> <td rowspan="2">MP5090</td> <td>202603D</td> </tr> <tr> <td>3/4/2028</td> </tr> <tr> <td rowspan="2">0.080 DGS</td> <td rowspan="2">N/A</td> <td>23825080A1</td> </tr> <tr> <td>10/5/2027</td> </tr> </tbody> </table>		Simulator	Serial #	Lot#/Exp	0.050	MP6291	202603A	3/3/2028	0.080	MP5089	202603B	3/3/2028	0.200	MP5090	202603D	3/4/2028	0.080 DGS	N/A	23825080A1	10/5/2027	Maintenance By: _____ Date: _____ ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Tank Sensor Tare ☐ Breath Tube Replacement ☐ Other:										
Simulator	Serial #	Lot#/Exp																															
0.050	MP6291	202603A																															
		3/3/2028																															
0.080	MP5089	202603B																															
		3/3/2028																															
0.200	MP5090	202603D																															
		3/4/2028																															
0.080 DGS	N/A	23825080A1																															
		10/5/2027																															
Optical Bench Adjustment By: <u>SLH</u> Barometric Pressure Gauge: <u>1011</u> ID#: <u>34420</u>		Department Inspection By: <u>SLH</u> Barometric Pressure ID#: <u>28427</u>																															
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td>MP6294</td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td>MP6295</td> <td>25090</td> <td>03/11/2027</td> </tr> <tr> <td>0.100</td> <td>MP6296</td> <td>25282</td> <td>08/12/2027</td> </tr> <tr> <td>0.200</td> <td>MP6297</td> <td>25020</td> <td>01/14/2027</td> </tr> <tr> <td>0.300</td> <td>MP6298</td> <td>25150</td> <td>05/06/2027</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>AG510701</td> <td>04/17/2027</td> </tr> </tbody> </table>		Simulator	Serial #	Lot #	Expiration	0.000	MP6294	N/A	N/A	0.040	MP6295	25090	03/11/2027	0.100	MP6296	25282	08/12/2027	0.200	MP6297	25020	01/14/2027	0.300	MP6298	25150	05/06/2027	0.080 DGS	N/A	AG510701	04/17/2027	Gauge: <u>1012</u> Instrument: <u>1014</u> Mouth Alcohol Solution Lot #: <u>2025-D</u> Exp: <u>09/25/2027</u> Acetone Stock Solution Lot #: <u>2025-B</u> Exp: <u>09/22/2027</u>			
Simulator	Serial #	Lot #	Expiration																														
0.000	MP6294	N/A	N/A																														
0.040	MP6295	25090	03/11/2027																														
0.100	MP6296	25282	08/12/2027																														
0.200	MP6297	25020	01/14/2027																														
0.300	MP6298	25150	05/06/2027																														
0.080 DGS	N/A	AG510701	04/17/2027																														
☒ Post Optical Bench Adjustment Stability Checks		<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td>MP6291</td> <td>202603A</td> <td>03/03/2028</td> </tr> <tr> <td>0.080</td> <td>MP6292</td> <td>202603B</td> <td>03/03/2028</td> </tr> <tr> <td>0.200</td> <td>MP6293</td> <td>202603D</td> <td>03/04/2028</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>23825080A1</td> <td>10/05/2027</td> </tr> </tbody> </table>				Simulator	Serial #	Lot #	Expiration	0.050	MP6291	202603A	03/03/2028	0.080	MP6292	202603B	03/03/2028	0.200	MP6293	202603D	03/04/2028	0.080 DGS	N/A	23825080A1	10/05/2027								
Simulator	Serial #	Lot #	Expiration																														
0.050	MP6291	202603A	03/03/2028																														
0.080	MP6292	202603B	03/03/2028																														
0.200	MP6293	202603D	03/04/2028																														
0.080 DGS	N/A	23825080A1	10/05/2027																														
Gauge ID #: <u>28427</u> Gauge: <u>1009</u> Instrument: <u>1010</u>		Attachments ☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☒ Optical Bench Adjustment																															
Notes/Suggested Service: Tech Review: Added initials in the inspection portion. SLH 6/19/2026		☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use																															
		<table style="width: 100%;"> <tr> <td style="width: 50%; text-align: center;"> Digitally signed by Taylor Gutschow Date: 2026.06.19 14:10:57 -04'00' </td> <td style="width: 50%; text-align: center;"> Digitally signed by Kaitlyn Spearin Date: 2026.06.22 10:48:14 -04'00' </td> </tr> </table>				Digitally signed by Taylor Gutschow Date: 2026.06.19 14:10:57 -04'00'	Digitally signed by Kaitlyn Spearin Date: 2026.06.22 10:48:14 -04'00'																										
Digitally signed by Taylor Gutschow Date: 2026.06.19 14:10:57 -04'00'	Digitally signed by Kaitlyn Spearin Date: 2026.06.22 10:48:14 -04'00'																																
		<table style="width: 100%;"> <tr> <td style="width: 50%; text-align: center;"> Taylor Gutschow Tech Review </td> <td style="width: 50%; text-align: center;"> Kaitlyn Spearin Admin Review </td> </tr> </table>				Taylor Gutschow Tech Review	Kaitlyn Spearin Admin Review																										
Taylor Gutschow Tech Review	Kaitlyn Spearin Admin Review																																

Stability Checks

0.050 g/210L 0.047 to 0.053 g/210L	0.080 g/210L 0.077 to 0.083 g/210L	0.200 g/210L 0.194 to 0.206 g/210L	DGS 0.080 g/210L 0.077 to 0.083 g/210L ≤0.003 g/210L of Wet																																																																																																																																																
ORANGE COUNTY SO Intoxilyzer - Alconol Analyzer Model 8000 SN 80-001420 06/11/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>12:47</td></tr> <tr><td>Control Test</td><td>0.049</td><td>12:48</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:48</td></tr> <tr><td>Control Test</td><td>0.049</td><td>12:49</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:50</td></tr> <tr><td>Control Test</td><td>0.049</td><td>12:50</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:51</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0490</td><td></td></tr> <tr><td>Std Dev</td><td>0.0000</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr> </tbody> </table> <i>Isabella Hughesworth</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	12:47	Control Test	0.049	12:48	Air Blank	0.000	12:48	Control Test	0.049	12:49	Air Blank	0.000	12:50	Control Test	0.049	12:50	Air Blank	0.000	12:51	Control Test Stats			Average	0.0490		Std Dev	0.0000		Rel Std Dev(%)	0.0000		ORANGE COUNTY SO Intoxilyzer - Alconol Analyzer Model 8000 SN 80-001420 06/11/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>12:59</td></tr> <tr><td>Control Test</td><td>0.077</td><td>12:59</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:00</td></tr> <tr><td>Control Test</td><td>0.077</td><td>13:00</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:01</td></tr> <tr><td>Control Test</td><td>0.077</td><td>13:02</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:02</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0770</td><td></td></tr> <tr><td>Std Dev</td><td>0.0000</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr> </tbody> </table> <i>Isabella Hughesworth</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	12:59	Control Test	0.077	12:59	Air Blank	0.000	13:00	Control Test	0.077	13:00	Air Blank	0.000	13:01	Control Test	0.077	13:02	Air Blank	0.000	13:02	Control Test Stats			Average	0.0770		Std Dev	0.0000		Rel Std Dev(%)	0.0000		ORANGE COUNTY SO Intoxilyzer - Alconol Analyzer Model 8000 SN 80-001420 06/11/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>14:17</td></tr> <tr><td>Control Test</td><td>0.196</td><td>14:18</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:18</td></tr> <tr><td>Control Test</td><td>0.196</td><td>14:19</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:20</td></tr> <tr><td>Control Test</td><td>0.195</td><td>14:20</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:21</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.1957</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.2951</td><td></td></tr> </tbody> </table> <i>Isabella Hughesworth</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	14:17	Control Test	0.196	14:18	Air Blank	0.000	14:18	Control Test	0.196	14:19	Air Blank	0.000	14:20	Control Test	0.195	14:20	Air Blank	0.000	14:21	Control Test Stats			Average	0.1957		Std Dev	0.0006		Rel Std Dev(%)	0.2951		ORANGE COUNTY SO Intoxilyzer - Alconol Analyzer Model 8000 SN 80-001420 06/11/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>12:01</td></tr> <tr><td>Control Test</td><td>0.081</td><td>12:01</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:02</td></tr> <tr><td>Control Test</td><td>0.081</td><td>12:02</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:02</td></tr> <tr><td>Control Test</td><td>0.081</td><td>12:03</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:03</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0810</td><td></td></tr> <tr><td>Std Dev</td><td>0.0000</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr> </tbody> </table> <i>Isabella Hughesworth</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	12:01	Control Test	0.081	12:01	Air Blank	0.000	12:02	Control Test	0.081	12:02	Air Blank	0.000	12:02	Control Test	0.081	12:03	Air Blank	0.000	12:03	Control Test Stats			Average	0.0810		Std Dev	0.0000		Rel Std Dev(%)	0.0000	
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	12:47																																																																																																																																																	
Control Test	0.049	12:48																																																																																																																																																	
Air Blank	0.000	12:48																																																																																																																																																	
Control Test	0.049	12:49																																																																																																																																																	
Air Blank	0.000	12:50																																																																																																																																																	
Control Test	0.049	12:50																																																																																																																																																	
Air Blank	0.000	12:51																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0490																																																																																																																																																		
Std Dev	0.0000																																																																																																																																																		
Rel Std Dev(%)	0.0000																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	12:59																																																																																																																																																	
Control Test	0.077	12:59																																																																																																																																																	
Air Blank	0.000	13:00																																																																																																																																																	
Control Test	0.077	13:00																																																																																																																																																	
Air Blank	0.000	13:01																																																																																																																																																	
Control Test	0.077	13:02																																																																																																																																																	
Air Blank	0.000	13:02																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0770																																																																																																																																																		
Std Dev	0.0000																																																																																																																																																		
Rel Std Dev(%)	0.0000																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	14:17																																																																																																																																																	
Control Test	0.196	14:18																																																																																																																																																	
Air Blank	0.000	14:18																																																																																																																																																	
Control Test	0.196	14:19																																																																																																																																																	
Air Blank	0.000	14:20																																																																																																																																																	
Control Test	0.195	14:20																																																																																																																																																	
Air Blank	0.000	14:21																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.1957																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	0.2951																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	12:01																																																																																																																																																	
Control Test	0.081	12:01																																																																																																																																																	
Air Blank	0.000	12:02																																																																																																																																																	
Control Test	0.081	12:02																																																																																																																																																	
Air Blank	0.000	12:02																																																																																																																																																	
Control Test	0.081	12:03																																																																																																																																																	
Air Blank	0.000	12:03																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0810																																																																																																																																																		
Std Dev	0.0000																																																																																																																																																		
Rel Std Dev(%)	0.0000																																																																																																																																																		

ORANGE COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000
06/11/2026

SN 80-001420 ✓
15:30:48

SMA

Auto Calibration

pg 1 of 2

<<<<< 3um >>>>> <<<<< 9um >>>>>

Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1

Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	0.0600	(-0.0130)	0.1120	(-0.0130)
Sample #2	0.0720	(0.0200)	0.1020	(0.0050)
Sample #3	0.0720	(0.0400)	0.1100	(0.0020)
Sample #4	0.1020	(0.0240)	0.1260	(0.0040)
Avg % Abs	0.0820	(0.0280)	0.1127	(0.0037)
STD DEV	0.0173	(0.0106)	0.0122	(0.0015)
REL STD DEV	21.123	(37.796)	10.846	(41.660)

Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1

Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	0.8430	(-0.0240)	1.5360	(-0.0010)
Sample #2	0.8350	(-0.0020)	1.5300	(-0.0070)
Sample #3	0.8280	(-0.0030)	1.5560	(-0.0180)
Sample #4	0.8300	(0.0110)	1.5360	(-0.0110)
Avg % Abs	0.8310	(0.0020)	1.5407	(-0.0120)
STD DEV	0.0036	(0.0078)	0.0136	(0.0056)
REL STD DEV	0.434	(390.513)	0.884	(46.398)

Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1

Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	1.9460	(-0.0190)	3.6280	(-0.0150)
Sample #2	1.9410	(-0.0060)	3.6200	(-0.0070)
Sample #3	1.9280	(-0.0170)	3.6250	(-0.0120)
Sample #4	1.9600	(-0.0230)	3.6360	(-0.0070)
Avg % Abs	1.9430	(-0.0153)	3.6270	(-0.0087)
STD DEV	0.0161	(0.0086)	0.0082	(0.0029)
REL STD DEV	0.828	(56.228)	0.226	(33.309)

Solution = 0.200 g/210L or 0.9524 mg/l, Samples = 4, Discarded = 1

Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	3.7260	(-0.0060)	6.9800	(-0.0110)
Sample #2	3.7490	(-0.0030)	6.9780	(-0.0090)
Sample #3	3.7310	(-0.0140)	6.9400	(-0.0100)
Sample #4	3.7560	(-0.0150)	6.9550	(0.0120)
Avg % Abs	3.7453	(-0.0107)	6.9577	(-0.0023)
STD DEV	0.0129	(0.0067)	0.0191	(0.0124)
REL STD DEV	0.344	(62.422)	0.275	(532.419)

Solution = 0.300 g/210L or 1.4286 mg/l, Samples = 4, Discarded = 1

Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	5.5330	(-0.0130)	10.1480	(-0.0220)
Sample #2	5.5210	(-0.0090)	10.1600	(-0.0270)
Sample #3	5.5330	(0.0000)	10.1380	(-0.0250)
Sample #4	5.5200	(0.0160)	10.1320	(-0.0160)
Avg % Abs	5.5247	(0.0023)	10.1433	(-0.0227)
STD DEV	0.0072	(0.0127)	0.0147	(0.0059)
REL STD DEV	0.131	(542.669)	0.145	(25.851)

ORANGE COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001420 ✓
06/11/2026 15:30:48 *SW*

Auto Calibration

pg 2 of 2

<<<< 3um >>>>

Zero Order Coef -212.45
First Order Coef 2534.93
Second Order Coef 16.31

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	-0.000	0.0001
0.040	0.040	-0.0000
0.100	0.100	-0.0003
0.200	0.200	0.0003
0.300	0.300	-0.0001

<<<< 9um >>>>

-149.32
1316.93
10.47

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	-0.000	0.0000
0.040	0.040	0.0000
0.100	0.100	-0.0001
0.200	0.200	0.0001
0.300	0.300	-0.0000

<<<< 3um >>>>

Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1
Sample

Sample #1	3498.00
Sample #2	3369.00
Sample #3	3382.00
Sample #4	3395.00
Avg	3382.0000
STD DEV	13.0000
REL STD DEV	0.384
H2O adjust (mg/l*10k)	427

<<<< 9um >>>>

3530.00
3486.00
3480.00
3489.00
3485.0000
4.5826
0.131
324

Barometric Pressure = 1010

*****CALIBRATION SUCCESSFUL***** ✓

Post-Stability Checks

80-00 **1420**

0.050 g/210L 0.047 to 0.053 g/210L	0.080 g/210L 0.077 to 0.083 g/210L	0.200 g/210L 0.194 to 0.206 g/210L	DGS 0.080 g/210L 0.077 to 0.083 g/210L ≤0.003 g/210L of Wet																																																																																																																																																
ORANGE COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001420 06/11/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>16:42</td></tr> <tr><td>Control Test</td><td>0.050</td><td>16:42</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:43</td></tr> <tr><td>Control Test</td><td>0.050</td><td>16:44</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:44</td></tr> <tr><td>Control Test</td><td>0.050</td><td>16:45</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:45</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0500</td><td></td></tr> <tr><td>Std Dev</td><td>0.0000</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr> </tbody> </table> <i>Shiglot</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	16:42	Control Test	0.050	16:42	Air Blank	0.000	16:43	Control Test	0.050	16:44	Air Blank	0.000	16:44	Control Test	0.050	16:45	Air Blank	0.000	16:45	Control Test Stats			Average	0.0500		Std Dev	0.0000		Rel Std Dev(%)	0.0000		ORANGE COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-101420 06/11/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>16:37</td></tr> <tr><td>Control Test</td><td>0.078</td><td>16:38</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:38</td></tr> <tr><td>Control Test</td><td>0.079</td><td>16:39</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:39</td></tr> <tr><td>Control Test</td><td>0.078</td><td>16:40</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:40</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0783</td><td></td></tr> <tr><td>Std Dev</td><td>0.0005</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.7370</td><td></td></tr> </tbody> </table> <i>Shiglot</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	16:37	Control Test	0.078	16:38	Air Blank	0.000	16:38	Control Test	0.079	16:39	Air Blank	0.000	16:39	Control Test	0.078	16:40	Air Blank	0.000	16:40	Control Test Stats			Average	0.0783		Std Dev	0.0005		Rel Std Dev(%)	0.7370		ORANGE COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001420 06/11/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>16:31</td></tr> <tr><td>Control Test</td><td>0.199</td><td>16:32</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:32</td></tr> <tr><td>Control Test</td><td>0.200</td><td>16:33</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:34</td></tr> <tr><td>Control Test</td><td>0.200</td><td>16:34</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:35</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.1997</td><td></td></tr> <tr><td>Std Dev</td><td>0.0005</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.2692</td><td></td></tr> </tbody> </table> <i>Shiglot</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	16:31	Control Test	0.199	16:32	Air Blank	0.000	16:32	Control Test	0.200	16:33	Air Blank	0.000	16:34	Control Test	0.200	16:34	Air Blank	0.000	16:35	Control Test Stats			Average	0.1997		Std Dev	0.0005		Rel Std Dev(%)	0.2692		ORANGE COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001420 06/11/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>16:51</td></tr> <tr><td>Control Test</td><td>0.080</td><td>16:51</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:52</td></tr> <tr><td>Control Test</td><td>0.080</td><td>16:52</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:53</td></tr> <tr><td>Control Test</td><td>0.080</td><td>16:53</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:53</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0800</td><td></td></tr> <tr><td>Std Dev</td><td>0.0000</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr> </tbody> </table> <i>Shiglot</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	16:51	Control Test	0.080	16:51	Air Blank	0.000	16:52	Control Test	0.080	16:52	Air Blank	0.000	16:53	Control Test	0.080	16:53	Air Blank	0.000	16:53	Control Test Stats			Average	0.0800		Std Dev	0.0000		Rel Std Dev(%)	0.0000	
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	16:42																																																																																																																																																	
Control Test	0.050	16:42																																																																																																																																																	
Air Blank	0.000	16:43																																																																																																																																																	
Control Test	0.050	16:44																																																																																																																																																	
Air Blank	0.000	16:44																																																																																																																																																	
Control Test	0.050	16:45																																																																																																																																																	
Air Blank	0.000	16:45																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0500																																																																																																																																																		
Std Dev	0.0000																																																																																																																																																		
Rel Std Dev(%)	0.0000																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	16:37																																																																																																																																																	
Control Test	0.078	16:38																																																																																																																																																	
Air Blank	0.000	16:38																																																																																																																																																	
Control Test	0.079	16:39																																																																																																																																																	
Air Blank	0.000	16:39																																																																																																																																																	
Control Test	0.078	16:40																																																																																																																																																	
Air Blank	0.000	16:40																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0783																																																																																																																																																		
Std Dev	0.0005																																																																																																																																																		
Rel Std Dev(%)	0.7370																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	16:31																																																																																																																																																	
Control Test	0.199	16:32																																																																																																																																																	
Air Blank	0.000	16:32																																																																																																																																																	
Control Test	0.200	16:33																																																																																																																																																	
Air Blank	0.000	16:34																																																																																																																																																	
Control Test	0.200	16:34																																																																																																																																																	
Air Blank	0.000	16:35																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.1997																																																																																																																																																		
Std Dev	0.0005																																																																																																																																																		
Rel Std Dev(%)	0.2692																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	16:51																																																																																																																																																	
Control Test	0.080	16:51																																																																																																																																																	
Air Blank	0.000	16:52																																																																																																																																																	
Control Test	0.080	16:52																																																																																																																																																	
Air Blank	0.000	16:53																																																																																																																																																	
Control Test	0.080	16:53																																																																																																																																																	
Air Blank	0.000	16:53																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0800																																																																																																																																																		
Std Dev	0.0000																																																																																																																																																		
Rel Std Dev(%)	0.0000																																																																																																																																																		

Verified instrument was able to read gas pressure. ☒

Florida Department of Law Enforcement

Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: ORANGE COUNTY SO
Time of Inspection: 13:47

Date of Inspection: 06/12/2026

Serial Number: 80-001420
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202603A Exp: 03/03/2028	0.08g/210L Test (g/210L) Lot#:202603B Exp: 03/03/2028	0.20g/210L Test (g/210L) Lot#:202603D Exp: 03/04/2028	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:23825080A1 Exp: 10/05/2027
0.000	0.049	0.078	0.197	0.080
0.000	0.050	0.078	0.197	0.080
0.000	0.049	0.078	0.198	0.080
0.000	0.050	0.078	0.197	0.080
0.000	0.050	0.078	0.198	0.080
0.000	0.049	0.078	0.197	0.080
0.000	0.049	0.078	0.198	0.080
0.000	0.049	0.079	0.197	0.080
0.000	0.049	0.078	0.197	0.079
0.000	0.050	0.078	0.198	0.079


Standard Deviations	0.0005	0.0003	0.0005	0.0004
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0004 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.


LEANDRA HIGGINBOTHAM
Signature and Printed Name

06/12/2026
Date



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
2331 Phillips Road
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-001420, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-001420</u>	UNCERTAINTY* $\pm$	
Owning Agency:	<u>ORANGE COUNTY SO</u>	0.050 g/ 210 L	0.004
Calibration Date:	<u>06/12/2026</u>	0.080 g/ 210 L	0.004
Calibration Time:	<u>13:47</u>	0.200 g/ 210 L	0.008
		0.080 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

06/12/2026

Date


LEANDRA HIGGINBOTHAM,
Department Inspector

FDLE/ATP Form 69 January 2026
Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality

Page 1 of 1