

Alcohol Testing Program - Instrument Processing Sheet

Agency: Haines City PD

Instrument Serial Number: 80-001007

Date In: 5/21/2026

DI Completion Date: 5/29/26

☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake By: IAS Date: 5/21/2026		Quality Checks By: IAS Date: 5/22/2026		Flow Adjustment By: SLH																																															
☒ Annual ☐ Dropped Off ☐ Registration ☐ Return from CMI / EE ☐ Training Instrument Visual Inspection ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/Accessories ☐ Power Cord ☐ Printer Cable ☐ Static Bag ☐ 12V DC Cable Notes:		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value: 161 ☒ Flow Verification (L/s) Flow Column #: ATP105 32 mm .132 (.139-.169) 36 mm .152 (.156-.190) 53 mm .234 (.228-.278) 103 mm .507 (.447-.547) ☒ Barometric Pressure Check Gauge ID #: 28427 Gauge: 1011 Instrument: 1015 ☒ Stability Checks <table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot#/Exp</th></tr></thead><tbody><tr><td>0.050</td><td>MP5088</td><td>202603A 3/3/2028</td></tr><tr><td>0.080</td><td>MP5089</td><td>202603B 3/3/2028</td></tr><tr><td>0.200</td><td>MP5090</td><td>202603D 3/4/2028</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>23825080A1 10/5/2027</td></tr></tbody></table>		Simulator	Serial #	Lot#/Exp	0.050	MP5088	202603A 3/3/2028	0.080	MP5089	202603B 3/3/2028	0.200	MP5090	202603D 3/4/2028	0.080 DGS	N/A	23825080A1 10/5/2027	Flow Column #: ATP102 ☒ 5L/min – 17mm ☒ 15L/min – 53mm ☒ 30L/min – 103mm ☒ R-Value: 166 ☒ Post Adjustment Verification (L/S) Flow Column #: ATP105 32 mm 0.156 (.139-.169) 36 mm 0.175 (.156-.190) 53 mm 0.246 (.228-.278) 103 mm 0.503 (.447-.547)																																
Simulator	Serial #	Lot#/Exp																																																	
0.050	MP5088	202603A 3/3/2028																																																	
0.080	MP5089	202603B 3/3/2028																																																	
0.200	MP5090	202603D 3/4/2028																																																	
0.080 DGS	N/A	23825080A1 10/5/2027																																																	
Optical Bench Adjustment By: _____		Department Inspection By: KTS																																																	
Barometric Pressure Gauge: _____ ID#: _____		Barometric Pressure ID#: 28427																																																	
<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.000</td><td></td><td>N/A</td><td>N/A</td></tr><tr><td>0.040</td><td></td><td></td><td></td></tr><tr><td>0.100</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.300</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table> ☐ Post Optical Bench Adjustment Stability Checks		Simulator	Serial #	Lot #	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			<table border="1"><thead><tr><th colspan="2">Gauge: 1010 Instrument: 1013</th></tr><tr><td>Mouth Alcohol Solution Lot #:</td><td>2025-D Exp: 9/25/27</td></tr><tr><td>Acetone Stock Solution Lot #:</td><td>2025-B Exp: 9/22/27</td></tr><tr><th>Simulator</th><th>Serial Number</th></tr><tr><td>0.000</td><td>MP5086</td></tr><tr><td>Interferent</td><td>MP5087</td></tr><tr><td>0.050</td><td>MP5088</td></tr><tr><td>0.080</td><td>MP5089</td></tr><tr><td>0.200</td><td>MP5090</td></tr></thead></table>				Gauge: 1010 Instrument: 1013		Mouth Alcohol Solution Lot #:	2025-D Exp: 9/25/27	Acetone Stock Solution Lot #:	2025-B Exp: 9/22/27	Simulator	Serial Number	0.000	MP5086	Interferent	MP5087	0.050	MP5088	0.080	MP5089	0.200	MP5090
Simulator	Serial #	Lot #	Expiration																																																
0.000		N/A	N/A																																																
0.040																																																			
0.100																																																			
0.200																																																			
0.300																																																			
0.080 DGS	N/A																																																		
Gauge: 1010 Instrument: 1013																																																			
Mouth Alcohol Solution Lot #:	2025-D Exp: 9/25/27																																																		
Acetone Stock Solution Lot #:	2025-B Exp: 9/22/27																																																		
Simulator	Serial Number																																																		
0.000	MP5086																																																		
Interferent	MP5087																																																		
0.050	MP5088																																																		
0.080	MP5089																																																		
0.200	MP5090																																																		
<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.050</td><td></td><td></td><td></td></tr><tr><td>0.080</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table> Gauge ID #: _____ Gauge: _____ Instrument: _____		Simulator	Serial #	Lot #	Expiration	0.050				0.080				0.200				0.080 DGS	N/A			<table border="1"><thead><tr><th colspan="2">Attachments</th></tr><tr><td>☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☐ Optical Bench Adjustment</td><td>☐ Post-Stability Checks ☒ Flow Adjustment ☐ Form 40 ☐ Other:</td></tr></thead></table>				Attachments		☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☐ Optical Bench Adjustment	☐ Post-Stability Checks ☒ Flow Adjustment ☐ Form 40 ☐ Other:																						
Simulator	Serial #	Lot #	Expiration																																																
0.050																																																			
0.080																																																			
0.200																																																			
0.080 DGS	N/A																																																		
Attachments																																																			
☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☐ Optical Bench Adjustment	☐ Post-Stability Checks ☒ Flow Adjustment ☐ Form 40 ☐ Other:																																																		
Notes/Suggested Service:		☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use Taylor Gutschow Digitally signed by Taylor Gutschow Date: 2026.06.04 16:07:23 -0400 Wen-Chi Pierson Digitally signed by Wen-Chi Pierson Date: 2026.06.05 08:55:42 -0400 Tech ReviewAdmin Review																																																	

80-001007

no equipment or
user error so performed
adjustment. *sut*

sut

HAINES CITY PD
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001007
05/22/2026
Software: 8100.27

Flow Rate Calibration*****

1: Rate (Liters/min) = 5

SQRT(Diff) = 6.707

2: Rate (Liters/min) = 15

SQRT(Diff) = 11.574

3: Rate (Liters/min) = 30

SQRT(Diff) = 21.094

Dependent Data Scale Factor = 100000 L/min

Independent Data Scale Factor = 256

Rounded Slope = 670

Rounded Intercept = -584771

Correlation = 0.99758

Stability Checks

0.050 g/210L	0.080 g/210L	0.200 g/210L	DGS 0.080 g/210L																																																																																																																																				
0.047 to 0.053 g/210L ✓	0.077 to 0.083 g/210L ✓	0.194 to 0.206 g/210L ✓	0.077 to 0.083 g/210L ✓ ≤0.003 g/210L of Wet ✓																																																																																																																																				
HAINES CITY PD Intoxilyzer - Alcolon Analyzer Model 8000 SN 80-001007 05/22/2026 Software: 8100.27 Test g/210L Time <table border="1"> <tr><td>Air Blank</td><td>0.000</td><td>15:59</td></tr> <tr><td>Control Test</td><td>0.049</td><td>15:59</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:00</td></tr> <tr><td>Control Test</td><td>0.050</td><td>16:01</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:01</td></tr> <tr><td>Control Test</td><td>0.049</td><td>16:02</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:02</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0493</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>1.1703</td><td></td></tr> </table>	Air Blank	0.000	15:59	Control Test	0.049	15:59	Air Blank	0.000	16:00	Control Test	0.050	16:01	Air Blank	0.000	16:01	Control Test	0.049	16:02	Air Blank	0.000	16:02	Control Test Stats			Average	0.0493		Std Dev	0.0006		Rel Std Dev(%)	1.1703		HAINES CITY PD Intoxilyzer - Alcolon Analyzer Model 8000 SN 80-001007 05/22/2026 Software: 8100.27 Test g/210L Time <table border="1"> <tr><td>Air Blank</td><td>0.000</td><td>16:06</td></tr> <tr><td>Control Test</td><td>0.079</td><td>16:07</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:07</td></tr> <tr><td>Control Test</td><td>0.078</td><td>16:08</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:08</td></tr> <tr><td>Control Test</td><td>0.078</td><td>16:09</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:10</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0783</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.7370</td><td></td></tr> </table>	Air Blank	0.000	16:06	Control Test	0.079	16:07	Air Blank	0.000	16:07	Control Test	0.078	16:08	Air Blank	0.000	16:08	Control Test	0.078	16:09	Air Blank	0.000	16:10	Control Test Stats			Average	0.0783		Std Dev	0.0006		Rel Std Dev(%)	0.7370		HAINES CITY PD Intoxilyzer - Alcolon Analyzer Model 8000 SN 80-001007 05/22/2026 Software: 8100.27 Test g/210L Time <table border="1"> <tr><td>Air Blank</td><td>0.000</td><td>16:16</td></tr> <tr><td>Control Test</td><td>0.198</td><td>16:17</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:17</td></tr> <tr><td>Control Test</td><td>0.197</td><td>16:18</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:19</td></tr> <tr><td>Control Test</td><td>0.197</td><td>16:19</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>16:20</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.1973</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.2926</td><td></td></tr> </table>	Air Blank	0.000	16:16	Control Test	0.198	16:17	Air Blank	0.000	16:17	Control Test	0.197	16:18	Air Blank	0.000	16:19	Control Test	0.197	16:19	Air Blank	0.000	16:20	Control Test Stats			Average	0.1973		Std Dev	0.0006		Rel Std Dev(%)	0.2926		HAINES CITY PD Intoxilyzer - Alcolon Analyzer Model 8000 SN 80-001007 05/22/2026 Software: 8100.27 Test g/210L Time <table border="1"> <tr><td>Air Blank</td><td>0.000</td><td>15:28</td></tr> <tr><td>Control Test</td><td>0.079</td><td>15:29</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:29</td></tr> <tr><td>Control Test</td><td>0.080</td><td>15:29</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:30</td></tr> <tr><td>Control Test</td><td>0.079</td><td>15:30</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:31</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0793</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.7277</td><td></td></tr> </table>	Air Blank	0.000	15:28	Control Test	0.079	15:29	Air Blank	0.000	15:29	Control Test	0.080	15:29	Air Blank	0.000	15:30	Control Test	0.079	15:30	Air Blank	0.000	15:31	Control Test Stats			Average	0.0793		Std Dev	0.0006		Rel Std Dev(%)	0.7277	
Air Blank	0.000	15:59																																																																																																																																					
Control Test	0.049	15:59																																																																																																																																					
Air Blank	0.000	16:00																																																																																																																																					
Control Test	0.050	16:01																																																																																																																																					
Air Blank	0.000	16:01																																																																																																																																					
Control Test	0.049	16:02																																																																																																																																					
Air Blank	0.000	16:02																																																																																																																																					
Control Test Stats																																																																																																																																							
Average	0.0493																																																																																																																																						
Std Dev	0.0006																																																																																																																																						
Rel Std Dev(%)	1.1703																																																																																																																																						
Air Blank	0.000	16:06																																																																																																																																					
Control Test	0.079	16:07																																																																																																																																					
Air Blank	0.000	16:07																																																																																																																																					
Control Test	0.078	16:08																																																																																																																																					
Air Blank	0.000	16:08																																																																																																																																					
Control Test	0.078	16:09																																																																																																																																					
Air Blank	0.000	16:10																																																																																																																																					
Control Test Stats																																																																																																																																							
Average	0.0783																																																																																																																																						
Std Dev	0.0006																																																																																																																																						
Rel Std Dev(%)	0.7370																																																																																																																																						
Air Blank	0.000	16:16																																																																																																																																					
Control Test	0.198	16:17																																																																																																																																					
Air Blank	0.000	16:17																																																																																																																																					
Control Test	0.197	16:18																																																																																																																																					
Air Blank	0.000	16:19																																																																																																																																					
Control Test	0.197	16:19																																																																																																																																					
Air Blank	0.000	16:20																																																																																																																																					
Control Test Stats																																																																																																																																							
Average	0.1973																																																																																																																																						
Std Dev	0.0006																																																																																																																																						
Rel Std Dev(%)	0.2926																																																																																																																																						
Air Blank	0.000	15:28																																																																																																																																					
Control Test	0.079	15:29																																																																																																																																					
Air Blank	0.000	15:29																																																																																																																																					
Control Test	0.080	15:29																																																																																																																																					
Air Blank	0.000	15:30																																																																																																																																					
Control Test	0.079	15:30																																																																																																																																					
Air Blank	0.000	15:31																																																																																																																																					
Control Test Stats																																																																																																																																							
Average	0.0793																																																																																																																																						
Std Dev	0.0006																																																																																																																																						
Rel Std Dev(%)	0.7277																																																																																																																																						
<i>LaBella Smykowski</i> Operator's Signature	<i>LaBella Smykowski</i> Operator's Signature	<i>LaBella Smykowski</i> Operator's Signature	<i>LaBella Smykowski</i> Operator's Signature																																																																																																																																				

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: HAINES CITY PD
Time of Inspection: 13:06

Date of Inspection: 05/29/2026

Serial Number: 80-001007
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202603A Exp: 03/03/2028	0.08g/210L Test (g/210L) Lot#:202603B Exp: 03/03/2028	0.20g/210L Test (g/210L) Lot#:202406N Exp: 06/20/2026	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:23825080A1 Exp: 10/05/2027
0.000	0.049	0.079	0.199	0.078
0.000	0.049	0.079	0.199	0.079
0.000	0.049	0.079	0.199	0.079
0.000	0.049	0.079	0.199	0.079
0.000	0.050	0.079	0.199	0.079
0.000	0.050	0.079	0.199	0.080
0.000	0.050	0.079	0.199	0.080
0.000	0.050	0.080	0.199	0.079
0.000	0.050	0.080	0.199	0.080
0.000	0.050	0.079	0.199	0.079


Standard Deviations	0.0005	0.0004	0.0000	0.0006
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0003 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



KATIE T SPEARIN
Signature and Printed Name

05/29/2026
Date



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
2331 Phillips Road
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-001007, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-001007</u>	UNCERTAINTY* $\pm$	
Owning Agency:	<u>HAINES CITY PD</u>	0.050 g/ 210 L	0.004
Calibration Date:	<u>05/29/2026</u>	0.080 g/ 210 L	0.004
Calibration Time:	<u>13:06</u>	0.200 g/ 210 L	0.008
		0.080 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration. *Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards. This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

FDLE/ATP Form 69 January 2026

Issuing Authority: Alcohol Testing Program

Date

05/29/2026


KATIE T SPEARIN,
Department Inspector

Service • Integrity • Respect • Quality