

Alcohol Testing Program - Instrument Processing Sheet

Agency: OSCEOLA COUNTY SO Instrument Serial Number: 80-004587  
 Date In: 1/8/2026 DI Completion Date: N/A  Ship  P/U  H/D  CMI  EE

<b>Intake By:</b> <u>KTS</u> <b>Date:</b> <u>1/8/2026</u>	<b>Quality Checks By:</b> <u>SLH</u> <b>Date:</b> <u>1/13/2026</u>	<b>Flow Adjustment By:</b> _____															
<input checked="" type="checkbox"/> Annual <input type="checkbox"/> Dropped Off <input type="checkbox"/> Registration <input type="checkbox"/> Return from CMI / EE <input type="checkbox"/> Training Instrument Visual Inspection <input checked="" type="checkbox"/> Case <input checked="" type="checkbox"/> Handle <input checked="" type="checkbox"/> Keyboard <input checked="" type="checkbox"/> Dry Gas Shelf <input checked="" type="checkbox"/> Feet <input checked="" type="checkbox"/> Breath Tube <input checked="" type="checkbox"/> Ports <input checked="" type="checkbox"/> Screws Tight Other Equipment/Accessories <input checked="" type="checkbox"/> Power Cord <input type="checkbox"/> Printer Cable <input checked="" type="checkbox"/> Static Bag <input type="checkbox"/> 12V DC Cable Notes:	<input checked="" type="checkbox"/> Breath Tube Screen <input checked="" type="checkbox"/> Replace External O-Rings <input checked="" type="checkbox"/> Instrument Set Up Verified <input checked="" type="checkbox"/> R-Value: <u>228 *</u> <input checked="" type="checkbox"/> Flow Verification (L/s) Flow Column #: <u>ATP 103 *</u> 32 mm <u>0.160</u> (.139-.169) 36 mm <u>0.175</u> (.156-.190) 53 mm <u>0.242</u> (.228-.278) 103 mm <u>0.515</u> (.447-.547) <input checked="" type="checkbox"/> Barometric Pressure Check Gauge ID #: <u>28427</u> Gauge: <u>1018</u> Instrument: <u>1017</u> <input checked="" type="checkbox"/> Stability Checks	Flow Column #: _____ <input type="checkbox"/> 5L/min – 17mm <input type="checkbox"/> 15L/min – 53mm <input type="checkbox"/> 30L/min – 103mm <input type="checkbox"/> R-Value: _____ <input type="checkbox"/> Post Adjustment Verification (L/S) Flow Column #: _____ 32 mm _____ (.139-.169) 36 mm _____ (.156-.190) 53 mm _____ (.228-.278) 103 mm _____ (.447-.547)															
	<table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot#/Exp</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td>MP5088</td> <td>202406K 06/19/2026</td> </tr> <tr> <td>0.080</td> <td>MP5089</td> <td>202406L 06/19/2026</td> </tr> <tr> <td>0.200</td> <td>MP5090</td> <td>202406N 06/20/2026</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>AG510701 04/17/2027</td> </tr> </tbody> </table>	Simulator	Serial #	Lot#/Exp	0.050	MP5088	202406K 06/19/2026	0.080	MP5089	202406L 06/19/2026	0.200	MP5090	202406N 06/20/2026	0.080 DGS	N/A	AG510701 04/17/2027	<b>Maintenance By:</b> _____ <b>Date:</b> _____ <input type="checkbox"/> Battery Replacement <input type="checkbox"/> Dry Gas Regulator Replacement <input type="checkbox"/> Tank Sensor Tare <input type="checkbox"/> Breath Tube Replacement <input type="checkbox"/> Other:
Simulator	Serial #	Lot#/Exp															
0.050	MP5088	202406K 06/19/2026															
0.080	MP5089	202406L 06/19/2026															
0.200	MP5090	202406N 06/20/2026															
0.080 DGS	N/A	AG510701 04/17/2027															

<b>Optical Bench Adjustment</b> By: <u>WKP/SLH</u>	<b>Department Inspection</b> By: _____																																								
Barometric Pressure Gauge: <u>1028/1024</u> ID#: <u>34419/28662</u>	Barometric Pressure ID#: _____																																								
<table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td>MP6294/6294</td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td>MP6295/6295</td> <td>25090</td> <td>03/11/2027</td> </tr> <tr> <td>0.100</td> <td>MP6296</td> <td>25282</td> <td>08/12/2027</td> </tr> <tr> <td>0.200</td> <td>MP6297</td> <td>25020</td> <td>01/14/2027</td> </tr> <tr> <td>0.300</td> <td>MP6298</td> <td>25150</td> <td>5/6/2027</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>28424080A3</td> <td>11/5/2026</td> </tr> </tbody> </table>	Simulator	Serial #	Lot #	Expiration	0.000	MP6294/6294	N/A	N/A	0.040	MP6295/6295	25090	03/11/2027	0.100	MP6296	25282	08/12/2027	0.200	MP6297	25020	01/14/2027	0.300	MP6298	25150	5/6/2027	0.080 DGS	N/A	28424080A3	11/5/2026	Gauge: _____ Instrument: _____ Mouth Alcohol Solution Lot #: _____ Exp: _____ Acetone Stock Solution Lot #: _____ Exp: _____ <table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td></td> </tr> <tr> <td>Interferent</td> <td></td> </tr> <tr> <td>0.050</td> <td></td> </tr> <tr> <td>0.080</td> <td></td> </tr> <tr> <td>0.200</td> <td></td> </tr> </tbody> </table>	Simulator	Serial Number	0.000		Interferent		0.050		0.080		0.200	
Simulator	Serial #	Lot #	Expiration																																						
0.000	MP6294/6294	N/A	N/A																																						
0.040	MP6295/6295	25090	03/11/2027																																						
0.100	MP6296	25282	08/12/2027																																						
0.200	MP6297	25020	01/14/2027																																						
0.300	MP6298	25150	5/6/2027																																						
0.080 DGS	N/A	28424080A3	11/5/2026																																						
Simulator	Serial Number																																								
0.000																																									
Interferent																																									
0.050																																									
0.080																																									
0.200																																									
<input checked="" type="checkbox"/> Post Optical Bench Adjustment Stability Checks																																									
<table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td>MP6291</td> <td>202406K</td> <td>06/19/2026</td> </tr> <tr> <td>0.080</td> <td>MP6292</td> <td>202406L</td> <td>06/19/2026</td> </tr> <tr> <td>0.200</td> <td>N/A</td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>N/A</td> <td></td> </tr> </tbody> </table>	Simulator	Serial #	Lot #	Expiration	0.050	MP6291	202406K	06/19/2026	0.080	MP6292	202406L	06/19/2026	0.200	N/A			0.080 DGS	N/A	N/A		<b>Attachments</b> <input type="checkbox"/> Form 41 <input checked="" type="checkbox"/> Stability Checks <input type="checkbox"/> Calibration Certificate <input checked="" type="checkbox"/> Optical Bench Adjustment <input checked="" type="checkbox"/> Post-Stability Checks <input type="checkbox"/> Flow Adjustment <input type="checkbox"/> Form 40 <input type="checkbox"/> Other:																				
Simulator	Serial #	Lot #	Expiration																																						
0.050	MP6291	202406K	06/19/2026																																						
0.080	MP6292	202406L	06/19/2026																																						
0.200	N/A																																								
0.080 DGS	N/A	N/A																																							
Gauge ID #: <u>N/A</u> Gauge: _____ Instrument: _____																																									

Notes/Suggested Service: Noted INT for one of the results for the 0.080 wet control. Due to other controls outside tolerance, no additional troubleshooting performed and will complete optical bench adjustment, SLH 1/13/25 * performed flow portion of quality checks on 1/20/2026, SLH 1/20/2026  First Optical Bench Adjustment (OBA) performed by WKP with subsequent failure of the 0.05 and 0.08 wet post-stabilities conducted. Second OBA initiated by SLH but the 0.04 adjustment standard had repeated faults detected similar to 1st OBA-sent to repair. SLH 2/13/2026	<input type="checkbox"/> Instrument Complies with Chapter 11D-8, FAC <input checked="" type="checkbox"/> Instrument Does Not Comply with Chapter 11D-8, FAC  <input type="checkbox"/> Return to/Place into Evidentiary Use <input checked="" type="checkbox"/> Remain Out of Evidentiary Use  <input type="checkbox"/> Conduct an Agency Inspection Before Evidentiary Use
Digitally signed by Taylor Gutschow Date: 2026.02.17 17:23:26 -05'00'	Digitally signed by Shayla Platt Date: 2026.02.20 23:03:56 -05'00'
<b>Taylor Gutschow</b>	<b>Shayla Platt</b>
<b>Tech Review</b>	<b>Admin Review</b>

# Stability Checks

<b>0.050 g/210L</b> 0.047 to 0.053 g/210L	<b>0.080 g/210L</b> 0.077 to 0.083 g/210L	<b>0.200 g/210L</b> 0.194 to 0.206 g/210L	<b>DGS 0.080 g/210L</b> 0.077 to 0.083 g/210L SO.003 g/210L of Wet																																																																																																																																																
<p><b>*NOTED</b></p> <p>OSCEOLA COUNTY SO                      Intoxilyzer - Alcohol Analyzer                      Model: 8000 SN 80-004587                      01/13/2026                      Software: 8100.27</p> <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>13:45</td></tr> <tr><td>Control Test</td><td>0.058</td><td>13:46</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:46</td></tr> <tr><td>Control Test</td><td>0.051</td><td>13:47</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:48</td></tr> <tr><td>Control Test</td><td>0.056</td><td>13:48</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:49</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0550</td><td></td></tr> <tr><td>Std Dev</td><td>0.0036</td><td></td></tr> <tr><td>Rel. Std Dev(%)</td><td>6.5555</td><td></td></tr> </tbody> </table>	Test	g/210L	Time	Air Blank	0.000	13:45	Control Test	0.058	13:46	Air Blank	0.000	13:46	Control Test	0.051	13:47	Air Blank	0.000	13:48	Control Test	0.056	13:48	Air Blank	0.000	13:49	Control Test Stats			Average	0.0550		Std Dev	0.0036		Rel. Std Dev(%)	6.5555		<p><b>*NOTED</b></p> <p>OSCEOLA COUNTY SO                      Intoxilyzer - Alcohol Analyzer                      Model: 8000 SN 80-104587                      01/13/2026                      Software: 8100.27</p> <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>13:39</td></tr> <tr><td>Control Test</td><td>0.084</td><td>13:40</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:41</td></tr> <tr><td>Control Test</td><td>INT*</td><td>13:41</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:42</td></tr> <tr><td>Control Test</td><td>0.084</td><td>13:43</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:43</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0853</td><td></td></tr> <tr><td>Std Dev</td><td>0.0023</td><td></td></tr> <tr><td>Rel. Std Dev(%)</td><td>2.7063</td><td></td></tr> </tbody> </table> <p>*Interferent Detect</p>	Test	g/210L	Time	Air Blank	0.000	13:39	Control Test	0.084	13:40	Air Blank	0.000	13:41	Control Test	INT*	13:41	Air Blank	0.000	13:42	Control Test	0.084	13:43	Air Blank	0.000	13:43	Control Test Stats			Average	0.0853		Std Dev	0.0023		Rel. Std Dev(%)	2.7063		<p><b>*NOTED</b></p> <p>OSCEOLA COUNTY SO                      Intoxilyzer - Alcohol Analyzer                      Model: 8000 SN 80-004587                      01/13/2026                      Software: 8100.27</p> <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>13:32</td></tr> <tr><td>Control Test</td><td>0.200</td><td>13:33</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:34</td></tr> <tr><td>Control Test</td><td>0.207</td><td>13:34</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:35</td></tr> <tr><td>Control Test</td><td>0.205</td><td>13:36</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:36</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.2040</td><td></td></tr> <tr><td>Std Dev</td><td>0.0135</td><td></td></tr> <tr><td>Rel. Std Dev(%)</td><td>1.7674</td><td></td></tr> </tbody> </table>	Test	g/210L	Time	Air Blank	0.000	13:32	Control Test	0.200	13:33	Air Blank	0.000	13:34	Control Test	0.207	13:34	Air Blank	0.000	13:35	Control Test	0.205	13:36	Air Blank	0.000	13:36	Control Test Stats			Average	0.2040		Std Dev	0.0135		Rel. Std Dev(%)	1.7674		<p><b>*NOTED</b></p> <p>DGS</p> <p>OSCEOLA COUNTY SO                      Intoxilyzer - Alcohol Analyzer                      Model: 8000 SN 80-004587                      01/13/2026                      Software: 8100.27</p> <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>13:18</td></tr> <tr><td>Control Test</td><td>0.084</td><td>13:19</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:19</td></tr> <tr><td>Control Test</td><td>0.086</td><td>13:20</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:20</td></tr> <tr><td>Control Test</td><td>0.084</td><td>13:20</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:21</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0847</td><td></td></tr> <tr><td>Std Dev</td><td>0.0012</td><td></td></tr> <tr><td>Rel. Std Dev(%)</td><td>1.3638</td><td></td></tr> </tbody> </table>	Test	g/210L	Time	Air Blank	0.000	13:18	Control Test	0.084	13:19	Air Blank	0.000	13:19	Control Test	0.086	13:20	Air Blank	0.000	13:20	Control Test	0.084	13:20	Air Blank	0.000	13:21	Control Test Stats			Average	0.0847		Std Dev	0.0012		Rel. Std Dev(%)	1.3638	
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	13:45																																																																																																																																																	
Control Test	0.058	13:46																																																																																																																																																	
Air Blank	0.000	13:46																																																																																																																																																	
Control Test	0.051	13:47																																																																																																																																																	
Air Blank	0.000	13:48																																																																																																																																																	
Control Test	0.056	13:48																																																																																																																																																	
Air Blank	0.000	13:49																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0550																																																																																																																																																		
Std Dev	0.0036																																																																																																																																																		
Rel. Std Dev(%)	6.5555																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	13:39																																																																																																																																																	
Control Test	0.084	13:40																																																																																																																																																	
Air Blank	0.000	13:41																																																																																																																																																	
Control Test	INT*	13:41																																																																																																																																																	
Air Blank	0.000	13:42																																																																																																																																																	
Control Test	0.084	13:43																																																																																																																																																	
Air Blank	0.000	13:43																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0853																																																																																																																																																		
Std Dev	0.0023																																																																																																																																																		
Rel. Std Dev(%)	2.7063																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	13:32																																																																																																																																																	
Control Test	0.200	13:33																																																																																																																																																	
Air Blank	0.000	13:34																																																																																																																																																	
Control Test	0.207	13:34																																																																																																																																																	
Air Blank	0.000	13:35																																																																																																																																																	
Control Test	0.205	13:36																																																																																																																																																	
Air Blank	0.000	13:36																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.2040																																																																																																																																																		
Std Dev	0.0135																																																																																																																																																		
Rel. Std Dev(%)	1.7674																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	13:18																																																																																																																																																	
Control Test	0.084	13:19																																																																																																																																																	
Air Blank	0.000	13:19																																																																																																																																																	
Control Test	0.086	13:20																																																																																																																																																	
Air Blank	0.000	13:20																																																																																																																																																	
Control Test	0.084	13:20																																																																																																																																																	
Air Blank	0.000	13:21																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0847																																																																																																																																																		
Std Dev	0.0012																																																																																																																																																		
Rel. Std Dev(%)	1.3638																																																																																																																																																		
<p><i>Shejlat</i>                      Operator's Signature</p>	<p><i>Shejlat</i>                      Operator's Signature</p>	<p><i>Shejlat</i>                      Operator's Signature</p>	<p><i>Shejlat</i>                      Operator's Signature</p>																																																																																																																																																

Needs a optical bench Adjustment. see

OSCEOLA COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-004587  
01/13/2026  
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	13:04
Control Test	0.201	13:05
Air Blank	0.000	13:06
Control Test	0.194	13:06
Air Blank	0.000	13:07
Control Test	0.190	13:08
Air Blank	0.000	13:08
Control Test Stats		
Average	0.1950	
Std Dev	0.0056	
Rel Std Dev(%)	2.8553	

- noted

Checked for leaks, retightened  
and changed o-rings on SIM  
prior to re-run.  
Sut

1/13/2026

  
Operator's Signature

OSCEOLA COUNTY SO  
Intoxilizer - Alcohol Analyzer  
Model 8000  
01/27/2026

SN 80-014587  
13:13:56

Auto Calibration  
Max Power Res Value = 86  
Auto Range Res Value = 58

SoI Value = 0.000 g/210L \*\*\*  
Fit Value = 0.0000 mg/l %%%  
Samples Taken = 4, Discarded = 1  
Sum Io = 12546, Sum Io = 13154

Sample % Abs (% Abs Ref)  
Sample #1 = 0.1760 (-0.0410)  
Sample #2 = 0.0600 (0.0370)  
Sample #3 = 0.1170 (0.0310)  
Sample #4 = 0.1620 (0.0580)  
Avg % Abs = 0.1130 (0.0421)  
STD DEV = 0.0511 (0.0142)  
REL STD DEV = 45.237 (33.756)

Sample % Abs (% Abs Ref)  
Sample #1 = 0.3200 (-0.0270)  
Sample #2 = 0.2620 (0.0191)  
Sample #3 = 0.2880 (0.0150)  
Sample #4 = 0.2520 (0.0210)  
Avg % Abs = 0.2673 (0.0150)  
STD DEV = 0.0186 (0.0087)  
REL STD DEV = 6.951 (58.119)

SoI Value = 0.040 g/210L \*\*\*  
Fit Value = 0.1905 mg/l %%%  
Samples Taken = 4, Discarded = 1  
Sum Io = 12542, Sum Io = 13164

Sample % Abs (% Abs Ref)  
Sample #1 = 0.7860 (-0.0070)  
Sample #2 = 0.7290 (0.0120)  
Sample #3 = 0.7960 (0.0210)  
Sample #4 = 0.7140 (0.0191)  
Avg % Abs = 0.7430 (0.01413)  
STD DEV = 0.0380 (0.0432)  
REL STD DEV = 5.113 (104.631)

Sample % Abs (% Abs Ref)  
Sample #1 = 1.6640 (-0.0080)  
Sample #2 = 1.5710 (0.0540)  
Sample #3 = 1.5510 (0.0340)  
Sample #4 = 1.5960 (0.0340)  
Avg % Abs = 1.5860 (0.0407)  
STD DEV = 0.1132 (0.0115)  
REL STD DEV = 0.834 (28.394)

SoI Value = 0.040 g/210L \*\*\*  
Fit Value = 0.1905 mg/l %%%  
Samples Taken = 4, Discarded = 1  
Sum Io = 12535, Sum Io = 13163

Sample % Abs (% Abs Ref)  
Sample #1 = 0.7580 (-0.0150)  
Sample #2 = 0.6430 (-0.0600)  
Sample #3 = 0.8500 (-0.0930)  
Sample #4 = 0.7580 (-0.0310)  
Avg % Abs = 0.8170 (-0.0613)  
STD DEV = 0.0512 (0.0311)  
REL STD DEV = 6.269 (50.579)

Sample % Abs (% Abs Ref)  
Sample #1 = 1.6220 (-0.0180)  
Sample #2 = 1.9950 (0.0150)  
Sample #3 = 1.9880 (0.0110)  
Sample #4 = 1.5620 (0.0370)  
Avg % Abs = 1.5817 (0.0210)  
STD DEV = 0.0174 (0.0140)  
REL STD DEV = 1.099 (56.567)

SoI Value = 0.040 g/210L \*\*\*  
Fit Value = 0.1905 mg/l %%%  
Samples Taken = 4, Discarded = 1  
Sum Io = 12552, Sum Io = 13165

Sample % Abs (% Abs Ref)  
Sample #1 = 0.9270 (0.0260)  
Sample #2 = 0.7460 (0.1400)  
Sample #3 = 0.7800 (0.1350)  
Sample #4 = 0.7670 (0.1250)  
Avg % Abs = 0.7643 (0.1327)  
STD DEV = 0.0172 (0.0075)  
REL STD DEV = 2.245 (5.657)

Sample % Abs (% Abs Ref)  
Sample #1 = 1.6280 (-0.0110)  
Sample #2 = 1.5870 (0.0530)  
Sample #3 = 1.5600 (0.0500)  
Sample #4 = 1.5660 (0.0590)  
Avg % Abs = 1.5710 (0.0540)  
STD DEV = 0.0142 (0.0046)  
REL STD DEV = 0.902 (8.486)

SoI Value = 0.100 g/210L \*\*\*  
Fit Value = 0.4762 mg/l %%%  
Samples Taken = 4, Discarded = 1  
Sum Io = 12538, Sum Io = 13163

Sample % Abs (% Abs Ref)  
Sample #1 = 1.9380 (-0.0380)  
Sample #2 = 1.8620 (0.0440)  
Sample #3 = 1.8800 (0.0580)  
Sample #4 = 1.7950 (0.0840)  
Avg % Abs = 1.8457 (0.0620)  
STD DEV = 0.0448 (0.0203)  
REL STD DEV = 2.427 (32.738)

Sample % Abs (% Abs Ref)  
Sample #1 = 3.6300 (-0.0100)  
Sample #2 = 3.5190 (0.0660)  
Sample #3 = 3.4980 (0.0790)  
Sample #4 = 3.5190 (0.0740)  
Avg % Abs = 3.5120 (0.0730)  
STD DEV = 0.0121 (0.0066)  
REL STD DEV = 0.345 (8.983)

SoI Value = 0.200 g/210L \*\*\*  
Fit Value = 0.9524 mg/l %%%  
Samples Taken = 4, Discarded = 1  
Sum Io = 12534, Sum Io = 13163

Sample % Abs (% Abs Ref)  
Sample #1 = 3.5440 (0.0080)  
Sample #2 = 3.4840 (0.0630)  
Sample #3 = 3.3830 (0.0710)  
Sample #4 = 3.3930 (0.0710)  
Avg % Abs = 3.4200 (0.0683)  
STD DEV = 0.0557 (0.0046)  
REL STD DEV = 1.627 (6.759)

Sample % Abs (% Abs Ref)  
Sample #1 = 6.800 (-0.0090)  
Sample #2 = 6.6620 (0.1370)  
Sample #3 = 6.5800 (0.1630)  
Sample #4 = 6.6030 (0.1650)  
Avg % Abs = 6.6150 (0.1550)  
STD DEV = 0.0423 (0.0156)  
REL STD DEV = 0.639 (10.078)

SoI Value = 0.300 g/210L \*\*\*  
Fit Value = 1.4286 mg/l %%%  
Samples Taken = 4, Discarded = 1  
Sum Io = 12539, Sum Io = 13160

Sample % Abs (% Abs Ref)  
Sample #1 = 5.2090 (-0.0100)  
Sample #2 = 5.1400 (0.0830)  
Sample #3 = 5.1510 (0.0800)  
Sample #4 = 5.0400 (0.0880)  
Avg % Abs = 5.1103 (0.0837)  
STD DEV = 0.0612 (0.0040)  
REL STD DEV = 1.197 (4.830)

Sample % Abs (% Abs Ref)  
Sample #1 = 9.8060 (-0.0190)  
Sample #2 = 9.6270 (0.1690)  
Sample #3 = 9.5990 (0.2100)  
Sample #4 = 9.5090 (0.2490)  
Avg % Abs = 9.5783 (0.2053)  
STD DEV = 0.0617 (0.0400)  
REL STD DEV = 0.644 (19.110)

\*\*\*\*\* AUTO CAL DATA \*\*\*\*\*  
<<<<< CHANNEL 1 >>>>>

SoI Val = 0.0000 mg/l or 0.000 g/210L  
% Abs = 0.113  
Std Dev = 0.05 Rel Std Dev = 45.24  
SoI Val = 0.1915 mg/l or 0.040 g/210L  
% Abs = 0.764  
Std Dev = 0.02 Rel Std Dev = 2.24  
SoI Val = 0.4762 mg/l or 0.100 g/210L  
% Abs = 1.846  
Std Dev = 0.04 Rel Std Dev = 2.43  
SoI Val = 0.9524 mg/l or 0.200 g/210L  
% Abs = 3.420  
Std Dev = 0.06 Rel Std Dev = 1.63  
SoI Val = 1.4286 mg/l or 0.300 g/210L  
% Abs = 5.110  
Std Dev = 0.06 Rel Std Dev = 1.20  
Zero Order Coef = -314.47  
First Order Coef = 2807.01  
Second Order Coef = 11.00  
Standard Deviation = 97.995741

<<<<< CHANNEL 2 >>>>>

SoI Val = 0.3000 mg/l or 0.000 g/210L  
% Abs = 0.267  
Std Dev = 0.02 Rel Std Dev = 6.95  
SoI Val = 0.1915 mg/l or 0.040 g/210L  
% Abs = 0.571  
Std Dev = 0.01 Rel Std Dev = 0.90  
SoI Val = 0.4762 mg/l or 0.100 g/210L  
% Abs = 3.512  
Std Dev = 0.01 Rel Std Dev = 0.35  
SoI Val = 0.9524 mg/l or 0.200 g/210L  
% Abs = 6.615  
Std Dev = 0.04 Rel Std Dev = 0.64  
SoI Val = 1.4286 mg/l or 0.300 g/210L  
% Abs = 9.578  
Std Dev = 0.06 Rel Std Dev = 0.64  
Zero Order Coef = -372.23  
First Order Coef = 1423.06  
Second Order Coef = 11.16  
Standard Deviation = 8.939223

Solution Stats Quadratic Fit Chan 1

Act	Fit	Residual
g/210L	g/210L	g/210L
0.000	0.000	-0.0001
0.040	0.139	0.0014
0.100	0.103	-0.0030
0.200	0.198	0.0023
0.300	0.301	-0.0007

Optical  
Bench  
calibration  
rejustment  
wkp 1/27/2026  
Page 3 of 2

SUA

Solution Stats Quadratic Fit Chan 2

Act	Fit	Residual
g/2.0L	g/2.0L	g/2.0L
0.000	0.000	-0.0002
0.040	0.040	0.0003
0.100	0.100	-0.0003
0.200	0.200	-0.0001
0.300	0.300	0.0001

Sol Value = 0.080 g/2.0L \*\*\*  
 Fit Value = 0.3810 mg/l %2%  
 Samples Taken = 4, Discard = 0  
 \*\*\*\*\* CHANNEL 1 \*\*\*\*\*  
 Sample #1 = 3252.00  
 Sample #2 = 3230.00  
 Sample #3 = 3352.00  
 Sample #4 = 3440.00  
 Average Result = 3340.6667  
 STD DEV = 105.4577  
 REL STD DEV = 3.157

\*\*\*\*\* CHANNEL 2 \*\*\*\*\*  
 Sample #1 = 3186.00  
 Sample #2 = 3211.00  
 Sample #3 = 3223.00  
 Sample #4 = 3231.00  
 Average Result = 3221.6667  
 STD DEV = 10.0664  
 REL STD DEV = 0.312

\*\*\*\*\*  
 Dry Gas H2O Adjust Results \*\*\*\*\*  
 Barometric Pressure = 1026  
 3 um H2O Adjust (mg/l\*10,000) = 469  
 9 um H2O Adjust (mg/l\*10,000) = 588  
 \*\*\*\*\* AUTO CAL PASS \*\*\*\*\*

optical Bench  
 calibration  
 adjustment wk P 1/27/2016  
 Page 2 of 2

SUA

80-004587

SMS

# Post-Cal Stability Checks

0.05g/210L 0.047 to 0.053	0.08g/210L 0.077 to 0.083	0.20g/210L 0.194 to 0.206	DGS 0.08g/210L 0.077 to 0.083																																																																								
<p>#2</p> <p>OSCEOLA COUNTY 50 Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-004587 01/27/2026 Software: 8100.27</p> <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.010</td><td>14:44</td></tr> <tr><td>Control Test</td><td>0.052</td><td>14:44</td></tr> <tr><td>Air Blank</td><td>0.010</td><td>14:45</td></tr> <tr><td>Control Test</td><td>0.052</td><td>14:46</td></tr> <tr><td>Air Blank</td><td>0.010</td><td>14:46</td></tr> <tr><td>Control Test</td><td>0.052</td><td>14:47</td></tr> <tr><td>Air Blank</td><td>0.010</td><td>14:48</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0520</td><td></td></tr> <tr><td>Std Dev</td><td>0.0000</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr> </tbody> </table> <p>Operator's Signature: <i>[Signature]</i></p>	Test	g/210L	Time	Air Blank	0.010	14:44	Control Test	0.052	14:44	Air Blank	0.010	14:45	Control Test	0.052	14:46	Air Blank	0.010	14:46	Control Test	0.052	14:47	Air Blank	0.010	14:48	Control Test Stats			Average	0.0520		Std Dev	0.0000		Rel Std Dev(%)	0.0000		<p>OSCEOLA COUNTY 50 Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-004587 01/27/2026 Software: 8100.27</p> <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>14:50</td></tr> <tr><td>Control Test</td><td>0.083</td><td>14:51</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:51</td></tr> <tr><td>Control Test</td><td>0.083</td><td>14:52</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:53</td></tr> <tr><td>Control Test</td><td>INT*</td><td>14:53</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:54</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0827</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.6984</td><td></td></tr> </tbody> </table> <p>*Interferent Detect</p> <p>Operator's Signature: <i>[Signature]</i></p>	Test	g/210L	Time	Air Blank	0.000	14:50	Control Test	0.083	14:51	Air Blank	0.000	14:51	Control Test	0.083	14:52	Air Blank	0.000	14:53	Control Test	INT*	14:53	Air Blank	0.000	14:54	Control Test Stats			Average	0.0827		Std Dev	0.0006		Rel Std Dev(%)	0.6984		<p>ceased post stability test</p>	<p>and optical bench calibrations SMS</p>
Test	g/210L	Time																																																																									
Air Blank	0.010	14:44																																																																									
Control Test	0.052	14:44																																																																									
Air Blank	0.010	14:45																																																																									
Control Test	0.052	14:46																																																																									
Air Blank	0.010	14:46																																																																									
Control Test	0.052	14:47																																																																									
Air Blank	0.010	14:48																																																																									
Control Test Stats																																																																											
Average	0.0520																																																																										
Std Dev	0.0000																																																																										
Rel Std Dev(%)	0.0000																																																																										
Test	g/210L	Time																																																																									
Air Blank	0.000	14:50																																																																									
Control Test	0.083	14:51																																																																									
Air Blank	0.000	14:51																																																																									
Control Test	0.083	14:52																																																																									
Air Blank	0.000	14:53																																																																									
Control Test	INT*	14:53																																																																									
Air Blank	0.000	14:54																																																																									
Control Test Stats																																																																											
Average	0.0827																																																																										
Std Dev	0.0006																																																																										
Rel Std Dev(%)	0.6984																																																																										


#1

- post calibration adjustment  
Stability  
SUA

OSCEOLA COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-004587  
01/27/2026  
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	14:24
Control Test	0.054	14:25
Air Blank	0.000	14:25
Control Test	0.052	14:26
Air Blank	0.000	14:27
Control Test	0.051	14:27
Air Blank	0.000	14:28
Control Test Stats		
Average	0.0523	
Std Dev	0.0015	
Rel Std Dev(%)	2.9188	

Changed 0-mg and placed  
SIM <sup>top</sup> cap <sub>SUA</sub> off and back  
on tight to rerun.  
SUA

  
Operator's Signature

1/27/2026 SAT  
2nd optical  
Bench adjustment

0.04 g/210L  
solution had  
faults detected  
on 1st and 2nd  
adjustment.

Prior to rerun  
I replaced  
SIM top to  
S/N MP 6289.

Fault occurred  
again.

OSCEOLA COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-004587  
01/27/2026 16:09:30

Auto Calibration  
Max Power Res Value = 87  
Auto Range Res Value = 59

Sol Value = 0.000 g/210L \*\*\*  
Fit value = 0.0000 mg/l %%%  
Samples Taken = 4, Discarded = 1  
Sum Io = 12661, Sum Io = 13224

<<<< CHANNEL 1 >>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 0.2610 (-0.0560)  
Sample #2 = 0.1460 (0.0850)  
Sample #3 = 0.1280 (0.0530)  
Sample #4 = 0.1110 (0.1040)  
Avg % Abs = 0.1283 (0.0807)  
STD DEV = 0.0175 (0.0258)  
REL STD DEV = 13.638 (31.952)

<<<< CHANNEL 2 >>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 0.3290 (-0.0190)  
Sample #2 = 0.2800 (0.0320)  
Sample #3 = 0.2390 (0.0630)  
Sample #4 = 0.2840 (0.0410)  
Avg % Abs = 0.2677 (0.0453)  
STD DEV = 0.0249 (0.0159)  
REL STD DEV = 9.305 (35.179)

Sol Value = 0.040 g/210L \*\*\*  
Fit value = 0.1905 mg/l %%%  
Samples Taken = 4, Discarded = 1  
Sum Io = 12666, Sum Io = 13218

<<<< CHANNEL 1 >>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 0.8770 (-0.0120)  
Sample #2 = 0.9130 (-0.0900)  
Sample #3 = 0.8010 (0.0400)  
Sample #4 = 0.7600 (0.0650)  
Avg % Abs = 0.8247 (0.0050)  
STD DEV = 0.0792 (0.0832)  
REL STD DEV = 9.604 (1664.333)

<<<< CHANNEL 2 >>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 1.6800 (-0.0040)  
Sample #2 = 1.5800 (0.0630)  
Sample #3 = 1.6090 (0.0620)  
Sample #4 = 1.5820 (0.0820)  
Avg % Abs = 1.5930 (0.0690)  
STD DEV = 0.0142 (0.0113)  
REL STD DEV = 0.890 (16.333)

Sol Value = 0.040 g/210L \*\*\*  
Fit value = 0.1905 mg/l %%%  
Samples Taken = 4, Discarded = 1  
Sum Io = 12677, Sum Io = 13214

<<<< CHANNEL 1 >>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 0.7380 (0.0390)  
Sample #2 = 0.9160 (0.0380)  
Sample #3 = 0.8570 (0.0260)  
Sample #4 = 0.8410 (0.0830)  
Avg % Abs = 0.8713 (0.0490)  
STD DEV = 0.0395 (0.0300)  
REL STD DEV = 4.533 (61.326)

<<<< CHANNEL 2 >>>>  
Sample % Abs (% Abs Ref)  
Sample #1 = 1.6670 (0.0070)  
Sample #2 = 1.6220 (0.0370)  
Sample #3 = 1.6290 (0.0350)  
Sample #4 = 1.6100 (0.0410)  
Avg % Abs = 1.6203 (0.0377)  
STD DEV = 0.0096 (0.0031)  
REL STD DEV = 0.593 (8.111)

Fault detected  
2x for 0.04 g/210L  
-Failed. SAT

Return Material Authorization

**Ship to:**

CMI, Inc.

Enforcement Electronics

Shipment to repair facility authorized by: Daniel Lyons on 2/2/2026

Items Returned:  Instrument  Supplies  Other Describe: \_\_\_\_\_

Instrument Model: Intoxilyzer 8000 Serial Number: 80-004587

Bill To Address:  Osceola County Sheriffs Office 2601 East Irlo Bronson Memorial Hwy Kissimmee, FL 34744	Ship to Address:  FDLE Offsite Mail Facility c/o FDLE HQ 813 B Lake Bradford Rd Tallahassee, FL 32304
--	--

Reason for Return: Optical bench calibration (x2) unsuccessful.
--

I require an estimate <b><u>BEFORE</u></b> any repairs will be authorized and/or conducted Please contact: <u>Daniel Lyons</u> Phone #: <u>407-799-1711</u> Email: <u>Daniel.lyons@osceolasheriff.org</u> ATP Contact Name: <u>LeAndra Higginbotham</u> ATP Email: <u>leandrahigginbotham@fdle.state.fl.us</u>
---