

Alcohol Testing Program - Instrument Processing Sheet

Agency: Charlotte County Sheriff's Office Instrument Serial Number: 80-001363
 Date In: 4/22/2026 DI Completion Date: N/A ☐ Ship ☐ P/U ☐ H/D ☒ CMI ☐ EE

Intake By: <u>TDG</u> Date: <u>4/23/2026</u>		Quality Checks By: <u>TDG</u> Date: <u>4/23/2026</u>		Flow Adjustment By: <u>TDG</u>																
☒ Annual ☐ Dropped Off ☐ Registration ☐ Return from CMI / EE ☐ Training Instrument Visual Inspection ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/Accessories ☐ Power Cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes: Agency reports receiving Ambient and Purge Fails.		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value: <u>200</u> ☒ Flow Verification (L/s) Flow Column #: <u>ATP106</u> 32 mm 0.140 (.139-.169) 36 mm 0.156 (.156-.190) 53 mm 0.222* (.228-.278) 103 mm 0.496 (.447-.547) ☒ Barometric Pressure Check Gauge ID #: <u>26932</u> Gauge: <u>1018</u> Instrument: <u>1017</u> ☒ Stability Checks <table border="1" style="width: 100%; border-collapse: collapse; margin-top: 5px;"> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot#/Exp</th> </tr> <tr> <td>0.050</td> <td>MP6286</td> <td>202406K 6/19/2026</td> </tr> <tr> <td>0.080</td> <td>MP6287</td> <td>202406L 6/19/2026</td> </tr> <tr> <td>0.200</td> <td>MP6288</td> <td>202406N 6/20/2026</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>AG526603 9/23/2027</td> </tr> </table>		Simulator	Serial #	Lot#/Exp	0.050	MP6286	202406K 6/19/2026	0.080	MP6287	202406L 6/19/2026	0.200	MP6288	202406N 6/20/2026	0.080 DGS	N/A	AG526603 9/23/2027	Flow Column #: <u>ATP101 (x2)</u> ☒ 5L/min – 17mm ☒ 15L/min – 53mm ☒ 30L/min – 103mm ☒ R-Value: <u>200 / 200</u> ☒ Post Adjustment Verification (L/S) Flow Column #: <u>ATP106 (x2)</u> 32 mm 0.132* / 0.128* (.139-.169) 36 mm 0.148* / 0.144* (.156-.190) 53 mm 0.230 / 0.222* (.228-.278) 103 mm 0.503 / 0.500 (.447-.547)	
Simulator	Serial #	Lot#/Exp																		
0.050	MP6286	202406K 6/19/2026																		
0.080	MP6287	202406L 6/19/2026																		
0.200	MP6288	202406N 6/20/2026																		
0.080 DGS	N/A	AG526603 9/23/2027																		
Optical Bench Adjustment By: _____		Department Inspection By: _____																		
Barometric Pressure Gauge: _____ ID#: _____		Barometric Pressure ID#: _____																		
Simulator	Serial #	Lot #	Expiration	Gauge: _____ Instrument: _____																
0.000		N/A	N/A	Mouth Alcohol Solution Lot #: _____ Exp: _____																
0.040				Acetone Stock Solution Lot #: _____ Exp: _____																
0.100				Simulator	Serial Number															
0.200				0.000																
0.300				Interferent																
0.080 DGS	N/A			0.050																
☐ Post Optical Bench Adjustment Stability Checks				0.080																
Simulator	Serial #	Lot #	Expiration	0.200																
0.050				Attachments																
0.080				☐ Form 41 ☒ Stability Checks ☐ Calibration Certificate ☐ Optical Bench Adjustment ☐ Post-Stability Checks ☒ Flow Adjustment ☐ Form 40 ☒ Other: Return Material Authorization																
0.200																				
0.080 DGS	N/A																			
Gauge ID #: _____ Gauge: _____ Instrument: _____																				
Notes/Suggested Service: *Outside nominal. (TDG 4/23/26) Could not get flow adjustment verification to pass. Instrument had a flow adjustment in December 2025. Sending to CMI. Compliance with Chapter 11D-8, FAC was not determined, but the instrument passed Stability Checks. Could not replicate the exception messages described. (TDG 5/15/26)				☐ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☐ Return to/Place into Evidentiary Use ☒ Remain Out of Evidentiary Use ☐ Conduct an Agency Inspection Before Evidentiary Use																
LeAndra Higginbotham Digitally signed by LeAndra Higginbotham Date: 2026.05.18 16:41:15 -04'00'				Kaitlyn Spearin Digitally signed by Kaitlyn Spearin Date: 2026.05.19 06:41:34 -04'00'																
Tech Review				Admin Review																

Flow Calibration Adjustment(s)

Performed by TDG

#1

CHARLOTTE COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001363
04/23/2026
Software: 8100.27

Flow Rate Calibration*****

1: Rate (Liters/min) = 5
SQRT(Diff) = 7.613
2: Rate (Liters/min) = 15
SQRT(Diff) = 11.914
3: Rate (Liters/min) = 30
SQRT(Diff) = 20.590
Dependent Data Scale Factor = 100000 L/min
Independent Data Scale Factor = 256
Rounded Slope = 741
Rounded Intercept = -871334
Correlation = 0.99701

#2

CHARLOTTE COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001363
04/23/2026
Software: 8100.27

Flow Rate Calibration*****

1: Rate (Liters/min) = 5
SQRT(Diff) = 5.383
2: Rate (Liters/min) = 15
SQRT(Diff) = 10.484
3: Rate (Liters/min) = 30
SQRT(Diff) = 19.898
Dependent Data Scale Factor = 100000 L/min
Independent Data Scale Factor = 256
Rounded Slope = 666
Rounded Intercept = -367437
Correlation = 0.99849

Stability Checks

0.05g/210L 0.047 to 0.053	0.08g/210L 0.077 to 0.083	0.20g/210L 0.194 to 0.206	DGS 0.08g/210L 0.077 to 0.083 ≤ 0.003 of Wet																																																																																																																																																												
CHARLOTTE COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001363 04/23/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>12:46</td></tr> <tr><td>Control Test</td><td>0.048</td><td>12:46</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:47</td></tr> <tr><td>Control Test</td><td>0.048</td><td>12:48</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:48</td></tr> <tr><td>Control Test</td><td>0.048</td><td>12:49</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:49</td></tr> <tr><td>Control Test</td><td>0.048</td><td>12:49</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0480</td><td></td></tr> <tr><td>Std Dev</td><td>0.0000</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr> </tbody> </table> Operator's Signature 	Test	g/210L	Time	Air Blank	0.000	12:46	Control Test	0.048	12:46	Air Blank	0.000	12:47	Control Test	0.048	12:48	Air Blank	0.000	12:48	Control Test	0.048	12:49	Air Blank	0.000	12:49	Control Test	0.048	12:49	Control Test Stats			Average	0.0480		Std Dev	0.0000		Rel Std Dev(%)	0.0000		CHARLOTTE COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001363 04/23/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>12:53</td></tr> <tr><td>Control Test</td><td>0.078</td><td>12:53</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:54</td></tr> <tr><td>Control Test</td><td>0.079</td><td>12:55</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:55</td></tr> <tr><td>Control Test</td><td>0.079</td><td>12:56</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:56</td></tr> <tr><td>Control Test</td><td>0.079</td><td>12:56</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0787</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.7339</td><td></td></tr> </tbody> </table> Operator's Signature 	Test	g/210L	Time	Air Blank	0.000	12:53	Control Test	0.078	12:53	Air Blank	0.000	12:54	Control Test	0.079	12:55	Air Blank	0.000	12:55	Control Test	0.079	12:56	Air Blank	0.000	12:56	Control Test	0.079	12:56	Control Test Stats			Average	0.0787		Std Dev	0.0006		Rel Std Dev(%)	0.7339		CHARLOTTE COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001363 04/23/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>13:15</td></tr> <tr><td>Control Test</td><td>0.200</td><td>13:16</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:16</td></tr> <tr><td>Control Test</td><td>0.199</td><td>13:17</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:17</td></tr> <tr><td>Control Test</td><td>0.199</td><td>13:18</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:19</td></tr> <tr><td>Control Test</td><td>0.199</td><td>13:19</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.1993</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.2896</td><td></td></tr> </tbody> </table> Operator's Signature 	Test	g/210L	Time	Air Blank	0.000	13:15	Control Test	0.200	13:16	Air Blank	0.000	13:16	Control Test	0.199	13:17	Air Blank	0.000	13:17	Control Test	0.199	13:18	Air Blank	0.000	13:19	Control Test	0.199	13:19	Control Test Stats			Average	0.1993		Std Dev	0.0006		Rel Std Dev(%)	0.2896		CHARLOTTE COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001363 04/23/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>13:06</td></tr> <tr><td>Control Test</td><td>0.079</td><td>13:06</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:07</td></tr> <tr><td>Control Test</td><td>0.079</td><td>13:07</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:08</td></tr> <tr><td>Control Test</td><td>0.078</td><td>13:08</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:08</td></tr> <tr><td>Control Test</td><td>0.078</td><td>13:08</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0787</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.7339</td><td></td></tr> </tbody> </table> Operator's Signature 	Test	g/210L	Time	Air Blank	0.000	13:06	Control Test	0.079	13:06	Air Blank	0.000	13:07	Control Test	0.079	13:07	Air Blank	0.000	13:08	Control Test	0.078	13:08	Air Blank	0.000	13:08	Control Test	0.078	13:08	Control Test Stats			Average	0.0787		Std Dev	0.0006		Rel Std Dev(%)	0.7339	
Test	g/210L	Time																																																																																																																																																													
Air Blank	0.000	12:46																																																																																																																																																													
Control Test	0.048	12:46																																																																																																																																																													
Air Blank	0.000	12:47																																																																																																																																																													
Control Test	0.048	12:48																																																																																																																																																													
Air Blank	0.000	12:48																																																																																																																																																													
Control Test	0.048	12:49																																																																																																																																																													
Air Blank	0.000	12:49																																																																																																																																																													
Control Test	0.048	12:49																																																																																																																																																													
Control Test Stats																																																																																																																																																															
Average	0.0480																																																																																																																																																														
Std Dev	0.0000																																																																																																																																																														
Rel Std Dev(%)	0.0000																																																																																																																																																														
Test	g/210L	Time																																																																																																																																																													
Air Blank	0.000	12:53																																																																																																																																																													
Control Test	0.078	12:53																																																																																																																																																													
Air Blank	0.000	12:54																																																																																																																																																													
Control Test	0.079	12:55																																																																																																																																																													
Air Blank	0.000	12:55																																																																																																																																																													
Control Test	0.079	12:56																																																																																																																																																													
Air Blank	0.000	12:56																																																																																																																																																													
Control Test	0.079	12:56																																																																																																																																																													
Control Test Stats																																																																																																																																																															
Average	0.0787																																																																																																																																																														
Std Dev	0.0006																																																																																																																																																														
Rel Std Dev(%)	0.7339																																																																																																																																																														
Test	g/210L	Time																																																																																																																																																													
Air Blank	0.000	13:15																																																																																																																																																													
Control Test	0.200	13:16																																																																																																																																																													
Air Blank	0.000	13:16																																																																																																																																																													
Control Test	0.199	13:17																																																																																																																																																													
Air Blank	0.000	13:17																																																																																																																																																													
Control Test	0.199	13:18																																																																																																																																																													
Air Blank	0.000	13:19																																																																																																																																																													
Control Test	0.199	13:19																																																																																																																																																													
Control Test Stats																																																																																																																																																															
Average	0.1993																																																																																																																																																														
Std Dev	0.0006																																																																																																																																																														
Rel Std Dev(%)	0.2896																																																																																																																																																														
Test	g/210L	Time																																																																																																																																																													
Air Blank	0.000	13:06																																																																																																																																																													
Control Test	0.079	13:06																																																																																																																																																													
Air Blank	0.000	13:07																																																																																																																																																													
Control Test	0.079	13:07																																																																																																																																																													
Air Blank	0.000	13:08																																																																																																																																																													
Control Test	0.078	13:08																																																																																																																																																													
Air Blank	0.000	13:08																																																																																																																																																													
Control Test	0.078	13:08																																																																																																																																																													
Control Test Stats																																																																																																																																																															
Average	0.0787																																																																																																																																																														
Std Dev	0.0006																																																																																																																																																														
Rel Std Dev(%)	0.7339																																																																																																																																																														

Return Material Authorization

Ship to:☒ CMI, Inc.☐ Enforcement ElectronicsShipment to repair facility authorized by: Vanessa Chandler on 5/12/2026Items Returned: ☒ Instrument ☐ Supplies ☐ Other Describe: _____Instrument Model: Intoxilyzer 8000 Serial Number: 80-001363

Bill To Address:

Charlotte County Sheriff's Office
Attn: Vanessa Chandler

Ship to Address:

Florida Department of Law Enforcement
Fort Myers Regional Operations Center
Attn: Taylor Gutschow
4700 Terminal Drive, Suite 1
Fort Myers, FL 33907

Reason for Return:

Instrument will not pass post flow adjustment verification. Agency reports the instrument gives Ambient and Purge Fail messages.

I require an estimate **BEFORE** any repairs will be authorized and/or conductedPlease contact: Vanessa Chandler Phone #: 941-833-6446Email: vchandler@ccsofl.netATP Contact Name: Taylor Gutschow ATP Email: TaylorGutschow@fdle.state.fl.us