

RECEIVED

# Florida Department of Law Enforcement Alcohol Testing Program

JUN 20 2016

FDLE  
Alcohol Testing Program

## DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: DELAND P.D.

Time of Inspection: 08:23

Date of Inspection: 06/17/2016

Serial Number: 80-001114

Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:201507A Exp: 07/14/2017	0.08g/210L Test (g/210L) Lot#:201601F Exp: 01/26/2018	0.20g/210L Test (g/210L) Lot#:201505A Exp: 05/12/2017	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG605301 Exp: 02/22/2018
0.000	0.046	0.080	0.201	0.079
0.000	0.046	0.081	0.202	0.077
0.000	0.046	0.081	0.203	0.077
0.000	0.046	0.081	0.204	0.077
0.000	0.045	0.081	0.203	0.077
0.000	0.046	0.080	0.203	0.077
0.000	0.046	0.081	0.203	0.077
0.000	0.046	0.081	0.203	0.077
0.000	0.046	0.081	0.202	0.077
0.000	0.046	0.081	0.203	0.077

Standard Deviations	0.0003	0.0004	0.0008	0.0006
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0005 Number of Simulators Used: 5

Remarks:

The above instrument complies (  ) does not comply (  ) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

JAKE L SHANAHAN

Signature and Printed Name

06/17/2016

Date

80-001114  
 Stability Checks  
 6/17/16

OSM  
 SK

DELAND P.O.  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-001114  
 06/17/2016  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	06:44
Control Test	0.078	06:45
Air Blank	0.000	06:45
Control Test	0.078	06:45
Air Blank	0.000	06:46
Control Test	0.078	06:46
Air Blank	0.000	06:47
Control Test Stats		
Average	0.0780	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

DGS

Operator's Signature

DELAND P.O.  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-001114  
 06/17/2016  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	06:53
Control Test	0.047	06:54
Air Blank	0.000	06:54
Control Test	0.048	06:55
Air Blank	0.000	06:55
Control Test	0.048	06:56
Air Blank	0.000	06:56
Control Test Stats		
Average	0.0477	
Std Dev	0.0006	
Rel Std Dev(%)	1.2112	

Operator's Signature

DELAND P.O.  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-001114  
 06/17/2016  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	06:48
Control Test	0.078	06:49
Air Blank	0.000	06:50
Control Test	0.078	06:50
Air Blank	0.000	06:51
Control Test	0.079	06:51
Air Blank	0.000	06:52
Control Test Stats		
Average	0.0783	
Std Dev	0.0006	
Rel Std Dev(%)	0.7370	

Operator's Signature

DELAND P.O.  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-001114  
 06/17/2016  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	06:58
Control Test	0.198	06:59
Air Blank	0.000	06:59
Control Test	0.199	07:00
Air Blank	0.000	07:00
Control Test	0.200	07:01
Air Blank	0.000	07:02
Control Test Stats		
Average	0.1990	
Std Dev	0.0010	
Rel Std Dev(%)	0.5025	

Operator's Signature