

Stability Tests - Jefferson CSO # 80-000956 9/4/16

JEFFERSON COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000956
09/06/2016
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:41
Control Test	0.048	11:42
Air Blank	0.000	11:42
Control Test	0.048	11:43
Air Blank	0.000	11:44
Control Test	0.048	11:44
Air Blank	0.000	11:45
Control Test Stats		
Average	0.0480	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

ppm

JEFFERSON COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000956
09/06/2016
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:36
Control Test	0.078	11:37
Air Blank	0.000	11:38
Control Test	0.077	11:38
Air Blank	0.000	11:39
Control Test	0.078	11:39
Air Blank	0.000	11:40
Control Test Stats		
Average	0.0777	
Std Dev	0.0006	
Rel Std Dev(%)	0.7434	

JEFFERSON COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000956
09/06/2016
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:38
Control Test	0.196	12:39
Air Blank	0.000	12:39
Control Test	0.195	12:40
Air Blank	0.000	12:41
Control Test	0.195	12:41
Air Blank	0.000	12:42
Control Test Stats		
Average	0.1953	
Std Dev	0.0006	
Rel Std Dev(%)	0.2956	

JEFFERSON COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000956
09/06/2016
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:33
Control Test	0.081	12:34
Air Blank	0.000	12:34
Control Test	0.080	12:35
Air Blank	0.000	12:35
Control Test	0.079	12:35
Air Blank	0.000	12:36
Control Test Stats		
Average	0.0800	
Std Dev	0.0010	
Rel Std Dev(%)	1.2500	

DGS

TSK

[Signature]

Operator's Signature

[Signature]

Operator's Signature

[Signature]

Operator's Signature

[Signature]

Operator's Signature