



Alcohol Testing Program

INSTRUMENT PROCESSING SHEET

Agency MARTIN County SO S/N 80-001189

Date In 9/13/16 Date Out 9/21/16 Ship P/U H/D CMI EE

Intake	Performed By <u>SP</u>	Quality Checks	Performed By <u>SP</u>	Flow Calibration	Performed By <u>SP</u>																																											
<input type="checkbox"/> Registration <input checked="" type="checkbox"/> Annual <input type="checkbox"/> Return from CMI <input type="checkbox"/> Return from Enforcement Electronics <input type="checkbox"/> Other _____ Visual Inspection: <u>OK</u> Case <u>OK</u> Handle <u>OK</u> Dry Gas Holder <u>OK</u> Feet <u>OK</u> Keyboard/Plug <u>OK</u> Back/Plugs <u>OK</u> Screws tight <u>OK</u> Breath Hose Other Equipment: <input type="checkbox"/> Power cord <input type="checkbox"/> Printer Cable <input type="checkbox"/> Other _____ Notes: _____ _____ _____		<input checked="" type="checkbox"/> Breath Tube Screen <input checked="" type="checkbox"/> Replace O-Rings <input checked="" type="checkbox"/> Instrument Set Up Verified <input checked="" type="checkbox"/> R-Value <u>179</u> <input checked="" type="checkbox"/> Flow Verification (L/s) Flow Column # <u>ATP102</u> 32mm <u>1104</u> (.139 - .169) 36mm <u>175</u> (.156 - .190) 53mm <u>240</u> (.228 - .278) 103mm <u>515</u> (.447 - .547) <input checked="" type="checkbox"/> Barometric Pressure Check Gauge ID # <u>28427</u> <input checked="" type="checkbox"/> Stability Checks		<input checked="" type="checkbox"/> Flow Calibration N/A <input type="checkbox"/> Flow Calibration Complete Flow Column # _____ <input type="checkbox"/> 5L/min - 17mm <input type="checkbox"/> 15L/min - 50mm <input type="checkbox"/> 30L/min - 103mm <input type="checkbox"/> R-Value _____ <input type="checkbox"/> Post Calibration Verification (L/s) Flow Column # _____ 32mm _____ (.139 - .169) 36mm _____ (.156 - .190) 53mm _____ (.228 - .278) 103mm _____ (.447 - .547)																																												
		<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td>0.05</td> <td>SD1018</td> <td>201507A 7-14-17</td> </tr> <tr> <td>0.08</td> <td>SD1011</td> <td>201601F 1-26-18</td> </tr> <tr> <td>0.20</td> <td>SD1025</td> <td>201604C 4-5-18</td> </tr> <tr> <td>0.08 DGS</td> <td>N/A</td> <td>AG612405 5-3-18</td> </tr> </tbody> </table>	Simulator	Serial #	Lot #/Exp	0.05	SD1018	201507A 7-14-17	0.08	SD1011	201601F 1-26-18	0.20	SD1025	201604C 4-5-18	0.08 DGS	N/A	AG612405 5-3-18		<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td>G4444</td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td>G2882</td> <td>16101</td> <td>2-2-18</td> </tr> <tr> <td>0.100</td> <td>G2078</td> <td>16001</td> <td>5-8-18</td> </tr> <tr> <td>0.200</td> <td>G2408</td> <td>16103</td> <td>6-14-18</td> </tr> <tr> <td>0.400</td> <td>G5358</td> <td>16102</td> <td>3-22-18</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>15615080A2</td> <td>7-5-17</td> </tr> </tbody> </table>	Simulator	Serial Number	Lot Number	Expiration	0.000	G4444	N/A	N/A	0.040	G2882	16101	2-2-18	0.100	G2078	16001	5-8-18	0.200	G2408	16103	6-14-18	0.400	G5358	16102	3-22-18	0.080 DGS	N/A	15615080A2	7-5-17	
Simulator	Serial #	Lot #/Exp																																														
0.05	SD1018	201507A 7-14-17																																														
0.08	SD1011	201601F 1-26-18																																														
0.20	SD1025	201604C 4-5-18																																														
0.08 DGS	N/A	AG612405 5-3-18																																														
Simulator	Serial Number	Lot Number	Expiration																																													
0.000	G4444	N/A	N/A																																													
0.040	G2882	16101	2-2-18																																													
0.100	G2078	16001	5-8-18																																													
0.200	G2408	16103	6-14-18																																													
0.400	G5358	16102	3-22-18																																													
0.080 DGS	N/A	15615080A2	7-5-17																																													
				<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.00</td> <td>SD1021 198P</td> </tr> <tr> <td>Interferent</td> <td>SD1021</td> </tr> <tr> <td>0.05</td> <td>SD1018</td> </tr> <tr> <td>0.08</td> <td>SD1011</td> </tr> <tr> <td>0.20</td> <td>SD1025</td> </tr> </tbody> </table>	Simulator	Serial Number	0.00	SD1021 198P	Interferent	SD1021	0.05	SD1018	0.08	SD1011	0.20	SD1025																																
Simulator	Serial Number																																															
0.00	SD1021 198P																																															
Interferent	SD1021																																															
0.05	SD1018																																															
0.08	SD1011																																															
0.20	SD1025																																															

RECEIVED
SEP 23 2016
FDLE
Alcohol Testing Program

Optical Bench Calibration	Performed By <u>SP</u>	Department Inspection	Performed By <u>SP</u>																				
<input type="checkbox"/> Optical Bench Calibration N/A <input checked="" type="checkbox"/> Optical Bench Calibration Complete Barometric Pressure Gauge <u>1015</u> ID # <u>26932</u>		<input checked="" type="checkbox"/> Barometric Pressure <u>1017</u> Gauge ID# <u>28427</u> <u>1014</u> Instrument Mouth Alcohol Solution Lot # <u>2015-A</u> Acetone Stock Solution Lot # <u>2016-B</u>																					
<input checked="" type="checkbox"/> Post Calibration Stability Checks		<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.05</td> <td>SD1018</td> <td>201507A</td> <td>7-14-17</td> </tr> <tr> <td>0.08</td> <td>SD1011</td> <td>201601F</td> <td>1-26-18</td> </tr> <tr> <td>0.20</td> <td>SD1025</td> <td>201604C</td> <td>4-5-18</td> </tr> <tr> <td>0.08 DGS</td> <td>N/A</td> <td>AG612405</td> <td>5-3-18</td> </tr> </tbody> </table>	Simulator	Serial Number	Lot Number	Expiration	0.05	SD1018	201507A	7-14-17	0.08	SD1011	201601F	1-26-18	0.20	SD1025	201604C	4-5-18	0.08 DGS	N/A	AG612405	5-3-18	
Simulator	Serial Number	Lot Number	Expiration																				
0.05	SD1018	201507A	7-14-17																				
0.08	SD1011	201601F	1-26-18																				
0.20	SD1025	201604C	4-5-18																				
0.08 DGS	N/A	AG612405	5-3-18																				
		<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.00</td> <td>SD1021 198P</td> </tr> <tr> <td>Interferent</td> <td>SD1021</td> </tr> <tr> <td>0.05</td> <td>SD1018</td> </tr> <tr> <td>0.08</td> <td>SD1011</td> </tr> <tr> <td>0.20</td> <td>SD1025</td> </tr> </tbody> </table>	Simulator	Serial Number	0.00	SD1021 198P	Interferent	SD1021	0.05	SD1018	0.08	SD1011	0.20	SD1025									
Simulator	Serial Number																						
0.00	SD1021 198P																						
Interferent	SD1021																						
0.05	SD1018																						
0.08	SD1011																						
0.20	SD1025																						

Notes: OPTICAL BENCH PERFORMED TO BRING VALUES CLOSER TO NOMINAL GAUGE: 1016 INSTRUMENT: 1003 SP

DC - KWB

- Instrument Complies with Chapter 11D-8, FAC
- Instrument Does Not Comply with Chapter 11D-8, FAC
- Return to/Place into Evidentiary Use
- Remain Out of Evidentiary Use
- Conduct an Agency Inspection Before Evidentiary Use

Brett Hubbard

Quality Control Review

9/23/16

Date

STABILITY CHECKS - INSTRUMENT # 80-001189 - MARTIN COUNTY SO
 9/20/16 SP

DES

MARTIN COUNTY SO
 Intoxilizer - Alcohol Analyzer
 Model 8000 SN 80-001189
 09/20/2016
 Software: 8100.27

MARTIN COUNTY SO
 Intoxilizer - Alcohol Analyzer
 Model 8000 SN 80-001189
 09/20/2016
 Software: 8100.27

MARTIN COUNTY SO
 Intoxilizer - Alcohol Analyzer
 Model 8000 SN 80-001189
 09/20/2016
 Software: 8100.27

MARTIN COUNTY SO
 Intoxilizer - Alcohol Analyzer
 Model 8000 SN 80-001189
 09/20/2016
 Software: 8100.27

MARTIN COUNTY SO
 Intoxilizer - Alcohol Analyzer
 Model 8000 SN 80-001189
 09/20/2016
 Software: 8100.27

MARTIN COUNTY SO
 Intoxilizer - Alcohol Analyzer
 Model 8000 SN 80-001189
 09/20/2016
 Software: 8100.27

MARTIN COUNTY SO
 Intoxilizer - Alcohol Analyzer
 Model 8000 SN 80-001189
 09/20/2016
 Software: 8100.27

Test	9/21/16	Time
Air Blank	0.050	10:15
Control Test	0.049	10:16
Air Blank	0.050	10:16
Control Test	0.049	10:17
Air Blank	0.050	10:17
Control Test	0.050	10:18
Air Blank	0.050	10:19
Control Test Stats		
Average	0.0433	
Std Dev	0.0065	
Rel Std Dev(%)	1.1703	

Test	9/21/16	Time
Air Blank	0.080	10:21
Control Test	0.079	10:22
Air Blank	0.080	10:22
Control Test	0.080	10:23
Air Blank	0.080	10:23
Control Test	0.079	10:24
Air Blank	0.080	10:25
Control Test Stats		
Average	0.0793	
Std Dev	0.0086	
Rel Std Dev(%)	0.7277	

Test	9/21/16	Time
Air Blank	0.000	10:28
Control Test	0.197	10:28
Air Blank	0.000	10:29
Control Test	0.197	10:30
Air Blank	0.000	10:30
Control Test	0.197	10:31
Air Blank	0.000	10:32
Control Test Stats		
Average	0.1970	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

Test	9/21/16	Time
Air Blank	0.000	10:34
Control Test	0.080	10:34
Air Blank	0.000	10:34
Control Test	0.080	10:35
Air Blank	0.000	10:35
Control Test	0.081	10:36
Air Blank	0.000	10:36
Control Test Stats		
Average	0.0803	
Std Dev	0.0006	
Rel Std Dev(%)	0.7187	

Test	9/21/16	Time
Air Blank	0.000	10:28
Control Test	0.197	10:28
Air Blank	0.000	10:29
Control Test	0.197	10:30
Air Blank	0.000	10:30
Control Test	0.197	10:31
Air Blank	0.000	10:32
Control Test Stats		
Average	0.1970	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

Test	9/21/16	Time
Air Blank	0.000	10:28
Control Test	0.197	10:28
Air Blank	0.000	10:29
Control Test	0.197	10:30
Air Blank	0.000	10:30
Control Test	0.197	10:31
Air Blank	0.000	10:32
Control Test Stats		
Average	0.1970	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

Test	9/21/16	Time
Air Blank	0.000	10:28
Control Test	0.197	10:28
Air Blank	0.000	10:29
Control Test	0.197	10:30
Air Blank	0.000	10:30
Control Test	0.197	10:31
Air Blank	0.000	10:32
Control Test Stats		
Average	0.1970	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

Test	9/21/16	Time
Air Blank	0.000	10:28
Control Test	0.197	10:28
Air Blank	0.000	10:29
Control Test	0.197	10:30
Air Blank	0.000	10:30
Control Test	0.197	10:31
Air Blank	0.000	10:32
Control Test Stats		
Average	0.1970	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

Test	9/21/16	Time
Air Blank	0.000	10:34
Control Test	0.080	10:34
Air Blank	0.000	10:34
Control Test	0.080	10:35
Air Blank	0.000	10:35
Control Test	0.081	10:36
Air Blank	0.000	10:36
Control Test Stats		
Average	0.0803	
Std Dev	0.0006	
Rel Std Dev(%)	0.7187	

SP

Operator's Signature

DES

SP

Operator's Signature

SP

Operator's Signature

SP

Operator's Signature

SK

MARTIN COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8100 SN 80-001189
19/21/2016 11:31:01

Auto Calibration
Max Power Res value = 23
Auto Range Res Value = 7

Sol Value = 0.000 g/210L ***
Fit value = 0.0000 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum lo = 12350, Sum hi = 14030
<<<< CHANNEL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.1970 (-0.0030)
Sample #2 = 0.1670 (0.0590)
Sample #3 = 0.1840 (0.1070)
Sample #4 = 0.1270 (0.1720)
Avg % Abs = 0.1593 (0.1127)
STD DEV = 0.0293 (0.0567)
REL STD DEV = 18.366 (50.337)

<<<< CHANNEL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.1420 (-0.0070)
Sample #2 = 0.1410 (0.0010)
Sample #3 = 0.1310 (0.0270)
Sample #4 = 0.1310 (0.0340)
Avg % Abs = 0.1343 (0.0207)
STD DEV = 0.0058 (0.0174)
REL STD DEV = 4.298 (84.134)

Sol Value = 0.040 g/210L ***
Fit value = 0.1905 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum lo = 12323, Sum hi = 14023
<<<< CHANNEL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.8990 (-0.0130)
Sample #2 = 0.9010 (0.0290)
Sample #3 = 0.8660 (0.0360)
Sample #4 = 0.8430 (0.0820)
Avg % Abs = 0.8700 (0.0490)
STD DEV = 0.0292 (0.0288)
REL STD DEV = 3.357 (58.760)

<<<< CHANNEL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 1.6050 (-0.0080)
Sample #2 = 1.5710 (0.0400)
Sample #3 = 1.6110 (0.0120)
Sample #4 = 1.5860 (0.0360)
Avg % Abs = 1.5893 (0.0293)
STD DEV = 0.0202 (0.0151)
REL STD DEV = 1.271 (51.626)

Sol Value = 0.100 g/210L ***
Fit value = 0.4762 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum lo = 12312, Sum hi = 14018

<<<< CHANNEL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 1.9980 (-0.0090)
Sample #2 = 1.9870 (0.0350)
Sample #3 = 1.9920 (0.0240)
Sample #4 = 2.0020 (0.0360)
Avg % Abs = 1.9937 (0.0317)
STD DEV = 0.0076 (0.0067)
REL STD DEV = 0.383 (21.026)

<<<< CHANNEL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 3.7910 (-0.0070)
Sample #2 = 3.7750 (0.0340)
Sample #3 = 3.7850 (0.0350)
Sample #4 = 3.7870 (0.0360)
Avg % Abs = 3.7823 (0.0350)
STD DEV = 0.0064 (0.0010)
REL STD DEV = 0.170 (2.857)

Sol Value = 0.200 g/210L ***
Fit value = 0.9524 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum lo = 12307, Sum hi = 14014

<<<< CHANNEL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 3.8110 (-0.0090)
Sample #2 = 3.8030 (0.0560)
Sample #3 = 3.8220 (0.0530)
Sample #4 = 3.7890 (0.0910)
Avg % Abs = 3.8047 (0.0667)
STD DEV = 0.0166 (0.0211)
REL STD DEV = 0.435 (31.690)

<<<< CHANNEL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 7.2730 (-0.0240)
Sample #2 = 7.2210 (0.0930)
Sample #3 = 7.2370 (0.0880)
Sample #4 = 7.2330 (0.1040)
Avg % Abs = 7.2303 (0.0950)
STD DEV = 0.0083 (0.0082)
REL STD DEV = 0.115 (8.616)

Sol Value = 0.400 g/210L ***
Fit value = 1.9048 mg/l %%%
Samples Taken = 4, Discarded = 1
Sum lo = 12294, Sum hi = 14005

<<<< CHANNEL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 7.1650 (-0.0220)
Sample #2 = 7.1390 (0.0810)
Sample #3 = 7.1350 (0.0850)
Sample #4 = 7.1150 (0.1100)
Avg % Abs = 7.1297 (0.0920)
STD DEV = 0.0129 (0.0157)
REL STD DEV = 0.180 (17.083)

OPTICAL BENCH
CALIBRATION
INSTRUMENT #
80-001189

9/21/16

MARTIN COUNTY
SO

SP

10

OPTICAL BENCH CALIBRATION - INSTRUMENT # 80-001189-
 MARTIN COUNTY SO 9/21/16
 SP

<<<< CHANNEL 2 >>>>
 Sample % Abs (% Abs Ref)
 Sample #1 = 13.4910 (-0.0110)
 Sample #2 = 13.4150 (0.1379)
 Sample #3 = 13.4030 (0.1493)
 Sample #4 = 13.3950 (0.1710)
 Avg % Abs = 13.4043 (0.1523)
 STD DEV = 0.0101 (0.0172)
 REL STD DEV = 0.075 (1.319)

Solution Stats Quadratic Fit Chan 1		
Act	Fit	Residual
g/210L	g/210L	g/210L
0.000	0.001	-0.0012
0.040	0.039	0.0013
0.100	0.099	0.0005
0.200	0.201	-0.0010
0.400	0.400	0.0002

***** AUTO CAL DATA *****
 <<<< CHANNEL 1 >>>>
 Sol Val = 0.0000 mg/l or 0.000 g/210L
 % Abs = 0.159
 Std Dev = 0.03 Rel Std Dev = 18.37
 Sol Val = 0.1905 mg/l or 0.040 g/210L
 % Abs = 0.870
 Std Dev = 0.03 Rel Std Dev = 3.36
 Sol Val = 0.4762 mg/l or 0.100 g/210L
 % Abs = 1.994
 Std Dev = 0.01 Rel Std Dev = 0.38
 Sol Val = 0.9524 mg/l or 0.200 g/210L
 % Abs = 3.805
 Std Dev = 0.02 Rel Std Dev = 0.44
 Sol Val = 1.9048 mg/l or 0.400 g/210L
 % Abs = 7.130
 Std Dev = 0.01 Rel Std Dev = 0.18
 Zero Order Coef = -337.92
 First Order Coef = 2474.86
 Second Order Coef = 34.04
 Standard Deviation = 51.468342

Solution Stats Quadratic Fit Chan 2		
Act	Fit	Residual
g/210L	g/210L	g/210L
0.000	0.001	-0.0007
0.040	0.039	0.0008
0.100	0.100	0.0004
0.200	0.201	-0.0006
0.400	0.400	0.0002

Sol Value = 0.000 g/210L ***
 Fit value = 0.0010 mg/l ****
 Samples Taken = 4, Discarded = 1
 ***** CHANNEL 1 *****
 Sample #1 = 2902.00
 Sample #2 = 2844.00
 Sample #3 = 2837.00
 Sample #4 = 2776.00
 Average Result = 2819.0000
 STD DEV = 37.4032
 REL STD DEV = 1.327

<<<< CHANNEL 2 >>>>
 Sol Val = 0.0000 mg/l or 0.000 g/210L
 % Abs = 0.134
 Std Dev = 0.01 Rel Std Dev = 4.30
 Sol Val = 0.1905 mg/l or 0.040 g/210L
 % Abs = 1.589
 Std Dev = 0.02 Rel Std Dev = 1.27
 Sol Val = 0.4762 mg/l or 0.100 g/210L
 % Abs = 3.782
 Std Dev = 0.01 Rel Std Dev = 0.17
 Sol Val = 0.9524 mg/l or 0.200 g/210L
 % Abs = 7.230
 Std Dev = 0.01 Rel Std Dev = 0.12
 Sol Val = 1.9048 mg/l or 0.400 g/210L
 % Abs = 13.404
 Std Dev = 0.01 Rel Std Dev = 0.08
 Zero Order Coef = -131.70
 First Order Coef = 1233.33
 Second Order Coef = 14.69
 Standard Deviation = 31.640590

***** CHANNEL 2 *****
 Sample #1 = 3345.00
 Sample #2 = 3354.00
 Sample #3 = 3338.00
 Sample #4 = 3359.00
 Average Result = 3351.3333
 STD DEV = 10.9697
 REL STD DEV = 0.327

 Dry Gas H2O Adjust Results *****
 Barometric Pressure = 1014
 3 um H2O Adjust (mg/l*10,000) = 990
 9 um H2O Adjust (mg/l*10,000) = 459
 **** AUTO CAL PASS

BK

WAB

POST CALIBRATION STABILITY CHECKS - INSTRUMENT # 80-001189 - MARTIN COUNTY SO

9/21/16

SP

DOS

MARTIN COUNTY SO
Intoxilyzer - Alcotest Analyzer
Model 8000
SN 80-001189
9/21/2016
Software: 8100.27

Test	9/21/16	Time
Air Blank	0.000	11:08
Control Test	0.049	11:09
Air Blank	0.000	11:09
Control Test	0.050	11:10
Air Blank	0.000	11:11
Control Test	0.050	11:11
Air Blank	0.000	11:12
Control Test Stats		
Average	0.0497	
Std Dev	0.0005	
Rel Std Dev(%)	1.1625	

SP

Operator's Signature

DOS

MARTIN COUNTY SO
Intoxilyzer - Alcotest Analyzer
Model 8000
SN 80-001185
9/21/2016
Software: 8100.27

Test	9/21/16	Time
Air Blank	0.000	11:13
Control Test	0.079	11:14
Air Blank	0.000	11:14
Control Test	0.080	11:15
Air Blank	0.000	11:16
Control Test	0.079	11:16
Air Blank	0.000	11:17
Control Test Stats		
Average	0.0793	
Std Dev	0.0006	
Rel Std Dev(%)	0.7277	

SP

Operator's Signature

MARTIN COUNTY SO
Intoxilyzer - Alcotest Analyzer
Model 8000
SN 80-001155
9/21/2016
Software: 8100.27

Test	9/21/16	Time
Air Blank	0.000	11:19
Control Test	0.199	11:19
Air Blank	0.000	11:20
Control Test	0.199	11:21
Air Blank	0.000	11:21
Control Test	0.198	11:22
Air Blank	0.000	11:23
Control Test Stats		
Average	0.1987	
Std Dev	0.0006	
Rel Std Dev(%)	0.2905	

SP

Operator's Signature

MARTIN COUNTY SO
Intoxilyzer - Alcotest Analyzer
Model 8000
SN 80-001185
9/21/2016
Software: 8100.27

Test	9/21/16	Time
Air Blank	0.000	11:23
Control Test	0.091	11:24
Air Blank	0.000	11:24
Control Test	0.091	11:24
Air Blank	0.000	11:25
Control Test	0.091	11:25
Air Blank	0.000	11:26
Control Test Stats		
Average	0.0910	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

SP

Operator's Signature

ASK