

STABILITY CHECKS - INSTRUMENT # 80-005146 - PATRICK AFB - 12/14/16 SP

PATRICK AFB
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005146
12/14/2016
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:12
Control Test	0.050	10:13
Air Blank	0.000	10:13
Control Test	0.049	10:14
Air Blank	0.000	10:14
Control Test	0.049	10:15
Air Blank	0.000	10:16
Control Test Stats		
Average	0.0493	
Std Dev	0.0006	
Rel. Std Dev(%)	1.1703	

SP

Operator's Signature

PATRICK AFB
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005146
12/14/2016
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:17
Control Test	0.081	10:17
Air Blank	0.000	10:18
Control Test	0.082	10:19
Air Blank	0.000	10:19
Control Test	0.081	10:20
Air Blank	0.000	10:20
Control Test Stats		
Average	0.0813	
Std Dev	0.0006	
Rel. Std Dev(%)	0.7099	

SP

Operator's Signature

PATRICK AFB
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005146
12/14/2016
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:06
Control Test	0.201	10:07
Air Blank	0.000	10:07
Control Test	0.201	10:08
Air Blank	0.000	10:08
Control Test	0.202	10:09
Air Blank	0.000	10:10
Control Test Stats		
Average	0.2013	
Std Dev	0.0006	
Rel. Std Dev(%)	0.2868	

SP

Operator's Signature

PATRICK AFB
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005146
12/14/2016
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:02
Control Test	0.079	10:03
Air Blank	0.000	10:03
Control Test	0.079	10:03
Air Blank	0.000	10:04
Control Test	0.079	10:04
Air Blank	0.000	10:05
Control Test Stats		
Average	0.0790	
Std Dev	0.0000	
Rel. Std Dev(%)	0.0000	

DGS

13K

SP

Operator's Signature

ADS