



Alcohol Testing Program

INSTRUMENT PROCESSING SHEET

Agency

FWC - Pensacola

S/N

80-000802

Date In

5/2/17

Date Out

5/30/17

☒ Ship☐ P/U☐ H/D☐ CMI☐ EE

Intake

Performed By

JP

- ☐ Registration
☒ Annual
☐ Return from CMI
☒ Return from Enforcement
 Electronics
☐ Other

Visual Inspection:

- ☒ Case ☒ Handle
☒ Dry Gas Holder ☒ Feet
☒ Keyboard/Plug ☒ Back/Plugs
☒ Screws tight ☒ Breath Hose

Other Equipment:

- ☐ Power cord
☐ Printer Cable
☒ Other State Bag

Notes:

Quality Checks

Performed By

JP

- ☒ Lab Temp °C 20.6
☒ Breath Tube Screen
☒ Replace External O-Rings
☒ Instrument Set Up Verified
☒ R-Value 213
☒ Flow Verification (L/s)
 Flow Column # ATP103
 32mm 160 (.139 - .169)
 36mm 179 (.156 - .190)
 53mm 232 (.228 - .278)
 103mm 503 (.447 - .547)

☒ Barometric Pressure CheckGauge ID # 26932☒ Stability Checks

Simulator	Serial #	Lot #/Exp
0.05	<u>DR2035</u>	<u>201603D</u> <u>3/8/18</u>
0.08	<u>SD3964</u>	<u>201611B</u> <u>11/15/18</u>
0.20	<u>G11621</u>	<u>201604C</u> <u>4/5/18</u>
0.08 DGS	N/A	<u>AG626604</u> <u>9/22/18</u>

Flow Calibration

Performed By

- ☒ Flow Calibration N/A
☐ Flow Calibration Complete
 Flow Column # ATP103
☐ 5L/min - 17mm
☐ 15L/min - 53mm
☐ 30L/min - 103mm
☐ R-Value
☐ Post Calibration Verification (L/s)
 Flow Column #
 32mm (.139 - .169)
 36mm (.156 - .190)
 53mm (.228 - .278)
 103mm (.447 - .547)

Maintenance

Performed By

- ☐ Battery Replacement
☐ Dry Gas Regulator Replacement
☐ Breath Tube Replacement
☐ Other

Quality Checks Cont. Performed By

Simulator Temperatures °C

External Digital Therm. ID#: 300505

- ☒ 34°C +- .2 Serial #: DR2035
☒ 34°C +- .2 Serial #: SD3964
☒ 34°C +- .2 Serial #: G11621

Calibration Adjustment

Performed By

- ☒ Calibration Adjustment N/A
☐ Calibration Adjustment Complete

Barometric Pressure Gauge ID #

Simulator	Serial Number	Lot Number	Expiration
0.000		N/A	N/A
0.040			
0.100			
0.200			
0.300			
0.080 DGS	N/A		

☐ Post Calibration Adjustment Stability Checks

Simulator	Serial Number	Lot Number	Expiration
0.05			
0.08			
0.20			
0.08 DGS	N/A		

Department Inspection

Performed By

- ☒ Barometric Pressure 1018 Gauge
 ID# 26932 1016 Instrument

Mouth Alcohol Solution Lot # 2016-AAcetone Stock Solution Lot # 2016-B

Simulator	Serial Number
0.00	<u>G2879</u>
Interferent	<u>G8144</u>
0.05	<u>DR2035</u>
0.08	<u>SD3964</u>
0.20	<u>G11621</u>

Attachments

- ☒ Form 41
☒ Pre-Stability Tests
☐ Flow Calibration
☐ Calibration Adjustment
☐ Post-Stability Tests
☐ Other

Notes/Suggested Service:

QA/DC/DB☒ Instrument Complies with Chapter 11D-8, FAC☐ Instrument Does Not Comply with Chapter 11D-8, FAC☒ Return to/Place into Evidentiary Use☐ Remain Out of Evidentiary Use☒ Conduct an Agency Inspection Before Evidentiary Use

Quality Control Review

Date



Alcohol Testing Program

INSTRUMENT PROCESSING SHEET

Agency FWC - PensacolaS/N 80-000892Date In 5/2/17

Date Out _____

☐ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake	Performed By <u>DP</u>	Quality Checks	Performed By _____	Flow Calibration	Performed By _____															
☐ Registration ☒ Annual ☐ Return from CMI ☒ Return from Enforcement Electronics ☐ Other _____		☐ Breath Tube Screen ☐ Replace O-Rings ☐ Instrument Set Up Verified ☐ R-Value _____ ☐ Flow Verification (L/s) Flow Column # _____ 32mm _____ (.139 - .169) 36mm _____ (.156 - .190) 53mm _____ (.228 - .278) 103mm _____ (.447 - .547)		☐ Flow Calibration N/A ☐ Flow Calibration Complete Flow Column # _____ ☐ 5L/min - 17mm ☐ 15L/min - 53mm ☐ 30L/min - 103mm ☐ R-Value _____ ☐ Post Calibration Verification (s) Flow Column # _____ 32mm _____ (.139 - .169) 36mm _____ (.156 - .190) 53mm _____ (.228 - .278) 103mm _____ (.447 - .547)																
Visual Inspection: ☒ Case ☒ Handle ☒ Dry Gas Holder ☒ Feet ☒ Keyboard/Plug ☒ Back/Plugs ☒ Screws tight ☒ Breath Hose		☐ Barometric Pressure Check Gauge ID # _____ ☐ Stability Checks																		
Other Equipment: ☐ Power cord ☐ Printer Cable ☒ Other <u>Static Bag</u>		<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td>0.05</td> <td></td> <td></td> </tr> <tr> <td>0.08</td> <td></td> <td></td> </tr> <tr> <td>0.20</td> <td></td> <td></td> </tr> <tr> <td>0.08 DGS</td> <td>N/A</td> <td></td> </tr> </tbody> </table>	Simulator	Serial #	Lot #/Exp	0.05			0.08			0.20			0.08 DGS	N/A			☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other _____	
Simulator	Serial #	Lot #/Exp																		
0.05																				
0.08																				
0.20																				
0.08 DGS	N/A																			
Notes: _____				Suggested Service _____ _____ _____																

Optical Bench Calibration				Performed By _____
☐ Optical Bench Calibration N/A ☐ Optical Bench Calibration Complete				
Barometric Pressure Gauge _____ ID # _____				
Simulator	Serial Number	Lot Number	Expiration	
0.000		N/A	N/A	
0.040				
0.100				
0.200				
0.400				
0.080 DGS	N/A			
☐ Post Calibration Stability Checks				
Simulator	Serial Number	Lot Number	Expiration	
0.05				
0.08				
0.20				
0.08 DGS	N/A			

Notes: QA/QC & RMB

Department Inspection		Performed By _____
☐ Barometric Pressure ID# _____ Gauge _____ Instrument _____		
Mouth Alcohol Solution Lot # _____		
Acetone Stock Solution Lot # _____		
Simulator	Serial Number	
0.00		
Interferent		
0.05		
0.08		
0.20		
Attachments ☐ Form 41 ☐ Pre-Stability Tests ☐ Flow Calibration ☐ Optical Bench Cal ☐ Post-Stability Tests ☐ Other _____		

- ☐ Instrument Complies with Chapter 11D-8, FAC
☐ Instrument Does Not Comply with Chapter 11D-8, FAC
☐ Return to/Place into Evidentiary Use
☐ Remain Out of Evidentiary Use
☐ Conduct an Agency Inspection Before Evidentiary Use

Quality Control Review

Date

Revised November 2012

FWC
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000892
05/30/2017
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:58
Control Test	0.049	10:59
Air Blank	0.000	11:00
Control Test	0.049	11:00
Air Blank	0.000	11:01
Control Test	0.048	11:02
Air Blank	0.000	11:02
Control Test Stats		
Average	0.0487	
Std Dev	0.0006	
Rel Std Dev(%)	1.1863	

P Murphy
Operator's Signature

FWC
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000892
05/30/2017
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:13
Control Test	0.080	11:14
Air Blank	0.000	11:15
Control Test	0.079	11:15
Air Blank	0.000	11:16
Control Test	0.080	11:16
Air Blank	0.000	11:17
Control Test Stats		
Average	0.0797	
Std Dev	0.0006	
Rel Std Dev(%)	0.7247	

P Murphy
Operator's Signature

FWC
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000892
05/30/2017
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:08
Control Test	0.201	11:08
Air Blank	0.000	11:09
Control Test	0.200	11:10
Air Blank	0.000	11:10
Control Test	0.200	11:11
Air Blank	0.000	11:12
Control Test Stats		
Average	0.2003	
Std Dev	0.0006	
Rel Std Dev(%)	0.2882	

P Murphy
Operator's Signature

FWC
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000892
05/30/2017
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:19
Control Test	0.078	11:19
Air Blank	0.000	11:20
Control Test	0.079	11:20
Air Blank	0.000	11:21
Control Test	0.079	11:21
Air Blank	0.000	11:22
Control Test Stats		
Average	0.0787	
Std Dev	0.0006	
Rel Std Dev(%)	0.7339	

DGS

P Murphy
Operator's Signature

DMS

30



Alcohol Testing Program

INSTRUMENT PROCESSING SHEET

Agency FWC-PensacolaS/N 80-000892Date In 3/23/17Date Out 3/29/17☐ Ship ☐ P/U ☐ H/D ☐ CMI ☒ EE

Intake ☐ Registration ☒ Annual ☐ Return from CMI ☒ Return from Enforcement Electronics ☐ Other _____ Visual Inspection: <u>OK</u> Case <u>OK</u> Handle <u>OK</u> Dry Gas Holder <u>OK</u> Feet <u>OK</u> Keyboard/Plug <u>OK</u> Back/Plugs <u>OK</u> Screws tight <u>OK</u> Breath Hose Other Equipment: ☒ Power cord ☐ Printer Cable ☒ Other <u>Static Bag</u> Notes: _____ _____ _____ _____	Quality Checks Performed By <u>PSM</u> ☒ Breath Tube Screen ☒ Replace O-Rings ☒ Instrument Set Up Verified ☒ R-Value <u>216</u> ☒ Flow Verification (L/s) Flow Column # <u>ATP102</u> 32mm <u>.160</u> (.139 - .169) 36mm <u>.175</u> (.156 - .190) 53mm <u>.234</u> (.228 - .278) 103mm <u>.495</u> (.447 - .547) ☒ Barometric Pressure Check Gauge ID # <u>28421</u> ☒ Stability Checks <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td>0.05</td> <td><u>DR2035</u></td> <td><u>201603D</u> <u>3/8/18</u></td> </tr> <tr> <td>0.08</td> <td><u>SD3964</u></td> <td><u>201601F</u> <u>1/26/18</u></td> </tr> <tr> <td>0.20</td> <td><u>SD1025</u></td> <td><u>201604C</u> <u>4/5/18</u></td> </tr> <tr> <td>0.08 DGS</td> <td>N/A</td> <td><u>AG626604</u> <u>9/22/18</u></td> </tr> </tbody> </table>	Simulator	Serial #	Lot #/Exp	0.05	<u>DR2035</u>	<u>201603D</u> <u>3/8/18</u>	0.08	<u>SD3964</u>	<u>201601F</u> <u>1/26/18</u>	0.20	<u>SD1025</u>	<u>201604C</u> <u>4/5/18</u>	0.08 DGS	N/A	<u>AG626604</u> <u>9/22/18</u>	Flow Calibration Performed By <u>PSM</u> ☒ Flow Calibration N/A ☐ Flow Calibration Complete Flow Column # _____ ☐ 5L/min - 17mm ☐ 15L/min - 53mm ☐ 30L/min - 103mm ☐ R-Value _____ ☐ Post Calibration Verification (L/s) Flow Column # _____ 32mm _____ (.139 - .169) 36mm _____ (.156 - .190) 53mm _____ (.228 - .278) 103mm _____ (.447 - .547)
Simulator	Serial #	Lot #/Exp															
0.05	<u>DR2035</u>	<u>201603D</u> <u>3/8/18</u>															
0.08	<u>SD3964</u>	<u>201601F</u> <u>1/26/18</u>															
0.20	<u>SD1025</u>	<u>201604C</u> <u>4/5/18</u>															
0.08 DGS	N/A	<u>AG626604</u> <u>9/22/18</u>															
		Maintenance Performed By _____ ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other _____ Suggested Service _____ _____ _____															

Optical Bench Calibration Performed By _____
☒ Optical Bench Calibration N/A
☐ Optical Bench Calibration Complete
 Barometric Pressure Gauge _____ ID # _____

Simulator	Serial Number	Lot Number	Expiration
0.000		N/A	N/A
0.040			
0.100			
0.200			
0.400			
0.080 DGS	N/A		

☐ Post Calibration Stability Checks

Simulator	Serial Number	Lot Number	Expiration
0.05			
0.08			
0.20			
0.08 DGS	N/A		

Notes: SENDING BACK TO ENFORCEMENT
ELECTRONICS FOR OPTICAL BENCH REPAIR /
CALIBRATION ADJUST
4/4/17
Scott Kirkland

Quality Control Review

Department Inspection

Performed By PSM

☒ Barometric Pressure 1014 Gauge
 ID# 28421 1013 Instrument

Mouth Alcohol Solution Lot # 2016-AAcetone Stock Solution Lot # 2016-B

Simulator	Serial Number
0.00	<u>SD1019</u>
Interferent	<u>SD1021</u>
0.05	<u>DR2035</u>
0.08	<u>SD3964</u>
0.20	<u>SD1025</u>

Attachments

- | | |
|---|---|
| ☒ Form 41 | ☐ Optical Bench Cal |
| ☒ Pre-Stability Tests | ☐ Post-Stability Tests |
| ☐ Flow Calibration | ☐ Other _____ |

- ☒ Instrument Complies with Chapter 11D-8, FAC
☐ Instrument Does Not Comply with Chapter 11D-8, FAC
☐ Return to/Place into Evidentiary Use
☒ Remain Out of Evidentiary Use
☐ Conduct an Agency Inspection Before Evidentiary Use

Date 4/4/17

FWC
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000892
03/29/2017
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:41
Control Test	0.049	08:42
Air Blank	0.000	08:42
Control Test	0.049	08:43
Air Blank	0.000	08:44
Control Test	0.049	08:44
Air Blank	0.000	08:45
Control Test Stats		
Average	0.0490	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

P Murphy
Operator's Signature

FWC
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000892
03/29/2017
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:36
Control Test	0.079	08:37
Air Blank	0.000	08:37
Control Test	0.079	08:38
Air Blank	0.000	08:39
Control Test	0.081	08:39
Air Blank	0.000	08:40
Control Test Stats		
Average	0.0797	
Std Dev	0.0012	
Rel Std Dev(%)	1.4494	

P Murphy
Operator's Signature

FWC
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000892
03/29/2017
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:51
Control Test	0.199	08:51
Air Blank	0.000	08:52
Control Test	0.200	08:53
Air Blank	0.000	08:53
Control Test	0.198	08:54
Air Blank	0.000	08:54
Control Test Stats		
Average	0.1990	
Std Dev	0.0010	
Rel Std Dev(%)	0.5025	

P Murphy
Operator's Signature

FWC
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000892
03/29/2017
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:56
Control Test	0.077	08:56
Air Blank	0.000	08:57
Control Test	0.077	08:57
Air Blank	0.000	08:58
Control Test	0.077	08:58
Air Blank	0.000	08:58
Control Test Stats		
Average	0.0770	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

DGS

P Murphy
Operator's Signature

BK



Alcohol Testing Program

INSTRUMENT PROCESSING SHEET

Agency FWC-PensacolaS/N 80-000892Date In 1/11/17Date Out 1/26/17☐ Ship ☐ P/U ☐ H/D ☐ CMI ☒ EE

Intake Performed By <u>DP</u> ☐ Registration ☒ Annual ☐ Return from CMI ☒ Return from Enforcement Electronics ☐ Other _____ Visual Inspection: <u>OK</u> Case <u>OK</u> Handle <u>OK</u> Dry Gas Holder <u>OK</u> Feet <u>OK</u> Keyboard/Plug <u>OK</u> Back/Plugs <u>OK</u> Screws tight <u>OK</u> Breath Hose Other Equipment: ☒ Power cord ☐ Printer Cable ☒ Other <u>Static Bag</u> Notes: _____ _____ _____ _____	Quality Checks Performed By _____ ☒ Breath Tube Screen ☒ Replace O-Rings ☒ Instrument Set Up Verified ☒ R-Value <u>212</u> ☐ Flow Verification (L/s) Flow Column # _____ 32mm _____ (.139 - .169) 36mm _____ (.156 - .190) 53mm _____ (.228 - .278) 103mm _____ (.447 - .547) ☐ Barometric Pressure Check Gauge ID # _____ ☐ Stability Checks <table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td>0.05</td> <td></td> <td></td> </tr> <tr> <td>0.08</td> <td></td> <td></td> </tr> <tr> <td>0.20</td> <td></td> <td></td> </tr> <tr> <td>0.08 DGS</td> <td>N/A</td> <td></td> </tr> </tbody> </table>	Simulator	Serial #	Lot #/Exp	0.05			0.08			0.20			0.08 DGS	N/A		Flow Calibration Performed By _____ ☐ Flow Calibration N/A ☐ Flow Calibration Complete Flow Column # _____ ☐ 5L/min - 17mm ☐ 15L/min - 53mm ☐ 30L/min - 103mm ☐ R-Value _____ ☐ Post Calibration Verification (L/s) Flow Column # _____ 32mm _____ (.139 - .169) 36mm _____ (.156 - .190) 53mm _____ (.228 - .278) 103mm _____ (.447 - .547) Maintenance Performed By _____ ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other _____ Suggested Service _____ _____ _____
Simulator	Serial #	Lot #/Exp															
0.05																	
0.08																	
0.20																	
0.08 DGS	N/A																

Optical Bench Calibration Performed By _____ ☐ Optical Bench Calibration N/A ☐ Optical Bench Calibration Complete Barometric Pressure Gauge _____ ID # _____ <table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td></td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.100</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.400</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table> ☐ Post Calibration Stability Checks <table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.05</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.08</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.20</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.08 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table>	Simulator	Serial Number	Lot Number	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.400				0.080 DGS	N/A			Simulator	Serial Number	Lot Number	Expiration	0.05				0.08				0.20				0.08 DGS	N/A			Department Inspection Performed By _____ ☐ Barometric Pressure _____ Gauge ID# _____ Instrument Mouth Alcohol Solution Lot # _____ Acetone Stock Solution Lot # _____ <table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.00</td> <td></td> </tr> <tr> <td>Interferent</td> <td></td> </tr> <tr> <td>0.05</td> <td></td> </tr> <tr> <td>0.08</td> <td></td> </tr> <tr> <td>0.20</td> <td></td> </tr> </tbody> </table> Attachments ☐ Form 41 ☐ Pre-Stability Tests ☐ Flow Calibration ☐ Optical Bench Cal ☐ Post-Stability Tests ☐ Other _____	Simulator	Serial Number	0.00		Interferent		0.05		0.08		0.20	
Simulator	Serial Number	Lot Number	Expiration																																																										
0.000		N/A	N/A																																																										
0.040																																																													
0.100																																																													
0.200																																																													
0.400																																																													
0.080 DGS	N/A																																																												
Simulator	Serial Number	Lot Number	Expiration																																																										
0.05																																																													
0.08																																																													
0.20																																																													
0.08 DGS	N/A																																																												
Simulator	Serial Number																																																												
0.00																																																													
Interferent																																																													
0.05																																																													
0.08																																																													
0.20																																																													

 Notes: Flow measurement NOT FUNCTIONAL
RETURNING FOR WARRANTY REPAIR

- ☐
- Instrument Complies with Chapter 11D-8, FAC
-
- ☐
- Instrument Does Not Comply with Chapter 11D-8, FAC
-
- ☐
- Return to/Place into Evidentiary Use
-
- ☒
- Remain Out of Evidentiary Use
-
- ☐
- Conduct an Agency Inspection Before Evidentiary Use

Quality Control Review

Date

Revised November 2012