



<b>TYPE OF TEST</b>	<b>SERIAL NUMBER</b>	<b>AGENCY</b>	<b>DATE</b>	<b>PERFORMED BY</b>
Stabilities	80-001056	Davie Police Department	8/1/2017	<i>DELL</i>

<b>0.05g/210L</b>	<b>0.08g/210L</b>	<b>0.20g/210L</b>	<b>DGS 0.08g/210L</b>
<b>SN: SD3967 Temp: 34.05</b>	<b>SN: SD3968 Temp: 34.08c</b>	<b>SN: SD3969 Temp: 34.96c</b>	<b>Lot AG626605</b>
<b>0.047 to 0.053</b> <input checked="" type="checkbox"/>	<b>0.077 to 0.083</b> <input checked="" type="checkbox"/>	<b>0.194 to 0.206</b> <input checked="" type="checkbox"/>	<b>0.077 to 0.083</b> <input checked="" type="checkbox"/>

DAVIE PD Intoxilyzer - Alcohol Analyzer Model 8000 08/01/2017 Software: 8100.27	DAVIE PD Intoxilyzer - Alcohol Analyzer Model 8000 08/01/2017 Software: 8100.27	DAVIE PD Intoxilyzer - Alcohol Analyzer Model 8000 08/01/2017 Software: 8100.27	DAVIE PD Intoxilyzer - Alcohol Analyzer Model 8000 08/01/2017 Software: 8100.27
Test	Test	Test	Test
Air Blank 0.000 06:03 Control Test 0.050 06:04 Air Blank 0.000 06:04 Control Test 0.051 06:05 Air Blank 0.000 06:06 Control Test 0.051 06:06 Air Blank 0.000 06:07 Control Test Stats Average 0.0507 Std Dev 0.0006 Rel Std Dev(%) 1.1395	Air Blank 0.000 06:09 Control Test 0.081 06:10 Air Blank 0.000 06:10 Control Test 0.080 06:11 Air Blank 0.000 06:12 Control Test 0.081 06:12 Air Blank 0.000 06:13 Control Test Stats Average 0.0807 Std Dev 0.0006 Rel Std Dev(%) 0.7157	Air Blank 0.000 06:14 Control Test 0.198 06:15 Air Blank 0.000 06:15 Control Test 0.196 06:16 Air Blank 0.000 06:17 Control Test 0.196 06:17 Air Blank 0.000 06:18 Control Test Stats Average 0.1967 Std Dev 0.0012 Rel Std Dev(%) 0.5871	Air Blank 0.000 06:19 Control Test 0.000 06:19 Air Blank 0.000 06:20 Control Test 0.081 06:20 Air Blank 0.000 06:21 Control Test 0.082 06:21 Air Blank 0.000 06:22 Control Test Stats Average 0.0810 Std Dev 0.0010 Rel Std Dev(%) 1.2346
Time	Time	Time	Time
----- 06:03 06:04 06:04 06:05 06:06 06:06 06:07 -----	----- 06:09 06:10 06:10 06:11 06:12 06:12 06:13 -----	----- 06:14 06:15 06:15 06:16 06:17 06:17 06:18 -----	----- 06:19 06:19 06:20 06:20 06:21 06:21 06:22 -----
Operator's Signature: <i>DELL</i>	Operator's Signature: <i>DELL</i>	Operator's Signature: <i>DELL</i>	Operator's Signature: <i>DELL</i>

*8/1/17*  
*DELL*

*DELL*