



Alcohol Testing Program

INSTRUMENT PROCESSING SHEET

Agency

Fruitland Park

S/N

80-001350

Date In

5/3/17

Date Out

5/4/17

☒ Ship☐ P/U☐ H/D☐ CMI☐ EE

Intake

Performed By

RDS

Quality Checks

Performed By

RDS

Flow Calibration

Performed By

RDS

- ☐ Registration
- ☒ Annual
- ☐ Return from CMI
- ☒ Return from Enforcement
- Electronics
- ☐ Other

Visual Inspection:

OK Case OK Handle

OK Dry Gas Holder OK Feet

OK Keyboard/Plug OK Back/Plugs

OK Screws tight OK Breath Hose

Other Equipment:

- ☐ Power cord
- ☐ Printer Cable
- ☒ Other Static Bag

Notes:

- ☒ Breath Tube Screen
- ☒ Replace O-Rings
- ☒ Instrument Set Up Verified
- ☒ R-Value 161
- ☒ Flow Verification (L/s)
- Flow Column # ATP103
- 32mm 0.156 (.139 - .169)
- 36mm 0.175 (.156 - .190)
- 53mm 0.238 (.228 - .278)
- 103mm 0.500 (.447 - .547)

☒ Barometric Pressure Check

Gauge ID # 26932

Stability Checks

Simulator	Serial #	Lot #/Exp
0.05	DR2035	201603D 3/8/18
0.08	SD3964	201611B 11/5/18
0.20	G11621	201604C 4/5/18
0.08 DGS	N/A	AG626605 9/22/18

- ☒ Flow Calibration N/A
- ☐ Flow Calibration Complete
- Flow Column #
- ☐ 5L/min - 17mm
- ☐ 15L/min - 53mm
- ☐ 30L/min - 103mm

☐ R-Value

☐ Post Calibration Verification (L/s)

Flow Column #

32mm (.139 - .169)

36mm (.156 - .190)

53mm (.228 - .278)

103mm (.447 - .547)

Maintenance

Performed By

- ☐ Battery Replacement
- ☐ Dry Gas Regulator Replacement
- ☐ Breath Tube Replacement
- ☐ Other

Suggested Service

Optical Bench Calibration

Performed By

- ☒ Optical Bench Calibration N/A
- ☐ Optical Bench Calibration Complete

Barometric Pressure Gauge ID #

Simulator	Serial Number	Lot Number	Expiration
0.000		N/A	N/A
0.040			
0.100			
0.200			
0.400			
0.080 DGS	N/A		

Post Calibration Stability Checks

Simulator	Serial Number	Lot Number	Expiration
0.05			
0.08			
0.20			
0.08 DGS	N/A		

Notes: Please update menu level 2 password to something unique. RDS

QA/QC OK SP

Department Inspection

Performed By

RDS

- ☒ Barometric Pressure
- ID# 28421
- 1007 Gauge
- 1007 Instrument

Mouth Alcohol Solution Lot # 2016-C

Acetone Stock Solution Lot # 2017-A

Simulator	Serial Number
0.00	G2879
Interferent	G8144
0.05	DR2035
0.08	SD3964
0.20	G11621

Attachments

- ☒ Form 41
- ☒ Pre-Stability Tests
- ☐ Flow Calibration
- ☐ Optical Bench Cal
- ☐ Post-Stability Tests
- ☐ Other

- ☒ Instrument Complies with Chapter 11D-8, FAC
- ☐ Instrument Does Not Comply with Chapter 11D-8, FAC
- ☒ Return to/Place into Evidentiary Use
- ☐ Remain Out of Evidentiary Use
- ☒ Conduct an Agency Inspection Before Evidentiary Use

Quality Control Review

Date

Revised November 2012

INSTRUMENT PROCESSING SHEET

Agency Fruitland Park PD

S/N 80-001350

Date In 2/9/17

Date Out 2/14/17

☐ Ship ☐ P/U ☐ H/D ☐ CMI ☒ EE

Intake ☐ Registration ☒ Annual ☐ Return from CMI ☐ Return from Enforcement Electronics ☐ Other _____ Visual Inspection: ☒ Case ☒ Handle ☒ Dry Gas Holder ☒ Feet ☒ Keyboard/Plug ☒ Back/Plugs ☒ Screws tight ☒ Breath Hose Other Equipment: ☐ Power cord ☐ Printer Cable ☒ Other <u>Static Bag</u> Notes: _____ _____ _____ _____	Quality Checks ☐ Breath Tube Screen ☐ Replace O-Rings ☐ Instrument Set Up Verified ☐ R-Value _____ ☐ Flow Verification (L/s) Flow Column # _____ 32mm _____ (.139 - .169) 36mm _____ (.156 - .190) 53mm _____ (.228 - .278) 103mm _____ (.447 - .547) ☐ Barometric Pressure Check Gauge ID # _____ ☐ Stability Checks <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td>0.05</td> <td></td> <td></td> </tr> <tr> <td>0.08</td> <td></td> <td></td> </tr> <tr> <td>0.20</td> <td></td> <td></td> </tr> <tr> <td>0.08 DGS</td> <td>N/A</td> <td></td> </tr> </tbody> </table>	Simulator	Serial #	Lot #/Exp	0.05			0.08			0.20			0.08 DGS	N/A		Flow Calibration ☐ Flow Calibration N/A ☐ Flow Calibration Complete Flow Column # _____ ☐ 5L/min - 4.7mm ☐ 15L/min - 53mm ☐ 30L/min - 103mm ☐ R-Value _____ ☐ Post Calibration Verification (L/s) Flow Column # _____ 32mm _____ (.139 - .169) 36mm _____ (.156 - .190) 53mm _____ (.228 - .278) 103mm _____ (.447 - .547) Maintenance ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other _____ Suggested Service <u>Power Issue? DNS</u> _____ _____
Simulator	Serial #	Lot #/Exp															
0.05																	
0.08																	
0.20																	
0.08 DGS	N/A																

Optical Bench Calibration ☐ Optical Bench Calibration N/A ☐ Optical Bench Calibration Complete Barometric Pressure Gauge _____ ID # _____ <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td></td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.100</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.400</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table> ☐ Post Calibration Stability Checks <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.05</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.08</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.20</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.08 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table> Notes: <u>Instrument would not start-up. Suspect toggle switch in the back. Sending to EE for repair. DNS</u> <u>SP</u>	Simulator	Serial Number	Lot Number	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.400				0.080 DGS	N/A			Simulator	Serial Number	Lot Number	Expiration	0.05				0.08				0.20				0.08 DGS	N/A			Department Inspection ☐ Barometric Pressure _____ Gauge ID# _____ Instrument Mouth Alcohol Solution Lot # _____ Acetone Stock Solution Lot # _____ <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.00</td> <td></td> </tr> <tr> <td>Interferent</td> <td></td> </tr> <tr> <td>0.05</td> <td></td> </tr> <tr> <td>0.08</td> <td></td> </tr> <tr> <td>0.20</td> <td></td> </tr> </tbody> </table> Attachments ☐ Form 41 ☐ Pre-Stability Tests ☐ Flow Calibration ☐ Optical Bench Cal ☐ Post-Stability Tests ☐ Other _____ ☒ Instrument Complies with Chapter 11D-8, FAC ☒ Instrument Does Not Comply with Chapter 11D-8, FAC ☐ Return to/Place into Evidentiary Use ☒ Remain Out of Evidentiary Use ☐ Conduct an Agency Inspection Before Evidentiary Use	Simulator	Serial Number	0.00		Interferent		0.05		0.08		0.20	
Simulator	Serial Number	Lot Number	Expiration																																																										
0.000		N/A	N/A																																																										
0.040																																																													
0.100																																																													
0.200																																																													
0.400																																																													
0.080 DGS	N/A																																																												
Simulator	Serial Number	Lot Number	Expiration																																																										
0.05																																																													
0.08																																																													
0.20																																																													
0.08 DGS	N/A																																																												
Simulator	Serial Number																																																												
0.00																																																													
Interferent																																																													
0.05																																																													
0.08																																																													
0.20																																																													

Quality Control Review

Date

Stability Checks #80-001350 Fruitland Park P.D. 5/4/17


 FRUITLAND PARK P.D.
 Intoxilyzer - Alcohol Analyzer
 Model 8000
 05/04/2017
 Software: 8100.27
 SN 80-001350

Test	g/210L	Time
Air Blank	0.000	09:16
Control Test	0.080	09:17
Air Blank	0.000	09:17
Control Test	0.078	09:17
Air Blank	0.000	09:18
Control Test	0.079	09:18
Air Blank	0.000	09:19
Control Test Stats		
Average	0.0790	
Std Dev	0.0010	
Rel Std Dev(%)	1.2658	

FRUITLAND PARK P.D.
 Intoxilyzer - Alcohol Analyzer
 Model 8000
 05/04/2017
 Software: 8100.27
 SN 80-001350

Test	g/210L	Time
Air Blank	0.000	09:31
Control Test	0.198	09:31
Air Blank	0.000	09:32
Control Test	0.197	09:33
Air Blank	0.000	09:33
Control Test	0.197	09:34
Air Blank	0.000	09:35
Control Test Stats		
Average	0.1973	
Std Dev	0.0006	
Rel Std Dev(%)	0.2926	

FRUITLAND PARK P.D.
 Intoxilyzer - Alcohol Analyzer
 Model 8000
 05/04/2017
 Software: 8100.27
 SN 80-001350

Test	g/210L	Time
Air Blank	0.000	09:26
Control Test	0.080	09:26
Air Blank	0.000	09:27
Control Test	0.080	09:28
Air Blank	0.000	09:28
Control Test	0.080	09:29
Air Blank	0.000	09:29
Control Test Stats		
Average	0.0800	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

FRUITLAND PARK P.D.
 Intoxilyzer - Alcohol Analyzer
 Model 8000
 05/04/2017
 Software: 8100.27
 SN 80-001350

Test	g/210L	Time
Air Blank	0.000	09:21
Control Test	0.050	09:22
Air Blank	0.000	09:22
Control Test	0.049	09:23
Air Blank	0.000	09:23
Control Test	0.049	09:24
Air Blank	0.000	09:25
Control Test Stats		
Average	0.0493	
Std Dev	0.0006	
Rel Std Dev(%)	1.1703	


 Operator's Signature


 Operator's Signature


 Operator's Signature


 Operator's Signature

SP

JP