







## INSTRUMENT PROCESSING SHEET

Agency Charlotte CSOS/N 80-000946Florida Department of  
Law EnforcementDate In 11/28/2023 DI Completion Date 12/05/2023☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake	By TDG	Quality Checks	By TDG	Date	12/01/2023	Flow Calibration	By	Date																																							
<input checked="" type="checkbox"/> Annual <input type="checkbox"/> Registration <input type="checkbox"/> Return from CMI / EE  Visual Inspection: <input checked="" type="checkbox"/> Case <input checked="" type="checkbox"/> Handle <input checked="" type="checkbox"/> Keyboard <input checked="" type="checkbox"/> Dry Gas Shelf <input checked="" type="checkbox"/> Feet <input checked="" type="checkbox"/> Breath Tube <input checked="" type="checkbox"/> Ports <input checked="" type="checkbox"/> Screws Tight  Other Equipment/ Accessories: <input type="checkbox"/> Power cord <input type="checkbox"/> Printer Cable <input checked="" type="checkbox"/> Static Bag <input type="checkbox"/> 12V DC Cable  Notes: <u>All reports the instrument detected interferents and gave purge fails.</u>		<input checked="" type="checkbox"/> Breath Tube Screen <input checked="" type="checkbox"/> Replace External O-Rings <input checked="" type="checkbox"/> Instrument Set Up Verified <input checked="" type="checkbox"/> R-Value <u>155</u> <input checked="" type="checkbox"/> Flow Verification (L/s) Flow Column # <u>ATP104</u> 32 mm <u>0.144</u> (.139 - .169) 36 mm <u>0.167</u> (.156 - .190) 53 mm <u>0.234</u> (.228 - .278) 103 mm <u>0.500</u> (.447 - .547) <input checked="" type="checkbox"/> Barometric Pressure Check Gauge ID # <u>26932</u> <input checked="" type="checkbox"/> Stability Checks				Flow Column # _____ <input type="checkbox"/> 5L/min - 17mm <input type="checkbox"/> 15L/min - 53mm <input type="checkbox"/> 30L/min - 103mm <input type="checkbox"/> R-Value _____ <input type="checkbox"/> Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)																																									
		<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #/Exp</th></tr></thead><tbody><tr><td>0.050</td><td>MP5094</td><td>202303K 03/29/2025</td></tr><tr><td>0.080</td><td>MP5095</td><td>202303L 03/29/2025</td></tr><tr><td>0.200</td><td>MP5096</td><td>202304C 04/05/2025</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>AG223802 08/26/2024</td></tr></tbody></table>		Simulator	Serial #	Lot #/Exp	0.050	MP5094	202303K 03/29/2025	0.080	MP5095	202303L 03/29/2025	0.200	MP5096	202304C 04/05/2025	0.080 DGS	N/A	AG223802 08/26/2024					<b>Maintenance</b> By TDG <u>  </u> <input type="checkbox"/> Battery Replacement <input type="checkbox"/> Dry Gas Regulator Replacement <input type="checkbox"/> Breath Tube Replacement <input checked="" type="checkbox"/> Other <u>Tared tank sensor on 12/1.</u>																								
Simulator	Serial #	Lot #/Exp																																													
0.050	MP5094	202303K 03/29/2025																																													
0.080	MP5095	202303L 03/29/2025																																													
0.200	MP5096	202304C 04/05/2025																																													
0.080 DGS	N/A	AG223802 08/26/2024																																													
<b>Calibration Adjustment</b> By <u>  </u> Barometric Pressure Gauge _____ ID # _____						<b>Department Inspection</b> By TDG <u>  </u> Barometric Pressure ID# <u>26932</u> Gauge <u>1018</u> Instrument <u>1019</u> Mouth Alcohol Solution Lot # <u>2023-A</u> Acetone Stock Solution Lot # <u>2022-B</u>																																									
<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.000</td><td></td><td>N/A</td><td>N/A</td></tr><tr><td>0.040</td><td></td><td></td><td></td></tr><tr><td>0.100</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.300</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table>		Simulator	Serial #	Lot #	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A							<table border="1"><thead><tr><th>Simulator</th><th>Serial Number</th></tr></thead><tbody><tr><td>0.000</td><td>MP5092</td></tr><tr><td>Interferent</td><td>MP5093</td></tr><tr><td>0.050</td><td>MP5094</td></tr><tr><td>0.080</td><td>MP5095</td></tr><tr><td>0.200</td><td>MP5096</td></tr></tbody></table>		Simulator	Serial Number	0.000	MP5092	Interferent	MP5093	0.050	MP5094	0.080	MP5095	0.200	MP5096
Simulator	Serial #	Lot #	Expiration																																												
0.000		N/A	N/A																																												
0.040																																															
0.100																																															
0.200																																															
0.300																																															
0.080 DGS	N/A																																														
Simulator	Serial Number																																														
0.000	MP5092																																														
Interferent	MP5093																																														
0.050	MP5094																																														
0.080	MP5095																																														
0.200	MP5096																																														
<input type="checkbox"/> Post Calibration Adjustment Stability Checks						<b>Attachments</b> <input checked="" type="checkbox"/> Form 41 <input type="checkbox"/> Post-Stability Checks <input checked="" type="checkbox"/> Stability Checks <input type="checkbox"/> Flow Calibration <u>TDG</u> <input checked="" type="checkbox"/> Calibration Certificate <input type="checkbox"/> Form 40 <u>12/18/2023</u> <input type="checkbox"/> Calibration Adjustment <input checked="" type="checkbox"/> Other <u>Extra Stabilities and AR Stabilities</u>																																									
<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.050</td><td></td><td></td><td></td></tr><tr><td>0.080</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table>		Simulator	Serial #	Lot #	Expiration	0.050				0.080				0.200				0.080 DGS	N/A							<input checked="" type="checkbox"/> Instrument Complies with Chapter 11D-8, FAC <input type="checkbox"/> Instrument Does Not Comply with Chapter 11D-8, FAC <input checked="" type="checkbox"/> Return to/Place into Evidentiary Use <input type="checkbox"/> Remain Out of Evidentiary Use <input checked="" type="checkbox"/> Conduct an Agency Inspection Before Evidentiary Use																					
Simulator	Serial #	Lot #	Expiration																																												
0.050																																															
0.080																																															
0.200																																															
0.080 DGS	N/A																																														
Notes/Suggested Service: <u>Checked breath tube screen and replaced o-rings on 11/29. After DGS Test failed during Stability Checks, the room was aired out. Ran additional DGS Stabilities the same day and a different day (see Extra Stabilities) using the same lot/expiration date. All passed. (TDG)</u>						Benjamin Siddoway <small>Digitally signed by Benjamin Siddoway Date: 2023.12.05 09:40:45 -0500</small> Phil Nicodemo <small>Digitally signed by Phil Nicodemo Date: 2023.12.19 14:35:47 -0500</small>																																									
<u>Added plastic caps for return to agency. (TDG)</u>						Tech Review / Date _____ Admin Review / Date _____																																									

Admin Review: Conducted 3 additional DGS stabilities (AR Stabilities) to conform to ATP TSOP. Used lot #01923080A3, expiration date 02/05/2025. Attached results. Updated Attachments on IPS. (TDG 12/18/2023)

# Stability Checks

0.05g/210L 0.047 to 0.053	0.08g/210L 0.077 to 0.083	0.20g/210L 0.194 to 0.206	DGS 0.08g/210L 0.077 to 0.083 $\leq 0.003$ of Wet																																																																																																																																																
<div>✓</div> <div>CHARLOTTE COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-000946 12/01/2023 Software: 8100.27</div> <table><tr><th>Test</th><th>g/210L</th><th>Time</th></tr><tr><td>Air Blank</td><td>0.000</td><td>12:30</td></tr><tr><td>Control Test</td><td>0.050</td><td>12:31</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:31</td></tr><tr><td>Control Test</td><td>0.050</td><td>12:32</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:32</td></tr><tr><td>Control Test</td><td>0.049</td><td>12:33</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:34</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0497</td><td></td></tr><tr><td>Std Dev</td><td>0.0006</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>1.1625</td><td></td></tr></table> <div>Operator's Signature </div>	Test	g/210L	Time	Air Blank	0.000	12:30	Control Test	0.050	12:31	Air Blank	0.000	12:31	Control Test	0.050	12:32	Air Blank	0.000	12:32	Control Test	0.049	12:33	Air Blank	0.000	12:34	Control Test Stats			Average	0.0497		Std Dev	0.0006		Rel Std Dev(%)	1.1625		<div>✓</div> <div>CHARLOTTE COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-000946 12/01/2023 Software: 8100.27</div> <table><tr><th>Test</th><th>g/210L</th><th>Time</th></tr><tr><td>Air Blank</td><td>0.000</td><td>12:37</td></tr><tr><td>Control Test</td><td>0.078</td><td>12:37</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:38</td></tr><tr><td>Control Test</td><td>0.078</td><td>12:39</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:39</td></tr><tr><td>Control Test</td><td>0.079</td><td>12:40</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:40</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0783</td><td></td></tr><tr><td>Std Dev</td><td>0.0006</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.7370</td><td></td></tr></table> <div>Operator's Signature </div>	Test	g/210L	Time	Air Blank	0.000	12:37	Control Test	0.078	12:37	Air Blank	0.000	12:38	Control Test	0.078	12:39	Air Blank	0.000	12:39	Control Test	0.079	12:40	Air Blank	0.000	12:40	Control Test Stats			Average	0.0783		Std Dev	0.0006		Rel Std Dev(%)	0.7370		<div>✓</div> <div>CHARLOTTE COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-000946 12/01/2023 Software: 8100.27</div> <table><tr><th>Test</th><th>g/210L</th><th>Time</th></tr><tr><td>Air Blank</td><td>0.000</td><td>12:44</td></tr><tr><td>Control Test</td><td>0.197</td><td>12:45</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:45</td></tr><tr><td>Control Test</td><td>0.197</td><td>12:46</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:46</td></tr><tr><td>Control Test</td><td>0.197</td><td>12:47</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:48</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.1970</td><td></td></tr><tr><td>Std Dev</td><td>0.0000</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr></table> <div>Operator's Signature </div>	Test	g/210L	Time	Air Blank	0.000	12:44	Control Test	0.197	12:45	Air Blank	0.000	12:45	Control Test	0.197	12:46	Air Blank	0.000	12:46	Control Test	0.197	12:47	Air Blank	0.000	12:48	Control Test Stats			Average	0.1970		Std Dev	0.0000		Rel Std Dev(%)	0.0000		<div>✗</div> <div>CHARLOTTE COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-000946 12/01/2023 Software: 8100.27</div> <table><tr><th>Test</th><th>g/210L</th><th>Time</th></tr><tr><td>Air Blank</td><td>0.000</td><td>12:44</td></tr><tr><td>Control Test</td><td>0.197</td><td>12:45</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:45</td></tr><tr><td>Control Test</td><td>0.197</td><td>12:46</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:46</td></tr><tr><td>Control Test</td><td>0.197</td><td>12:47</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:48</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.1970</td><td></td></tr><tr><td>Std Dev</td><td>0.0000</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr></table> <div>Operator's Signature </div>	Test	g/210L	Time	Air Blank	0.000	12:44	Control Test	0.197	12:45	Air Blank	0.000	12:45	Control Test	0.197	12:46	Air Blank	0.000	12:46	Control Test	0.197	12:47	Air Blank	0.000	12:48	Control Test Stats			Average	0.1970		Std Dev	0.0000		Rel Std Dev(%)	0.0000	
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	12:30																																																																																																																																																	
Control Test	0.050	12:31																																																																																																																																																	
Air Blank	0.000	12:31																																																																																																																																																	
Control Test	0.050	12:32																																																																																																																																																	
Air Blank	0.000	12:32																																																																																																																																																	
Control Test	0.049	12:33																																																																																																																																																	
Air Blank	0.000	12:34																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0497																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	1.1625																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	12:37																																																																																																																																																	
Control Test	0.078	12:37																																																																																																																																																	
Air Blank	0.000	12:38																																																																																																																																																	
Control Test	0.078	12:39																																																																																																																																																	
Air Blank	0.000	12:39																																																																																																																																																	
Control Test	0.079	12:40																																																																																																																																																	
Air Blank	0.000	12:40																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0783																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	0.7370																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	12:44																																																																																																																																																	
Control Test	0.197	12:45																																																																																																																																																	
Air Blank	0.000	12:45																																																																																																																																																	
Control Test	0.197	12:46																																																																																																																																																	
Air Blank	0.000	12:46																																																																																																																																																	
Control Test	0.197	12:47																																																																																																																																																	
Air Blank	0.000	12:48																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.1970																																																																																																																																																		
Std Dev	0.0000																																																																																																																																																		
Rel Std Dev(%)	0.0000																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	12:44																																																																																																																																																	
Control Test	0.197	12:45																																																																																																																																																	
Air Blank	0.000	12:45																																																																																																																																																	
Control Test	0.197	12:46																																																																																																																																																	
Air Blank	0.000	12:46																																																																																																																																																	
Control Test	0.197	12:47																																																																																																																																																	
Air Blank	0.000	12:48																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.1970																																																																																																																																																		
Std Dev	0.0000																																																																																																																																																		
Rel Std Dev(%)	0.0000																																																																																																																																																		



065  
CHARLOTTE COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-000946  
12/01/2023  
Software: 8100.27

Extra  
Stability  
Tests

065  
CHARLOTTE COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-000946  
12/05/2023  
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:50
Control Test	0.078	12:50
Air Blank	0.000	12:50
Control Test	0.077	12:51
Air Blank	0.000	12:51
Control Test	0.078	12:52
Air Blank	0.000	12:52
Control Test	0.077	12:52
Air Blank	0.000	12:53
Control Test	0.078	12:53
Air Blank	0.000	12:54
Control Test	0.078	12:54
Air Blank	0.000	12:54
Control Test	0.078	12:55
Air Blank	0.000	12:55
Control Test	0.078	12:56
Air Blank	0.000	12:56
Control Test	0.078	12:56
Air Blank	0.000	12:57
Control Test	0.078	12:57
Air Blank	0.000	12:58
Control Test	0.078	12:58
Air Blank	0.000	12:59
Control Test	0.078	12:59
Air Blank	0.000	12:59
Control Test	0.077	13:00
Air Blank	0.000	13:00
Control Test	0.078	13:00
Air Blank	0.000	13:01
Control Test	0.078	13:01
Air Blank	0.000	13:02
Control Test	0.078	13:02
Air Blank	0.000	13:03
Control Test	0.078	13:03
Air Blank	0.000	13:03
Control Test	0.078	13:04
Air Blank	0.000	13:04
Control Test	0.078	13:04
Air Blank	0.000	13:05
Control Test	0.078	13:05
Air Blank	0.000	13:06
Control Test Stats		
Average	0.0778	
Std Dev	0.0004	
Rel Std Dev(%)	0.4706	

Test	g/210L	Time
Air Blank	0.000	09:25
Control Test	0.080	09:26
Air Blank	0.000	09:26
Control Test	0.080	09:26
Air Blank	0.000	09:27
Control Test	0.079	09:27
Air Blank	0.000	09:28
Control Test	0.079	09:28
Air Blank	0.000	09:28
Control Test	0.079	09:29
Air Blank	0.000	09:29
Control Test	0.080	09:30
Air Blank	0.000	09:30
Control Test	0.080	09:30
Air Blank	0.000	09:31
Control Test	0.079	09:31
Air Blank	0.000	09:32
Control Test	0.080	09:32
Air Blank	0.000	09:33
Control Test	0.081	09:33
Air Blank	0.000	09:33
Control Test	0.080	09:34
Air Blank	0.000	09:34
Control Test	0.079	09:34
Air Blank	0.000	09:35
Control Test	0.080	09:35
Air Blank	0.000	09:36
Control Test	0.080	09:36
Air Blank	0.000	09:37
Control Test	0.080	09:37
Air Blank	0.000	09:37
Control Test	0.079	09:38
Air Blank	0.000	09:38
Control Test	0.080	09:38
Air Blank	0.000	09:39
Control Test	0.079	09:39
Air Blank	0.000	09:40
Control Test	0.080	09:40
Air Blank	0.000	09:41
Control Test	0.079	09:41
Air Blank	0.000	09:41
Control Test Stats		
Average	0.0797	
Std Dev	0.0006	
Rel Std Dev(%)	0.7372	

MG  
Operator's Signature

MG  
Operator's Signature

# Florida Department of Law Enforcement Alcohol Testing Program

## DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: CHARLOTTE COUNTY SO  
Time of Inspection: 15:22

Date of Inspection: 12/05/2023

Serial Number: 80-000946  
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202303K Exp: 03/29/2025	0.08g/210L Test (g/210L) Lot#:202303L Exp: 03/29/2025	0.20g/210L Test (g/210L) Lot#:202304C Exp: 04/05/2025	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:01923080A3 Exp: 02/05/2025
0.000	0.050	0.079	0.197	0.080
0.000	0.050	0.080	0.197	0.080
0.000	0.050	0.079	0.197	0.079
0.000	0.050	0.080	0.197	0.079
0.000	0.050	0.079	0.197	0.079
0.000	0.050	0.079	0.197	0.078
0.000	0.050	0.080	0.197	0.079
0.000	0.050	0.079	0.197	0.079
0.000	0.050	0.079	0.197	0.079
0.000	0.050	0.079	0.196	0.079


Standard Deviations	0.0000	0.0004	0.0003	0.0005
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0003 Number of Simulators Used: 5

Remarks:

The above instrument complies ( ☒ ) does not comply ( ☐ ) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



TAYLOR D GUTSCHOW

Signature and Printed Name

12/05/2023  
Date





# Calibration Certificate

Florida Department of Law Enforcement  
Alcohol Testing Program  
4700 Terminal Drive, Suite 1  
Ft. Myers, FL 33907

This is to certify the calibration of Intoxilyzer 8000 serial number 80-000946, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-000946</u>	UNCERTAINTY* $\pm$
Owning Agency:	<u>CHARLOTTE COUNTY SO</u>	0.050 g/ 210 L 0.004
Calibration Date:	<u>12/05/2023</u>	0.080 g/ 210 L 0.004
Calibration Time:	<u>15:22</u>	0.200 g/ 210 L 0.007
		0.080 g/ 210 L Dry Gas Control 0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within  $\pm 0.005$  or 5%, whichever is greater, of the target alcohol concentration.

\*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ( $k=3$ ).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

## TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards. This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

12/05/2023

Date

TAYLOR D GUTSCHOW,

Department Inspector

FDLE/ATP Form 69 December 2021

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality

AR: Stabilities  
MG

TDG 12/18/2023

CHARLOTTE COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-000946  
12/18/2023  
Software: 8100.27

0.65

Test	g/210L	Time
Air Blank	0.000	14:50
Control Test	0.080	14:50
Air Blank	0.000	14:51
Control Test	0.079	14:51
Air Blank	0.000	14:52
Control Test	0.079	14:52
Air Blank	0.000	14:53
Control Test Stats		
Average	0.0793	
Std Dev	0.0006	
Rel Std Dev(%)	0.7277	

MG

Operator's Signature





## INSTRUMENT PROCESSING SHEET

Agency Charlotte CSOS/N 80-000946Florida Department of  
Law EnforcementDate In 07/13/2023 DI Completion Date 07/17/2023☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake	By TDG	Quality Checks	By TDG	Date <u>07/14/2023</u>	Flow Calibration	By	Date																																								
<input checked="" type="checkbox"/> Annual <input type="checkbox"/> Registration <input type="checkbox"/> Return from CMI / EE  Visual Inspection: <input checked="" type="checkbox"/> Case <input checked="" type="checkbox"/> Handle <input checked="" type="checkbox"/> Keyboard <input checked="" type="checkbox"/> Dry Gas Shelf <input checked="" type="checkbox"/> Feet <input checked="" type="checkbox"/> Breath Tube <input checked="" type="checkbox"/> Ports <input checked="" type="checkbox"/> Screws Tight  Other Equipment/ Accessories: <input type="checkbox"/> Power cord <input type="checkbox"/> Printer Cable <input checked="" type="checkbox"/> Static Bag <input type="checkbox"/> 12V DC Cable  Notes: <u>AI reports Purge Fail and</u> <u>Int Det messages.</u>		<input checked="" type="checkbox"/> Breath Tube Screen <input checked="" type="checkbox"/> Replace External O-Rings <input checked="" type="checkbox"/> Instrument Set Up Verified <input checked="" type="checkbox"/> R-Value <u>154</u> <input checked="" type="checkbox"/> Flow Verification (L/s) Flow Column # <u>ATP104</u> 32 mm <u>0.148</u> (.139 - .169) 36 mm <u>0.167</u> (.156 - .190) 53 mm <u>0.238</u> (.228 - .278) 103 mm <u>0.511</u> (.447 - .547) <input checked="" type="checkbox"/> Barometric Pressure Check Gauge ID # <u>26932</u> <input checked="" type="checkbox"/> Stability Checks			Flow Column # _____ <input type="checkbox"/> 5L/min - 17mm <input type="checkbox"/> 15L/min - 53mm <input type="checkbox"/> 30L/min - 103mm <input type="checkbox"/> R-Value _____ <input type="checkbox"/> Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)																																										
		<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #/Exp</th></tr></thead><tbody><tr><td>0.050</td><td>MP5094</td><td>202201C 01/11/2024</td></tr><tr><td>0.080</td><td>MP5095</td><td>202201D 01/18/2024</td></tr><tr><td>0.200</td><td>MP5096</td><td>202201E 01/18/2024</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>AG223802 08/26/2024</td></tr></tbody></table>	Simulator	Serial #	Lot #/Exp	0.050	MP5094	202201C 01/11/2024	0.080	MP5095	202201D 01/18/2024	0.200	MP5096	202201E 01/18/2024	0.080 DGS	N/A	AG223802 08/26/2024			<b>Maintenance</b> By _____ <input type="checkbox"/> Battery Replacement <input type="checkbox"/> Dry Gas Regulator Replacement <input type="checkbox"/> Breath Tube Replacement <input type="checkbox"/> Other _____																											
Simulator	Serial #	Lot #/Exp																																													
0.050	MP5094	202201C 01/11/2024																																													
0.080	MP5095	202201D 01/18/2024																																													
0.200	MP5096	202201E 01/18/2024																																													
0.080 DGS	N/A	AG223802 08/26/2024																																													
<b>Calibration Adjustment</b> By _____		<b>Department Inspection</b> By <u>TDG</u>																																													
Barometric Pressure Gauge _____ ID # _____		Barometric Pressure ID# <u>26932</u>																																													
<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.000</td><td></td><td>N/A</td><td>N/A</td></tr><tr><td>0.040</td><td></td><td></td><td></td></tr><tr><td>0.100</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.300</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table>		Simulator	Serial #	Lot #	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			Gauge <u>1017</u> Instrument <u>1018</u> Mouth Alcohol Solution Lot # <u>2021-D</u> Acetone Stock Solution Lot # <u>2022-B</u>		<table border="1"><thead><tr><th>Simulator</th><th>Serial Number</th></tr></thead><tbody><tr><td>0.000</td><td>MP5092</td></tr><tr><td>Interferent</td><td>MP5093</td></tr><tr><td>0.050</td><td>MP5094</td></tr><tr><td>0.080</td><td>MP5095</td></tr><tr><td>0.200</td><td>MP5096</td></tr></tbody></table>		Simulator	Serial Number	0.000	MP5092	Interferent	MP5093	0.050	MP5094	0.080	MP5095	0.200	MP5096		
Simulator	Serial #	Lot #	Expiration																																												
0.000		N/A	N/A																																												
0.040																																															
0.100																																															
0.200																																															
0.300																																															
0.080 DGS	N/A																																														
Simulator	Serial Number																																														
0.000	MP5092																																														
Interferent	MP5093																																														
0.050	MP5094																																														
0.080	MP5095																																														
0.200	MP5096																																														
<input type="checkbox"/> Post Calibration Adjustment Stability Checks		<b>Attachments</b>																																													
<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.050</td><td></td><td></td><td></td></tr><tr><td>0.080</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table>		Simulator	Serial #	Lot #	Expiration	0.050				0.080				0.200				0.080 DGS	N/A			<input checked="" type="checkbox"/> Form 41 <input checked="" type="checkbox"/> Stability Checks <input checked="" type="checkbox"/> Calibration Certificate <input type="checkbox"/> Calibration Adjustment		<input type="checkbox"/> Post-Stability Checks <input type="checkbox"/> Flow Calibration <input checked="" type="checkbox"/> Form 40 <input type="checkbox"/> Other _____																							
Simulator	Serial #	Lot #	Expiration																																												
0.050																																															
0.080																																															
0.200																																															
0.080 DGS	N/A																																														
Notes/Suggested Service: <u>Added plastic caps for return to</u> <u>agency. Unable to duplicate the messages reported by</u> <u>the Agency Inspector; if they persist after the instrument</u> <u>is returned, contact the Department Inspector. (TDG)</u>		<input checked="" type="checkbox"/> Instrument Complies with Chapter 11D-8, FAC <input type="checkbox"/> Instrument Does Not Comply with Chapter 11D-8, FAC <input checked="" type="checkbox"/> Return to/Place into Evidentiary Use <input type="checkbox"/> Remain Out of Evidentiary Use <input checked="" type="checkbox"/> Conduct an Agency Inspection Before Evidentiary Use																																													
		Israel Soto <small>Digitally signed by Israel Soto Date: 2023.07.18 07:53:52 -04'00'</small>		Phil Nicodemo <small>Digitally signed by Phil Nicodemo Date: 2023.07.21 13:47:09 -04'00'</small>																																											
		Tech Review / Date		Admin Review / Date																																											



# Florida Department of Law Enforcement Alcohol Testing Program

## AGENCY INSPECTION REPORT - INTOXILYZER 8000

Agency: CHARLOTTE COUNTY SO  
Time of Inspection: 13:02

Date of Inspection: 07/14/2023

Serial Number: 80-000946  
Software: 8100.27

Check or Test	YES	NO
Date and/or Time Adjusted		No
Diagnostic Check (Pre-Inspection): OK		No
Alcohol Free Subject Test: 0.000		No
Mouth Alcohol Test: Slope Not Met		No
Interferent Detect Test: Interferent Detect		No
Diagnostic Check (Post-Inspection): OK		No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#: Exp:	0.08g/210L Test (g/210L) Lot#: Exp:	0.20g/210L Test (g/210L) Lot#: Exp:	0.08 g/210L Dry Gas Std Test (g/210L) Lot#: Exp:

Number of Simulators Used: \_\_\_\_\_

Remarks:

AI NOT CONDUCTED. BYPASSED TO BRING OUT OF DISABLED MODE.

*Not determined* *MG*  
*7/14/2023*

The above instrument complies ( ☒ ) does not comply ( ) with Chapter 11D-8, FAC.

I certify that I hold a valid Florida Department of Law Enforcement Agency Inspector Permit and that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

*Taylor D Gutschow*

Signature and Printed Name TAYLOR D GUTSCHOW

07/14/2023  
Date



Type of Test	Serial Number	Agency	Date	Performed By
Stabilities	80-000946	Charlotte CSO	7/14/2023	TDG MG

0.05g/210L	0.08g/210L	0.20g/210L	DGS 0.08g/210L																																																																																																																																				
0.047 to 0.053	0.077 to 0.083	0.194 to 0.206	0.077 to 0.083																																																																																																																																				
✓	✓	✓	✓																																																																																																																																				
<div>CHARLOTTE COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 07/14/2023 Software: 8100.27</div> <div>Testg/210LTime</div> <table><tr><td>Air Blank</td><td>0.000</td><td>13:41</td></tr><tr><td>Control Test</td><td>0.049</td><td>13:41</td></tr><tr><td>Air Blank</td><td>0.000</td><td>13:42</td></tr><tr><td>Control Test</td><td>0.048</td><td>13:43</td></tr><tr><td>Air Blank</td><td>0.000</td><td>13:43</td></tr><tr><td>Control Test</td><td>0.049</td><td>13:44</td></tr><tr><td>Air Blank</td><td>0.000</td><td>13:44</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0487</td><td></td></tr><tr><td>Std Dev</td><td>0.0006</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>1.1863</td><td></td></tr></table>	Air Blank	0.000	13:41	Control Test	0.049	13:41	Air Blank	0.000	13:42	Control Test	0.048	13:43	Air Blank	0.000	13:43	Control Test	0.049	13:44	Air Blank	0.000	13:44	Control Test Stats			Average	0.0487		Std Dev	0.0006		Rel Std Dev(%)	1.1863		<div>CHARLOTTE COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 07/14/2023 Software: 8100.27</div> <div>Testg/210LTime</div> <table><tr><td>Air Blank</td><td>0.000</td><td>13:48</td></tr><tr><td>Control Test</td><td>0.077</td><td>13:48</td></tr><tr><td>Air Blank</td><td>0.000</td><td>13:48</td></tr><tr><td>Control Test</td><td>0.077</td><td>13:49</td></tr><tr><td>Air Blank</td><td>0.000</td><td>13:50</td></tr><tr><td>Control Test</td><td>0.078</td><td>13:50</td></tr><tr><td>Air Blank</td><td>0.000</td><td>13:51</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0773</td><td></td></tr><tr><td>Std Dev</td><td>0.0006</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.7466</td><td></td></tr></table>	Air Blank	0.000	13:48	Control Test	0.077	13:48	Air Blank	0.000	13:48	Control Test	0.077	13:49	Air Blank	0.000	13:50	Control Test	0.078	13:50	Air Blank	0.000	13:51	Control Test Stats			Average	0.0773		Std Dev	0.0006		Rel Std Dev(%)	0.7466		<div>CHARLOTTE COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 07/14/2023 Software: 8100.27</div> <div>Testg/210LTime</div> <table><tr><td>Air Blank</td><td>0.000</td><td>13:55</td></tr><tr><td>Control Test</td><td>0.198</td><td>13:56</td></tr><tr><td>Air Blank</td><td>0.000</td><td>13:56</td></tr><tr><td>Control Test</td><td>0.198</td><td>13:57</td></tr><tr><td>Air Blank</td><td>0.000</td><td>13:57</td></tr><tr><td>Control Test</td><td>0.198</td><td>13:58</td></tr><tr><td>Air Blank</td><td>0.000</td><td>13:58</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.1980</td><td></td></tr><tr><td>Std Dev</td><td>0.0000</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr></table>	Air Blank	0.000	13:55	Control Test	0.198	13:56	Air Blank	0.000	13:56	Control Test	0.198	13:57	Air Blank	0.000	13:57	Control Test	0.198	13:58	Air Blank	0.000	13:58	Control Test Stats			Average	0.1980		Std Dev	0.0000		Rel Std Dev(%)	0.0000		<div>068</div> <div>CHARLOTTE COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 07/14/2023 Software: 8100.27</div> <div>Testg/210LTime</div> <table><tr><td>Air Blank</td><td>0.000</td><td>14:00</td></tr><tr><td>Control Test</td><td>0.078</td><td>14:00</td></tr><tr><td>Air Blank</td><td>0.000</td><td>14:01</td></tr><tr><td>Control Test</td><td>0.077</td><td>14:01</td></tr><tr><td>Air Blank</td><td>0.000</td><td>14:01</td></tr><tr><td>Control Test</td><td>0.078</td><td>14:02</td></tr><tr><td>Air Blank</td><td>0.000</td><td>14:02</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0777</td><td></td></tr><tr><td>Std Dev</td><td>0.0006</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.7434</td><td></td></tr></table>	Air Blank	0.000	14:00	Control Test	0.078	14:00	Air Blank	0.000	14:01	Control Test	0.077	14:01	Air Blank	0.000	14:01	Control Test	0.078	14:02	Air Blank	0.000	14:02	Control Test Stats			Average	0.0777		Std Dev	0.0006		Rel Std Dev(%)	0.7434	
Air Blank	0.000	13:41																																																																																																																																					
Control Test	0.049	13:41																																																																																																																																					
Air Blank	0.000	13:42																																																																																																																																					
Control Test	0.048	13:43																																																																																																																																					
Air Blank	0.000	13:43																																																																																																																																					
Control Test	0.049	13:44																																																																																																																																					
Air Blank	0.000	13:44																																																																																																																																					
Control Test Stats																																																																																																																																							
Average	0.0487																																																																																																																																						
Std Dev	0.0006																																																																																																																																						
Rel Std Dev(%)	1.1863																																																																																																																																						
Air Blank	0.000	13:48																																																																																																																																					
Control Test	0.077	13:48																																																																																																																																					
Air Blank	0.000	13:48																																																																																																																																					
Control Test	0.077	13:49																																																																																																																																					
Air Blank	0.000	13:50																																																																																																																																					
Control Test	0.078	13:50																																																																																																																																					
Air Blank	0.000	13:51																																																																																																																																					
Control Test Stats																																																																																																																																							
Average	0.0773																																																																																																																																						
Std Dev	0.0006																																																																																																																																						
Rel Std Dev(%)	0.7466																																																																																																																																						
Air Blank	0.000	13:55																																																																																																																																					
Control Test	0.198	13:56																																																																																																																																					
Air Blank	0.000	13:56																																																																																																																																					
Control Test	0.198	13:57																																																																																																																																					
Air Blank	0.000	13:57																																																																																																																																					
Control Test	0.198	13:58																																																																																																																																					
Air Blank	0.000	13:58																																																																																																																																					
Control Test Stats																																																																																																																																							
Average	0.1980																																																																																																																																						
Std Dev	0.0000																																																																																																																																						
Rel Std Dev(%)	0.0000																																																																																																																																						
Air Blank	0.000	14:00																																																																																																																																					
Control Test	0.078	14:00																																																																																																																																					
Air Blank	0.000	14:01																																																																																																																																					
Control Test	0.077	14:01																																																																																																																																					
Air Blank	0.000	14:01																																																																																																																																					
Control Test	0.078	14:02																																																																																																																																					
Air Blank	0.000	14:02																																																																																																																																					
Control Test Stats																																																																																																																																							
Average	0.0777																																																																																																																																						
Std Dev	0.0006																																																																																																																																						
Rel Std Dev(%)	0.7434																																																																																																																																						
<div>Operator's Signature</div> <div>MG</div>	<div>Operator's Signature</div> <div>MG</div>	<div>Operator's Signature</div> <div>MG</div>	<div>Operator's Signature</div> <div>MG</div>																																																																																																																																				

Comments:

# Florida Department of Law Enforcement

## Alcohol Testing Program

### DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: CHARLOTTE COUNTY SO  
Time of Inspection: 13:08

Date of Inspection: 07/17/2023

Serial Number: 80-000946  
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202201C Exp: 01/11/2024	0.08g/210L Test (g/210L) Lot#:202201D Exp: 01/18/2024	0.20g/210L Test (g/210L) Lot#:202201E Exp: 01/18/2024	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG223802 Exp: 08/26/2024
0.000	0.049	0.077	0.198	0.078
0.000	0.049	0.078	0.198	0.077
0.000	0.049	0.077	0.198	0.077
0.000	0.049	0.077	0.198	0.077
0.000	0.049	0.077	0.198	0.077
0.000	0.048	0.077	0.197	0.078
0.000	0.049	0.077	0.198	0.077
0.000	0.048	0.077	0.197	0.077
0.000	0.049	0.077	0.197	0.077
0.000	0.049	0.077	0.198	0.077

Standard Deviations	0.0004	0.0003	0.0004	0.0004
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0003 Number of Simulators Used: 5

Remarks:

The above instrument complies ( ☒ ) does not comply ( ☐ ) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



TAYLOR D GUTSCHOW

Signature and Printed Name

07/17/2023  
Date





# Calibration Certificate

Florida Department of Law Enforcement  
Alcohol Testing Program  
4700 Terminal Drive, Suite 1  
Ft. Myers, FL 33907

This is to certify the calibration of Intoxilyzer 8000 serial number 80-000946, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-000946</u>	UNCERTAINTY* $\pm$
Owning Agency:	<u>CHARLOTTE COUNTY SO</u>	0.050 g/ 210 L 0.004
Calibration Date:	<u>07/17/2023</u>	0.080 g/ 210 L 0.004
Calibration Time:	<u>13:08</u>	0.200 g/ 210 L 0.007
		0.080 g/ 210 L Dry Gas Control 0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within  $\pm 0.005$  or 5%, whichever is greater, of the target alcohol concentration.

\*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ( $k=3$ ).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

## TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards. This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

  
TAYLOR D GUTSCHOW,  
Department Inspector

07/17/2023  
Date

FDLE/ATP Form 69 December 2021  
Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality