



## INSTRUMENT PROCESSING SHEET

Agency Collier CSOS/N 80-001651Florida Department of  
Law EnforcementDate In 03/13/2023DI Completion Date 04/28/2023☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake	By TDG	Quality Checks	By TDG	Date <u>04/11/2023</u>	Flow Calibration	By	Date																																																													
<input checked="" type="checkbox"/> Annual <input type="checkbox"/> Registration <input checked="" type="checkbox"/> Return from CMI / EE  Visual Inspection: <input checked="" type="checkbox"/> Case <input checked="" type="checkbox"/> Handle <input checked="" type="checkbox"/> Keyboard <input checked="" type="checkbox"/> Dry Gas Shelf <input checked="" type="checkbox"/> Feet <input checked="" type="checkbox"/> Breath Tube <input checked="" type="checkbox"/> Ports <input checked="" type="checkbox"/> Screws Tight  Other Equipment/ Accessories: <input type="checkbox"/> Power cord <input type="checkbox"/> Printer Cable <input checked="" type="checkbox"/> Static Bag <input type="checkbox"/> 12V DC Cable  Notes: _____ _____ _____ _____ _____ _____ _____ _____ _____ _____		<input checked="" type="checkbox"/> Breath Tube Screen <input checked="" type="checkbox"/> Replace External O-Rings <input checked="" type="checkbox"/> Instrument Set Up Verified <input checked="" type="checkbox"/> R-Value <u>114</u> <input checked="" type="checkbox"/> Flow Verification (L/s) Flow Column # <u>ATP104</u> 32 mm <u>0.152</u> (.139 - .169) 36 mm <u>0.164</u> (.156 - .190) 53 mm <u>0.238</u> (.228 - .278) 103 mm <u>0.511</u> (.447 - .547) <input checked="" type="checkbox"/> Barometric Pressure Check Gauge ID # <u>28663</u> <input checked="" type="checkbox"/> Stability Checks			Flow Column # _____ <input type="checkbox"/> 5L/min - 17mm <input type="checkbox"/> 15L/min - 53mm <input type="checkbox"/> 30L/min - 103mm <input type="checkbox"/> R-Value _____ <input type="checkbox"/> Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)																																																															
				<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #/Exp</th></tr></thead><tbody><tr><td>0.050</td><td>MP6286</td><td>202201C 01/11/2024</td></tr><tr><td>0.080</td><td>MP4864</td><td>202201D 01/18/2024</td></tr><tr><td>0.200</td><td>MP6288</td><td>202201E 01/18/2024</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>AG223802 08/26/2024</td></tr></tbody></table>	Simulator	Serial #	Lot #/Exp	0.050	MP6286	202201C 01/11/2024	0.080	MP4864	202201D 01/18/2024	0.200	MP6288	202201E 01/18/2024	0.080 DGS	N/A	AG223802 08/26/2024	<b>Maintenance</b> By _____ <input type="checkbox"/> Battery Replacement <input type="checkbox"/> Dry Gas Regulator Replacement <input type="checkbox"/> Breath Tube Replacement <input type="checkbox"/> Other _____ _____ _____ _____ _____ _____																																																
Simulator	Serial #	Lot #/Exp																																																																		
0.050	MP6286	202201C 01/11/2024																																																																		
0.080	MP4864	202201D 01/18/2024																																																																		
0.200	MP6288	202201E 01/18/2024																																																																		
0.080 DGS	N/A	AG223802 08/26/2024																																																																		
<b>Calibration Adjustment</b> By TDG Barometric Pressure Gauge <u>1015</u> ID # <u>28199</u> <table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.000</td><td>MP4864</td><td>N/A</td><td>N/A</td></tr><tr><td>0.040</td><td>MP5098</td><td>21410</td><td>09/30/2023</td></tr><tr><td>0.100</td><td>MP5099</td><td>22310</td><td>08/11/2024</td></tr><tr><td>0.200</td><td>MP5100</td><td>22050</td><td>02/07/2024</td></tr><tr><td>0.300</td><td>MP5101</td><td>22220</td><td>06/15/2024</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>AG115904</td><td>06/08/2023</td></tr></tbody></table> <input checked="" type="checkbox"/> Post Calibration Adjustment Stability Checks <table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.050</td><td>MP5094</td><td>202201C</td><td>01/11/2024</td></tr><tr><td>0.080</td><td>MP5095</td><td>202201D</td><td>01/18/2024</td></tr><tr><td>0.200</td><td>MP5096</td><td>202201E</td><td>01/18/2024</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>AG223802</td><td>08/26/2024</td></tr></tbody></table>				Simulator	Serial #	Lot #	Expiration	0.000	MP4864	N/A	N/A	0.040	MP5098	21410	09/30/2023	0.100	MP5099	22310	08/11/2024	0.200	MP5100	22050	02/07/2024	0.300	MP5101	22220	06/15/2024	0.080 DGS	N/A	AG115904	06/08/2023	Simulator	Serial #	Lot #	Expiration	0.050	MP5094	202201C	01/11/2024	0.080	MP5095	202201D	01/18/2024	0.200	MP5096	202201E	01/18/2024	0.080 DGS	N/A	AG223802	08/26/2024	<b>Department Inspection</b> By TDG Barometric Pressure ID# <u>26932</u> Gauge <u>1012</u> Instrument <u>1014</u> Mouth Alcohol Solution Lot # <u>2021-D</u> Acetone Stock Solution Lot # <u>2022-B</u> <table border="1"><thead><tr><th>Simulator</th><th>Serial Number</th></tr></thead><tbody><tr><td>0.000</td><td>MP5092</td></tr><tr><td>Interferent</td><td>MP5093</td></tr><tr><td>0.050</td><td>MP5094</td></tr><tr><td>0.080</td><td>MP5095</td></tr><tr><td>0.200</td><td>MP5096</td></tr></tbody></table> <b>Attachments</b> <table border="1"><tbody><tr><td><input checked="" type="checkbox"/> Form 41 <input checked="" type="checkbox"/> Stability Checks <input checked="" type="checkbox"/> Calibration Certificate <input checked="" type="checkbox"/> Calibration Adjustment</td><td><input checked="" type="checkbox"/> Post-Stability Checks <input type="checkbox"/> Flow Calibration <input checked="" type="checkbox"/> Form 40 (x2) TDG <input type="checkbox"/> Other <u>5/5/2023</u></td></tr></tbody></table>			Simulator	Serial Number	0.000	MP5092	Interferent	MP5093	0.050	MP5094	0.080	MP5095	0.200	MP5096	<input checked="" type="checkbox"/> Form 41 <input checked="" type="checkbox"/> Stability Checks <input checked="" type="checkbox"/> Calibration Certificate <input checked="" type="checkbox"/> Calibration Adjustment	<input checked="" type="checkbox"/> Post-Stability Checks <input type="checkbox"/> Flow Calibration <input checked="" type="checkbox"/> Form 40 (x2) TDG <input type="checkbox"/> Other <u>5/5/2023</u>
Simulator	Serial #	Lot #	Expiration																																																																	
0.000	MP4864	N/A	N/A																																																																	
0.040	MP5098	21410	09/30/2023																																																																	
0.100	MP5099	22310	08/11/2024																																																																	
0.200	MP5100	22050	02/07/2024																																																																	
0.300	MP5101	22220	06/15/2024																																																																	
0.080 DGS	N/A	AG115904	06/08/2023																																																																	
Simulator	Serial #	Lot #	Expiration																																																																	
0.050	MP5094	202201C	01/11/2024																																																																	
0.080	MP5095	202201D	01/18/2024																																																																	
0.200	MP5096	202201E	01/18/2024																																																																	
0.080 DGS	N/A	AG223802	08/26/2024																																																																	
Simulator	Serial Number																																																																			
0.000	MP5092																																																																			
Interferent	MP5093																																																																			
0.050	MP5094																																																																			
0.080	MP5095																																																																			
0.200	MP5096																																																																			
<input checked="" type="checkbox"/> Form 41 <input checked="" type="checkbox"/> Stability Checks <input checked="" type="checkbox"/> Calibration Certificate <input checked="" type="checkbox"/> Calibration Adjustment	<input checked="" type="checkbox"/> Post-Stability Checks <input type="checkbox"/> Flow Calibration <input checked="" type="checkbox"/> Form 40 (x2) TDG <input type="checkbox"/> Other <u>5/5/2023</u>																																																																			
Notes/Suggested Service: <u>Intox locked up while bypassing the AI to bring out of Disabled Mode. Cycled the power, logged back in, and successfully bypassed AI. Told the Agency Inspector to contact his Department Inspector if similar issues are experienced in the field. (TDG)</u>  <u>Agency name was blank in Settings; updated after Form 40 was printed. Form 40 lacks the agency name. (TDG)</u>				<input checked="" type="checkbox"/> Instrument Complies with Chapter 11D-8, FAC <input type="checkbox"/> Instrument Does Not Comply with Chapter 11D-8, FAC <input checked="" type="checkbox"/> Return to/Place Into Evidentiary Use <input type="checkbox"/> Remain Out of Evidentiary Use <input checked="" type="checkbox"/> Conduct an Agency Inspection Before Evidentiary Use																																																																
Israel Soto <small>Digitally signed by Israel Soto Date: 2023.05.01 14:01:12 -0400</small>				Phil Nicodemo <small>Digitally signed by Phil Nicodemo Date: 2023.05.10 10:20:48 -0400</small>																																																																
Tech Review / Date				Admin Review / Date																																																																

Admin Review: Printed test Form 40 to ensure agency name would appear after being updated. Added second Form 40 to Attachments. Added to packet. (TDG 5/5/2023)

# Florida Department of Law Enforcement Alcohol Testing Program

## AGENCY INSPECTION REPORT - INTOXILYZER 8000

Agency:

Time of Inspection: 10:35

Date of Inspection: 04/11/2023

Serial Number: 80-001651

Software: 8100.27

Check or Test	YES	NO
Date and/or Time Adjusted		No
Diagnostic Check (Pre-Inspection): OK		No
Alcohol Free Subject Test: 0.000		No
Mouth Alcohol Test: Slope Not Met		No
Interferent Detect Test: Interferent Detect		No
Diagnostic Check (Post-Inspection): OK		No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#: Exp:	0.08g/210L Test (g/210L) Lot#: Exp:	0.20g/210L Test (g/210L) Lot#: Exp:	0.08 g/210L Dry Gas Std Test (g/210L) Lot#: Exp:

Number of Simulators Used: \_\_\_\_\_

Remarks:

AI NOT CONDUCTED. BYPASSED TO OPERATE INSTRUMENT.

The above instrument complies ( Not determined ) does not comply ( ) with Chapter 11D-8, FAC.

I certify that I hold a valid Florida Department of Law Enforcement Agency Inspector Permit and that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

*Taylor D Gutschow*

TAYLOR D GUTSCHOW

Signature and Printed Name

04/11/2023  
Date



# Florida Department of Law Enforcement Alcohol Testing Program

## AGENCY INSPECTION REPORT - INTOXILYZER 8000

Agency: COLLIER COUNTY SO  
Time of Inspection: 15:52

Date of Inspection: 05/04/2023

Serial Number: 80-001651  
Software: 8100.27

Check or Test	YES	NO
Date and/or Time Adjusted		No
Diagnostic Check (Pre-Inspection): OK		No
Alcohol Free Subject Test: 0.000		No
Mouth Alcohol Test: Slope Not Met		No
Interferent Detect Test: Interferent Detect		No
Diagnostic Check (Post-Inspection): OK		No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#: Exp:	0.08g/210L Test (g/210L) Lot#: Exp:	0.20g/210L Test (g/210L) Lot#: Exp:	0.08 g/210L Dry Gas Std Test (g/210L) Lot#: Exp:

Number of Simulators Used: \_\_\_\_\_

Remarks:  
TEST

*Not determined*

*76  
5/4/2023*

The above instrument complies ( ☒ ) does not comply ( ☐ ) with Chapter 11D-8, FAC.

I certify that I hold a valid Florida Department of Law Enforcement Agency Inspector Permit and that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

*Taylor D Gutschow*

TAYLOR D GUTSCHOW

Signature and Printed Name

05/04/2023  
Date

Type of Test	Serial Number	Agency	Date	Performed By
Stabilities	80-001651	Collier CSO	04/11/2023	TDG MG


0.05g/210L	0.08g/210L	0.20g/210L	0.08g/210L	0.08g/210L
0.047 to 0.053	0.077 to 0.083	0.194 to 0.206	0.077 to 0.083	0.077 to 0.083
✓	✓	✓	✗	✓

COLLIER COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000  
04/11/2023  
Software: 8100.27  
SN 80-001651

Test g/210L Time


Air Blank	0.000	11:39
Control Test	0.048	11:40
Air Blank	0.000	11:40
Control Test	0.048	11:41
Air Blank	0.000	11:41
Control Test	0.049	11:42
Air Blank	0.000	11:43
Control Test	0.0483	
Average	0.0006	
Std Dev	1.1945	
Rel Std Dev(%)		

Operator's Signature 

COLLIER COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000  
04/11/2023  
Software: 8100.27  
SN 80-001651

Test g/210L Time

Air Blank	0.000	11:48
Control Test	0.079	11:48
Air Blank	0.000	11:49
Control Test	0.078	11:50
Air Blank	0.000	11:50
Control Test	0.077	11:51
Air Blank	0.000	11:52
Control Test	0.0780	
Average	0.0010	
Std Dev	1.2821	
Rel Std Dev(%)		


Operator's Signature 

COLLIER COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000  
04/11/2023  
Software: 8100.27  
SN 80-001651

Test g/210L Time


Air Blank	0.000	11:55
Control Test	0.200	11:55
Air Blank	0.000	11:56
Control Test	0.199	11:57
Air Blank	0.000	11:57
Control Test	0.198	11:58
Air Blank	0.000	11:59
Control Test	0.1990	
Average	0.0010	
Std Dev	0.5025	
Rel Std Dev(%)		

Operator's Signature 

COLLIER COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000  
04/11/2023  
Software: 8100.27  
SN 80-001651

Test g/210L Time

Air Blank	0.000	12:00
Control Test	0.076	12:01
Air Blank	0.000	12:01
Control Test	0.077	12:01
Air Blank	0.000	12:02
Control Test	0.077	12:02
Air Blank	0.000	12:03
Control Test	0.0767	
Average	0.0006	
Std Dev	0.7531	
Rel Std Dev(%)		

Operator's Signature 

**Comments:** An optical cal adjust will be conducted. MG 4/11/2023



COLLIER COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-001651  
04/27/2023 10:07:24

Optical Cal  
ML 4/27/23

Auto Calibration

pg 1 of 2

<<<< 3um >>>> <<<< 9um >>>>

Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1

Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	0.1130	(-0.0230)	0.1170	(-0.0070)
Sample #2	0.0790	(0.0250)	0.1130	(0.0100)
Sample #3	0.0740	(0.0320)	0.1130	(-0.0150)
Sample #4	0.0940	(0.0440)	0.1140	(0.0040)
Avg % Abs	0.0823	(0.0337)	0.1133	(-0.0003)
STD DEV	0.0104	(0.0096)	0.0006	(0.0131)
REL STD DEV	12.642	(28.542)	0.509	(3915.354)

Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1

Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	0.8560	(-0.0240)	1.5650	(-0.0020)
Sample #2	0.8690	(-0.0070)	1.5660	(0.0120)
Sample #3	0.8630	(-0.0040)	1.5470	(0.0210)
Sample #4	0.8640	(0.0250)	1.5550	(0.0100)
Avg % Abs	0.8653	(0.0047)	1.5560	(0.0143)
STD DEV	0.0032	(0.0177)	0.0095	(0.0059)
REL STD DEV	0.371	(378.706)	0.613	(40.880)

Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1

Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	2.0040	(-0.0090)	3.6920	(-0.0030)
Sample #2	1.9890	(0.0060)	3.6500	(0.0360)
Sample #3	1.9620	(0.0330)	3.6720	(0.0250)
Sample #4	1.9790	(0.0450)	3.6610	(0.0360)
Avg % Abs	1.9767	(0.0280)	3.6610	(0.0323)
STD DEV	0.0137	(0.0200)	0.0110	(0.0064)
REL STD DEV	0.691	(71.339)	0.300	(19.642)

Solution = 0.200 g/210L or 0.9524 mg/l, Samples = 4, Discarded = 1

Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	3.8230	(-0.0070)	7.1120	(-0.0230)
Sample #2	3.8080	(0.0320)	7.0380	(0.0320)
Sample #3	3.8040	(0.0440)	7.0250	(0.0600)
Sample #4	3.7870	(0.0470)	7.0080	(0.0640)
Avg % Abs	3.7997	(0.0410)	7.0237	(0.0520)
STD DEV	0.0112	(0.0079)	0.0150	(0.0174)
REL STD DEV	0.293	(19.359)	0.214	(33.530)

Solution = 0.300 g/210L or 1.4286 mg/l, Samples = 4, Discarded = 1

Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	5.6250	(-0.0140)	10.3140	(-0.0200)
Sample #2	5.5540	(0.0410)	10.1840	(0.0940)
Sample #3	5.5690	(0.0420)	10.1800	(0.1030)
Sample #4	5.5490	(0.0720)	10.1460	(0.1320)
Avg % Abs	5.5573	(0.0517)	10.1700	(0.1097)
STD DEV	0.0104	(0.0176)	0.0209	(0.0199)
REL STD DEV	0.187	(34.096)	0.205	(18.107)

COLLIER COUNTY SO  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-001651  
04/27/2023 10:07:24

Optical Cal  
MLG 4/27/23

Auto Calibration

pg 2 of 2

<<<< 3um >>>>  
-----  
Zero Order Coef -217.40  
First Order Coef 2460.13  
Second Order Coef 26.99  
-----

<<<< 9um >>>>  
-----  
-138.20  
1289.03  
12.65  
-----

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	-0.000	0.0003
0.040	0.041	-0.0006
0.100	0.100	0.0002
0.200	0.200	0.0001
0.300	0.300	-0.0000

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	0.000	-0.0002
0.040	0.040	0.0001
0.100	0.100	0.0002
0.200	0.200	-0.0003
0.300	0.300	0.0001

<<<< 3um >>>>  
-----

<<<< 9um >>>>  
-----





Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1

Sample	3um	9um
Sample #1	3057.00	3226.00
Sample #2	3105.00	3264.00
Sample #3	3145.00	3265.00
Sample #4	3121.00	3291.00
Avg	3123.6667	3273.3333
STD DEV	20.1329	15.3080
REL STD DEV	0.645	0.468
H2O adjust (mg/l*10k)	686	536

Barometric Pressure = 1014

\*\*\*\*\*CALIBRATION SUCCESSFUL\*\*\*\*\*

Type of Test	Serial Number	Agency	Date	Performed By
Stabilities (Post-Call)	80-001651	Collier CO	04/27/2023	TDG MLC

0.05g/210L	0.08g/210L	0.20g/210L	DGS 0.08g/210L
0.047 to 0.053	0.077 to 0.083	0.194 to 0.206	0.077 to 0.083
COLLIER COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 04/27/2023 Software: 8100.27	COLLIER COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 04/27/2023 Software: 8100.27	COLLIER COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 04/27/2023 Software: 8100.27	COLLIER COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 04/27/2023 Software: 8100.27
<p>Test: g/210L Time</p> <p>Air Blank 0.000 12:13</p> <p>Control Test 0.049 12:14</p> <p>Air Blank 0.000 12:14</p> <p>Control Test 0.048 12:15</p> <p>Air Blank 0.000 12:15</p> <p>Control Test 0.049 12:16</p> <p>Air Blank 0.000 12:17</p> <p>Control Test Stats</p> <p>Average 0.0487</p> <p>Std Dev 0.0006</p> <p>Rel Std Dev(%) 1.1863</p>	<p>Test: g/210L Time</p> <p>Air Blank 0.000 12:21</p> <p>Control Test 0.078 12:22</p> <p>Air Blank 0.000 12:22</p> <p>Control Test 0.078 12:23</p> <p>Air Blank 0.000 12:23</p> <p>Control Test 0.078 12:24</p> <p>Air Blank 0.000 12:24</p> <p>Control Test Stats</p> <p>Average 0.0780</p> <p>Std Dev 0.0000</p> <p>Rel Std Dev(%) 0.0000</p>	<p>Test: g/210L Time</p> <p>Air Blank 0.000 12:28</p> <p>Control Test 0.202 12:28</p> <p>Air Blank 0.000 12:29</p> <p>Control Test 0.200 12:30</p> <p>Air Blank 0.000 12:30</p> <p>Control Test 0.199 12:31</p> <p>Air Blank 0.000 12:31</p> <p>Control Test Stats</p> <p>Average 0.2003</p> <p>Std Dev 0.0015</p> <p>Rel Std Dev(%) 0.7625</p>	<p>Test: g/210L Time</p> <p>Air Blank 0.000 12:04</p> <p>Control Test 0.080 12:04</p> <p>Air Blank 0.000 12:04</p> <p>Control Test 0.080 12:05</p> <p>Air Blank 0.000 12:05</p> <p>Control Test 0.081 12:06</p> <p>Air Blank 0.000 12:06</p> <p>Control Test Stats</p> <p>Average 0.0803</p> <p>Std Dev 0.0006</p> <p>Rel Std Dev(%) 0.7187</p>
Operator's Signature 	Operator's Signature 	Operator's Signature 	Operator's Signature 

Comments:



# Florida Department of Law Enforcement Alcohol Testing Program

## DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: COLLIER COUNTY SO  
Time of Inspection: 13:00

Date of Inspection: 04/28/2023

Serial Number: 80-001651  
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202201C Exp: 01/11/2024	0.08g/210L Test (g/210L) Lot#:202201D Exp: 01/18/2024	0.20g/210L Test (g/210L) Lot#:202201E Exp: 01/18/2024	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG223802 Exp: 08/26/2024
0.000	0.049	0.079	0.200	0.080
0.000	0.049	0.079	0.200	0.080
0.000	0.049	0.079	0.200	0.080
0.000	0.049	0.079	0.200	0.081
0.000	0.049	0.079	0.200	0.080
0.000	0.049	0.079	0.200	0.081
0.000	0.049	0.079	0.201	0.081
0.000	0.050	0.079	0.200	0.081
0.000	0.049	0.079	0.200	0.081
0.000	0.049	0.079	0.200	0.081


Standard Deviations	0.0003	0.0000	0.0003	0.0005
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0002 Number of Simulators Used: 5

Remarks:

The above instrument complies ( ☒ ) does not comply ( ☐ ) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

 TAYLOR D GUTSCHOW  
Signature and Printed Name

04/28/2023  
Date





# Calibration Certificate

Florida Department of Law Enforcement  
Alcohol Testing Program  
4700 Terminal Drive, Suite 1  
Ft. Myers, FL 33907

This is to certify the calibration of Intoxilyzer 8000 serial number 80-001651, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-001651</u>	UNCERTAINTY * $\pm$
Owning Agency:	<u>COLLIER COUNTY SO</u>	0.050 g/ 210 L 0.004
Calibration Date:	<u>04/28/2023</u>	0.080 g/ 210 L 0.004
Calibration Time:	<u>13:00</u>	0.200 g/ 210 L 0.007
		0.080 g/ 210 L Dry Gas Control 0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within  $\pm 0.005$  or 5%, whichever is greater, of the target alcohol concentration.

\*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence (k=3).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

## TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards. This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

  
TAYLOR D GUTSCHOW,  
Department Inspector

04/28/2023

Date

FDLE/ATP Form 69 December 2021

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality