

INSTRUMENT PROCESSING SHEET

Agency Fair Harbor Beach PD S/N 80-006931
 Date In 10/18/17 Date Out 10/19/17 Ship P/U H/D CMI EE

Intake Performed By <u>SP</u> <input checked="" type="checkbox"/> Registration <input checked="" type="checkbox"/> Annual <input checked="" type="checkbox"/> Return from CMI <input type="checkbox"/> Return from Enforcement Electronics <input type="checkbox"/> Other _____ Visual Inspection: <input checked="" type="checkbox"/> Case <u>ca</u> <input checked="" type="checkbox"/> Handle <u>ca</u> <input checked="" type="checkbox"/> Dry Gas Holder <u>ca</u> <input checked="" type="checkbox"/> Feet <u>ca</u> <input checked="" type="checkbox"/> Keyboard/Plug <u>ca</u> <input checked="" type="checkbox"/> Back/Plugs <u>ca</u> <input checked="" type="checkbox"/> Screws tight <u>ca</u> <input checked="" type="checkbox"/> Breath Hose <u>ca</u> Other Equipment: <input checked="" type="checkbox"/> Power cord <input type="checkbox"/> Printer Cable <input checked="" type="checkbox"/> Other <u>Static Bag</u> Notes: _____ _____ _____	Quality Checks Performed By <u>SP</u> <input checked="" type="checkbox"/> Lab Temp °C <u>22.1</u> <input checked="" type="checkbox"/> Breath Tube Screen <input checked="" type="checkbox"/> Replace External O-Rings <input checked="" type="checkbox"/> Instrument Set Up Verified <input checked="" type="checkbox"/> R-Value <u>150</u> <input checked="" type="checkbox"/> Flow Verification (L/s) Flow Column # <u>ATP103</u> 32mm <u>.160</u> (.139 - .169) 36mm <u>.171</u> (.156 - .190) 53mm <u>.234</u> (.228 - .278) 103mm <u>.480</u> (.447 - .547) <input checked="" type="checkbox"/> Barometric Pressure Check Gauge ID # <u>28421</u> <input checked="" type="checkbox"/> Stability Checks <table border="1" style="width:100%; border-collapse: collapse; margin-top: 5px;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td>0.05</td> <td><u>G2835</u></td> <td><u>201707D</u> <u>7/25/19</u></td> </tr> <tr> <td>0.08</td> <td><u>DR1279</u></td> <td><u>201707E</u> <u>7/25/19</u></td> </tr> <tr> <td>0.20</td> <td><u>SD1025</u></td> <td><u>201707C</u> <u>7/24/19</u></td> </tr> <tr> <td>0.08 DGS</td> <td>N/A</td> <td><u>AG626605</u> <u>9/22/18</u></td> </tr> </tbody> </table>	Simulator	Serial #	Lot #/Exp	0.05	<u>G2835</u>	<u>201707D</u> <u>7/25/19</u>	0.08	<u>DR1279</u>	<u>201707E</u> <u>7/25/19</u>	0.20	<u>SD1025</u>	<u>201707C</u> <u>7/24/19</u>	0.08 DGS	N/A	<u>AG626605</u> <u>9/22/18</u>	Flow Calibration Performed By _____ <input checked="" type="checkbox"/> Flow Calibration N/A <input type="checkbox"/> Flow Calibration Complete Flow Column # _____ <input type="checkbox"/> 5L/min - 17mm <input type="checkbox"/> 15L/min - 453mm <input type="checkbox"/> 30L/min - 103mm <input type="checkbox"/> R-Value _____ <input type="checkbox"/> Post Calibration Verification (L/s) Program Flow Column # _____ 32mm _____ (.139 - .169) 36mm _____ (.156 - .190) 53mm _____ (.228 - .278) 103mm _____ (.447 - .547) Maintenance Performed By _____ <input type="checkbox"/> Battery Replacement <input type="checkbox"/> Dry Gas Regulator Replacement <input type="checkbox"/> Breath Tube Replacement <input type="checkbox"/> Other _____ Quality Checks Cont. Performed By <u>SP</u> Simulator Temperatures °C External Digital Therm. ID#: <u>300504</u> <input checked="" type="checkbox"/> 34°C +/- .2 Serial #: <u>G2835</u> <input checked="" type="checkbox"/> 34°C +/- .2 Serial #: <u>DR1279</u> <input checked="" type="checkbox"/> 34°C +/- .2 Serial #: <u>SD1025</u>
Simulator	Serial #	Lot #/Exp															
0.05	<u>G2835</u>	<u>201707D</u> <u>7/25/19</u>															
0.08	<u>DR1279</u>	<u>201707E</u> <u>7/25/19</u>															
0.20	<u>SD1025</u>	<u>201707C</u> <u>7/24/19</u>															
0.08 DGS	N/A	<u>AG626605</u> <u>9/22/18</u>															

RECEIVED
OCT 23 2017
FDLE
Alcohol Testing Program

Calibration Adjustment Performed By _____ <input checked="" type="checkbox"/> Calibration Adjustment N/A <input type="checkbox"/> Calibration Adjustment Complete Barometric Pressure Gauge _____ ID # _____ <table border="1" style="width:100%; border-collapse: collapse; margin-top: 5px;"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td></td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.100</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.300</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table> <input type="checkbox"/> Post Calibration Adjustment Stability Checks <table border="1" style="width:100%; border-collapse: collapse; margin-top: 5px;"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.05</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.08</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.20</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.08 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table>	Simulator	Serial Number	Lot Number	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			Simulator	Serial Number	Lot Number	Expiration	0.05				0.08				0.20				0.08 DGS	N/A			Department Inspection Performed By <u>SP</u> <input checked="" type="checkbox"/> Barometric Pressure <u>1020</u> Gauge ID# <u>28421</u> <u>1020</u> Instrument Mouth Alcohol Solution Lot # <u>2016-C</u> Acetone Stock Solution Lot # <u>2017-A</u> <table border="1" style="width:100%; border-collapse: collapse; margin-top: 5px;"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.00</td> <td><u>G2880</u></td> </tr> <tr> <td>Interferent</td> <td><u>G2840</u></td> </tr> <tr> <td>0.05</td> <td><u>G2835</u></td> </tr> <tr> <td>0.08</td> <td><u>DR1279</u></td> </tr> <tr> <td>0.20</td> <td><u>SD1025</u></td> </tr> </tbody> </table>	Simulator	Serial Number	0.00	<u>G2880</u>	Interferent	<u>G2840</u>	0.05	<u>G2835</u>	0.08	<u>DR1279</u>	0.20	<u>SD1025</u>
Simulator	Serial Number	Lot Number	Expiration																																																										
0.000		N/A	N/A																																																										
0.040																																																													
0.100																																																													
0.200																																																													
0.300																																																													
0.080 DGS	N/A																																																												
Simulator	Serial Number	Lot Number	Expiration																																																										
0.05																																																													
0.08																																																													
0.20																																																													
0.08 DGS	N/A																																																												
Simulator	Serial Number																																																												
0.00	<u>G2880</u>																																																												
Interferent	<u>G2840</u>																																																												
0.05	<u>G2835</u>																																																												
0.08	<u>DR1279</u>																																																												
0.20	<u>SD1025</u>																																																												

Attachments	
<input checked="" type="checkbox"/> Form 41 <input checked="" type="checkbox"/> Pre-Stability Tests <input type="checkbox"/> Flow Calibration	<input type="checkbox"/> Calibration Adjustment <input type="checkbox"/> Post-Stability Tests <input type="checkbox"/> Other _____

Notes/Suggested Service:
QA/QC ✓ SP

JJ Dahan 10/23/17
 Quality Control Review Date

<input checked="" type="checkbox"/> Instrument Complies with Chapter 11D-8, FAC <input type="checkbox"/> Instrument Does Not Comply with Chapter 11D-8, FAC <input checked="" type="checkbox"/> Return to/Place into Evidentiary Use <input type="checkbox"/> Remain Out of Evidentiary Use <input checked="" type="checkbox"/> Conduct an Agency Inspection Before Evidentiary Use

FT WALTON BEACH PD
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-006931
 10/19/2017
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:50
Control Test	0.049	10:51
Air Blank	0.000	10:51
Control Test	0.049	10:52
Air Blank	0.000	10:52
Control Test	0.049	10:53
Air Blank	0.000	10:54
Control Test Stats		
Average	0.0490	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

P Murphy

 Operator's Signature

FT WALTON BEACH PD
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-006931
 10/19/2017
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:45
Control Test	0.077	10:46
Air Blank	0.000	10:46
Control Test	0.078	10:47
Air Blank	0.000	10:48
Control Test	0.079	10:48
Air Blank	0.000	10:49
Control Test Stats		
Average	0.0780	
Std Dev	0.0010	
Rel Std Dev(%)	1.2821	

P Murphy

 Operator's Signature

FT WALTON BEACH PD
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-006931
 10/19/2017
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:00
Control Test	0.198	11:01
Air Blank	0.000	11:01
Control Test	0.197	11:02
Air Blank	0.000	11:02
Control Test	0.197	11:03
Air Blank	0.000	11:04
Control Test Stats		
Average	0.1973	
Std Dev	0.0006	
Rel Std Dev(%)	0.2926	

P Murphy

 Operator's Signature

FT WALTON BEACH PD
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-006931
 10/19/2017
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:10
Control Test	0.080	11:10
Air Blank	0.000	11:11
Control Test	0.080	11:11
Air Blank	0.000	11:12
Control Test	0.080	11:12
Air Blank	0.000	11:13
Control Test Stats		
Average	0.0800	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

DGS

P Murphy

 Operator's Signature

SP

*10/23/17
 JO*