



INSTRUMENT PROCESSING SHEET

Agency FFWCC

S/N 80-003196

Florida Department of Law Enforcement

Date In 03/18/2019

DI Completion Date 3/20/19

Ship P/U H/D CMI EE

<b>Intake</b> Performed By <u>JD</u> <input checked="" type="checkbox"/> Annual <input type="checkbox"/> Registration <input type="checkbox"/> Return from CMI / EE Visual Inspection: <input checked="" type="checkbox"/> Case <input checked="" type="checkbox"/> Handle <input checked="" type="checkbox"/> Keyboard <input checked="" type="checkbox"/> Dry Gas Shelf <input checked="" type="checkbox"/> Feet <input checked="" type="checkbox"/> Breath Tube <input checked="" type="checkbox"/> Ports <input checked="" type="checkbox"/> Screws Tight Other Equipment/ Accessories: <input checked="" type="checkbox"/> Power cord <input checked="" type="checkbox"/> Printer Cable <input type="checkbox"/> Static Bag <input checked="" type="checkbox"/> 12V DC Cable Notes: <u>Pelican Case</u>		<b>Quality Checks</b> Performed By <u>BK</u> <input checked="" type="checkbox"/> Breath Tube Screen <input checked="" type="checkbox"/> Replace External O-Rings <input checked="" type="checkbox"/> Instrument Set Up Verified <input checked="" type="checkbox"/> R-Value <u>126</u> <input checked="" type="checkbox"/> Flow Verification (L/s) Flow Column # <u>ATP102</u> 32 mm <u>.144</u> (.139 - .169) 36 mm <u>.156</u> (.156 - .190) 53 mm <u>.234</u> (.228 - .278) 103 mm <u>.507</u> (.447 - .547) <input checked="" type="checkbox"/> Barometric Pressure Check Gauge ID # <u>30793</u> <input checked="" type="checkbox"/> Stability Checks		<b>Flow Calibration</b> Performed By _____ Flow Column # _____ <input type="checkbox"/> 5L/min - 17mm <input type="checkbox"/> 15L/min - 53mm <input type="checkbox"/> 30L/min - 103mm <input type="checkbox"/> R-Value _____ <input type="checkbox"/> Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)																																							
<b>Final Release Date</b> <b>FDLE</b>  <b>MAR 21 2019</b>  <b>Alcohol Testing Program</b>		<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td><u>SD1021</u></td> <td><u>201707D</u> <u>07/25/2019</u></td> </tr> <tr> <td>0.080</td> <td><u>DR1275</u></td> <td><u>201707E</u> <u>07/25/2019</u></td> </tr> <tr> <td>0.200</td> <td><u>SD1011</u></td> <td><u>201707C</u> <u>07/24/2019</u></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td><u>AG805701</u> <u>2/26/20</u></td> </tr> </tbody> </table>		Simulator	Serial #	Lot #/Exp	0.050	<u>SD1021</u>	<u>201707D</u> <u>07/25/2019</u>	0.080	<u>DR1275</u>	<u>201707E</u> <u>07/25/2019</u>	0.200	<u>SD1011</u>	<u>201707C</u> <u>07/24/2019</u>	0.080 DGS	N/A	<u>AG805701</u> <u>2/26/20</u>	<b>Maintenance</b> Performed By _____ <input type="checkbox"/> Battery Replacement <input type="checkbox"/> Dry Gas Regulator Replacement <input type="checkbox"/> Breath Tube Replacement <input type="checkbox"/> Other _____																								
Simulator	Serial #	Lot #/Exp																																									
0.050	<u>SD1021</u>	<u>201707D</u> <u>07/25/2019</u>																																									
0.080	<u>DR1275</u>	<u>201707E</u> <u>07/25/2019</u>																																									
0.200	<u>SD1011</u>	<u>201707C</u> <u>07/24/2019</u>																																									
0.080 DGS	N/A	<u>AG805701</u> <u>2/26/20</u>																																									
<b>Calibration Adjustment</b> Performed By _____ Barometric Pressure Gauge _____ ID # _____		<b>Department Inspection</b> Performed By <u>BK</u> Barometric Pressure ID# <u>28421/28421</u> Gauge <u>1022/1021</u> Instrument <u>1020/1019</u> Mouth Alcohol Solution Lot # <u>2018-B/2018-B</u> Acetone Stock Solution Lot # <u>2018-A/2018-A</u>																																									
<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td></td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.100</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.300</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table>		Simulator	Serial Number	Lot Number	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td><u>G11621/G11621</u></td> </tr> <tr> <td>Interferent</td> <td><u>DR3855/DR3855</u></td> </tr> <tr> <td>0.050</td> <td><u>SD1021/SD1021</u></td> </tr> <tr> <td>0.080</td> <td><u>DR1275/DR1275</u></td> </tr> <tr> <td>0.200</td> <td><u>SD1011/SD1011</u></td> </tr> </tbody> </table>		Simulator	Serial Number	0.000	<u>G11621/G11621</u>	Interferent	<u>DR3855/DR3855</u>	0.050	<u>SD1021/SD1021</u>	0.080	<u>DR1275/DR1275</u>	0.200	<u>SD1011/SD1011</u>
Simulator	Serial Number	Lot Number	Expiration																																								
0.000		N/A	N/A																																								
0.040																																											
0.100																																											
0.200																																											
0.300																																											
0.080 DGS	N/A																																										
Simulator	Serial Number																																										
0.000	<u>G11621/G11621</u>																																										
Interferent	<u>DR3855/DR3855</u>																																										
0.050	<u>SD1021/SD1021</u>																																										
0.080	<u>DR1275/DR1275</u>																																										
0.200	<u>SD1011/SD1011</u>																																										
<input type="checkbox"/> Post Calibration Adjustment Stability Checks		<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> <th>Lot Number</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table>		Simulator	Serial Number	Lot Number	Expiration	0.050				0.080				0.200				0.080 DGS	N/A																						
Simulator	Serial Number	Lot Number	Expiration																																								
0.050																																											
0.080																																											
0.200																																											
0.080 DGS	N/A																																										
Notes/Suggested Service: <u>Performed additional diagnostic after bumped power cord to confirm. BK</u>		<table border="1"> <tr> <td> <input checked="" type="checkbox"/> Form 41 <u>X2</u>  <input checked="" type="checkbox"/> Stability Checks  <input checked="" type="checkbox"/> Calibration Certificate  <input type="checkbox"/> Calibration Adjustment         </td> <td> <input type="checkbox"/> Post-Stability Checks  <input type="checkbox"/> Flow Calibration  <input type="checkbox"/> Form 40  <input checked="" type="checkbox"/> Other <u>Diagnostic</u> </td> </tr> </table>		<input checked="" type="checkbox"/> Form 41 <u>X2</u> <input checked="" type="checkbox"/> Stability Checks <input checked="" type="checkbox"/> Calibration Certificate <input type="checkbox"/> Calibration Adjustment	<input type="checkbox"/> Post-Stability Checks <input type="checkbox"/> Flow Calibration <input type="checkbox"/> Form 40 <input checked="" type="checkbox"/> Other <u>Diagnostic</u>																																						
<input checked="" type="checkbox"/> Form 41 <u>X2</u> <input checked="" type="checkbox"/> Stability Checks <input checked="" type="checkbox"/> Calibration Certificate <input type="checkbox"/> Calibration Adjustment	<input type="checkbox"/> Post-Stability Checks <input type="checkbox"/> Flow Calibration <input type="checkbox"/> Form 40 <input checked="" type="checkbox"/> Other <u>Diagnostic</u>																																										
		<input checked="" type="checkbox"/> Instrument Complies with Chapter 11D-8, FAC <input type="checkbox"/> Instrument Does Not Comply with Chapter 11D-8, FAC <input checked="" type="checkbox"/> Return to/Place into Evidentiary Use <input type="checkbox"/> Remain Out of Evidentiary Use <input checked="" type="checkbox"/> Conduct an Agency Inspection Before Evidentiary Use																																									
		Tech Review / Date <u>Pagan 3/20/19</u> Admin Review / Date <u>J. Deham 3/21/19</u>																																									

# Florida Department of Law Enforcement Alcohol Testing Program

## DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: FWC  
Time of Inspection: 14:03

Date of Inspection: 03/20/2019

Serial Number: 80-003196  
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:201707D Exp: 07/25/2019	0.08g/210L Test (g/210L) Lot#:201707E Exp: 07/25/2019	0.20g/210L Test (g/210L) Lot#:201707C Exp: 07/24/2019	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG805701 Exp: 02/26/2020
0.000	0.048	0.080	0.195	0.082
0.000	0.049	0.081	0.197	0.082
0.000	0.050	0.081	0.198	0.082
0.000	0.050	0.081	0.198	0.082
0.000	0.049	0.081	0.198	0.083
0.000	0.049	0.080	0.199	0.082
0.000	0.049	0.081	0.198	0.082
0.000	0.049	0.081	0.198	0.082
0.000	0.050	0.079	0.198	0.082
0.000	0.050	0.080	0.197	0.082

Standard Deviations	0.0006	0.0007	0.0010	0.0003
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0006 Number of Simulators Used: 5

Remarks:

Pgan

The above instrument complies (  ) does not comply (  ) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

BRETT H KIRKLAND  
 \_\_\_\_\_  
 Signature and Printed Name

03/20/2019  
 Date

3/21/19  
JD

# Florida Department of Law Enforcement Alcohol Testing Program

## DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: FWC  
Time of Inspection: 11:54

Date of Inspection: 03/20/2019

Serial Number: 80-003196  
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK		No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:201707D Exp: 07/25/2019	0.08g/210L Test (g/210L) Lot#:201707E Exp: 07/25/2019	0.20g/210L Test (g/210L) Lot#:201707C Exp: 07/24/2019	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG805701 Exp: 02/26/2020
0.000	0.049	0.080	0.195	0.079
0.000	0.049	0.080	0.198	0.079
0.000	0.049	0.081	0.198	0.078
0.000	0.049	0.080	0.199	0.079
0.000	0.049	0.081	0.198	0.079
0.000	0.049	0.080	0.197	0.078
0.000	0.049	0.081	0.197	0.078
0.000	0.049	0.080	0.197	0.079
0.000	0.049	0.080	0.197	0.078
0.000	0.049	0.080	0.196	0.079

Standard Deviations	0.0000	0.0004	0.0011	0.0005
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0005 Number of Simulators Used: 5

**Remarks:**

BUMPED POWER CORD DURING DIAGNOSTIC CHECK WILL REPEAT INSPECTION. Non-compliance: BUMPED POWER CORD DURING DIAGNOSTIC CHECK

*ADM*

The above instrument complies ( ) does not comply ( X ) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

*Brett Kirkland*

BRETT H KIRKLAND

Signature and Printed Name

03/20/2019  
Date

*3/21/19  
JA*

FWC  
Intoxilyzer - Alcohol Analyzer  
Model: 8000 SN 80-003196  
03/20/2019  
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:43
Control Test	0.078	08:43
Air Blank	0.000	08:43
Control Test	0.079	08:44
Air Blank	0.000	08:44
Control Test	0.078	08:44
Air Blank	0.000	08:45
Control Test Stats		
Average	0.0783	
Std Dev	0.0006	
Rel Std Dev(%)	0.7370	

065

YBK  
Operator's Signature

WBB

80-003196 check  
stability  
YBK 3/20/19

FWC  
Intoxilyzer - Alcohol Analyzer  
Model: 8000 SN 80-003196  
03/20/2019  
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:02
Control Test	0.196	09:02
Air Blank	0.000	09:03
Control Test	0.197	09:04
Air Blank	0.000	09:04
Control Test	0.198	09:05
Air Blank	0.000	09:05
Control Test Stats		
Average	0.1970	
Std Dev	0.0010	
Rel Std Dev(%)	0.5076	

YBK  
Operator's Signature

FWC  
Intoxilyzer - Alcohol Analyzer  
Model: 8000 SN 80-003196  
03/20/2019  
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:54
Control Test	0.080	08:55
Air Blank	0.000	08:56
Control Test	0.081	08:56
Air Blank	0.000	08:57
Control Test	0.081	08:57
Air Blank	0.000	08:58
Control Test Stats		
Average	0.0807	
Std Dev	0.0006	
Rel Std Dev(%)	0.7157	

YBK  
Operator's Signature

FWC  
Intoxilyzer - Alcohol Analyzer  
Model: 8000 SN 80-003196  
03/20/2019  
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:07
Control Test	0.049	09:08
Air Blank	0.000	09:08
Control Test	0.049	09:09
Air Blank	0.000	09:09
Control Test	0.049	09:10
Air Blank	0.000	09:10
Control Test Stats		
Average	0.0490	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

YBK  
Operator's Signature

3/21/19  
J2



# Calibration Certificate

Florida Department of Law Enforcement  
Alcohol Testing Program  
2729 Fort Knox Blvd.  
Bldg. 2, Suite 1300  
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-003196, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-003196</u>	UNCERTAINTY * ±	
Owning Agency:	<u>FWC</u>	0.050 g/ 210 L	0.004
Calibration Date:	<u>03/20/2019</u>	0.080 g/ 210 L	0.004
Calibration Time:	<u>14:03</u>	0.200 g/ 210 L	0.007
		0.080 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.  
Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.  
\*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence (k=3).

### TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Thermometer temperatures are checked with NIST traceable Eutechnics 4400 digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the uses of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

*WSP*

03/20/2019 Date  
*Brett Kirkland*  
BRETT H KIRKLAND,  
Department Inspector

FDLE/ATP Form 69 July 2018  
Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality

*3/21/19  
J2*

FAC  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-003196  
03/20/2019  
Software: 8100.27

-----  
: DIAGNOSTICS :  
-----

Voltage/Current Test	OK
RAM Test	OK
EEPROM Checksum Test	OK
Real Time Clock Test	OK
DSP Test	OK
Analytical Stability Test	OK
Internal Printer Test	OK
Modem Test	OK
Temperature Regulation Test	OK

BK

OPM

3/21/19  
JD