



INSTRUMENT PROCESSING SHEET

Agency Florida Highway PatrolS/N 80-007441Florida Department of
Law EnforcementDate In 3/23/2021 DI Completion Date 4/2/2021☐ Ship ☒ P/U ☐ H/D ☐ CMI ☐ EE

Intake	By <u>RAW</u>	Quality Checks	By <u>IS</u>	Date <u>04/02/2021</u>	Flow Calibration	By	Date																														
☒ Annual ☒ Registration ☐ Return from CMI / EE Visual Inspection: ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/ Accessories: ☒ Power cord ☐ Printer Cable ☒ Static Bag ☒ 12V DC Cable Notes: <u>New Registration</u>		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value <u>185</u> ☒ Flow Verification (L/s) Flow Column # <u>ATP-102</u> 32 mm <u>0.144</u> (.139 - .169) 36 mm <u>0.160</u> (.156 - .190) 53 mm <u>0.230</u> (.228 - .278) 103 mm <u>0.496</u> (.447 - .547) ☒ Barometric Pressure Check Gauge ID # <u>30793</u> ☒ Stability Checks			Flow Column # _____ ☐ 5L/min - 17mm ☐ 15L/min - 53mm ☐ 30L/min - 103mm ☐ R-Value _____ ☐ Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)																																
		<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #/Exp</th></tr></thead><tbody><tr><td>0.050</td><td>MP5088</td><td>202010A 10-05-2022</td></tr><tr><td>0.080</td><td>MP5089</td><td>202010B 10-05-2022</td></tr><tr><td>0.200</td><td>MP5090</td><td>202010D 10-06-2022</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>AG011102 04-20-2022</td></tr></tbody></table>			Simulator	Serial #	Lot #/Exp	0.050	MP5088	202010A 10-05-2022	0.080	MP5089	202010B 10-05-2022	0.200	MP5090	202010D 10-06-2022	0.080 DGS	N/A	AG011102 04-20-2022	<table border="1"><thead><tr><th colspan="2">Maintenance</th><th>By</th></tr></thead><tbody><tr><td>☐ Battery Replacement</td><td></td><td></td></tr><tr><td>☐ Dry Gas Regulator Replacement</td><td></td><td></td></tr><tr><td>☐ Breath Tube Replacement</td><td></td><td></td></tr><tr><td>☐ Other</td><td></td><td></td></tr></tbody></table>			Maintenance		By	☐ Battery Replacement			☐ Dry Gas Regulator Replacement			☐ Breath Tube Replacement			☐ Other		
Simulator	Serial #	Lot #/Exp																																			
0.050	MP5088	202010A 10-05-2022																																			
0.080	MP5089	202010B 10-05-2022																																			
0.200	MP5090	202010D 10-06-2022																																			
0.080 DGS	N/A	AG011102 04-20-2022																																			
Maintenance		By																																			
☐ Battery Replacement																																					
☐ Dry Gas Regulator Replacement																																					
☐ Breath Tube Replacement																																					
☐ Other																																					
		<table border="1"><thead><tr><th colspan="2">DI Temp. Checks</th><th>By <u>IS</u></th></tr></thead><tbody><tr><td>☒ Lab Temp °C</td><td><u>23.70</u></td><td></td></tr><tr><td>External Digital Therm. ID#:</td><td><u>300505</u></td><td></td></tr><tr><td>☒ 34°C +/-2</td><td>Serial #: <u>MP5088</u></td><td></td></tr><tr><td>☒ 34°C +/-2</td><td>Serial #: <u>MP5089</u></td><td></td></tr><tr><td>☒ 34°C +/-2</td><td>Serial #: <u>MP5090</u></td><td></td></tr></tbody></table>			DI Temp. Checks		By <u>IS</u>	☒ Lab Temp °C	<u>23.70</u>		External Digital Therm. ID#:	<u>300505</u>		☒ 34°C +/-2	Serial #: <u>MP5088</u>		☒ 34°C +/-2	Serial #: <u>MP5089</u>		☒ 34°C +/-2	Serial #: <u>MP5090</u>																
DI Temp. Checks		By <u>IS</u>																																			
☒ Lab Temp °C	<u>23.70</u>																																				
External Digital Therm. ID#:	<u>300505</u>																																				
☒ 34°C +/-2	Serial #: <u>MP5088</u>																																				
☒ 34°C +/-2	Serial #: <u>MP5089</u>																																				
☒ 34°C +/-2	Serial #: <u>MP5090</u>																																				
Calibration Adjustment		Department Inspection																																			
By _____		By <u>IS</u>																																			
Barometric Pressure Gauge _____ ID # _____		Barometric Pressure ID# <u>30793</u>																																			
<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.000</td><td></td><td>N/A</td><td>N/A</td></tr><tr><td>0.040</td><td></td><td></td><td></td></tr><tr><td>0.100</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.300</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table>		Simulator	Serial #	Lot #	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			Gauge <u>1027</u> Instrument <u>1027</u> Mouth Alcohol Solution Lot # <u>2021-A</u> Acetone Stock Solution Lot # <u>2020-A</u>							
Simulator	Serial #	Lot #	Expiration																																		
0.000		N/A	N/A																																		
0.040																																					
0.100																																					
0.200																																					
0.300																																					
0.080 DGS	N/A																																				
☐ Post Calibration Adjustment Stability Checks		<table border="1"><thead><tr><th>Simulator</th><th>Serial Number</th></tr></thead><tbody><tr><td>0.000</td><td>MP5086</td></tr><tr><td>Interferent</td><td>MP5087</td></tr><tr><td>0.050</td><td>MP5088</td></tr><tr><td>0.080</td><td>MP5089</td></tr><tr><td>0.200</td><td>MP5090</td></tr></tbody></table>			Simulator	Serial Number	0.000	MP5086	Interferent	MP5087	0.050	MP5088	0.080	MP5089	0.200	MP5090																					
Simulator	Serial Number																																				
0.000	MP5086																																				
Interferent	MP5087																																				
0.050	MP5088																																				
0.080	MP5089																																				
0.200	MP5090																																				
<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.050</td><td></td><td></td><td></td></tr><tr><td>0.080</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table>		Simulator	Serial #	Lot #	Expiration	0.050				0.080				0.200				0.080 DGS	N/A			<table border="1"><thead><tr><th colspan="2">Attachments</th></tr></thead><tbody><tr><td>☒ Form 41</td><td>☐ Post-Stability Checks</td></tr><tr><td>☒ Stability Checks</td><td>☐ Flow Calibration</td></tr><tr><td>☒ Calibration Certificate</td><td>☐ Form 40</td></tr><tr><td>☐ Calibration Adjustment</td><td>☒ Other <u>Form 47</u></td></tr></tbody></table>			Attachments		☒ Form 41	☐ Post-Stability Checks	☒ Stability Checks	☐ Flow Calibration	☒ Calibration Certificate	☐ Form 40	☐ Calibration Adjustment	☒ Other <u>Form 47</u>			
Simulator	Serial #	Lot #	Expiration																																		
0.050																																					
0.080																																					
0.200																																					
0.080 DGS	N/A																																				
Attachments																																					
☒ Form 41	☐ Post-Stability Checks																																				
☒ Stability Checks	☐ Flow Calibration																																				
☒ Calibration Certificate	☐ Form 40																																				
☐ Calibration Adjustment	☒ Other <u>Form 47</u>																																				
Notes/Suggested Service: <u>On alcohol free/mouth alcohol test, accidentally did mouth alcohol test first instead of alcohol free test. IS</u>		☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use																																			
		2021.04.14 David Eliezer Reyes-Rivera Digitally signed by David Eliezer Reyes-Rivera Date: 2021.04.05 09:32:18 -04'00'																																			
		Tech Review / Date _____ Admin Review / Date _____																																			

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: FHP

Time of Inspection: 12:13

Date of Inspection: 04/02/2021

Serial Number: 80-007441

Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202010A Exp: 10/05/2022	0.08g/210L Test (g/210L) Lot#:202010B Exp: 10/05/2022	0.20g/210L Test (g/210L) Lot#:202010D Exp: 10/06/2022	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG011102 Exp: 04/20/2022
0.000	0.050	0.081	0.202	0.081
0.000	0.050	0.081	0.201	0.081
0.000	0.050	0.081	0.202	0.080
0.000	0.050	0.081	0.202	0.080
0.000	0.050	0.081	0.202	0.080
0.000	0.050	0.081	0.202	0.080
0.000	0.050	0.081	0.202	0.081
0.000	0.050	0.081	0.201	0.081
0.000	0.050	0.081	0.202	0.080
0.000	0.050	0.081	0.202	0.080
Standard Deviations	0.0000	0.0000	0.0004	0.0005

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0002 Number of Simulators Used: 5

Remarks:

A F / M A: DID MA FIRST INSTEAD OF AF.

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



ISRAEL SOTO

Signature and Printed Name

04/02/2021
Date

Stability Checks

FHP
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-007441
04/02/2021
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:29
Control Test	0.050	09:29
Air Blank	0.000	09:30
Control Test	0.050	09:31
Air Blank	0.000	09:31
Control Test	0.050	09:32
Air Blank	0.000	09:32
Control Test Stats		
Average	0.0500	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

Operator's Signature

FHP
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-007441
04/02/2021
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:34
Control Test	0.081	09:34
Air Blank	0.000	09:35
Control Test	0.080	09:36
Air Blank	0.000	09:36
Control Test	0.080	09:37
Air Blank	0.000	09:38
Control Test Stats		
Average	0.0803	
Std Dev	0.0006	
Rel Std Dev(%)	0.7187	

Operator's Signature

FHP
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-007441
04/02/2021
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:40
Control Test	0.202	09:41
Air Blank	0.000	09:41
Control Test	0.200	09:42
Air Blank	0.000	09:42
Control Test	0.200	09:43
Air Blank	0.000	09:44
Control Test Stats		
Average	0.2007	
Std Dev	0.0012	
Rel Std Dev(%)	0.5754	

Operator's Signature

FHP
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-007441
04/02/2021
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:48
Control Test	0.081	09:48
Air Blank	0.000	09:49
Control Test	0.081	09:49
Air Blank	0.000	09:50
Control Test	0.080	09:50
Air Blank	0.000	09:51
Control Test Stats		
Average	0.0807	
Std Dev	0.0006	
Rel Std Dev(%)	0.7157	

Operator's Signature

wet

Dry



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
2729 Fort Knox Blvd.
Bldg. 2, Suite 1300
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-007441, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-007441</u>	UNCERTAINTY* $\pm$	
Owning Agency:	<u>FHP</u>	0.050 g/ 210 L	0.005
Calibration Date:	<u>04/02/2021</u>	0.080 g/ 210 L	0.004
Calibration Time:	<u>12:13</u>	0.200 g/ 210 L	0.007
		0.080 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.
*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Thermometer temperatures are checked with NIST traceable Eutechnics 4400 digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the uses of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

 Digitally signed by Israel Soto
Date: 2021.04.05 07:56:20
-04'00'

Israel Soto

04/02/2021

Date

ISRAEL SOTO,

Department Inspector

FDLE/ATP Form 69 January 2021

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality

Page 1 of 1



Florida Department of
Law Enforcement

REQUEST FOR REGISTRATION

MAKE AND MODEL OF INSTRUMENT: Intoxilyzer 8000

SERIAL NUMBER: 80-007441

OWNING AGENCY: Florida Highway Patrol

DATE OF DEPARTMENT INSPECTION: 04-02-2021

AGENCY INSPECTOR: Susan Barge

ADDRESS: 6030 CR 2321

CITY, STATE, ZIP: Panama City, FL. 32404

TELEPHONE NUMBER: 352-620-4701

FAX NUMBER: _____

EMAIL ADDRESS (if available): Susanbarga@flhsmv.gov

For Program Office Use Only:

- 3/24/21
- ☒ Registration Issued
 - ☒ Instrument Added to Evidentiary Instrument Database
 - ☒ Instrument Added to Monthly Statistics Database
 - ☒ Contact Information Added to Instrument Database