



INSTRUMENT PROCESSING SHEET

Agency Okaloosa County SOS/N 80-000763Florida Department of
Law EnforcementDate In 12-6-2023DI Completion Date 12-6-2023☐ Ship ☒ P/U ☐ H/D ☐ CMI ☐ EE

Intake	By PN	Quality Checks	By IS	Date <u>12-6-2023</u>	Flow Calibration	By IS	Date <u>12-6-2023</u>																																													
☒ Annual ☐ Registration ☐ Return from CMI / EE Visual Inspection: ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/ Accessories: ☐ Power cord ☐ Printer Cable ☐ Static Bag ☐ 12V DC Cable Notes: _____ _____ _____ _____ _____ _____ _____ _____ _____ _____		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value <u>204</u> ☒ Flow Verification (L/s) Flow Column # <u>ATP-105</u> 32 mm <u>0.171</u> (.139 - .169) 36 mm <u>0.191</u> (.156 - .190) 53 mm <u>0.261</u> (.228 - .278) 103 mm <u>0.535</u> (.447 - .547) ☒ Barometric Pressure Check Gauge ID # <u>28421</u> ☒ Stability Checks			☒ Flow Column # <u>ATP-102</u> ☒ 5L/min - 17mm ☒ 15L/min - 53mm ☒ 30L/min - 103mm ☒ R-Value <u>206</u> ☒ Post Calibration Verification (L/s) Flow Column # <u>ATP-105</u> 32 mm <u>0.148</u> (.139 - .169) 36 mm <u>0.164</u> (.156 - .190) 53 mm <u>0.234</u> (.228 - .278) 103 mm <u>0.523</u> (.447 - .547)																																															
		<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #/Exp</th></tr></thead><tbody><tr><td>0.050</td><td>MP6291</td><td>202303K 03-29-2025</td></tr><tr><td>0.080</td><td>MP6292</td><td>202303L 03-29-2025</td></tr><tr><td>0.200</td><td>MP6293</td><td>202304C 04-05-2025</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>AG229803 10-25-2024</td></tr></tbody></table>			Simulator	Serial #	Lot #/Exp	0.050	MP6291	202303K 03-29-2025	0.080	MP6292	202303L 03-29-2025	0.200	MP6293	202304C 04-05-2025	0.080 DGS	N/A	AG229803 10-25-2024	<table border="1"><thead><tr><th colspan="2">Maintenance</th><th>By _____</th></tr></thead><tbody><tr><td>☐ Battery Replacement</td><td colspan="2"></td></tr><tr><td>☐ Dry Gas Regulator Replacement</td><td colspan="2"></td></tr><tr><td>☐ Breath Tube Replacement</td><td colspan="2"></td></tr><tr><td>☐ Other _____</td><td colspan="2"></td></tr><tr><td>_____</td><td colspan="2"></td></tr><tr><td>_____</td><td colspan="2"></td></tr><tr><td>_____</td><td colspan="2"></td></tr><tr><td>_____</td><td colspan="2"></td></tr><tr><td>_____</td><td colspan="2"></td></tr></tbody></table>			Maintenance		By _____	☐ Battery Replacement			☐ Dry Gas Regulator Replacement			☐ Breath Tube Replacement			☐ Other _____			_____			_____			_____			_____			_____		
Simulator	Serial #	Lot #/Exp																																																		
0.050	MP6291	202303K 03-29-2025																																																		
0.080	MP6292	202303L 03-29-2025																																																		
0.200	MP6293	202304C 04-05-2025																																																		
0.080 DGS	N/A	AG229803 10-25-2024																																																		
Maintenance		By _____																																																		
☐ Battery Replacement																																																				
☐ Dry Gas Regulator Replacement																																																				
☐ Breath Tube Replacement																																																				
☐ Other _____																																																				

Calibration Adjustment By IS _____				Department Inspection By IS _____																																																
Barometric Pressure Gauge <u>1023</u> ID # <u>28663</u>				Barometric Pressure ID# <u>28427</u>																																																
<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.000</td><td>MP5091</td><td>N/A</td><td>N/A</td></tr><tr><td>0.040</td><td>MP5082</td><td>22460</td><td>12-28-2024</td></tr><tr><td>0.100</td><td>MP5086</td><td>22430</td><td>11-30-2024</td></tr><tr><td>0.200</td><td>MP5084</td><td>22400</td><td>10-12-2024</td></tr><tr><td>0.300</td><td>MP5085</td><td>23070</td><td>03-06-2025</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>06723080A5</td><td>04-05-2025</td></tr></tbody></table>				Simulator	Serial #	Lot #	Expiration	0.000	MP5091	N/A	N/A	0.040	MP5082	22460	12-28-2024	0.100	MP5086	22430	11-30-2024	0.200	MP5084	22400	10-12-2024	0.300	MP5085	23070	03-06-2025	0.080 DGS	N/A	06723080A5	04-05-2025	Gauge <u>1018</u> Instrument <u>1021</u> Mouth Alcohol Solution Lot # <u>2023-A</u> Acetone Stock Solution Lot # <u>2022-B</u>																				
Simulator	Serial #	Lot #	Expiration																																																	
0.000	MP5091	N/A	N/A																																																	
0.040	MP5082	22460	12-28-2024																																																	
0.100	MP5086	22430	11-30-2024																																																	
0.200	MP5084	22400	10-12-2024																																																	
0.300	MP5085	23070	03-06-2025																																																	
0.080 DGS	N/A	06723080A5	04-05-2025																																																	
☒ Post Calibration Adjustment Stability Checks				<table border="1"><thead><tr><th>Simulator</th><th>Serial Number</th></tr></thead><tbody><tr><td>0.000</td><td>MP6289</td></tr><tr><td>Interferent</td><td>MP6290</td></tr><tr><td>0.050</td><td>MP6291</td></tr><tr><td>0.080</td><td>MP6292</td></tr><tr><td>0.200</td><td>MP6293</td></tr></tbody></table>				Simulator	Serial Number	0.000	MP6289	Interferent	MP6290	0.050	MP6291	0.080	MP6292	0.200	MP6293																																	
Simulator	Serial Number																																																			
0.000	MP6289																																																			
Interferent	MP6290																																																			
0.050	MP6291																																																			
0.080	MP6292																																																			
0.200	MP6293																																																			
<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.050</td><td>MP6291</td><td>202303K</td><td>03-29-2025</td></tr><tr><td>0.080</td><td>MP6292</td><td>202303L</td><td>03-29-2025</td></tr><tr><td>0.200</td><td>MP6293</td><td>202304C</td><td>04-05-2025</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>AG229803</td><td>10-25-2024</td></tr></tbody></table>				Simulator	Serial #	Lot #	Expiration	0.050	MP6291	202303K	03-29-2025	0.080	MP6292	202303L	03-29-2025	0.200	MP6293	202304C	04-05-2025	0.080 DGS	N/A	AG229803	10-25-2024	Attachments ☒ Form 41 ☒ Post-Stability Checks ☒ Stability Checks ☒ Flow Calibration ☒ Calibration Certificate ☒ Form 40 IS 12-6-2023 ☒ Calibration Adjustment ☐ Other _____																												
Simulator	Serial #	Lot #	Expiration																																																	
0.050	MP6291	202303K	03-29-2025																																																	
0.080	MP6292	202303L	03-29-2025																																																	
0.200	MP6293	202304C	04-05-2025																																																	
0.080 DGS	N/A	AG229803	10-25-2024																																																	
Notes/Suggested Service: <u>Agency Inspector indicated November Agency Inspection had Control Outside Tolerance and repeated Alcohol Free Test. Optical Bench Calibration Adjustment performed. Internal printer not printing full printouts, external printer used to print documents. IS</u> <u>Tech Review: Checked Form 40 under Attachments. IS 12-6-2023</u> _____ _____				☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use																																																
Benjamin Siddoway Digitally signed by Benjamin Siddoway Date: 2023.12.06 15:18:07 -05'00'				Phil Nicodemo Digitally signed by Phil Nicodemo Date: 2023.12.06 15:44:42 -05'00'																																																
Tech Review / Date _____				Admin Review / Date _____																																																

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: OKALOOSA COUNTY SO
Time of Inspection: 13:55

Date of Inspection: 12/06/2023

Serial Number: 80-000763
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202303K Exp: 03/29/2025	0.08g/210L Test (g/210L) Lot#:202303L Exp: 03/29/2025	0.20g/210L Test (g/210L) Lot#:202304C Exp: 04/05/2025	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG229803 Exp: 10/25/2024
0.000	0.049	0.079	0.198	0.078
0.000	0.049	0.078	0.198	0.078
0.000	0.049	0.078	0.198	0.078
0.000	0.049	0.079	0.198	0.078
0.000	0.049	0.078	0.198	0.078
0.000	0.049	0.078	0.198	0.078
0.000	0.049	0.079	0.198	0.078
0.000	0.049	0.079	0.198	0.078
0.000	0.050	0.079	0.198	0.078
0.000	0.049	0.079	0.198	0.078
Standard Deviations	0.0003	0.0005	0.0000	0.0000

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0002 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Israel Soto

ISRAEL SOTO

Signature and Printed Name

12/06/2023
Date

Stability checks

OKALOOSA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000763
12/06/2023
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:06
Control Test	0.051	09:07
Air Blank	0.000	09:07
Control Test	0.050	09:08
Air Blank	0.000	09:08
Control Test	0.051	09:09
Air Blank	0.000	09:10
Control Test Stats		
Average	0.0507	
Std Dev	0.0006	
Rel Std Dev(%)	1.1395	

OKALOOSA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000763
12/06/2023
Software: 8100.27

Test g/210L Time

.050 #1

.050 #2

Internal Printer
did not print
stability check
results. Connected
external printer
and continued
with stability checks

✓

✓

Operator's Signature

stability checks

OKALOOSA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000763
12/06/2023
Software: 8100.27

Test	g/210L	Time
-----	-----	-----
Air Blank	0.000	09:10
Control Test	0.080	09:11
Air Blank	0.000	09:12
Control Test	0.080	09:12
Air Blank	0.000	09:13
Control Test	0.080	09:14
Air Blank	0.000	09:14
Control Test Stats		
Average	0.0800	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

wet



Operator's Signature

Stability checks

OKALOOSA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000763
12/06/2023
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:15
Control Test	0.199	09:16
Air Blank	0.000	09:17
Control Test	0.199	09:17
Air Blank	0.000	09:18
Control Test	0.199	09:18
Air Blank	0.000	09:19
Control Test Stats		
Average	0.1990	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	



Operator's Signature

Stability Checks

OKALOOSA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000763
12/06/2023
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:20
Control Test	0.080	09:21
Air Blank	0.000	09:21
Control Test	0.080	09:21
Air Blank	0.000	09:22
Control Test	0.081	09:22
Air Blank	0.000	09:23
Control Test Stats		
Average	0.0803	
Std Dev	0.0006	
Rel Std Dev(%)	0.7187	

Dry



Operator's Signature



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
2331 Phillips Road.
Suite B1032
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-000763, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-000763</u>	UNCERTAINTY* $\pm$	
Owning Agency:	<u>OKALOOSA COUNTY SO</u>	0.050 g/ 210 L	0.004
Calibration Date:	<u>12/06/2023</u>	0.080 g/ 210 L	0.004
Calibration Time:	<u>13:55</u>	0.200 g/ 210 L	0.007
		0.080 g/ 210 L Dry Gas Control	0.005

All results are reported in g/210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

Israel Soto
Soto
Digitally signed by Israel Soto
Date: 2023.12.06 15:11:44 -05'00'

12/06/2023

Date

ISRAEL SOTO,
Department Inspector

FDLE/ATP Form 69 March 2022
Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality

OKALOOSA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000763
12/06/2023 09:27:33

Optical Bench
calibration Adjustment

✓

Auto Calibration

pg 1 of 2

<<<<< 3um >>>>>			<<<<< 9um >>>>>		

Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1					
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)	
Sample #1	0.0630	(-0.0080)	0.1390	(0.0020)	
Sample #2	0.0870	(0.0000)	0.1480	(0.0190)	
Sample #3	0.1100	(-0.0020)	0.1640	(0.0170)	
Sample #4	0.1410	(0.0020)	0.1790	(0.0100)	
Avg % Abs	0.1127	(0.0000)	0.1637	(0.0153)	
STD DEV	0.0271	(0.0020)	0.0155	(0.0047)	
REL STD DEV	24.052	(0.000)	9.472	(30.821)	

Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1					
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)	
Sample #1	0.8010	(-0.0130)	1.4990	(0.0020)	
Sample #2	0.8370	(-0.0330)	1.5280	(-0.0210)	
Sample #3	0.8060	(-0.0150)	1.4930	(-0.0020)	
Sample #4	0.7910	(0.0000)	1.5040	(0.0000)	
Avg % Abs	0.8113	(-0.0160)	1.5083	(-0.0077)	
STD DEV	0.0235	(0.0165)	0.0179	(0.0116)	
REL STD DEV	2.891	(103.267)	1.187	(151.177)	

Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1					
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)	
Sample #1	1.8800	(-0.0340)	3.5210	(-0.0170)	
Sample #2	1.8620	(0.0050)	3.4980	(0.0030)	
Sample #3	1.8840	(-0.0060)	3.5130	(-0.0090)	
Sample #4	1.8680	(0.0260)	3.5000	(0.0160)	
Avg % Abs	1.8713	(0.0083)	3.5037	(0.0033)	
STD DEV	0.0114	(0.0163)	0.0081	(0.0125)	
REL STD DEV	0.608	(195.100)	0.232	(375.100)	

Solution = 0.200 g/210L or 0.9524 mg/l, Samples = 4, Discarded = 1					
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)	
Sample #1	3.5630	(-0.0150)	6.6010	(-0.0090)	
Sample #2	3.5520	(-0.0010)	6.6410	(-0.0050)	
Sample #3	3.5900	(0.0090)	6.6800	(-0.0090)	
Sample #4	3.6040	(0.0110)	6.6750	(-0.0120)	
Avg % Abs	3.5820	(0.0063)	6.6653	(-0.0087)	
STD DEV	0.0269	(0.0064)	0.0212	(0.0035)	
REL STD DEV	0.751	(101.512)	0.318	(40.522)	

Solution = 0.300 g/210L or 1.4286 mg/l, Samples = 4, Discarded = 1					
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)	
Sample #1	5.2920	(-0.0100)	9.6900	(-0.0170)	
Sample #2	5.2580	(0.0220)	9.6960	(0.0040)	
Sample #3	5.2680	(0.0250)	9.6880	(0.0230)	
Sample #4	5.2710	(0.0570)	9.6980	(0.0280)	
Avg % Abs	5.2657	(0.0347)	9.6940	(0.0183)	
STD DEV	0.0068	(0.0194)	0.0053	(0.0127)	
REL STD DEV	0.129	(55.960)	0.055	(69.067)	

OKALOOSA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000763
12/06/2023 09:27:33

Optical Bench Calibration Adjustment

Auto Calibration

✓

pg 2 of 2

<<<< 3um >>>>

Zero Order Coef -294.94
First Order Coef 2676.76
Second Order Coef 17.55

<<<< 9um >>>>

-221.87
1384.68
11.54

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	0.000	-0.0001
0.040	0.040	0.0003
0.100	0.100	-0.0003
0.200	0.200	0.0001
0.300	0.300	-0.0000

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	0.000	-0.0001
0.040	0.040	0.0002
0.100	0.100	-0.0002
0.200	0.200	0.0001
0.300	0.300	0.0000

<<<< 3um >>>>

<<<< 9um >>>>

Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1

Sample

Sample #1	3109.00
Sample #2	3147.00
Sample #3	3100.00
Sample #4	2985.00
Avg	3077.3333
STD DEV	83.3447
REL STD DEV	2.708
H2O adjust (mg/l*10k)	732

3313.00
3324.00
3273.00
3270.00
3289.0000
30.3480
0.923
520

Barometric Pressure = 1022

*****CALIBRATION SUCCESSFUL*****

Post Stability Checks

OKALOOSA COUNTY SO

Intoxilyzer - Alcohol Analyzer

Model 8000

SN 80-000763

12/06/2023

Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:14
Control Test	0.050	11:14
Air Blank	0.000	11:15
Control Test	0.049	11:15
Air Blank	0.000	11:16
Control Test	0.049	11:17
Air Blank	0.000	11:17
Control Test Stats		
Average	0.0493	
Std Dev	0.0006	
Rel Std Dev(%)	1.1703	



Operator's Signature

Post Stability Checks

OKALOOSA COUNTY SO

Intoxilyzer - Alcohol Analyzer

Model 8000

SN 80-000763

12/06/2023

Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:19
Control Test	0.080	11:20
Air Blank	0.000	11:20
Control Test	0.079	11:21
Air Blank	0.000	11:22
Control Test	0.079	11:22
Air Blank	0.000	11:23
Control Test Stats		
Average	0.0793	
Std Dev	0.0006	
Rel Std Dev(%)	0.7277	

wet



Operator's Signature

Post stability checks

OKALOOSA COUNTY SO

Intoxilyzer - Alcohol Analyzer

Model 8000

SN 80-000763

12/06/2023

Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:25
Control Test	0.198	11:26
Air Blank	0.000	11:26
Control Test	0.197	11:27
Air Blank	0.000	11:27
Control Test	0.197	11:28
Air Blank	0.000	11:29
Control Test Stats		
Average	0.1973	
Std Dev	0.0006	
Rel Std Dev(%)	0.2926	



Operator's Signature

Post Stability Checks

OKALOOSA COUNTY SO

Intoxilyzer - Alcohol Analyzer

Model 8000

SN 80-000763

12/06/2023

Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:30
Control Test	0.079	11:31
Air Blank	0.000	11:31
Control Test	0.080	11:31
Air Blank	0.000	11:32
Control Test	0.079	11:32
Air Blank	0.000	11:33
Control Test Stats		
Average	0.0793	
Std Dev	0.0006	
Rel Std Dev(%)	0.7277	

Dry



Operator's Signature

Flow
calibration
Adjustment

✓

OKALOOSA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000763
12/06/2023
Software: 8100.27

Flow Rate Calibration*****

1: Rate (Liters/min) = 5

SQRT(Diff)) = 7.414

2: Rate (Liters/min) = 15

SQRT(Diff)) = 11.746

3: Rate (Liters/min) = 30

SQRT(Diff)) = 20.711

Dependent Data Scale Factor = 100000 L/min

Independent Data Scale Factor = 256

Rounded Slope = 722

Rounded Intercept = -790713

Correlation = 0.99651

Florida Department of Law Enforcement Alcohol Testing Program

AGENCY INSPECTION REPORT - INTOXILYZER 8000

Agency: OKALOOSA COUNTY SO
Time of Inspection: 08:15

Date of Inspection: 12/06/2023

Serial Number: 80-000763
Software: 8100.27

Check or Test	YES	NO
Date and/or Time Adjusted		No
Diagnostic Check (Pre-Inspection): OK		No
Alcohol Free Subject Test: 0.000		No
Mouth Alcohol Test: Slope Not Met		No
Interferent Detect Test: Interferent Detect		No
Diagnostic Check (Post-Inspection): OK		No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#: Exp:	0.08g/210L Test (g/210L) Lot#: Exp:	0.20g/210L Test (g/210L) Lot#: Exp:	0.08 g/210L Dry Gas Std Test (g/210L) Lot#: Exp:

Number of Simulators Used: _____

Remarks:

BYPASSED AI TO OPERATE INSTRUMENT. COMPLIANCE NOT DETERMINED

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I hold a valid Florida Department of Law Enforcement Agency Inspector Permit and that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



PHIL NICODEMO

Signature and Printed Name

12/06/2023
Date



Agency Okaloosa County SO

S/N 80-000763

Date In 11-7-2023 DI Completion Date 11-9-2023

☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake		By IS		Quality Checks		By IS		Date 11-9-2023		Flow Calibration		By		Date																																
☒ Annual ☐ Registration ☐ Return from CMI / EE Visual Inspection: ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/ Accessories: ☐ Power cord ☐ Printer Cable ☐ Static Bag ☐ 12V DC Cable Notes: _____ _____ _____ _____ _____ _____ _____ _____ _____ _____				☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value <u>206</u> ☒ Flow Verification (L/s) Flow Column # <u>ATP-105</u> 32 mm <u>0.167</u> (.139 - .169) 36 mm <u>0.183</u> (.156 - .190) 53 mm <u>0.253</u> (.228 - .278) 103 mm <u>0.535</u> (.447 - .547) ☒ Barometric Pressure Check Gauge ID # <u>30793</u> ☒ Stability Checks				Flow Column # _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value _____ ☐ Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)																																						
				<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td rowspan="2">0.050</td> <td rowspan="2">MP6291</td> <td>202303K</td> </tr> <tr> <td>03-29-2025</td> </tr> <tr> <td rowspan="2">0.080</td> <td rowspan="2">MP6292</td> <td>202303L</td> </tr> <tr> <td>03-29-2025</td> </tr> <tr> <td rowspan="2">0.200</td> <td rowspan="2">MP6293</td> <td>202304C</td> </tr> <tr> <td>04-05-2025</td> </tr> <tr> <td rowspan="2">0.080 DGS</td> <td rowspan="2">N/A</td> <td>AG229803</td> </tr> <tr> <td>10-25-2024</td> </tr> </tbody> </table>		Simulator	Serial #	Lot #/Exp	0.050	MP6291	202303K	03-29-2025	0.080	MP6292	202303L	03-29-2025	0.200	MP6293	202304C	04-05-2025	0.080 DGS	N/A	AG229803	10-25-2024			<table border="1"> <thead> <tr> <th>Maintenance</th> <th>By</th> </tr> </thead> <tbody> <tr> <td>☐ Battery Replacement</td> <td></td> </tr> <tr> <td>☐ Dry Gas Regulator Replacement</td> <td></td> </tr> <tr> <td>☐ Breath Tube Replacement</td> <td></td> </tr> <tr> <td>☐ Other _____</td> <td></td> </tr> <tr> <td>_____</td> <td></td> </tr> <tr> <td>_____</td> <td></td> </tr> <tr> <td>_____</td> <td></td> </tr> <tr> <td>_____</td> <td></td> </tr> </tbody> </table>		Maintenance	By	☐ Battery Replacement		☐ Dry Gas Regulator Replacement		☐ Breath Tube Replacement		☐ Other _____		_____		_____		_____		_____	
Simulator	Serial #	Lot #/Exp																																												
0.050	MP6291	202303K																																												
		03-29-2025																																												
0.080	MP6292	202303L																																												
		03-29-2025																																												
0.200	MP6293	202304C																																												
		04-05-2025																																												
0.080 DGS	N/A	AG229803																																												
		10-25-2024																																												
Maintenance	By																																													
☐ Battery Replacement																																														
☐ Dry Gas Regulator Replacement																																														
☐ Breath Tube Replacement																																														
☐ Other _____																																														

Calibration Adjustment				By _____	
Barometric Pressure Gauge _____		ID # _____			
Simulator	Serial #	Lot #	Expiration		
0.000		N/A	N/A		
0.040					
0.100					
0.200					
0.300					
0.080 DGS	N/A				
☐ Post Calibration Adjustment Stability Checks					
Simulator	Serial #	Lot #	Expiration		
0.050					
0.080					
0.200					
0.080 DGS	N/A				
Notes/Suggested Service: Internal printer not printing. External printer used to print Stability Checks. IS					

Department Inspection		By IS _____	
Barometric Pressure ID# <u>30793</u>			
Gauge <u>1018</u>	Instrument <u>1017</u>		
Mouth Alcohol Solution Lot # <u>2023-A</u>			
Acetone Stock Solution Lot # <u>2022-B</u>			
Simulator	Serial Number		
0.000	MP6289		
Interferent	MP6290		
0.050	MP6291		
0.080	MP6292		
0.200	MP6293		
Attachments			
☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☐ Calibration Adjustment		☐ Post-Stability Checks ☐ Flow Calibration ☒ Form 40 ☐ Other _____	
☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC			
☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use			
☒ Conduct an Agency Inspection Before Evidentiary Use			
Taylor Gutschow Digitally signed by Taylor Gutschow Date: 2023.11.09 14:01:15 -05'00'		Phil Nicodemo Digitally signed by Phil Nicodemo Date: 2023.11.15 10:28:56 -05'00'	
Tech Review / Date _____		Admin Review / Date _____	

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: OKALOOSA COUNTY SO
Time of Inspection: 10:02

Date of Inspection: 11/09/2023

Serial Number: 80-000763
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202303K Exp: 03/29/2025	0.08g/210L Test (g/210L) Lot#:202303L Exp: 03/29/2025	0.20g/210L Test (g/210L) Lot#:202304C Exp: 04/05/2025	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG229803 Exp: 10/25/2024
0.000	0.050	0.080	0.199	0.079
0.000	0.050	0.080	0.200	0.080
0.000	0.050	0.080	0.199	0.080
0.000	0.050	0.080	0.199	0.079
0.000	0.050	0.080	0.199	0.080
0.000	0.051	0.081	0.199	0.079
0.000	0.050	0.080	0.199	0.079
0.000	0.051	0.080	0.199	0.079
0.000	0.050	0.080	0.199	0.079
0.000	0.051	0.081	0.199	0.080
Standard Deviations	0.0004	0.0004	0.0003	0.0005

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0004 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Israel Soto

ISRAEL SOTO

Signature and Printed Name

11/09/2023
Date

stability checks

OKALOOSA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000763
11/09/2023
Software: 8100.27

Test	g/210L	Time

Air Blank	0.000	07:57
Control Test	0.050	07:58
Air Blank	0.000	07:58
Control Test	0.050	07:59
Air Blank	0.000	08:00
Control Test	0.051	08:00
Air Blank	0.000	08:01
Control Test Stats		
Average	0.0503	
Std Dev	0.0006	
Rel Std Dev(%)	1.1471	



Operator's Signature

stability checks

OKALOOSA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000763
11/09/2023
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:02
Control Test	0.081	08:02
Air Blank	0.000	08:03
Control Test	0.081	08:04
Air Blank	0.000	08:04
Control Test	0.080	08:05
Air Blank	0.000	08:05
Control Test Stats		
Average	0.0807	
Std Dev	0.0006	
Rel Std Dev(%)	0.7157	

wet



Operator's Signature

stability checks

OKALOOSA COUNTY SO

Intoxilyzer - Alcohol Analyzer

Model 8000

SN 80-000763

11/09/2023

Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:07
Control Test	0.199	08:08
Air Blank	0.000	08:09
Control Test	0.199	08:09
Air Blank	0.000	08:10
Control Test	0.199	08:10
Air Blank	0.000	08:11
Control Test Stats		
Average	0.1990	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	



Operator's Signature

stability checks

OKALOOSA COUNTY SO

Intoxilyzer - Alcohol Analyzer

Model 8000

SN 80-000763

11/09/2023

Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:20
Control Test	0.079	08:20
Air Blank	0.000	08:21
Control Test	0.079	08:21
Air Blank	0.000	08:22
Control Test	0.080	08:22
Air Blank	0.000	08:22
Control Test Stats		
Average	0.0793	
Std Dev	0.0006	
Rel Std Dev(%)	0.7277	

Dry



Operator's Signature



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
2331 Phillips Road.
Suite B1032
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-000763 , manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-000763</u>	UNCERTAINTY* $\pm$
Owning Agency:	<u>OKALOOSA COUNTY SO</u>	0.050 g/ 210 L 0.004
Calibration Date:	<u>11/09/2023</u>	0.080 g/ 210 L 0.004
Calibration Time:	<u>10:02</u>	0.200 g/ 210 L 0.007
		0.080 g/ 210 L Dry Gas Control 0.005

All results are reported in g/210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.
This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

Israel Soto
Soto
Digitally signed by Israel Soto
Date: 2023.11.09 13:08:31 -05'00'

11/09/2023

Date

ISRAEL SOTO,
Department Inspector

FDLE/ATP Form 69 March 2022
Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality

Florida Department of Law Enforcement Alcohol Testing Program

AGENCY INSPECTION REPORT - INTOXILYZER 8000

Agency: OKALOOSA COUNTY SO

Time of Inspection: 07:19

Date of Inspection: 11/09/2023

Serial Number: 80-000763

Software: 8100.27

Check or Test	YES	NO
Date and/or Time Adjusted		No
Diagnostic Check (Pre-Inspection): OK		No
Alcohol Free Subject Test: 0.000		No
Mouth Alcohol Test: Slope Not Met		No
Interferent Detect Test: Interferent Detect		No
Diagnostic Check (Post-Inspection): OK		No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#: Exp:	0.08g/210L Test (g/210L) Lot#: Exp:	0.20g/210L Test (g/210L) Lot#: Exp:	0.08 g/210L Dry Gas Std Test (g/210L) Lot#: Exp:

Number of Simulators Used: _____

Remarks:

BYPASSED AI TO OPERATE INSTRUMENT

N/A Not a compliance check

The above instrument complies (☒) does not comply () with Chapter 11D-8, FAC.

I certify that I hold a valid Florida Department of Law Enforcement Agency Inspector Permit and that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Israel Soto

ISRAEL SOTO

Signature and Printed Name

11/09/2023
Date