



INSTRUMENT PROCESSING SHEET

Agency Gulfport Police DepartmentS/N 80-001061Florida Department of
Law EnforcementDate In 06/19/2025 DI Completion Date 07/17/2025☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake	By <u>TDG</u>	Date <u>07/11/2025</u>	Quality Checks	By <u>TDG</u>	Date <u>07/15/2025</u>	Flow Calibration	By _____	Date _____																																								
☒ Annual ☐ Registration ☐ Return from CMI / EE Visual Inspection: ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/ Accessories: ☒ Power cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes: <u>Agency Inspector Michael Marotta reports a Modem Fail.</u>			☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value <u>230</u> ☒ Flow Verification (L/s) Flow Column # <u>ATP104</u> 32 mm <u>0.156</u> (.139 - .169) 36 mm <u>0.175</u> (.156 - .190) 53 mm <u>0.253</u> (.228 - .278) 103 mm <u>0.496</u> (.447 - .547) ☒ Barometric Pressure Check Gauge ID # <u>28199</u> ☒ Stability Checks			☐ Flow Column # _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value _____ ☐ Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)																																										
			<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #/Exp</th></tr></thead><tbody><tr><td>0.050</td><td>MP6286</td><td>202406K 06/19/2026</td></tr><tr><td>0.080</td><td>MP6287</td><td>202406L 06/19/2026</td></tr><tr><td>0.200</td><td>MP6288</td><td>202406N 06/20/2026</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>AG429602 10/22/2026</td></tr></tbody></table>			Simulator	Serial #	Lot #/Exp	0.050	MP6286	202406K 06/19/2026	0.080	MP6287	202406L 06/19/2026	0.200	MP6288	202406N 06/20/2026	0.080 DGS	N/A	AG429602 10/22/2026	<table border="1"><thead><tr><th colspan="2">Maintenance By _____ Date _____</th></tr><tr><td>☐ Battery Replacement</td><td></td></tr><tr><td>☐ Dry Gas Regulator Replacement</td><td></td></tr><tr><td>☐ Breath Tube Replacement</td><td></td></tr><tr><td>☐ Other _____</td><td></td></tr></thead></table>			Maintenance By _____ Date _____		☐ Battery Replacement		☐ Dry Gas Regulator Replacement		☐ Breath Tube Replacement		☐ Other _____																
Simulator	Serial #	Lot #/Exp																																														
0.050	MP6286	202406K 06/19/2026																																														
0.080	MP6287	202406L 06/19/2026																																														
0.200	MP6288	202406N 06/20/2026																																														
0.080 DGS	N/A	AG429602 10/22/2026																																														
Maintenance By _____ Date _____																																																
☐ Battery Replacement																																																
☐ Dry Gas Regulator Replacement																																																
☐ Breath Tube Replacement																																																
☐ Other _____																																																
Calibration Adjustment By _____			Department Inspection By <u>TDG</u>																																													
Barometric Pressure Gauge _____ ID # _____			Barometric Pressure ID# <u>28199</u> Gauge <u>1019</u> Instrument <u>1017</u> Mouth Alcohol Solution Lot # <u>2024-A</u> Acetone Stock Solution Lot # <u>2023-B</u>																																													
<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.000</td><td></td><td>N/A</td><td>N/A</td></tr><tr><td>0.040</td><td></td><td></td><td></td></tr><tr><td>0.100</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.300</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table>			Simulator	Serial #	Lot #	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			<table border="1"><thead><tr><th>Simulator</th><th>Serial Number</th></tr></thead><tbody><tr><td>0.000</td><td>MP6284</td></tr><tr><td>Interferent</td><td>MP6285</td></tr><tr><td>0.050</td><td>MP6286</td></tr><tr><td>0.080</td><td>MP6287</td></tr><tr><td>0.200</td><td>MP6288</td></tr></tbody></table>						Simulator	Serial Number	0.000	MP6284	Interferent	MP6285	0.050	MP6286	0.080	MP6287	0.200	MP6288
Simulator	Serial #	Lot #	Expiration																																													
0.000		N/A	N/A																																													
0.040																																																
0.100																																																
0.200																																																
0.300																																																
0.080 DGS	N/A																																															
Simulator	Serial Number																																															
0.000	MP6284																																															
Interferent	MP6285																																															
0.050	MP6286																																															
0.080	MP6287																																															
0.200	MP6288																																															
☐ Post Calibration Adjustment Stability Checks			<table border="1"><thead><tr><th colspan="2">Attachments</th></tr></thead><tbody><tr><td>☒ Form 41</td><td>☐ Post-Stability Checks</td></tr><tr><td>☒ Stability Checks</td><td>☐ Flow Calibration</td></tr><tr><td>☒ Calibration Certificate</td><td>☐ Form 40</td></tr><tr><td>☐ Calibration Adjustment</td><td>☐ Other _____</td></tr></tbody></table>						Attachments		☒ Form 41	☐ Post-Stability Checks	☒ Stability Checks	☐ Flow Calibration	☒ Calibration Certificate	☐ Form 40	☐ Calibration Adjustment	☐ Other _____																														
Attachments																																																
☒ Form 41	☐ Post-Stability Checks																																															
☒ Stability Checks	☐ Flow Calibration																																															
☒ Calibration Certificate	☐ Form 40																																															
☐ Calibration Adjustment	☐ Other _____																																															
Notes/Suggested Service: <u>Could not duplicate Modem Fail described by agency. (TDG 7/17/25)</u>			☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use																																													
			Digitally signed by Shayla Platt Platt Date: 2025.07.19 11:51:58 -04'00'																																													
			Phil Nicodemo Nicodemo Date: 2025.07.18 12:04:33 -04'00'																																													
			Tech Review / Date _____ Admin Review / Date _____																																													

Stability Checks

0.05g/210L	0.08g/210L	0.20g/210L	DGS 0.08g/210L
0.047 to 0.053	0.077 to 0.083	0.194 to 0.206	0.077 to 0.083
✓ SUPPORT EQ Intoxalyzer - Alcohol Analyzer Model: 8000 07/15/2025 SN: 80-00016 Software: 8100.27 Test Air Blank 0.000 14:33 Control Test 0.050 14:34 Air Blank 0.000 14:35 Control Test 0.050 14:35 Air Blank 0.000 14:36 Control Test 0.051 14:37 Air Blank 0.000 14:37 Control Test 0.053 14:37 Average 0.0503 Std Dev 0.0006 Rel. Std Dev(%) 1.1471	✓ SUPPORT EQ Intoxalyzer - Alcohol Analyzer Model: 8000 07/15/2025 SN: 80-00016 Software: 8100.27 Test Air Blank 0.000 14:39 Control Test 0.083 14:40 Air Blank 0.000 14:41 Control Test 0.082 14:41 Air Blank 0.000 14:42 Control Test 0.082 14:43 Air Blank 0.000 14:43 Control Test 0.083 14:43 Average 0.0823 Std Dev 0.0006 Rel. Std Dev(%) 0.7012	✓ SUPPORT EQ Intoxalyzer - Alcohol Analyzer Model: 8000 07/15/2025 SN: 80-00016 Software: 8100.27 Test Air Blank 0.000 14:51 Control Test 0.203 14:51 Air Blank 0.000 14:52 Control Test 0.203 14:52 Air Blank 0.000 14:53 Control Test 0.202 14:54 Air Blank 0.000 14:54 Control Test 0.2027 14:54 Average 0.2027 Std Dev 0.0006 Rel. Std Dev(%) 0.2849	✓ SUPPORT EQ Intoxalyzer - Alcohol Analyzer Model: 8000 07/15/2025 SN: 80-00016 Software: 8100.27 Test Air Blank 0.000 14:57 Control Test 0.080 14:57 Air Blank 0.000 14:58 Control Test 0.081 14:58 Air Blank 0.000 14:59 Control Test 0.080 14:59 Air Blank 0.000 14:59 Control Test 0.080 14:59 Control Test Status 0.0803 Average 0.0803 Std Dev 0.0006 Rel. Std Dev(%) 0.7187
Operator's Signature	Operator's Signature	Operator's Signature	Operator's Signature

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: GULFPORT PD

Time of Inspection: 13:28

Date of Inspection: 07/17/2025

Serial Number: 80-001061

Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202406K Exp: 06/19/2026	0.08g/210L Test (g/210L) Lot#:202406L Exp: 06/19/2026	0.20g/210L Test (g/210L) Lot#:202406N Exp: 06/20/2026	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG429602 Exp: 10/22/2026
0.000	0.051	0.082	0.203	0.080
0.000	0.050	0.082	0.204	0.080
0.000	0.050	0.081	0.203	0.080
0.000	0.051	0.081	0.203	0.080
0.000	0.050	0.082	0.204	0.080
0.000	0.051	0.082	0.204	0.080
0.000	0.051	0.082	0.204	0.080
0.000	0.051	0.082	0.204	0.079
0.000	0.051	0.082	0.204	0.080
0.000	0.051	0.083	0.204	0.080

Standard Deviations	0.0004	0.0005	0.0004	0.0003
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0004 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.


TAYLOR D GUTSCHOW
Signature and Printed Name

07/17/2025
Date



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
4700 Terminal Drive, Suite 1
Ft. Myers, FL 33907

This is to certify the calibration of Intoxilyzer 8000 serial number 80-001061, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-001061</u>	UNCERTAINTY* $\pm$	
Owning Agency:	<u>GULFPORT PD</u>	0.050 g/ 210 L	0.004
Calibration Date:	<u>07/17/2025</u>	0.080 g/ 210 L	0.004
Calibration Time:	<u>13:28</u>	0.200 g/ 210 L	0.007
		0.080 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.
*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.
This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

Taylor
Gutschow
Digitally signed by Taylor
Gutschow
Date: 2025.07.17 15:16:55
-04'00'

07/17/2025

Date

TAYLOR D GUTSCHOW,
Department Inspector

FDLE/ATP Form 69 December 2021

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality