



INSTRUMENT PROCESSING SHEET

Agency Volusia County SOS/N 80-001132Florida Department of
Law EnforcementDate In 2/3/2025DI Completion Date N/A☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake By <u>ALL</u> Date <u>02/03/2025</u> ☒ Annual ☐ Registration ☐ Return from CMI / EE Visual Inspection: ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/ Accessories: ☐ Power cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes: _____ _____ _____ _____ _____ _____ _____ _____ _____ _____	Quality Checks By <u>DA</u> Date <u>02/04/2025</u> ☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value <u>201</u> ☒ Flow Verification (L/s) Flow Column # <u>ATP102</u> 32 mm <u>0.140</u> (.139 - .169) 36 mm <u>0.156</u> (.156 - .190) 53 mm <u>0.226</u> (.228 - .278) 103 mm <u>0.492</u> (.447 - .547) ☒ Barometric Pressure Check Gauge ID # <u>30793</u> ☒ Stability Checks <table border="1" style="width:100%; border-collapse: collapse;"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #/Exp</th></tr></thead><tbody><tr><td>0.050</td><td>MP5088</td><td>202406K 06/19/2026</td></tr><tr><td>0.080</td><td>MP5089</td><td>202406L 06/19/2026</td></tr><tr><td>0.200</td><td>MP5090</td><td>202406N 06/20/2026</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>AG429602 10/22/2026</td></tr></tbody></table>	Simulator	Serial #	Lot #/Exp	0.050	MP5088	202406K 06/19/2026	0.080	MP5089	202406L 06/19/2026	0.200	MP5090	202406N 06/20/2026	0.080 DGS	N/A	AG429602 10/22/2026	Flow Calibration By _____ Date _____ Flow Column # _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value _____ ☐ Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547) Maintenance By _____ Date _____ ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other _____ _____ _____ _____ _____ _____																																																					
Simulator	Serial #	Lot #/Exp																																																																				
0.050	MP5088	202406K 06/19/2026																																																																				
0.080	MP5089	202406L 06/19/2026																																																																				
0.200	MP5090	202406N 06/20/2026																																																																				
0.080 DGS	N/A	AG429602 10/22/2026																																																																				
Calibration Adjustment By _____ Barometric Pressure Gauge _____ ID # _____ <table border="1" style="width:100%; border-collapse: collapse;"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.000</td><td></td><td>N/A</td><td>N/A</td></tr><tr><td>0.040</td><td></td><td></td><td></td></tr><tr><td>0.100</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.300</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table> ☐ Post Calibration Adjustment Stability Checks <table border="1" style="width:100%; border-collapse: collapse;"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.050</td><td></td><td></td><td></td></tr><tr><td>0.080</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table>		Simulator	Serial #	Lot #	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			Simulator	Serial #	Lot #	Expiration	0.050				0.080				0.200				0.080 DGS	N/A			Department Inspection By _____ Barometric Pressure ID# _____ Gauge _____ Instrument _____ Mouth Alcohol Solution Lot # _____ Acetone Stock Solution Lot # _____ <table border="1" style="width:100%; border-collapse: collapse;"><thead><tr><th>Simulator</th><th>Serial Number</th></tr></thead><tbody><tr><td>0.000</td><td></td></tr><tr><td>Interferent</td><td></td></tr><tr><td>0.050</td><td></td></tr><tr><td>0.080</td><td></td></tr><tr><td>0.200</td><td></td></tr></tbody></table> Attachments <table style="width:100%;"><tr><td>☐ Form 41</td><td>☐ Post-Stability Checks</td></tr><tr><td>☒ Stability Checks</td><td>☐ Flow Calibration</td></tr><tr><td>☐ Calibration Certificate</td><td>☐ Form 40</td></tr><tr><td>☐ Calibration Adjustment</td><td>☐ Other _____</td></tr></table> ☐ Instrument Complies with Chapter 11D-8, FAC ☒ Instrument Does Not Comply with Chapter 11D-8, FAC ☐ Return to/Place into Evidentiary Use ☒ Remain Out of Evidentiary Use ☐ Conduct an Agency Inspection Before Evidentiary Use Taylor Gutschow Digitally signed by Taylor Gutschow Date: 2025.07.28 13:54:15 +0400 Tech Review / Date _____ Admin Review / Date _____	Simulator	Serial Number	0.000		Interferent		0.050		0.080		0.200		☐ Form 41	☐ Post-Stability Checks	☒ Stability Checks	☐ Flow Calibration	☐ Calibration Certificate	☐ Form 40	☐ Calibration Adjustment	☐ Other _____
Simulator	Serial #	Lot #	Expiration																																																																			
0.000		N/A	N/A																																																																			
0.040																																																																						
0.100																																																																						
0.200																																																																						
0.300																																																																						
0.080 DGS	N/A																																																																					
Simulator	Serial #	Lot #	Expiration																																																																			
0.050																																																																						
0.080																																																																						
0.200																																																																						
0.080 DGS	N/A																																																																					
Simulator	Serial Number																																																																					
0.000																																																																						
Interferent																																																																						
0.050																																																																						
0.080																																																																						
0.200																																																																						
☐ Form 41	☐ Post-Stability Checks																																																																					
☒ Stability Checks	☐ Flow Calibration																																																																					
☐ Calibration Certificate	☐ Form 40																																																																					
☐ Calibration Adjustment	☐ Other _____																																																																					
Notes/Suggested Service: <u>Ran stabilities and DGS produced nominal values. However, upon running the simulators, all of the ethanol control solutions yielded zeroes. Instrument needs to be sent to repair. DA 02/04/2025</u> <u>No response on repairs. Returning to agency in disabled mode. SP 7/26/25</u> Admin Review: Flow measurement outside range. SP 7/30/25																																																																						

Stability Checks

80-001132 DA

2/4/25

LOUISIANA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001132
02/04/2025
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:54
Control Test	0.000	09:55
Air Blank	0.000	09:56
Control Test	0.000	09:56
Air Blank	0.000	09:57
Control Test	0.000	09:57
Air Blank	0.000	09:58
Control Test	0.000	09:58
Control Test Stats		
Average	0.0000	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

0.05g/210L

DA

Operator's Signature

LOUISIANA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001132
02/04/2025
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:00
Control Test	0.000	10:01
Air Blank	0.000	10:01
Control Test	0.000	10:02
Air Blank	0.000	10:02
Control Test	0.000	10:03
Air Blank	0.000	10:04
Control Test	0.000	10:04
Control Test Stats		
Average	0.0000	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

0.05g/210L

DA

Operator's Signature

LOUISIANA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001132
02/04/2025
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:20
Control Test	0.000	10:20
Air Blank	0.000	10:21
Control Test	0.000	10:21
Air Blank	0.000	10:22
Control Test	0.000	10:23
Air Blank	0.000	10:23
Control Test	0.000	10:23
Control Test Stats		
Average	0.0000	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

wet

0.08g/210L

DA

Operator's Signature

LOUISIANA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001132
02/04/2025
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	10:10
Control Test	0.000	10:11
Air Blank	0.000	10:11
Control Test	0.000	10:12
Air Blank	0.000	10:12
Control Test	0.000	10:13
Air Blank	0.000	10:13
Control Test	0.000	10:13
Control Test Stats		
Average	0.0000	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

0.200g/210L

DA

Operator's Signature

LOUISIANA COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001132
02/04/2025
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:35
Control Test	0.078	09:35
Air Blank	0.000	09:36
Control Test	0.077	09:36
Air Blank	0.000	09:36
Control Test	0.077	09:37
Air Blank	0.000	09:37
Control Test	0.000	09:37
Control Test Stats		
Average	0.0773	
Std Dev	0.0006	
Rel Std Dev(%)	0.7466	

DGS

DA

Operator's Signature

Repeated 0.05g/210L test to ensure correct stability test was selected. Checked all seals for leaks and ensured proper connection with each stability test and all results from wet bath concentrations were zeroes.

DA 2/4/25