

Alcohol Testing Program - Instrument Processing Sheet

Agency: Niceville PD

Instrument Serial Number: 80-001310

Date In: 8/3/2026

DI Completion Date: 8/6/2026

☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake By: IAS Date: 8/3/2026																				
☒ Annual ☐ Dropped Off ☐ Registration ☒ Return from CMI / EE ☐ Training Instrument Visual Inspection ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/Accessories ☒ Power Cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes:			Quality Checks By: WKP Date: 8/6/2026 ☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value: 122 _____ ☒ Flow Verification (L/S) Flow Column #: ATP 105 32 mm .152 (.139-.169) 36 mm .167 (.156-.190) 53 mm .242 (.228-.278) 103 mm .507 (.447-.547) ☒ Barometric Pressure Check Gauge ID #: 28427 Gauge: 1017 Instrument: 1019 ☒ Stability Checks <table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot#/Exp</th></tr></thead><tbody><tr><td>0.050</td><td>MP6291</td><td>202603A 3/3/2028</td></tr><tr><td>0.080</td><td>MP6292</td><td>202603B 3/3/2028</td></tr><tr><td>0.200</td><td>MP6293</td><td>202603D 3/4/2028</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>23725080A1 10/5/2027</td></tr></tbody></table>			Simulator	Serial #	Lot#/Exp	0.050	MP6291	202603A 3/3/2028	0.080	MP6292	202603B 3/3/2028	0.200	MP6293	202603D 3/4/2028	0.080 DGS	N/A	23725080A1 10/5/2027
Simulator	Serial #	Lot#/Exp																		
0.050	MP6291	202603A 3/3/2028																		
0.080	MP6292	202603B 3/3/2028																		
0.200	MP6293	202603D 3/4/2028																		
0.080 DGS	N/A	23725080A1 10/5/2027																		
Flow Adjustment By: _____ Flow Column #: _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value: _____ ☐ Post Adjustment Verification (L/S) Flow Column #: _____ 32 mm _____(.139-.169) 36 mm _____(.156-.190) 53 mm _____(.228-.278) 103 mm _____(.447-.547)																				
Maintenance By: _____ Date: _____ ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Tank Sensor Tare ☐ Breath Tube Replacement ☐ Other:																				
Optical Bench Adjustment By: _____			Department Inspection By: WKP																	
Barometric Pressure Gauge: _____ ID#: _____			Barometric Pressure ID#: 28427																	
Simulator	Serial #	Lot #	Expiration	Gauge: 1017	Instrument: 1019															
0.000		N/A	N/A	Mouth Alcohol Solution Lot #: 2025-D	Exp: 9/25/2027															
0.040				Acetone Stock Solution Lot #: 2025-B	Exp: 9/22/2027															
0.100				Simulator	Serial Number															
0.200				0.000	MP6289															
0.300				Interferent	MP6290															
0.080 DGS	N/A			0.050	MP6291															
☐ Post Optical Bench Adjustment Stability Checks				0.080	MP6292															
Simulator	Serial #	Lot #	Expiration	0.200	MP6293															
0.050				Attachments																
0.080				☒ Form 41																
0.200				☒ Stability Checks																
0.080 DGS	N/A			☒ Calibration Certificate																
				☐ Optical Bench Adjustment																
Gauge ID #: _____				☐ Post-Stability Checks																
Gauge: _____ Instrument: _____				☐ Flow Adjustment																
				☐ Form 40																
				☒ Other: Calculation Form																
Notes/Suggested Service:				☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use																
				LeAndra Higginbotham Tech Review Digitally signed by LeAndra Higginbotham Date: 2026.08.13 10:58:57 -04'00'																
				Kaitlyn Spearin Admin Review Digitally signed by Kaitlyn Spearin Date: 2026.08.13 11:59:19 -04'00'																

Stability Checks

0.05g/210L 0.047 to 0.053 ✓	0.08g/210L 0.077 to 0.083 ✓	0.20g/210L 0.194 to 0.206 ✓	DGS 0.08g/210L 0.077 to 0.083 ✓ ≤0.003 of Wet ✓																																																																																																																																																
NICEVILLE PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001310 08/06/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>09:38</td></tr> <tr><td>Control Test</td><td>0.050</td><td>09:38</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:39</td></tr> <tr><td>Control Test</td><td>0.049</td><td>09:39</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:40</td></tr> <tr><td>Control Test</td><td>0.049</td><td>09:41</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:41</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0493</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>1.1703</td><td></td></tr> </tbody> </table> <i>menp</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	09:38	Control Test	0.050	09:38	Air Blank	0.000	09:39	Control Test	0.049	09:39	Air Blank	0.000	09:40	Control Test	0.049	09:41	Air Blank	0.000	09:41	Control Test Stats			Average	0.0493		Std Dev	0.0006		Rel Std Dev(%)	1.1703		NICEVILLE PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001310 08/06/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>09:48</td></tr> <tr><td>Control Test</td><td>0.078</td><td>09:48</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:49</td></tr> <tr><td>Control Test</td><td>0.079</td><td>09:50</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:50</td></tr> <tr><td>Control Test</td><td>0.078</td><td>09:51</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:51</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0780</td><td></td></tr> <tr><td>Std Dev</td><td>0.0000</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr> </tbody> </table> <i>menp</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	09:48	Control Test	0.078	09:48	Air Blank	0.000	09:49	Control Test	0.079	09:50	Air Blank	0.000	09:50	Control Test	0.078	09:51	Air Blank	0.000	09:51	Control Test Stats			Average	0.0780		Std Dev	0.0000		Rel Std Dev(%)	0.0000		NICEVILLE PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001310 08/06/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>09:53</td></tr> <tr><td>Control Test</td><td>0.201</td><td>09:53</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:54</td></tr> <tr><td>Control Test</td><td>0.201</td><td>09:55</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:55</td></tr> <tr><td>Control Test</td><td>0.200</td><td>09:56</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:56</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.2007</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.2877</td><td></td></tr> </tbody> </table> <i>menp</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	09:53	Control Test	0.201	09:53	Air Blank	0.000	09:54	Control Test	0.201	09:55	Air Blank	0.000	09:55	Control Test	0.200	09:56	Air Blank	0.000	09:56	Control Test Stats			Average	0.2007		Std Dev	0.0006		Rel Std Dev(%)	0.2877		NICEVILLE PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001310 08/06/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>10:03</td></tr> <tr><td>Control Test</td><td>0.078</td><td>10:03</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>10:04</td></tr> <tr><td>Control Test</td><td>0.079</td><td>10:04</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>10:05</td></tr> <tr><td>Control Test</td><td>0.079</td><td>10:05</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>10:06</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0787</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.7339</td><td></td></tr> </tbody> </table> <i>menp</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	10:03	Control Test	0.078	10:03	Air Blank	0.000	10:04	Control Test	0.079	10:04	Air Blank	0.000	10:05	Control Test	0.079	10:05	Air Blank	0.000	10:06	Control Test Stats			Average	0.0787		Std Dev	0.0006		Rel Std Dev(%)	0.7339	
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	09:38																																																																																																																																																	
Control Test	0.050	09:38																																																																																																																																																	
Air Blank	0.000	09:39																																																																																																																																																	
Control Test	0.049	09:39																																																																																																																																																	
Air Blank	0.000	09:40																																																																																																																																																	
Control Test	0.049	09:41																																																																																																																																																	
Air Blank	0.000	09:41																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0493																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	1.1703																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	09:48																																																																																																																																																	
Control Test	0.078	09:48																																																																																																																																																	
Air Blank	0.000	09:49																																																																																																																																																	
Control Test	0.079	09:50																																																																																																																																																	
Air Blank	0.000	09:50																																																																																																																																																	
Control Test	0.078	09:51																																																																																																																																																	
Air Blank	0.000	09:51																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0780																																																																																																																																																		
Std Dev	0.0000																																																																																																																																																		
Rel Std Dev(%)	0.0000																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	09:53																																																																																																																																																	
Control Test	0.201	09:53																																																																																																																																																	
Air Blank	0.000	09:54																																																																																																																																																	
Control Test	0.201	09:55																																																																																																																																																	
Air Blank	0.000	09:55																																																																																																																																																	
Control Test	0.200	09:56																																																																																																																																																	
Air Blank	0.000	09:56																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.2007																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	0.2877																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	10:03																																																																																																																																																	
Control Test	0.078	10:03																																																																																																																																																	
Air Blank	0.000	10:04																																																																																																																																																	
Control Test	0.079	10:04																																																																																																																																																	
Air Blank	0.000	10:05																																																																																																																																																	
Control Test	0.079	10:05																																																																																																																																																	
Air Blank	0.000	10:06																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0787																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	0.7339																																																																																																																																																		

Florida Department of Law Enforcement

Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: NICEVILLE PD
Time of Inspection: 13:12

Date of Inspection: 08/06/2026

Serial Number: 80-001310
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#: 202603A Exp: 03/03/2028	0.08g/210L Test (g/210L) Lot#: 202603B Exp: 03/03/2028	0.20g/210L Test (g/210L) Lot#: 202603D Exp: 03/04/2028	0.08 g/210L Dry Gas Std Test (g/210L) Lot#: 23725080A1 Exp: 10/05/2027
0.000	0.047	0.078	0.199	0.078
0.000	0.048	0.078	0.199	0.078
0.000	0.049	0.078	0.199	0.077
0.000	0.048	0.078	0.199	0.077
0.000	0.049	0.078	0.199	0.077
0.000	0.049	0.077	0.199	0.077
0.000	0.049	0.078	0.198	0.077
0.000	0.049	0.077	0.198	0.077
0.000	0.049	0.078	0.198	0.077
0.000	0.049	0.078	0.199	0.077

Standard Deviations	0.0006	0.0004	0.0004	0.0004
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0004 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



WEN-CHI K PIERSON

Signature and Printed Name

08/06/2026
Date

Florida Department of Law Enforcement
Alcohol Testing Program - Average Calculations for Calibration Certificate

Calculations By: WKP

Serial Number: 80-001310

Date: 8/7/2026

Solution:

0.05 g/210 L Test: 0.486 / 10 = 0.049 g/210L
(SUM) (AVERAGE)

0.08 g/210L Test: 0.778 / 10 = 0.078 g/210L
(SUM) (AVERAGE)

0.20 g/210L Test: 1.987 / 10 = 0.199 g/210L
(SUM) (AVERAGE)

0.08 g/210L Dry Gas Std Test: 0.772 / 10 = 0.077 g/210L
(SUM) (AVERAGE)

Comments:

Average = (Sum of observations) / (Total number of observations)

All results are reported in g/210L



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
2331 Phillips Road
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-001310, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-001310</u>	Average of $y \pm$ Uncertainty*	
Owning Agency:	<u>NICEVILLE PD</u>	0.049 g/ 210 L	0.004
Calibration Date:	<u>08/06/2026</u>	0.078 g/ 210 L	0.004
Calibration Time:	<u>13:12</u>	0.199 g/ 210 L	0.008
		0.077 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

08/06/2026

Date

Wen-Chi
Pierson

Digitally signed by Wen-Chi Pierson
Date: 2026.08.07 09:22:10
-04'00'

WEN-CHI K PIERSON,
Department Inspector

FDLE/ATP Form 69 August 2026

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality

Page 1 of 1