

Alcohol Testing Program - Instrument Processing Sheet

Agency: Baker County SO

Instrument Serial Number: 80-001287

Date In: 6/3/2026

DI Completion Date: 6/9/2026

☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake By: SLH Date: 6/3/2026				Quality Checks By: WKP Date: 6/5/2026				Flow Adjustment By: _____			
☒ Annual ☐ Dropped Off ☐ Registration ☐ Return from CMI / EE ☐ Training Instrument Visual Inspection ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/Accessories ☒ Power Cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes:				☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value: 188 ☒ Flow Verification (L/s) Flow Column #: ATP 105 32 mm 0.152 (.139-.169) 36 mm 0.167 (.156-.190) 53 mm 0.246 (.228-.278) 103 mm 0.527 (.447-.547) ☒ Barometric Pressure Check Gauge ID #: 28662 Gauge: 1017 Instrument: 1012 ☒ Stability Checks				Flow Column #: _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value: _____ ☐ Post Adjustment Verification (L/S) Flow Column #: _____ 32 mm _____ (.139-.169) 36 mm _____ (.156-.190) 53 mm _____ (.228-.278) 103 mm _____ (.447-.547)			
				Simulator		Serial #		Lot#/Exp		Maintenance By: WKP Date: 6/5/2026	
				0.050		MP5088		202603A 3/3/2028		☐ Battery Replacement	
				0.080		MP5089		202603B 3/3/2028		☒ Dry Gas Regulator Replacement	
				0.200		MP5090		202603D 3/4/2028		☒ Tank Sensor Tare	
0.080 DGS		N/A		23825080A1 10/5/2027		☐ Breath Tube Replacement					
Optical Bench Adjustment By: _____				Department Inspection By: WKP							
Barometric Pressure Gauge: _____ ID#: _____				Barometric Pressure ID#: 34419							
Simulator		Serial #		Lot #		Expiration		Gauge: 1018 Instrument: 1014			
0.000				N/A		N/A		Mouth Alcohol Solution Lot #: 2025-D Exp: 9/25/2027			
0.040								Acetone Stock Solution Lot #: 2025-B Exp: 9/22/2027			
0.100								Simulator			
0.200								Serial Number			
0.300								0.000			
0.080 DGS		N/A						MP5086			
☐ Post Optical Bench Adjustment Stability Checks								Interferent			
Simulator		Serial #		Lot #		Expiration		MP5087			
0.050								0.050			
0.080								MP5088			
0.200								0.080			
0.080 DGS		N/A						MP5089			
Gauge ID #: _____								0.200			
Gauge: _____ Instrument: _____								MP5090			
Notes/Suggested Service: Admin Correction: Added instrument correspondence relating to the dry gas regulator replacement to the respective instrument correspondence folder. WKP 6/15/2026				Attachments							
				☒ Form 41				☐ Post-Stability Checks			
				☒ Stability Checks				☐ Flow Adjustment			
				☒ Calibration Certificate				☐ Form 40			
				☐ Optical Bench Adjustment				☐ Other:			
				☒ Instrument Complies with Chapter 11D-8, FAC							
				☐ Instrument Does Not Comply with Chapter 11D-8, FAC							
				☒ Return to/Place into Evidentiary Use							
				☐ Remain Out of Evidentiary Use							
				☒ Conduct an Agency Inspection Before Evidentiary Use							
				LeAndra Higginbotham				Kaitlyn Spearin			
				Digitally signed by LeAndra Higginbotham Date: 2026.06.14 14:13:58 -04'00'				Digitally signed by Kaitlyn Spearin Date: 2026.06.15 11:11:37 -04'00'			
				Tech Review				Admin Review			

Stability Checks

0.05g/210L 0.047 to 0.053	0.08g/210L 0.077 to 0.083	0.20g/210L 0.194 to 0.206	DGS 0.08g/210L 0.077 to 0.083 ≤0.003 of Wet																																																																																																																																																
BAKER COUNTY SO Intoxilyzer - Alconol Analyzer Model 8000 SN 80-001287 06/05/2026 Software: 8.00.27 <table><thead><tr><th>Test</th><th>g/210L</th><th>Time</th></tr></thead><tbody><tr><td>Air Blank</td><td>0.000</td><td>12:14</td></tr><tr><td>Control Test</td><td>0.052</td><td>12:14</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:15</td></tr><tr><td>Control Test</td><td>0.051</td><td>12:15</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:16</td></tr><tr><td>Control Test</td><td>0.051</td><td>12:17</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:17</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0513</td><td></td></tr><tr><td>Std Dev</td><td>0.0006</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>1.1247</td><td></td></tr></tbody></table> <i>menpi</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	12:14	Control Test	0.052	12:14	Air Blank	0.000	12:15	Control Test	0.051	12:15	Air Blank	0.000	12:16	Control Test	0.051	12:17	Air Blank	0.000	12:17	Control Test Stats			Average	0.0513		Std Dev	0.0006		Rel Std Dev(%)	1.1247		BAKER COUNTY SO Intoxilyzer - Alconol Analyzer Model 8000 SN 80-001287 06/05/2026 Software: 8.00.27 <table><thead><tr><th>Test</th><th>g/210L</th><th>Time</th></tr></thead><tbody><tr><td>Air Blank</td><td>0.000</td><td>12:30</td></tr><tr><td>Control Test</td><td>0.079</td><td>12:31</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:31</td></tr><tr><td>Control Test</td><td>0.079</td><td>12:32</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:33</td></tr><tr><td>Control Test</td><td>0.079</td><td>12:33</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:34</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0790</td><td></td></tr><tr><td>Std Dev</td><td>0.0000</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr></tbody></table> <i>menpi</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	12:30	Control Test	0.079	12:31	Air Blank	0.000	12:31	Control Test	0.079	12:32	Air Blank	0.000	12:33	Control Test	0.079	12:33	Air Blank	0.000	12:34	Control Test Stats			Average	0.0790		Std Dev	0.0000		Rel Std Dev(%)	0.0000		BAKER COUNTY SO Intoxilyzer - Alconol Analyzer Model 8000 SN 80-001287 06/05/2026 Software: 8.00.27 <table><thead><tr><th>Test</th><th>g/210L</th><th>Time</th></tr></thead><tbody><tr><td>Air Blank</td><td>0.000</td><td>12:25</td></tr><tr><td>Control Test</td><td>0.200</td><td>12:25</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:26</td></tr><tr><td>Control Test</td><td>0.198</td><td>12:27</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:27</td></tr><tr><td>Control Test</td><td>0.198</td><td>12:28</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:29</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.1987</td><td></td></tr><tr><td>Std Dev</td><td>0.0112</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.5812</td><td></td></tr></tbody></table> <i>menpi</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	12:25	Control Test	0.200	12:25	Air Blank	0.000	12:26	Control Test	0.198	12:27	Air Blank	0.000	12:27	Control Test	0.198	12:28	Air Blank	0.000	12:29	Control Test Stats			Average	0.1987		Std Dev	0.0112		Rel Std Dev(%)	0.5812		BAKER COUNTY SO Intoxilyzer - Alconol Analyzer Model 8000 SN 80-001287 06/05/2026 Software: 8.00.27 <table><thead><tr><th>Test</th><th>g/210L</th><th>Time</th></tr></thead><tbody><tr><td>Air Blank</td><td>0.000</td><td>12:09</td></tr><tr><td>Control Test</td><td>0.080</td><td>12:09</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:10</td></tr><tr><td>Control Test</td><td>0.080</td><td>12:10</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:11</td></tr><tr><td>Control Test</td><td>0.080</td><td>12:11</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:11</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0803</td><td></td></tr><tr><td>Std Dev</td><td>0.0006</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.7187</td><td></td></tr></tbody></table> <i>menpi</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	12:09	Control Test	0.080	12:09	Air Blank	0.000	12:10	Control Test	0.080	12:10	Air Blank	0.000	12:11	Control Test	0.080	12:11	Air Blank	0.000	12:11	Control Test Stats			Average	0.0803		Std Dev	0.0006		Rel Std Dev(%)	0.7187	
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	12:14																																																																																																																																																	
Control Test	0.052	12:14																																																																																																																																																	
Air Blank	0.000	12:15																																																																																																																																																	
Control Test	0.051	12:15																																																																																																																																																	
Air Blank	0.000	12:16																																																																																																																																																	
Control Test	0.051	12:17																																																																																																																																																	
Air Blank	0.000	12:17																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0513																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	1.1247																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	12:30																																																																																																																																																	
Control Test	0.079	12:31																																																																																																																																																	
Air Blank	0.000	12:31																																																																																																																																																	
Control Test	0.079	12:32																																																																																																																																																	
Air Blank	0.000	12:33																																																																																																																																																	
Control Test	0.079	12:33																																																																																																																																																	
Air Blank	0.000	12:34																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0790																																																																																																																																																		
Std Dev	0.0000																																																																																																																																																		
Rel Std Dev(%)	0.0000																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	12:25																																																																																																																																																	
Control Test	0.200	12:25																																																																																																																																																	
Air Blank	0.000	12:26																																																																																																																																																	
Control Test	0.198	12:27																																																																																																																																																	
Air Blank	0.000	12:27																																																																																																																																																	
Control Test	0.198	12:28																																																																																																																																																	
Air Blank	0.000	12:29																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.1987																																																																																																																																																		
Std Dev	0.0112																																																																																																																																																		
Rel Std Dev(%)	0.5812																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	12:09																																																																																																																																																	
Control Test	0.080	12:09																																																																																																																																																	
Air Blank	0.000	12:10																																																																																																																																																	
Control Test	0.080	12:10																																																																																																																																																	
Air Blank	0.000	12:11																																																																																																																																																	
Control Test	0.080	12:11																																																																																																																																																	
Air Blank	0.000	12:11																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0803																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	0.7187																																																																																																																																																		

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: BAKER COUNTY SO
Time of Inspection: 15:05

Date of Inspection: 06/09/2026

Serial Number: 80-001287
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202603A Exp: 03/03/2028	0.08g/210L Test (g/210L) Lot#:202603B Exp: 03/03/2028	0.20g/210L Test (g/210L) Lot#:202603D Exp: 03/04/2028	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:23825080A1 Exp: 10/05/2027
0.000	0.049	0.080	0.198	0.079
0.000	0.049	0.080	0.199	0.079
0.000	0.049	0.080	0.199	0.079
0.000	0.050	0.079	0.199	0.078
0.000	0.049	0.080	0.198	0.079
0.000	0.050	0.080	0.199	0.078
0.000	0.050	0.080	0.199	0.078
0.000	0.050	0.080	0.199	0.078
0.000	0.050	0.079	0.199	0.079
0.000	0.050	0.080	0.199	0.078

Standard Deviations	0.0005	0.0004	0.0004	0.0005
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0004 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Wen-Chi K Pierson WEN-CHI K PIERSON
Signature and Printed Name

06/09/2026
Date



Florida Department of Law Enforcement
Alcohol Testing Program
2331 Phillips Road
Tallahassee, FL 32308

Calibration Certificate

This is to certify the calibration of Intoxilyzer 8000 serial number 80-001287, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-001287</u>	UNCERTAINTY* $\pm$
Owning Agency:	<u>BAKER COUNTY SO</u>	0.050 g/ 210 L 0.004
Calibration Date:	<u>06/09/2026</u>	0.080 g/ 210 L 0.004
Calibration Time:	<u>15:05</u>	0.200 g/ 210 L 0.008
		0.080 g/ 210 L Dry Gas Control 0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

Wen-Chi

Pierson

Digitally signed by
Wen-Chi Pierson
Date: 2026.06.10
09:14:21 -04'00'

06/09/2026

Date

WEN-CHI K PIERSON,

Department Inspector

FDLE/ATP Form 69 January 2026

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality